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
EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 - 20099 SESTO SAN GIOVANNI (MILANO) - ITALY
PRODUCT	HP1gG CODE: HP.G.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

WE HEREBY DECLARE THAT THE ABOVE MENTIONED PRODUCT MEETS THE PROVISIONS OF THE COUNCIL DIRECTIVE 98/79/EC FOR IN VITRO DIAGNOSTIC DEVICES.

ISO CERTIFICATE(S)	UNI CEI EN ISO 13485-N° 50 100 5931/B RELEASED BY CERTIFICATION BODY TÜV Italia S.r.l.
--------------------	---

PLACE & DATE OF FIRST ISSUE	MILANO - MARCH 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) - MAY 2018
SIGNATURE Legal Representative Dr. ssa Fiorenza Scozzesi	 DIA.PRO DIAGNOSTIC BIOPROBES S.R.L.

Rev: 05/2018

DIA.PRO Diagnostic Bioprobes S.r.l.
Sede legale e lab.: Via G. Carducci, 27 - 20099 Sesto S. Giovanni (MI) - Italia
Tel. +39 02 27007161/6450 • Fax +39 02 44386771 • <http://www.diapro.it> • E-mail: info@diapro.it
Capitale sociale €50.000,00 I.V. - P.IVA: 11924660159 - Reg. Imp. 11924660159 - REA 1509959

HBsAg one

Version ULTRA

Fourth generation Enzyme Immunoassay (ELISA) for the determination of Hepatitis B surface Antigen or HBsAg in human serum and plasma

- for "in vitro" diagnostic use only -

DIA.PRO
Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy
Phone +39 02 27007161
Fax +39 02 2607726
e-mail: info@diagpro.it

REF SAGI/ULTRA CE
96/192/480/960 Tessi

HBsAg One version ULTRA

A. INTENDED USE
Fourth generation Enzyme Immunoassay (ELISA) for the one-step determination of Hepatitis B surface Antigen or HBsAg in human plasma and sera.
The kit is intended for the screening of blood units, is able to detect HBsAg mutants and finds application in the follow-up of HBV-infected patients.
For "in vitro" diagnostic use only.

B. INTRODUCTION
The World Health Organization (WHO) defines Hepatitis B Virus infection as follows:

Hepatitis B is one of the major diseases of mankind and is a serious global public health problem. Hepatitis B virus (HBV) causes hepatitis A, B, C, D, and E. All of these viruses can cause an acute disease with symptoms lasting several weeks including yellowing of the skin and eyes (jaundice); dark urine; extreme fatigue; nausea; vomiting and abdominal pain. It can take several months to a year to feel fit again. HBsAg is a virus that can cause chronic infection in which the patient never gets rid of the virus and many years later develops carcinoma of the liver or liver cancer.

HBV is the most serious type of viral hepatitis and the only type causing chronic hepatitis for which a vaccine is available. Hepatitis B virus is transmitted by contact with blood or body fluids of an infected person in the same way as human immunodeficiency virus (HIV), the virus that causes AIDS. However, HBV is 50 to 100 times more infectious than HIV. About 10% of people getting infected with HBV are: (a) perinatally from mother to baby at the time of birth; (b) by sexual transmission; (c) unsafe injections and transfusions; (d) sexual contact.

Worldwide, most infections occur from infected mother to child, from child to child contact in household settings, and from reuse of un-sterilized needles and syringes in many developing countries, almost all children become infected with the virus. In many industrialized countries (e.g. Western Europe and North America), the pattern of transmission is transmission between courtesans, mother-to-child and child-to-child (hand-to-hand, hepatitis B vaccination, shared needles, and intravenous drug use). However, the majority of infections in these countries are acquired during young adulthood for sexual activity, and injecting drug use. In addition, hepatitis B virus is the major infectious occupational hazard of health workers, and most health care workers have received hepatitis B vaccine.

Hepatitis B virus is not spread by contaminated food or water, and cannot be spread casually in the workplace. High rates of chronic HBV infection are seen in India, China, and several countries in Southeastern Europe. In the Middle East and Indian sub-continent, about 5% of the population are infected. Infection is less common in Western Europe and North America, where less than 1% are chronically infected.

Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 40% to 50% of children infected between 1 to 4 years of age develop chronic infection. The risk of death from HBV-related liver disease is about 1% per person who become chronically infected during childhood. However, patients with chronic hepatitis B have not yet been infected in many countries where 10% to 15% of children have not yet been infected. In many countries where the rate of chronic infection has been reduced, it has been found that the rate of chronic infection has been reduced in both the general population and in groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programs.

Hepatitis B vaccine has an outstanding record of safety, and effectiveness. Since 1982, over one billion doses of hepatitis B vaccine have been used worldwide. The vaccine is given as a series of three intramuscular doses. Studies have shown that the vaccine is 95% effective in preventing children and adults from developing chronic hepatitis B. The vaccine is also effective in preventing chronic hepatitis B in 15% of children who have not yet been infected. In many countries where the rate of chronic infection has been reduced, it has been found that the rate of chronic infection has been reduced in both the general population and in groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programs.

Hepatitis B surface Antigen or HBsAg is the most important protein of the envelope of Hepatitis B Virus, responsible for acute and chronic viral hepatitis.

The surface antigen contains the determinant "a", common to all the known viral subtypes, immunologically distinguished by two distinct subgroups (ay and ad).
The ability to detect HBsAg with high sensitive immunoassays in the last years has led to an understanding of its distribution and epidemiology worldwide and to radically decrease the risk of infection in transfusion.

C. PRINCIPLE OF THE TEST
A mix of mouse monoclonal antibodies specific to the determinants "a", "d" and "y" of HBsAg is fixed to the surface of microwells. Patient's serum/plasma is added to the microwell together with a second mix of mouse monoclonal antibodies, conjugated with Horseradish Peroxidase (HRP) and directed against a different epitope of the determinant "a" and against "preS".

The specific immunocomplex, formed in the presence of HBsAg in the sample, is captured by the solid phase.
At the end of the one-step incubation, microwells are washed to remove unbound serum proteins and HRP conjugate.
The chromogen/substrate is then added and, in the presence of captured HBsAg immunocomplex, the colorless substrate is hydrolyzed by the bound HRP conjugate to a colored end-product. After blocking the enzymatic reaction, its optical density is measured by an ELISA reader.

The color intensity is proportional to the amount of HBsAg present in the sample.
The version ULTRA is particularly suitable for automated screenings and is able to detect "s" mutants.

D. COMPONENTS

The standard configuration contains reagents to perform 192 tests and is made of the following components:

- 1. Microparticle MICROPOLATE**
n° 2 - 12 strips of 8 breakable wells coated with anti HBsAg affinity purified mouse monoclonal antibodies, specific to "a", "y" and "d" determinants, and sealed into a bag with desiccant.
- 2. Negative Control [CONTROL -]**
1x4 (mL/vial). Ready to use control. It contains goat serum, 10 mM phosphate buffer, pH 7.4±0.1, 0.09% Na-azide and 0.1% Kathon GC as preservatives. The negative control is pale yellow color coded.
- 3. Positive Control [CONTROL +]**
1x4 (mL/vial). Ready to use control. It contains goat serum, non infectious recombinant HBsAg, 10 mM phosphate buffer, pH 7.4±0.1, 0.02% gentamicin sulphate and 0.1% Kathon GC as preservatives. The positive control is color coded green.
- 4. Calibrator [CAL...]**
n° 2 vials. Lyophilized calibrator. To be dissolved with EIA grade water as reported in the label. Contains (total-bovine serum) non infectious recombinant HBsAg at 0.5 (U/ml) (2% v/v) International standard for HBsAg, NIBSC code 00/598), 10 mM phosphate buffer, pH 7.4±0.1, 0.02% gentamicin sulphate and 0.1% Kathon GC as preservatives.
- 5. Wash buffer concentrate [WASHBUF 20X]**
2x60 (mL/bottle). 20X concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer, pH 7.3±0.1, 0.2% Tween 20 and 0.1% Kathon GC.

REFERENCES

1. Aach R.D., Gristham J.W., Pecker S.W., Detection of Australia antigen by radioimmunoassay, Proc.Natl.Acad.Sci.USA, 68:1595, 1971.
2. Blumberg B.S., Srinock A.I., London W.T., Hepatitis and leukemia: their relation to Australia antigen. Bull.N.Y.Acad.Med., 44:1580, 1968.
3. Bonadio A., Davis M., Malheja R., The use of enzyme-linked immunosorbent assay for screening hybridoma antibodies against hepatitis B surface antigen. J.Immunol.Meth., 49:1, 1982.
4. Sedgwick C.H., Barrer D.T., Enzyme immunoassay for hepatitis B surface antigen. J.Clin.Med., 61: 305, 1977.
5. Fazakas S., De St.Groth, Scheldens D., Immunological and serological antibodies: strategy and tactics. J.Immunol.Meth., 35: 1, 1980.
6. Reesink H.W., et al., Comparison of six 3rd generation tests for the detection of HBeAg. Vox Sang., 38:51, 1980.
7. Rook G.A.W., Chromogens for the enzyme-linked immunosorbent assay (ELISA) using horseradish peroxidase. Labr.Rev., 52: 281, 1979.
8. Schroder J., Monoclonal antibodies: a new tool for research and immunodiagnosis. Med.Biol., 58:281, 1981.
9. Coleman P.C., Chen Y.C., Mutschauer K., Immunoassay detection of hepatitis B surface antigen mutants. J.Med.Virol., 1992:59(1):9-24

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia Pro Diagnostic Bioprobes S.r.l.
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

CE
0318

HBS AB

A. INTENDED USE
Enzyme ImmunoAssay (EUSA) for both the quantitative and qualitative determination of antibodies to the Surface Antigen of Hepatitis B Virus in human plasma and sera.
For "in vitro" diagnostic use only.

B. INTRODUCTION
The World Health Organization (WHO) defines Hepatitis B Virus infection as follows:

"Hepatitis B is one of the major diseases of mankind and is a serious global public health problem. Hepatitis means inflammation of the liver, and the most common cause is infection with one of 5 viruses, called hepatitis A,B,C,D, and E. All of these viruses can cause an acute disease with symptoms lasting several weeks including yellowing of the skin and eyes (jaundice); dark urine; extreme fatigue; nausea; vomiting and abdominal pain. It can take several months to a year to feel fit again. Hepatitis B virus can cause chronic infection in which the patient never gets rid of the virus and many years later develops cirrhosis of the liver or liver cancer.

HBV is the most serious type of viral hepatitis and the only type causing chronic hepatitis for which a vaccine is available. Hepatitis B virus is transmitted by contact with blood or body fluids of an infected person in the same way as human immunodeficiency virus (HIV), the virus that causes AIDS. However, HBV is 50 to 100 times more infectious than HIV. The main ways of getting infected with HBV are: (a) perinatal (from mother to baby at the birth); (b) child to child transmission; (c) unsafe injections and transfusions; (d) sexual contact.

Worldwide, most infections occur from infected mother to child, from child to child contact, in household settings and from reuse of unsterilized needles and syringes. In many developing countries, almost all children become infected with the virus. In many industrialized countries (e.g. Western Europe and North America), the pattern of transmission is different. In these countries, mother-to-infant and child-to-child transmission accounted for up to one third of chronic infections before childhood hepatitis B vaccination programmes were implemented. However, the majority of infections in these countries are acquired during young adulthood by sexual activity, and injecting drug use. In addition, hepatitis B virus is the major infectious occupational hazard of health workers, and most health care workers have received hepatitis B vaccine.

Hepatitis B virus is not spread by contaminated food or water, and cannot be spread casually in the workplace. High rates of chronic HBV infection are also found in the southern parts of Eastern and Central Europe, in the Middle East and Indian sub-continent, about 5% are chronically infected. Infection is less common in Western Europe and North America, where less than 1% are chronically infected.

Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 30% to 50% of children infected between 1 to 4 years of age develop chronic infection. The risk of death from HBV-related liver cancer or

cirrhosis is approximately 25% for persons who become chronically infected during childhood.

Chronic hepatitis B in some patients is treated with drugs called *interferon* or *lamivudine*, which can help some patients. Patients with cirrhosis are sometimes given liver transplants, with varying success. It is preferable to prevent this disease with vaccine than to try and cure it.

Hepatitis B vaccine has an outstanding record of safety and effectiveness. Since 1992, over one billion doses of hepatitis B vaccine have been used worldwide. The vaccine is given as a series of three intramuscular doses. Studies have shown that the vaccine is 95% effective in preventing children and adults from developing chronic infection if they have not yet been infected. In many countries where 8% to 15% of children used to become chronically infected with HBV, the rate of chronic infection has been reduced to less than 1% in immunized groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programmes."

Hepatitis B surface Antigen (HBsAg) is the major structural polypeptide of the envelope of the Hepatitis B Virus (HBV). This antigen is composed mainly of the type common determinant "s" and the type specific determinants "dr" and "y", present only on the specific serotypes.
Upon infection, a strong immunological response develops firstly against the type specific determinants and in a second time against the "s" determinant.

Anti "s" antibodies are however recognised to be most effective in the neutralisation of the virus, protecting the patient from other infections and leading it to convalescence.
The detection of HBsAb has become important for the follow up of patients infected by HBV and the monitoring of recipients upon vaccination with synthetic and natural HBsAg.

C. PRINCIPLE OF THE TEST

Microplates are coated with a preparation of highly purified HBsAg that in the first incubation with sample specifically captures anti HBsAg antibodies to the solid phase.
After washing, captured antibodies are detected by an HBsAg, labeled with peroxidase (HRP), that specifically binds the second available binding site of these antibodies.
The enzyme specifically bound to wells, by acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of HBsAb in the sample and can be detected by an ELISA reader.
The amount of antibodies may be quantified by means of a standard curve calibrated against the WHO reference preparation.
Samples are pre-treated in the well with an individual diluent able to block interference present in vaccinated individuals.

D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

1. Microplate: MICROPLATE

8x12 microwell strips coated with purified heat-inactivated HBsAg of both serotypes (ad and ay) from human origin and sealed into a bag with desiccant.
Allow the microplate to reach room temperature before opening. Reseal the microplate in the bag with desiccant and store at 4°C.

HBSAb

Enzyme Immunoassay for qualitative/quantitative determination of antibodies to Hepatitis B surface Antigen in human serum and plasma

- for "in vitro" diagnostic use only -

DIA.PRO
Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy
Phone: +39 02 27007161
Fax: +39 02 26007726
e-mail: info@diagpro.it

REFERENCES

1. Engvall E. et al., *J Immunochimistry*, 8, 871-874, 1971.
2. Engvall E. et al., *J Immunol*, 109, 129-135, 1971.
3. Remington J.S. and Klein J.O. In "Infectious diseases of the fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
4. Volk W.A. In "Essential of Medical Microbiology", 2nd ed., pp 729, G.B. Lippincott Company, Philadelphia, New York, S.J.oss, Toronto
5. Strydom D.R. et al., *AntiInMed*, 63 : 838, 1975.
6. Barser L.F., Dodd R.J., Sandler S.G. In "Viral Hepatitis: Laboratory and Clinical Science", F. Denhardt, J. Denhardt eds., M. Dekker Inc., New York, 215-230, 1983.
7. Cassatt T., *BritMed Bull.*, 28 : 156, 1972
8. Larber J.J. et al., *J Immunol*, 106 : 1056, 1971
9. Mushawar H.K. et al., *Ann J Clin Pathol*, 76 : 773, 1981.
10. Awood G.R., *Immunol Today*, 5 : 185, 1984
11. Aach H.D., *Lancet* 7874 : 190-193, 1974.
12. Jng W. et al., *J Hepatol*, 9 : 201-207, 1988
13. F. Covari et al., *Boll. Ist. Sieroter. Milan*, 63 : 14-18, 1984
14. Davidson et al., *J Natl Cancer Inst.*, 59 : 1451-1467, 1977
15. Espockey et al., *J Natl Cancer Inst.*, 59 : 1451-1467, 1977
16. S.H. Bell et al., *N.E.J. Med.*, 315 : 209-214, 1986
17. H. Hoonagale et al., *Hepatology*, 7 : 758-763, 1987
18. C. Johnson, *J Gen Virol*, 67 : 1215-1235
19. W. Jng et al., *J Hepatol*, 6 : 201-207, 1988
20. P. Richter et al., *Nephtologie*, 7 : 114-117, 1986
21. W. Stenius et al., *N.E.J. Med.*, 303 : 833-836, 1980
22. P. Toffler et al., *Nature*, 317 : 469-495, 1985
23. A.J. Zuckerman et al., in "Hepatitis Viruses of Man" Academic Press, London, 1979

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia, Pro Diagnostice Bioprobes SH
Via G. Carducci n° 27 - Sesto San Giovanni (MI) - Italy

HDV Ab

A. INTENDED USE
Competitive Enzyme Immunoassay (EUSA) for the qualitative determination of antibodies to Hepatitis Delta Virus or HDV in human plasma and sera with a "two-steps" methodology.

The kit is used for the follow-up of patients infected by HDV. For "in vitro" diagnostic use only.

B. INTRODUCTION

The Hepatitis Delta Virus or HDV is a RNA defective virus composed of a core presenting the delta-specific antigen, encapsulated by HBsAg, that requires the helper function of HBV to support its replication. Infection by HDV occurs in the presence of acute or chronic HBV infection. When acute delta and acute HBV simultaneously occur, the illness becomes severe and clinical and biochemical features may be indistinguishable from those of HBV infection alone. In contrast, a patient with chronic HBV infection can support HDV replication indefinitely, usually with a less severe illness appearing as a clinical exacerbation.

The determination of HDV specific serological markers (HDV Ag, HDV Ab, HDV IgM and HDV IgG) represents in these cases an important tool for the clinician for the classification of the etiological agent, for the follow up of infected patients and their treatment. The detection of HDV total antibodies allows the classification of the illness and the monitoring of the seroconversion event.

C. PRINCIPLE OF THE TEST

Anti-HDV antibodies, if present in the sample, compete with a virus-specific polyclonal IgG, labeled with peroxidase (HRP), for a fixed amount of rec-HDV coated on the microplate. The test is carried out with a two steps incubation competitive system. First the sample is added to the plate and specific anti HDV antibodies bind to the adsorbed antigen. After washing, an enzyme conjugated antibody to HDV is added and binds to the free portion of the antigen coated. After washing a chromogen/substrate mixture is dispensed. The concentration of the bound enzyme on the solid phase becomes inversely proportional to the amount of anti-HDV antibodies in the sample and its activity is detected by the added chromogen/substrate. The concentration of HDV-specific antibodies in the sample is determined by means of a cut-off value that allows for the semi quantitative detection of anti-HDV antibodies.

D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

1. Microplate: **MICROPLATE**

8x12 microwell strips coated with recombinant HDV-specific antigen and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Negative Control: **CONTROL -**

1x2 microwell. Ready to use. Contains goat serum proteins, 100 mM Tris-HCl buffer pH 7.4 +/-0.1, 0.3% Sodium Azide and 0.1% Kathon GC as preservatives. The negative control is colour coded pale yellow.

3. Positive Control: **CONTROL +**

1x2 microwell. Ready to use. Contains goat serum proteins, high titer anti HDV antibodies, 100 mM Tris-HCl buffer pH 7.4 +/-0.1, 0.09% Sodium Azide and 0.1% Kathon GC as preservatives. The positive control is colour coded green.

4. Calibrator: **CAL**
...
1 vial. Lyophilised. To be dissolved with EIA grade water as reported in the label. Contains bovine serum proteins, low titer human antibodies to HDV, 0.2 mg/ml gentamicin sulphate and 0.1% Kathon GC as preservatives.
Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

5. Wash buffer concentrate: **WASHBUF 20X**
1x60ml/bottle. 20X concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

6. Enzyme conjugate: **CONJ**
1x1 microwell. Ready-to-use solution. Contains 5% bovine serum albumin, 10 mM Tris buffer pH 6.8 +/-0.1, Horseradish peroxidase conjugated antibody to HDV in presence of 0.2 mg/ml gentamicin sulphate and 0.1% Kathon GC as preservatives. The component is colour coded red.

7. Chromogen/substrate: **SUBS TMB**
1x1 microwell. Contains a 50 mM citrate-phosphate buffered solution at pH 5.5-3.8, 4% DMSO, 0.02% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide of H₂O₂.
Note: To be stored protected from light as sensitive to strong illumination.

8. Sulphuric Acid: **H2SO4 0.31M**
1x1 microwell. Contains 0.3 M H₂SO₄ solution. Attention: Irritant (H315, H319, P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).
Plate sealers n° 2

Instructions for Use n° 1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes in the range 10-1000 µl and disposable plastic tips.
2. EIA grade water (double distilled or deionized, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minutes range or higher.
4. Absorbent paper tissues.
5. Calibrated EUSA microplate thermostat incubator (dry or wet) set at +37°C.
6. Calibrated EUSA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated EUSA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor, responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, face-free gloves and glasses. The use of any sharp (needles) or cutting blades devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication, "Biosafety in Microbiological and Biomedical Laboratories", ed 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.

DIA.PRO
Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy

Phone +39 02 27007161
Fax +39 02 26007726
e-mail: info@diagno.it

REF DAB CE
96 Tests

6. Pipette 100 µl TBH₂O₂ mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at room temperature for 20 min.

Important note: Do not expose to strong direct light as a high background might be generated.

7. Pipette 100 µl Sulphinic Acid into all the wells using the same pipetting sequence as in step n° 6 to stop the enzymatic reaction. Addition of the stop solution will turn the negative control and negative samples from blue to yellow.

8. Measure the colour intensity of the solution in each well, as described in section 1.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, strongly recommended), blanking the instrument on A1.

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self-oxidation of the chromogen can occur leading to a higher background.
- The use of the Calibrator, a low positive control, is not mandatory for the assay as the CAL does not enter into the cut-off calculation. The CAL may be used as a low tier positive control when a laboratory internal quality verification is required by the management. When used for such purpose, dispense 100 µl of it, possibly in duplicate.

Control/Calibrator	100 µl
Samples	100 µl
1st Incubation	60 min
Temperature	+37°C
Washing step	4-5 cycles
Enzyme Conjugate	100 µl
2nd Incubation	60 min
Temperature	+37°C
Washing step	4-5 cycles
TMB/H ₂ O ₂ mix	100 µl
3rd Incubation	20 min
Temperature	r.t.
Sulphinic Acid	100 µl
Reading OD	450nm & 620nm

An example of dispensation scheme (including CAL) is reported in the table below.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	S1	S9										

Legend: BLK = Blank NC = Negative Control
CAL = Calibrator PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A check is performed on the negative and positive controls any time, and on the Calibrator in addition when the kit OD_{450nm} or Co/S values have been matched in the analysis. Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.100 OD _{450nm} value
Negative Control (NC)	> 1.000 OD _{450nm} after blanking If lower carefully control the washing procedure and decrease the number of cycles or the soaking time
Positive Control (PC)	OD _{450nm} > 1.500 Coefficient of variation < 30%
Calibrator (CAL)	PC < OD _{450nm} < NC-PC/5

If the results of the test match the requirements stated above, proceed to the next section. If they don't, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD _{450nm}	1. that the Chromogen/Substrate solution has not been contaminated during the assay 2. that the washing procedure and the washer settings are as validated in the pre qualification after blanking
Negative Control (NC) < 1.000 OD _{450nm}	1. that the proper washing solution has been used and the washer has been primed with it before use; 2. that no mistake has been done in the assay (pipetting, dispensation, of positive control) 3. that no contamination of positive control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 4. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate; 5. that the washer needles are not blocked or partially obstructed.
Calibrator OD _{450nm} outside the range	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of Calibrator) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred
Positive Control OD _{450nm} > NC/10	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution (ex.: dispensation of negative control instead of Positive Control) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred

If any of the above problems have occurred, report the problem to the supervisor for further actions.

P. RESULTS

The results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-off} = (\text{NC} + \text{PC}) / 5$$

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station,

ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Results are interpreted as ratio between the cut-off value and the sample OD_{450nm} or Co/S. Results are interpreted according to the following table:

Co/S	Interpretation
< 0.9	Negative
0.9 – 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HDV.

Any patient showing an equivocal result should be re-tested on a second sample taken 1-2 weeks after the initial sample. A positive result is indicative of HDV infection and therefore the patient should be treated accordingly.

Important notes:

- Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
- When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
- Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

An example of calculation is reported below.

The following data must not be used instead of real figures obtained by the user.

Negative Control: 2.100 – 2.200 – 2.000 OD_{450nm}
Mean Value: 2.100 OD_{450nm}
Higher than 1.000 – Accepted

Positive Control: 0.100 OD_{450nm}
Lower than NC/10 – Accepted

Cut-Off = (2.100 + 0.100) / 5 = 0.440

Calibrator: 0.300-0.280 OD_{450nm}
Mean Value: 0.280 OD_{450nm}
Within the range PC ≤ OD_{450nm} < (NC+PC)/5 – Accepted

Sample 1: 0.020 OD_{450nm}
Sample 2: 1.900 OD_{450nm}
Sample 1 Co/S > 1.1 positive
Sample 2 Co/S < 0.9 negative

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

1. LIMIT OF DETECTION

In accordance of an international standard, the sensitivity of the assay has been calculated by means of the product named Accudata™ 121 supplied by Boston Biomedica Inc. – USA.

The table below reports the OD_{450nm} shown by this preparation when diluted in Fetal Calf Serum to prepare a limiting dilution curve, in three different lots.

DAB/CE Lot#	Lot# 102			Lot# 0101			Lot# 0403		
	Co/S	OD ₄₅₀ nm	Value	Co/S	OD ₄₅₀ nm	Value	Co/S	OD ₄₅₀ nm	Value
Accudata™ # 127	1x	0.171	3.0	0.153	2.9	0.156	2.8	0.197	2.7
	2x	0.197	2.7	0.176	2.6	0.179	2.5	0.200	2.4
	4x	0.200	2.2	0.220	2.1	0.202	2.2	0.258	1.7
	8x	0.298	1.7	0.285	1.6	0.271	1.6	0.412	1.2
	16x	0.412	1.2	0.405	1.1	0.402	1.1	0.450	0.9
	32x	0.450	0.9	0.482	0.8	0.482	0.8	0.550	0.6
	64x	0.550	0.6	0.582	0.5	0.582	0.5	0.650	0.4
	128x	0.650	0.4	0.682	0.3	0.682	0.3	0.750	0.2
CTRL (d)	2.464	0.00000	2.261	0.00000	2.114	0.00000			

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY

The diagnostic performances were evaluated in a clinical trial conducted by the Department of Gastro-hepatology, Prof. M. Rizzetto, S. Giovanni Battista Hospital, Torino, Italy, on more than 400 samples against a reference kit. Negative, positive and potentially interfering samples were examined in the trial. Both plasma, derived with different standard techniques of preparation (chrae, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed. Results are briefly reported in the tables below:

Sensitivity > 98 %
Specificity > 98 %

3. PRECISION
The mean values obtained from a study, conducted on two samples of different anti-HDV antibody reactivity, examined in 16 replicates in three separate runs for three lots of product, is reported below:

DAB/CE lot#1102				
Negative Control (N = 16)				
Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.242	2.428	2.433	2.401
Std Deviation	0.113	0.168	0.122	0.114
CV %	4.8	4.4	5.0	4.7
Calibrator (N = 16)				
Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.298	0.289	0.285	0.291
Std Deviation	0.023	0.027	0.025	0.025
CV %	7.7	9.3	9.1	8.7
Co/S				
Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.209	2.237	2.245	2.230
Std Deviation	0.105	0.108	0.108	0.107
CV %	4.7	4.8	4.8	4.8

DAB/CE lot#1103

Negative Control (N = 16)				
Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.209	2.237	2.245	2.230
Std Deviation	0.105	0.108	0.108	0.107
CV %	4.7	4.8	4.8	4.8

HBcAb

A. INTENDED USE
Competitive Enzyme Immunoassay (EUSA) for the determination of antibodies to Hepatitis B core Antigen in human plasma and sera. The kit is intended for the screening of blood units and the follow up of HBV-infected patients. For *in vitro* diagnostic use only.

B. INTRODUCTION
The World Health Organization (WHO) defines Hepatitis B as follows:

"Hepatitis B is one of the major diseases of mankind and is a serious global public health problem. Hepatitis means inflammation of the liver, and the most common cause is infection with one of 5 viruses, called hepatitis A, B, C, D, and E. All of these viruses can cause an acute disease with symptoms lasting several weeks, including yellowing of the skin and eyes (jaundice), dark urine, extreme fatigue, nausea, vomiting and abdominal pain. It can take several months to a year to feel fit again. Hepatitis B virus can cause chronic infection in which the patient never feels fit of the virus and many years later develops cirrhosis of the liver or liver cancer.

HBV is the most serious type of viral hepatitis and the only type causing chronic hepatitis for which a vaccine is available. Hepatitis B virus is transmitted by contact with blood or body fluids of an infected person in the same way as human immunodeficiency virus (HIV), the virus that causes AIDS. However, HBV is 50 to 100 times more infectious than HIV. The main ways of getting infected with HBV are (a) perinatal (from mother to baby at the birth); (b) child-to-child transmission; (c) unsafe injections and transfusions; (d) sexual contact.

Worldwide, most infections occur from infected mother to child, from child to child contact in household settings, and from reuse of unsterilized needles and syringes. In many developing countries, almost all children become infected with the virus. In many industrialized countries (e.g. Western Europe and North America), the pattern of transmission is different. In these countries, mother-to-infant and child-to-child transmission accounted for up to one third of chronic infections before childhood hepatitis B vaccination programmes were implemented. However, the majority of infections in these countries are acquired during young adulthood by sexual activity and injecting drug use. In addition, hepatitis B virus is the major infectious occupational hazard of health workers, and most health care workers have received hepatitis B vaccine.

Hepatitis B virus is not spread by contaminated food or water, and cannot be spread casually in the workplace. High rates of chronic HBV infection are also found in the southern parts of Eastern and Central Europe. In the Middle East and Indian subcontinent, about 5% are chronically infected. Infection is less common in Western Europe and North America, where less than 1% are chronically infected.

Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 30% to 50% of children infected between 1 to 4 years of age develop chronic

infection. The risk of death from HBV-related liver cancer or cirrhosis is approximately 25% for persons who become chronically infected during childhood.

Chronic hepatitis B in some patients is treated with drugs called *interferon* or *lamivudine*, which can help some patients. Patients with cirrhosis are sometimes given liver transplants, with varying success. It is preferable to prevent this disease with vaccine than to try and cure it.

Hepatitis B vaccine has an outstanding record of safety and effectiveness. Since 1982, over one billion doses of hepatitis B vaccine have been used worldwide. The vaccine is given as a series of three intramuscular doses. Studies have shown that the vaccine is 85% effective in preventing children and adults from developing chronic infection if they have not yet been infected. In many countries where 8% to 15% of children used to become chronically infected with HBV, the rate of chronic infection has been reduced to less than 1% in immunized groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programmes.

Hepatitis B core Antigen (or HBcAg) is the major component of the core particles of HBV. HBcAg is composed of a single polypeptide of about 17 kD that is released upon disaggregating the core particles; the antigen contains at least one immunodominant epitope.

Upon primary infection with HBV, HBcAg antibodies are one of the first markers of HBV hepatitis seen in the serum of the patient, slightly later than HBeAg, the serum e antigen. Anti-HBcAg antibodies are produced usually at 4-6 weeks and their presence is detectable even years after infection. Isolated HBcAb, in absence of other HBV markers, have been observed in infected blood units, suggesting the use of this test for screening HBV. In addition of HBeAg, HBcAb has become important for the classification of the viral agent, together with the detection of the other markers of HBV infection, in sera and plasma.

C. PRINCIPLE OF THE TEST

The assay is based on the principle of competition where the antibodies in the sample compete with a monoclonal antibody for a fixed amount of antigen on the solid phase.

A purified recombinant HBcAg is added to the microwells. The patient's serum/plasma is added to the microwell together with an additive able to block interferences present in the sample.

In the second incubation after washing, a monoclonal antibody, conjugated with Horseradish Peroxidase (HRP) and specific for HBcAg is added and binds to the free re-HBcAg coated on the plastic.

After incubation, microwells are washed to remove any unbound conjugate and then the chromogenic substrate is added in the presence of peroxidase enzyme. The colorless substrate is hydrolyzed to a colored end-product.

The color intensity is inversely proportional to the amount of antibodies to HBcAg present in the sample.

D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

1. Microplate MICROWELL

96 microwell strips coated with recombinant HBcAg and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; resal unused strips in the bag with desiccant and store at 2-8°C.

HBcAb

Competitive Enzyme Immunoassay for the determination of antibodies to Hepatitis B core Antigen in human serum and plasma

- for "in vitro" diagnostic use only -

DIA.PRO
Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy
Phone +39 02 27007161
Fax +39 02 26687226
e-mail: info@dia.pro.it

REF. BCAB CE
96 Tests

2. Negative Control **CONTROL**
1x1.0ml/vial. Ready to use. Contains 5% bovine serum albumin, 10 mM phosphate buffer pH 7.4 +/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The negative control is pale yellow color coded.

3. Positive Control **CONTROL**
1x1.0ml/vial. Ready to use. Contains 5% bovine serum albumin, anti-HBc antibodies at a concentration of about 10 PEI U/ml, calibrated on PEI HBc Reference Material 82), 10 mM phosphate buffer pH 7.4 +/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The positive control is green color coded.

4. Calibrator **CAL**
n°1 vial. Lyophilised. To be dissolved with EIA grade water as reported in the label. Contains fetal bovine serum, human antibodies to HBcAg at a concentration of 2 PEI U/ml +/-10% (calibrated on PEI HBc Reference Material 82) and 0.1% Kathon GC as preservative.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

5. Wash buffer concentrate **WASHBUF 20X**
1x60ml/vial. 20x concentrated solution.
Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0 +/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

6. Enzyme Conjugate **CONJ**
1x15ml/vial. Ready-to-use solution. Contains 5% bovine serum albumin, 10 mM Tris buffer, pH 6.8 +/-0.1, Horseradish peroxidase conjugated mouse monoclonal antibody to HBcAg in presence of 0.3 mg/ml gentamicin sulphate and 0.1% Kathon GC as preservatives. The component is red colour coded.

7. Chromogen/Substrate **SUBS TMB**
1x15ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.6 +/-0.1, 0.03% tetra-methyl-benzidine (TMB), 0.02% hydrogen peroxide (H₂O₂) and 4% dimethylsulphoxide **Note: To be stored protected from light as sensitive to strong illumination.**

8. Specimen Diluent **DILSPF**
4x3ml/vial. 10 mM Tris buffered solution pH 8.0 +/-0.1 containing 0.1% Kathon GC for the pre-treatment of samples and controls in the plate, blocking interference.
Note: Use all the content of one vial before opening a second one. The reagent is sensitive to oxidation.

9. Sulphuric Acid **H₂SO₄ 0.3 M**
1x15ml/vial. Contains 0.3 M H₂SO₄ solution.
Attention: Irritant (Xi) R36/38; S22/63/30

10. Plate sealing foil n°2

11. Instruction manual n°1

E MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, chemical treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELSA microplate thermostatic incubator (90° or wet) set at +37°C.
6. Calibrated ELSA microwell reader with 450nm (reading) and with 620-630nm (blinking) filters.
7. Calibrated ELSA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out the type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, face-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in disposal procedures as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
6. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
7. Do not interchange components between different lots of the kits. It is recommended that components between two lots of the same lot should not be interchanged.
8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
10. Avoid sample contamination between kit reagents by using disposable tips and changing them between the use of each one.
11. Do not use the kit after the expiration date stated on external (primary container) and internal (vial) labels.
12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
13. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 15-18 hrs or heat inactivation by autoclaving at 121°C for 20 min.
15. Accidental spills have to be absorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
16. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
17. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and

disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND RECOMMENDATIONS

1. Blood is drawn aseptically by venopuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples, especially sodium azide as this chemical would affect the enzymatic activity of the conjugate.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
4. Haemolysed (red) and visibly hyperfibrinemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2/-8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.8µm filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6-reusers of the device and up to 6 months.

1. Microplates:

Allow the microplate to reach room temperature (about 1 h) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storage.

In this case, call Dia.P.R.'s customer service. Unused strips have to be placed back inside the aluminum pouch after the desiccant supplied, firmly zipped and stored at +2-8°C. After next opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Positive Control:

Ready to use. Mix well on vortex before use.

4. Calibrator:

Add the volume of ELSA grade water, reported on the label, to the lyophilised powder; let fully dissolve and then gently mix on vortex. The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.

5. Wash buffer concentrate:

The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-to-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles. **Note: Once diluted, the wash solution is stable for 1 week at +2/-8°C.**

6. Enzyme conjugate:
Ready to use. Mix well on vortex before use.
Avoid contamination of the liquid with oxidizing chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

7. Chromogen/Substrate:

Ready to use. Mix well on vortex before use.
Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

8. Specimen Diluent

Ready to use solution. Mix gently on vortex before use. Use all the content of one vial before opening a second one. The reagent is sensitive to oxidation.

9. Sulphuric Acid:

Ready to use. Mix well on vortex before use.
Attention: Irritant (Xi) R36/38; S22/63/30
Legends: R36/38 = Irritating to eyes and skin.
S 22/63/30 = In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

1. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a linearity of ±2%.
2. The ELSA incubator has to be set at +37°C (tolerance of ±0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELSA tests.
3. The ELSA washer is extremely important to the overall performance of the assay. The washer must be carefully validated and correctly optimized using the kit for routine laboratory tests. Usually, 4-9 washing cycles (aspiration) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly the "rinse" it is recommended to obtain a clear and negative kit control/calibrator and well appeared and negative and positive reference samples, and check to match the values reported below in the "Reagents Validation of the test" and "Assay Performance". Regular calibration of the reader delivered and maintenance (cleaning, disinfection and clearing of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The ELSA microplate reader has to be equipped with a reading filter of 450nm, and with a second filter (620-630nm, strongly recommended) for blinking purposes. Its standard performance should be: (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to 2.0; (c) linearity to ≥ 20% repeatability ± 1%; (d) blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.

- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values performed in the sections 'Validation of Test' and 'Assay Performances'. The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and attention) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run. Diarrho's customer service offers support to the user in the setting and checking of instruments used in combination with the kit. In order to assure full compliance with the requirements described, support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

- Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
- Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen (TM3) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no leakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Describe the Calibrator as described above and gently mix.
- Allow all the error components to reach room temperature (about 17°C) and then mix gently on vortex all liquid components.
- Set the ELISA incubator at -37°C and prepare the ELISA washer by heating with the diluted washing solution, according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the microplates are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be performed according to the procedure given below, taking care to maintain the same incubation time for all the samples being tested.

- Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
- Leave the A1 well empty for blanking purposes.
- Dispense 50 µl Specimen Diluent into all the control and sample wells.
- Pipette 50 µl of the Negative Control in triplicate, 50 µl of the Calibrator in duplicate and then 50 µl of the Positive Control in single. Then dispense 50 µl of each of the samples.
- Incubate the microplate for 60 min at +37°C.

Important note: Strips have to be sealed with the adhesive sealing foil, only when test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section L.3)
- Pipette 100 µl Enzyme Conjugate in all the wells, except A1; incubate the microplate for 60 min at +37°C.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

- When the second incubation is finished, wash the microwells as previously described (section L.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1 included.

Important note: Do not expose to strong direct light, as a high background might be generated.

- Incubate the microplate protected from light at room temperature (18-24°C) for 20 minutes. Wells dispensed with negative control and negative samples will turn from clear to blue (competitive method).
- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9 to stop the enzymatic reaction. Addition of the stop solution will turn the negative control and negative samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading and recommended), blanking the instrument on A1.

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false results.
- Reading has to be performed immediately after the addition of the Stop solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.

N. ASSAY SCHEME

Specimen Diluent	50 µl
Control/Calibrator and samples	50 µl
1° Incubation Temperature	60 min +37°C
Wash	n° 4-5
Enzyme Conjugate	100 µl
2° Incubation Temperature	60 min +37°C
Wash	n° 4-5
TMB/H2O2 mix	100 µl
3° Incubation Temperature	20 min r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

An example of dispensation scheme is reported below:

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	S1	S9										

Legend: BLK = Blank
 NC = Negative Control
 CAL = Calibrator
 PC = Positive Control
 S = Sample

O. INTERNAL QUALITY CONTROL

A check is performed on the control/calibrator any time the kit is used in order to verify whether the expected OD450nm or CoS values have been matched in the analysis. Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.050 OD450nm value
Negative Control (NC)	> 1.000 OD450nm after blanking
Calibrator (about 2 PEI U/ml)	CoS > 1
Positive Control	< 0.200 OD450nm

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.050 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative control (NC) OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that no reagent is expired; 3. that no mistake has been done in the assay procedure (dispensation of positive control treated at 20°C)
Calibrator CoS < 1	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the coating (dispensation of negative control treated at 20°C); 3. that the procedure has been correctly performed; 4. that no external contamination of the positive control has occurred
Positive Control > 0.200 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the coating (dispensation of negative control treated at 20°C); 3. that the procedure has been correctly performed; 4. that no external contamination of the positive control has occurred

P. RESULTS
 The results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-Off} = (\text{NC} + \text{PC}) / 5$$

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formula is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Results are interpreted according to the following table:

CoS	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HBV.
 Any patient showing an equivocal result should be re-tested on a second sample taken 1-2 weeks after the initial sample.
 The blood unit should not be transfused.
 A positive result is indicative of HBV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

- Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
- When test results are transferred from the laboratory to a data handling system, attention must be paid to avoid erroneous data transfer of viral hepatitis infection has to be taken by and reported to the patient by a suitably qualified medical doctor.

An example of calculation is reported below.

The following data must not be used instead of real figures obtained by the user:

- Negative Control: 2.000 – 2.200 – 2.000 OD450nm
- Mean Value: 2.100 OD450nm
- Higher than 1.000 – Accepted
- Positive Control: 0.100 OD450nm
- Lower than 0.200 – Accepted
- Cut-Off = $(2.100 + 0.100) / 5 = 0.440$
- Calibrator: 0.400-0.360 OD450nm
- Mean value: 0.380 OD450nm
- CoS > 1 – Accepted
- Sample 1: 0.028 OD450nm
- Sample 2: 1.890 OD450nm
- Sample 1 CoS > 1.1 positive
- Sample 2 CoS < 0.9 negative

If any of the above problems have occurred, report the problem to the supervisor for further actions.

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

1. LIMIT OF DETECTION:

The sensitivity of the assay has been calculated by means of the reference preparation for HcAb supplied by Paul Ehrlich Institute (PEI Hbc Reference Material 82). The assay shows a sensitivity of about 1,25 PEI U/ml. The table below reports the CoV's values shown by a standard diluted as suggested by the manufacturer to prepare a limiting dilution curve in Feal Cell Serum (FCS).

PEI U/ml	Lot 1001	Lot 0702	Lot 0702/22	Lot 1202
9	22,6	18,0	19,0	17,7
2,5	8,0	5,5	5,4	5,0
1,25	1,1	1,3	1,0	1,0
0,625	0,4	0,4	0,4	0,4

In addition Accurun 1 – series 3000 – supplied by Boston Biomedica Inc, USA, was tested to determine its CoV's value. Results are reported in the table below:

Accurun 1 – series 3000

Value	Lot 1001	Lot 0702	Lot 1202
CoV's	2,9	2,3	2,2

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY

The Performance Evaluation of the device was carried out in a trial conducted on more than total 6000 samples.

2.1 Diagnostic Specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. A total of more 5000 unselected donors, including 1-time donors, were examined. In a first study 2023 samples were tested against a US company as reference. A specificity of 99,5% was found. In a second study 1588 samples were examined against a European company. A specificity of 99,7% was found. In the last study 1563 samples were assayed against the same US company; a value of 99,8% was found.

In addition to the above population, 206 samples from hospitalized patients were tested against the European company. A value of 99,3% specificity was found. Moreover, diagnostic specificity was assessed by testing 164 potentially interfering specimens (other infectious diseases, pregnant women, hemodialyzed, iatrogenic, etc.) against the European company. A value of specificity of 100% was assessed. Finally, both human plasma, derived with different standard techniques of preparation (Citrate, EDTA and heparin), and human sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

2.2 Diagnostic Sensitivity

It defines as the probability of the assay of scoring positive in the presence of specific analyte. 373 possible specimens were tested against the European company; a diagnostic sensitivity of 99,7% was found.

3. PRECISION

The mean values obtained from a study conducted on three lots and on two samples of different anti-HbcAg reactivity, examined in 16 replicates in three separate runs is reported below:

BCAB CE lot # 1202

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	1,943	1,930	1,924	1,932
Std.Deviation	0,081	0,076	0,102	0,087
CV %	4,2	4,0	5,3	4,5

Calibrator (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0,142	0,147	0,146	0,145
Std.Deviation	0,014	0,017	0,018	0,016
CV %	9,8	11,4	12,1	11,1
CoV's	2,8	2,7	2,8	2,7

BCAB CE lot # 0702

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	2,163	2,110	2,106	2,128
Std.Deviation	0,167	0,095	0,139	0,111
CV %	4,9	4,2	6,6	5,2

Calibrator (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0,192	0,193	0,195	0,190
Std.Deviation	0,018	0,023	0,019	0,020
CV %	10,0	12,0	9,9	10,6
CoV's	2,5	2,2	2,3	2,3

BCAB CE lot # 0702/22

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	2,278	2,298	2,130	2,169
Std.Deviation	0,135	0,128	0,158	0,140
CV %	5,9	6,0	7,5	6,5

Calibrator (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0,193	0,190	0,199	0,194
Std.Deviation	0,023	0,023	0,027	0,025
CV %	12,1	12,3	13,6	12,6
CoV's	2,4	2,2	2,2	2,3

The variability shown in the tables did not result in sample misclassification.

S. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with subsequent alteration of the level of the analyte. This test is suitable only for testing single samples and not pooled ones. Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, serology, as well as other diagnostic data should be considered.

REFERENCES

1. Aach R.D., Gisham J.W., Parker S.W., Prociak/Acad.Sci.USA, 68: 1956, 1971.
2. Blumberg B.S., Srinak A.I., London W.T., Hepatitis and leukemia: their relation to Australia antigen. Bull.N.Y.Acad.Med., 44:1565, 1968.
3. Bonello A., Devo M., Mantega R., J Immunol.Meth., 49:1, 1982.
4. Caldwell C.W., Barpet J.T., Clin.Chim.Acta 81: 305, 1977.
5. Frazee S., De S.Geron, Schiedegger D., J Immunol.Meth., 53: 1, 1980.
6. Resnik N.V., et al., Volk.Sang., 39:61, 1990.
7. Rock G.A.V., Leprosy, 52:281, 1981.
8. Schroder J.D. Med.Sid., 36: 281, 1981.
9. Almeida J.D. et al., Lancet, II: 1225, 1971.
10. Hoornagle J.H. et al., Lancet, II: 899, 1973.
11. Hoenigle J.H. et al., NE.J.Med., 290: 1936, 1974.
12. Katsnelki J.H. et al., Clin.Pain., 31:537, 1978.
13. Szymanski J., et al., Ann.Epidemiol., 104: 236, 1976.
14. Gishman B., et al., J Immunol. Methods, 15(2):219-231, 2002.
15. Schriener RS and Kramps JA, Rev.Sci.Tech., 17(2):550-561, 1998.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Produced by
Dia.Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

HBe Ag&Ab

A. INTENDED USE
Enzyme immunoassay (ELISA) for the determination of Hepatitis B Virus e Antigen and Antibody in human plasma and sera. The kit is intended for the follow-up of acute infection and of chronic patients under therapy.
For *in vitro* diagnostic use only.

B. INTRODUCTION

Hepatitis B Virus (HBV) Antigen of HBeAg is known to be intimately associated with Hepatitis B Virus of HBV replication and the presence of infectious Dane particles in the blood. Recently, it has been found that HBeAg is a product of integrated degradation of Hepatitis B core Antigen or HBeAg, occurring in hepatocytes, whose expression is under the control of the promoter region of HBV genome.

If HBeAg is considered a specific marker of infectivity, the presence of anti HBeAg antibodies in blood is recognised to be a clinical sign of recovery from infection to convalescence. The determination of these two analyses in samples from HBV patients has become important for the diagnosis of the phase of illness and as a prognostic value in the follow up of infected patients.

C. PRINCIPLE OF THE TEST

HBeAg: If present in the sample, it is captured by a specific monoclonal antibody, in the 1st incubation.

In the 2nd incubation, after washing, a tracer, composed of a mix of two specific anti HBeAg monoclonal antibodies, labeled with peroxidase (HRP), is added to the microplate and binds to the captured HBeAg. The concentration of the bound enzyme on the solid phase is proportional to the amount of HBeAg in the sample and its activity is detected by adding the chromogen/substrate in the 3rd incubation. The presence of HBeAg in the sample is determined by means of a cut-off value that allows for the semiquantitative detection of the antigen.

HBeAb

Anti HBeAg antibodies, if present in the sample, compete with a recombinant HBeAg preparation for a fixed amount of an anti HBeAg antibody, coated on the microplate wells.

The competitive assay is carried out in two incubations, the first with the sample and reHBeAg, and the second with a tracer, composed of two anti HBeAg monoclonal antibodies, labeled with peroxidase (HRP). The concentration of the bound enzyme on the solid phase becomes inversely proportional to the amount of anti HBeAg antibodies in the sample and its activity is detected by adding the chromogen/substrate in the third incubation.

The concentration of HBeAg specific antibodies in the sample is determined by means of a cut-off value that allows for the semi quantitative detection of anti HBeAg antibodies.

D. COMPONENTS

The kit contains reagents for total 96 tests.

1. Microplate: [MICROPLATE]

n° 1 coated microplate
12 strips of 8 breakable wells coated with anti HBeAg specific monoclonal antibody, postcoated with bovine serum proteins and sealed into a bag with desiccant. Follow the microplate to reach room temperature before opening. (seal unused strips in the bag with desiccant and store at 2,8°C).

2. Negative Control: [CONTROL -]

1x20µl/vial. Ready to use control. It contains bovine serum, 0,09% sodium azide and 0,1% Kathon GC as preservatives.
The negative control is colorless.

3. Antigen Positive Control: [CONTROL + Ag]

1x10µl/vial. Ready to use control. It contains 2% bovine serum albumin, non infectious recombinant HBeAg, 100 mM Tris buffer pH 7.4+/-0.1, 0,09% sodium azide and 0,1% Kathon GC as preservatives.
The positive control is green color coded.

4. Antibody Positive Control: [CONTROL + Ab]
1x10µl/vial. Ready to use control. It contains 2% bovine serum albumin, human anti HBeAg positive plasma at about 10 PEI U/ml, 100 mM Tris buffer, pH 7.4+/-0.1, 0,09% sodium azide and 0,1% Kathon GC as preservatives. The label is red colored.
The positive control is yellow color coded.

5. Antigen Calibrator: [CALIB. ...m]

n° 1 vial. Lyophilized calibrator for HBeAg. To be dissolved with EIA grade water as reported in the label. It contains fetal bovine serum, non infectious recombinant HBeAg at 1 PEI U/ml +/-10%, 0,02% gentamicine sulphate and 0,1% Kathon GC as preservatives.

Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

6. Antibody Calibrator: [CALIB. ...m]

n° 1 vial. Lyophilized calibrator for anti HBeAg antibody. To be dissolved with EIA grade water as reported in the label. It contains fetal bovine serum, positive plasma at 0,25 PEI U/ml +/-10%, 0,02% gentamicine sulphate and 0,1% Kathon GC as preservatives. The label is red colored.
Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

7. Wash buffer concentrate: [WASH-BUF 20X]

1x60ml/bottle. 20x concentrated solution.
Once diluted, the wash solution contains 10 mM phosphate buffer pH 7,0+/-0,2, 0,05% Tween 20 and 0,1% Kathon GC.

8. Enzyme conjugate: [CONJ]

1x10µl/vial. Ready to use conjugate. It contains horseradish peroxidase conjugated with a mix of monoclonal antibodies to HBeAg, 10 mM Tris buffer pH 6,8+/-0,1, 2% BSA, 0,1% Kathon GC and 0,02% gentamicine sulphate as preservatives.
The reagent is red color coded.

9. HBe Antigen: [Ag-HBe]

1x10µl/vial. Ready to use reagent. It contains recombinant HBeAg, fetal bovine serum, buffered solution pH 8,0+/-0,1, 0,1% Kathon GC and 0,09% sodium azide as preservatives.
The reagent is blue color coded.

10. Chromogen/substrate: [SUBST TMB]

1x10µl/vial. Ready-to-use component. It contains a 50 mM citrate-phosphate buffered solution at pH 3,5-3,8, 4% dimethylsulphoxide, 0,03% tetra-methyl-benzidine or TMB and 0,02% hydrogen peroxide or H₂O₂.
Note: To be stored protected from light as sensitive to strong illumination.

11. Sulphuric Acid: [H₂SO₄ 0,3 M]

1x10µl/vial. It contains 0,3 M H₂SO₄ solution.
Attention: Irritant (H315, H319). P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

12. Plate sealing foil: n° 2

13. Package insert: n° 1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (15µl, 100µl and 500µl) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermosatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (tearing) and with 620-630nm (blinking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

HBe Ag&Ab

Enzyme Immunoassay (ELISA) for the determination of Hepatitis B Virus "e" Antigen and Antibody in human plasma and sera.

- for "in vitro" diagnostic use only -

DIA.PRO
Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy
Phone +39 02 27007161
Fax +39 02 26007726
e-mail: info@dia.pro.it

REF. HBE CE
96 Tests

F. WARNINGS AND PRECAUTIONS

- The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
- All the personnel involved in performing the assay have to wear protective laboratory clothes, face-free gloves and glasses. The use of any sharp (needles) or cutting (bladed) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
- The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the ChromogenSubstrate (MSB) from strong light and avoid vibration of the bench surface where the test is undertaken.
- Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
- Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
- Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
- Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
- Avoid cross-contamination between kit reagents by using disposable (primary container) and internal (vials) labels.
- Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what is reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
- Vials produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste generated from the original substances. In particular, liquid waste generated from samples has to be treated as potentially infective material and inactivated. Sterilizer or household bleach treatment with a 10% final concentration of household bleach for 15-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
- Accidental spills have to be decontaminated with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
- The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water.
- Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND RECOMMENDATIONS

- Blood is drawn aseptically, by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- Avoid any addition of the enzymatic activity of the conjugate, generating false negative results.
- Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results.
- Hemolysed and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- Sera and plasma can be stored at +2...-8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be freeze/thawed more than once as this may generate particles that could affect the test result.

H. PREPARATION OF COMPONENTS AND WARNINGS

- If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µm filters to clean up the sample for testing.
- A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.
- Microplate:** Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing. In this case, call Dia Pro's customer service. Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2...-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.
- Negative Control:** Ready to use. Mix well on vortex before use.
- Antigen Positive Control:** Ready to use. Mix well on vortex before use.
- Antibody Positive Control:** Ready to use. Mix well on vortex before use.
- Antigen Calibrator:** Add the volume of ELISA grade water, reported on the label, to the lyophilized powder, let fully dissolve and then gently mix on vortex. **Note: The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.**
- Antibody Calibrator:** Add the volume of ELISA grade water, reported on the label, to the lyophilized powder, let fully dissolve and then gently mix on vortex. **Note: The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.**
- Wash buffer concentrate:** The whole content of the 20x concentrated solution has to be diluted with distilled water up to 200 ml and mixed gently end-over-end before using. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles. **Note: Once diluted, the wash solution is stable for 1 week at +2...8°C.**
- Enzyme conjugate:** Ready to use. Mix well on vortex before use. Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.
- HBs Antigen:** Ready to use. Mix well on vortex before use. Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.
- ChromogenSubstrate:** Ready to use. Mix well on vortex before use. Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.
- Sulphuric Acid:** Ready to use. Mix well on vortex before use. Antigen: H115, H131B, P305, P306, P331+P333, P331+P313, P302+P352, P332+P313, P305+P351+P338, P331+P313, P302+P353.

Legend:
 Warning H statements:
 H315 – Causes skin irritation.
 H319 – Causes serious eye irritation.

I. PRECAUTIONARY P STATEMENTS:

P201 – Wear protective gloves/protective clothing/eye protection/face protection.

- P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.
- P332 + P313 – IF skin irritation occurs: Get medical attention.
- P305 + P351 + P338 – IF IN EYES: Rinse carefully with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P337 + P313 – If eye irritation persists: Get medical advice/attention.
- P362 + P363 – Take off contaminated clothing and wash it before reuse.

L. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

- Microplates have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a linearity of +2%.
- The ELISA incubator has to be set at +37°C (tolerance of +1-0,5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests: 4-5 washing cycles (aspiration + dispensation of 500µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
- Incubation times have a tolerance of ±5%.
- The ELISA reader has to be equipped with a reading filter of 450nm and with a second filter (520-650nm, strongly recommended) for blanking purposes. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceeds 20-50 runs per run.

Legend:
 Warning H statements:
 H315 – Causes skin irritation.
 H319 – Causes serious eye irritation.

M. ASSAY PROCEDURE

- Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
 - Leave the A1 well empty for blanking purposes.
 - Pipette 100 µl of the Negative Control, in triplicate, 100 µl of the Antigen Calibrator in duplicate and then 100 µl of the Antigen Positive Control in single.
 - Then dispense 100 µl of samples in the proper wells.
 - Check for the presence of samples in wells by naked eye (there is a marked colour difference between empty and full wells) or by reading at 450/620nm (samples show OD values higher than 0.100).
 - Incubate the microplate for 60 min at +37°C.
- Important note:** Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.
- When the first incubation is finished, wash the microwells as previously described (section I.3)
 - Dispense 100 µl Enzyme Conjugate in all wells, except for A1, used for blanking operations.
- Important note:** Be careful not to touch the inner surface of the well with the pipette tip and not to immerse the top of it into samples or controls. Contamination might occur.
- Check that the reagent has been dispensed properly and then incubate the microplate for 60 min at +37°C.
 - When the second incubation is finished, wash the microwells as previously described (section I.3)
 - Pipette 100 µl ChromogenSubstrate into all the wells. A1 included.
- Important note:** Do not expose to strong direct light as a high background might be generated.
- Incubate the microplate protected from light at room temperature (16-24°C) for 20 minutes. Wells dispensed with positive control and positive samples will turn from clear to blue.
 - Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 11. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.
 - Measure the color intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-650nm filter (background) subtraction, strongly recommended). Blanking the instrument on A1.

Legend:
 Warning H statements:
 H315 – Causes skin irritation.
 H319 – Causes serious eye irritation.

- Pipette 50 µl of the Negative Control in triplicate, 50 µl of the Antibody Calibrator in duplicate and then 50 µl of the Antibody Control in single.
- Then dispense 50 µl of samples in the proper wells.
- Check for the presence of samples in wells by naked eye (there is a marked color difference between empty and full wells) or by reading at 450/620nm (samples show OD values higher than 0,100).
- Dispense then 50 µl of HBe Antigen in all the wells, except for A1.
- Incubate the microplate for 60 min at +37°C.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section I.3)
- Finally proceed as described for the HBeAg assay from point 8 to the last one.

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.
- The Calibrator (CAL) does not affect the cut-off calculation and therefore the test results calculation. The Calibrator may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

HBe antigen test

Controls and calibrator	100 µl
Samples	100 µl
1st Incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme Conjugate	100 µl
2nd Incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
1MBH₂O₂ mix	100 µl
3rd Incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

Hbe antibody test

Controls and calibrator	50 µl
Samples	50 µl
Neutralising antigen	50 µl
1st Incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzymatic conjugate	100 µl
2nd Incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
1MBH₂O₂ mixture	100 µl
3rd Incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

An example of dispensation scheme is reported below:

	Microplate											
	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	PC	S9										
I	S1	S9										

Legend: BLK = Blank // NC = Negative Control

PC = Positive Control // CAL = Calibrators // S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as qualified.

Control that the following data are matched:

HBe Antigen	
Blank well	OD450nm
Check	< 0,100 OD450nm
Negative Control (NC)	< 0,150 OD450nm after blanking
Antigen Calibrator	S/Co > 2,0
Positive Control (PC)	> 1,500 OD450nm
HBe Antibody	
Blank well	OD450nm
Check	> 1,000 OD450nm
Negative Control (NC)	> 1,000 OD450nm after blanking
Antibody Calibrator	OD450nm < NC/1,5
Positive Control (PC)	OD450nm < NC/10

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, don't proceed any further and perform the following checks:

HBeAg	Check
Blank well	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Positive Control	2. that the washing procedure and the washer settings are as validated in the pre qualification study;
Negative Control	3. that no mistake has been done in the assay before use;
Antibody Calibrator	4. that no contamination of the positive control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate;
Enzyme Conjugate	5. that microplates have not become contaminated with positive samples or with the enzyme conjugate needles are not blocked or partially obstructed.

Calibrator	S/Co < 2
1. that the procedure has been correctly performed;	
2. that no mistake has occurred during the distribution (ex.: dispersion of negative control material);	
3. that the washing procedure and the washer settings are as validated in the pre qualification study;	
4. that no external contamination of the calibrator has occurred.	
Positive Control	< 1,500 OD450nm
1. that the procedure has been correctly performed; needles has occurred during the distribution of the control (dispersion of negative control instead of positive control);	
2. that the washing procedure and the washer settings are as validated in the pre qualification study;	
4. that no external contamination of the positive control has occurred.	

Hbe antibody

Blank well	OD450nm
Check	> 1,000 OD450nm
Negative Control (NC)	> 1,000 OD450nm after blanking
Antibody Calibrator	OD450nm < NC/1,5
Positive Control (PC)	OD450nm < NC/10

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, don't proceed any further and perform the following checks:

Calibrator	OD450nm > NC/1,5
1. that the procedure has been correctly performed;	
2. that no mistake has occurred during its distribution (ex.: dispersion of negative control material);	
3. that the washing procedure and the washer settings are as validated in the pre qualification study;	
4. that no external contamination of the calibrator has occurred.	
Positive Control	> 1,500 OD450nm
1. that the procedure has been correctly performed;	
2. that no mistake has occurred during the distribution of the control;	
3. that the washing procedure and the washer settings are as validated in the pre qualification study;	
4. that no external contamination of the positive control has occurred.	

If any of the above problems have occurred, report the problem to the supervisor for further actions.

P. CALCULATION OF THE CUT-OFF

The results are calculated by means of a cut-off value determined with the following formula:

$$NC + 0,100 = \text{Cut-Off (Co)}$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

HBeAg:

$$(NC + PC) / 3 = \text{Cut-Off (Co)}$$

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Results are interpreted as follows:

S/Co	Interpretation
< 0,9	Negative
0,9 - 1,1	Equivocal
> 1,1	Positive

HBeAg:

Co/S	Interpretation
< 0,9	Negative
0,9 - 1,1	Equivocal
> 1,1	Positive

HBeAb:

Note:
S = OD450nm of the sample
Co = cut-off value

An example of calculation for HBeAg assay is reported below.

The following data must not be used instead or real figures obtained by the user.

Negative Control:	0,020 - 0,030 - 0,025 OD450nm
Mean Value:	0,025 OD450nm
Lower than 6,153 - Accepted	
Positive Control:	2,489 OD450nm
Higher than 1,500 - Accepted	
Cut-Off = 0,025 + 0,100 = 0,125	
Calibrator:	0,520 - 0,540 OD450nm
Mean value:	0,530 OD450nm
S/Co higher than 2,0 - Accepted	S/Co = 4,2
Sample 1:	0,030 OD450nm
Sample 2:	1,800 OD450nm
Sample 1 S/Co < 0,9 = negative	
Sample 2 S/Co > 1,1 = positive	

An example of calculation for HBeAb is reported below.

The following data must not be used instead or real figures obtained by the user.

Negative Control:	2,100 - 2,200 - 2,000 OD450nm
Mean Value:	2,100 OD450nm
Higher than 1,000 - Accepted	
Positive Control:	0,100 OD450nm
Lower than NC/10 - Accepted	
Cut-Off = (2,100 + 0,100) / 3 = 0,733	
Calibrator:	0,720 - 0,760 OD450nm
Mean value:	0,740 OD450nm
OD450nm < NC/1,5 - Accepted	
Sample 1:	0,020 OD450nm
Sample 2:	1,900 OD450nm
Sample 1 Co/S > 1,1 = positive	
Sample 2 Co/S < 0,9 = negative	

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory director to reduce the risk of judgment errors and misinterpretations.
2. The identification of the clinical status of a HBV patient (acute, chronic, asymptomatic hepatitis) has to be done on the basis also of the other markers of HBV infection (HBeAg, HBeSAb, HBsAb, HBeIgm).
3. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCE CHARACTERISTICS

A) HBeAg

1. Limit of detection
The limit of detection of the assay has been calculated by means of the International Standard for HBeAg, supplied by Paul Elich Institute (PEI). The data obtained by examining the limit of detection on three lots is reported in the table below.

LotID	HBE/CE	PEI U/ml
0103	0.25	
0103Z2	0.25	
0303	0.25	

In addition the preparation Accunum # 51, produced by Boston Biomedica Inc., USA, has been tested, upon dilution in FCS. Results are reported for three lots of products.

BBI's Accunum 51 (S/Ce)

LotID	1 x	2 x	4 x	8 x	15x
HBE/CE	4.1	1.6	0.9	0.6	0.4
0103	4.1	1.7	0.9	0.6	0.4
0103Z2	4.0	1.6	0.9	0.5	0.4

2. Diagnostic Sensitivity

The diagnostic sensitivity has been tested on panels of samples classified positive by a US FDA approved kit. Positive samples were collected from different HBV pathologies (acute, chronic), bearing HBeAg reactivity. An overall value > 98% has been found in the study conducted on a total number of more than 200 samples. Moreover the Panel of Seroconversion code PHM 935B, produced by BBI, was examined. Data are reported below and compared with those reported by BBI for two other commercial products:

Sample ID	HBE/CE S/Ce	Abbott EIA S/Ce	Serim EIA S/Ce
21	0.2	0.2	0.2
22	3.7	4.3	6.3
23	1.9	3.2	3.1
24	1.1	2.4	1.5
25	1.0	2.1	1.2
26	0.6	1.7	0.7
27	0.2	0.8	0.3
28	0.2	0.6	0.2
29	0.2	0.4	0.2
30	0.2	0.3	0.2
31	0.1	0.3	0.2
32	0.1	0.3	0.2

Finally the Performance Panel code PHU 201, produced by BBI, was tested. Data are reported below and compared with those reported by BBI for an other commercial product.

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.555	0.573	0.562	0.563
Std Deviation	0.026	0.025	0.024	0.025
CV %	4.7	4.3	4.2	4.4
S/Ce	4.2	4.4	4.4	4.3

Member	PEI U/ml	HBE/CE	Serim EIA
1	3	3.3	7.0
2	6	17.5	21.8
3	26	30.1	37.1
4	31	29.4	23.5
5	1	1.1	2.2
6	2	2.3	6.0
7	35	35.1	24.6
8	38	29.2	31.9
9	4	16.6	10.8
10	1	3.3	0.2
11	1	0.2	1.2
12	<1	0.2	1.2
13	<1	0.9	1.4
14	<1	0.2	0.2
15	<1	0.4	0.1
16	<1	0.5	0.1
17	<1	0.3	0.2
18	<1	0.2	0.2
19	<1	0.2	0.1
20	<1	0.2	0.1
21	<1	0.3	1.0
22	<1	0.3	0.1
23	<1	0.2	0.2
24	<1	0.2	0.2
25	<1	0.3	0.2

HBE/CE lot # 0303

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.032	0.032	0.033	0.032
Std Deviation	0.002	0.002	0.003	0.002
CV %	7.4	6.2	7.3	7.8

3. Diagnostic Specificity

The diagnostic specificity has been determined on panels of negative samples from normal individuals and blood donors, classified negative with a FDA approved kit. Both plasma, derived with different standard techniques of preparation (citrate, EDTA, and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed. Frozen specimens have also been tested to check whether this interferes with the performance of the test. No interference was observed on clean and particle free samples. Samples derived from patients with different viral (HCV and HAV) and non viral pathologies of the liver that may interfere with the test were examined. No cross reaction was observed.

The Performance Evaluation study conducted in a qualified external reference center on more than 500 samples has provided a value > 98%.

4. Precision

It has been calculated on two samples examined in 16 replicate in three different runs on three lots. The values found were as follows: HBE/CE lot # 0103

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.032	0.032	0.033	0.032
Std Deviation	0.002	0.002	0.003	0.002
CV %	7.4	6.2	7.3	7.8

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.556	0.575	0.575	0.568
Std Deviation	0.027	0.029	0.028	0.028
CV %	4.7	5.3	4.9	4.9
S/Ce	4.4	4.4	4.4	4.4

HBE/CE lot # 0103Z2

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.033	0.031	0.030	0.032
Std Deviation	0.003	0.003	0.002	0.003
CV %	7.9	8.5	7.4	8.0

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.555	0.573	0.562	0.563
Std Deviation	0.026	0.025	0.024	0.025
CV %	4.7	4.3	4.2	4.4
S/Ce	4.2	4.4	4.4	4.3

HBE/CE lot # 0303

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.033	0.034	0.036	0.034
Std Deviation	0.003	0.003	0.004	0.003
CV %	9.7	9.8	9.2	9.8

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.573	0.572	0.564	0.569
Std Deviation	0.023	0.026	0.025	0.024
CV %	4.1	4.8	4.5	4.5
S/Ce	4.5	4.3	4.1	4.3

B) HBe Antibody

1. Limit of detection
The limit of detection of the assay has been calculated by means of the International Standard for HBeAb, supplied by Paul Elich Institute (PEI). The data obtained by examining the limit of detection on three lots is reported in the table below.

LotID	HBE/CE	PEI U/ml
0103	0.25	
0103Z2	0.25	
0303	0.25	

In addition the preparation Accunum # 52, produced by Boston Biomedica Inc., USA, has been tested, upon dilution in FCS. Results are reported for three lots of products.

Accunum 52 (CvS)

LotID	1 x	2 x	4 x	8 x	15x
HBE/CE	1.0	0.8	0.5	0.4	0.4
0103	1.0	0.8	0.6	0.5	0.4
0103Z2	1.0	0.8	0.6	0.5	0.4
0303	1.0	0.8	0.5	0.4	0.4

2. Diagnostic sensitivity

The diagnostic sensitivity has been tested on panels of samples classified positive for HBeAb by a US FDA approved kit. Positive samples were collected from different HBV pathologies bearing anti-HBeAg antibody reactivity. An overall value > 98% has been found in the study conducted on a total number of more than 200 samples. Moreover the Panel of Seroconversion code PHM 935B, produced by BBI, was examined. Data are reported below and compared with those reported by BBI for two other commercial products.

Sample ID	HBE/CE C/S	Abbott EIA C/S	Serim EIA C/S
21	0.4	0.4	0.5
22	0.4	0.4	0.5
23	0.4	0.5	0.5
24	0.4	0.5	0.6
25	0.4	0.5	0.5
26	0.5	0.6	0.6
27	0.5	0.6	0.6
28	0.7	0.8	0.7
29	0.6	0.9	0.7
30	0.6	1.0	0.9
31	1.0	1.3	1.1
32	1.0	1.2	1.0

Finally the Performance Panel code PHU 201, produced by BBI, was tested. Data are reported below and compared with those reported by BBI for an other commercial product.

PEI 0.25 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.560	0.561	0.576	0.564
Std Deviation	0.036	0.037	0.046	0.039
CV %	6.5	6.6	8.0	7.3
C/S	0.9	1.0	1.0	1.0

Member	PEI U/ml	HBE/CE	Serim EIA
1	-	0.3	0.5
2	-	0.2	0.6
3	-	0.2	0.4
4	-	0.2	0.4
5	-	0.2	0.5
6	-	0.2	0.5
7	-	1.9	0.6
8	-	1.8	0.6
9	-	0.3	0.5
10	-	0.4	0.9
11	-	0.4	0.9
12	-	0.4	0.9
13	2	4.4	9.1
14	1	3.8	2.9
15	<1	1.0	1.5
16	>50	4.3	120.9
17	<1	1.0	1.0
18	5	5.9	21.9
19	1	5.0	4.5
20	11	7.0	4.5
21	2	1.9	10.0
22	26	28.1	90.7
23	<1	0.3	0.5
24	<1	0.8	1.3
25	50	28.1	167.4

HBE/CE lot # 0303

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	2.444	2.420	2.471	2.458
Std Deviation	0.129	0.150	0.142	0.144
CV %	5.2	6.6	5.7	5.9

3. Diagnostic specificity

The clinical specificity has been determined as described before for HBeAg. The Performance Evaluation study conducted in a qualified external reference center on more than 500 samples has provided a value > 98%.

4. Precision

It has been calculated on two samples examined in 16 replicate in three different runs on three lots. The values found were as follows:

HBE/CE lot # 0103

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.601	0.590	0.778	0.649
Std Deviation	0.060	0.060	0.050	0.051
CV %	1.0	1.0	6.7	1.9

PEI 0.25 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.561	0.560	0.578	0.564
Std Deviation	0.036	0.037	0.046	0.039
CV %	6.5	6.6	8.0	7.3
C/S	1.0	1.0	1.0	1.0

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	2.315	2.351	2.413	2.363
Std Deviation	0.127	0.144	0.148	0.139
CV %	5.5	6.1	6.0	5.9

PEI 0.25 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.567	0.579	0.575	0.571
Std Deviation	0.041	0.050	0.046	0.046
CV %	5.4	6.3	5.8	5.8
C/S	1.0	1.0	1.0	1.0

HBE/CE lot # 0303Z2

Negative Control (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	2.334	2.415	2.437	2.395
Std Deviation	0.146	0.155	0.158	0.153
CV %	6.3	6.4	6.5	6.4

PEI 0.25 U/ml (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.560	0.561	0.576	0.564
Std Deviation	0.036	0.037	0.046	0.039
CV %	6.5	6.6	8.0	7.3
C/S	0.9	1.0	1.0	1.0

5. LIMITATIONS

Frozen samples containing fibrin particles or aggregates may generate false positive results.

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.

This test is suitable only for testing single samples and not pooled ones. Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.

REFERENCES

1. Engvall E. and Perlmann P., *J. Immunochimistry*, 8, 871-874, 1971
2. Engvall E. and Perlmann P., *J. Immunol.*, 109, 129-135, 1971
3. Remington J.S. and Klein J.O., In "Infectious diseases of the fetus and newborn infant", Sanders, Philadelphia, London, Toronto.
4. Volk W.A., In "Essential of Medical Microbiology", 2nd ed, pp 729, G.B.Lippincott Company, Philadelphia, New York, St. José, Toronto.
5. Strydom D.R., Bryan J.A. and Dixon R.E., *Ann.Intl.Med.*, 83, pp 838, 1975.
6. Barker L.F., Garely R.J., Lorenz D.E., *Viral Hepatitis*, 581-587, 1978.
7. Cossart Y., *Brd.Med.Bull.*, 28, pp 156, 1972
8. Lander J.J., Alter H. and Purcell R., *J. Immunol.*, 106, pp 1086, 1971
9. Mushawar I.K., Dienstag J.L., Polesky H.F., et al., *Ann.J.Clin.Pathol.*, 76, pp 773, 1981.
10. Ling C.M., Mushawar I.K., et al., *Infection and Immunity*, 24: 235, 1979.
11. Mushawar I.K., Overby L.R. et al., *J. Med. Virol.*, 2: 77, 1978
12. Adeniyi J., Prosser G.G. et al., *J. Med. Dis.*, 141: 293, 1980
13. Magnus L.O., Linnam A. et al., *J. Am. Med. Assoc.*, 231: 566, 1975
14. Krogman S., Overby L.R. et al., *NIEngl.J. Med.*, 300: 101, 1979

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

0318

HCV AB

HCV AB

**Version 4.0 Enzyme Immunoassay
for the determination of
anti Hepatitis C Virus antibody
in human serum and plasma**

- for "in vitro" diagnostic use only -

A. INTENDED USE
Version 4.0 Enzyme Immunoassay (EISA) for the detection of antibodies to Hepatitis C Virus in human plasma and sera. It is intended for the screening of blood units and the follow-up of HCV-infected patients. For "in vitro" diagnostic use only.

B. INTRODUCTION
The World Health Organization (WHO) define Hepatitis C infection as follows:

"Hepatitis C is a viral infection of the liver which had been referred to as parenterally transmitted "non A, non B hepatitis" until identification of the causative agent in 1989. The discovery and characterization of the hepatitis C virus (HCV) led to the understanding of its primary role in post-transfusion hepatitis and its tendency to induce persistent infection.

HCV is a major cause of acute hepatitis and chronic liver disease, including cirrhosis and liver cancer. Globally, an estimated 170 million persons are chronically infected with HCV and 3 to 4 million persons are newly infected each year. HCV is spread primarily by direct contact with human blood. The major causes of HCV infection worldwide are the use of unsterilized blood transfusions, and re-use of needles and syringes, that have not been adequately sterilized. No vaccine is currently available to prevent hepatitis C and treatment for chronic hepatitis C is too costly for most persons in developing countries to afford. Thus from a global perspective, the greatest impact on hepatitis C disease burden will likely be achieved by focusing efforts on reducing the risk of HCV transmission from nosocomial exposures (e.g. blood transfusions, unsafe injection practices) and high-risk behaviours (e.g. injection drug use).

Hepatitis C virus (HCV) is one of the viruses (A, B, C, D, and E), which together account for the vast majority of cases of viral hepatitis. It is an enveloped RNA virus in the *Flaviviridae* family which appears to have a narrow host range. Humans and chimpanzees are the only known species susceptible to infection, with both species developing similar disease.

An important feature of the virus is the relative mutability of its genome, which in turn is probably related to the high propensity (80%) of inducing chronic infection. HCV is clustered into several distinct genotypes which may be important in determining the severity of the disease and the response to treatment.

The incubation period of HCV infection before the onset of clinical symptoms ranges from 15 to 150 days. In acute infections, the most common symptoms are fatigue and jaundice; however, the majority of cases (between 60% and 70%), even those that develop chronic infection, are asymptomatic. About 80% of newly infected patients progress to develop chronic infection. Chronic develops in about 10% to 20% of persons with chronic infection, and liver cancer develops in 1% to 5% of persons with chronic infection over a period of 20 to 30 years. Most patients suffering from liver cancer, who do not have hepatitis B virus infection have evidence of HCV infection. The mechanisms by which HCV infection leads to liver cancer are still unclear. Hepatitis C also exacerbates the severity of underlying liver disease when it coexists with other hepatic conditions. In particular, liver disease progresses more rapidly among persons with

alcoholic liver disease and HCV infection. HCV is spread primarily by direct contact with human blood. Transmission through blood transfusions that are not screened for HCV infection, through the reuse of inadequately sterilized needles, syringes or other medical equipment, or through needle-sharing among drug users, is well documented. Sexual and perinatal transmission may also occur, although less frequently. Other modes of transmission such as social, cultural, and behavioural practices using percutaneous procedures (e.g. ear and body piercing, circumcision, tattooing) can occur if inadequately sterilized equipment is used. HCV is not spread by sneezing, hugging, coughing, food or water, sharing eating utensils or casual contact.

In both developed and developing countries, high risk groups include injecting drug users, recipients of unscreened blood, haemophiliacs, dialysis patients and persons with multiple sex partners who engage in unprotected sex. In developed countries, it is estimated that 90% of persons with chronic HCV infection are current and former injecting drug users and those with a history of transfusion of unscreened blood or blood products. In many developing countries, where unscreened blood and blood products are still being used, the major means of transmission are unsterilized injection equipment and unscreened blood transfusions. In addition, people who use traditional scarification and circumcision practices are at risk if they use or re-use unsterilized tools.

WHO estimates that about 170 million people, 3% of the world's population, are infected with HCV and are at risk of developing liver cirrhosis and/or liver cancer. The prevalence of HCV infection in some countries in Africa, the Eastern Mediterranean, South-East Asia and the Western Pacific (when prevalence data are available) is high compared to some countries in North America and Europe.

Diagnostic tests for HCV are used to prevent infection through screening of donor blood and plasma, to establish the clinical diagnosis and to make better decisions regarding medical management of a patient. Diagnostic tests commercially available today are based on Enzyme Immunoassay assays (EIA) for the detection of HCV specific antibodies. EIAs can detect more than 95% of chronically infected patients but can detect only 50% to 70% of acute infections. A recombinant immunoblot assay (RIBA) that identifies antibodies which react with individual HCV antigens is often used as a supplemental test for confirmation of a positive EIA result. Testing for HCV circulating by amplification tests RNA (e.g. polymerase chain reaction or PCR, branched DNA assay) is also being utilized for confirmation of serological results as well as for assessing the effectiveness of antiviral therapy. A positive result indicates the presence of acute infection and a potential for spread of the infection and/or the development of chronic liver disease.

Antiviral drugs such as interferon taken alone or in combination with ribavirin, can be used for the treatment of persons with chronic hepatitis C, but the cost of treatment is very high. Treatment with interferon alone is effective in about 10% to 20% of patients. Interferon combined with ribavirin is effective in about 30% to 50% of patients. Ribavirin does not appear to be effective when used alone.

There is no vaccine against HCV. Research is in progress but the high mutability of the HCV genome complicates vaccine development. Lack of knowledge of any protective immune response following HCV infection also impedes vaccine research. It is not known whether the immune system is able to eliminate the virus.

DIA.PRO
Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy
Phone +39 02 27007161
Fax +39 02 26887226.
e-mail: info@dipro.it

REF CVAB CE
96.192.480.960 TS&S

Some studies, however, have shown the presence of virus neutralizing antibodies in patients with HIV infection. In the absence of a vaccine, all precautions to prevent infection must be taken including (a) screening and testing of blood and organ donors; (b) Virus inactivation of plasma derived products; (c) implementation and maintenance of infection control practices in health care settings; including appropriate sterilization of medical and dental equipment; (d) promotion of behaviour change among the general public and health care workers to reduce overuse of injections and to use safe injection practices; and (e) Risk reduction counselling for persons with high-risk drug and sexual practices."

The genome encodes for structural components, a nucleocapsid protein and two envelope glycoproteins, and functional constituents involved in the virus replication and protein processing. The nucleocapsid-encoding region seems to be the most conservative among the isolates obtained all over the world.

C. PRINCIPLE OF THE TEST

Microplates are coated with HCV-specific antigens derived from core and its regions encoding for conservative and immunodominant antigenic determinants (Core peptide, nonstructural NS3, NS4 and NS5 peptides). The solid phase is first treated with the diluted sample and HCV antibodies are captured, by the antigens. After washing out all the other components of the sample, in the presence of a peroxidase substrate, the HCV antibodies are detected by the addition of polyoxonal specific anti HgDdM antibodies, labeled with peroxidase (HRP). The enzyme captured on the solid phase acting on the substrate (peroxidase mixture) generates an optical signal that is proportional to the amount of anti HCV antibodies present in the sample. A cut-off value (at optical densities) is interpreted into HCV antibody negative and positive results.

D. COMPONENTS

Code CVAB/CE contains reagents for 192 tests.

- 1. Microplate MICROPLATE n° 2 microplates
- 12 strips of 8 microwells coated with Core peptide, recombinant NS3, NS4 and NS5 peptides. Plates are sealed in a bag with desiccant.

2. Negative Control CONTROL

1x4 Dmivial. Ready to use control. It contains 1% goat serum proteins, 10 mM Na-acetate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. The negative control is olive green colour coded.

3. Positive Control CONTROL

1x4 Dmivial. Ready to use control. It contains 1% goat serum proteins, human antibodies positive to HCV, 10 mM Na-acetate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. The Positive Control is blue colour coded.

4. Calibrator CAL

n° 2 vials. Lyophilized calibrator. To be dissolved with the volume of ELA grade water reported on the label. It contains bovine serum proteins, human antibodies to HCV whose antigen is adsorbed on the MBSC Working Standard code 99488-002-141, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.3 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

Number of tests	36	480	960
Code	CVAB/CE/36	CVAB/CE/480	CVAB/CE/960
1. Microplate	n° 1	n° 5	n° 10
2. Negative Control	1x4 Dmivial	1x4 Dmivial	1x4 Dmivial
3. Positive Control	1x2 Dmivial	1x2 Dmivial	1x2 Dmivial
4. Calibrator	n° 1 vial	n° 5 vials	n° 10 vials
5. Wash buffer core	1560ml/bottle	540ml/bottles	4x150ml/bottles
6. Enz. Conjugate	1x150ml/bottle	2x40ml/bottles	4x40ml/bottles
7. Chromog/Subst	1x150ml/bottle	1x40ml/bottle	4x40ml/bottles
8. Assay Diluent	1x150ml/bottle	2x40ml/bottles	4x40ml/bottles
9. Sulphuric Acid	1x150ml/bottle	2x40ml/bottles	4x40ml/bottles
10. Sample Diluent	1x150ml/bottle	2x40ml/bottles	4x40ml/bottles
11. Plate seal foil	n° 2	n° 1	n° 1
12. Pack. insert	n° 1	n° 1	n° 1

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

- 5. Wash buffer concentrate WASHBUF 20X 2x60ml/bottle. 20X concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0 +/-0.2, 0.05% Tween 20 and 0.05% Kathon GC.

- 6. Enzyme Conjugate CONJ 2x150ml/bottle. Ready to use and pink/red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG and IgM, 5% BSA, 10 mM Tris buffer pH 6.8 +/-0.1, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

- 7. Chromogen/Substrate SUBS TMB 2x150ml/bottle. Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H2O2.

- 8. Assay Diluent DILUS 1x150ml/bottle. 10 mM Tris buffered solution pH 8.0 +/-0.1 containing 0.1% Kathon GC for the pre-treatment of samples and controls in the plate, blocking interference.

- 9. Sulphuric Acid H2SO4 O.3 M 1x320ml/bottle. It contains 0.3 M H2SO4 solution. Attention: Irritant (H315); H319; P280; P302+P352; P332+P313; P305+P351+P338; P337+P313; P362+P363

- 10. Sample Diluent DILSFE 2x50ml/bottle. It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

- 11. Plate sealing foils n° 4 bluish green in the presence of sample.

- 12. Package insert n° 1 Important note: Only upon specific request, Dia.Pro can supply reagents for 96, 480, 960 tests, as reported below:

E. MATERIALS REQUIRED BUT NOT PROVIDED

- 1. Calibrated Micropipettes (20ul and 10ul) and disposable plastic tips.
- 2. ELA grade water (distilled or deionised, charcoal treated) to remove oxidizing chemicals used as disinfectants).
- 3. Timer with 60 minute range or higher.
- 4. Absorbent paper tissues.
- 5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of +37°C.
- 6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
- 7. Calibrated ELISA microplate washer.
- 8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

- 1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
- 2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
- 3. All the personnel involved in performing the assay have to wear protective laboratory clothes, lab-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in bioassay procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Bioassay in Microbiological and Biomedical Laboratories", ed 1994.
- 4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
- 5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate from strong light and avoid vibration of the bench surface where the test is undertaken.
- 6. Upon receipt, store the kit at 2, 8°C into a temperature controlled refrigerator or cold room.
- 7. Do not interchange components between different lots of the Kits. It is recommended that components between two kits of the same lot should not be interchanged.
- 8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
- 9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
- 10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
- 11. Do not use the kit after the expiration date stated on the external container and internal vials/bottles.
- 12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's Laboratory Safety Manual: "Microbiological and Biomedical Laboratories", ed. 1994.
- 13. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into autoclaved workstations, in order to avoid contamination.
- 14. Waste management during the use of the kit has to be described in compliance with national directives and laws concerning laboratory waste of chemical, biological and medical substances. In particular, liquid waste generated from sample washing procedure, from residuals of control and from samples has to be treated as potentially infective material and inactivated

- 15. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be housed in proper containers designated for laboratory/hospital waste.
- 16. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
- 17. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND RECOMMENDATIONS

- 1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- 2. Avoid any addition of preservatives to samples, especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
- 3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
- 4. Haemolysed (red) and visibly hyperfibrinemic (milky) samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- 5. Sera and plasma can be stored at +2°, -8°C for up to seven days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- 6. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.3µm filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-use of the device and up to 6 months.

1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before operating the container. Check that the desiccant is not turned to dark green, indicating a defect of manufacturing. In this case call Dia.Pro's customer service. Unused strips have to be placed back into the aluminium pouch, in presence of desiccant supplied, firmly zipped and stored at +2°, -8°C.

When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infective, even if HCV, eventually present in the control, has been chemically inactivated.

4. Calibrator:

Dissolve 6 ELA grade water reported on this label. Mix well on vortex before use.

Handle this component as potentially infective, even if HIV, eventually present in the control, has been chemically inactivated.

Note: When dissolved the Calibrator is not stable. Store in aliquots at -20°C.

5. Wash buffer concentrate:
 The 20x concentrated solution has to be diluted with EA grade water up to 1200 ml and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution.
 In the preparation avoid foaming as the presence of bubbles could give origin to a bad washing efficiency.
Note: Once diluted, the wash solution is stable for 1 week at 4-8°C.

6. Enzyme conjugate:
 Ready to use. Mix well on vortex before use.
 Be careful not to contaminate the liquid with oxidizing chemicals, alcohol, detergents or microbeads.
 If this component has to be transferred use only plastic, possibly sterile disposable containers.

7. Chromogen/Substrate:
 Ready to use. Mix well on vortex before use.
 Be careful not to contaminate the liquid with oxidizing chemicals, alcohol, detergents or microbeads.
 Do not expose to strong illumination, oxidizing agents and metallic surfaces.
 If this component has to be transferred use only plastic, possibly sterile disposable containers.

8. Assay Diluent:
 Ready to use. Mix well on vortex before use.

9. Sulphuric Acid:
 Ready to use. Mix well on vortex before use.
 Attention: Iritant (H315; H319; P280; P302+P352; P332+P313; P505+P561+P538; P333+P313; P562+P563).
 Precautionary statements:
 P280 – Wear protective gloves/protective clothing/eye protection/face protection.
 P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.
 P303 + P361 + P353 – IF SKIN irritation occurs: Get medical advice/attention.
 P308 + P313 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P337 + P313 – If eye irritation persists: Get medical advice/attention.
 P362 + P363 – Take off contaminated clothing and wash it before reuse.

10. Sample Diluent:
 Ready to use. Mix well on vortex before use.

1. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submerged to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the assay. They should also be regularly maintained in order to ensure a precision of 1% and a tightness of +4.2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations. Provided that the instrument is validated for the incubation of ELISA tests. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully

validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the sections "Validation of Test" and "Assay Performances". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

4. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performance should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to 2.0; (c) linearity to 2.0; (d) repeatability 2-1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and aspiration) has to be validated and correctly set. Particular attention has to be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended for blood screening when the number of samples to be tested exceeds 20-30 units per run.

7. When using automatic devices, in case the vial holder of the instrument does not fit with the vials supplied in the kit, transfer the solution into appropriate containers and label them with the same label peeled out from the original vial. This operation is important in order to avoid mismatching content of vials when transferring them. When the vials are over, return the secondary labeled containers to 2,8°C, if they are not immediately used.

8. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to best conformance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transparent pipette, no spillage of liquid is present inside the box. Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Calibrator, as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.

6. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.

7. Check that the ELISA reader has been turned on at least 20 minutes before reading.
 8. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
 9. Check that the microplates are set to the required volume.
 10. Check that all the other equipment is available and ready to use.
 11. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

Automated assay:
 In case the test is carried out automatically with an ELISA system, we suggest to make the instrument aspirate 200 µl Sample Diluent and then 10 µl sample.
 All the mixture is then carefully dispensed directly into the appropriate sample well of the microplate. Before the next sample is aspirated, needles have to be duly washed to avoid any cross-contamination among samples.
 Do not dilute control/calibrator as they are ready to use.
 Dispense 200 µl control/calibrator in the appropriate control/calibrator wells.

Important Note: Visually monitor that samples have been diluted and dispensed into appropriate wells. This is simply achieved by checking that the colour of dispensed samples has turned to dark bluish-green while the colour of the negative control has remained olive green.

For the next operations follow the operative instructions reported below for the Manual Assay.
 It is strongly recommended to check that the time lag between the dispensation of the first and the last sample will be calculated by the instrument and taken into consideration by delaying the first washing operation accordingly.

Manual assay:

1. Place the required number of Microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µl of Negative Control in triplicate, 200 µl Calibrator in duplicate and 200 µl Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µl of Sample Diluent (DILSPE) to all the sample wells; then dispense 10 µl sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µl Assay Diluent (DILAS) into all the control/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for 45 min at +37°C.

Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

6. Wash the microplate with an automatic washer by delivering and aspirating 350µl/well of diluted washing solution (see paragraph 3.5.1 section 1.3).
 7. Pipette 100µl Enzyme Conjugate to each well, except the 1st blanking well and cover the reader. Check that the pink/red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for 45 min at +37°C.
9. Wash microwells as in step 5.
10. Pipette 100µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at room temperature (18-24°C) for 15 minutes.

Important note: Do not expose to strong direct illumination. High background might be generated.

11. Pipette 100µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of acid will turn the positive control and positive samples from blue to yellow/brown.
12. Measure the colour intensity of the solution in each well, as described in section 5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

Important notes:

1. If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.
3. Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.
4. The Calibrator (CAL) does not affect the cut-off calculation and therefore the test results calculation. The Calibrator may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

Method	Operations
Controls & Calibrator	200 µl
Samples	200µl dil.+10µl
Assay Diluent (DILAS)	50 µl
1 st incubation	45 min
Temperature	+37°C
Wash step	4-5-cycles
Enzyme conjugate	100µl
2 nd incubation	45 min
Temperature	+37°C
Wash step	4-5-cycles
1A1B1+2D2	100 µl
3 rd incubation	15 min
Temperature	RT
Sulphuric Acid	100 µl
Reading OD	450nm

An example of dispensation scheme is reported below:

Microplate												
A	1	2	3	4	5	6	7	8	9	10	11	12
B	BLK	S2										
C	NC	S3										
D	NC	S4										
E	CALL	S6										
F	CALL	S7										
G	PC	S8										
H	ST	S9										

BLK = Blank NC = Negative Control S = Sample
 CALL = Calibrator PC = Positive Control

O. INTERNAL QUALITY CONTROL
 A check is carried out on the controls and the calibrator any time the kit is used in order to verify whether their OD450nm values are as expected and reported in the table below.

Check	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	< 0.050 mean OD450nm value after blanking
Calibrator	S/Co > 1.1
Positive Control	> 1.000 OD450nm value

If the results of the test match the requirements stated above, proceed to the next section.
 If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromopropylsulfate solution has not anti-contaminated during the assay.
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre-qualification study. 2. that the proper washing solution has been used and the washer has been primed with it before the assay. 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control).
Calibrator	1. that the procedure has been correctly executed. 2. that no mistake has been done in its distribution (dispensation of negative control instead of calibrator) or the washer settings are as validated in the pre-qualification study. 4. that no external contamination of the calibrator has occurred.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly executed. 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of positive control, in this case, the negative control will have an OD450nm value > 0.150, too). 3. that the washing procedure and the washer settings are as validated in the pre-qualification study. 4. that no external contamination of the positive control has occurred.

Should these problems happen, after checking, report any residual problem to the supervisor for further actions.

P. CALCULATION OF THE CUT-OFF
 The tests results are calculated by means of a cut-off value determined with the following formula on the mean OD450nm value of the Negative Control (NC):

$$NC + 0.350 = \text{Cut-Off (Co)}$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Important note: When the calculation of results is done by the operating system of an ELISA automated work station be sure that the proper formulation is used to calculate the cut-off value and generate the right interpretations of results.

Q. INTERPRETATION OF RESULTS
 Test results are interpreted as ratio of the sample OD450nm and the Cut-Off value (as S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HCV or that the blood unit may be transfused.
 Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.
 A positive result is indicative of HCV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

- Important notes:**
1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.
 2. Any positive result should be confirmed by an alternative method capable to detect IgG and IgM antibodies (confirmation test) before a diagnosis of viral hepatitis is formulated.
 3. As proved in the Performance Evaluation of the product, the assay is able to detect seroconversion to anti HCV core antibodies earlier than some other commercial kits. Therefore a positive result, not confirmed with these commercial kits, does not have to be ruled out as a false positive result! The sample has to be anyway submitted to a confirmation test (suggested upon request by DiaPro srl, code CC0491).
 4. As long as the assay is able to detect also IgM antibodies some discrepant results with other commercial products for the detection of anti HCV antibodies - lacking and high conjugate in the formulation of the enzyme tracer and therefore missing light reactivity - may be present. The real positivity of the sample for antibodies to HCV should be then confirmed by examining also light reactivity, important for the diagnosis of HCV infection.
 5. When test results are transmitted from the laboratory to an Informatics centre, attention has to be done to avoid wrong data transfer.
 6. Diagnosis of viral hepatitis infection has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below:

The following data must not be used instead of real figures obtained by the user:

Negative Control: 0.019 - 0.020 - 0.021 OD450nm
 Mean Value: 0.020 OD450nm
 Lower than 0.050 - Accepted
 Positive Control: 2.189 OD450nm
 Higher than 1.000 - Accepted
 Cut-Off = 0.020 + 0.350 = 0.370

Calibrator: 0.550 - 0.530 OD450nm S/Co = 1.4
 Mean value: 0.540 OD450nm
 S/Co higher than 1.1 - Accepted
 Sample 1: 0.070 OD450nm
 Sample 2: 1.690 OD450nm
 Sample 1 S/Co < 0.9 = negative
 Sample 2 S/Co > 1.1 = positive

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

1. LIMIT OF DETECTION

The limit of detection of the assay has been calculated by means of the British Working Standard for anti-HCV, NIBSC code 98/588-003-W1. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined.

Dilution	Lot #1	Lot #2
Factor	S/Co	S/Co
1 X	2.0	2.0
2 X	1.1	1.2
4 X	0.7	0.8
8 X	0.5	0.5
Negative plasma	0.3	0.3

In addition the sample coded Accurn 1 - series 3000 - supplied by Boston Biomedical Inc., USA, has been evaluated. In lot# showing the results below:

CVABCE Lot ID	Accurn 1 Series	S/Co
1201	3000	1.5
6002	3000	1.5
1202	3000	1.9

In addition, 7 samples, tested positive for HCV Ab with Ortho HCV 3.0 S/NA, code 830320, lot# EKED68-1, were diluted in HCV Ab negative plasma (CVABCE, lot# 1202, and Ortho. The following table reports the data obtained:

Sample n°	Limit Dilution	CVABCE S/Co	Ortho 3.0 S/Co
1	256 X	1.9	1.3
2	256 X	2.4	0.7
3	256 X	1.9	1.0
4	128 X	2.5	3.2
5	85 X	3.3	1.4
6	128 X	2.2	0.8
7	135 X	3.2	2.2

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY
 The Performance Evaluation of the device was carried out in a trial conducted on more than total 5000 samples.

2.1 Diagnostic specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. In addition to the first study, where a total of 5043 unselected blood donors (including 14 time donors), 210 hospitalized patients and 162 potentially interfering specimens (other infectious diseases, E.coli antibody positive, patients affected by non viral hepatitis diseases, dialysis patients, pregnant women, hemolyzed, lipemic, etc.), were examined, the diagnostic specificity was recently assessed by testing a total of 2876 negative blood donors on six different lots. A value of specificity of 100% was found.
 No false reactivity due to the method of specimen preparation has been observed. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well, to check for interferences due to collection and storage.
 No interference was observed.

2.2 Diagnostic Sensitivity

It is defined as the probability of the assay of scoring positive in the presence of specific analyte.
 The diagnostic sensitivity has been assessed externally on a total number of 359 specimens; a diagnostic sensitivity of 100% was found. Internally more than other 50 positive samples were tested, providing a value of diagnostic sensitivity of again 100%.
 Positive samples, from indications carried out by different genotypes of HCV were tested as well.
 Furthermore, most of seroconversion panels available from Boston Biomedical Inc., USA, (PHV) and Zepionetix, USA (HCV) have been studied.
 Results are reported below for some of them.

Panel	N° samples	DiProCo	Ortho**
PHV/901	11	9	9
PHV/904	7	2	4
PHV/905	9	9	4
PHV/906	7	7	7
PHV/907	7	3	7
PHV/908	13	10	8
PHV/909	3	2	2
PHV/910	5	3	3
PHV/911	5	3	3
PHV/912	4	1	2
PHV/914	4	2	1
PHV/915	4	2	5
PHV/916	4	3	0
PHV/917	10	4	3
PHV/918	10	6	6
PHV/919	8	2	0
PHV/920	3	3	3
HCV/10039	10	6	6
HCV/6272	5	2	0
HCV/10155	9	6	7
HCV/10155	9	5	4

Note: ** Positive samples detected
 ** HCV v.3.0

Finally the Product has been tested on the panel EFS Ac HCV, lot n° 0106.03.22/01/14, supplied by the Etablissement Français Du Sang (EFS), France, with the following results:

EFS Panel Ac HCV

Sample	1st run S/Co	2nd run S/Co	3rd run S/Co	Results
HCV 1	2.2	2.4	2.6	positive
HCV 2	1.6	2.0	2.1	positive
HCV 3	1.5	1.7	1.6	positive
HCV 4	5.2	6.5	5.5	positive
HCV 5	1.6	1.8	1.6	positive
HCV 6	0.4	0.4	0.4	negative

The variability shown in the tables above did not result in sample misclassification.

S. LIMITATIONS

Repeatable false positive results, not confirmed by RIBA or similar confirmation techniques, were assessed as less than 0.1% of the normal population.
Frozen samples containing fibrin particles or aggregates after thawing have been observed to generate some false results.

REFERENCES

3. PRECISION:
It has been calculated on two samples, one negative and one low positive, examined in 16 replicates in three separate runs. Results are reported as follows:

Lot # 1202

Negative Sample (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.094	0.099	0.096	0.096
Std.Deviation	0.008	0.007	0.008	0.007
CV %	8.7	6.6	7.9	7.7

Cal # 2 - 7K (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.396	0.403	0.418	0.406
Std.Deviation	0.023	0.029	0.027	0.026
CV %	5.9	7.1	6.4	6.5
S/Co	1.1	1.1	1.2	1.1

Lot # 0602

Negative Sample (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.097	0.096	0.094	0.096
Std.Deviation	0.009	0.010	0.006	0.009
CV %	8.9	10.1	6.4	9.1

Cal # 2 - 7K (N = 16)

Mean values	1st run	2nd run	3rd run	Average value
OD 450nm	0.400	0.395	0.393	0.396
Std.Deviation	0.021	0.025	0.028	0.024
CV %	5.4	6.2	6.6	6.1
S/Co	1.2	1.2	1.1	1.2

Lot # 0602Z

Negative Sample (N = 16)

Mean values	1st run	2nd run	3rd run	Average
OD 450nm	0.087	0.091	0.088	0.089
Std.Deviation	0.009	0.007	0.008	0.008
CV %	10.0	8.2	8.5	8.9

Cal # 2 - 7K (N = 16)

Mean values	1st run	2nd run	3rd run	Average
OD 450nm	0.386	0.390	0.391	0.389
Std.Deviation	0.023	0.021	0.023	0.022
CV %	6.0	5.3	5.8	5.7
S/Co	1.1	1.2	1.2	1.2

1. CDC. Public Health Service interagency guidelines for screening donors of blood plasma or tissues and serum for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(NC):1-17.

2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:652-65.

3. McCullien GM, Alter MJ, Meyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Garin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*. Edizon Minerva Medica, Turin, 1997; 267-70.

4. D'Adda MC. Chronic liver disease and cirrhosis. In Eastaugh JF, ed. *Diseases of the Liver*. London: Chapman and Hall, 1996; 11-22.

5. US Department of Health and Human Services. Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447. 615-45.

5. Alter MJ, Hadler SC, Juszon RU, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;263:2323-26.

6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.

7. Alter HJ, Purcell RH, Holland PV, Fornsone SM, Morrow AG, Montenegro Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:839-41.

8. Scott LB, Wright EC, Zimmerman HJ, McCullien GM, VA Cooperative Studies Group. Seroprevalence of hepatitis B, hepatitis C, and transfusion-associated hepatitis. *Am J Med Sci* 1975;270:355-62.

9. Fornsone SM, Kasikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.

10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW, Hagan JF. A cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:389-92.

11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major hepatitis C virus antigen. *Am J Clin Pathol* 1991;94:104-16.

12. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major hepatitis C virus antigen. *Am J Clin Pathol* 1991;94:104-16.

13. Aach RG, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first and second generation assays. *N Engl J Med* 1991;325:1535-39.

14. Alter MJ, Margolis HS, Krawczynski K, Juszon RN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.

15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.

16. Vanara JG, Nelson KE, Madoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1026-31.

17. Zelle JB, Jain S, Krawczynski K, et al. Seroprevalence of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:50-5.

18. Fingerhuth MJ, Jasinsku DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.

19. Garfin RS, Vatoro D, Galia N, Doherty MC, Nelson KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency virus, and human T-lymphotropic viruses. *Am J Public Health* 1993;83:655-61.

20. Bessler DG, Alter HJ, Deinberg JL, Forberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.

21. Tricoli CL, Hollinger FB, Hooks WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.

22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;61:283-9.

23. Tokars JL, Miller EP, Alter MJ, Adjuino ML. National surveillance of dialysis-associated diseases in the United States, 1995. *ASAIO J* 1998;44:99-107.

24. Gerson DH, Chairakis E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:86-71.

25. Weinstock HS, Bolen G, Reingold AL, Pollen LB. Hepatitis C virus infection among injection drug users: a clinic for sexually transmitted diseases. *JAMA* 1995;273:392-4.

26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ, Hagan JF. Hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.

27. Burkholder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect Dis* 1994;269:253-5.

28. Thomas DL, Zaitman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore—an analysis of 309 sex partners. *J Infect Dis* 1995;171:768-75.

29. Thomas DL, Factor SH, Klein GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.

30. Cooper BW, Kusell A, Titon RC, Goodman R, Lovitz RC. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;15:92-5.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

0318
 Manufacturer:
 Dia Pro, Diagnostic Bioprobes Srl,
 Via G. Carducci n° 27 - Sesto San Giovanni (MI) - Italy

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Este certificado ampara todas las marcas de estos productos incluidas por el fabricante en su declaración de conformidad. / This certificate covers all trademarks of these products included by the manufacturer in his declaration of conformity.
Este certificado ampara a todas las marcas de estos productos incluidas por el fabricante en su declaración de conformidad. / This certificate covers all trademarks of these products included by the manufacturer in his declaration of conformity.