4 L x 1: 105-012291-00

1 L x 4: 105-012292-00

Instruction for Use

Hematology Analyzer Reagent (M-6LH LYSE)

English (1~4)

German (5~9)

Spanish (10~13)

French (14~18)

Italian (19~22)

Turkish (23~26)

Portuguese (27~31)

Russian (32~36)

Russian-Russia (37~41)

Find IFUs in more languages on our website:

<https://www.mindray.com/en/resources-center/ivd-edocs>

© 2012-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.
The issue date of this manual is: 2024.09

Intellectual Property Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

mindray

MINDRAY

are the trademarks, registered or otherwise, of Mindray in China and other countries. All other trademarks that appear in this manual are used only for informational or editorial purposes. They are the property of their respective owners.

Responsibility on the Manufacturer Party

Contents of this manual are subject to change without prior notice.

All information contained in this manual is believed to be correct. Mindray shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this manual.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- 1) Use of this product is conducted as instructed in this operation manual by Mindray authorized personnel or skilled/trained clinical professionals;
- 2) The product is stored and transported as instructed in the operation manual and the labeling;
- 3) The product is used only on the specified models.

Warranty

Free service will be provided to:

- Products that are covered by the warranty terms of Mindray.

Pay service will be provided to:

- 1) Expired products;
- 2) Products within the warranty period, but are damaged by improper use or human error, stored in improper environment, or damaged by force majeure or other reason not caused by the product itself.

Return Procedure

Follow the procedure below to return products to Mindray if needed.

Contact Mindray customer service department, inform them the product name and lot No. (the name and lot No. are labeled on the outer package of the product). If the name and lot No. are not legible, your return request will not be accepted. Please specify the product name, model, lot No. and the return reason.

M-6LH LYSE

PRODUCT

M-6LH LYSE

MODEL

M-6LH

SPECIFICATIONS

4Lx1, 1Lx4

INTENDED USE

The M-6LH LYSE formulated to measure the hemoglobin-related parameters.

PRINCIPLE

The M-6LH LYSE is formulated to lyse red blood cells, release hemoglobin in red blood cells and transform it into methemoglobin, which enables the measurement of hemoglobin-related parameters by using the colorimetric method.

ACTIVE INGREDIENTS

Surfactant.....≤0.5%

Tris Buffer.....≤0.5%

STORAGE AND EXPIRATION DATE

This product shall be stored in ventilated room free from corrosive gas at 2°C~30°C with the humidity lower than 90%. The shelf life is 2 years.

The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

APPLICABLE INSTRUMENT

This product applies to BC-6000 / BC-6000Plus / BC-6100 / BC-6100Plus / BC-6200 / BC-6600Plus / BC-6700Plus / BC-6800Plus / BC-7500[B] CRP / BC-7500[N] CRP / BC-7500[R] CRP / BC-7500[NR] CRP / BC-7300[B] CRP / BC-7300[N] CRP / BC-7300[R] CRP / BC-7300[NR] CRP / BC-7500[B] CS / BC-7500[N] CS / BC-7500[R] CS / BC-7500[NR] CS / BC-700[B] / BC-700[R] / BC-720[R] / BC-760[B] / BC-760[R] / BC-780[R] / BC-700[B] CS / BC-700[R] CS / BC-760[B] CS / BC-760[R] CS / BC-6800 / BC-6600 / BC-780[B] CS / BC-780[R] CS / BC-770[B] CS / BC-770[R] CS / BC-7600 / BC-7800 / BC-7900 Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

SAMPLE REQUIREMENT

Fresh human whole blood anticoagulated. Do not use contaminated samples.

MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

INSTRUCTIONS FOR USE

- 1) Restore the product to usage temperature;
- 2) Open the external packing and insert the pickup tube into the reagent container based on the matching color of the reagent cap and the connector of analyzer cap assembly;
- 3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
- 4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis. See the operator's manual of the analyzer for details.

CUT-OFF VALUE/REFERENCE INTERVALS

Not applicable.

RESULT ELABORATION

Not applicable.

LIMITATIONS

Not applicable.

PRODUCT SPECIFICATIONS

1. Appearance: Transparent liquid without deposits, suspended grains or flocks.
2. Blank count results: The blank count results of the reagent tested on the Mindray Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1 Blank Count Requirements

Parameter	Blank Count Requirements
HGB	≤1g/L

PRECAUTIONS

1. For professional use in vitro diagnosis only.
2. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
3. Do not use the reagent if it is frozen.
4. If the blank count is abnormal after the reagent is transported, place it still for 24

- hours at room temperature before using it.
- If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
 - Dispose of the waste, residual and contaminated packing based on local regulations.
 - The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
 - Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
 - Disposal of waste liquid and materials should be in accordance with local guidelines.
 - The Material Safety Data Sheet (SDS) is available upon request
 - All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
 - Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

GRAPHICAL SYMBOLS

		
Batch code	Use-by date	Temperature limit
		
HUMIDITY LIMITATION	<i>In vitro</i> diagnostic medical device	Manufacturer
		
CATALOGUE NUMBER	Authorized representative in the European Community	Consult instructions for use

European Conformity

Unique Device Identifier

REFERENCES

Not applicable.

INTELLECTUAL PROPERTY STATEMENT

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

COMPANY CONTACT

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
E-mail Address: service@mindray.com
Tel: +86 755 81888998
Fax: +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

Approval Date of the Operator's Manual

2024.07

Revision bar indicates update to previous version.

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