FOR INFORMATION ONLY.
WHEN PERFORMING
THE ASSAY ALWAYS REFER
TO PACKAGE INSERT
SUPPLIED
WITH THE KIT

RO

CanAg CA242 EIA

REF

101-10

IVD

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EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIIFICAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Kölblik kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunautoti iki/Izletot līdz/Houdbaar
tot/Brukes innen/Użyty przed/
Prazo de validade/Expiră la/
Использовать до/Använd före/
Použite ně do/ Uporabno do/
Upotrebljivo do/Son Kullanma Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzzszám/Codice del lotto/
Partijos kodas/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Dاتا на производство/Datum výroby/
Produktionsdato/Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/
Date de fabrication/Datum proizvodnje/
Gyártási idő/Data di produzione/
Pagaminimo data/Ražošanas datums/
Productiedatum/Fremstillingsdato/
Data produkcji/Data de fabrico/Data fabricației/Дата производства/
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

IVD

In Vitro Diagnostic Medical Device/
Медицински уред за диагностика
ин vitro/Diagnostický zdravotnícký
prostředek in vitro/Medicinsk udstyr til
in vitro-diagnostik/In-vitro-Diagnostikum/
Ιατροτεχνολογικό προϊόν για διάγνωση
In Vitro/Dispositivo médico para
diagnóstico in vitro/In vitro diagnostiline
meditsiiniseade/Dispositif médical de
diagnostic in vitro/Diagnostički medicinski
uređaj In Vitro/In vitro orvosdiagnostikai
eszköz/Dispositivo medico per test
diagnostici in vitro/In Vitro Diagnostinė
Medicinos Priemonė/Mediciniska ierice
in vitro diagnostikai/In vitro-diagnostisch
medisch instrument/In vitro diagnostisk
medisinsk utstyr/Wyrób medyczny do
diagnostyki in vitro/Dispositivo Médico
de Diagnóstico In Vitro/Dispozitiv medical
pentru diagnostic in vitro/Только для
диагностики In Vitro/Endast för in
vitro-diagnostik/ Zdravotnicka pomůcka na
diagnostiku in vitro/In vitro diagnostični
pripomoček/Diagnostički medicinski
uređaj In Vitro/<96> testleri için yeterlilik
içerir

REF

Catalogue number/Каталожен номер/
Katalogové číslo/Katalognummer/
Bestellnummer/Αριθμός καταλόγου/
Número de catálogo/Katalogi number/
Numéro de catalogue/Kataloški broj/
Katalógusszám/Numero di catalogo/
Katalogo numeris/Numurs katalogā/
Catalogusnummer/Katalognummer/
Numer katalogowy/Número do catálogo/
Număr de catalog/Номер по каталогу/
Produktnummer/Katalogové číslo/
Kataloška številka/Kataloški broj/
Katalog numarası

Temperature limitation/
Температурни граници/
Teplotní omezení/
Temperaturbegrænsning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limites de temperatura/
Temperatuuri piirang/
Limite de température/
Temperaturno ograničenje/
Hőmérsékletre vonatkozó korlátozás/
Limiti di temperatura/
Temperatūriniai apribojimai/
Temperatūras ierobežojums/
Temperatuurbeperking/
Temperaturbegrensninger/
Temperatury graniczne/
Limite de temperatura/
Limite de temperatură/
Температурный режим/
Temperaturbegrænsning/
Teplotné obmedzenie
Omejitve temperature/
Temperaturno ograničenje/
Sıcaklık sınırlaması/

Contains sufficient for <96> tests/Съдържа
достатъчно количество за тестове
<96>/Lze použít pro <96> testů/Ineholder
tilstrækkeligt/Inhalt ausreichend für <96>
Prüfungen/Περιεχόμενο επαρκές για
«96» εξετάσεις/Contenido suficiente para
<96> ensayos/Kogusest piisab <96> testi
läbiviimiseks/Contenu suffisant pour "96"
tests/Sadržj dovoljno za <96> testova/A
doboz tartalma <96> vizsgálat elvégzéséhez
elegendő/Contenuto sufficiente per "96"
saggi/Turinys skirtas atlikti <96> tyrimus/
Saturis pietiekams <96> testiem/Inhoud
voldoende voor "96" testen/til "96" test/
Tilstrækkelig innhold for <96> prøver/
Wystarczy na wykonanie <96> testów/
Conteúdo suficiente para "96" ensaios/
Conținut suficient pentru 96 de teste/
Содержит достаточные количества для
«96» определений/Innehåller tillräckligt
till "96" antal tester/Obsah postačuje na
tento počet testov: <96>/Vsebinsa zadostuje
za <96> testov/Sadržina dovoljna za <96>
testova/<96> testleri için yeterlilik içerir

Consult Instructions for Use/
Прочетете инструкцията за
употреба/Konzultujte s návodom
k použití/Se brugsanvisning/Siehe
Gebrauchsanweisung/Συμβουλευτείτε
τις Οδηγίες σχετικά με τη χρήση/
Consulte las instrucciones de uso/
Vt kasutusjuhendit/Consulter le mode
d'emploi/Pročítajte upute za uporabu/
Olvassa el a használati utasítást/
Consultare le istruzioni per l'uso/Dél
naudojimo žiūrėkite instrukcijas/Izlasiet
lietošanas instrukciju/Raadpleeg de
instructies voor gebruik/Les instruksene
før bruk/Sprawdzić w instrukcji użycia/
Consulte as Instruções de Utilização/
Consultajti instrukcijunile de utilizare/
Обратитесь к инструкции по
применению/Se brugsanvisning/
Prečítajte si návod na používanie/
Pročítajte uputstvo za upotrebu/
Kullanım Talimatlarına Bakınız

CONT

Contents of kit/Съдържание на набора/
Obsah soupravy/Kitkets indhold/Inhalt
des Kits/Περιεχόμενα του kit/Contenido
del kit/Komplekt sisaldab/Contenu du
kit/Sadržaj opreme/A készlet tartalma/
Contenuto del kit/Rinkinio turinys/
Komplekta saturs/Inhoud van de set/
Settets innhold/Zawartość zestawu/
Conteúdo do kit/Conținutul setului/
Компоненты набора/Kit innehåll/
Obsah súpravy/Vsebina kompleta/Sadržaj
opreme/Kitin ičindėkiler

Biological risks/Биологическа
опасност/Biologická rizika/Biologisk
fare/Biologische Gefahren/Βιολογικοί
κίνδυνοι/Riesgos biológicos/
Bioloigilised ohud/Risques biologiques/
Biološkli rizici/Biologiai kockázatok/Rischi
biologici/Biologinis pavojus/Biologiskais
risks/Biologische risico's/Biologische
risikoer/Zagrozenie biologiczne/Riscos
biológicos/ Biologisk risk/Pericole
biologice/Биологическая опасность/
Biologicky rizikové/Biologické riziká/
Biološkli rizici/Biolojik riskler

ORIG HUM

Human/C човешки производ/Lidské/
Humanit/Human/δείγματα αναφοράς/
Humano/Inimpăritolu/Humaine/Ljudskog
porjekla/Humán/Origine Umana/
Žmogaus kilmės/Cilvēku izcelsmes/
Human/Menneske/Ludzka/Humano/
Origine umana/Человеческого
происхождения/Human/Ludské/
Humanega izvora/Ljudskog porekla/İnsan

ORIG MOU

From mouse/C миши производ/Мыši/
Fra mus/Maus/από ποττίκι/de ratón/
Hiirtelt/De souris/Mišijeg porjekla/
Egérbőli/Murino/Pelės kilmės/No peles/
Van muizen/Fra mus/Mysia/Do rato/De
la șoareci/Мышиного происхождения/
Från mus/Myšie/Mišjega izvora/Mišijeg
porekla/Faređen

ORIG BOV

Bovine/C говежди производ/
Hovězí/Bovin/Rind/από βοοειδή/
Bovino/Veistelt/Bovine/Rogate stoke/
Szarvasmarha/Bovino/Jaučio/No
liellopa/Bovien/Bovin/Wolowy/Bovino/
Origine bovină/крупного рогатого
скота/Från ko/Hovädzie/Govejega
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтваряне с/
Rozfedte pomocí/Rekonstitueres med/
Rekonstituieren mit/Ανασύσταση με/
Reconstituir con/Lahjendamine/
Reconstituer avec/Rekonstituiraite s/
Feloldáshoz/Ricostituire con/Atkurti,
ištirpdant su/Atšķaidīt ar/Reconstituie
met/Rekonstitueres med/Odtworzyć
za pomocą/Reconstituir com/A
se reconstitui cu/Растворить в/
Rekonstituera med/Rozriedte pomocou/
Rekonstituiraite z/s/Ponovno formiranje
sa/Yeniden oluşturuur

Manufacturer/Производитель/Výrobce/
Producent/Hersteller/Κατασκευαστής/
Fabricante/Tootja/Fabricant/Proizvođač/
Gyártó/Fabbicante/Gamintojas/
Ražotājs/Fabrikant/Produsent/
Producent/Fabricante/Producător/
Производитель/Tilverkare/ Výrobca/
Izdelovalec/Proizvođač/Üretici

INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE

Please visit our website www.fdi.com/ifu to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website www.fdi.com. If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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- V, Roselli M, Guadagni F. (2001) Comparative analysis of CA 242 and CA 19-9 serum tumor markers in colorectal cancer patients. A longitudinal evaluation. *Anticancer Res.* Mar-Apr;21(2B):1263-70.
8. Carpelan-Holmström, M., Haglund, K., Lundin, J., Järvinen, H and Roberts, P., (1996) Pre-operative serum levels of CA 242 and CEA predict outcome in colorectal cancer. *Eur. J. Cancer* 32(7), 1156-1161.
 9. Engarås B. (2003) Individual cutoff levels of carcinoembryonic antigen and CA 242 indicate recurrence of colorectal cancer with high sensitivity. *Dis Colon Rectum.* Mar;46(3):313-21.
 10. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
 11. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
 12. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
 13. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
 14. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

CanAg[®] is a registered trademark of Fujirebio Diagnostics AB

Fujirebio Diagnostics AB

Elof Lindälvs gata 13

SE-414 58 Göteborg

Sweden

Phone + 46 31 85 70 30

Fax + 46 31 85 70 40

info@fdab.com

www.fdab.com

Kit Contents

ITEM	QUANTITY
Microplate, streptavidin coated, 96 wells	1 plate, 12 x 8 breakable
CALIBRATORS	
CA242 0 U/mL	1 x 0.75 mL
CA242 15 U/mL	1 x 0.75 mL
CA242 50 U/mL	1 x 0.75 mL
CA242 100 U/mL	1 x 0.75 mL
CA242 150 U/mL	1 x 0.75 mL
CONTROLS	
CA242 Control 1	1 x 0.75 mL
CA242 Control 2	1 x 0.75 mL
Biotin Anti-CA242 monoclonal antibody	1 x 15 mL
Tracer, HRP Anti-CA242	1 x 0.75 mL
TMB HRP-Substrate	1 x 12 mL
Stop Solution	1 x 15 mL
Wash Concentrate	1 x 50 mL

LITERATURE REFERENCES

- Duffy MJ, van Dalen A., Haglund C., et al. Clinical utility of biochemical markers in colorectal cancer: European Group on tumour Markers (EGTM) guidelines. *Eur. J. Cancer* 2003; 39:718-737.
- Carpelan-Holmström M., Haglund C., Lundin J., Järvinen H and Roberts P., (1996) Pre-operative serum levels of CA 242 and CEA predict outcome in colorectal cancer. *Eur. J. Cancer* 32(7), 1156-1161.
- Carpelan-Holmström M., Haglund C., Lundin J., et al. Independent prognostic value of preoperative serum markers CA242, specific tissue polypeptide antigen and human chorionic gonadotropin beta, but not of carcinoembryonic antigen or tissue polypeptide antigen in colorectal cancer. *Br. J. Cancer* 1996; 74:925-929.
- Von Kleist, Hesse, Kananeeh, (1996) Comparative Evaluation of Four Tumor Markers, CA 242, CA 19-9, TPA, and CEA in Carcinomas of the Colon, *Anticancer Research* 16: 2325-2332.
- Carpelan-Holmström M., Haglund C., Kuusela P., et al. Preoperative serum levels of CEA and CA242 in colorectal cancer. *Br. J. Cancer* 1995; 71:868-872.
- Hall NR., Finan PJ., Stephenson MB., et al. The role of CA242 and CEA in surveillance following curative resection for colorectal cancer. *Br. J. Cancer* 1994; 70:549-553.
- Nilsson O., Johansson C., Glimelius B., Persson B., Nørgaard-Pedersen B., Andrén-Sandberg Å., and Lindholm L. (1992): Sensitivity and specificity of CA242 in gastro-intestinal cancer. A comparison with CEA, CA50 and CA19-9. *Br J Cancer* 65, 215-221.
- Carpelan-Holmström M., Haglund C., et al. Serum CA 242 and CEA detect different patients with recurrent colorectal cancer. *Anticancer Res* 1996, 16, 981-986.
- Kawa, S., Tokoo, M., Hasebe, O., Hayashi, K., Imai, H., Oguchi, H., Kiyosawa, K., Furuta, S., and Homma, T., (1994) Comparative study of CA242 and CA19-9 for the diagnosis of pancreatic cancer, *Br. J. Cancer*, 70, 481-486.
- Lundin J., Roberts PJ., et al. (1995) Prognostic Significance of Serum CA 242 in Pancreatic Cancer. A Comparison with CA 19-9. *Anticancer Research* 15: 2181-2186.
- Dahlén U., Karlsson B., Lindholm L., Nilsson O., (1993) Development of an enzyme immuno-assay for determination of the tumour associated antigen CA242, *J. Tumor Marker Oncology* 8, 3, p 111.


Fujirebio Diagnostics Inc.
201 Great Valley Parkway
Malvern, PA 19355-1307
USA
+1 610 240 3800 phone
+1 610 240 3803 fax
info@fdi.com e-mail
www.fdi.com website

Fujirebio Diagnostics AB
Majnabbeterminalen
SE-414 55 Göteborg
Sweden
+46 31 85 70 30 phone
+46 31 85 70 40 fax
info@fdab.com e-mail
www.fdab.com website

PRODUCT INFORMATION

CanAg CA242 EIA

The tumor marker CA242

...is defined by the monoclonal antibody C242. Serum levels of CA242 are low in healthy subjects and subjects with benign diseases, while elevated levels are commonly found in serum from patients with gastro-intestinal cancer. The CA242 marker may be used as an aid in the management of patients with known or suspected gastro-intestinal carcinomas as a complement to existing clinical and laboratory methods.

USA:
For Research Use Only
– not for use in diagnostic procedures

Colorectal cancer

Colorectal cancer is one of the most common malignancies in the world. The main indication for CA242 is for prognosis and therapy monitoring of patients with colorectal carcinoma. Patients with high preoperative serum concentration of CA242 have a poorer clinical outcome than patients with low concentration⁽¹⁻³⁾.

The most widely used tumor marker in colorectal cancer, CEA, lacks sensitivity and specificity in early stage disease⁽¹⁾. CEA and CA242 are expressed independently of each other and a combination of both markers have shown to improve the diagnostic sensitivity compared to CEA separate^(4,5,6). This increased sensitivity is most pronounced in Dukes A-C patients^(5,7). CA242 was also more sensitive for detection of lung metastases than CEA, which showed to be more sensitive in liver metastases in local recurrences⁽⁸⁾. Based on these findings, CA242 may complement CEA for preoperative and postoperative measurements for a better patient management.


Figure 1. Sensitivity of CEA and CA242 and the combination of CEA/242 in colorectal cancer. Carpelan-Holmström et al Br J cancer 1995, 71:868-872.


Figure 2. Sensitivity of CA242 and CA19-9 in the pancreas cancer diagnostic. Kawa et al. Br J Cancer 1994; 70:481-486


Pancreatic cancer

Pancreatic cancer is a frequent neoplasia in industrialized countries. Despite the development of new diagnostic techniques, patients with pancreatic cancer continue to have an extremely poor prognosis. The major reason is that pancreatic carcinoma becomes symptomatic at a late stage, when curative intervention is limited.

There is a need for sensitive and specific biomarkers for the diagnosis and management of patients with pancreatic cancer. The sensitivity of CA242 at 90% specificity ranges from 70-80% in patients with confirmed pancreatic cancer⁽⁹⁾. Increasing levels may indicate progression or poor treatment response, while decreasing concentrations indicate a positive response to therapy. Preoperative serum level of CA242 has been shown to be an independent prognostic factor of overall survival⁽¹⁰⁾.


Assay procedure

Below is an illustration of the CanAg CA242 EIA two-step assay procedure. Results are available within four hours.

CanAg CA242 EIA

The CanAg CA242 EIA is a solid-phase immunoassay based on the sandwich technique⁽¹¹⁾.

SPECIFICATIONS

Results within:	4 hours
Detection limit:	≤ 1 U/mL
Measuring range:	1–150 µg/L
Sample volume:	25 µL
Hook effect:	No hook up to 150 000 µg/L
Stability:	18 months at 2-8° C
Calibrator range:	0–150 U/mL
Incubation temp:	20–28° C
Recovery:	± 10%
Detection:	620 nm or 405 nm

ORDERING INFORMATION

Prod. No. 101-10
CanAg CA242 EIA
For 96 determinations