FCE-CONT-5-v5 (12/2018)

Français - FR

USAGE PRÉVU

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English - EN

INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of ELITech quantitative tests listed in the value sheet.

COMPOSITION

ELITROL I and ELITROL II are lyophilized control sera prepared from human serum.

MATERIALS REQUIRED BUT NOT PROVIDED

- General Laboratory equipment.

WARNINGS AND PRECAUTIONS

- These controls are for professional *in vitro* diagnostic use only.
- These controls are prepared from human material tested by FDA approved methods and found to be non reactive for HbsAg, anti-HCV antibody and anti-HIV1&2 antibodies. However as no test method can rule out the potential risk of infection with absolute certainty, handle cautiously as potentially infectious.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.
- Carefully open the vial, avoiding the loss of lyophilizate.
- Add in exactly 5 mL of distilled/deionized water.
- Carefully close the vial and dissolve the contents completely by occasional gentle swirling during 30 minutes.
- Do not shake to avoid the formation of foam.

Important : To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.

PRODUCT DETERIORATION

- The control may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particles would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the control if the damages of packaging might have an effect on the product performance (leakages).

STABILITY AND STORAGE

Prior to reconstitution :

- Stored at 2-8 °C and protected from light.
- Do not use after expiration date indicated on the vial label.

After reconstitution :

- Stability of the components:
Between 2 and 8 °C : 5 days
Between -25 and -15 °C : 4 weeks (when frozen once)
- Exceptions :
- Stability of total bilirubin (when stored protected from light) :
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 5 days (when frozen once)
- Stability of direct bilirubin (when stored protected from light) :
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 5 days (when frozen once)
- Stability of alkaline phosphatase:
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 4 weeks (when frozen once)

Note :

- Store control sera tightly capped and protected from light after reconstitution.
- Mix content prior to use.
- Frozen specimen should not be thawed at 37°C, and repeated freezing and thawing should be avoided, as well as violent mixing, which may denature ferritin.

PROCEDURE

The control sera are handled as samples according to the instructions use of each reagent. It is recommended to use these sera together with samples, at least once a day and after each calibration.

CONTROL VALUES

The concentrations and activities of the components are lot-specific. The information of traceability values are indicated in the data sheet enclosed in the kit. Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values.

For users with Selectra Touch ProSoftware, use the values included in the barcode available on the value sheet.

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TO WHOM TO BE CONCERNED

We, Seppim S.A.S., manufacturers of Elitech Clinical Systems reagents, having our factory at Zone Industrielle, 61500 Sées - France, confirm that our clinical reagents have been validated on Vital Scientific equipment. As such available Elitech Clinical Systems reagent applications for Vital Scientific instruments are CE-IVD compliant.

Reagents, other than Elitech Clinical Systems reagents, are not validated on Vital Scientific equipments, and we also can't know the impact of other reagents on Vital Scientific equipments.

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22 mai 2012

Signed on behalf of the manufacturer
Valérie GOURDON
Regulatory Affairs Manager
COMPANY SEPPIM S.A.S

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English - EN

INTENDED USE

ELITech Clinical Systems ELICAL 2 is a multi-parametric calibrator for in vitro diagnostic use in the calibration of quantitative ELITech Clinical Systems methods.

DESCRIPTION

ELITICAL 2 is a lyophilized calibrator prepared from human serum. The concentrations and activities of the calibrator components have been adjusted to ensure optimal calibration of the clinical chemistry reagents.

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- General Laboratory equipment.

  WARNINGS AND PRECAUTIONS

- The calibrator is for professional in vitro diagnostic use only.
- The calibrator is prepared from human material tested by FDA or CE approved methods and found to be negative for HbsAg, anti-HCV antibody and anti-HIV1&2 antibodies. However as no test method can rule out the potential risk of infection with absolute certainty, handle cautiously as potentially infectious.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local and legal requirements.

PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.
- Add in exactly 3 mL of distilled/deionized water.
- Carefully close the vial and dissolve the contents completely by occasional gentle swirling within 30 minutes.
- Do not shake to avoid the formation of foam.

PRODUCT DETERIORATION

The calibrator may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particles would indicate deterioration.

- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

DAMAGED PACKAGING

Do not use the calibrator if the damages of packaging might have an effect on the product performance (leakages).

STABILITY AND STORAGE

Prior to reconstitution :

- Stored at 2-8   C and protected from light.
- Do not use after expiration date indicated on the vial label.

After reconstitution :

- Stability of the components:
 Between 15 and 25   C : 8 hours
 Between 2 and 8   C : 2 days
 Between -25 and -15   C : 4 weeks (when frozen once)

Exceptions:

- Stability of direct bilirubin (when stored protected from light):
 Between 15 and 25   C : 3 hours
 Between 2 and 8   C : 8 hours
 Between -25 and -15   C : 2 weeks (when frozen once)

- Stability of total bilirubin (when stored protected from light):
 Between 15 and 25   C : 6 hours
 Between 2 and 8   C : 1 day
 Between -25 and -15   C : 2 weeks (When frozen once)

Note: Store vials tightly capped and protected from light.

PROCEDURE

To calibrate with ELICAL 2 follow the procedure described in the instructions for use of the corresponding reagent.

CALIBRATION VALUES

The concentrations and activities of the components are lot-specific. The information of traceability values are indicated in the data sheet enclosed in the kit.

For users with Selectra Touch Pro Software, use the values included in the barcode available on the value sheet.

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AUTHORIZATION LETTER

TO WHOM IT MAY CONCERN:

MEDICA CORPORATION, having facilities at 5 Oak Park Drive, Bedford, MA 01730, USA, do hereby authorize the company:

GBG-MLD SRL
65 Tighina Site
Office 607
Chisinau, MD-2001
Republic of Moldova

to be our **DISTRIBUTOR** for the **EasyLyte®**, **EasyElectrolytes™**, **EasyBloodGas™** and **EasyStat®** analyzers as well as associated reagents and consumables in **Moldova**.


GBG-MLD SRL is authorized by **MEDICA CORPORATION** to enter tenders and quote for all aforementioned products.

GBG-MLD SRL is authorized by **MEDICA CORPORATION** to present offers on our behalf to tenders placed by the government and other institutions for Medica products and consumables

GBG-MLD SRL responsibilities include sales of the **EasyLyte®**, **EasyElectrolytes™**, **EasyBloodGas™** and **EasyStat®** analyzers and providing service as well as maintaining a supply of reagents and replacement parts.

GBG-MLD SRL is also authorized to provide warranty service for the Medica **EasyLyte®**, **EasyElectrolytes™**, **EasyBloodGas™** and **EasyStat®** analyzers.

This authorization is effective immediately and is valid until December 31, 2022, unless revoked earlier in writing by Medica Corporation.

David Hagopian
VP, Sales & Marketing
MEDICA CORPORATION

2/4/2020
Date

Medica Corporation
5 Oak Park Drive
Bedford, Massachusetts 01730
Tel: 781 275 4892
Fax: 781 275 2731
www.medicacorp.com

Declaration of Conformity

Product Name:

EasyLyte and accessories per attachment

Model/Type:

EasyLyte Na/K, Na/K/Cl, Na/K/Li, Na/K/Cl/Li,
Na/K/Ca/pH, Na/K/Cl/Ca/Li
EasyElectrolytes Na/K/Cl, Na/K/Li

Manufacturer

Medica Corporation
5 Oak Park Drive, Bedford, Massachusetts, 01730, USA

Representative

EC REP Emurgo Europe, Prinsessegracht 20,
2514 AP The Hague, The Netherlands
Tel: +31 70 345 8570
Fax: +31 70 346 7299

Means of Conformity

Medica Corporation declares that the products listed are covered by Annex III of Directive 98/79/EC. These products are self-certified since they are for professional use only and are not listed on Annex II, List A or Annex II, List B of Directive 98/79/EC. In addition, they are in conformity with the Annex I, "Essential Requirements" and provisions of council Directive 98/79/EC for In-Vitro Diagnostic Medical Devices; Directive 2011/65/EU Restriction of Hazardous Substance in Electrical and Electronic Equipment, and the corresponding national laws of the Member States.

Place and Date: Bedford, Massachusetts, USA, September 27, 2018

Signature:

Name: Photios Makris, Ph.D.
Title: VP, Regulatory Affairs

EasyLyte Accessories		
Catalog No.	Accessory	EDMA Code
2004	EasyLyte Na/K Analyzer	21 07 11 02
2014	EasyLyte Plus Na/K/Cl Analyzer	21 07 11 02
2015	EasyLyte Lithium Na/K/Li Analyzer	21 07 11 02
2016	EasyLyte Calcium Na/K/Ca/pH Analyzer	21 07 11 02
2021	EasyLyte Na/K/Cl/Li Analyzer	21 07 11 02
2030	EasyLyte EXPAND Analyzer, Na/K/Cl/Ca-Li	21 07 11 02
2070	EasyLyte EasySampler	21 07 11 02
2101	EasyLyte K+ Electrode	11 04 01 06
2102	EasyLyte Na+ Electrode	11 04 01 07
2113	EasyLyte Cl- Electrode	11 04 01 03
2106	EasyLyte Li+ Electrode	11 04 01 04
2150	EasyLyte Ca++ Electrode	11 04 01 02
2151	EasyLyte pH Electrode	11 70 31 02
2152	EasyLyte Disposable Reference Electrode	11 04 04 01
2103	EasyLyte Reference Electrode	11 04 04 01
2258	EasyLyte Membrane Assembly	21 07 11 02
2120	EasyLyte Na/K 800 ml Solutions Pack	11 04 04 02
2121	EasyLyte Na/K/Cl 800mL Solutions Pack	11 04 04 02
2122	EasyLyte Na/K/Li 800mL Solutions Pack	11 04 04 02
2123	EasyLyte Na/K/Ca/pH 800mL Solutions Pack	11 04 04 02
2028	EasyLyte Na/K/Cl/Li 400mL Solution Pack	11 04 04 02
2109	EasyLyte Na/K 400mL Solutions Pack	11 04 04 02
2112	EasyLyte Na/K/Cl 400mL Solutions Pack	11 04 04 02
2115	EasyLyte Na/K/Li 400mL Solutions Pack	11 04 04 02
2114	EasyLyte Na/K/Ca/pH 400mL Solutions Pack	11 04 04 02
2026	EasyLyte Na/K/Cl/Li 800mL Solution Pack	11 04 04 02
2124	EasyLyte Na/K/Cl/Ca-Li 800mL Solutions Pack	11 04 04 02
2814	EasyQC Bi-Level Quality Control Kit	11 50 02 04
2815	EasyQC Tri-Level Quality Control Kit	11 50 02 04
2843	EasyLyte Quality Control Sample Cups (60)	21 07 11 02
2118	Daily Cleaning Solution Kit	11 01 01 27
2598	EasyLyte Daily Cleaner Cup	21 07 11 02
2108	EasyLyte Solutions Valve	21 07 11 02
2107	EasyLyte Sample Probe	21 07 11 02
2257	EasyLyte Sample Detector	21 07 11 02

EasyLyte Accessories, continued			Catalog No.	Accessory	EDMA Code	EDMA Code
			4002	EasyElectrolyte Na/K/Cl Analyzer		21 07 11 02
			4003	EasyElectrolyte Na/K/Li Analyzer		21 07 11 02
			4102	Reagent Module, Na/K/Cl		11 04 04 02
			4103	Reagent Module, Na/K/Li		11 04 04 02
			7205	EasyElectrolyte/EasyStat Na+ Electrode		11 04 01 07
			7206	EasyElectrolyte/EasyStat K+ Electrode		11 04 01 06
			4203	EasyElectrolyte Cl- Electrode		11 04 01 03
			4204	EasyElectrolyte Li+ Electrode		11 04 01 04
			6204	EasyElectrolyte/EasyStat/EasyBloodGas Reference Electrode		11 04 04 01
			4207	EasyElectrolyte Spacer Electrode		11 04 01 90
			4301	EasyElectrolyte Troubleshooting Kit		21 07 11 02
			2118	Daily Cleaning Solution Kit		11 01 01 27
			4402	EasyStat/EasyBloodGas/EasyElectrolyte Red Test Dye Solution		11 30 01 11
			4403	EasyElectrolyte Urine Diluent		11 04 04 90
			2814	Bi-Level Quality Control Kit		11 50 02 04
			2815	Tri-Level Quality Control Kit		11 50 02 04
			4405	EasyElectrolyte Na/K/Cl Demonstration Kit		21 07 11 02
			4406	EasyElectrolyte Na/K/Li Demonstration Kit		21 07 11 02
			4404	EasyElectrolyte Capillary Tube Kit		21 07 11 02
			4306	EasyElectrolyte Sampler		21 07 11 02
			6504	EasyBloodGas/EasyElectrolyte Pump Tube		21 07 11 02
			6505	EasyStat/EasyBloodGas/EasyElectrolyte Printer Paper		21 07 11 02
			4506	EasyElectrolyte Sensor Module		21 07 11 02
			4507	EasyElectrolyte Valve Module		21 07 11 02
			4508	EasyStat/EasyBloodGas/EasyElectrolyte Compression Plate		21 07 11 02
			7302	Probe Wipers		21 07 11 02
			4522	EasyElectrolyte Daily Cleaner Sample Cups		21 07 11 02
			4539	EasyElectrolyte Sensor Module, Li+		21 07 11 02
			6537	EasyElectrolyte/EasyStat/EasyBloodGas Serial Cable, 9-pin		21 07 11 02
			6520	EasyElectrolyte/EasyStat/EasyBloodGas Barcode Reader Kit		21 07 11 02

CERTIFICATE OF REGISTRATION

This is to certify that the quality management system of:

Medica Corporation

Main Site: 5 Oak Park Drive

Bedford, Massachusetts 01730 United States

has been assessed by Intertek as conforming to the requirements of:

ISO 13485:2016

The quality management system is applicable to:

The Design, Development, Manufacture, Service, Distribution of in-vitro diagnostic medical devices, in-vitro diagnostic test kits, in-vitro diagnostic reagents, in-vitro diagnostic analyzers/software used in the diagnosis and management of cancer, immune status, disease status, autoimmune status, cardiac markers, protein metabolism, endocrine disorders, blood analytes, urinalysis, blood gases.

Certificate Number:

0082581-01

Initial Certification Date:

2009-04-17

Certificate Issue Date:

2019-01-01

Certificate Expiry Date:

2021-04-16

A handwritten signature in black ink, appearing to read 'Calin Moldovean'.

Calin Moldovean

President

Intertek Testing Services NA Ltd.,
1829, 32nd avenue, Lachine, QC, H8T 3J1,
Canada

EasyLyte EasyBloodGas EasyStat

Training Certificate

This is to certify that

Seaver, Sergio
Of Global Biomedical Corp

has completed training for the operation and service of the
EasyLyte, EasyBloodGas, and EasyStat analyzers.

MEDICA

Date November 25, 2004

Randall Rollins
Signed: Randall Rollins
Technical Service Manager

Declaration of Conformity

helena
Biosciences Europe

HL-7-0664DC DOI 2015/08 (1)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5267L	Thromboplastin L	55983

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 06 Aug 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0135DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.


Product Code	Description	GMDN Classification Code
5183	Routine Control SA	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0137DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5186	Routine Control N	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0512DC DOI 2015/08 (5)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5556	Clauss Fibrinogen 50	55997

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed

Date: 12 Aug 2015

Tel +44 (0)191 482 8440

Fax +44 (0)191 482 8442

info@helena-biosciences.com

www.helena-biosciences.com

Helena Biosciences Europe

Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,

United Kingdom

Clauss Fibrinogen 50

REF 5556
REF 5556H

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-1771P 2012/03 (7)

Clause Fibrinogen 50 Instructions for use

INTENDED USE
For the quantitative determination of fibrinogen, based on the Clauss method, in human citrated plasma on Helena Biosciences Europe AC-4, C-1, C-2, C-4 Coagulation Systems.

Clauss Fibrinogen 50 is intended for the quantitative determination of fibrinogen in human plasma. Clauss' developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of citrated plasma after the addition of Thrombin (500 NIH units/mL). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use. Decreased levels can be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenemia (no detectable fibrinogen), hypofibrinogenemia (<1 mg/mL), and dysfibrinogenemia (abnormal fibrinogen molecule).

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HBsAg), HIV-1 and 2 antibody and HCV antibody, however they should be handled with the same precaution as a human plasma sample.

COMPOSITION

Component	Content	Description	Preparation
Thrombin 50 NIH/mL	5 x 4 mL (REF 5556) 5 x 2 mL (REF 5556H)	Contains approximately 200 (REF 5556) / 100 (REF 5556H) units bovine Thrombin with stabilisers.	Reconstitute each vial with: 4 mL (REF 5556) 2 mL (REF 5556H) of distilled water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake.
Fibrinogen Calibrator	2 x 1 mL (REF 5556) 1 x 1 mL (REF 5556H)	Contains 1 mL of lyophilized normal human plasma assayed for fibrinogen levels.	Reconstitute each vial with exactly 1 mL of distilled water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
Owren's Buffer	2 x 25 mL (REF 5556) 1 x 25 mL (REF 5556H)	Contains 25 mL of buffer which contains butyrol and sodium chloride as preservative.	The buffer is ready for use as packaged.

Each kit contains instructions for use.
Each kit contains kit specific reference values inset.

ITEMS REQUIRED BUT NOT PROVIDED

Refer to the Instrument Operator Manual and/or application reagents contained in this kit are also available as separate items.
REF 5374 Clauss Fibrinogen (Thrombin only)
REF 5378 Clauss Fibrinogen (Thrombin only)
REF 5379 Fibrinogen Calibrator
REF 5375 Owren's Buffer

STORAGE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.
Thrombin 50 NIH/mL Once reconstituted, the reagent is stable for 8 hours at 15 - 30°C, 1 week at 2 - 8°C or 1 month at -20°C.
Fibrinogen Calibrator Once reconstituted, the reagent is stable for 4 hours at 2 - 8°C.
Owren's Buffer Store at 2 - 8°C once opened.

SAMPLE COLLECTION AND PREPARATION

Paste or siliconized glass should be used throughout. Blood (10 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 10 minutes. Plasma should be kept at 2 - 8°C or 18 - 24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

PROCEDURE

Refer to the appropriate Instrument Operator Manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Expected values for fibrinogen in healthy adults are 150-350 mg/dL (1.5-3.5 g/L).

LIMITATIONS

Hepatic levels >0.6 units/mL, and fibrinolytic degradation products >100 mg/mL may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, to assay using an appropriate dilution to bring values into the standard range.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supply the following controls available for use with this product:

REF 5301	Specialty Assayed Control N
REF 5302	Specialty Assayed Control A
REF 5166	Routine Control N
REF 5167	Routine Control A
REF 5163	Routine Control SA

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data. The Clauss Fibrinogen assay is designed to give a linear calibration from 1.5 - 6.5 g/L.

Reproducibility			Inter-assay precision		
Intra-assay precision			Inter-assay precision		
Fibrinogen (g/L)	n	CV (%)	Fibrinogen (g/L)	n	CV (%)
1.24	5	4.88	1.02	10	6.17
3.01	5	3.58	3.62	10	3.75

REFERENCES

- Clauss A (1957) Gerinnungsphysiologische Schnellmethode zur Bestimmung des Fibrinogens, *Acta Haematol*, 17:237-240
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci.* 8:145-182
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI H21-A5
- Scully RE et al. (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, 302(37):37-48
- Okuno T, Sekirno V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, 36(9):199-201

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1. Clauss A (1957) Gerinnungsphysiologische Schnellmethode zur Bestimmung des Fibrinogens, *Acta Haematol*; 17:237-246
2. Shaw TS (1977) Assays for Fibrinogen and Its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci*; 8:145-192
3. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays: Approved Guideline, 5th edn. CLSI H21-A5
4. Scully PE et al. (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, 302(3):37-48
5. Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*; 36(8):186-201

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
Product Code	Description	GMDN Classification Code
5392	Thrombin Time	55987

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Instructions for use

INTENDED PURPOSE

The Thrombin Time kit is intended for carrying out clot based haemostasis assays.

The Thrombin Time reagent contains bovine thrombin. The reagent can be used manually or on semi-automated and automated instruments. The test is commonly applied to detect various causes of interference with normal blood coagulation. Prolongation of the thrombin clotting time can be taken as a qualitative indication of abnormal fibrinogen levels (high or low), or the presence of interfering substances such as Fibrin Degradation Products (FDPs) or heparin. Quantitative evaluation of the possible causes of prolonged thrombin clotting time should be performed as follows: studies, such as aPTT or chromogenic assay for heparin, Clausen thrombin, FDP determinations, heparin neutralization by protamine sulphate or polybrene, normal plasma mixing studies or reptilase assay to distinguish between hypofibrinogenemia and FDP effects.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in vitro diagnostic use only - DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the risk to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Composant	Contient	Description	Préparation
Thrombin Time	10 x 1 mL (REF 5392H) 10 x 2 mL (REF 5392) 10 x 5 mL (REF 5377)	Each vial contains a lyophilized preparation of bovine thrombin with buffers and stabilizers. The reconstituted reagent contains ± 10 NIH units/ml of thrombin. The reagent should be a white lyophilized plug.	Reconstitute the reagent with the recommended volume of purified water: 1 mL - REF 5392H 2 mL - REF 5392 5 mL - REF 5377 Allow to stand for 5 minutes then mix gently by inversion and transfer to a plastic tube. Reconstituted reagent should be a clear, colourless solution.

ITEMS REQUIRED BUT NOT PROVIDED

Any high quality mechanical, opto-mechanical or photo-optical coagulation instrument capable of performing thrombin clotting time testing may be used.

STORAGE, SHELF LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thrombin Reconstituted reagent should be stored at -20°C and is stable for 14 days or at -80°C for 1 month. Reconstituted reagent can also be stored on board a Sysmex CA1500 photo-optical instrument for up to 7 days.

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Plastic or siliconized glass should be used throughout. Blood (or plasma) should be collected into 0.2% or 5% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at -2°C or -18°C to -24°C . Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at -37°C prior to testing. Do not keep at -37°C for more than 5 minutes.

PROCEDURE

- Collect and prepare the blood sample according to the directions outlined in "Sample Collection And Preparation".
- Reconstitute the control plasma according to the package insert provided with each control.
- Prepare the reagent for use in the procedure according to the reconstitution instructions in the "Composition" section.
- Perform all tests in duplicate. Calculate the mean clotting time of the duplicate determinations to the nearest 0.1 seconds. Individual values should be within $\pm 5\%$ of the mean value.

Manual Method

- Pipette 0.2 mL of the patient plasma or control plasma into a reaction tube.
- Inoculate at 37°C for 3 minutes.
- Pipette 0.1 mL of Thrombin Time into the reaction tube containing patient or control plasma whilst simultaneously starting a timer.
- Record the time taken for clot formation to the nearest 0.1 second.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application notes.

INTERPRETATION OF RESULTS

The results of the Thrombin Time test should be reported to the nearest 0.1 seconds. The normal range (usually the mean ± 2 standard deviations) should be established by each laboratory. Results outside the normal range should be considered abnormal and follow-up testing performed.

LIMITATIONS

Expected values for the Thrombin Time test will vary from one laboratory to another depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of sample storage are all very important. Plasma sample collection and storage conditions should be standardized and carefully controlled. Unexpected results should be confirmed by additional testing.

As well as the cause of prolonged thrombin times indicated, a report has suggested that many systemic amyloidosis patients with bleeding complications may have a circulating inhibitor which prolongs the Thrombin Time test. Also, therapeutic levels of heparin may entirely abolish clotting in the Thrombin Time test, although neutralisation with protamine sulphate or polybrene should correct the Thrombin Time.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. Controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with the product:

REF 5186 Routine Control N
REF 5187 Routine Control A
REF 5188 Routine Control SA

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives:

Reproductibility	Intra-assay precision	Inter-assay precision
Sample n	CV (%)	CV (%)
Normal 10	12.5	1.29

BIBLIOGRAPHY

- Laposata et al. The Clinical Haemostasis Handbook, Yearbook Medical Publishers Inc., p219, 1989
- Thompson AR and Barker LA. Manual of Haemostasis and Thrombosis, 3rd Ed., F.A. Davis Co., p62, 1983
- DeMott WR. Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432-433, 1990
- Clinical and Laboratory Standards Institute (2006) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn, CLSI H21-A5
- Gastineau DA et al. (1991) Inhibitor of the Thrombin Time in Systemic Amyloidosis: A Common Coagulation Abnormality, Blood, 77:2637-40

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1. Lupatelli A: *The Clinical Hemostasis Handbook*. Yearbook Medical Publishers Inc. #219, 1965
2. Thompson AR and Harker LA: *Manual of Hemostasis and Thrombosis*, 3rd Ed., F.A. Davis Co., Philadelphia, 1983
3. DeMott WR: *Laboratory Test Handbook*, 2nd Ed., Jacobs D.S. et al (Eds., Lexi-Comp Inc., #432-433, 1990
4. Clinical and Laboratory Standards Institute (2006) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays. Approved Guideline, 5th edn. CLSI: H21-A5
5. DeMott WR, et al: *Effect of Aspirin and Tissue Inhibitor of Metalloproteinase on Systemic Amyloidosis: A Common Gastrointestinal Abnormality*. *Blood* 77:2637-40

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Helena Biosciences Europe

Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442

Info info@helena-biosciences.com
www www.helena-biosciences.com

Helena Laboratories (UK) Ltd trading as Helena Biosciences Europe. Registered in England: 1790207. ISO 13485

Training certificate

helena
Biosciences Europe

This is to certify that

Sergiu Sorocovici

from

IM Global Biomarketing Group

has received training on the following:

Electrophoresis products: SAS-1/2, V8

Haemostasis products: C-series, AC-4, AggRAM and reagents

Service training: AC-4

Signed:

Date: 31st October - 4th November 2011

Tel +44 (0)191 482 8440 info@helena-biosciences.com Queensway South, Team Valley Trading Estate,
Fax +44 (0)191 482 8442 techsupport@helena-biosciences.com Gateshead, Tyne and Wear NE11 1 OSD, United Kingdom

www.helena-biosciences.com

Instrument Training

ELITechGroup
VITAL

SCIENTIFIC

Vital Scientific BV hereby declares that the participant has attended a four days seminar for service engineers and the participant is now a certified engineer for the declared instruments.

Participant: Mr. A. Legun

Company: Global Biomarketing Group-Moldova SRL
Moldova

Instrument: Vitalab: XL Series
E Series
Junior Series
Dry ISE
Micro Series
ProXS


Date of training: April 20th – April 23rd, 2010

System Support Manager:

Jan Oostendorp

System Support Engineer:

Frank v.d. Korput

To: Whom it May Concern

Regulatory status of parts & accessories

As mentioned on the current Declarations of Conformity of our Clinical Chemistry Analyzers also the accessories conform to the provisions of the EU Directive on In Vitro Diagnostic Medical Devices (98/79/EC). This applies to the parts and accessories as mentioned in the attached list.

'IVD accessory' means an article which, whilst not being an IVD medical device, is intended specifically by its manufacturer to be used together with an IVD device to enable that IVD device to be used in accordance with its intended purpose.

ELITechGroup B.V.

Adriaan P. Intveld
Manager Quality Assurance & Regulatory Affairs

Part number	Description	IVD medical device	IVD accessory	general laboratory use	spare part	supporting part
1540-001	Anti-Slip sheet					✓
2206-007	Cooling Liquid (1 L)					✓
3062-021	Sample cup (1000 pcs)		✓			
3062-033	Sample tube 6 ml (500 pcs)					✓
3062-040	Water container 10 L					✓
3062-041	Water container 5 L					✓
3066-155	Syringe 100 µl		✓			
3066-156	Syringe 1 ml		✓			
3069-040	Keyboard Dust cover					✓
3069-047	Keyboard Dust cover					✓
3070-518	Cap holder					✓
3070-538	Cap rotor left					✓
3070-539	Cap rotor right					✓
3201-002	Dichromate 8 Abs (25ml)		✓			
3365-192	USB Stick					✓
3374-003	Main cable (USA)					✓
3374-039	Pumpunit cable		✓			
3374-066	Mains cable					✓
3374-097	Serial Null-modem cable					✓
3374-286	USB Extension cable					✓
4804-038	Reagent identification Disc					✓
6001-826	Diluted Waste container		✓			
6001-827	Concentrated Waste container		✓			
6001-860	Water container		✓			
6001-861	Tube assay (analyser)		✓			
6001-872	Tube assay (cooling unit)		✓			
6002-102	Assorter unit				✓	
6002-386	System software on CD		✓			
6002-706	Reaction Rotor set (3 pcs)		✓			
6002-726	System Disc		✓			
6002-817	Bottle 30 ml (20 pcs)		✓			
6002-818	Bottle 15 ml (20 pcs)		✓			
6002-904	Water container 5 L		✓			
6002-910	Assorter unit				✓	
6002-913	External tubing		✓			
6003-074	System software on USB stick		✓			
6003-444	Diluted Waste Container 5 L		✓			
6003-466	Keyboard Support option					✓
6003-797	CW Waste Container 2L		✓			
6003-808	Assorter unit				✓	

Declaration of conformity

Avantor Performance Materials Poland S.A. who is an established manufacturer of reagents and products for diagnostic in vitro located at:

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Herewith declares the following:

Reagents mentioned in attached list are labeled with J.T.Baker label, comply with the In Vitro Diagnostic Medical Devices Directive 98/79/EC and the requirements of ISO 13485 Standard.

This declaration is the basis for CE marking of the In Vitro Diagnostic Medical Devices.

The products are not part of List A and List B of Annex II of the IVD Directive 98/79/EC but are subject to self registration.

This declaration is valid for all the IVD medical devices described above and which are placed on the market by ourselves on or after the date hereof and which bear the CE marking.

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Date: March 18, 2019

Letter of Authorization

Avantor Performance Materials Poland S.A., reg. No. 0000010108 who is an established manufacturer of Hematology- Reagents, Stains, Controls and Calibrators and products for Histology located at:

Sowińskiego 11
44-101 Gliwice
Poland

herewith confirms that:

I.M Global Biomarketing Group Moldova S.R.L
Republic of Moldova
MD-2001, Chisinau
Tighina str. 65, 607 office
Tel (373 22) 549 120, 549 121
Fax (373,22) 547 373

is authorized to act as our distributor for our hematology/histology reagents and controls (Products) in Moldova

We declare that we will supply the Products for the needs of tenders.
We declare that we will supply the Products for tenders with warranty as per the Avantor General Conditions of Sale.

Furthermore I.M Global Biomarketing Group is duly entitled to:

- Register, promote, offer, negotiate prices and sell our Products in Moldova;
- carry out the required product training of the medical and technical personnel who will use these products.

The product specialists of I.M Global Biomarketing Group have been duly trained and are qualified for providing all services in regards to consulting, sales, maintenance and training.

In all the above activities I.M Global Biomarketing Group is acting in its own name and on its own account.

This authorization letter is valid until about 1 year after date.

Avantor Performance Materials S.A.
Poland

A handwritten signature in black ink, appearing to read 'H van den Berg', with a stylized flourish underneath.

H van den Berg,
Marketing Product Manager Diagnostics

BUREAU VERITAS
Certification

Certificate

Awarded to

Avantor Performance Materials Poland S.A.

ul. Sowlińskiego 11, 44-101 GLIWICE
POLAND

Bureau Veritas Certification certify that the Management System of the above organisation has been audited and found to be in accordance with the requirements of the management system standards detailed below

STANDARD

ISO 9001:2015

SCOPE OF SUPPLY

SALES OF CHEMICAL SERVICES AND CHEMICAL PRODUCTS
INCLUDING FINE CHEMICALS, ENNOBLED CHEMICALS, HIGH
PURITY SOLVENTS, CHEMICAL SERVICES.

PRODUCTION AND TESTING OF CHEMICAL PRODUCTS INCLUDING
FINE CHEMICALS, ENNOBLED CHEMICALS AND HIGH PURITY SOLVENTS.

Certification Cycle Start Date: **15 September 2018**

*Subject to the continued satisfactory operation of the organisation's Management System,
this certificate is valid until:* **14 September 2021**

To check this certificate validity please call: +48 22 549 04 00

*Further clarification regarding the scope of this certificate and the applicability of the management system requirements
may be obtained by consulting the organisation.*

Issue Date: 29 June 2018

Certificate Number: **PL008875/P**

AC 081
QMS

Certificate of Completion

This is to certify

Mr. Alexei Legun

Has successfully completed

The technical maintenance training course

On

Fully Automatic Blood Cell Counter

PCE-210

Particle(Blood Cell)Counter

PCE-170/PCE-170N

Hemoglobin meter

Hb-20N

March 24, 2005

Hiroshi Shimosaka

President

ERMA INC.

Diluid* Erma

Intended use

Diluid* Erma is a specially filtered, non-sterile blood diluting fluid for use in cell counting and sizing.

The reagent is designed for automated instrumentation, capable to monitor a three-part WBC differential, based on the aperture impedance principle and electronically adjusted to operate at an osmolality of 330 ± 20 mOsm/kg. Diluid* Erma should be used in combination with CyMet* ERMA III Diff and Lyzerglobin* PCE.

Summary and principle

The reagent is used to dilute whole blood prior to counting and sizing of RBC, PLT and WBC. Content of the reagent maintains stability of RBC, PLT and WBC during counting.

Content: Diluid* Erma is water based and contains:

NaCl, Na₂SO₄, procaine HCl and preservatives in an inorganic buffer compound.

Warning and precautions

Harmful if swallowed. Avoid contact with eyes, skin and clothing.

Storage and stability: Diluid* Erma is stable for three years at 18-30°C.

Indications of deterioration

There are no visible deteriorations; otherwise the reagent shouldn't give optimised performance. The reagent can be used through out shelf life, giving optimised performance. No guarantees to reagent performance are given after the expiry date.

Instructions for use

Diluid* Erma should be used undiluted according to instrument manufacturer's instructions for use and should be connected as listed in the Operators manual. Reagent may be used with Hypochlorite 0.5% or Proclean* as a cleaning agent. Furthermore reagent may be used with next lysing reagents: with CyMet* ERMA III Diff and Lyzerglobin* PCE.

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ProClean*

Intended use

ProClean* is a specially filtered, non-sterile cleaning fluid for use in cleaning of cell counters.

The product is designed for semi-automated and automated instrumentation, capable to clean blood diluting parts of the instrument.

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Content

ProClean* is water based and contains:

Proteolytic enzym, poly-oxy-ethylene-alkyl-alcohol, NaCl, Na₂SO₄ and preservatives in an inorganic buffer compound. ProClean Contains a purple inert dye.

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Storage and stability

ProClean* is stable for two years at 18-30°C.

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Instructions for use

ProClean* should be used undiluted according to instrument manufacturer's instructions for use and should be connected as listed in the Operators manual. Reagent may be used with all kinds of Diluids* and CyMet's*.

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Avantor™ Performance Materials
Teugseweg 20 – 7418 AM Deventer – The Netherlands
Tel: +31 (0)570 687500
The devices as mentioned in this sheet comply with the
in Vitro Diagnostic Medical Device Directive 98/79/EG

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REF 3460.0500 CyMet * Erma III Diff 500 ml flacon HDPE

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Hypochlorite 0.5%

Intended use

Hypochlorite 0.5% is a specially filtered, non-sterile cleaning fluid for use in cleaning of cell counters.

The reagent is designed for semi-automated and automated instrumentation, capable to clean blood diluting parts of the instrument.

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Content

Hypochlorite 0.5% is water based and contains:

Sodium hypochlorite (0.5% active chlorine) and poly-oxy-ethylene-alkyl-alcohol.

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Storage and stability

Hypochlorite 0.5 % is stable for one year at 18-30°C.

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Instructions for use

Hypochlorite 0.5% should be used undiluted according to instrument manufacturer's instructions for use and should be connected as listed in the Operators manual.

Reagent may be used with all kinds of Diluids* and CyMet's*.

Pack size

REF 3917.1000 Hypochlorite 0.5% 1 liter bottle

REF 3917.5000 Hypochlorite 0.5% 5 liter bottle

* Trademark of Avantor™ Performance Materials - Deventer – The Netherlands

Avantor™ Performance Materials
Teugseweg 20 – 7418 AM Deventer – The Netherlands
Tel: +31 (0)570 687500
The devices as mentioned in this sheet comply with the
In Vitro Diagnostic Medical Device Directive 98/79/EG

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3 Diff Control pentru Coulter etc (cu opțiune III Diff.) series

373	low	2.5 ml	4	382	low	4.5 ml	4
4				0			
373	norm	2.5 ml	4	382	Norm	4.5 ml	4
5				1			
373	high	2.5 ml	4	382	High	4.5 ml	4
6				2			

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The reagent is designed for automated instrumentation, capable to monitor a three-part WBC differential, based on the aperture impedance principle and electronically adjusted to operate at an osmolality of 330 ± 20 mOsm/kg.

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Content: Diluid* III Diff is water based and contains:

NaCl, Na_2SO_4 , procaine HCl and preservatives in an inorganic buffer compound.

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Storage and stability: Diluid* III Diff is stable for three years at 18-30°C.

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Instructions for use

Diluid* Erma should be used undiluted according to instrument manufacturer's instructions for use and should be connected as listed in the Operators manual. Reagent may be used with Hypochlorite 0.5% or Proclean* as a cleaning agent. Furthermore reagent may be used with next lysing reagents: with Cymet BS and Detectoterge BS.

* Trademark of Avantor™ Performance Materials - Deventer –

DetectoTerge* BS

Intended use

DetectoTerge* BS is a specially filtered, non-sterile cleaning fluid for use in cleaning of cell counters. The product is designed for automated instrumentation, capable to clean blood diluting parts of the instrument. Intended to be used in vitro for the examination of specimens derived from the human body.

In daily routine DetectoTerge* BS is connected to the instrument as an on-line cleaner.

Between every sample the blood aspiration path is cleaned with DetectoTerge* BS. DetectoTerge* BS should be used on the BeneSphera 3 part differential hematology analyzer.

Summary and principle

The reagent is used to clean blood diluting parts prior to remove cell fragments from the instrument.

The solution is isotonic and contains surface active components to remove remaining blood cells in the fluid path. It also contains a harmless green dye visualising the reagent in the tubing of the instrument.

Content

DetectoTerge* BS is water based and contains:

Poly-oxy-ethylene-alkyl-alcohol, NaCl, Na₂SO₄ and preservatives in an inorganic buffer compound.

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Storage and stability

DetectoTerge* BS is stable for three years at 18-30°C.

Indications of deterioration

There are no visible deteriorations; otherwise the reagent shouldn't give optimised performance. The reagent can be used throughout shelf life, giving optimised performance. No guarantees to reagent performance are given after the expiry date.

Instructions for use

DetectoTerge* BS should be used undiluted according to instrument manufacturer's instructions for use and should be connected as listed in the Operators manual.

Reagent may be used with Diluid* III Diff and CyMet* BS3 CN Free.


Pack size

REF 2970.0900PE

DetectoTerge* BS

900 millilitres HDPE bottle

* Trademark of Avantor Performance Materials B.V. – Deventer – Holland

Avantor™ Performance Materials
Tegelseweg 20 – 7418 AM Deventer – The Netherlands
Tel. +31 (0)570 687500
The devices as mentioned in this sheet comply with the
In Vitro Diagnostic Medical Device Directive 98/79/EG

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Dia.Pro
Diagnostic
Bio**Probes**

Letter of Authorization

We, "Dia.Pro Diagnostic Bioprobes S.r.l." located at Via G. Carducci, Nr. 27 – Sesto San Giovanni (Milan) 20099, Italy, authorize

GLOBAL BIOMARKETING GROUP – MOLDOVA SRL
Str. Tighina 65, Oficiu 607
MD-2001 CHISINAU
REP. MOLDOVA


as our **exclusive distributor for the territory of the Republic of Moldova**, to participate in various tenders with **Dia.Pro ELISA** products.

We, Dia.Pro Diagnostic Bioprobes S.r.l shall supply our distributor GLOBAL BIOMARKETING GROUP – MOLDOVA SRL with all products in strict compliance with the existing "Distribution Agreement" rev.0117 valid until 31-Dec-2020, with possibility of renewal upon agreement between both parties for an additional period.

Dia.Pro Diagnostic Bioprobes S.r.l will grant the supply of all awarded tenders until their natural expiry, of which a documental proof has to be provided to Dia.Pro by the distributor GLOBAL BIOMARKETING GROUP – MOLDOVA SRL.

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Este certificado ampara todas las marcas de estos productos incluidas por el fabricante en su declaración de conformidad. / This certificate covers all trademarks of these products included by the manufacturer in his declaration of conformity.
Este certificado ampara todas las marcas de estos productos incluidas por el fabricante en su declaración de conformidad. / *This certificate covers all trademarks of these products included by the manufacturer in his declaration of conformity.*
**DIAGNOSTIC BIOPROBES S.r.l.**