



**WUHU RUIJIN MEDICAL INSTRUMENT
& DEVICE CO.,LTD**

Address/Tel/Fax: No.33, EastWanchun Road, Wuhu Economic
& Technology Development Zone, Wuhu Anhui China;
Tel: 86-553-2672511 Fax: 86-553-2672510
Website: www.whruijin.com



***Medical electric drill -
canulate drill RJ-PD (ND-2511)***

***Holds K-wires or be used in
intramedullary interlocking nail surgery
or arthroscopic surgery.***

» Feature

01 *Variable speed, convenient for use*

02 *Imported motor*

03 *High coaxiality, low noise*

04 *Can be autoclavable up to 135°C*



Composition »

*One handpiece
One charger
Two batteries
Two battery asptic transfer rings
One drill chuck key
One instruction manual
One Aluminum case*

» Main Technical Data

| | | |
|----------------------------|--------------------------|------------------------------------------------------------|
| Hand Piece | Output Power | 150W |
| | Unit rise in temperature | ≤25°C |
| Charger | Input Voltage | AC 110-220V , 50-60Hz |
| Battery | Voltage | 14.4V |
| | Battery capacity | 1800mAh |
| | Battery feature | Environmental friendly nickel-metal hydrogen battery |
| | Charging time | Two hours |
| Non-loading speed | 0-680/0-1000rpm | |
| Non-loading noise | ≤65dB (A) | |
| Axial run out | ≤0.15mm | |
| Rated torque | 5.5N.m/4N.m | |
| Impulsive torque | 10.5N.m/8.5N.m | |
| Drill chuck clamping range | 0.8-7.9mm | |
| Canulated diameter | ≤4.2mm | |

» Product details



» Accessories

Drill chuck & Key



Charger




Battery of 14. 4v



Aspetic battery transfer rings



Aluminum case packing



Sterilization case is additional cost.



*Can connect various flexible reamer,
k wire adapter and acetabular reamer,
changable during the process of surgery.*

Basic steps of sterilization, cleaning and maintenance of power tools.

- 1.Disaccessibly: All the consumable parts (saw blades, drill bits, etc) and battery should be removed before cleaning.*
- 2.Cleaning: Please wash and rinse with flowing water. Note: use soft brush to clean, then wipe with soft dry cloth. All instruments are prohibited from soak cleaning.*
- 3.Maintenance: Maintain by medical lubricant every time when cleaning and wiping instruments. Apply lubricant on every mobilizable part.
Note: Applying lubricant and rotate parts to lubricate to the full.*
- 4.Assembly: Assembly and functional test. Note: After cleaning and applying lubricants, install all parts one by one and conduct function test and then check its intactness.*
- 5.Sterilization: Forevacuum high temperature for 132 centigrade and high pressure within 20 minutes; natural cooling time is more than 40 minutes. Expoxymethane and formaldehyde are prohibited. Instruments must be kept dry before sterilization.*