#### - Precisão

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4. Schumann, G., *et al., Clin Chem Lab Med.*, (2002), **40**, 718.

5. Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples.* (2002). WHO/DIL/ LAB/99.1 Rev.2..

6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.

7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCC Press, (1997).

8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCC Press, (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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#### English - EN

#### INTENDED USE

ALT/GPT 4+1 SL is an *in vitro* diagnostic reagent intended for the quantitative determination of alanine aminotransferase (ALT) in human serum and plasma samples on analyzers or semi-automatic analyzers. This *in vitro* diagnostic device is for professional use only.

#### CLINICAL SIGNIFICANCE <sup>(1-3)</sup>

Alanine aminotransferase (ALT) is a transaminase also known as glutamate pyruvate transaminase (GPT). ALT catalyzes the transfer of the amino group of L-alanine to α-ketoglutarate to give L-glutamate. The highest levels are found in the liver and the kidneys. ALT levels are markedly increased in acute hepatitis (viral or toxic), and to a lesser extent in chronic hepatitis, cirrhosis, icterus, liver carcinoma or following various drugs. ALT levels can also be increased in heart diseases. ALT is more liver specific than AST (aspartate aminotransferase). Measurement of AST and ALT may assist in differentiating hepatitis from extrahepatic cell damage. ALT measurement is indicated to screen patients at risk of liver diseases and to help diagnose or monitor liver diseases.

#### LIMITATION OF USE

The quantitative assay of alanine aminotransferase (ALT) alone cannot be used to diagnose a disease or a specific pathology.

The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

#### METHOD & PRINCIPLE <sup>(4)</sup>

IFCC method without pyridoxal phosphate (P-5'-P). Kinetic.

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## BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Panteghini, M. & Blais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
- Dufour, D.R., *The Pancreas: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 651 and appendix.
- Winn-Deen, E.S., *et al. Clin. Chem.*, (1988), **34**, 2005
- Guder W.G, *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples*, (2002). WHO/DIL/LAB/99.1 Rev.2
- Schumann, G., *et al. Clin. Chem. Lab Med.*, (2006), **44**, 1146.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCPress, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCPress, (1995).

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
**- Special Programming Instructions :**  
see § INSTALLATION AND USE

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## English - EN

### INTENDED USE

AMYLASE SL is an *in vitro* diagnostic reagent intended for the quantitative determination of amylase in human serum and plasma samples on analyzers. This *in vitro* diagnostic device is for professional use only.

### CLINICAL SIGNIFICANCE <sup>(1-2)</sup>

$\alpha$ -Amylase is an enzyme mainly from pancreatic or salivary origin that hydrolyses 1,4- $\alpha$ -glycosidic bonds, thus helping for starch digestion. An increase of serum amylase can be observed in acute pancreatitis mainly, as well as in many other situations such as other pancreatic or intra-abdominal pathologies, and salivary gland diseases. During acute pancreatitis, serum amylase level starts rising after a few hours and can be increased up to 6-fold, peaking after 12 to 72h, and returning to normal level within about 4 days. In clinical practice, measurement of amylase activity is mainly used to help diagnose and monitor acute pancreatitis, and may also help diagnose and follow-up other pancreatic diseases.


### LIMITATION OF USE

The quantitative assay of amylase alone cannot be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

### METHOD & PRINCIPLE <sup>(3)</sup>

CNP-G<sub>3</sub> / Standardization against IFCC method – Kinetic

CNP-G<sub>3</sub> substrate is hydrolyzed by catalytic action of  $\alpha$ -amylase to produce CNP (2-chloro-4-nitrophenol) which absorbs at 405 nm.

CNP-G<sub>3</sub> = 2-Chloro-4-nitrophenyl- $\alpha$ -maltotrioside  
CNP-G<sub>2</sub> = 2-Chloro-4-nitrophenyl- $\alpha$ -maltoside

The rate of increase in absorbance is measured at 405 nm and is directly proportional to the activity  $\alpha$ -amylase in the sample.

### COMPOSITION

#### Reagent: R

Good's buffer, pH 6.15

CNP-G<sub>3</sub> 2.27 mmol/L

Sodium azide < 0.1 % (w/w)

Also contains calcium and chloride salts, as well as potassium thiocyanate for optimal performance.

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- Normal saline solution (NaCl 9 g/L).

- Analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

### PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.

- Reagent R contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Saliva and sweat contain amylase. It is therefore recommended to wear gloves and a mask to avoid the contamination of the reagent.

- Contact with acids liberates very toxic gas.

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## BIBLIOGRAPHIE/BIBLIOGRAPHY

### BIBLIOGRAFÍA/BIBLIOGRAFIA

- Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R. Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
- Dufour, D.R. *The liver: Function and Chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, A.H.B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 172.
- Howard, C.R., *Breastmilk & Breastfeeding Jaundice*, Pediatric Clinical Advisor: Instant Diagnosis and Treatment, 2<sup>nd</sup> Ed., Garfunkel, L.C., Kaczorowski, J.M. & Christy, C., (Mosby Inc. eds.), (2007), 82.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), 59, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press (1995).

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- Special Programming instructions : see § INSTALLATION AND USE

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## BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFÍA/BIBLIOGRAFIA

1. Wu, A. H. B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 306
2. Sanhai, W.R. *et al.*, *Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 677 and appendix.
3. Apple, F. S. & Jaffe, A. S., *Cardiovascular Disease, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 614.
4. Schumann, G., *et al.*, *Clin Chem Lab Med.*, (2002), **40**, 635.
5. Neumeier, D., *et al.*, *Clin Chim Acta.*, (1976), **73**, 445.
6. Klein, G., *et al.*, *Clin Chem.*, (2001), **47**, Suppl. A30.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCPress, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCPress, (1995).

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
- **Special Programming Instructions :** see § INSTALLATION AND USE

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
G-6-P : D-Glucose-6-Phosphate  
G-6-PDH : Glucose-6-Phosphate Déshydrogénase

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#### BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Kaplan, J.M., First, M.R., *Renal Function. Clinical Chemistry: Theory, Analysis, Correlation*. 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 567 and appendix.
- Delaney, M.P., Price, C.P., Lamb, E.J., *Kidney Function and Disease. Tietz Fundamentals of Clinical Chemistry*. 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 631.
- Wu, A. H. B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 316.
- Vasiliades, J., Clin. Chem., (1976), **22**, 1664.
- Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/ LAB/99.1 Rev.2.
- Ceriotti, F. *et al., Reference Intervals for Serum Creatinine Concentrations: Assessment of Available Data for Global Application, Clin\_Chem.*, (2008), **54**, 559.
- Young D.S., **Effects of drugs on clinical laboratory tests**, 4<sup>th</sup> Ed., AACCC Press (1995).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., **Effects of preanalytical variables on clinical laboratory tests**, 2<sup>nd</sup> Ed., AACCC Press (1997).

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### English - EN

#### INTENDED USE

CREATININE JAFFE is an *in vitro* diagnostic reagent intended for the quantitative determination of creatinine in human serum samples on analyzers or semi-automatic analyzers.

The standard is intended for the calibration of the reagent. These *in vitro* diagnostic devices are for professional use only.

#### CLINICAL SIGNIFICANCE <sup>(1-3)</sup>

Creatinine is the waste product of creatine, a compound present mainly in muscle tissues. It is usually produced a relatively constant rate and is freely filtered by the glomerulus, which makes it an excellent marker of renal function, with serum creatinine concentration increasing in case of acute or chronic renal insufficiency.

In practice, creatinine measurement in blood performed to help diagnose kidney diseases, to screen for early kidney failure for patients at risk, to monitor treatments against these pathologies or to monitor renal function during treatments that may impair this function.

#### LIMITATION OF USE

The quantitative assay of creatinine alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

#### METHOD & PRINCIPLE <sup>(4)</sup>

Jaffe – Kinetic.

In alkaline medium, creatinine reacts with picrate to form a red complex. The rate of increase of absorbance measured at 505 nm is proportional to creatinine concentration.

#### COMPOSITION

**Reagent 1: R1**  
Sodium hydroxide 0.31 mol/L  
Also contains phosphate  
**Reagent 2: R2**  
Picric acid 8.7 mmol/L  
**Standard: Std**  
Creatinine 2 mg/dL  
177 µmol/L

#### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

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## PRECAUTIONS FOR USE AND WARNINGS

- The reagent R1 is classified as hazardous :

- Consult Safety Data Sheet (SDS) for a proper handling.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent/standard vials from different kits.

### STABILITY

**Store at 2-25 °C and protect from light. Do not freeze.**

Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

**On board stability :**

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### PREPARATION

The reagent and standard are ready to use.

### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

### SAMPLES

#### Specimen <sup>(5)</sup>

- Serum.

- Using any other specimen type should be validated by the laboratory.

#### Warnings and precautions

Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

#### Storage and stability <sup>(5)</sup>

- 7 days at 20-25 °C.

- 7 days at 4-8 °C.

- 3 months at - 20 °C.

### REFERENCE VALUES <sup>(3,6)</sup>

<i>Serum</i>	mg/dL	µmol/L
Adults		
Men <span> </span> :	0.8 - 1.3	71 - 115
Women <span> </span> :	0.6 - 1.2	53 - 106

***Note :*** *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

### ❖INSTALLATION AND USE

**For use on Selectra Pro analyzers :**

- Consult operator manual.

- **Special Programming Instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

### PROCEDURE

<b>Manual Procedure</b>	
<i>Wavelength<span> </span>:</i>	<i>505 nm</i>
<i>Optical path<span> </span>:</i>	<i>1 cm</i>
<i>Sample/reagent ratio<span> </span>:</i>	<i>1:10</i>
<i>Temperature<span> </span>:</i>	<i>37<span> </span>°C</i>
<i>Read against distilled water.</i>	

<b>Working reagent (1 volume of R1 + 1 volume of R2)</b>	1000 µL
<b>Sample or Standard</b>	100 µL

*Mix and read the absorbance (A1) 10 seconds after the sample or standard addition. Exactly 2 minutes after the first reading, take second reading (A2).*

#### Automatic Procedure

These reagents may be used in several automatic analyzers. For Selectra Analyzers, validated applications are available on request.

### CALCULATION

(A2 - A1) Sample x n = standard concentration
(A2 - A1) Standard

**Conversion factor :** mg/dL x 88.40 = µmol/L
mg/dL x 10 = mg/L

### CALIBRATION

Creatinine Standard 2 mg/dL and ELICAL 2 are traceable to ID-MS (Isotope Dilution - Mass Spectrometry) reference method.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.
Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

### PERFORMANCES

Performances were obtained on Selectra Pro, following CLSI technical recommendations, under controlled environmental conditions.

#### - Measuring range

0.20 - 15.00 mg/dL (18 - 1326 µmol/L)

Do not report results outside the measuring range.

**Limit of Detection (LoD) and Limit of Quantification (LoQ)**

**LoD** = 0.05 mg/dL (4 µmol/L)

**LoQ** = 0.20 mg/dL (18 µmol/L)

#### - Precision

Imprecision data has been obtained on 2 Selectra Pro analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		<b>Mean</b>	<b>Within-run</b>	<b>Total</b>
	<b>n</b>	<b>mg/dL</b>	<b>µmol/L</b>	<b>CV (%)</b>
<b>Level 1</b>	80	0.95	84	1.1
<b>Level 2</b>	80	1.51	133	1.1
<b>Level 3</b>	80	4.98	440	0.9

#### - Correlation

A comparative study has been performed between CREATININE JAFFE reagent on a Selectra ProM analyzer and a similar commercially available system on 108 human serum samples.

The sample concentrations ranged from 0.16 to 14.80 mg/dL (14 - 1308 µmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.999

Linear regression: y = 0.939x + 0.14 mg/dL (12 µmol/L).

#### - Limitations/Interferences

- Antibiotics containing cephalosporin lead to significant false-positive values.<sup>(7)</sup>

- The presence of ketone bodies can cause artificially high results.<sup>(1)</sup>

- Studies have been performed to determine the level of interference from different compounds. The following creatinine levels were tested: 1.50 mg/dL and 5.00 mg/dL.

No significant interference is defined by a recovery ±10% of the initial value.

**Triglycerides:** No significant interference up to 2 100 mg/dL (23.7 mmol/L).

**Unconjugated bilirubin:** No significant interference up to 15.0 mg/dL (257 µmol/L).

**Conjugated bilirubin:** No significant interference up to 5.9 mg/dL (101 µmol/L).

**Hemoglobin:** No significant interference up to 500 mg/dL.

**Glucose :** No significant interference up to 500 mg/dL (27.8 mmol/L).

**Creatine:** No significant interference up to 10.0 mg/dL (762 µmol/L).

**Uric acid** : No significant interference up to 20.0 mg/dL (1190 µmol/L).

**L-dopa :** No significant interference up to 20.0 mg/dL.

**Áscorbic Ácido:** No significant interference up to 20 mg/dL.

**Acetaminophen :** No significant interference up to 30 mg/dL.

**Acetylsalicylic acid :** No significant interference up to 200 mg/dL.

**Methyl dopa :** No significant interference up to 20 mg/dL.

**Pyruvate :** No significant interference up to 3.0 mg/dL (341 µmol/L)

**Proteins :** No significant interference between 3.0 g/dL and 11.0 g/dL.

- Do not use icteric samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.<sup>(8)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(7,9)</sup>

#### - On board stability/Calibration frequency

**On Board Stability:** 2 days

**Calibration frequency:** At the end of onboard stability period,change the bottle with a new one filled with fresh reagent.

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

*These performances have been obtained using Selectra Pro analyzer. Results may vary if a different instrument or a manual procedure is used.*

*❖The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.*

### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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## English - EN

### INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems quantitative tests listed in the value sheet.  
This *in vitro* diagnostic device is for professional use only.

### COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.  
- Sodium azide < 0.1 % (w/w)  
- Concentrations for each analyte to test are lot-specific.

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### TRACEABILITY

The traceability is indicated in the value sheet enclosed in the kit.

### PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved or cleared in compliance with European Directive 98/79/EC, Annex II, List A. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.  
- Take precautions when handling broken glass vials as sharp edges can injure the user.  
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.  
- Take normal precautions and adhere to good laboratory practice.  
- Use clean or single use laboratory equipment only to avoid contaminations.  
- Consult Safety Data Sheet (SDS) for a proper handling.

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### PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.  
- Add in exactly 3 mL of distilled or deionized water.  
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

### PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.  
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).  
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

### STABILITY

#### Prior to reconstitution :

- Stored at 2-8 °C and protected from light.  
- Do not use after expiration date indicated on the vial label.

#### After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.

- Stability of the components:  
Between 15 and 25 °C : 8 hours  
Between 2 and 8 °C : 2 days  
Between -25 and -15 °C : 4 weeks (when frozen once)

#### Exceptions:

- Stability of total bilirubin (when stored protected from light):  
Between 15 and 25 °C : 6 hours  
Between 2 and 8 °C : 1 day  
Between -25 and -15 °C : 2 weeks (When frozen once)  
- Stability of direct bilirubin (when stored protected from light):  
Between 15 and 25 °C : 3 hours  
Between 2 and 8 °C : 8 hours  
Between -25 and -15 °C : 2 weeks (when frozen once)

**Note:** Store vials tightly capped and protected from light.

### PROCEDURE

To use ELICAL 2, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

### LIMITATIONS

Using ELICAL 2 has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet).

Using any other system should be validated by the laboratory.

### VALUES

The concentrations are indicated in the value sheet enclosed in the kit.  
For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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### English - EN

### INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of ELITech Clinical Systems quantitative tests listed in the value sheet. These *in vitro* diagnostic devices are for professional use only.

### COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.  
- Concentrations for each analyte to test are lot-specific.

### MATERIALS REQUIRED BUT NOT PROVIDED

- ELITech Clinical Systems reagents listed in the value sheet.  
- General Laboratory equipment.

### TRACEABILITY

The traceability is indicated on the instructions for use of the ELITech Clinical Systems reagent used in combination with the recommended calibrator.

### PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.  
- Take precautions when handling broken glass vials as sharp edges can injure the user.  
- Take normal precautions and adhere to good laboratory practice.  
- Use clean or single use laboratory equipment only to avoid contaminations.  
- Consult Safety Data Sheet (SDS) for a proper handling.

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### PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.  
- Carefully open the vial, avoiding the loss of lyophilizate.  
- Add exactly 5 mL of distilled or deionized water.  
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

*Important: To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.*

### PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.  
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).  
- Damage to the container may impact on product

performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

### STABILITY

#### Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.  
- Do not use after the expiry date stated on the vial label.

#### After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.  
- Stability of the components:  
Between 2 and 8 °C : 5 days  
Between -25 and -15 °C : 4 weeks (when frozen once)

#### Exceptions:

- Stability of total bilirubin (when stored protected from light):  
Between 2 and 8 °C : 24 hours  
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of direct bilirubin (when stored protected from light):  
Between 2 and 8 °C : 8 hours  
Between -25 and -15 °C : 5 days (when frozen once)

#### Stability of alkaline phosphatase:

Between 2 and 8 °C : 24hours  
Between -25 and -15 °C : 4 weeks (when frozen once)

#### Note:

- Store the vials tightly capped and protected from light after reconstitution.  
- Mix content thoroughly prior to use.  
- Thawing frozen specimen at 37°C, and vigorous mixing (vortex), should be avoided as this may denature ferritin.

### PROCEDURE

To use ELITROL I & II, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

### LIMITATIONS

Using ELITROL I & II has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet). Using any other system should be validated by the laboratory.

### VALUES

The concentrations are indicated in the value sheet enclosed in the kit. Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values. For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

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### ASSISTANCE TECHNIQUE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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- Symbols used are defined on ISO-15223-1 standard, except those presented below.

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PIT-SLSY-EN-V1 (05/2022)

### INTENDED USE

This solution is intended for washing probes and cuvettes of Selectra Pro Series analyzers.

This accessory for an *in vitro* diagnostic device is for professional use only.

### COMPOSITION

Sodium azide < 0.1 % (w/w)  
Also contains surfactants

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### PRECAUTIONS OF USE AND WARNINGS

- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Solution contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of this solution always flush with copious amounts of water to prevent azide buildup.
- Consult Safety Data Sheet (SDS) for a proper handling.

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Do not use after expiration date indicated on the vial labels.

Diluted solution stability: 2 weeks.

Open vial stability: until product's expiry date.

### PREPARATION

Solution needs to be diluted 1/400 in distilled or deionized water, as indicated in the analyzer's user manual.

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### INSTALLATION AND USE

Consult Selectra Pro operator manual corresponding.

### LIMITATIONS

The use of SYSTEM SOLUTION has been validated for the Selectra Pro Series analyzers.

Using a different analyzer should be validated by the laboratory.

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	Content
	Solution
	Modification from previous version
	European Conformity

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PIT-SOL-EN-v1 (05/2022)

### INTENDED USE

These solutions are intended for washing Selectra Pro Series analyzers.  
These accessories for an *in vitro* diagnostic device are for professional use only.

### COMPOSITION

**Acid Solution:** Acidic solution (pH 1).  
**System Cleaning Solution:** Alkaline sodium hypochlorite solution.

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- Distilled or deionized water, as indicated in the analyzer's user manual.
- General Laboratory equipment (e.g. pipette).
- Selectra Pro Family analyzers and accessories.
- Do not use materials that are not required as indicated above.

### PRECAUTIONS OF USE AND WARNINGS

- **Acid solution** and **System Cleaning Solution** are classified as hazardous

- Consult Safety Data Sheet (SDS) for a proper handling.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.

### STABILITY

**Store at 15-25 °C and protect from light. Do not freeze.**

Do not use after expiration date indicated on the vial labels.  
Open vial stability for both solutions : until product's expiry date.

On board-stability:

**Acid Solution** (pure) : 2 weeks  
**System Cleaning Solution** (diluted) : 2 weeks

### PREPARATION

**Acid Solution** is ready to use.  
**System Cleaning Solution** needs to be diluted 1/10 in distilled or deionized water, as indicated in the analyzer's user manual.

### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate product deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the container may impact on product performance. Do not use the solution if there is physical evidence of deterioration (e.g. leakages or punctured container).

### INSTALLATION AND USE

**Special programming instructions:** Verify if the Probe and Test incompatibilities settings in the Selectra Pro instrument are set as indicated in the Special Programming Instructions table on page 2 of this instructions for use.

Consult the Selectra Pro operator manual on how to program these Probe and Test incompatibilities.

### LIMITATIONS

The use of ACID SOLUTION and SYSTEM CLEANING SOLUTION has been validated for the Selectra Pro analyzer.  
Using a different analyzer should be validated by the laboratory.

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### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements.

By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

### TECHNICAL ASSISTANCE

Contact your local distributor or ELITech Clinical Systems SAS (ccsupport@elitechgroup.com).

### SYMBOLS

Symbols used are defined in ISO 15223-1 standard, except those presented below :

<b>CONT</b>	Content
<b>SOL</b>	Solution
	Modification from previous version
<b>CE</b>	European Conformity

## SPECIAL PROGRAMMING INSTRUCTIONS TABLE

Program the Probe & Test incompatibilities in the Selectra Pro instrument as indicated in the table below

### TOUCHPRO SOFTWARE

#### Probe Incompatibilities

TEST	Type	Reagent
CHOLESTEROL HDL 2G SL	Link	Acid Solution
HDL CHOLESTEROL	Link	Acid Solution
VITAMIN D	Link	Acid Solution
TOTAL PROTEIN PLUS	Incompatibility	MAGNESIUM XYLIDYL
CHOLESTEROL SL	Incompatibility	MAGNESIUM XB

#### Test Incompatibilities

TEST	Type	Reagent
TRIGLYCERIDES SL	Link	Acid Solution*
VITAMIN D	Link	Acid Solution

\*To be programmed only if LIPASE reagent is installed on-board

LIPASE and TRIGLYCERIDES must be run in separate routines.

### OTHER SOFTWARE

#### Needle Incompatibility

CHOLESTEROL HDL 2G SL	Link <<	HCl
HDL CHOLESTEROL	Link <<	HCl
VITAMIN D	Link <<	HCl
TOTAL PROTEIN PLUS	Incompatibility :	MAGNESIUM XYLIDYL
CHOLESTEROL SL	Incompatibility :	MAGNESIUM XB

#### Cuvette Incompatibility

TRIGLYCERIDES SL	Link <<	HCl*
VITAMIN D	Link <<	HCl

\*To be programmed only if LIPASE reagent is installed on-board

LIPASE and TRIGLYCERIDES must be run in separate routines.