

Declaration of Conformity

EUROIMMUN Medizinische Labordiagnostika AG
Seekamp 31, D-23560 Lübeck, Germany

declare under our sole responsibility that the ELISA products

Anti-Borrelia ELISA (IgG)

EI 2132-9601 G

Anti-Borrelia ELISA (IgM)

EI 2132-9601 M

CSF: Anti-Borrelia ELISA (IgM)

EI 2132-9601-L M

(product name, order no)

meet the demands of

Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998

Lübeck, 25.04.2008

(Place and date of issue)

Wolfgang Schlumberger, PhD
- Member of the Board -

Susanne Aleksandrowicz
- Member of the Board -



Anti-Borrelia plus VISE ELISA (IgG)
Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EL 2132-9601-2 G	Borrelia burgdorferi, Borrelia afzelii, Borrelia garinii, Borrelia VISE Antigen	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative or quantitative *in vitro* assay for human antibodies of the IgG class against Borrelia antigens in serum or plasma for the diagnosis of infection with Borrelia and associated diseases: Erythema chronicum migrans, Lyme-disease, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculitis and neuro-borreliosis.

Application: Clinical diagnosis of borreliosis can be achieved by determination of Borrelia-specific antibodies of classes IgG and IgM. The Anti-Borrelia plus VISE ELISA (IgG) is based on an optimised lysate mixture of the most relevant human pathogenic Borrelia strains and also contains recombinant VISE. Due to its complete antigen spectrum, the ELISA offers a very high sensitivity and is therefore ideally suited for antibody screening. Samples with positive or borderline ELISA results should always be further investigated using immunoblot.

Principle of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with a mix of whole antigen extracts of Borrelia burgdorferi sensu stricto, Borrelia afzelii, Borrelia garinii and recombinant VISE of Borrelia burgdorferi. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	[STRIPS]
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	[CAL 1]
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	[CAL 2]
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	[CAL 3]
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	[POS CONTROL]
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	[NEG CONTROL]
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	[CONJUGATE]
8. Sample buffer ready for use	light blue	1 x 100 ml	[SAMPLE BUFFER]
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	[WASH BUFFER 10x]
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	[SUBSTRATE]
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	[STOP SOLUTION]
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
[LOT] Lot description In vitro diagnostic medical device	CE		Storage temperature Unopened, usable until

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C), approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- Coated wells: Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bags).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.

- Calibrators and controls: Ready for use. The reagents must be mixed thoroughly before use.

- Enzyme conjugate: Ready for use. The enzyme conjugate must be mixed thoroughly before use.

- Sample buffer: Ready for use.

- Wash buffer: The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).

For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.

- Chromogen/substrate solution: Ready for use. Close the bottle immediately after use, as the contents are sensitive to light ☼. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.

- Stop solution: Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBSAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples must be incubated within one working day.

Sample dilution: Patient samples are diluted 1:101 with sample buffer.
Example: Add 10 µl of sample to 1.0 ml sample buffer and mix well by vortexing (sample cassettes are not suitable for mixing).

NOTE: The calibrators and controls are pre-diluted and ready for use, do not dilute them.



Incubation

For semiquantitative analysis: incubate calibrator 2 along with the positive and negative controls and patient samples. For quantitative analysis: incubate calibrators 1, 2 and 3 along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation:
(1st step)

Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.
Incubate for 30 minutes at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values. Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells.
Incubate for 30 minutes at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation:
(3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells.
Incubate for 15 minutes at room temperature (+18°C to +25°C), protect from direct sunlight.

Stopping:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm within 30 minutes of adding the stop solution. Prior to measuring slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices:

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer 1, Analyzer 1-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry. Automated test performance using other fully automated, open system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the semiquantitative analysis of 24 patient sample (P 1 to P 24).
The pipetting protocol for microtiter strips 7-10 is an example for the quantitative analysis of 24 patient sample (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.
The wells can be broken off individually from the strips. Therefore, the number of tests performed can be matched to the number of samples, minimising reagent wastage.
Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of the calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

- Ratio <0.8: negative
- Ratio ≥0.8 to <1.1: borderline
- Ratio ≥1.1: positive



Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as ">200 RU/ml". It is recommended that the sample be retested at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4. The upper limit of the normal range of non-infected persons (cut-off value) recommended by EUROIMMUN is 20 relative units (RU)/ml. EUROIMMUN recommends interpreting results as follows:

<16 RU/ml: negative
 ≥16 to <22 RU/ml: borderline
 ≥22 RU/ml: positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends to retest the samples.

A negative serological result does not exclude an infection. Particularly in the early phase of an infection, antibodies may not yet be present or are only present in such small quantities that they are not detectable. In case of a borderline result, a secure evaluation is not possible. If there is a clinical suspicion and a negative test result, we recommend clarification by means of other diagnostic methods and/or the serological investigation of a follow-up sample. A positive result indicates that there has been contact with the pathogen. In the determination of pathogen-specific IgM antibodies, polyclonal stimulation of the immune system or antibody persistence may affect the diagnostic relevance of positive findings. Significant titer increases (exceeding factor 2) and/or seroconversion in a follow-up sample taken after 7 to 10 days can indicate an acute infection. To investigate titer changes, sample and follow-up sample should be incubated in adjacent wells of the ELISA microplate within the same test run. For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against Borrelia, the calibration is performed in relative units (RU)/ml.

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.



The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is provided by particularly suitable Borrelia strains (Borrelia burgdorferi sensu stricto, Borrelia atzeili, Borrelia garinii) as well as recombinant VlsE (variable major protein-like sequence, expressed) of Borrelia burgdorferi. VlsE is a newly characterised surface protein of Borrelia which is expressed exclusively in vivo and which contains conserved and highly immunogenic epitopes. The cultured borrelia have been solubilised by using sodium dodecyl sulphate. The used antigen mixture contains all relevant proteins.

Linearity: The linearity of the Anti-Borrelia plus VlsE ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R² for all sera was > 0.95. The Anti-Borrelia plus VlsE ELISA (IgG) is linear at least in the tested concentration range (8 RU/ml to 124 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-Borrelia plus VlsE ELISA (IgG) is 0.7 RU/ml.

Cross reactivity: Cross reactivity of the Anti-Borrelia plus VlsE ELISA (IgG) was evaluated in a study performed on 256 patient sera with antibodies against Treponema pallidum and 16 patient sera with antibodies against Leptospira. 95 sera of the Anti-Treponema pallidum-positive sera and 1 serum of the Anti-Leptospira-positive sera presented reactivity. Therefore, the cross reactivity amounts to 37.1% and 6.3%, respectively. Cross reactions with antibodies against other spirochaetes cannot be excluded. The Anti-Borrelia plus VlsE ELISA is designed as a screening test which should be followed by a confirmatory test (Western- or line blot). Further sera from patients with different possible influencing factors were investigated. An overview of these results can be found in the following table.

Possible influencing factors	n	Anti-Borrelia plus VlsE ELISA (IgG) positive
Acute EBV infection	33	0%
Anti-HSV 1	12	0%
Anti-CMV	12	0%
Anti-VZV	12	0%
Anti-Adenovirus	12	0%
Anti-RSV	12	0%
Anti-Parainfluenza	12	0%
Anti-Influenza-A virus	12	0%
Anti-Influenza-B virus	12	0%
Anti-Mycoplasma pneumoniae	12	0%
Anti-Measles virus	12	0%
Anti-Mumps virus	12	0%
Anti-Rubella virus	12	0%
Anti-Toxoplasma gondii	11	0%
Anti-Chlamydia pneumoniae	11	0%
Anti-Helicobacter pylori	12	0%
ANA + DNS (AAb)	36	2.8%
Rheuma factor	37	0%
Neurological diseases	54	3.7%



Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 40 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The Intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 8 different test runs.

Intra-assay variation, n = 20		
Serum	Mean value (RU/ml)	CV (%)
1	52	3.5
2	97	3.5
3	128	3.6

Inter-assay variation, n = 4 x 6		
Serum	Mean value (RU/ml)	CV (%)
1	55	3.6
2	103	3.7
3	138	3.8

Specificity and sensitivity: Sera of 165 patients with suspected borreliosis were analysed using the EUROIMMUN Anti-Borrelia plus VISE ELISA (IgG) and the EUROIMMUN Anti-Borrelia EUROLINE-Westernblot (IgG) as a reference method. The test showed a specificity of 90.2% and a sensitivity of 100%.

ELISA	EUROIMMUN EUROLINE-Westernblot	
	positive	borderline
positive	60	11
borderline	0	1
negative	0	0
n = 165		negative
		9
		1
		83

Borderline results are not included by the calculation of specificity and sensitivity.

Clinical study: 138 patients with clinically characterised borreliosis in different disease stages were screened with the Anti-Borrelia plus VISE-ELISA (IgG) and Anti-Borrelia-ELISA (IgM). Sensitivities between 89 to 96% were found.

Clinic	EUROIMMUN Anti-Borrelia plus VISE-ELISA (IgG)		IgG and/or IgM positive
	IgG positive	IgM positive	
erythema migrans, n = 97	60	68	86 (88.7%)
	22	5	25 (96.2%)
	13	4	14 (93.3%)
mono- or poly- arthritis, n = 26			
neuroborreliosis, n = 15			
n = 138			

All ELISA negative results were confirmed as negative by an Anti-Borrelia Westernblot.

Reference range: The levels of the anti-Borrelia antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 20 RU/ml, 5% of the blood donors were anti-Borrelia positive (IgG), which reflects the known percentage of infections in adults.



Clinical significance

The history of Lyme disease, a contagious condition caused by *Borrelia burgdorferi* and transmitted to man by ticks, offers infectious a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis. Fifty years later, Häuser noted the affection was transmitted by ticks. Independently, Atzeius, then Lipschütz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemia of arthritis occurred in children in Lyme, Connecticut, USA. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemia, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that ticks transmitted "reptonema-like spirochaetes", which were later identified as the causal agent: *Borrelia burgdorferi*.

Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaeaceae family. This family includes five genera: *borrelia*, *spirochaetes*, *crispspira*, *treponema* and *leptospira*. Among about 30 tick species feeding on humans, *Ixodes ricinus* is the most frequent tick species biting humans in Europe. It is the vector of *Borrelia burgdorferi*, which causes Lyme disease, of anaplasmata phagocytophilum and the tick-borne encephalitis virus. *I. ricinus* ticks pass through three developmental stages: larvae, nymphs and adults (females and males). The density of this tick species can be very high, reaching in some places more than 300 ticks/100 m². Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100,000 individuals caused by the species *B. burgdorferi sensu stricto*, *B. atzeius*, and *B. garinii* (in rare cases also by *B. spielmanni* and *B. bavariensis*). The human seroprevalence rate of antibodies against *Borrelia burgdorferi* in the normal population of Germany and other Central European countries is about 8% (for IgG); in highly endemic areas even considerably higher with ranges from 18% to 52%. In East Asia, e.g. in an endemic area in China, the seroprevalence amounts to 26%. Persons working in the area of forestry display anti-Borrelia antibodies in about 40%, hunters in more than 50% of cases.

B. Borrelia burgdorferi infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a *Borrelia burgdorferi* infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting. Lymphadenopathies are observed in a few cases (lymphadenitis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular borreliosis and borreliosis with hepatitis, with regional lymphadenitis and mixtures of these. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against *Borrelia burgdorferi* can be detected serologically in 50% to 90% of patients during stage I. Humans produce highly specific antibodies against the outer surface protein C (OspC) shortly after infection with *Borrelia burgdorferi*. The main antigen VISE, which is exclusively expressed *in vivo*, is the most sensitive antigen for IgG antibody detection, whereas OspC is the most sensitive antigen for IgM antibody detection.

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears on one side and after several weeks on the opposite side. Also arthritis particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found. Less frequently are cardiological manifestations such as myocarditis and pericarditis. Antibodies against *Borrelia burgdorferi* can be detected in 50% to 90% of patients in stage II. In the early phase of this stage mainly IgM antibodies are found but in the late phase often only IgG antibodies occur. However, IgM antibodies can persist for a long time. By additionally determining antibodies against VISE the serological hit rate can be increased by 20% to 30%. VISE displays the highest sensitivity of all antigens tested for IgG detection.



Stage III: The typical manifestations of a *Borrelia burgdorferi* infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibody titers are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable.

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing. An antibody index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in the blood/CSF barrier function.

Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi*: ELISA, indirect immunofluorescence (IIFT), passive haemagglutination and immunoblot. Generally, ELISA or IIFT are used as screening test for preliminary characterisation of the serum sample. According to guidelines in the USA and in Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIFT) as a first step, which, if reactive, is followed by an immunoblot. In fresh infections it is recommended performing ELISA/IIFT and immunoblot in parallel, since some weak reactions become detectable earlier in the blot than in the screening test. VISE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection should be included in the test. By additionally determining antibodies against VISE the serological hit rate can be increased by 20% compared to whole extract Westernblots. Of all recombinant antigens tested, VISE possesses the highest sensitivity for the detection of a *Borrelia* infection.

Literature references

- Burgdorfer W, Barbour AG, Hayes SF, Benach JL, Granwaldt E, Davis JP. Lyme disease – a tick-borne Spirochetosis? *Science* 216 (1982) 1317-1319.
- Christova I. Enzyme-linked immunosorbent assay, immunofluorescent assay, and recombinant immunoblotting in the serodiagnosis of early Lyme borreliosis. *Int J Immunopathol Pharmacol* 16 (2003) 261-268.
- EUROIMMUN AG, Komorowski L, Janssen A, Stöcker W, Probst C. Dimerisation of recombinant OspC leads to an antigen with enhanced potential for active vaccination. 12. International Conference on Lyme Borreliosis and other Tick-Borne Diseases, Ljubljana, Slovenia (2010).
- EUROIMMUN AG, Komorowski L, Probst C, Janssen A, Stöcker W. Polypeptide und Verfahren zur spezifischen Detektion von Antikörpern bei Patienten mit einer Borrelien-Infektion, Europäische Patentanmeldung EP2199303 (2008).
- EUROIMMUN AG, Meyer W, Janssen A, Schepfer T, Schlumberger W, Stöcker W. Lyme borreliosis: Prevalence of antibodies against non-proteinic (lipid) antigens. *Int J Med Microbiol* 298S2 (Suppl. 45) (2008) 39.
- EUROIMMUN AG, Probst C, Ott A, Schepfer T, Meyer W, Stöcker W, Komorowski L. N-Terminal disulfide-bridging of *Borrelia* outer surface protein C increases its diagnostic and vaccine potentials. *Ticks and Tick-Borne Diseases* 3 (2012) 1-7.
- EUROIMMUN AG, Meyer W, Schepfer T, Schlumberger W, Stöcker W. Antibodies Against the Newly Identified Major Antigen Vise: A Milestone in the Serological Diagnosis of Lyme Borreliosis. Poster zum EUREGIO Congress of Clinical Chemistry and Laboratory Medicine, Aachen (2003).
- EUROIMMUN AG, Steinhagen K, Schlumberger W, Stöcker W. Nachweis einer spezifischen Intrastrukturellen Antikörpersynthese mit modernen ELISA-Testsystemen: Hohe Trefferrate bei Multipler Sklerose und Neuroborreliose. *J Lab Med* 25 (2001) 135-149.
- EUROIMMUN AG, Meyer W, Janssen A, Kühn B, Schepfer T, Suer W, Komorowski L, Probst C, Schlumberger W, Stöcker W. EUROLINE *Borrelia* RN-AT: Marriage of authentic recombinant and native *Borrelia* antigens in a line immunosassay format. *International Journal of Medical Microbiology* 298S2 (Suppl. 45): 39-40 (2008).
- Goettler G, Schulte-Spechtel U, Hillermann R, Liegl G, Wilske B, Fingerle V. Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of Vise and DbpA homologues. *J Clin Microbiol* 43 (2005) 3602-3609.
- Hauser U, Lehnert G, Lobenzanzer R, Wilske B. Interpretation Criteria for Standardized WesternBlots for Three European Species of *Borrelia burgdorferi* Sensu Lato. *J Clin Microbiol* 35 (1997) 1433-1444.
- Johnson R, Schmid GP, Hyde FW, Steigenwalt AG, Brenner DJ. *Borrelia burgdorferi* sp. Nov.: Etiologic agent of Lyme disease. *Int J Syst Bacteriol* 34 (1984) 496-497.
- Oehme R, Hartelt K, Backe H, Brockmann S, Kimmig P. Foci of tick-borne diseases in south-west Germany. *Int J Med Microbiol* 291 (2002) 22-29.
- Robert-Koch-Institut. Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg. *Epidemiologisches Bulletin* 14 (1998) 93-97.
- Schulte-Spechtel U, Lehnert G, Liegl G, Fingerle V, Heimert C, Johnson BJ, Wilske B. Significant improvement of the recombinant *Borrelia*-specific immunoglobulin G immunoblot test by addition of Vise and a DbpA homologue derived from *Borrelia garinii* for diagnosis of early neuroborreliosis. *J Clin Microbiol* 41 (2003) 1299-1303.
- Tilton RC, Ryan RW. The laboratory diagnosis of Lyme disease. *J Clin Immunossay* 16 (1993) 208-214.
- Wilske B, Fingerle V, Schulte-Spechtel U. Microbiological and serological diagnosis of Lyme borreliosis. *FEMS Immunol Med Microbiol* 49 (2007) 13-21.
- Wilske B. Epidemiology and diagnosis of Lyme borreliosis. *Ann Med* 37 (2005) 568-579.
- Wilske B, Zöller L, Brade V, Effert H, Göbel UB, Stanek G, Pfister HW, Miq 12/2000 Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik: Lyme-Borreliose. Urban & Fischer Verlag München Jena (2000).
- Zhang JR, Hardham JM, Barbour AG, Norris SJ. Antigenic Variation in Lyme Disease *Borreliae* by Promiscuous Recombination of VMP-like Sequence Cassettes. *Cell* 89 (1997) 275-285.



Anti-Borrelia ELISA (IgM) Test instruction


ORDER NO	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2132-9601 M	Borrelia burgdorferi, Borrelia atzeltii, Borrelia garinii	IgM	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative or quantitative *in vitro* assay for human antibodies of the IgM class against Borrelia antigens in serum or plasma for the diagnosis of Infection with Borrelia. Associated diseases: Erythema chronicum migrans, lymphadenitis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculitis, neuroborreliosis.

Application: Clinical diagnosis of borreliosis can be achieved by determination of Borrelia-specific antibodies of classes IgG and IgM. The Anti-Borrelia ELISA (IgM) is based on an optimised lysate mixture of the most relevant human pathogenic Borrelia strains. Specific cultivation methods help to ensure a large proportion of OspC in the test. Due to its complete antigen spectrum, the ELISA offers a very high sensitivity and is therefore ideally suited for antibody screening. Samples with positive or borderline ELISA results should always be further investigated using immunoblot.

Principle of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with antigen extracts of Borrelia burgdorferi sensu stricto, Borrelia atzeltii and Borrelia garinii. In the first reaction step, diluted patient samples are incubated with the wells. In the case of positive samples, specific IgM antibodies (also IgA and IgG) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgM (enzyme conjugate), a catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	STRIPS
2. Calibrator 1, 200 RU/ml (IgM, human), ready for use	dark red	1 x 2,0 ml	CAL 1
3. Calibrator 2, 20 RU/ml (IgM, human), ready for use	red	1 x 2,0 ml	CAL 2
4. Calibrator 3, 2 RU/ml (IgM, human), ready for use	light red	1 x 2,0 ml	CAL 3
5. Positive control, (IgM, human), ready for use	blue	1 x 2,0 ml	POS CONTROL
6. Negative control, (IgM, human), ready for use	green	1 x 2,0 ml	NEG CONTROL
7. Enzyme conjugate peroxidase-labelled anti-human IgM (goat), ready for use	red	1 x 12 ml	CONJUGATE
8. Sample buffer containing IgG/RF absorbent (anti-human IgG anti-body preparation obtained from goat), ready for use	green	1 x 100 ml	SAMPLE BUFFER
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	WASH BUFFER 10X
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	SUBSTRATE
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	STOP SOLUTION
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
LOT In vitro diagnostic medical device	 Storage temperature Unopened usable until		

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

Coated wells: Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature (+18°C to +25°C) to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bags). Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.

Calibrators and controls: Ready for use. The reagents must be mixed thoroughly before use.

Enzyme conjugate: Ready for use. The enzyme conjugate must be mixed thoroughly before use.

Sample buffer: Ready for use. The green coloured sample buffer contains IgG/RF absorbent. Serum or plasma samples diluted with this sample buffer are only to be used for the determination of IgM antibodies.

Wash buffer: The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).

For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water. The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.

Chromogen/substrate solution: Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.

Stop solution: Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.





Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Introduction: Before the determination of specific antibodies of class IgM, antibodies of class IgG should be removed from the patient sample. This procedure must be carried out in order to prevent any rheumatoid factors from reacting with specifically bound IgG, which would lead to false positive IgM test results, and to prevent that specific IgG displace IgM from the antigen, which would lead to false IgM negative test results.

Functional principle: The sample buffer (green coloured) contains an anti-human antibody preparation from goat. IgG of a serum or plasma sample is bound with high specificity by these antibodies and precipitated. If the sample also contains rheumatoid factors, these will be absorbed by the IgG/anti-human IgG complex.

Separation properties:

- All IgG subclasses are bound and precipitated by the anti-human IgG antibodies.
- Human serum IgG in concentrations of up to 15 mg per ml are removed (average serum IgG concentration in adults: 12 mg per ml).
- Rheumatoid factors are also removed.
- The recovery rate of the IgM fraction is almost 100%.

Performance: The patient samples for analysis are diluted 1:101 with green coloured sample buffer. For example, add 10 µl sample to 1.0 ml sample buffer and mix well by vortexing. Sample pipettes are not suitable for mixing. Incubate the mixture for at least 10 minutes at room temperature (+18°C to +25°C). Subsequently, it can be pipetted into the microplate wells according to the pipetting protocol.

Notes:

- Antibodies of the class IgG should not be analysed with this mixture.
- It is possible to check the efficacy of the IgG/RF absorbent for an individual patient sample by performing an IgG test in parallel to the IgM test using the mixture. If the IgG test is negative, the IgM result can be considered as reliable.
- The calibrators and controls are ready for use, do not dilute them.

For semiquantitative analysis: incubate calibrator 2 along with the positive and negative controls and patient samples. For quantitative analysis: incubate calibrators 1, 2 and 3 along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation:
(1st step)

Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for 30 minutes at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g. less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values. Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgM) into each of the microplate wells. Incubate for 30 minutes at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation:
(3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for 15 minutes at room temperature (+18°C to +25°C), protect from direct sunlight.

Stopping:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm within 30 minutes of adding the stop solution. Prior to measuring, carefully shake the microplate to ensure a homogeneous distribution of the solution.



Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.
Automated test performance using other fully automated, open system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

A	C 2	P 6	P 14	P 22										
B	pos	P 7	P 15	P 23										
C	neg	P 8	P 16	P 24										
D	P 1	P 9	P 17											
E	P 2	P 10	P 18											
F	P 3	P 11	P 19											
G	P 4	P 12	P 20											
H	P 5	P 13	P 21											

The pipetting protocol for microtiter strips 1-4 is an example of the semiquantitative determination of antibodies in 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example of the quantitative determination of antibodies in 24 patient samples (P 1 to P 24).

Calibrators (C 1 to C 3), positive (pos.) and negative (neg.) control as well as the patient samples have been incubated in one well each. The reliability of the ELISA test can be improved by duplicate determinations of each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimises reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio < 0.8	negative
Ratio ≥ 0.8 to < 1.1	borderline
Ratio ≥ 1.1	positive



Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibrators against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as ">200 RU/ml". It is recommended that the sample be retested at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4. The upper limit of the reference range of non-infected persons (cut-off value) recommended by EUROIMMUN is 20 relative units RU/ml. EUROIMMUN recommends interpreting results as follows:

<16 RU/ml:	negative
≥16 to <22 RU/ml:	borderline
≥22 RU/ml:	positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends to retest the samples.

A negative serological result does not exclude an infection. Particularly in the early phase of an infection, antibodies may not yet be present or are only present in such small quantities that they are not detectable. In case of a borderline result, a secure evaluation is not possible. If there is a clinical suspicion and a negative test result, we recommend clarification by means of other diagnostic methods and/or the serological investigation of a follow-up sample. A positive result indicates that there has been contact with the pathogen. In the determination of pathogen-specific IgM antibodies, polyclonal stimulation of the immune system or antibody persistence may affect the diagnostic relevance of positive findings. Significant titer increases (exceeding factor 2) and/or seroconversion in a follow-up sample taken after 7-10 days can indicate an acute infection. To investigate titer changes, sample and follow-up sample should be incubated in adjacent wells of the ELISA microplate within the same test run. For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against Borrelia, the calibration is performed in relative units (RU/ml).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio values determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.



The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extraction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is provided by particularly suitable Borrelia strains (Borrelia burgdorferi sensu stricto, Borrelia atzeili and Borrelia garinii). The cultured bacteria were solubilised using sodium dodecyl sulphate. The used antigen mixture contains all relevant proteins.

Linearity: The linearity of the Anti-Borrelia ELISA (IgM) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95. The Anti-Borrelia ELISA (IgM) is linear at least in the tested concentration range (3 RU/ml to 177 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-Borrelia ELISA (IgM) is 1.3 RU/ml.

Cross reactivity: Cross reactivity of the Anti-Borrelia ELISA (IgM) was evaluated in a study performed on 265 patient sera with antibodies against Treponema pallidum and 18 patient sera with antibodies against Leptospira. 26 sera of the Anti-Treponema pallidum-positive sera and 3 sera of the Anti-Leptospira-positive sera presented reactivity. Therefore, the cross reactivity amounts to 9.9% and 16.7%, respectively. Cross reactions with antibodies against other spirochaetes cannot be excluded. The Anti-Borrelia ELISA is designed as a screening test which should be followed by a confirmatory test (Western- or line blot). Further sera from patients with different possible influencing factors were investigated. An overview of these results can be found in the following table.

Possible influencing factors	n	Anti-Borrelia ELISA (IgM) positive
Acute EBV infection	27	0%
Anti-CMV	19	0%
Anti-Measles virus	13	0%
Anti-Mumps virus	12	0%
Anti-Toxoplasma gondii	14	0%
Anti-VZV	15	0%
Anti-Rubella virus	10	0%
Anti-HSV	5	0%
Anti-HAV	7	0%
Anti-HBV	8	0%
Anti-HAMA	4	0%
ANA + DNS (AAB)	45	0%
Rheuma factor	39	0%
Neurological diseases	54	1.9%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.



Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

Intra-assay variation, n = 20		
Serum	Mean value (RU/ml)	CV (%)
1	75	4.3
2	103	2.9
3	115	3.4

Inter-assay variation, n = 4 x 6		
Serum	Mean value (RU/ml)	CV (%)
1	78	6.9
2	109	6.2
3	122	4.7

Specificity and sensitivity: Sera from 150 patients with suspected borreliosis were analysed using the EUROIMMUN Anti-Borrelia ELISA (IgM) and the EUROIMMUN Anti-Borrelia burgdorferi Westernblot (IgM) as a reference method. The test showed a specificity of 96.4% and a sensitivity of 100%.

ELISA	EUROIMMUN Westernblot	
	positive	negative
n = 150	34	4
	0	106

Clinical study: 138 patients with clinically characterised borreliosis in different disease stages were screened with the Anti-Borrelia plus VisE-ELISA (IgG) and Anti-Borrelia-ELISA (IgM). Sensitivities between 89-96% were found.

Clinic	EUROIMMUN Anti-Borrelia plus VisE-ELISA (IgG), EUROIMMUN Anti-Borrelia-ELISA (IgM)		
	IgG positive	IgM positive	IgG und/or IgM positive
n = 138	60	68	86 (88.7%)
	22	5	25 (96.2%)
erythema migrans, mono- or poly-arthritis, neuroborreliosis	n = 97	n = 26	n = 15
	60	22	13

All ELISA negative results were confirmed as negative by an Anti-Borrelia Westernblot.

Reference range: The levels of the anti-Borrelia antibodies (IgM) were analysed with this EUROIMMUN ELISA in 500 healthy blood donors. With a cut-off of 20 RU/ml, 1.6% of the blood donors were anti-Borrelia positive (IgM).

Clinical significance

The history of Lyme disease, a contagious condition caused by Borrelia burgdorferi and transmitted to man by ticks, offers infectiousologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis. Fifty years later, Hauser noted the affection was transmitted by ticks. Independently, Atzeili, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemia of arthritis occurred in children in Lyme, Connecticut, USA. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemia, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that ticks transmitted "Treponema-like spirochaetes", which were later authenticated as the causal agent: Borrelia burgdorferi.



Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaceae family. This family includes five genera: borrelia, spirochaetes, cristispira, teponema and leptospira. Among about 30 tick species feeding on humans, Ixodes ricinus is the most frequent tick species biting humans in Europe. It is the vector of *Borrelia burgdorferi*, which causes Lyme disease, of *Anaplasma phagocytophilum* and the tick-borne encephalitis virus. I. ricinus ticks pass through three developmental stages: larvae, nymphs and adults (females and males). The density of this tick species can be very high, reaching in some places more than 300 ticks/100 m². Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere. In Europe, especially in Central Europe and Scandinavia, there are up to 155 cases per 100,000 individuals caused by the species *B. burgdorferi sensu stricto*, *B. afzelii*, and *B. garinii* (in rare cases also by *B. spielmanii* and *B. bavariensis*). The human seroprevalence rate of antibodies against *Borrelia burgdorferi* in the normal population of Germany and other Central European countries is about 8% (for IgG); in highly endemic areas even considerably higher with ranges from 18% to 52%. In East Asia, e.g. in an endemic area in China, the seroprevalence amounts to 26%. Persons working in the area of forestry display anti-*Borrelia* antibodies in about 40%, hunters in more than 50% of cases.

A *Borrelia burgdorferi* infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a *Borrelia burgdorferi* infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting. Lymphadenopathies are observed in a few cases (lymphadenosis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgia, influenza-like, cardiovascular borreliosis and borreliosis with hepatitis, with regional lymphadenitis and mixtures of these. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against *Borrelia burgdorferi* can be detected serologically in 50% to 90% of patients during stage I. Humans produce highly specific antibodies against the outer surface protein C (OspC) shortly after infection with *Borrelia burgdorferi*. The main antigen VlsE, which is exclusively expressed *in vivo*, is the most sensitive antigen for IgG antibody detection, whereas OspC is the most sensitive antigen for IgM antibody detection.

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears on one side and after several weeks on the opposite side. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found. Less frequently are cardiologic manifestations such as myocarditis and pericarditis. Antibodies against *Borrelia burgdorferi* can be detected in 50% to 90% of patients in stage II. In the early phase of this stage mainly IgM antibodies are found but in the late phase often only IgG antibodies occur. However, IgM antibodies can persist for a long time. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% to 30%. VlsE displays the highest sensitivity of all antigens tested for IgG detection.

Stage III: The typical manifestations of a *Borrelia burgdorferi* infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibody titres are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable.

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing. An antibody index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in the blood/CSF barrier function.



Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi*: ELISA, indirect immunofluorescence (IIFT), passive haemagglutination and immunoblot. Generally, ELISA or IIFT are used as screening test for preliminary characterisation of the serum sample. According to guidelines in the USA and in Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIFT) as a first step, which, if reactive, is followed by an immunoblot. In fresh infections it is recommended performing ELISA/IIFT and immuno-blot in parallel, since some weak reactions become detectable earlier in the blot than in the screening test. VlsE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection should be included in the test. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% compared to whole extract Westernblots. Of all recombinant antigens tested, VlsE possesses the highest sensitivity for the detection of a *Borrelia* infection.

Literature references

1. Burgdorfer W, Barbour AG, Hayes SF, Benach JL, Grunwaldt E, Davis JP. Lyme disease - a tick-borne Spirochetosis? Science 216 (1982) 1317-1319.
2. Christova I. Enzyme-linked immunosorbent assay, immunofluorescent assay, and recombinant immunoblotting in the serodiagnosis of early Lyme borreliosis. Int J Immunopathol Pharmacol 16 (2003) 261-268.
3. EUROIMMUN AG, Komorowski L, Janssen A, Stöcker W, Probst C. Dimerisation of recombinant OspC leads to an antigen with enhanced potential for active vaccination. 12. International Conference on Lyme Borreliosis and other Tick-Borne Diseases, Lubljana, Slovenia (2010).
4. EUROIMMUN AG, Komorowski L, Probst C, Janssen A, Stöcker W. Poly peptide und Verfahren zur spezifischen Detektion von Antikörpern bei Patienten mit einer Borrelien-Infektion. Europäische Patentanmeldung EP2199303 (2009).
5. EUROIMMUN AG, Meyer W, Janssen A, Schaper T, Schlumberger W, Stöcker W. Lyme borreliosis: Prevalence of antibodies against non-proteinic (lipid) antigens. Int J Med Microbiol 298S2 (Suppl. 45) (2008) 39.
6. EUROIMMUN AG, Probst C, Ott A, Schaper T, Meyer W, Stöcker W, Komorowski L. N-Terminal disulfide-bridging of *Borrelia* outer surface protein C increases its diagnostic and vaccine potentials. Ticks and Tick-Borne Diseases 3 (2012) 1-7.
7. EUROIMMUN AG, Meyer W, Schaper T, Schlumberger W, Stöcker W. Antibodies Against the Newly Identified Major Antigen VlsE: A milestone in the Serological Diagnosis of Lyme Borreliosis. Poster zum EUREGIO, Congress of Clinical Chemistry and Laboratory Medicine, Aachen (2003).
8. EUROIMMUN AG, Steinhagen K, Schlumberger W, Stöcker W. Nachweis einer spezifischen intrathekalen Antikörpersynthese mit modernen ELISA-Testsystemen. Hohe Trefferrate bei Multipler Sklerose und Neuroborreliose. J Lab Med 25 (2001) 135-149.
9. EUROIMMUN AG, Meyer W, Janssen A, Kühn B, Schaper T, Suer W, Komorowski L, Probst C, Schlumberger W, Stöcker W. EUROLINE *Borrelia* RN-AT: Marriage of authentic recombinant and native *borrelia* antigens in a line immunoassay format. International Journal of Medical Microbiology 298S2 (Suppl. 45): 39-40 (2008).
10. Goettner G, Schulte-Spechtel U, Hillenmann R, Liegl G, Wislke B, Fingerle V. Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of VlsE and DbpA homologues. J Clin Microbiol 43 (2005) 3602-3609.
11. Hauser U, Lehnert G, Lobentanzer R, Wislke B. Interpretation Criteria for Standardized WesternBlots for Three European Species of *Borrelia burgdorferi*. Sensu Latio J Clin Microbiol 35 (1997) 1433-1444.



12. Johnson R, Schmid GP, Hyde FW, Steigenwalt AG, Brenner DJ. *Borrelia burgdorferi* sp. Nov.: Etiologic agent of Lyme disease. *Int J Syst Bacteriol* 34 (1984) 496-497.
13. Oehme R, Hartert K, Backe H, Brockmann S, Kimmig P. Foci of tick-borne diseases in south-west Germany. *Int J Med Microbiol* 291 (2002) 22-29.
14. Robert-Koch-Institut. Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg. *Epidemiologisches Bulletin* 14 (1998) 93-97.
15. Schulte-Spechtel U, Lehner G, Liegl G, Fingerle V, Heimert C, Johnson BJ, Wilske B. Significant improvement of the recombinant Borrelia-specific immunoglobulin G immuno-blot test by addition of VISE and a DbpA homologue derived from *Borrelia garinii* for diagnosis of early neuroborreliosis. *J Clin Microbiol* 41 (2003) 1299-1303.
16. Tilton RC, Ryan RW. The laboratory diagnosis of Lyme disease. *J Clin Immunoassay* 16 (1993) 208-214.
17. Wilske B, Fingerle V, Schulte-Spechtel U. Microbiological and serological diagnosis of Lyme borreliosis. *FEMS Immunol Med Microbiol* 49 (2007) 13-21.
18. Wilske B. Epidemiology and diagnosis of Lyme borreliosis. *Ann Med* 37 (2005) 568-579.
19. Wilske B, Zöller L, Brade V, Effert H, Göbel UB, Stanek G, Pfister HW. MIQ 12/2000 Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik: Lyme-Borreliose. Urban & Fischer Verlag München Jena (2000).
20. Zheng JR, Hardham JM, Barbour AG, Norris SJ. Antigenic Variation in Lyme Disease *Borreliae* by Promiscuous Recombination of VMP-like Sequence Cassettes. *Cell* 89 (1997) 275-285.



Declaration of Conformity

EUROIMMUN Medizinische Labordiagnostika AG
Seekamp 31, D-23560 Lübeck, Germany

declare under our sole responsibility that the WESTERBLOT products

Anti-Borrelia burgdorferi-WESTERNBLOT (IgG)	DY 2132-#### G
Anti-Borrelia burgdorferi-WESTERNBLOT (IgM)	DY 2132-#### M

(product name, order no)

meet the demands of

Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998 and its transpositions in national laws which apply to it.

Conformity assessment procedure: Annex III

Lübeck, 04.10.2016
(Place and date of issue)

Dr. Wolfgang Schlumberger
- Member of the Board -

Susanne Aleksandrowicz
- Member of the Board -




Anti-Borrelia burgdorferi-WESTERNBLOT (IgG)
Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
DY 2132-3001 G	Borrelia burgdorferi	IgG	Antigen coated	30 x 01 (30)
DY 2132-24001 G	(complete antigens)	IgG	membrane strips	240 x 01 (240)

Indication: The Westernblot test kit provides a qualitative in vitro assay for human autoantibodies of the immunoglobulin class IgG against Borrelia in serum or plasma for the diagnosis of Lyme borreliosis and associated diseases (Erythema chronicum migrans, lymphadenitis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculoneuritis and neuroborreliosis).

Principles of the test: The test kit contains test strips with electrophoretically separated antigen extracts of Borrelia burgdorferi. The blot strips will be blocked and incubated in the first reaction step with diluted patient samples. In the case of positive samples, specific antibodies of the class IgG (and IgA, IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Format	Format	Symbol
1. Test strips Single strips with electrophoretically separated Borrelia burgdorferi antigens	30 x 1	240 x 1	[STRIPS]
2. Evaluation matrix with control strip Test strip incubated with a positive control serum	1 pattern	1 pattern per test strip lot	---
3. Enzyme conjugate Alkaline phosphatase-labelled anti-human IgG (goat), 10x concentrate	2 x 3 ml	16 x 3 ml	[CONJUGATE 10x]
4. Universal buffer 10x concentrate	1 x 100 ml	8 x 100 ml	[BUFFER 10x]
5. Substrate solution Nitroblue tetrazolium chloride/5-Bromo-4-chloro-3-indolylphosphate (NBT/BCIP), ready for use	1 x 50 ml	8 x 50 ml	[SUBSTRATE]
6. Adhesive foil	---	8 sheets	---
7. Test instruction	1 booklet	1 booklet	---
LOT			
IVD			
In vitro diagnostic medical device			

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Undiluted patient samples and incubated blot strips should be handled as infectious waste. Other reagents do not need to be collected separately, unless stated otherwise in official regulations.

The following components are not provided in the test kits but can be ordered at EUROIMMUN under the respective order numbers. Performance of the test requires an incubation tray:

ZD 9895-0130 Incubation tray with 30 channels (black)

ZD 9898-0144 Incubation tray with 44 channels (black, for the EUROBlotOne and EUROBlotCamera system)

For the creation of work protocols and the evaluation of incubated test strips using EUROlineScan green paper and adhesive foil are required:

ZD 9880-0101 Green paper (1 sheet)

ZD 9885-0116 Adhesive foil for approx. 16 test strips

ZD 9885-0130 Adhesive foil for approx. 30 test strips

if a visual evaluation is to be performed in individual cases, the required evaluation protocol can be ordered under: ZD 2132-0101 Evaluation protocol visual Anti-Borrelia burgdorferi-WESTERNBLOT.

Preparation and stability of the reagents

Note: The bag containing the blot strips is printed with a number in addition to the test kit lot number. This number refers to the strip batch and is also printed on the corresponding evaluation template. These two numbers must match to ensure correct evaluation of test results.

All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. Unopened, reagents are stable until the indicated expiry date when stored at +2°C to +8°C. After initial opening, reagents are stable for 12 months or until the expiry date, if earlier, unless stated otherwise in the instructions. Opened reagents must also be stored at +2°C to +8°C and protected from contamination.

- Coated test strips: Ready for use. Open the packing with the test strips only when the strips have reached room temperature to prevent condensation on the strips. After removal of the strips the packing should be sealed tightly and stored at +2°C to +8°C. To ensure correct evaluation of results, the lot number on the bag must match the lot number on the strips as well as on the evaluation matrix.

- Enzyme conjugate: The enzyme conjugate is supplied as a 10x concentrate. For the preparation of the ready for use enzyme conjugate the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with ready for use diluted universal buffer. For 1 test strip dilute 0.15 ml anti-human enzyme conjugate with 1.35 ml ready for use diluted universal buffer. The ready for use diluted enzyme conjugate should be used at the same working day.

- Universal buffer: The universal buffer is supplied as a 10x concentrate. For the preparation of the ready for use universal buffer the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with deionised or distilled water. For the incubation of 1 test strip 1.5 ml buffer concentrate should be diluted with 13.5 ml deionised or distilled water. The ready for use diluted universal buffer should be used at the same working day.

- Substrate solution: Ready for use. Close bottle immediately after use, as the contents are sensitive to light.

Warning: The control of human origin has tested negative for HBsAg, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.



Preparation and stability of the serum or plasma samples

Sample material: Human serum or EDTA, heparin or citrate plasma

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: The patient samples for analysis are diluted 1:51 in ready for use diluted universal buffer. For example, add 30 µl of sample to 1.5 ml ready for use diluted universal buffer and mix well by vortexing. Sample pipettes are not suitable for mixing.

Incubation

Blocking: According to the number of serum samples to be tested fill each channel of the incubation tray with 1.5 ml ready for use diluted universal buffer and a blot strip. Remove the required amount of blot strips from the packing using a pair of tweezers. The number on the test strip should be visible. Incubate for 15 minutes at room temperature (+18°C to +25°C) on a rocking shaker. Afterwards aspirate off all the liquid.

Sample incubation:
(1st step) Fill each channel with 1.5 ml of the diluted serum samples. Incubate at room temperature (+18°C to +25°C) for 30 minutes on a rocking shaker.

Wash: Aspirate off the liquid from each channel and wash 3 x 5 minutes each with 1.5 ml working strength universal buffer on a rocking shaker.

Conjugate incubation:
(2nd step) Pipette 1.5 ml diluted enzyme conjugate (alkaline phosphatase-conjugated anti-human IgG) into each channel. Incubate for 30 minutes at room temperature (+18°C to +25°C) on a rocking shaker.

Wash: Aspirate off the liquid from each channel. Wash as described above.

Substrate incubation:
(3rd step) Pipette 1.5 ml substrate solution into the channels of the incubation tray. Incubate for 10 minutes at room temperature (+18°C to +25°C) on a rocking shaker.


Stopping: Aspirate off the liquid from each channel and wash each strip 3 x 1 minute with deionised or distilled water.

For automated incubation with the EUROBlotMaster select the program Euro02 Inf WB30.

For automated incubation with the EUROBlotOne select the program Euro 01/02.

Anti-Borrelia burgdorferi-WESTERNBLOT (IgG)

Incubation protocol



Anti-Borrelia burgdorferi-WESTERNBLOT (IgG)
Evaluation and interpretation

Handling: For evaluation of incubated test strips we generally recommend using the EUROlineScan software. After stopping the reaction using deionised or distilled water, place the incubated test strips onto the adhesive foil of the green work protocol using a pair of tweezers. The position of the test strips can be corrected while they are wet. As soon as all test strips have been placed onto the protocol, they should be pressed hard using filter paper and left to air-dry. After they have dried, the test strips will be stuck to the adhesive foil. The dry test strips are then scanned using a flatbed scanner (EUROIMMUN AG) and evaluated with EUROlineScan. Alternatively, imaging and evaluation is possible directly from the incubation trays (EUROBioCamera and EUROBioOne). For general information about the EUROlineScan program please refer to the EUROlineScan user manual (EUROIMMUN AG). The code for entering the Test in EUROlineScan is B.b._WB_IgG.

If a visual evaluation must be performed in exceptional cases, hold the evaluation matrix next to the stuck-on blot strips and position it so that the black band above the number on the blot strips lines up with the alignment bar of the evaluation matrix. The lot number on the evaluation matrix must match the lot number on the blot strips. Clearly recognisable bands on the blot strips which concur with the labelled bands on the evaluation matrix are noted in the evaluation protocol.

Antigens: The antigen source for the EUROIMMUN Anti-Borrelia burgdorferi-WESTERNBLOT is provided by a particularly suitable Borrelia burgdorferi strain. The cultured Borrelia have been solubilised using sodium dodecyl sulphate followed by a separation of the solubilised protein using discontinuous polyacrylamide gel electrophoresis according to molecular mass and transfer of the separated proteins to nitrocellulose. From each test kit, 2 control test strips have been removed and incubated with a reference serum. One of these stained strips is included in the kit, the other remains with EUROIMMUN for documentation purposes. Diagnostically relevant antigens have been characterised with monoclonal reference antibodies from the German National Reference Laboratory for Borreliae.

Specificity of the antigens on the test strips:

Band	Antigen	Specificity
83 kDa	Membrane-vesical protein, p83	Degradation product of p100, high specificity
75 kDa	Heat shock protein, p75	Unspecific
62 kDa	Heat shock protein, p62	Unspecific
57/59 kDa	p57 and p59	Unspecific
50 kDa	p50	Unspecific
47 kDa	p47	Probably genus specific
43 kDa	p43	Unspecific
41 kDa	Flagellin, p41	Genus specific, cross reactivity to other spirochaetaeae and bacteria having flagella
39 kDa	Bmp A, p39	High specificity
34 kDa	Osp B, p34	Outer surface protein B, high specificity
32 kDa	p32	Unspecific
31 kDa	Osp A, p31	Outer surface protein A, high specificity
29 kDa	p29	Probably specific, poorly investigated
28 kDa	p28	Unspecific
25 kDa	Osp C, p25	Outer surface protein C, high specificity
21/22 kDa	p21/22	High specificity
18 kDa	p18	Probably specific
17 kDa	p17	Poorly investigated

In the lower part of the test strip there is a conjugate control membrane chip (IgA, IgG and IgM). Below the conjugate control, there is a membrane chip with a control band (Control).

Attention: A correctly performed determination of antibodies of class IgG against the antigens described above is indicated by a positive reaction of the control band and a positive reaction of the IgG band.

If one of these bands only shows a very weak reaction or none at all, the result is not valid.

Specificity of the antigens: Borrelia burgdorferi antigens can generally be divided into three categories.

Category	Antigens
1	Cross-reacting and undefined antigens with the molecular mass 17 kDa, 28 kDa, 32 kDa, 43 kDa, 47 kDa, 50 kDa, 57 kDa, 59 kDa, 62 kDa and 75 kDa.
2	A genus-specific antigen with a molecular mass of 41 kDa (flagellin).
3	Species-specific and highly specific antigens with the molecular mass 18 kDa, 21/22 kDa, 25 kDa, 31 kDa, 34 kDa, 39 kDa and 83 kDa.

IgG-class antibodies against Borrelia burgdorferi

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration, as negative sera sometimes produce weak signals in individual bands. Based on experience, the results of the Borrelia burgdorferi WESTERNBLOT test can be divided into negative, borderline and positive results.

Result	Characteristics
Negative	No bands, or weak intensities of some antigens from categories 1 and 2.
Borderline	A distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol) and several distinctive signals from categories 1 and 2. It is recommended that a new sample be taken and the test repeated after a few weeks.
Positive	More than one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol). In addition, and particularly in the case of patients in the late stage of the disease, numerous bands from categories 1 and 2 can be observed.

IgM-class antibodies against Borrelia burgdorferi

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration. The results of the Borrelia burgdorferi WESTERNBLOT test can be divided into negative, borderline and positive results.

In the early phase of a Borrelia infection IgM antibodies are typically directed against Osp C (p25). IgM antibodies against other specific Borrelia antigens are not considered definitive indicators of a fresh Borrelia infection.

IgM antibodies against flagellin (p41) can represent the initial response of the body to Borrelia burgdorferi. However, an unspecific reaction cannot be excluded, since it is known that antibodies directed against other microorganisms cross react with Borrelia burgdorferi flagellin (p41). For this reason, a single band at the position of flagellin (p41) in IgM detection should not be considered proof of a fresh infection with Borrelia burgdorferi. If only the flagellin band (p41) reacts positively, the test should be repeated several weeks later with a fresh blood sample.

In the serological investigation of a Borrelia infection the determination of antibodies of class IgM often yields unclear results. IgM antibodies can sometimes be found in serum years after an infection or following antibiotic treatment. Therefore, the detection of IgM antibodies does not necessarily indicate a fresh infection. A negative IgM result does not exclude a fresh infection. With a second infection, only antibodies of class IgG and not IgM can be formed.

In the late stage of borreliosis a positive IgM result does not provide any additional information, due to the antibody persistence mentioned above. The cause of such false-positive IgM results often remains unclear. They are observed, for example, in infectious mononucleosis, herpes virus infections and various autoimmune diseases.

For the diagnosis of a fresh Borrelia infection, a positive IgM result should be confirmed with a positive IgG result using a fresh blood sample 3 to 6 weeks later.

Result	Characteristics
Negative	No antigen bands recognisable or weak intensities of some bands of category 1
Borderline	One antigen band of category 2 (flagellin, p41) or a weak band of category 3. It is recommended that a fresh sample be taken and the test repeated after a few weeks. Antibodies against Osp C are characteristic for a fresh infection.
Positive	At least one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol). Antibodies against Osp C are characteristic for a fresh infection.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.



Test characteristics

Measurement range: The Westernblot is a qualitative method. No measurement range is provided.

Inter- and intra-assay variation: The inter-assay variation was determined by multiple analyses of characterised samples over several days. The intra-assay variation was determined by multiple analyses of characterised samples on one day. In every case, the intensity of the bands was within the specified range. This Westernblot displays excellent inter- and intra-assay reproducibility.

Interference: Haemolytic, lipaemic and icteric sera showed no effect on the analytical results.

Prevalence: Sera from 156 clinically characterised patients and 517 healthy blood donors were investigated with the EUROIMMUN Anti-Borrelia burgdorferi-Westernblot.

Clinically characterised sera	n	Prevalence		
		IgG	IgM	IgM/IgG
Erythema migrans	108	64%	67%	84%
Neuroborreliosis	32	85%	44%	88%
Arthritis	10	100%	30%	100%
Acrodermatitis	6	83%	50%	100%
Serum samples	n	Prevalence		
Healthy blood donors *	517	IgG 5%	IgM 4%	

*Medical University of Luebeck

The prevalence of anti-Borrelia antibodies from the healthy blood donor samples agrees with the values from literature [14].

Cross reactivity: The quality of the antigen used (whole antigen, SDS extract) and the antigen source (Borrelia burgdorferi sensu stricto) ensure high specificity of the Westernblot. The determination of cross reactivity is not necessary with the Westernblot, since specific reactions and unspecific or cross reactions can be differentiated directly with this test system.

Clinical significance

The history of Lyme disease, a contagious condition caused by Borrelia burgdorferi and transmitted to humans by ticks, offers infectiousologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis. Fifty years later, Hausser noted the affection was transmitted by ticks. Independently, Atzeius, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemia of arthritis occurred in children in Lyme, Connecticut, USA. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemia, Steele and co-workers defined "Lyme disease" as a neurological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that ticks transmitted "treponema-like spirochaetes", which were later authenticated as the causal agent: Borrelia burgdorferi.



Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaceae family. This family includes five genera: Borrelia, spirochetes, cristispira, treponema and leptospira. Among about 30 tick species feeding on humans, Ixodes ricinus is the most frequent tick species biting humans in Europe. It is the vector of Borrelia burgdorferi, which causes Lyme disease, of anaplasma phagocytophilum and the tick-borne encephalitis virus. I. ricinus ticks pass through three developmental stages: larvae, nymphs and adults (females and males). The density of this tick species can be very high, reaching in some places more than 300 ticks/100 m². Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100,000 individuals caused by the species B. burgdorferi sensu stricto, B. atzeii, and B. garinii. The human seroprevalence rate of antibodies against Borrelia burgdorferi in the normal population of Germany and other Central European countries is about 8% (for IgG), in highly endemic areas even considerably higher with ranges from 18% to 52%. In East Asia, e.g. in an endemic area in China, the seroprevalence amounts to 26%. Persons working in the area of forestry display anti-Borrelia antibodies in about 40%, hunters in more than 50% of cases.

A Borrelia burgdorferi infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a Borrelia burgdorferi infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting. Lymphadenopathies are observed in a few cases (lymphadenosus cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular borreliosis and borreliosis with hepatitis, with regional lymphadenitis and mixtures of these. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against Borrelia burgdorferi can be detected serologically in 50% to 90% of patients during stage I. Humans produce highly specific antibodies against the outer surface protein C (OspC) shortly after infection with Borrelia burgdorferi. The main antigen Vlse, which is exclusively expressed in vivo, is the most sensitive antigen for IgG antibody detection, whereas OspC is the most sensitive antigen for IgM antibody detection.

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears at one side and after several weeks on the opposite side. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found. Less frequently are cardiological manifestations such as myocarditis and pericarditis. Antibodies against Borrelia burgdorferi can be detected in 50% to 90% of patients in stage II. In the early phase of this stage mainly IgM antibodies are found but in the late phase often only IgG antibodies occur. However, IgM antibodies can persist for a long time. By additionally determining antibodies against Vlse the serological hit rate can be increased by 20% to 30%. Vlse displays the highest sensitivity of all antigens tested for IgG detection.

Stage III: The typical manifestations of a Borrelia burgdorferi infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibody titres are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable.

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing. An antibody index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in the blood/CSF barrier function.



Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi*. ELISA – indirect immunofluorescence (iIFIT) – passive haemagglutination and immunoblot. Generally, ELISA or iIFIT are used as screening test for preliminary characterisation of the serum sample. According to guidelines in the USA and in Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (iIFIT) as a first step, which, if reactive, is followed by an immunoblot. In fresh infections it's recommended performing ELISA/iIFIT and immunoblot in parallel, since some weak reactions become detectable earlier in the blot than in the screening test. ViSE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection should be included in the test. By additionally determining antibodies against ViSE the serological hit rate can be increased by 20% compared to whole extract Westernblots. Of all recombinant antigens tested, ViSE possesses the highest sensitivity for the detection of a *Borrelia* infection.

Literature references

- Burgdorfer W, Barbour AG, Hayes SF, Benach JL, Gunwaldt E, Davis JP. Lyme disease – a tick-borne Spirochetosis? *Science* 216 (1982) 1317-1319.
- Christova I. Enzyme-linked immunosorbent assay, immunofluorescent assay, and recombinant immunoblotting in the serodiagnosis of early Lyme borreliosis. *Int J Immuno-pathol Pharmacol* 16 (2003) 261-268.
- EUROIMMUN AG, Komorowski L, Janssen A, Stöcker W, Probst C. Dimersation of recombinant OspC leads to an antigen with enhanced potential for active vaccination. 12. International Conference on Lyme Borreliosis and other Tick-Borne Diseases, Ljubljana, Slovenia (2010).
- EUROIMMUN AG, Komorowski L, Probst C, Janssen A, Stöcker W. Poly peptide und Verfahren zur spezifischen Detektion von Antikörpern bei Patienten mit einer Borrelieninfektion. Europäische Patentanmeldung EP2199303 (2008).
- EUROIMMUN AG, Meyer W, Janssen A, Schepfer T, Schlumberger W, Stöcker W. Lyme borreliosis: Prevalence of antibodies against non-proteinic (lipid) antigens. *Int J Med Microbiol* 298S2 (Suppl. 45) (2008) 39.
- EUROIMMUN AG, Probst C, Ott A, Schepfer T, Meyer W, Stöcker W, Komorowski L. N-Terminal disulfide-bridging of *Borrelia* outer surface protein C increases its diagnostic and vaccine potentials. *Ticks and Tick-Borne Diseases* 3 (2012) 1-7.
- EUROIMMUN AG, Meyer W, Schepfer T, Schlumberger W, Stöcker W. Antibodies Against the Newly Identified Major Antigen ViSE: A milestone in the Serological Diagnosis of Lyme Borreliosis. Poster zum EUREGIO Congress of Clinical Chemistry and Laboratory Medicine, Aachen (2003).
- EUROIMMUN AG, Steinhagen K, Schlumberger W, Stöcker W. Nachweis einer spezifischen intrathekalen Antikörpersynthese mit modernen ELISA-Testsystemen: Hohe Trefferrate bei Multiplex Sklerose und Neuroborreliose. *J Lab Med* 25 (2001) 135-149.
- EUROIMMUN AG, Meyer W, Janssen A, Kühn B, Schepfer T, Suer W, Komorowski L, Probst C, Schlumberger W, Stöcker W. EUROLINE *Borrelia* RN-AT: Marriage of authentic recombinant and native *borrelia* antigens in a line immunoassay format. *International Journal of Medical Microbiology* 298S2 (Suppl. 45): 39-40 (2008).
- Goethner G, Schulte-Spechtel U, Hillermann R, Liegl G, Wilske B, Fingerle V. Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of ViSE and DbpA homologues. *J Clin Microbiol* 43 (2005) 3602-3609.
- Hauser U, Lehnert G, Lobentanzer R, Wilske B. Interpretation Criteria for Standardized WesternBlots for Three European Species of *Borrelia burgdorferi* Sensu Lato. *J Clin Microbiol* 35 (1997) 1433-1444.
- Johnson R, Schmid GP, Hyde FW, Steigewalt AG, Brenner DJ. *Borrelia burgdorferi* sp. Nov. Etiologic agent of Lyme disease. *Int J Syst Bacteriol* 34 (1984) 496-497.
- Oehme R, Hartelt K, Backe H, Brockmann S, Kirming P. Foci of tick-borne diseases in southwest Germany. *Int J Med Microbiol* 291 (2002) 22-29.
- Robert-Koch-Institut. Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg. *Epidemiologisches Bulletin* 14 (1998) 93-97.
- Schulte-Spechtel U, Lehnert G, Liegl G, Fingerle V, Heimerl C, Johnson BJ, Wilske B. Significant improvement of the recombinant *Borrelia*-specific immunoglobulin G immunoblot test by addition of ViSE and a DbpA homologue derived from *Borrelia garinii* for diagnosis of early neuroborreliosis. *J Clin Microbiol* 41 (2003) 1299-1303.
- Tilton RC, Ryan RW. The laboratory diagnosis of Lyme disease. *J Clin Immunassay* 16 (1993) 208-214.
- Wilske B, Fingerle V, Schulte-Spechtel U. Microbiological and serological diagnosis of Lyme borreliosis. *FEMS Immunol Med Microbiol* 49 (2007) 13-21.
- Wilske B. Epidemiology and diagnosis of Lyme borreliosis. *Ann Med* 37 (2005) 588-579.
- Wilske B, Zöller L, Brade V, Eifert H, Göbel UB, Stanek G, Pfister HW, MQ 12/2000 Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik: Lyme-Borreliose. Urban & Fischer Verlag München Jena (2000).
- Zhang JR, Hardham JM, Barbour AG, Norris SJ. Antigenic Variation in Lyme Disease *Borrelia* by Promiscuous Recombination of VMP-like Sequence Cassettes. *Cell* 89 (1997) 275-285.





Anti-Borrelia burgdorferi-WESTERNBLOT (IgM) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
DY 2132-3001 M	Borrelia burgdorferi (complete antigens)	IgM	Antigen coated membrane strips	30 x 01 (30)
DY 2132-24001 M				240 x 01 (240)

Indication: The Westernblot test kit provides a qualitative in vitro assay for human autoantibodies of the immunoglobulin class IgM against Borrelia in serum or plasma for the diagnosis of Lyme borreliosis. Associated diseases: Erythema chronicum migrans, lymphadenitis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculoneuritis and neuroborreliosis.

Principles of the test: The test kit contains test strips with electrophoretically separated antigen extracts of Borrelia burgdorferi. The blot strips will be blocked and incubated in the first reaction step with diluted patient samples. In the case of positive samples, specific antibodies of the class IgM (and IgA, IgG) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgM (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Format	Format	Symbol
1. Test strips Single strips with electrophoretically separated Borrelia burgdorferi antigens	30 x 1	240 x 1	STRIPS
2. Evaluation matrix with control strip Test strip incubated with a positive control serum	1 pattern	1 pattern per test strip lot	
3. Enzyme conjugate Alkaline phosphatase-labelled anti-human IgM (goat), 10x concentrate	2 x 3 ml	16 x 3 ml	CONJUGATE 10X
4. Universal buffer 10x concentrate	1 x 100 ml	8 x 100 ml	BUFFER 10X
5. Substrate solution Nitroblue tetrazolium chloride/5-Bromo-4-chloro-3-indolylphosphate (NB/TB/CLIP), ready for use	1 x 50 ml	8 x 50 ml	SUBSTRATE
6. Adhesive foil	---	8 sheets	
7. Test instruction	1 booklet	1 booklet	
LOT			
INVD			

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Undiluted patient samples and incubated blot strips should be handled as infectious waste. Other reagents do not need to be collected separately, unless stated otherwise in official regulations.

The following components are not provided in the test kits but can be ordered at EUROIMMUN under the respective order numbers.

Performance of the test requires an incubation tray:
ZD 9895-0130 Incubation tray with 30 channels (black)
ZD 9898-0144 Incubation tray with 44 channels (black) for the EURObioOne and EURObioCamera system)

For the creation of work protocols and the evaluation of incubated test strips using EUROlineScan green paper and adhesive foil are required:

ZD 9880-0101 Green paper (1 sheet)

ZD 9885-0116 Adhesive foil for approx. 16 test strips

ZD 9885-0130 Adhesive foil for approx. 30 test strips

If a visual evaluation is to be performed in individual cases, the required evaluation protocol can be ordered under: ZD 2132-0101 Evaluation protocol visual Anti-Borrelia-burgdorferi-WESTERNBLOT

Preparation and stability of the reagents

Note: The bag containing the blot strips is printed with a number in addition to the test kit lot number. This number refers to the strip batch and is also printed on the corresponding evaluation template. These two numbers must match to ensure correct evaluation of test results.

All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. Unopened, reagents are stable until the indicated expiry date when stored at +2°C to +8°C. After initial opening, reagents are stable for 12 months or until the expiry date, if earlier, unless stated otherwise in the instructions. Opened reagents must also be stored at +2°C to +8°C and protected from contamination.

- **Coated test strips:** Ready for use. Open the packing with the test strips only when the strips have reached room temperature to prevent condensation on the strips. After removal of the strips the packing should be sealed tightly and stored at +2°C to +8°C. To ensure correct evaluation of results, the lot number on the bag must match the lot number on the strips as well as on the evaluation matrix.

- **Enzyme conjugate:** The enzyme conjugate is supplied as a 10x concentrate. For the preparation of the ready for use enzyme conjugate the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with ready for use diluted universal buffer. For 1 test strip dilute 0.15 ml anti-human IgG concentrate with 1.35 ml ready for use diluted universal buffer. The ready for use diluted enzyme conjugate should be used at the same working day.

- **Universal buffer:** The universal buffer is supplied as a 10x concentrate. For the preparation of the ready for use universal buffer the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with deionised or distilled water. For the incubation of 1 test strip 1.5 ml buffer concentrate should be diluted with 13.5 ml deionised or distilled water. The ready for use diluted universal buffer should be used at the same working day.

- **Substrate solution:** Ready for use. Close bottle immediately after use, as the contents are sensitive to light.

Warning: The control of human origin has tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.



Preparation and stability of the serum or plasma samples

Anti-Borrelia-burgdorferi WESTERNBLOT (IgM)

Sample material: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: The patient samples for analysis are diluted 1:51 in ready for use diluted universal buffer. For example, add 30 µl of sample to 1.5 ml ready for use diluted universal buffer and mix well by vortexing. Sample pipettes are not suitable for mixing.

Incubation protocol

Incubation

Blocking: According to the number of serum samples to be tested fill each channel of the incubation tray with 1.5 ml ready for use diluted universal buffer and a blot strip. Remove the required amount of blot strips from the packing using a pair of tweezers. The number on the test strip should be visible. Incubate for 15 minutes at room temperature (+18°C to +25°C) on a rocking shaker. Afterwards aspirate off all the liquid.

Sample incubation: (1st step) Fill each channel with 1.5 ml of the diluted serum samples and incubate at room temperature (+18°C to +25°C) for 30 minutes on a rocking shaker.

Wash: Aspirate off the liquid from each channel and wash 3 x 5 minutes each with 1.5 ml working strength universal buffer on a rocking shaker.

Conjugate incubation: (2nd step) Pipette 1.5 ml diluted enzyme conjugate (alkaline phosphatase-conjugated anti-human IgM) into each channel and incubate for 30 minutes at room temperature (+18°C to +25°C) on a rocking shaker.


Wash: Aspirate off the liquid from each channel. Wash as described above.

Substrate incubation: (3rd step) Pipette 1.5 ml substrate solution into the channels of the incubation tray. Incubate for 10 minutes at room temperature (+18°C to +25°C) on a rocking shaker.

Stopping: Aspirate off the liquid from each channel and wash each strip 3 x 1 minute with deionised or distilled water.

For automated incubation with the EUROBlotMaster select the program Euro02 Int WB30.

For automated incubation with the EUROBlotOne select the program Euro 01/02.





Evaluation and Interpretation of the results of the Anti-Borrelia-burgdorferi WESTERNBLOT (IgM)

Handling: For evaluation of incubated test strips we generally recommend using the EUROlineScan software. After stopping the reaction using deionised or distilled water, place the incubated test strips onto the adhesive foil of the green work protocol using a pair of tweezers. The position of the test strips can be corrected while they are wet. As soon as all test strips have been placed onto the protocol, they should be pressed hard using filter paper and left to air-dry. After they have dried, the test strips will be stuck to the adhesive foil. The dry test strips are then scanned using a flatbed scanner (EUROIMMUN AG) and evaluated with EUROlineScan. Alternatively, imaging and evaluation is possible directly from the incubation trays (EUROBioCamera and EUROBioOne). For general information about the EUROlineScan program please refer to the EUROlineScan user manual (EUROIMMUN AG). The code for entering the Test in EUROlineScan is B_b_WB_IgM.

If a visual evaluation must be performed in exceptional cases, hold the evaluation matrix next to the stuck-on blot strips and position it so that the black band above the number on the blot strips lines up with the alignment bar of the evaluation matrix. The lot number on the evaluation matrix must match the lot number on the blot strips. Clearly recognisable bands on the blot strips which concur with the labeled bands on the evaluation matrix are noted in the evaluation protocol.

Antigens: The antigen source for the EUROIMMUN Anti-Borrelia-burgdorferi WESTERNBLOT is provided by a particularly suitable *Borrelia burgdorferi* strain [2, 3]. The cultured borrelia have been solubilised using sodium dodecyl sulphate followed by a separation of the solubilised protein using discontinuous polyacrylamide gel electrophoresis according to molecular mass and transfer of the separated proteins to nitrocellulose. From each test kit, 2 control test strips have been removed and incubated with a reference serum. One of these stained strips is included in the kit, the other remains with EUROIMMUN for documentation purposes. Diagnostically relevant antigens have been characterised with monoclonal reference antibodies from the German National Reference Laboratory for Borreliae [39].

Band	Antigen	Specificity
83 kDa	Membrane-vesical protein, p 83	Degradation product of p 100, high specificity.
75 kDa	Heat shock protein, p 75	Unspecific.
62 kDa	Heat shock protein, p 62	Unspecific.
57/59 kDa	p 57 and p 59	Unspecific.
50 kDa	p 50	Unspecific.
47 kDa	p 47	Probably genus specific.
43 kDa	p 43	Unspecific.
41 kDa	Flagellin, p 41	Genus specific, cross reactivity to other spirochaetae and bacteria having flagella
39 kDa	Bmp A, p 39	High specificity.
34 kDa	Osp B, p 34	Outer surface protein B, high specificity.
32 kDa	p 32	Unspecific.
31 kDa	Osp A, p 31	Outer surface protein A, high specificity.
29 kDa	p 29	Probably specific, poorly investigated.
28 kDa	p 28	Unspecific.
25 kDa	Osp C, p 25	Outer surface protein C, high specificity.
21/22 kDa	p 21/22	High specificity.
18 kDa	p 18	Probably specific.
17 kDa	p 17	Poorly investigated.

In the lower part of the test strip there is a conjugate control membrane chip (IgA, IgG and IgM). Below the conjugate control, there is a membrane chip with a control band (Control).

Attention: A correctly performed determination of antibodies of class IgM against the antigens described above is indicated by a positive reaction of the control band and a positive reaction of the IgM band.

If one of these bands only shows a very weak reaction or none at all, the result is not valid.

Specificity of the antigens: *Borrelia burgdorferi* antigens can generally be divided into three categories [40].

Category	Antigens
1	Cross-reacting and undefined antigens with the molecular mass 17 kDa, 28 kDa, 32 kDa, 43 kDa, 47 kDa, 50 kDa, 57 kDa, 59 kDa, 62 kDa and 75 kDa.
2	A genus-specific antigen with a molecular mass of 41 kDa (flagellin).
3	Species-specific and highly specific antigens with the molecular mass 18 kDa, 21/22 kDa, 25 kDa, 31 kDa, 34 kDa, 39 kDa and 83 kDa.

IgG class antibodies against *Borrelia burgdorferi*:

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration, as negative sera sometimes produce weak signals in individual bands. Based on experience, the results of the *Borrelia burgdorferi* WESTERNBLOT test can be divided into negative, borderline and positive results.

Result	Characteristics
Negative	No bands, or weak intensities of some antigens from categories 1 and 2.
Borderline	A distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol) and several distinctive signals from categories 1 and 2. It is recommended that a new sample be taken and the test repeated after a few weeks.
Positive	More than one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol). In addition, and particularly in the case of patients in the late stage of the disease, numerous bands from categories 1 and 2 can be observed.

IgM class antibodies against *Borrelia burgdorferi*:

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration. The results of the *Borrelia burgdorferi* WESTERNBLOT test can be divided into negative, borderline and positive results.

In the early phase of a *Borrelia burgdorferi* infection IgM antibodies are typically directed against Osp C (p 25), IgM antibodies against other specific *Borrelia burgdorferi* antigens are not considered definitive indicators of a fresh *Borrelia burgdorferi* infection.

IgM antibodies against flagellin (p 41) can represent the initial response of the body to *Borrelia burgdorferi*. However, an unspecific reaction cannot be excluded, since it is known that antibodies directed against other microorganisms cross react with *Borrelia burgdorferi* flagellin (p 41). For this reason, a single band at the position of flagellin (p 41) in IgM detection should not be considered proof of a fresh infection with *Borrelia burgdorferi*. If only the flagellin band (p 41) reacts positively, the test should be repeated several weeks later with a fresh blood sample.

In the serological investigation of a *Borrelia burgdorferi* infection the determination of antibodies of class IgM often yields unclear results. IgM antibodies can sometimes be found in serum years after an infection or following antibiotic treatment. Therefore, the detection of IgM antibodies does not necessarily indicate a fresh infection. A negative IgM result does not exclude a fresh infection. With a second infection, only antibodies of class IgG and not IgM can be formed.

In the late stage of borreliosis a positive IgM result does not provide any additional information, due to the antibody persistence mentioned above. The cause of such false-positive IgM results often remains unclear. They are observed, for example, in infectious mononucleosis, herpes virus infections and various autoimmune diseases [39].

For the diagnosis of a fresh *Borrelia burgdorferi* infection, a positive IgM result should be confirmed with a positive IgG result using a fresh blood sample 3-6 weeks later.

Result	Characteristics
Negative	No antigen bands recognisable or weak intensities of some bands of category 1
Borderline	One antigen band of category 2 (flagellin, p 41) or a weak band of category 3. It is recommended that a fresh sample be taken and the test repeated after a few weeks. Antibodies against Osp C are characteristic for a fresh infection.
Positive	At least one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol). Antibodies against Osp C are characteristic for a fresh infection.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.



Test characteristics

Measurement range: The Westernblot is a qualitative method. No measurement range is provided

Inter- and intra-assay variation: The inter-assay variation was determined by multiple analyses of characterised samples over several days. The intra-assay variation was determined by multiple analyses of characterised samples on one day. In every case, the intensity of the bands was within the specified range. This Westernblot displays excellent inter- and intra-assay reproducibility.

Interference: Haemolytic, lipaemic and icteric sera showed no effect on the analytical results.

Prevalence: Sera from 156 clinically characterised patients and 517 healthy blood donors were investigated with the EUROIMMUN Anti-Borrelia burgdorferi-WESTERNBLOT.

Clinically characterised sera	n	IgG	IgM	IgM/IgG
Erythema migrans	108	64%	67%	84%
Neuroborreliosis	32	85%	44%	88%
Arthritis	10	100%	30%	100%
Acrodermatitis	6	83%	50%	100%

Serum samples	n	IgG	IgM	Prevalence
Healthy blood donors *	517	5%	4%	4%

*Medical University of Luebeck

The prevalence of anti-Borrelia antibodies from the healthy blood donor samples agrees with the values from literature [7].

Cross reactivity: The quality of the antigen used (whole antigen, SDS extract) and the antigen source (Borrelia burgdorferi sensu stricto) ensure high specificity of the Westernblot. The determination of cross reactivity is not necessary with the Westernblot, since specific reactions and unspecific or cross reactions can be differentiated directly with this test system.

Clinical significance

The history of Lyme disease, a contagious condition caused by Borrelia burgdorferi, transmitted to humans by ticks, offers infectiologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis [1]. Fifty years later Hauser noted the affection was transmitted by ticks [1]. Independently, Atzeilus, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemic of arthritis occurred in children in Lyme, Connecticut, USA [1]. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemic, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that tick bites transmitted a Spirochaeta which was later authenticated as the causal agent: Borrelia burgdorferi [1, 2, 3].

Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaeaceae family. This family includes five genera: borrelia, spirochaetes, cristispira, treponema and leptospira. Among about 30 tick species feeding on humans, Ixodes ricinus is the most frequent tick species biting humans in Europe. It is the vector of Borrelia burgdorferi which causes Lyme disease, and of the tick-borne encephalitis virus and males) [8]. The density of this tick species may be very high, reaching in some places more than 300 ticks/100 m² [8]. Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere [9, 10, 11]. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100 000 individuals caused by the species B. burgdorferi sensu stricto, B. atzeili, and B. garinii [4, 6, 10, 11, 12]. The human seroprevalence rates of antibodies against Borrelia burgdorferi in the normal population of Germany and other Central European countries ranges from 18% to 52% [13]. In East Asia, e.g. in China, there are 26% on the average [14]. Persons working in the area of forestry display anti-borrelia antibodies in about 40%, hunters in more than 50% of cases [7, 15].

A Borrelia burgdorferi infection can manifest itself in the areas of dermatology, neurology, ophthalmology, pneumatology and internal medicine [7, 16]. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a Borrelia burgdorferi infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner [17]. The erythema is accompanied by influenza-like general symptoms with fever, shaking chills, headaches and vomiting [18]. Lymphadenopathies are observed in a few cases (lymphadenosis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular, hepatitis, regional lymphadenitis, and mixed [19, 20, 21]. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against Borrelia burgdorferi can be detected serologically in 50% to 90% of patients in Stage I. The prevalence of specific IgG antibodies is considerably lower [22, 23]. However, serological tests often provide negative results in this stage of the illness [3, 19]. In order to overcome this problem a new test was created including the newly identified VISE (variable major protein-like sequence, expressed), which can be considered as the major antigen for Borrelia serology. Over 85% of IgG-positive sera can be identified at a glance by assessing the VISE band in incubated Westernblot strips. VISE allows detection of antibodies against all Borrelia species. The risk of a false negative reaction due to species differences is ten times lower [24]. Humans also produce highly specific borrelia-specific antibodies against outer surface protein C (OspC), shortly after infection with Borrelia burgdorferi sensu stricto: immunoglobulin M (IgM) OspC and immunoglobulin G (IgG) OspC borrelial antibodies [25]. VISE is the most sensitive antigen for IgG antibody detection, OspC for IgM antibody detection [3, 26, 27].

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve parestheses, Bannwarth's lymphocytic meningoradiculoneuritis [20, 28, 29]. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears at one side and after several weeks on the opposite side [20]. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found [30]. Rarer are cardiological manifestations, such as myocarditis and pericarditis. Antibodies against Borrelia burgdorferi can be detected in 50% to 90% of patients in stage II [24]. In the early phase of this stage, mainly IgM antibodies are found, but in the late phase, often only IgG antibodies occur. However, IgM antibodies can persist for a long time [3, 30]. By additionally determining antibodies against VISE the serological hit rate can be increased by 20%-30% [31]. VISE displays the highest sensitivity of all antigens tested for IgG detection [26].

Stage III: The typical manifestations of a Borrelia burgdorferi infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibodies are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable [3, 24].



For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing [21, 29, 32]. An Antibody Index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in blood/CSF barrier function [33].

Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi* [34, 35, 36, 37, 38, 39]: ELISA, indirect immunofluorescence, passive haemagglutination and immunoblot [26]. Many investigators employ an ELISA test in parallel with an indirect immunofluorescence test for preliminary characterisation of the serum sample [26, 34]. According to guidelines of the USA and Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIF) as first step. If reactive, followed by Westernblot [3, 24, 40]. In fresh infections it is recommended performing ELISA/IIF and Westernblot in parallel, since some weak reactions become detectable earlier in Westernblot than in the screening tests [12, 24]. VISE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection, should be included in the test [3]. By additionally determining antibodies against VISE the serological hit rate can be increased by 20% compared to whole extract Westernblots and by 30% compared to recombinant antigen Westernblots [24, 27, 31]. Of all recombinant antigens tested, VISE possesses the highest sensitivity for the detection of a *Borrelia* infection [31, 32, 33].

Literature references

1. Belach S. Lyme disease. *Presse Med* 14 (1995) 81-87.
2. Burgdorfer W, Barbour AG, Hayes SF, Benach JL, Gunwaldt E, Davis JP. Lyme disease - a tick-borne Spirochetosis? *Science* 216 (1982) 1317-1319.
3. Johnson R, Schmid GP, Hyde FW, Steigerwalt AG, Brenner DJ. *Borrelia burgdorferi* sp. nov.: Etiologic agent of Lyme disease. *Int J Syst Bacteriol* 34 (1984) 496-497.
4. Wilske B. Epidemiology and diagnosis of Lyme borreliosis. *Ann Med* 37 (2005) 568-579.
5. Derdakova M, Lencakova D. Association of genetic variability within the *Borrelia burgdorferi* sensu lato with the ecology, epidemiology of Lyme borreliosis in Europe. *Ann Agric Environ Med* 12 (2005) 165-172.
6. Steffen I, Hirsch HH. Diagnostic tests of Lyme borreliosis. *Ther Umsch* 62 (2005) 737-744.
7. Robert-Koch-Institut. Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg. *Epidemiologisches Bulletin* 14 (1998) 93-97.
8. Gem L. The biology of the Ixodes ricinus tick. *Ther Umsch* 62 (2005) 707-712.
9. Holzer BR. Tick borne diseases. *Ther Umsch* 62 (2005) 757-763.
10. Bazovska S, Machacova E, Spalekova M, Kontrosova S. Reported incidence of Lyme disease in Slovakia and antibodies to *B. burgdorferi* antigens detected in healthy population. *Bratisl Lek Listy* 106 (2005) 270-273.
11. Egger M. Lyme borreliosis - an overview. *Ther Umsch* 62 (2005) 731-736.
12. Hauser U, Lehner G, Lobentanz R, Wilske B. Interpretation Criteria for Standardized Western Blots for Three European Species of *Borrelia burgdorferi* Sensu Lato. *J Clin Microbiol* 35 (1997) 1433-1444.
13. Oehme R, Harrelt K, Backe H, Brockmann S, Kimmig P. Foci of tick-borne diseases in southwest Germany. *Int J Med Microbiol* 291 (2002) 22-29.
14. Ai CX, Zhang WF, Zhao JH. Sero-epidemiology of Lyme disease in an endemic area in China. *Microbiol Immunol* 38 (1994) 505-509.
15. Cetin E, Soitoudoh M, Auer H, Stanek G. Paradigm Burgendland: risk of *Borrelia burgdorferi* sensu lato infection indicated by variable seroprevalence rates in hunters. *Wien Klin Wochenschr* 118 (2006) 677-681.
16. Escudero-Nieto R, Guerrero-Espejo A. Diseases produced by *Borrelia*. *Enferm Infecc Microbiol Clin* 23 (2005) 232-240.
17. Vorobeva NN, Sumilvaia ON. Clinical manifestations in acute Lyme diseases. *Med Parazitol (Mosk)* 4 (2003) 3-7.
18. Roscoe C, Epperly T. Tick-borne relapsing fever. *Am Fam Physician* 15 (2005) 2039-2044.
19. Coulter P, Lema C, Flayhart D, Umhardt AS, Aucott JN, Auwaerter PG, Dumlner JS. Two-year evaluation of *Borrelia burgdorferi* culture and supplemental tests for definitive diagnosis of Lyme disease. *J Clin Microbiol* 43 (2005) 5080-5084.
20. Miodzikowska-Albrecht J, Zarowski M, Stenborn B, Winczewska-Wiktor A, Gurdia B, Wigowska-Sowinska J. Bilateral facial nerve palsy in the course of neuroborreliosis in children-dynamics, laboratory tests and treatment. *Rocz Akad Med Białymst* 50 Suppl 1 (2005) 64-69.
21. Christen HJ, Hanefeld F, Effert H, Thomssen R. Epidemiology and clinical manifestations of neuroborreliosis. *Acta Paediatr Suppl* 386 (1993) 1-75.
22. Bykasin AV, Godfrey HP, Carbonaro CA, Wormser GP, Agüero-Rosenfeld ME, Cabello FC. *Borrelia burgdorferi* BmpA, BmpB, and BmpD proteins are expressed in human infection and contribute to P39 immunoblot reactivity in patients with Lyme disease. *Clin Diagn Lab Immunol* 12 (2005) 935-940.
23. Christova I. Enzyme-linked immunosorbent assay, immunofluorescent assay, and recombinant immunoblotting in the serodiagnosis of early Lyme borreliosis. *Int J Immunopathol Pharmacol* 16 (2003) 261-268.
24. Wilske B, Fingerle V, Schulte-Spechtel U. Microbiological and serological diagnosis of Lyme borreliosis. *FEMS Immunol Med Microbiol* 49 (2007) 13-21.
25. Lovich SD, Jobe DA, Schnell RF, Callister SM. Borrelial OSpC antibodies specific for a highly conserved epitope are immunodominant in human Lyme disease and do not occur in mice or hamsters. *Clin Diagn Lab Immunol* 12 (2005) 746-751.
26. Goettner G, Schulte-Spechtel U, Hillermann R, Liegl G, Wilske B, Fingerle V. Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of VISE and DbpA homologues. *J Clin Microbiol* 43 (2005) 3602-3609.
27. Zhang JR, Hardham JM, Barbour AG, Norris SJ. Antigenic Variation in Lyme Disease Borreliae by Promiscuous Recombination of VMP-like Sequence Cassettes. *Cell* 89 (1997) 275-285.
28. Brooks CS, Vuppala SR, Jett AM, Akins DR. Identification of *Borrelia burgdorferi* outer surface proteins. *Infect Immun* 74 (2006) 296-304.
29. Steinhoff AP, Smith MJ, Shah SS, Coffin SE. Neuroborreliosis with progression from pseudotumor cerebri to aseptic meningitis. *Pediatr Infect Dis J* 25 (2006) 91-92.
30. Shi C, Wolfe J, Russell JQ, Forner K, Hardin N, Anguita J, Budd RC. Fas ligand deficiency impairs host inflammatory response. *Infect Immun* 74 (2006) 1156-1160.
31. Schulte-Spechtel U, Lehner G, Liegl G, Fingerle V, Heimerl C, Johnson BJ, Wilske B. Significant improvement of the recombinant *Borrelia*-specific immunoglobulin G immunoblot test by addition of VISE and a DbpA homologue derived from 2003. 1299-1303.
32. Schulte-Spechtel U, Lehner G, Liegl G, Fingerle V, Heimerl C, Johnson B, Wilske B. Significant improvement of the recombinant *Borrelia*-IgG immunoblot for serodiagnosis of early neuroborreliosis. *Int J Med Microbiol* 293 (2004) 158-160.
33. Reiber H, Lange P. Virus-spezifische Antikörper in Liquor und Serum. ELISA-Analytik und Auswertung mittels Antikörper-Index- und Quotientendiagramm. *LabMed* 15 (1997) 204-207.
34. Jovicic VJ, Grego EM, Lako BL, Risovic BM, Lepsanovic ZA, Stajkovic NT. Improved serodiagnosis of early Lyme borreliosis: immunoblot with local *Borrelia* atzei strain. *APMIS* 111 (2003) 1053-1059.



35. Komorowski L, Stöcker W (Erfinder). Serum-freies Medium für die Kultivierung von Borrelien und dessen Benutzung. EUROIMMUN Medizinische Labordiagnostika AG. Patent. Aktenzeichen: DE 102006045252.6-41 (2006).
36. Stöcker W, Fauer H, Krause C, Barth E, Martinez A. Verfahren zur Optimierung der automatischen Fluoreszenzkernung in der Immundiagnostik. Deutsche Patentanmeldung DE 10 2006 027 516.0 (2006).
37. Meyer W, Scheper T, Lehmann H. EUROIMMUN Medizinische Labordiagnostika AG. Selbstklebende Blotmembranen. Eingetragenes deutsches Gebrauchsmuster DE 202 15 269.5 (2003).
38. Morini M, Rateike M, Müller M. EUROIMMUN Medizinische Labordiagnostika AG. Lichtquelle für ein Aurichthfluoreszenzzykroskop. Eingetragenes deutsches Gebrauchsmuster DE 20 2004 010 121U (2004).
39. Wilske B, Zöller L, Brade V, Eifert H, Göbel UB, Stanek G, Pfister HW. MIQ 12/2000 Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik. Lyme-Borreliose. Urban & Fischer Verlag München Jena (2000).
40. Tilton RC, Ryan RW. The laboratory diagnosis of Lyme disease. J Clin Immunocassay 16 (1993) 208-214.
41. Ma B, Christen B, Leung D, Vigo-Pelfrey C. Serodiagnosis of Lyme Borreliosis by Western Immunoblot: Reactivity of varus significant antibodies against Borrelia burg-dorferi. J Clin Microbiol; Feb. 370-376 (1992).
42. Zöller L, Burkard S, Schäfer H. Validity of Western Immunoblot band patterns in the serodiagnosis of Lyme Borreliosis. J Clin Microbiol; Jan; 174-182 (1991).
43. Moskophidis M, Lütner B. Wertigkeit des Immunoblots in der Serodiagnostik der Lyme-Borreliose. LAB-med; 19: 231-237 (1995).

Anti-West Nile Virus ELISA (IgG) Test instruction

ORDER NO	ANTIBODIES AGAINST West Nile virus	IG CLASS	SUBSTRATE	FORMAT
EI 2662-9601 G		IgG	Ag-coated microplate wells	96 x 01 (96)

Principle of the test: The ELISA test kit provides a semiquantitative or quantitative *in vitro* assay for human antibodies of the IgG class against West Nile virus in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with West Nile virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG (also IgA and IgM) antibodies will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens, 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	[STRIPS]
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	[CAL1]
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	[CAL2]
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	[CAL3]
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	[POS CONTROL]
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	[NEG CONTROL]
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	[CONJUGATE]
8. Sample buffer ready for use	light blue	1 x 100 ml	[SAMPLE BUFFER]
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	[WASH BUFFER 10X]
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	[SUBSTRATE]
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	[STOP SOLUTION]
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
14. Protective foil	---	2 pieces	
LOT In vitro diagnostics			Storage temperature Unopened usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

EUROIMMUN

Medizinische
Laboragnostika
AG



Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below. The thermostat adjusted ELISA incubator must be set at 37°C +/- 1°C.

Coated wells: Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.

- Calibrators and controls: Ready for use. The reagents must be mixed thoroughly before use.
- Enzyme conjugate: Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- Sample buffer: Ready for use.
- Wash buffer: The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength diluted wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- Chromogen/substrate solution: Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- Stop solution: Ready for use.

Warning: The controls used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted 1:101 in sample buffer. For example, dilute 10 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.



Incubation

For semiquantitative analysis incubate calibrator 2 along with the positive and negative controls and patient samples. For quantitative analysis incubate calibrators 1, 2 and 3 along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step)

Transfer 100 µl of the calibrator, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.

For manual processing of microplate wells, cover the finished test plate with the protective foil. When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate 60 minutes at 37°C ± 1°C.

Washing:

Manual: Remove the protective foil and empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Remove the protective foil and empty the wells and subsequently wash 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g. less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells.

When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate 30 minutes at room temperature (+18°C to +25°C)

Empty the wells. Wash as described above.

Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for 15 minutes at room temperature (+18°C to +25°C) (protect from direct sunlight).

Stopping the reaction:
Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm within 30 minutes of adding the stop solution. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry. Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the semiquantitative analysis of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the quantitative analysis of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimizes reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of the calibrator 2. Calculate the ratio according the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$


EUROIMMUN recommends interpreting results as follows:

$$\begin{array}{ll} \text{Ratio } < 0.8: & \text{negative} \\ \text{Ratio } > 0.8 \text{ to } < 1.1: & \text{borderline} \\ \text{Ratio } \geq 1.1: & \text{positive} \end{array}$$



In cases of borderline test results, an additional patient sample should be taken 7 days later and re-tested in parallel with the first patient-sample. The results of both samples allow proper evaluation of titer changes.

Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction of a serum sample lies above the value of calibrator 1 (200 RU/ml). The result should be given as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range of non-infected persons (cut-off value) recommended by EUROIMMUN is 20 relative units (RU)/ml. EUROIMMUN recommends interpreting results as follows:

<16 RU/ml:	negative
≥16 to <22 RU/ml:	borderline
≥22 RU/ml:	positive

Evaluation information: For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class IgG) can be helpful. Diagnosis can be secured by the determination of the titer change in two serum samples taken at an interval of at least 7 days and analysed in parallel.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.



Test characteristics

Calibration: As no quantified international reference serum exists for antibodies against West Nile virus, the calibration is performed in relative units (RU).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibration sera are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is a recombinant, detergent-extracted glycoprotein E of WNV from the membrane fraction of human cells.

Linearity: The linearity of the Anti-West Nile Virus ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95. The Anti-West Nile Virus ELISA (IgG) is linear at least in the tested concentration range (10 RU/ml to 160 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free samples plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-West Nile Virus ELISA (IgG) is 0.4 RU/ml.

Cross reactivity: Cross reactivities to other flaviviruses cannot be excluded. They were recognized with anti-TBE positive and anti-Dengue virus positive samples.

Antibodies against	n	Anti-West Nile Virus ELISA (IgG)
Adenovirus	12	0%
Chlamydia pneumoniae	12	0%
CMV	12	0%
EBV-CA	12	0%
Helicobacter pylori	12	0%
HSV-1	12	0%
Influenza virus A	12	0%
Influenza virus B	12	0%
Measles virus	12	0%
Mumps virus	12	0%
Mycoplasma pneumoniae	12	0%
Parainfluenza virus Pool	12	0%
RSV	12	0%
Rubella virus	12	0%
Toxoplasma gondii	9	0%
VZV	12	0%



Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

Intra-assay variation, n = 20		
Serum	Mean value (RU/ml)	CV (%)
1	101	6.1
2	104	3.3
3	149	2.7

Inter-assay variation, n = 4 x 6		
Serum	Mean value (RU/ml)	CV (%)
1	103	5.1
2	120	4.9
3	171	4.2

Sensitivity and specificity: Samples from 295 patients (origin: Europe) were investigated using the EUROIMMUN Anti-West Nile Virus ELISA and a neutralization test (NT) (performed by RKI, Berlin) as a reference method. The specificity was 96.9%, with a sensitivity of 99.5%. Values for 4 of the samples were borderline and were not included in the calculation.

Reference range: The levels of anti-West Nile virus antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 20 RU/ml, 1.0% of the blood donors were anti-West Nile virus positive (IgG).

Clinical significance

West Nile virus (WNV) is an enveloped single-stranded RNA virus of the Flaviviridae family [1]. This family comprises around 100 virus types that are presently categorized into the three known species Flavivirus, Pestivirus and Hepacivirus [1, 2, 3, 4, 5]. West Nile virus received its name in 1937 when it was first isolated from a blood sample of an elderly woman living in the West Nile district in Uganda, who had fever of unknown cause accompanied by neurological disorders [6]. Further isolates were achieved only in 1951 from the sera of children with weak, unspecific symptoms, namely in Egypt where the virus is endemically distributed. At that time mice and embryonated hen's eggs were used for virus detection [1].

WNV is present not only in tropical areas, but also in moderate climate regions [2, 3, 4, 5]. Significant epidemics were observed in 1951/52 and 1957 in Israel and 1974 and 1983/84 in South Africa [1]. In the Algeria (1994), Rumania (1996/97), the Czech Republic (1997), the Democratic Republic of Congo (1998), Russia, North America (1999) and Israel (2000) [2, 3, 7, 8, 9, 10, 11, 12].

In the USA 149 infections with 18 cases of death were recorded from 1999 to 2001. In 2002 this number rose to 4156 infections and 284 deaths. In 2003 to 9558 infections and 262 deaths [7, 8, 10]. Currently the virus has been detected in seven Canadian provinces, in 48 USA states and in Mexico, as well as in Puerto Rico, the Dominican Republic, Jamaica, Guadeloupe and El Salvador [3, 9, 11].

Since 1958, when antibodies against the WNV were first detected in the sera of two Albanians, repeated outbreaks of West Nile fever have occurred in Southern and Eastern Europe and meanwhile also in Central and West Europe [8, 9, 10, 12, 14, 15]. Its emergence and rapid spread is credited to world climate change, long-distance travel and globalization of economic trade [9, 12, 13, 15, 16, 17]. Consequent monitoring of West Nile activity by controlling sera of exposed persons is essential [15, 17, 18, 19]. Seroprevalence studies in endemic regions have shown an infection spread of up to 40% [14]. West Nile virus is therefore the Flavivirus with the largest distribution area [2, 3, 4, 5, 9].



WNV is transmitted by a number of mosquitoes. In the Mediterranean region and in Africa mosquitoes of the *Culex univittatus* complex species are the main arthropod hosts, while in North America WNV could be detected in 37 mosquito species, with *Culex pipiens* being the main vector [1, 2, 3, 4, 7, 9]. In India *Culex vishnui* and in France *Culex modestus* were identified as the main vectors [1]. In total WNV could be found in more than 40 mosquito and in several tick species [3, 4, 5].

Birds represent the vertebrate reservoir [1, 20, 21]. Alone in the USA WNV has been found in more than 162 species of birds. Many of them showed clinical symptoms and thousands of birds died after contracting a natural infection [7, 8, 10, 20, 21]. Birds that survive develop lifelong immunity. Acting as coincidental hosts mammals can also become infected when bitten by an infected mosquito [10, 13, 14]. Transmission has additionally been documented via breast milk, bone marrow transplants, liver and heart transplants, blood transfusions, lab accidents such as open wounds during handling of infected brain tissue as well as transplacental transmission [2, 3, 4, 5, 8, 22, 23, 24, 25, 26, 27, 28]. Other than humans, mostly only horses became ill after an infection [2, 13, 21]. As well as WNV infection via a mosquito bite, a second natural infection source is possible in animals, namely via feeding on infected prey [3, 5]. An experimental infection of cats was successfully achieved by feeding them infected mice [20, 21].

70% to 80% of the humans infected with WNV showed no symptoms [8]. In the remaining 20%-30%, signs of sudden flu-like symptoms appear after an incubation period of 2-6 days with fever ranging from 38.5 to 40°C lasting for 3-5 days, nausea, shivering, head and back aches, joint and muscle pain and other unspecific symptoms such as loss of appetite, dizziness, vomiting, diarrhoea, coughing and a sore throat [1, 2, 3, 4, 5].

Typical for epidemics occurring fever are exanthema on the breast, back and upper extremities and general lymph node swelling [1, 3, 26]. Severe clinical cases of WNV infections are characterized by myocarditis, pancreatitis and hepatitis and since 1996, also neurological disorders, as WNV is now capable of crossing the blood-brain barrier [3]. The neurological symptoms begin after a short prodrome phase approximately 1-7 days after infection and become manifest in the form of encephalitis and meningoencephalitis accompanied by stiffness, spasms and shivering as the result of damage done to the basal ganglia [3, 8, 29, 30].

Another widespread symptom is general muscle weakness similar to the Guillain-Barre syndrome and also polio-like paralysis [29, 30]. Approximately 4%-14% of the hospitalized patient cases are fatal [29]. High risk factors are old age and a weak immune system [1, 3, 29].

An infection with WNV during pregnancy can cause miscarriages, congenital meningitis, birth defects in approx. 10% of the cases and in an additional 10% of newborns growth disturbances [3, 8, 27, 28].

The diagnosis of WNV can be performed by virus detection or by detection of specific antibodies [1, 4, 17, 18, 19]. As virus isolation from serum or cerebrospinal fluid or virus detection using Real-Time Reverse Transcription-Polymerase Chain Reaction (RT-PCR) is usually unsuccessful due to short viraemia and low virus titers, the detection of specific WNV antibodies using ELISA and IFA has gained importance [11, 17, 18, 19, 31, 32, 33, 34, 35, 36, 37, 38].

Specific IgM antibodies in serum can be determined using ELISA or IIFT [19]. Antibodies of class IgM are detectable in serum from the second day after initial symptoms of the illness occur. A four-fold increase in titer of the respective class of antibody is considered proof of a WNV infection.

If the IgM test is negative, even though the symptoms indicate a WNV infection, a second serum sample should be taken and tested for IgM antibodies a few days later. A combination of ELISA and IIFT provides close to 100% reliability [17, 18, 31, 32, 33, 34]. Anti-WNV IgM antibodies persist for 2 to 3 months, often for more than a year [4, 17, 18, 19, 22].

Antibodies of class IgG are detectable approx. 2 days after the appearance of IgM antibodies [11, 19, 32, 33, 34, 40]. Two to four weeks after a positive IgM result the infection can be confirmed and its severity and prognosis evaluated using a qualitative and quantitative test for the detection of specific WNV IgG antibodies in the patient serum [34].



For the reliable differentiation between acute and past infections the detection of low-avidity IgG antibodies gives evidence for a primary or an acute WNV infection, while high-avidity antibodies indicate a past or reactivated WNV infection [39, 40, 41]. EUROIMMUN offers additional test systems for determination of IgG avidity in both ELISA and IIFT formats. The detection of low-avidity antibodies using ELISA and IIFT in parallel is possible for WNV as it is for *Toxoplasma gondii*, rubella virus, EBV-EA, EBV-CA and Corona virus [39, 40, 41].

As the degree of similarity within the Flavivirus family is high antibody cross reactions can occur [19, 31, 42]. Therefore samples that are positive for specific IgM and/or IgG antibodies against WNV should be titrated and investigated on all relevant Flavivirus IIFT substrates for cross reactions. By comparing the titer strengths the initial result can be confirmed or disproved by the second detection and an infection with another Flavivirus identified as the source of illness [17, 32, 34].

To supplement and extend the current Anti-West Nile Virus ELISA and Anti-West Nile Virus IIFT (each IgG or IgM or avidity) BIOCHIP Mosaic and Profiles for the detection (IIFT) of infections with Flaviviruses and the BIOCHIP Mosaic Fever Profile 1: South-East Asia have been developed. With these tests specific antibodies (IgG and IgM) against several infectious agents can be investigated simultaneously [36, 37, 38, 43, 44, 45]. These supplementary tests allow similar or ambiguous disease symptoms and potential cross reactions to be clarified and differential diagnostic issues to be addressed [1].

A specific antiviral therapy for WNV encephalitis is not available at present [1, 2, 3, 4, 5, 14, 35]. Intensive medical care is the only possibility to positively influence the illness. Eradication of WNV is impossible due to the natural bird-mosquito cycle [1]. A vaccine with formalin inactivated WNV is available for horses [2, 10]. Therefore public education, individual precautionary measures and protection against insect bites are essential contributions to preventing WNV infections [2, 3, 4, 5, 14, 16].

Literature references

- Falke D, Flavi-Viren. In: Hahn H, Falke D, Kaufmann SHE, Ullmann U. *Medizinische Mikrobiologie und Infektiologie*. Springer-Verlag; 5. Auflage (2004) 523-529.
- Campbell GL, Marfin AA, Lanciotti RS, Gubler DJ. West Nile Virus. *Lancet Infect Dis* 2 (2002) 519-529.
- Nosal B, Pellizzari R. West Nile virus. *CMAJ* 168 (2003) 1443-1444.
- Gernwehr BP, Lillbridge KM, Higgs S, Mason PW, Aronson JF, Campbell GA, Barrett ADT. West Nile virus: where are we now? *The Lancet Infect Dis* 4 (2004) 547-551.
- Stock I. West Nile virus. An unusual flavivirus with increased importance. [Article in German] *Med Monatsschr Pharm* 27 (2004) 112-120.
- Smithburn KC, Hughes TP, Burke AW et al. A neurotropic virus isolated from the blood of a native of Uganda. *Am J Trop Med Hyg* 20 (1940) 471-472.
- Kilpatrick AM, Kramer LD, Jones MJ, Marra PP, Daszak P. West Nile virus epidemics in North America are driven by shifts in mosquito feeding behavior. *PLoS Bio* 4 (2006) 82.
- Zeller HG, Schuffenecker I. West Nile virus: an overview of its spreads in Europe and the Mediterranean basin in contrast to its spread in the Americas. *Eur J Clin Microbiol Infect Dis* 23 (2004) 147-156.
- Rappole JH, Compton BW, Leimgruber P, Robertson J, King DJ, Renner SC. Modeling movement of West Nile virus in the Western Hemisphere. *Vector Borne Zoonotic Dis* 6 (2006) 128-139.
- Hubalek Z. European experience with the West Nile virus ecology and epidemiology: could it be relevant for the New World? *Viral Immunol* 13 (2000) 415-426.
- Tilley PA, Wallie R, Chow A, Jayaraman GC, Fonseca K, Drebot MA, Prekasitis J, Fox J. Clinical utility of commercial enzyme immunoassays during the inaugural season of West Nile virus activity, Alberta, Canada. *J Clin Microbiol* 43 (2005) 4691-4695.
- Hubalek Z, Halouzka J. West Nile-Fever – a re-emerging mosquito-borne viral disease in Europe. *Emerg Infect Dis* 5 (1999) 643-650.
- Pauli G. West Nile virus. Prevalence and significance as a zoonotic pathogen. [Article in German] *Bundesgesundheitsblatt Gesundheitsforschung Gesundheitsschutz* 47 (2004) 655-660.
- Robert-Koch-Institut. *West-Nil-Fieber*. *Epidemiol Bull* 13 (2003) 96-99.
- Epslein PR. West Nile virus and the climate. *J Urban Health* 78 (2001) 367-371.
- Gubler DJ. The global emergence/resurgence of arboviral diseases as public health problems. *Arch Med Res* 33 (2002) 330-342.
- Koraka P, Zeller H, Niedrig M, Osterhaus AD, Groen J. Reactivity of serum from patients with flavivirus infection measured by immunofluorescence assay and ELISA. *Microbes Infect* 12 (2002) 1209-1215.
- Shi PY, Wong SJ. Serologic diagnosis of West-Nile virus infection. *Expert Rev Mol Diagn* 6 (2003) 733-741.
- Tardai G, Rula S, Chitu V, Rossi C, Tsai TF, Cernescu C. Evaluation of immunoglobulin M (IgM) and IgG enzyme immunoassays in serologic diagnosis of West Nile virus infection. *J Clin Microbiol* 38 (2000) 2232-2239.
- Rappole JH, Hubalek Z. Migratory birds and West Nile virus. *J Appl Microbiol* 94 (2003) 47-58.
- Kilpatrick AM, Daszak P, Jones MJ, Marra PP, Kramer LD. Host heterogeneity dominates West Nile virus transmission. *Proc Biol Sci* 22 (2006) 2321-2333.
- Prince HE, Tobler LH, Lape-Nixon M, Foster GA, Stramer SL, Busch MP. Development and persistence of West Nile virus-specific immunoglobulin M (IgM), IgA, and IgG in viremic blood donors. *J Clin Microbiol* 43 (2005) 4316-4320.
- Iwamoto M, Jennigan DB, Glasch A, Treppka MJ, Blackmore CG, Hellinger WC, Pham SM, Zaki S, Lanciotti RS, Lance-Parker SE, DiazGranados CA, Winquist AG, Perrino CA, Wiersma S, Hillier KL, Goodman JL, Marfin AA, Chamberland ME, Petersen LR, West Nile Virus in Transplant Recipients Investigation Team. Transmission of West Nile virus from an organ donor to four transplant recipients. *N Engl J Med* 348 (2003) 2196-2203.
- Montgomery SP, Brown JA, Kluennert M, Smith TL, Crall N, Lanciotti RS, Macedo de Oliveira A, Bo T, Marfin AA. 2003 West Nile Virus Transmission-Associated Transmission Investigation Team. Transfusion-associated transmission of West Nile virus, United States 2003 through 2005. *Transfusion* 46 (2006) 2038-2046.
- Cameron C, Reeves J, Antonishyn N, Tilley P, Alport T, Towns D, Lane D, Saldanha J. West Nile virus in Canadian blood donors. *Transfusion* 45 (2005) 487-491.
- Anderson RG, Horn KB, Hoang MP, Gottlieb E, Benin B. Punctate exanthem of West Nile Virus infection: report of 3 cases. *J Am Acad Dermatol* 51 (2004) 820-823.
- O'Leary DR, Kuhn S, Kniss KL, Hinkley AF, Rasmussen SA, Page WJ, Kightlinger LK, Beecham BD, Miller TK, Neitzel DF, Michaels SR, Campbell GL, Lanciotti RS, Hayes EB. Birth outcomes following West Nile Virus infection of pregnant women in the United States: 2003-2004. *Pediatrics* 117 (2006) 537-545.
- Pastley JE, Hinkley AF, O'Leary DR, Kramer WC, Lanciotti RS, Campbell GL, Hayes EB. West Nile virus infection among pregnant women in a northern Colorado community, 2003 to 2004. *Pediatrics* 117 (2006) 814-820.
- Guarner J, Shieh WJ, Hunter S, Paddock CD, Morken T, Campbell GL, Marfin AA, Zaki SR. Clinicopathologic study and laboratory diagnosis of 23 cases with West Nile virus encephalomyelitis. *Hum Pathol* 35 (2004) 983-990.
- Bouffard JP, Riudavets MA, Holman R, Rushing EJ. Neuropathology of the brain and spinal cord in human West Nile virus infection. *Clin Neuropathol* 23 (2004) 59-61.






31. Hogrefe WR, Moore R, Lape-Nixon M, Wagner M, Prince HE. Performance of immunoglobulin G (IgG)- and IgM-enzyme-linked immunosorbent assays using a West Nile virus recombinant antigen (preM/E) for detection of West Nile virus- and other flavivirus-specific antibodies. *J Clin Microbiol* 42 (2004) 4641-4648.
32. Niedig M, Sonnenberg* K, Steinhagen* K, Paweska JT (*EUROIMMUN AG). Comparison of ELISA and immunoassays for measurement of IgG and IgM antibody to West Nile virus in human sera against virus neutralisation. *J Virol Meth* in press (2006).
33. Malen AK, Martins TB, Hill HR, Litwin CM. Evaluations of commercial West Nile virus immunoglobulin G (IgG) and IgM enzyme immunoassays show the value of continuous validation. *J Clin Microbiol* 42 (2004) 727-733.
34. Malen AK, Stipanovich PJ, Martins TB, Hill HR, Litwin CM. Detection of IgG and IgM to West Nile virus. Development of an immunofluorescence assay. *Am J Clin Pathol* 119 (2003) 508-515.
35. Sanchez MD, Pierson TC, McAllister D, Hanna SL, Purfer BA, Valentine LE, Murtadha MM, Hoxie JA, Doms RW. Characterization of neutralizing antibodies to West Nile virus. *Virology* 25 (2005) 70-82.
36. EUROIMMUN AG. Stöcker W, Wessel S, Morrin M, Müller M. Konstante Lichtquelle für die Fluoreszenzmikroskopie. Deutsche Patentanmeldung (Offenlegungsschrift) DE 10 2006 027 518.7 (2006).
37. Stöcker* W, Fauer* H, Krause* C, Barth E, Martinez A (*EUROIMMUN AG). Verfahren zur Optimierung der automatischen Fluoreszenzerkennung in der Immun Diagnostik. Deutsche Patentanmeldung (Offenlegungsschrift) DE 10 2006 027 516.0 und WO2007/140952 (2006).
38. EUROIMMUN AG. Stöcker W. Process for analyses to be carried out on immobilized biological tissue. USA-Patent Nr. US4647543 (1983/1987).
39. Lewett PN, Sonnenberg* K, Sidaway F, Shead S, Niedig M, Steinhagen* K, Horsman GB, Drebot MA (*EUROIMMUN AG). Use of Immunoglobulin G Avidity Assays for Differentiation of Primary from Previous Infections with West Nile Virus. *Journal of Clinical Microbiology* Vol 43 (2005) 5873-5875.
40. Fox JL, Hazell SL, Tobler LH, Busch MP. Immunoglobulin G avidity in differentiation between early and late antibody responses to West Nile virus. *Clin Vaccine Immunol* 13 (2006) 33-36.
41. Lewett PN, Sonnenberg* K, Shead S, Sidaway F, Niedig M, Kreher* P, Arp* B, Erdmann* S, Jahn A, Steinhagen* K, Drebot M, Horsman GB (*EUROIMMUN AG). IgG avidity determinations but not IgM tests are capable to differentiate recent from previous infections with West Nile Virus. Poster zum 21. Annual clinical virology symposium and annual meeting Pan American society for clinical virology, USA, Florida, Clearwater beach (2005).
42. Stasny K, Kiermayr S, Heinz FX. Cryptic properties of a cluster of dominant flavivirus cross-reactive antigenic sites. *J Virol* 80 (2006) 9557-9568.
43. Stöcker* W, Horbach E, Slenczka W, Schaper J, Böttger K, Herrmann C (*EUROIMMUN AG). Automatische Bestückung von Objektträgern mit „BIOCHIPs“: Eine neue Technologie für die indirekte Immunfluoreszenz. *Lab. med.* (1989) 171-172.
44. Sonnenberg* K, Hinrichs E, Müller-Kunert* E, Schlumberger* W, Stöcker* W (*EUROIMMUN AG). Einsatz der BIOCHIP-Technologie in der Serologie von Influenza-, EBV- und Legionella-Infektionen. Poster zum 3. Deutschen Kongress für Infektions- und Tropenmedizin, 15. bis 18. März 1995 in Berlin-Hohenschönhausen. Abstrakt im Kongressband (1995).
45. EUROIMMUN AG. Stöcker W. Die BIOCHIP-Technologie. Ein neuer Weg der Produktion Festphase-gebundener Bioreagenzien für das medizinische und biologische Laboratorium. Dokumentation zum Vorschlag des Rektors der Medizinischen Universität Lübeck für den Schmidt-Römhild-Technologiepreis 1989. Eigenverlag der Firma EUROIMMUN (1990).

Anti-West Nile Virus ELISA (IgM) Test Instruction

ORDER NO.	ANTIBODIES AGAINST West Nile Virus	IG CLASS IgM	SUBSTRATE Ag-coated microplate wells	FORMAT 96 x 01 (96)
EL 2662-9601 M				

Principle of the test: The ELISA test kit provides a semiquantitative in vitro assay for human antibodies of the IgM class against West Nile virus in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with West Nile virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgM (also IgA and IgG) antibodies will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgM (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens, 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	[STRIPS]
2. Calibrator (IgM, human), ready for use	dark red	1 x 2.0 ml	[CAL]
3. Positive control (IgM, human), ready for use	blue	1 x 2.0 ml	[POS CONTROL]
4. Negative control (IgM, human), ready for use	green	1 x 2.0 ml	[NEG CONTROL]
5. Enzyme conjugate peroxidase-labelled anti-human IgM (goat), ready for use	red	1 x 12 ml	[CONJUGATE]
6. Sample buffer buffer containing IgG/Rf-Absorbent (Anti-human IgG antibody preparation obtained from goat), ready for use	green	1 x 100 ml	[SAMPLE BUFFER]
7. Wash buffer 10x concentrate	colourless	1 x 100 ml	[WASH BUFFER 10x]
8. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	[SUBSTRATE]
9. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	[STOP SOLUTION]
10. Test instruction	---	1 booklet	
11. Quality control certificate	---	1 protocol	
12. Protective foil	---	3 pieces	
LOT	Storage temperature		
IVD	In vitro diagnostics  Unopened usable until		

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

EUROIMMUN

Medizinische
Labor Diagnostika
AG



Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below. The thermostat adjusted ELISA incubator must be set at 37°C +/- 1°C.

Coated wells: Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag). Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.

Calibrator and controls: Ready for use. The reagents must be mixed thoroughly before use.

Enzyme conjugate: Ready for use. The enzyme conjugate must be mixed thoroughly before use.

Sample buffer: Ready for use. The green coloured sample buffer contains IgG/Rf absorbent. Serum or plasma samples diluted with this sample buffer are only to be used for the determination of IgM antibodies.

Wash buffer: The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).

For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water. The working strength diluted wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.

Chromogen/substrate solution: Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.

Stop solution: Ready for use.

Warning: The control sera used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.



Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Introduction: Before the determination of specific antibodies of class IgM, antibodies of class IgG should be removed from the patient sample. This procedure must be carried out in order to prevent any rheumatoid factors from reacting with specifically bound IgG, which would lead to false positive IgM test results, and to prevent specific IgG displacing IgM from the antigen, which would lead to false IgM negative test results.

Functional principle: The sample buffer (green coloured) contains an anti-human antibody preparation from goat. IgG from a serum sample is bound with high specificity by these antibodies and precipitated. If the sample also contains rheumatoid factors, these will be absorbed by the IgG/anti-human IgG complex.

Separation properties:

- All IgG subclasses are bound and precipitated by the anti-human IgG antibodies.
- Human serum IgG in concentrations of up to 15 mg per ml are removed (average serum IgG concentration in adults: 12 mg per ml).
- Rheumatoid factors are also removed.
- The recovery rate of the IgM fraction is almost 100%.

Performance: The patient samples for analysis are diluted 1:101 with sample buffer. For example, add 10 µl serum to 1.0 ml sample buffer and mix well. Incubate the mixture for at least 10 minutes at room temperature. Subsequently, it can be pipetted into the microplate wells according to the pipetting protocol.

Notes:

- Antibodies of the class IgG should not be analyzed with this mixture.
- It is possible to check the efficacy of the IgG/Rf absorbent for an individual patient sample by performing an IgG test in parallel to the IgM test using the mixture. If the IgG test is negative, the IgM result can be considered as reliable.
- The calibrator and controls containing IgM antibodies are pre-diluted and ready for use, do not dilute them.



Incubator

(Partly) manual test performance

Sample incubation:
(1st step)

Transfer 100 µl of the calibrator, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.

For manual processing of microplate wells, cover the finished test plate with the protective foil. When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.
Incubate 60 minutes at 37°C ± 1°C.

Washing:

Manual: Remove the protective foil and empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Remove the protective foil and empty the wells and subsequently wash 3 times with 450 µl of working strength wash buffer (program setting e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.
Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgM) into each of the microplate wells.

When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate 30 minutes at room temperature (+18°C to +25°C).

Empty the wells. Wash as described above.

Substrate incubation:
(3. step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for 15 minutes at room temperature (+18°C to +25°C) (protect from direct sunlight).

Stopping the reaction:
(3rd step)

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm within 30 minutes of adding the stop solution. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the Analyzer I, Analyzer I-2P and the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry. Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C	P 6	P 14	P 22								
B	pos.	P 7	P 15	P 23								
C	neg.	P 8	P 16	P 24								
D	P 1	P 9	P 17									
E	P 2	P 10	P 18									
F	P 3	P 11	P 19									
G	P 4	P 12	P 20									
H	P 5	P 13	P 21									

The above pipetting protocol is an example of the semiquantitative analysis of antibodies in 24 patient samples (P 1 to P 24).

Calibrator (C), positive (pos.) and negative (neg.) control as well as the patient samples have been incubated in one well each. The reliability of the ELISA test can be improved by duplicate determinations of each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimizes reagent wastage. Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

The extinction value of the calibrator defines the upper limit of the reference range of non-infected persons (cut-off) recommended by EUROIMMUN. Values above the indicated cut-off are to be considered as positive, those below as negative.

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator. Use the following formula to calculate the ratio:

$$\frac{\text{Extinction of the control/patient sample}}{\text{Extinction of calibrators}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8: negative
Ratio ≥ 0.8 to <1.1: borderline
Ratio ≥ 1.1 : positive



Evaluation information: For duplicate determinations the mean of the two-values should be taken, if the two values deviate substantially from one another the sample should be retested.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class (IgG) can be helpful. Diagnosis can be secured by the determination of the titer change in two serum samples taken at an interval of at least 7 days and analysed in parallel.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

Test characteristics

Calibration: As no quantified international reference serum exists for antibodies of the IgM class against West Nile virus, results are provided in the form of ratios which are a relative measure for the concentration of antibodies.

For every group of tests performed, the extinction values of the calibrators and the ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the control sera are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is a recombinant, detergent-extracted glycoprotein E of WNV from the membrane fraction of human cells.

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-West Nile ELISA (IgM) is a ratio value of 0.03.

Cross reactivity: Cross reactivities to other flaviviruses cannot be excluded. They were recognized with anti-TBE positive and anti-Dengue virus positive samples.

Antibodies against	n	Anti-West Nile ELISA (IgM)
Borrelia burgdorferi	9	0%
CMV	8	0%
EBV-CA	9	0%
HSV-1/2	2	0%
Measles virus	10	0%
Mumps virus	9	0%
Parvovirus B 19	8	0%
Rubella virus	10	0%
Toxoplasma gondii	10	0%
VZV	4	0%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.



Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

Intra-assay variation, n = 20		
Serum	Mean value (Ratio)	CV (%)
1	1.9	5.9
2	2.0	2.2
3	2.4	4.0

Inter-assay variation, n = 4 x 6		
Serum	Mean value (Ratio)	CV (%)
1	1.9	3.9
2	2.1	9.4
3	2.4	6.8

Sensitivity and specificity: Study I: For the determination of the sensitivity 18 clinically and serologically precharacterised sera (Robert Koch Institute, Berlin, Germany) were tested with the EUROIMMUN Anti-West Nile Virus ELISA (IgM). The sensitivity amounted to 94.4%.

Study II: 99 patient sera, characterised as positive at the Saskatchewan Disease Control Laboratory (Canada) using several serological methods, were investigated with the EUROIMMUN Anti-West Nile Virus ELISA (IgM). The sensitivity was 93.7%. For the determination of the specificity 500 blood donor samples were investigated. The specificity was 99.8%.

Reference range: The levels of anti-West Nile virus antibodies (IgM) were analyzed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off ratio of 1.0, 0.2% of the blood donors were anti-West Nile virus positive (IgM).

Clinical significance

West Nile virus (WNV) is an enveloped single-stranded RNA virus of the Flaviviridae family [1]. This family comprises around 100 virus types that are presently categorized into the three known species Flavivirus, Pestivirus and Hepacivirus [1, 2, 3, 4, 5]. West Nile virus received its name in 1937 when it was first isolated from a blood sample of an elderly woman living in the West Nile district in Uganda, who had fever of unknown cause accompanied by neurological disorders [6]. Further isolates were achieved only in 1951 from the sera of children with weak, unspecific symptoms, namely in Egypt where the virus is endemically distributed. At that time mice and embryonated hen's eggs were used for virus detection [1].

WNV is present not only in tropical areas, but also in moderate climate regions [2, 3, 4, 5]. Significant epidemics were observed in 1951/52 and 1957 in Israel and 1974 and 1983/84 in South Africa [1]. In the mid 90's the virus changed its virulence causing an epidemic accumulation of WNV encephalitis in Algeria (1994), Rumania (1996/97), the Czech Republic (1997), the Democratic Republic of Congo (1998), Russia, North America (1999) and Israel (2000) [2, 3, 7, 8, 9, 10, 11, 12]. In the USA 149 infections with 18 cases of death were recorded from 1999 to 2001. In 2002 this number rose to 4156 infections and 284 deaths, in 2003 to 9858 infections and 262 deaths [7, 8, 10]. Currently the virus has been detected in seven Canadian provinces, in 48 USA states and in Mexico, as well as in Puerto Rico, the Dominican Republic, Jamaica, Guadeloupe and El Salvador [3, 9, 11].

Since 1958, when antibodies against the WNV were first detected in the sera of two Albanians, repeated outbreaks of West Nile fever have occurred in Southern and Eastern Europe and meanwhile also in Central and West Europe [8, 9, 10, 12, 14, 15]. Its emergence and rapid spread is credited to world climate change, long-distance travel and globalization of economic trade [9, 12, 13, 15, 16, 17]. Consequent monitoring of West Nile activity by controlling sera of exposed persons is essential [15, 17, 18, 19]. Seroprevalence studies in endemic regions have shown an infection spread of up to 40% [14]. West Nile virus is therefore the Flavivirus with the largest distribution area [2, 3, 4, 5, 9].



WNV is transmitted by a number of mosquitoes. In the Mediterranean region and in Africa mosquitoes of the *Culex univittatus* complex species are the main arthropod hosts, while in North America WNV could be detected in 37 mosquito species, with *Culex pipiens* being the main vector [1, 2, 3, 4, 7, 9]. In India *Culex vishnui* and in France *Culex modestus* were identified as the main vectors [1]. In total WNV could be found in more than 40 mosquito and in several tick species [3, 4, 5].

Birds represent the vertebrate reservoir [1, 20, 21]. Alone in the USA WNV has been found in more than 162 species of birds. Many of them showed clinical symptoms and thousands of birds died after contracting a natural infection [7, 8, 10, 20, 21]. Birds that survive develop lifelong immunity. Acting as coincidental hosts mammals can also become infected when bitten by an infected mosquito [10, 13, 14]. Transmission has additionally been documented via breast milk, bone marrow transplants, liver and heart transplants, blood transfusions, lab accidents such as open wounds during handling of infected brain tissue as well as transplacental transmission [2, 3, 4, 5, 8, 22, 23, 24, 25, 26, 27, 28]. Other than humans, mostly only horses became ill after an infection [2, 13, 21]. As well as WNV infection via a mosquito bite, a second natural infection source is possible in animals, namely via feeding on infected prey [3, 5]. An experimental infection of cats was successfully achieved by feeding them infected mice [20, 21].

70% to 80% of the humans infected with WNV showed no symptoms [8]. In the remaining 20%-30%, signs of sudden flu-like symptoms appear after an incubation period of 2-6 days with fever ranging from 38.5 to 40°C, lasting for 3-5 days, nausea, shivering, head and back aches, joint and muscle pain and other unspecific symptoms such as loss of appetite, dizziness, vomiting, diarrhoea, coughing and a sore throat [1, 2, 3, 4, 5].

Typical for epidemics occurring fever are exanthema on the breast, back and upper extremities and general lymph node swelling [1, 3, 26]. Severe clinical cases of WNV infections are characterized by myocarditis, pancreatitis and hepatitis and since 1996, also neurological disorders, as WNV is now capable of crossing the blood-brain barrier [3]. The neurological symptoms begin after a short prodrome phase approximately 1-7 days after infection and become manifest in the form of encephalitis and meningoencephalitis accompanied by stiffness, spasms and shivering as the result of damage done to the basal ganglia [3, 8, 29, 30].

Another widespread symptom is general muscle weakness similar to the Guillain-Barré syndrome and also polio-like paralysis [29, 30]. Approximately 4%-14% of the hospitalized patient cases are fatal [29]. High risk factors are old age and a weak immune system [1, 3, 29].

An infection with WNV during pregnancy can cause miscarriages, congenital meningitis, birth defects in approx. 10% of the cases and in an additional 10% of newborns growth disturbances [3, 8, 27, 28].

The diagnosis of WNV can be performed by virus detection or by detection of specific antibodies [1, 4, 17, 18, 19]. As virus isolation from serum or cerebrospinal fluid or virus detection using Real-Time Reverse Transcription-Polymerase Chain Reaction (RT-PCR) is usually unsuccessful due to short viraemia and low virus titers, the detection of specific WNV antibodies using ELISA and IFA has gained importance [11, 17, 18, 19, 31, 32, 33, 34, 35, 36, 37, 38].

Specific IgM antibodies in serum can be determined using ELISA or IIFT [19]. Antibodies of class IgM are detectable in serum from the second day after initial symptoms of the illness occur. A four-fold increase in titer of the respective class of antibody is considered proof of a WNV infection.

If the IgM test is negative, even though the symptoms indicate a WNV infection, a second serum sample should be taken and tested for IgM antibodies a few days later. A combination of ELISA and IIFT provides close to 100% reliability [17, 18, 31, 32, 33, 34]. Anti-WNV IgM antibodies persist for 2 to 3 months, often for more than a year [4, 17, 18, 19, 22].

Antibodies of class IgG are detectable approx. 2 days after the appearance of IgM antibodies [11, 19, 32, 33, 34, 40]. Two to four weeks after a positive IgM result the infection can be confirmed and its severity and prognosis evaluated using a qualitative and quantitative test for the detection of specific WNV IgG antibodies in the patient serum [34].

For the reliable differentiation between acute and past infections the detection of low-avidity IgG antibodies gives evidence for a primary or an acute WNV infection, while high-avidity antibodies indicate a past or reactivated WNV infection [39, 40, 41].



EUROIMMUN offers additional test systems for determination of IgG avidity in both ELISA and IIFT formats. The detection of low-avidity antibodies using ELISA and IIFT in parallel is possible for WNV as it is for *Toxoplasma gondii*, rubella virus, EBV-EA, EBV-CA and Corona virus [39, 40, 41].

As the degree of similarity within the Flavivirus family is high antibody cross reactions can occur [19, 31, 42]. Therefore samples that are positive for specific IgM and/or IgG antibodies against WNV should be titrated and investigated on all relevant Flavivirus IIFT substrates for cross reactions. By comparing the titer strengths the initial result can be confirmed or disproved by the second detection and an infection with another Flavivirus identified as the source of illness [17, 32, 34].

To supplement and extend the current Anti-West Nile Virus ELISA and Anti-West Nile Virus IIFT (each IgG or IgM or avidity) BIOCHIP Mosaics and Profiles for the detection (IIFT) of infections with Flaviviruses and the BIOCHIP Mosaic Fever Profile 1: South-East Asia have been developed with these tests specific antibodies (IgG and IgM) against several infectious agents can be investigated simultaneously [36, 37, 38, 43, 44, 45]. These supplementary tests allow similar or ambiguous disease symptoms and potential cross reactions to be clarified and differential diagnostic issues to be addressed [1].

A specific antiviral therapy for WNV encephalitis is not available at present [1, 2, 3, 4, 5, 14, 35]. Intensive medical care is the only possibility to positively influence the illness. Eradication of WNV is impossible due to the natural bird-mosquito cycle [1]. A vaccine with formalin inactivated WNV is only available for horses [2, 10]. Therefore public education, individual precautionary measures and protection against insect bites are essential contributions to preventing WNV infections [2, 3, 4, 5, 14, 16].

Literature references

- Falke D, Flavi-Viren. In: Hahn H, Falke D, Kaufmann SHE, Ullmann U. *Medizinische Mikrobiologie und Infektiologie*. Springer-Verlag; 5. Auflage (2004) 523-529.
- Campbell GL, Marfin AA, Lanciotti RS, Gubler DJ. West Nile Virus. *Lancet Infect Dis* 2 (2002) 519-529.
- Nosal B, Pellizzari R. West Nile virus. *CMAJ* 168 (2003) 1443-1444.
- Granwehr BP, Lillbridge KM, Higgs S, Mason PW, Aronson JF, Campbell GA, Barrett ADT. West Nile virus: where are we now? *The Lancet Infect Dis* 4 (2004) 547-551.
- Stoek I. West Nile virus. An unusual flavivirus with increased importance. [Article in German] *Med Monatsschr Pharm* 27 (2004) 112-120.
- Smithburn KC, Hughes TP, Burke AW et al. A neurotropic virus isolated from the blood of a native of Uganda. *Am J Trop Med Hyg* 20 (1940) 471-472.
- Kilpatrick AM, Kramer LD, Jones MJ, Marra PP, Daszak P. West Nile virus epidemics in North America are driven by shifts in mosquito feeding behavior. *PLoS Bio* 4 (2006) 82.
- Zeller HG, Schuffenecker I. West Nile virus: an overview of its spreads in Europe and the Mediterranean basin in contrast to its spread in the Americas. *Eur J Clin Microbiol Infect Dis* 23 (2004) 147-156.
- Rappole JH, Compton BW, Leimgruber P, Robertson J, King DJ, Renner SC. Modeling movement of West Nile virus in the Western Hemisphere. *Vector Borne Zoonotic Dis* 6 (2006) 128-139.
- Hubalek Z. European experience with the West Nile virus ecology and epidemiology: could it be relevant for the New World? *Viral Immunol* 13 (2000) 445-426.
- Tilley PA, Walle R, Chow A, Jayaraman GC, Fonseca K, Drebot MA, Preiksaitis J, Fox J. Clinical utility of commercial enzyme immunoassays during the inaugural season of West Nile virus activity, Alberta, Canada. *J Clin Microbiol* 43 (2005) 4691-4695.
- Hubalek Z, Halouzka J. West Nile Fever – a re-emerging mosquito-borne viral disease in Europe. *Emerg Infect Dis* 5 (1999) 643-650.
- Pauli G. West Nile virus: Prevalence and significance as a zoonotic pathogen. [Article in German] *Bundesgesundheitsblatt Gesundheitsforschung Gesundheitsschutz* 47 (2004) 653-660.
- Robert-Koch-Institut. *West-Nil-Fieber*. Epidemiol Bull 13 (2003) 96-99.
- Erstein PR. West Nile virus and the climate. *J Urban Health* 78 (2001) 367-371.
- Gubler DJ. The global emergence/resurgence of arboviral diseases as public health problems. *Arch Med Res* 33 (2002) 330-342.
- Koraka P, Zeller H, Niedrig M, Osterhaus AD, Groen J. Reactivity of serum from patients with flavivirus infection measured by immunofluorescence assay and ELISA. *Microbes Infect* 12 (2002) 1209-1215.
- Shi PY, Wong SJ. Serologic diagnosis of West-Nile virus infection. *Expert Rev Mol Diagn* 6 (2003) 733-741.
- Tardei G, Ruta S, Chiti V, Rossi C, Tsai TF, Cernescu C. Evaluation of immunoglobulin M (IgM) and IgG enzyme immunoassays in serologic diagnosis of West Nile virus infection. *J Clin Microbiol* 38 (2000) 2232-2239.
- Rappole JH, Hubalek Z. Migratory birds and West Nile virus. *J Appl Microbiol* 94 (2003) 47-58.
- Kilpatrick AM, Daszak P, Jones MJ, Marra PP, Kramer LD. Host heterogeneity dominates West Nile virus transmission. *Proc Biol Sci* 22 (2006) 2327-2333.
- Prince HE, Tobler LH, Lape-Nixon M, Foster GA, Stramer SL, Busch MP. Development and persistence of West Nile virus-specific immunoglobulin M (IgM), IgA, and IgG in viremic blood donors. *J Clin Microbiol* 43 (2005) 4316-4320.
- Iwanoto M, Jernigan DB, Guasch A, Trejka MJ, Blackmore CG, Hellinger WC, Pharm SM, Zaki S, Lanciotti RS, Lanco-Parker SE, DiazGranados CA, Winquist AG, Perino CA, Wiersma S, Hillier KL, Goodman JL, Marfin AA, Chalmazand ME, Petersen LR. West Nile Virus in Transplant Recipients Investigation Team. Transmission of West Nile Virus from an organ donor to four transplant recipients. *N Engl J Med* 348 (2003) 2196-2203.
- Montgomery SP, Brown JA, Kuehnert M, Smith TL, Crall N, Lanciotti RS, Macedo de Oliveira A, Boo T, Marfin AA. 2003 West Nile Virus Transfusion-Associated Transmission Investigation Team. Transfusion-associated transmission of West Nile virus, United States 2003 through 2005. *Transfusion* 46 (2006) 2038-2046.
- Cameron C, Reeves J, Antonishyn N, Tilley P, Apport T, Towns D, Lane D, Saldanha J. West Nile virus in Canadian blood donors. *Transfusion* 45 (2005) 487-491.
- Anderson RC, Horn KB, Hoang MP, Gottlieb E, Benmin B. Punctate exanthem of West Nile Virus infection: report of 3 cases. *J Am Acad Dermatol* 51 (2004) 820-823.
- O'Leary DR, Kuhn S, Kniss KL, Hinckley AF, Rasmussen SA, Pape WJ, Kightlinger LK, Beecham BD, Miller TK, Neizel DF, Michels SR, Campbell GL, Lanciotti RS, Hayes EB. Birth outcomes following West Nile Virus infection of pregnant women in the United States. 2003-2004. *Pediatrics* 117 (2006) 537-545.
- Paisley JE, Hinckley AF, O'Leary DR, Kramer WC, Lanciotti RS, Campbell GL, Hayes EB. West Nile virus infection among pregnant women in a northern Colorado community, 2003 to 2004. *Pediatrics* 117 (2006) 814-820.
- Gartner J, Shieh WJ, Hunter S, Paddock CD, Morken T, Campbell GL, Marfin AA, Zaki SR. Clinicopathologic study and laboratory diagnosis of 23 cases with West Nile virus encephalomyelitis. *Hum Pathol* 35 (2004) 983-990.
- Bouffard JP, Ruidavets MA, Holman R, Rushing EJ. Neuropathology of the brain and spinal cord in human West Nile virus infection. *Clin Neuropathol* 23 (2004) 59-61.





31. Hogrefe WR, Moore R, Lape-Nixon M, Wagner M, Prince HE. Performance of immunoglobulin G (IgG) and IgM enzyme-linked immunosorbent assays using a West Nile virus recombinant antigen (preM/E) for detection of West Nile virus- and other flavivirus-specific antibodies. *J Clin Microbiol* 42 (2004) 4641-4648.
32. Niedrig M, Sonnenberg* K, Steinhagen* K, Paweska JT (*EUROIMMUN AG). Comparison of ELISA and immunoassays for measurement of IgG and IgM antibody to West Nile virus in human sera against virus neutralisation. *J Virol Meth* in press (2006).
33. Malan AK, Martins TB, Hill HR, Litwin CM. Evaluations of commercial West Nile virus immunoglobulin G (IgG) and IgM enzyme immunoassays show the value of continuous validation. *J Clin Microbiol* 42 (2004) 727-733.
34. Malan AK, Stipanovich PJ, Martins TB, Hill HR, Litwin CM. Detection of IgG and IgM to West Nile virus. Development of an immunofluorescence assay. *Am J Clin Pathol* 119 (2003) 508-515.
35. Sanchez MD, Pierson TC, McAllister D, Hanna SL, Puffer BA, Valentine LE, Muradha MM, Hoxie JA, Doms RW. Characterization of neutralizing antibodies to West Nile virus. *Virology* 25 (2005) 70-82.
36. EUROIMMUN AG. Stöcker W, Wessel S, Morrin M, Müller M, Konstante Lichtquelle für die Fluoreszenzmikroskopie. Deutsche Patentanmeldung (Offenlegungsschrift) DE 10 2006 027 518.7 (2006).
37. Stöcker* W, Fauer* H, Krause* C, Barth E, Martinez A (*EUROIMMUN AG). Verfahren zur Optimierung der automatischen Fluoreszenzerkennung in der Immun Diagnostik. Deutsche Patentanmeldung (Offenlegungsschrift) DE 10 2006 027 516.0 und WO2007/40952 (2006).
38. EUROIMMUN AG. Stöcker W. Process for analyses to be carried out on immobilized biological tissue. USA-Patent Nr. US4647543 (1983/1987).
39. Levett PN, Sonnenberg* K, Sidaway F, Shead S, Niedrig M, Steinhagen* K, Horsman GB, Drebot MA (*EUROIMMUN AG). Use of Immunoglobulin G Avidity Assays for Differentiation of Primary from Previous Infections with West Nile Virus. *Journal of Clinical Microbiology* Vol 43 (2005) 5873-5875.
40. Fox JL, Hazell SL, Tobler LH, Busch MP. Immunoglobulin G avidity in differentiation between early and late antibody responses to West Nile virus. *Clin Vaccine Immunol* 13 (2006) 33-36.
41. Levett PN, Sonnenberg* K, Shead S, Sidaway F, Niedrig M, Kreher* P, Arp* B, Erdmann* S, Jahn A, Steinhagen* K, Drebot M, Horsman GB (*EUROIMMUN AG). IgG avidity determinations but not IgM tests are capable to differentiate recent from previous infections with West Nile Virus. Poster zum 21. Annual clinical virology symposium and annual meeting Pan American society for clinical virology. USA, Florida, Clearwater beach (2005).
42. Stiasny K, Kiermayr S, Heinz FX. Cryptic properties of a cluster of dominant flavivirus cross-reactive antigenic sites. *J Virol* 80 (2006) 9557-9566.
43. Stöcker* W, Horbach E, Stenczka W, Schaper J, Böttger K, Herrmann C (*EUROIMMUN AG). Automatische Bestückung von Objektträgern mit „BIOCHIPs“: Eine neue Technologie für die indirekte Immunfluoreszenz. *Lab. med.* (1989) 171-172.
44. Sonnenberg* K, Hinrichs E, Müller-Kunert* E, Schlumberger* W, Stöcker* W (*EUROIMMUN AG). Einsatz der BIOCHIP-Technologie in der Serologie von Influenza-, EBV- und Legionella-Infektionen. Poster zum 3. Deutschen Kongress für Infektions- und Tropenmedizin, 15. bis 18. März 1995 in Berlin-Hohenschönhausen. Abstrakt im Kongressband (1995).
45. EUROIMMUN AG. Stöcker W. Die BIOCHIP-Technologie: Ein neuer Weg der Produktion Festphase-gebundener Bioreagenzien für das medizinische und biologische Laboratorium. Dokumentation zum Vorschlag des Rektors der Medizinischen Universität Lübeck für den Schmidt-Römhild-Technologiepreis 1989. Eigenverlag der Firma EUROIMMUN (1990).



Declaration of Conformity

EUROIMMUN Medizinische Labordiagnostika AG
Seekamp 31, D-23560 Lübeck, Germany

declare under our sole responsibility that the ELISA products

Anti-West Nile Virus ELISA (IgG)
Anti-West Nile Virus ELISA (IgM)
Avidity: Anti-West Nile Virus ELISA (IgG)
EI 2662-9601 G
EI 2662-9601 M
EI 2662-9601-1 G

(product name, order no)

meet the demands of

Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998 and its transpositions in national laws which apply to it.
Conformity assessment procedure: Annex III

Lübeck, 13.10.2016
(Place and date of issue)

Susanne Aleksandrowicz
- Member of the Board -

Dr. Wolfgang Schlumberger
- Member of the Board -



Avidity determination of antibodies against Measles Viruses (IgG)
Test instruction


ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2610-9601-1 G	Measles viruses	IgG	Ag-coated microplate wells	96 x 01 (96)

Background

The differentiation between fresh and long-standing infections is one of the greatest challenges in serology. Until now this was based mainly on determination of specific antibodies of the immunoglobulin class IgM, which generally only appear initially. However, the detection of these antibodies is often unreliable and problematic due to interfering factors such as persistence of the IgM response, too weak or delayed IgM production, and unspecific IgM production through polyclonal B-cell stimulation.

In recent years additional determination of the antibody avidity has become an established method for identification of primary infections. The immune system reacts to an infection by first forming low-avidity antibodies. With continued disease duration, IgG that are more precisely adapted to the antigens are produced – the avidity increases. If high-avidity IgG are detectable in the serum, it can be assumed that the infection is at a late stage.

Contents of the test system: EI 2610-9601-1 G:

Component	Colour	Format	Symbol
1. Test kit Anti-Measles Viruses ELISA (IgG, order number EI 2610-9601 G)	---	---	---
2. Positive control HA	red	1 x 1,3 ml	POS CONTROL HA
3. Positive control LA	blue	1 x 1,3 ml	POS CONTROL LA
4. Urea solution for Anti-Measles ELISA, ready for use	yellow	1 x 12 ml	UREA
5. Phosphate buffer ready for use	light blue	1 x 12 ml	PBS BUFFER
6. Test instruction	---	1 booklet	---
LOT In-vitro determination	 Storage temperature Unopened usable until		

Storage and stability: The test kit has to be stored at a temperature between 2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents are to be disposed of according to official regulations.

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18 C to +25 C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2 C to +8 C and protected from contamination, unless stated otherwise below.

- Controls: Ready for use. The reagents must be mixed thoroughly before use.
- Urea solution: Ready for use. The urea solution included in this test system may only be used for the avidity determination of antibodies against Measles.
- Phosphate buffer: Ready for use.

Warning: The calibrators and controls used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.

Preparation and stability of the patient samples

Sample material: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8 C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted 1:101 sample buffer. For example: dilute 10 µl serum to 1,0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The controls are prediluted and ready for use, do not dilute them.



Incubation

Sample incubation: (1. step)

Transfer 100 µl of the controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.
Incubate for 30 minutes at room temperature (+18°C to +25°C).

Wash:

Manual: Empty the wells and subsequently wash 1 time using 300 µl of working strength wash buffer.
Automatic: Wash reagent wells 1 time with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer. Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Urea incubation: (2. step)

Pipette 200 µl of urea solution into each of the microplate wells of the first microtiter strip and 200 µl of phosphate buffer into each of the microplate wells of the second microtiter strip.
Incubate for 10 minutes at room temperature (+18°C to +25°C).

Wash:

Empty the wells. Wash as described above, but wash 3 times using working strength wash buffer for each wash.

Attention: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values.
Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Conjugate incubation: (3. step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells.
Incubate for 30 minutes at room temperature.

Wash:

Empty the wells. Wash as described above, but wash 3 times using working strength wash buffer for each wash.

Substrate incubation: (4. step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells.
Incubate for 15 minutes at room temperature (+18°C to +25°C) protect from direct sunlight.

Stopping the reaction:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm within 30 minutes of adding the stop solution. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



Pipetting protocol

1	2	3	4	5	6	7	8	9	10	11	12
A pos	pos	P 7	P 7	P 16	P 15						
HA	HA										
G pos	pos	P 8	P 8	P 16	P 16						
LA	LA										
C P 1	P 1	P 9	P 9	P 17	P 17						
D P 2	P 2	P 10	P 10	P 18	P 18						
E P 3	P 3	P 11	P 11								
F P 4	P 4	P 12	P 12								
G P 5	P 5	P 13	P 13								
H P 6	P 6	P 14	P 14								

The above pipetting protocol is an example of the avidity determination of IgG antibodies in 18 patient samples (P 1 to P 18).

Controls (pos HA and pos LA) as well as the patient samples have been incubated in duplicate in one well each of two different microtiter strips. The reagent wells of the microtiter strips 1, 3, 5 etc. are treated with urea solution after the incubation with patient samples, the reagent wells of the microtiter strips 2, 4, 6 etc. are treated with phosphate buffer.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimises reagent wastage.

Both positive controls with high-avidity and low-avidity antibodies serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

The presence of low-avidity antibodies in a patient's serum has been proved if the ELISA extinction value is considerably reduced by urea treatment. For an objective interpretation the relative avidity index (RAI) is calculated and expressed in percent using the extinction values with and without urea treatment.

$$\text{Extinction of the sample with urea treatment} \times 100 = \text{relative avidity index (RAI) in \%}$$

$$\text{Extinction of the sample without urea treatment}$$

The upper limit of the range of low-avidity antibodies (cut-off value) recommended by EUROIMMUN is 40% RAI. Values below the indicated cut-off are to be considered as an indication of low-avidity antibodies, values between 40% and 60% RAI as equivocal, values above 60% RAI as an indication of high-avidity antibodies. If a result is classified as equivocal, it is recommended to collect a second sample not less than 7 days later and to test it together with the first sample.

RAI < 40%: indication of low-avidity antibodies
RAI 40% - 60%: equivocal
RAI > 60%: indication of high-avidity antibodies

Reliable results in the the measurement of IgG antibody avidity can only be yielded if the patient sample contains a diagnostically significant concentration of specific antibodies.

Generally, the determination of the relative avidity index is not helpful in samples which have an O.D. of < 0.140 after incubation without urea treatment.



For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

Attention:

In some patients with an acute infection, very high titers of IgG antibodies can be found. Even though the specific IgG antibody population is in different maturation stages, both, high-avidity and low-avidity, the vacant antigen epitopes are predominantly occupied by high-avidity antibodies in high titer samples. The determination of the avidity of the whole specific IgG antibody population can lead to false high RAI values in results.



False high RAI values were found in some cases of acute infections when the extinction value of the IgG measurement without urea treatment was >1,200.

It is recommended for samples with extinction values of >1,200 to repeat the avidity determination with a higher sample dilution (e.g. 1:401). If low avidity of IgG antibodies is already found at extinction values of >1,200, no further testing is necessary.

Test characteristics

- A panel of 104 sera from patients with the following diseases was investigated:
- 10 patients with clinically and serologically diagnosed acute measles infection
- 14 patients with acute EBV infection and resulting polyclonal stimulation
- 29 healthy blood donors
- 44 patients with MMR vaccination more than 6 months ago
- 7 patients with SSPE after measles infection

In these panels, patients with acute infection showed an average RAI of 13%, whereas the RAI determined for the remaining patients without acute measles infection was > 80%.



Clinical significance

The measles virus (MV) is the most instantly recognisable member of Morbilliviruses, a group of viruses belonging to the Paramyxoviridae family [1]. No animal reservoir is known. The measles virus causes an acute febrile illness which occurs mainly in childhood and is very infectious [2, 3]. In 1999, measles still caused worldwide 873,000 deaths per year [1, 4, 5]. Today they are less frequent, because of vaccination, especially in the western hemisphere [6, 7]. However, measles epidemics are still observed in some countries [2, 3, 4, 5, 6, 8, 9]. Individuals acutely infected with the virus exhibit a wide range of clinical symptoms ranging from a characteristic mild self-limiting infection to death [1, 2, 8, 10].

MV infections are characterised by an incubation period of about 10 days, flu-like symptoms with fever, malaise, catarrh of the upper respiratory tract, cough, congestion and conjunctivitis. Soon afterwards the measles rash, a typical exanthema, appears first near the ears, then on the forehead, in the face and over the rest of the body [1, 5, 8, 11].

Complications arising from MV infections include secondary bacterial pneumonia, otitis media (approx. 1%), encephalitis (approx. 1%), myocarditis, miscarriage and a condition called subacute sclerosing pan-encephalitis (SSPE) [5, 12, 13]. Persistent MV infection of the otic capsule is an aetiological factor in otosclerosis [9, 14]. Anti-measles IgG for the serological diagnosis of otosclerotic hearing loss has a high specificity and sensitivity [15]. SSPE is a progressive, generally fatal brain disorder caused by chronic measles virus infection. It occurs about 7 to 10 years after the infection and generally kills within 3 years from the onset of the symptoms. The patients suffer from behavioural changes, cognitive deterioration, vision problems and eventually advanced neurological symptoms, such as severe spasms, and finally severe physical and mental impairment that leads to death [13]. Males are more commonly affected than females.

The risk of SSPE from measles was underestimated according to older data [12]. Actual papers put it at closer to 6.5 to 11 cases of SSPE per 100,000 measles infections; that means 7 to 13 times higher than the earlier estimates [1, 9, 12].



Women with acute measles infection during pregnancy and a negative result for measles-specific antibodies were observed e.g. in Japan, India, Thailand, Kenya and Brazil: 3 of 4 pregnancies ended in preterm delivery, spontaneous abortion or stillbirth; 2 of 4 neonates were found to have congenital measles with a positive result for IgM antibodies [5].

Antibodies against MV can be found in the serum of almost all patients during and after a measles infection. IgM antibodies develop soon after the onset of symptoms and can be measured using ELISA or indirect immunofluorescence tests [16, 20, 21, 22]. 50% of patients have IgM antibodies within three days, more than 90% within 10 days after occurrence of the rash [15, 17]. The Anti-Measles Virus IgM ELISA is more rapid and sensitive for the serological diagnosis of measles infections than other tests [15, 18, 19]. MV infections often cause an increase in heterologic antibodies. The statistically evaluated detection rate for antibodies is significantly higher for ELISA and IFT in comparison with e.g. neutralisation tests [16, 20]. IgG and IgM antibodies against MV are reliable markers to confirm suspected measles infections.

Measles myelitis or encephalitis can be verified by detecting antibodies against measles in the cerebrospinal fluid (CSF) [23, 24, 25, 26]. These specific antibodies are synthesised in the brain [24]. The CSF-serum quotient (L:SQ) allows to differentiate between a blood-derived and a pathological, brain-derived specific antibody fraction in CSF, taking into account individual changes in the blood/CSF barrier function [24, 25, 26, 27]. Therefore it is necessary to confirm the presence of antibodies against MV using ELISA both in CSF and in the serum. During measles myelitis or encephalitis an intrathecal synthesis of antibodies against MV in CSF takes place. Due to the fact that specific antibodies can pass the blood-cerebrospinal fluid barrier by diffusion from serum to CSF it is necessary to determine the relative CSF/serum quotient (CSQ_{rel}, synonym: antibody specificity index) [24, 25, 26]. The quotient is calculated from the amount of specific anti-measles virus IgG antibodies in total CSF IgG in proportion to the amount of specific IgG antibodies in total serum IgG. During conversion the CSF/serum quotient of the pathogen-specific IgG-antibody concentrations CSQ_{path-spec} (IgG) is put into relation to the CSF/serum quotient of the total IgG concentrations CSQ_{total} (IgG) [27]. A relative CSQ result above 1.5 indicates the production of specific antibodies in the central nervous system (CNS) and the involvement of the CNS in the disease [25, 26].

With respect to the severe complications known from measles infections, the Robert Koch Institute in Germany recommends vaccinating small children, with a first shot between the age of 11 to 14 months and a second between 15 and 23 months [2, 4, 10, 11]. Neutralisation activity and persistence of antibodies are induced in response to the immunisation [6, 15].

Life-long immunity is generally developed. However, antibody levels are 8 to 10 times lower in post-vaccination sera than in convalescent sera [6, 19, 28]. A passive immunisation with specific immunoglobulin concentrates is usually given to immunosuppressed seronegative individuals, such as tumour patients and recipients of transplants, as well as to seronegative pregnant women after exposure to the virus.

The European Regional Office of the WHO aims at eliminating measles from the region in the following years by area-wide vaccination campaigns [4, 9, 10]. This is expected to limit the number of apparent infections and especially of severe courses of the disease. For the diagnosis of the remaining cases of measles infection and of infections acquired outside Europe as well as for the clarification of atypical courses of the disease in partly immunised patients the antibody determination in serum and CSF will be of growing importance [3, 5, 8, 9].



Literature references

1. Rima BK, Duprex WP. Morbilliviruses and human disease. *J Pathol* 208 (2006) 199-214.
2. Hogg GG, Darlington RJ, Hogg KG, Lester R. Measles immunity and immunisation status in Australian children 1 to 4 years of age. *J Paediatr Child Health* 42 (2006) 165-169.
3. Karimi A, Arjomandi A, Alborzi A, Rasouli M, Kadivar MR, Obood B, Pourabbas B. Prevalence of measles antibody in children of different ages in Shiraz, Islamic Republic of Iran. *East Mediterr Health J* 10 (2004) 469-473.
4. Cutts FT, Steinglass R. Should measles be eradicated. *Br Med J* 316 (1998) 765-767.
5. de Barros EN, Silva EM. Epidemiologic surveillance of measles and rubella in Campinas (SP), Brazil: the reliability of the data. [Article in Portuguese]. *Rev Panam Salud Publica* 19 (2006) 172-178.
6. Gay N, Ramsay M, Cohen B, Hesketh L, Morgan-Capner P, Brown D, Miller E. The epidemiology of measles in England and Wales since the 1994 vaccination campaign. *Commun Dis Rep CDR Rev* 7 (1997) 17-21.
7. de Melker H, Pebody RG, Edmunds WJ, Levy-Bruhl D, Valle M, Rota MC, Salmazo S, van den Hof S, Berbers G, Salou P, Spaendonck MCV, Crovari P, Davidkin I, Gabutti G, Hesketh L, Morgan-Capner P, Plesner AM, Raux M, Tische A, Miller E. The seroepidemiology of measles in Western Europe. *Epidemiol Infect* 126 (2001) 249-259.
8. Ceylan A, Etem M, Korukluoglu G, Acemoglu H, Kara IH, Eten PG, Arslan C, Ay ME. An epidemic caused by measles virus type D6 in Turkey. *Turk J Pediatr* 47(2005) 309-315.
9. Ota MO, Moss WJ, Griffin DE. Emerging diseases: measles. *J Neurovirol* 11 (2005) 447-454.
10. Shann F. A little bit of measles does you good. *Br Med J* 219 (1999) 4-5.
11. CDC. Measles outbreaks still occur among school-age children and travelers. *MMWR* 46 (1997) 242-245.
12. Smith TC. Measles vaccine doesn't cause SSPE. *Aetiology* (2006).
13. Yokoyama T, Sakurai M, Aota Y, Wakabayashi Y, Ohyanishi K. An adult case of acute disseminated encephalomyelitis accompanied with measles infection. *Intern Med* 44 (2005) 1204-1205.
14. Singh MP, Ratho RK, Panda N, Mishra B. Otsclerosis – do we have a viral aetiology? *Nepal Med Coll J* 7 (2005) 129-130.
15. Tische A, Gerke E, Strauss J, Smelhausova M, Mirazova M. Immunoglobulin specific and conventional methods in the serodiagnosis of measles. [Article in German]. *Z Gesamte Hyg* 35 (1989) 367-369.
16. Roodhri F, Roustai MH, Mostafaei A, Salemanjandeh H, Foroshani RS, Sabahi F. Development of an enzyme-linked immunosorbent assay for immunoglobulin M antibodies against measles virus. *Clin Diagn Lab Immunol* 10 (2003) 439-442.
17. Rosser E, Miller H, McCulloch B, Sullivan L, Ward K. Comparison of immunofluorescence and enzyme immunoassay for detection of measles-specific immunoglobulin M antibody. *J Clin Microbiol* 29 (1991) 1069-1071.
18. Pedersen JR, Antonsdottir A, Ewald T, Mordhorst CH. Detection of measles IgM antibodies by enzyme-linked immunosorbent assay (ELISA). *Acta Patrol Microbiol Immunol Scand [B]* 90 (1982) 153-160.
19. Glikmann G, Petersen I, Mordhorst CH. Prevalence of IgG-antibodies to mumps and measles virus in non-vaccinated children. *Dan Med Bull* 35 (1988) 185-187.



20. Bouche FB, Brons NH, Houard S, Schneider F, Müller CP. Evaluation of hemagglutinin protein-specific immunoglobulin-M-for-diagnosis-of-measles-by-an-enzyme-linked-immunosorbent-assay based on recombinant protein produced in a high-efficiency mammalian expression system. *J Clin Microbiol* 36 (1998) 3509-3513.
21. EUROIMMUN AG. Stöcker W, Fauer H, Krause C, Barth E, Martinez A. Verfahren zur Optimierung der automatischen Fluoreszenzzerknennung in der Immun Diagnostik. Deutsche Patentanmeldung (Offenlegungsschrift) DE 10 2006 027 516 0 und WO2007/140952 (2006).
22. EUROIMMUN AG. Aktuelle Themen der Autoimmundiagnostik und der Infektions-Serologie. Wissenschaftliches Fortbildungsseminar mit Vorträgen von Prof. Dr. G. Wick, Prof. Dr. N. Sepp, Prof. Dr. F. Delsenhammer, Dr. med. W. Stöcker, Prof. Dr. Gerold Stanek, Prof. Dr. R. Würzner, A. Krapf, Dr. med. Dipl. oec. med. J. Brunner, Innsbruck (2007).
23. Dennis RH, Herb E. Immunological diagnosis in viral infections of the central nervous system: course of antibody titres against homo- and heterologous viruses. *Med Microbiol Immunol* 178 (1989) 255-268.
24. Reiber HQ, Lange P. Quantification of Virus-Specific Antibodies in Cerebrospinal Fluid and Serum: Sensitive and Specific Detection of Antibody Synthesis in Brain. *Clin Chem* 37 (1991) 1153-1160.
25. Reiber H, Peter JB. Cerebrospinal fluid analysis: disease-related data patterns and evaluation programs. *J Neurol Sci* 184 (2001) 101-122.
26. Reiber H, Lange P. Virus-spezifische Antikörper in Liquor und Serum. ELISA-Analytik und Auswertung mittels Antikörper-Index und Quotientendiagramm. *Lab Med* 15 (1991) 204-207.
27. Reiber H, Ungelohr S, Jacobi C. The intrathecal, polyspecific and oligoclonal immune response in multiple sclerosis. *Multiple Sclerosis* 4 (1998) 111-117.
28. Dine MS, Hutchins SS, Thomas A, Williams I, Bellini WJ, Redd SC. Persistence of vaccine-induced antibody to measles 26-33 years after vaccination. *J Infect Dis* 189 (2004) 123-130.





Anti-Measles Virus ELISA (IgG)

Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI2610-9601 G	Measles virus	IgG	Ag-coated microplate wells	96 x 01 (96)

Indication: measles

Principles of the test: The ELISA test kit provides a quantitative or semiquantitative in vitro assay for human antibodies of the IgG class against measles virus in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with measles virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	[STRIPS]
2. Calibrator 1 5000 IU/l (IgG, human), ready for use		1 x 2.0 ml	[CAL 1]
3. Calibrator 2 1000 IU/l (IgG, human), ready for use	red coloured in decreasing intensity	1 x 2.0 ml	[CAL 2]
4. Calibrator 3 250 IU/l (IgG, human), ready for use		1 x 2.0 ml	[CAL 3]
5. Calibrator 4 50 IU/l (IgG, human), ready for use		1 x 2.0 ml	[CAL 4]
6. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	[POS-CONTROL]
7. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	[NEG-CONTROL]
8. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	[CONJUGATE]
9. Sample buffer ready for use	light blue	1 x 100 ml	[SAMPLE BUFFER]
10. Wash buffer 10x concentrate	colourless	1 x 100 ml	[WASH BUFFER 10X]
11. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	[SUBSTRATE]
12. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	[STOP SOLUTION]
13. Test instruction		1 booklet	
14. Quality control certificate		1 protocol	
LOT Lot			Storage temperature
VD In vitro determination			<input type="checkbox"/> Unopened, usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C), approx. 30 minutes, before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

Coated wells: Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.

Calibrators and controls: Ready for use. The reagents must be mixed thoroughly before use.

Enzyme conjugate: Ready for use. The enzyme conjugate must be mixed thoroughly before use.

Sample buffer: Ready for use.

Wash buffer: The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).

For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.

Chromogen/substrate solution: Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The Chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.

Stop solution: Ready for use.

Warning: The controls and calibrators used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted 1:101 in sample buffer. For example: dilute 10 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.



Incubation

For semiquantitative analysis incubate calibrator 3 along with the positive and negative controls and patient samples. For quantitative analysis incubate calibrators 1 to 4 along with the positive and negative controls and patient samples.
(Partly) manual test performance

Sample incubation:
(1st step)

Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for 30 minutes at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles), too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labeled anti-human IgG) into each of the microplate wells. Incubate for 30 minutes at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation:
(3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for 15 minutes at room temperature (+18°C to +25°C) (protect from direct sunlight).

Stopping the reaction:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm within 30 minutes of adding the stop solution. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer 1 or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry.

Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.



Pipetting protocol

1	2	3	4	5	6	7	8	9	10	11	12
A	C 3	P 6	P 14	P 22		C 1	E 3	P 11	P 18		
B	pos	P 7	P 15	P 23		C 2	P 4	P 12	P 20		
C	neg	P 8	P 16	P 24		C 3	P 5	P 13	P 21		
D	P 1	P 9	P 17			C 4	P 6	P 14	P 22		
E	P 2	P 10	P 18			pos	P 7	P 15	P 23		
F	P 3	P 11	P 19			neg	P 8	P 16	P 24		
G	P 4	P 12	P 20				P 1	P 9	P 17		
H	P 5	P 13	P 21				P 2	P 10	P 18		

The pipetting protocol for microtiter strips 1-4 is an example for the semiquantitative analysis of 24 patient samples (P 1 to P 24).
The pipetting protocol for microtiter strips 7-10 is an example for the quantitative analysis of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 4), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The reagent wells are break-off format. Therefore, the number of tests performed can be matched to the number of samples, minimizing reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 3. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 3}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

$$\begin{array}{ll} \text{Ratio} < 0,8: & \text{negative} \\ \text{Ratio} \geq 0,8 \text{ to } < 1,1: & \text{borderline} \\ \text{Ratio} \geq 1,1: & \text{positive} \end{array}$$

In cases of borderline test results, an additional patient sample should be taken 7 days later and retested in parallel with the first patient sample. The results of both samples allow proper evaluation of titer changes.

Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 4 calibrators against the corresponding units (linear/non-linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction of a serum sample lies above the value of calibrator 1 (5000 IU/l), the result should be given as ">5000 IU/l". It is recommended that the sample be re-tested at a dilution of 1:400. The result in IU/l read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the reference range of non-infected persons (cut-off value) recommended by EUROIMMUN is 250 International Units (IU). EUROIMMUN recommends interpreting results as follows:

<200 IU/l:	negative
≥200 to <275 IU/l:	borderline
≥275 IU/l:	positive

Evaluation information: For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class IgG) can be helpful. Diagnosis can be secured by the determination of the titer change in two serum samples taken at an interval of at least 7 days and analysed in parallel.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

Test characteristics

Calibration: The controls of the Anti-Measles Virus ELISA (IgG) were calibrated using the 3rd international standard serum NIBSC 97/648 (anti-measles and anti-polio virus serum, National Institute for Biological Standards and Control, Hertfordshire, England, approved as international reference preparation by the WHO Expert Committee on Biological Standardization). The NIBSC 97/648 serum contains 3 International Units (IU) per ampoule by definition and was resuspended in a concentration of 3 IU/ml.

For every group of tests performed, the extinction values of the calibrators and the international units determined for the positive and negative controls must be within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.



The activity of the enzyme used is temperature-dependent and the extinction values may vary, if a thermostat is not used. The higher the room temperature during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is provided by inactivated cell lysates of Vero cells infected with the "Edmonston" strain of measles viruses.

Linearity: The linearity of the anti-measles viruses ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95. The Anti-Measles Virus ELISA (IgG) is linear at least in the tested concentration range (52 IU/l - 4865 IU/l).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-Measles Virus ELISA (IgG) is 8 IU/l.

Cross reactivity: The quality of the antigen used ensures a high specificity of the ELISA. Sera from patients with infections caused by various agents were investigated with the EUROIMMUN Anti-measles virus ELISA (IgG)

Antibodies against	n	Anti-Measles virus ELISA (IgG)
Adenovirus	8	0%
CMV	6	0%
EBV-CA	11	0%
HSV-1	3	0%
Influenza virus type A	5	0%
Influenza virus type B	11	0%
Mumps virus	4	0%
Mycoplasma pneumoniae	4	0%
Parainfluenza virus types 1-4	11	0%
RSV	9	0%
Rubella virus	6	0%
Toxoplasma	3	0%
VZV	5	0%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for hemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

Serum	Intra-assay variation, n = 20	
	Mean value (IU/l)	CV (%)
1	830	8.0
2	3410	6.6
3	3725	5.6

Serum	Inter-assay variation, n = 4 x 6	
	Mean value (IU/l)	CV (%)
1	796	11.6
2	3635	5.0
3	3946	6.8



Specificity and sensitivity: 112 clinically characterized patient samples (interlaboratory test samples from INSTANB, Germany; Labqualty, Finland and NEOAS, UK) were examined with the EUROIMMUN Anti-Measles Virus ELISA (IgG). The test showed a specificity and a sensitivity of 100% each.

	INSTANB / Labqualty / NEOAS (IgG)		
	positive	borderline	negative
EUROIMMUN	89	1	0
Anti-Measles Virus ELISA (IgG)	0	0	0
	0	0	22

n = 112

Reference range: The levels of anti-measles virus antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 250 IU/l, 94% of the blood donors were anti-measles virus positive (IgG), which reflects the known percentage of infections in adults.

Clinical significance

The measles virus (MV) is the most instantly recognizable member of Morbilliviruses, a group of viruses belonging to the Paramyxoviridae family [1]. No animal reservoir is known. The measles virus causes an acute feverish illness which occurs mainly in childhood and is very infectious [2, 3]. In 1999, measles still caused worldwide 873,000 deaths per year [1, 4, 5]. Today they are less frequent because of vaccination, especially in the western hemisphere [6, 7]. However, measles epidemics are still observed in some countries [2, 3, 4, 5, 8, 9]. Individuals acutely infected with the virus exhibit a wide range of clinical symptoms ranging from a characteristic mild self-limiting infection to death [1, 2, 8, 10].

MV infections are characterised by an incubation period of about 10 days. Flu-like symptoms with fever, malaise, catarrh of the upper respiratory tract, cough, congestion and conjunctivitis. Soon afterwards the measles rash, a typical exanthema, appears first near the ears, then on the forehead, in the face and over the rest of the body [1, 5, 8, 11].

Complications arising from MV infections include secondary bacterial pneumonia, otitis media (approx. 1%), encephalitis (approx. 1%), myocarditis, miscarriage and a condition called subacute sclerosing panencephalitis (SSPE) [5, 12, 13]. Persistent MV infection of the oitic capsule is an aetiological factor in otosclerosis [9, 14]. Anti-measles IgG for the serological diagnosis of otosclerotic hearing loss has a high specificity and sensitivity [15]. SSPE is a progressive, generally fatal brain disorder caused by chronic measles virus infection. It occurs about 7 to 10 years after the infection and generally kills within 3 years from the onset of the symptoms. The patients suffer from behavioural changes, cognitive deterioration, vision problems and eventually advanced neurological symptoms, such as severe spasms, and finally severe physical and mental impairment that leads to death [13]. Males are more commonly affected than females. The risk of SSPE from measles was underestimated according to older data [12]. Actual papers put it at closer to 6.5 to 11 cases of SSPE per 100,000 measles infections; that means 7 to 13 times higher than the earlier estimates [1, 9, 12].

Women with acute measles infection during pregnancy and a negative result for measles-specific antibodies were observed e.g. in Japan, India, Thailand, Kenya and Brazil: 3 of 4 pregnancies ended in preterm delivery, spontaneous abortion or stillbirth; 2 of 4 neonates were found to have congenital measles with a positive result for IgM antibodies [5].

Antibodies against MV can be found in the serum of almost all patients during and after a measles infection. IgM antibodies develop soon after the onset of symptoms and can be measured using ELISA or indirect immunofluorescence tests (IIFT) [16, 20, 21, 22]. 50% of patients have IgM antibodies within three days, more than 90% within 10 days after occurrence of the rash [15, 17]. The Anti-Measles Virus IgM ELISA is more rapid and sensitive for the serological diagnosis of measles infections than other tests [15, 18, 19].



MV infections often cause an increase in heterologous antibodies. The statistically evaluated detection rate for antibodies is significantly higher for ELISA and IIFT in comparison with e.g. neutralisation tests [16, 20]. IgG and IgM antibodies against MV are reliable markers to confirm suspected measles infections.

Measles myelitis or encephalitis can be verified by detecting antibodies against measles in the cerebrospinal fluid (CSF) [23, 24, 25, 26]. These specific antibodies are synthesised in the brain [24]. The CSF-serum quotient (LSQ) allows to differentiate between a blood-derived and a pathological, brain-derived specific antibody fraction in CSF, taking into account individual changes in the blood/CSF barrier function [24, 25, 26, 27]. Therefore it is necessary to confirm the presence of antibodies against MV using ELISA both in CSF and in the serum. During measles myelitis or encephalitis an irritable synthesis of antibodies against MV in CSF takes place. Due to the fact that specific antibodies can pass the blood-cerebrospinal fluid barrier by diffusion from serum to CSF it is necessary to determine the relative CSF/serum quotient (CSQrel., synonym: antibody specificity index) [24, 25, 26]. The quotient is calculated from the amount of specific anti-measles virus IgG antibodies in total CSF IgG in proportion to the amount of specific IgG antibodies in total serum IgG. During conversion the CSF/serum quotient of the pathogen-specific IgG-antibody concentrations CSQpath-spec. (IgG) is put into relation to the CSF/serum quotient of the total IgG concentrations CSQtotal (IgG) [27]. A relative CSQ result above 1.5 indicates the production of specific antibodies in the central nervous system (CNS) and the involvement of the CNS in the disease [25, 26].

With respect to the severe complications known from measles infections, the Robert Koch Institute in Germany recommends vaccinating small children, with a first shot between the age of 11 to 14 months and a second between 15 and 23 months [2, 4, 10, 11]. Neutralisation activity and persistence of antibodies are induced in response to the immunisation [6, 15].

Life-long immunity is generally developed. However, antibody levels are 8 to 10 times lower in post-vaccination sera than in convalescent sera [6, 19, 28]. A passive immunisation with specific immunoglobulin concentrates is usually given to immunosuppressed seronegative individuals, such as tumour patients and recipients of transplants, as well as to seronegative pregnant women after exposure to the virus.

The European Regional Office of the WHO aims at eliminating measles from the region in the following years by area-wide vaccination campaigns [4, 9, 10]. This is expected to limit the number of apparent infections and especially of severe courses of the disease. For the diagnosis of the remaining cases of measles infection and of infections acquired outside Europe as well as for the clarification of atypical courses of the disease in partly immunised patients the antibody determination in serum and CSF will be of growing importance [3, 5, 8, 9].



Literature references

1. Rima BK, Duprex WP. Morbilliviruses and human disease. *J Pathol* 208 (2006) 199-214.
2. Hogg GG, Darlington RJ, Hogg KG, Lester R. Measles immunity and immunisation status in Australian children 1 to 4 years of age. *J Paediatr Child Health* 42 (2006) 165-169.
3. Karimi A, Ajomandi A, Alborzi A, Rasouli M, Kadiyar MR, Obood B, Pourabbas B. Prevalence of measles antibody in children of different ages in Shiraz, Islamic Republic of Iran. *East Mediterr Health J* 10 (2004) 468-473.
4. Cutts FT, Steinglass R. Should measles be eradicated. *Br Med J* 316 (1998) 765-767.
5. de Barros EN, Silva EM. Epidemiologic surveillance of measles and rubella in Campinas (SP), Brazil: the reliability of the data. [Article in Portuguese] *Rev Panam Salud Publica* 19 (2006) 172-178.
6. Gay N, Ramsay M, Cohen B, Hesketh L, Morgan-Capner P, Brown D, Miller E. The epidemiology of measles in England and Wales since the 1994 vaccination campaign. *Commun Dis Rep CDR Rev* 7 (1997) 17-21.
7. de Melker H, Pebody RG, Edmunds WJ, Levy-Bruhl D, Valle M, Rota MC, Salmasso S, van den Hof S, Berbers G, Salou P, Speendonck MCV, Crovari P, Davidkin I, Gabutti G, Hesketh L, Morgan-Capner P, Plesner AM, Raux M, Tische A, Miller E. The seroepidemiology of measles in Western Europe. *Epidemiol Infect* 126 (2001) 249-259.
8. Ceylan A, Ertem M, Kouklouglu G, Acemoglu H, Kara IH, Ertan PG, Arslan C, Ay ME. An epidemic caused by measles virus type D6 in Turkey. *Turk J Pediatr* 47(2005) 309-315.
9. Ota MO, Moss WJ, Griffin DE. Emerging diseases: measles. *J Neurovrol* 11 (2005) 447-454.
10. Shann F. A little bit of measles does you good. *Br Med J* 219 (1999) 4-5.
11. CDC. Measles outbreaks still occur among school-age children and travelers. *MMWR* 46 (1997) 242-245.
12. Smith TC. Measles vaccine doesn't cause SSPE. *Aetiology* (2006).
13. Yokoyama T, Sakurai M, Aota Y, Wakabayashi Y, Ohyashiki K. An adult case of acute disseminated encephalomyelitis accompanied with measles infection. *Intern Med* 44 (2005) 1204-1205.
14. Singh MP, Ratho RK, Panda N, Mishra B. Orosclerosis - do we have a viral aetiology? *Nepal Med Coll J* 7 (2005) 129-130.
15. Tische A, Gerke E, Strauss J, Smelhausova M, Mrazova M. Immunoglobulin specific and conventional methods in the serodiagnosis of measles. [Article in German] *Z Gesamte Hyg* 35 (1989) 367-369.
16. Roodhart F, Rousstal MH, Mostatafe A, Soleimanjahi H, Foroshani RS, Sabahi F. Development of an enzyme-linked immunosorbent assay for immunoglobulin M antibodies against measles virus. *Clin Diagn Lab Immunol* 10 (2003) 439-442.
17. Rossier E, Miller H, McCulloch B, Sullivan L, Ward K. Comparison of immunofluorescence and enzyme immunoassay for detection of measles-specific immunoglobulin M antibody. *J Clin Microbiol* 29 (1991) 1059-1071.
18. Pedersen JR, Antonsdottir A, Ewald T, Mordhorst CH. Detection of measles IgM antibodies by enzyme-linked immunosorbent assay (ELISA). *Acta Pathol Microbiol Immunol Scand [B]* 90 (1982) 153-160.
19. Gilkman G, Petersen I, Mordhorst CH. Prevalence of IgG-antibodies to mumps and measles virus in non-vaccinated children. *Dan Med Bull* 35 (1988) 185-187.
20. Bouche FB, Bruns NH, Howard S, Schneider F, Müller CP. Evaluation of hemagglutinin protein-specific immunoglobulin M for diagnosis of measles by an enzyme-linked immunosorbent assay based on recombinant protein produced in a high-efficiency mammalian expression system. *J Clin Microbiol* 36 (1998) 3509-3513.
21. Stöcker* W, Fauer* H, Krause* C, Barth E, Martinez A (*EUROIMMUN AG). Verfahren zur Optimierung der automatischen Fluoreszenzkerennung in der Immun Diagnostik. Deutsche Patentanmeldung (Offenlegungsschrift) DE 10 2006 027 516.0 und WO2007/140952 (2006).
22. EUROIMMUN AG. Aktuelle Themen der Autoimmun Diagnostik und der Infektions-Serologie. Wissenschaftliches Fortbildungsseminar mit Vorträgen von Prof. Dr. G. Wick, Prof. Dr. N. Sepp, Prof. Dr. F. Dörsenhammer, Dr. med. W. Stöcker, Prof. Dr. Gerold Starek, Prof. Dr. R. Würzner, A. Kraft, Dr. med. Dipl. oec. med. J. Brunner, Innsbruck (2007).
23. Dennis RH, Herb E. Immunological diagnosis in viral infections of the central nervous system: course of antibody titres against homo- and heterologous viruses. *Med Microbiol Immunol* 178 (1989) 255-268.
24. Reiber HQ, Lange P. Quantification of Virus-Specific Antibodies in Cerebrospinal Fluid and Serum: Sensitive and Specific Detection of Antibody Synthesis in Brain. *Clin Chem* 37 (1991) 1153-1160.
25. Reiber H, Peter JB. Cerebrospinal fluid analysis: disease-related data patterns and evaluation programs. *J Neurol Sci* 184 (2001) 101-122.
26. Reiber H, Lange P. Virus-spezifische Antikörper in Liquor und Serum. ELISA-Analytik und Auswertung mittels Antikörper-Index und Quotientendiagramm. *Lab Med* 15 (1991) 204-207.
27. Reiber H, Ungefehr S, Jacobi C. The intrathecal, polyspecific and oligoclonal immune response in multiple sclerosis. *Multiple Sclerosis* 4 (1998) 1111-1117.
28. Dine MS, Hutchins SS, Thomas A, Williams I, Bellini WJ, Redd SC. Persistence of vaccine-induced antibody to measles 26-33 years after vaccination. *J Infect Dis* 189 (2004) 123-130.

