Revolution™ Ascend

User Manual

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
5818564-1EN
Revision 1
US English
© 2021 General Electric Company
All rights reserved.

SECRET DE AFACERI
CONFIDENTIAL

Table 8-15 Icons (Table continued)

Number	Icon	Function	Description
8		Fast Speed	When pressed in conjunction with Table Up/Down or Cradle In/Out, causes those functions to operate at a faster speed.
9		Internal Landmark	Press to establish the table's reference point when positioning the patient. This reference point is normally the anatomic reference point used when positioning the patient. For example, if the patient's anatomic reference
10		External Landmark	<p>point is the sternal notch, then the sternal notch would be centered to the laser alignment light.</p> <ul style="list-style-type: none"> • For Internal Landmark, the gantry displays a table location of 0 mm. This sets the zero point for which S and I scan locations are centered around. • For External Landmark, the gantry displays a table location approximately 240 mm from the internal landmark, depending on table characterization. <p>A landmark must be set before you click Confirm settings. At Done Scanning, the landmark is cleared. For scan setup details, see 11.3.3.1 Set up and position the patient on page 262 procedure.</p>
11		Laser Alignment Light	<p>Press to toggle all laser alignment lights on/off and to move the gantry components from the park or idle position to the alignment lights position. Alignment lights are used to establish landmark locations. Three alignment lights are displayed:</p> <ul style="list-style-type: none"> • Axial, divides anatomy into superior and inferior sections. • Sagittal, divides the anatomy into right and left sections. • Coronal, divides anatomy into anterior and posterior sections. <div style="border: 1px solid black; padding: 5px; width: fit-content; margin: 10px 0;">  CAUTION </div> <p> For patient safety, it is important to always have patients close their eyes anytime the laser alignment lights are on.</p> <p>The sagittal, coronal, and transverse alignment lights are within ± 1 mm of the gantry coordinate system (x, y and z-axis) centered at gantry iso center.</p> <p>NOTE</p> <p>The laser alignment light switch is provided as an alternative to beam attenuators.</p>

3. In the Scan Data Manager, highlight the exam, series, and scan you want to delete. Hold down the **Control** key while selecting to choose multiple records.
4. Right click and select **Delete Scan Data**

CT01	1483	
CT01	1482	
CT01	1481	Restore and Lock Scan Data
CT01	1480	Delete Scan Data
CT01	1479	
CT01	1478	

19.8.9 Add reconstruction

You can add a new reconstruction from the raw scan data.

1. In File Manager, highlight any exam to activate the control panel icons.
2. From the File Manager control panel, click the **Scan Data Manager** icon .
3. In the Scan Data Manager, right-click the exam, series, or group for which you want to add a new reconstruction.
4. Click **Open Exam Session**.
5. The Reconstruction and Post Processing window opens. The item you selected in the Scan Data Manager will be selected in the task list. Use the task list to create the reconstruction (as you would do after a scan).

The Exam Summary page opens on the Scan Display (as it does at the completion of a scan).

19.8.10 Save full or anonymized raw scan data

To save raw scan data, you can send it to a USB. You can send the full scan data, or an anonymized version of the data.

Sending anonymized data does not guaranty that data will be rendered anonymous in compliance with applicable data privacy laws. Review the scan data before sharing the data with third parties to ensure compliance with applicable privacy laws or regulations.

Raw data saved anonymously and then restored will reconstruct with patient identifying information removed.

If you select a scan file for saving, you cannot select it for another process while it is in the save queue.

NOTE

Scan must be idle when you start this procedure. Each cardiac scan file takes approximately 5 minutes to save. Scan times can be longer based on the length of the exposure.

Make sure the selected scan files do not exceed the space available on the USB.

1. Connect a USB device to the USB port on the front of the media tower.
2. In File Manager, highlight any exam to activate the control panel icons.
3. From the File Manager control panel, click the **Scan Data Manager** icon .

Table 23-2 Scanner Geometry (Table continued)

Parameter	Revolution™ Ascend	HiSpeed CT/i
Focal spot to det	949 mm	1100 mm
SFOV	500 mm	480 mm
Bore	745 mm	700 mm

23.3.2 Patient Scanning

Table 23-3 Table Scan Length

Scan	Scan Length (mm)		
	VT1700V Table	VT2000 Table	VT2000x Table
Scout Scans	1600 mm	1900 mm	1900 mm
Axial/Cardiac Scans	1730 mm	2000 mm	2000 mm
Helical Scans	1580 mm	1890 mm	1890 mm

NOTE

Helical scan range varies based on the helical pitch and gantry rotation speed selected.

23.3.3 Table Capacity

Standard Weight Capacity Table – Table maximum capacity is 227 kg (500 lbs)

Heavy Weight Capacity Table – Table maximum capacity is 306 kg (675 lbs)

23.4 Network

23.4.1 Purpose of Revolution™ Ascend Scanner Connection to Network

The Revolution™ Ascend Scanner is intended to be connected to a network in order to support the following functionality:

- DICOM services to retrieve images from other DICOM-compliant machines
- DICOM services to push images to other DICOM-compliant machines
- DICOM services to query for images on other DICOM-compliant machines
- DICOM services to confirm that images have been permanently stored on a DICOM-compliant machine
- DICOM services to get DICOM modality worklist information from a remote hospital or radiology department information system computer
- DICOM services to allow a Modality Performed Procedure Step to be communicated to the Hospital/ Radiology information system
- DICOM services to verify the remote DICOM system is connected properly to the scanner device
- Services to provide authentication and authorization against Enterprise directory servers

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5818565-1EN
Revision 1
US English
© 2021 General Electric Company
All rights reserved.

SECRET DE AFACERI
CONFIDENTIAL

Figure 4-25 Representative Images (Reference 21CFR 1020.33(d)(3)(ii))


Table 4-13

Number	Phantom	Technique
1	Helical Head (Typical Noise 0.43%)	120 kV, 140 mA, 1 sec., 5 mm, 0.516:1 Pitch
2	Helical Body (Typical Noise 1.4%)	120 kV, 140 mA, 1 sec., 5 mm, 0.516:1 Pitch
3	Axial Head (Typical Noise 0.42%)	120 kV, 135 mA, 2 sec., 5 mm
4	Axial Body (Typical Noise 1.22%)	120 kV, 135 mA, 2 sec., 5 mm

4.20 DOSIMETRY (Reference: IEC60601-2-44:2009+A1:2012 Clause203.5.2.4.1, IEC60601-2-44:2009+A1:2012 Clause203.109.1)

Dosimetry information is provided in terms of the CTDI and CTDI_w dose indices. Optionally CTDI_{vol} and its associated DLP (dose length product) is automatically computed and displayed on the patient Rx

menu to assist in managing patient dose. This section provides a brief description to help you better understand these dose reporting standards.

4.20.1 General Information

Dose is the amount of energy imparted by the X-ray beam at a given point in an exposed material (patient tissue, phantom, air, etc.) and is measured in units of mGy (milliGray). The old unit was the RAD, which equals 10 mGy. Dose is dependent on the energy absorption factors of the material and on the X-ray exposure. The X-ray exposure is measured in C/kg (coulombs per kilogram) and is dependent on the technique factors used for the scan. An absorbed dose of 1 mGy means that 1 Joule per gram of energy has been imparted. The dose is generally proportional to the exposure, which increases with increasing mA, kV and scan time and decreases with increasing patient size. The X-ray exposure to a point occurs from both direct X-ray from the tube and from scattered X-ray due to adjacent material exposure.

Patient biological risk is related to dose but is also highly dependent on the specific organs exposed and the age and gender of the patient. The effective dose is a way to characterize the stochastic risk to the patient population. The effective dose is the sum of the doses weighted in accordance with the specific radio-sensitivity of the particular organs or tissues exposed. Weighting values are published in ICRP 60 (International Committee on Radiation Protection, Publication 60). The effective dose is a whole body dose equivalent value that has been scaled to represent the dose of the exposed organs. Although we can accurately describe the X-ray exposure potential to a patient for a CT scan, we cannot easily determine the patient dose or risk in terms of effective dose. This is because each patient is anatomically unique and the specific details of his or her anatomy along with the source exposure must be processed using time-consuming monte-carlo computer programs (or other more approximate methods) to predict how radiation will be scattered and accumulated within various patient organs.

Since it is not possible to characterize the specific dose given to individual patients, the CT dose indices are provided to help make relative comparisons. These dose index values can be used to compare CT systems and to help select appropriate operating conditions for scanning. However, it is important to recognize that the dose reported by these indices is inversely proportional to phantom size (see [Figure 4-26 Relationship between dose and phantom size for head and body filters at 120kV. Similar curves are obtained for the 80, 100, and 140 kVs. on page 108](#)). This means that for the same scan technique (protocol), smaller phantoms (patients) will produce a higher absorbed dose than larger phantoms (patients) - see "Influence of phantom diameter, kVp and scan mode upon computed tomography dose index", Edward L. Nickoloff, Ajoy K. Dutta, and Zheng F. Lu, Medical Physics 30, 395 (2003). **Therefore, it is critical to remember that the body scan FOV's uses the 32cm CTDI phantom and all pediatric and head filter uses the 16cm CTDI phantom for dose reporting purposes (CTDI_{vol} display on Scan Rx Menu).** [Table 4-15 SFOV selection vs. CTDI phantom used for dose reporting \(IEC60601-2-44:2009+A1:2012 Clause203.112\) on page 108](#) indicates the phantom size used for each SFOV.

At this manual, dosimetry information was measured with the following instrument (Calibration Tolerance: $\pm 5\%$):

1. Radiation Meter RadCal 9015
2. Radiation meter Converter RadCal 9060
3. CT Probe RadCal 10x5 - 3CT (for CTDI)
4. Ion-chamber RadCal 10x5 - 1800 (for stray (scatter) radiation)

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
AW Server 3.2

Powering advanced applications across the healthcare enterprise

Continuous improvement in image acquisition capabilities are rapidly expanding your everyday clinical practice and providing new ways to improve patient care. As a result, a solution for manipulating and distributing access to volumetric images becomes as important as the acquisition modalities themselves.

AW Server is your platform supporting a broad portfolio of efficient and automated workflows that allow you to manage your time, equipment resources and costs while enhancing your team's ability to collaborate and provide diagnoses. The AW Server adds speed, efficiency, and diagnostic flexibility to your workflow. Its applications portfolio helps boost your diagnostic confidence as you analyze and evaluate exams from angiography to X-ray and almost everything in between. With its novel, easy to use, intuitive interfaces and reporting flexibility, the AW Server streamlines your workflow and helps make your entire department more productive.

Overview

AW Server 3.2 introduces an extensive portfolio of advanced applications that is virtualization-ready. With a local cache receiving DICOM images from your modalities and priors from your PACS, this streamlined advanced visualization workflow engine enables rapid preparation and communication of 3D results throughout the enterprise with access to innovative AW applications.

What's new

- Designed with VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic confidence and productivity.
- Simplified user experience with powerful support for advanced applications workflows.
- Virtualization-ready platform provides enterprise-wide access to advanced processing applications.
- Support for dual monitors expands screen space to hang more views.
- New hardware supports expanded number of slices - up to 160,000.
- PACS integration interface to streamline advanced visualization workflow for improved productivity.

Visit us:

www.gehealthcare.com/aw/platforms/aw-server/

Features

- Choice of IT-friendly rack mount chassis or virtualized deployment on VMWare® environment
- Receives DICOM® images from modalities and PACS
- Patient list for management of images on server cache.
- Advanced search for fast and easy search of a patient's exam history on PACS or any other DICOM-compliant device.
- **Enhanced Quick Filters of the Patient List filters studies by Modality, Date, End Review Status or Exam Description.**
- Offers pre-processing for automation and acceleration of workflows
- Powers advanced visualization and image processing applications
- Delivers diagnostic quality client software to PC and MAC®¹ clients and specific GE Healthcare CT consoles
- Dual monitor client support for symmetric displays²
- 2D Viewer for image display, manipulation, annotation, review.
- Integrated Filmer with enhanced flexibility to perform filming and data exporting tasks.
- End Review automates routine filming and networking tasks with just one click.
- Supports Saved State selection
- Provides common tools by 2D and 3D applications, e.g. Window/Level, Zoom, Pan, distance measurement, 2D ROI selection²
- Supports core AW Volume Viewer³ software

- Supports Advanced Applications² for automated post-processing vessel analysis and oncology with streamlined reporting
- Supports IT remote "lights-out" management
- 3rd Party PACS integration supports launching advanced applications directly from the PACS workflow.

Server Requirements

AW Server 3.2 may be purchased as a turnkey solution that includes off-the-shelf enterprise-class server hardware supporting concurrent processing of up to 160,000 slices.

Alternatively, AW Server 3.2 may be purchased for use with existing VMware installations for large enterprises that choose to centralize their infrastructure. AW Server 3.2 is also available as a virtual machine whose tight integration with GE Universal Viewer provides a single environment for 2D and advanced visualization radiology workflows.

Indication for Use

AW Server is a medical software system that allows multiple users to remotely access AW applications from compatible computers on a network. The system allows networking, selection, processing and filming of multimodality DICOM images. Both the client and server software are only for use with off the shelf hardware technology that

meets defined minimum specifications.

The device is not intended for diagnosis of mammography images. The device is not intended for diagnosis of lossy compressed images. For other images, trained physicians may use the images as a basis for diagnosis upon ensuring that monitor quality, ambient light conditions and image compression ratios are consistent with clinical application.

Notes:

- ¹ Requires purchase of Parallels® 10 and Windows® 7, Windows 8.1, Windows 10 32-bit or 64-bit software.
- ² Maximum resolution recommendations vary based on the network bandwidth and latency; maximum resolution might be limited by the specific clinical applications.
- ³ Requires purchase of appropriate application licenses.

Product Details

2D Viewer

The 2D Viewer is an application used to display, manipulate, annotate and review 2D images by a trained physician for diagnostic interpretations.

Display customization

Display customization allows the user to manage layouts to display data. Key benefits include:

- Allows modification of number of exams or series displayed to facilitate either single or multi-exam reviews.
- Standard layouts give you flexibility to tailor the image display from 1x1 to 8x8. The current layout is retained if the images are from the same modality and orientation.
- Annotation levels allow selection of the image information fields to display.
- Toolbar customization puts you in control of which buttons are displayed to meet your needs.

Study Navigation

- The navigator lets you assign a series to a view on the fly.
- Cine mode also provides temporal, spatial, or manual playback loops.
- With two exams or series loaded, Cine mode supports a side-by-side display format with synchronized playback loops for more efficient comparisons.

Image Review

- The initial image window and level setting is based on the DICOM header. Once displayed, several methods of adjusting image window and level are provided.
- **Provides routine image manipulation features:**
 - Flip/Rotate, Zoom, Pan, Magnifying Glass,
 - Inverse Video: inverts grayscale color map,
 - Display normal: lets you return the image to its default viewing parameters.

You can access the following features using a single mouse click directly on the image. Having these direct manipulation tools easily at hand gives you fast access with less distraction from your review task by eliminating the need to return to the graphical user interface controls.

- Several tools are available:
 - 2D distance
 - Angle

- Report cursor
 - Box ROI, elliptical ROI and free-hand ROI
- The program continuously updates statistics on the fly.

- The Annotation feature lets you highlight areas of interest by adding text and line/arrow pointers to any image feature.
- The Copy/Paste/Erase feature lets you copy/paste/remove any text or graphic placed on an image.
- Cross-reference indicates the position of the current slice over the scout or localizer or any other non-parallel series.
- Save lets you store a copy of the image as it appears on the screen for future review.
- The Key Images feature allows you to flag images and create key objects.

Batch Filming

The Print Series feature lets you automatically batch film an entire series with a single keystroke.

Filmer

The integrated Filmer enhances the efficiency of the review station and gives you greater filming and exporting flexibility. This feature supports two modes: the Mini Filmer mode and Full-Screen mode, which provides the ability for customization and film layout template creation.

Three key mechanisms give you flexibility:

1. Free format filming
2. DICOM structured reporting (SR)
3. Data export (HTML/PDF and JPEG, PNG, MPEG, or AVI)

With the Filmer you can easily extract significant images from any AW application (2D Viewer, Volume Viewer¹ etc.).

- You can film images individually by dragging and dropping to the on-screen Filmer, or by a single keystroke.
- Multiple image formatting lets you film multiple images in a single page frame
- With Film MID you can send multiple images to a single Filmer frame
- Batch filming is supported from applications that provide that capability (Print Series in the Viewer, Batch Film Protocols in Volume Viewer¹).

Communication between applications and the Filmer is accomplished in the Mini Filmer mode, which provides a minimized footprint. The Mini Filmer mode provides the following:

- Compatibility with Batch Filming from Volume Viewer¹

- Store and position images transferred from an application.
- Type of export (film, media, database).
- Rapid switching between Full-Screen Filmer and application.

The resultant electronic films become a quick summary of the patient study and radiological interpretation to be reviewed by clinicians and physicians. They can also be saved as independent files (DICOM SR and Secondary Captures) for teaching purposes.

The Filmer works on the What You See Is What You Get (WYSIWYG) approach and the electronic film can contain one or several pages with specific layouts for each page or all pages.

A flexible Edit Mode provides the ability to easily add, manipulate, format, or delete images from the film. Images can contain text and graphics from measurements and your annotations, and may be window/leveled, magnified, flipped, rotated, or cine. You can add additional annotation to the image in edit mode.

The Preview Mode displays the film as it will be printed or exported, taking into account the layout applied to each page and the compression level specified for non-DICOM exporting (JPEG/PNG and MPEG/AVI/QTVR).

Printing

- AW Server 3.2 includes printing from local client Windows application and Network DICOM Printing (B&W and color).

Volume Viewer Application

- The AW Server supports AW's 3D image analysis package with a broad range of multimodality image analysis tools including Multiplanar Reformat, Volume Rendering and Navigation.
- For detailed description of Volume Viewer features, please refer to the Volume Viewer Product Data Sheet.

Workflow Capabilities

- 3rd party PACS integration provides an interface to launch advanced applications directly from the PACS client in the same context¹.
- Ability to manually adjust image compression during user interaction including lossless to multiple lossy levels to improve performance for sluggish network situations. Static images will always be displayed at full resolution.

- Support for Saved State creation and selection.
- Common tools by 2D and 3D applications (Window/Level, Zoom, Pan, Distance measurement, 2D ROI selection).
- AW Server 3.2 provides local DICOM cache allowing rapid launching of 3D applications and temporary storage of derived imagery.
- Allows you to save the exam to the PACS database in an intermediate post processing state at any time and then to restore it, allowing multiple radiologists or technologists to contribute to post processing studies.

IT and Security Capabilities

- Active directory integration allows enterprise level user authentication. Directories supported include Microsoft® Active Directory® (MSAD) and other LDAP authentication services
- Provides a centralized user management system in case integration with enterprise directory is unavailable.
- Centralized web-based administration interface for configuration, diagnostics, usage and day-to-day server management.
- Provides group privileges for full access to exams or to enforce mandatory exam search by configurable criteria.
- Centralized Enterprise Audit trails capture user access to patient data using industry standard protocols: IHE ATNA compliant audit messages (RFC3881), TCP, BSD Syslog, and Reliable Syslog (Cooked profile Experimental) transport to third party audit repositories. Audit logs are kept on local storage and managed with log rotation.
- Display of annotations and patient protected information can be controlled for privacy
- Provides support for deployment of client software through Microsoft Windows® Installer (MSI), making it easy to manage enterprise client PCs.
- Features a stateful server firewall to track network connections and detect malicious intrusion attempts.
- Allows thin client access over the internet with either http port (80) or https port (443).³

Notes

¹ Requires purchase of separate application licenses.

² Please refer to Volume Viewer datasheet for details.

³ For security reasons, GE Healthcare always recommends VPN access to the server.

Performance

Performance and interactivity of AW Advanced Applications depends on the network bandwidth, latency and client workstation configuration. For optimal performance, a minimum network bandwidth of 40 Mbps with a latency of 2.7ms or lower is required. High bandwidth such as a hardwired local area network (LAN) permits faster interactions at lower compression ratios.

AW Server may also be used over WIFI or WAN/Internet although performance will depend on round-trip latency between the client workstation and the server. A minimum of 3 Mbps bandwidth per client with latency less than 50ms is recommended for reasonable performance when compression is used.

To optimize performance, GE Healthcare's "Smart Compression" technology applies a user-selectable level of compression to displayed images only while you page, browse or rotate the images. Images are displayed automatically at full fidelity once user interaction stops. A clear visual display on the images indicates when compression is being applied. The device is not intended to meet diagnostic quality criteria with lossy-compressed images.

The AW Server 3.2 supports the following compression levels:

- Lossless
- Lossy 15:1 (compression ratio is 15:1)
- Lossy 22:1 (compression ratio is 22:1)
- Lossy 33:1 (compression ratio is 33:1)

Turn-Key Server Specifications

AW Server 3.2 is deployed on one of two physical server configurations supporting three sizes: L, XL and XXL. The L size supports 3 active concurrent users¹ sharing 40,000 slices in memory. **The XL size supports 6 concurrent users sharing 80,000 slices and the XXL up to 12 users sharing 160,000 slices.** Turn-key server hardware specifications listed here are minimum server hardware capabilities. Specific server hardware is subject to change without prior notice. AW Server 3.2 software is also supported on some previously deployed AW Server systems although slice counts⁵ may be restricted by available memory.

Server Hardware L

Processor 2 Intel® Xeon® E5-2640 10-core CPUs
Memory 64GB RAM⁵ supporting 40,000 slices
Disk drive 300GB RAID -1 OS disk, 2.4TB RAID-6 storage disk
Network 4 1-Gbps NICs
Management Dedicated Embedded Lights Out Manager
Power Fully redundant power and cooling
OS GE HELIOS 6.8

Server Hardware XL, XXL

Processor 2 Intel Xeon E5-2697A 16-core CPUs
Memory 256GB RAM⁵ supporting 80,000 (XL) or 160,000 (XXL) slices
Disk drive 300GB RAID-1 OS disk, 7.2TB RAID-6 storage disk
Network 4 1-Gbps NIC
Additional NIC: 2 ports 10Gbps
Management Dedicated Embedded Lights Out Manager
Power Fully redundant power and cooling
OS GE HELIOS 6.8

Virtualized Server Specifications

AW Server 3.2 is available for deployment as a software only system on existing VMware installations. Customers must provide a VMWare environment for which CPU and RAM overcommit are not recommended. The use of hyperthreading is not recommended and if enabled may result in performance degradation for compute-intensive features of AW Server. Hard disk drives are required to store all VM data with thick provisioning.

- Intel Xeon CPU with SSE® 4.1 (45 nm or better)
- VMware ESXi® 5.0u1, 5.1, 5.5, 6.0

VM Configuration Requirements

Virtual machines supporting AW Server must be configured with the following resources:
L: 8 core vCPU, 24GB vRAM, 2 NICs, 70GB OS, 2TB vHDD (IOPS 800 or better). L system supports up to 3 concurrent users¹ with 8,000 shared slices
UV-L: Universal Viewer integration requires 8 core vCPU, 24GB vRAM, 2 NICs, 70GB OS, 2TB vHDD (IOPS 800 or better). This deployment supports up to 4 concurrent users¹ with 16,000 shared slices
XXL: 24 core vCPU, 64GB vRAM, 2 NICs, 70GB OS, 2TB vHDD (IOPS 800 or better). XXL system supports up to 12 concurrent users¹ with 40,000 shared slices and can also be tightly integrated to support Universal Viewer.

Server Licensing

AW Server 3.2 may be purchased either as a single physical installation or as a virtual installation. Application licensing is structured to allow easy addition of clinically specific workflow tools. Applications licenses are sold separately, allowing you to purchase only what you need to support your facility.

Client Requirements

AW Server delivers user-installable client software to enable access to applications and data on the server.

The minimum HW configuration is as follows:

Processor	2.33GHz Pentium®4 minimum, Pentium dual core processors recommended for optimal performance
Memory	1024 MB minimum
Disk drive	250MB free space available
Network card	100 Mbps minimum (1000 Mbps recommended)
Internet connection	Customer-provided IPSEC VPN, for internet/WAN operation
Mouse	Two or three-button mouse. Three button mouse suggested for best use of functions

Monitors

AW Server 3.2 supports identical dual monitor clients with a recommended combined resolution of 4MP (2x2MP), larger monitors may be used³. Advanced applications may take advantage of the dual screen support to allow larger or more numerous views.

Screen resolution

Minimum: 1024 x 768 with full color (24 bit)

Recommended: dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048)

Maximum: 6MP combined³

Supported client OS software

- Windows 7 SP1 32 and 64 bit,
- Windows 8.1 32 and 64bit
- Windows 10 32 and 64 bit
- Mac Parallels® (Mac OS X 10.10, Parallels 9, Windows 7 SP1 32bit)

User interface languages

- Chinese⁴
- Danish
- Dutch
- English
- Finnish
- French
- German
- Italian
- Japanese⁴
- Korean⁴
- Norwegian
- Portuguese
- Russian⁴
- Spanish
- Swedish

Keyboard layouts

- Danish
- Dutch
- Finnish
- French
- Canadian French
- German
- Italian
- Norwegian
- Portuguese
- Spanish
- Swedish
- United States English
- United Kingdom English


Standards

- DICOM Storage Service Class for RT, CT, MR, CR, X-ray (Angio and R&F), Digital X-ray (DX), MG, NM, PET, Key Image Notes (KIN), Structured Reporting (SR), Grayscale Softcopy Presentation State (GSPS), U/S, Secondary Capture, Secondary Capture Color DICOM Image Objects.
- DICOM Service Class User (SCU) for image send and Service Class Provider (SCP) for image receive
- DICOM Storage Commitment Service Class User (SCU)
- DICOM Print (Color and B&W)

Regulatory Compliance

This product complies with:

- European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC
- EN 980:2008 - Symbols for use in the labeling of medical devices
- EN 1041:2008 - Information supplied by the manufacturer with medical devices
- ISO 14971:2007 / EN ISO 14971:2012 -Medical devices - Application of risk management to medical devices
- IEC 62304:2006 / EN 62304:2006 +AC:2008 - Medical device software Software life-cycle processes
- IEC 62366:2007 / EN 62366:2008 - Medical Devices - Application of usability engineering to medical devices

Notes

- ¹ Number of supported clients depends on the server configuration, purchased licenses, network and client monitor resolution.
- ² Capacity calculated at the advertised value of the hard disk manufacturer, available disk space might vary.
- ³ Maximum resolution recommendations vary based on the network bandwidth and latency; maximum resolution might be limited by the specific clinical applications.
- ⁴ Limited to Windows clients only, not available on Linux clients.
- ⁵ Not all hardware platforms will support or be upgradeable to all RAM amounts or slice counts.

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems. Our "healthymagination" vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

Manufacturer address:
GE Hungary Kft.
Madarász Viktor utca 6.
1131 Budapest
Hungary

European registered place of business:
GE Medical Systems SCS
283 rue de la Minière
78530 Buc France

Copyright ©2016 General Electric Company. General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Timing and availability remain at GE's discretion and are subject to change and applicable regulatory approvals. Contact your GE representative for the most current information.

GE and the GE Monogram are trademarks of General Electric Company.

Microsoft, Windows are registered trademarks of Microsoft Corporation in the United States and/or other countries.

Mac and Mac OS are trademarks of Apple Inc., registered in the U.S. and other countries.

Intel, Core, Pentium and Xeon are trademarks of Intel Corporation in the United States and/or other countries.

Parallels is a registered trademark of Parallels Software International, Inc.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

All other trademarks are the property of their respective owners.

DOC1830515

GE Healthcare, a division of General Electric Company.

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
DOC1830515

GE Healthcare, o divizie a General Electric Company.

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Volume Viewer

GE Healthcare's premium 3D Advanced Visualization and Image Processing platform brings you powerful capabilities to help you visualize and interpret your multimodality imaging data with confidence and ease.

AutoLaunch

Preprocessing enabler

Imaging modalities provide a wealth of diagnostic information, but also present a challenge to your Radiology workflow as image volumes continue to increase while your time for reading and reporting decreases.

Volume Viewer provides you with a rich 3D image processing toolset aimed at creating and displaying the views you need with little user input and streamlining interpretation and reporting by providing the visualization tools you need with minimum clicks.

With an intuitive, modern user interface, Volume Viewer helps a wide diversity of users to learn and master the expanding portfolio of sophisticated tools and applications that it provides. This customizable user interface maximizes the real estate allocated for displaying images and provides access to the users' favorite tools directly from the image viewport. In addition, a customizable toolbar provides one click access to commonly used tools to facilitate a productive review.

Overview

Volume Viewer provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-ray, PET, PET/MR and PET/CT datasets. Volume Viewer also features a broad portfolio of high-performance analysis tools, automating routine tasks and helping to make 3D image processing a stress-free component of your routine workflow.

Volume Viewer is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- Customizable protocol-driven workflows supporting CT, MR, PET and Innova 3D XA modalities, designed to help enhance your clinical review productivity
- Compare multiple exams from different modalities, at loading² and on the fly²
- Adaptive layouts that take advantage of the dual monitor capabilities in landscape and portrait orientations
- Easy point-and-click segmentation tools for segmenting structures of interest
- CT Dual Energy and GSI⁴ support with dedicated workflows and tools
- Advanced 4D PET support with Q Freeze⁶ tool and 4D PET Dynamic workflow
- Interactive Summary Table collecting measurements as they are deposited on the images.
- Rich set of export tools including export for 3D printing³.

<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/volume-viewer>

SECRET DE AFACERI

CONFIDENTIAL

General Features

- **Unified user interface**, with one click access to your favorite tools and more space for displaying images.
- **Customizable protocol-driven workflows** designed to help enhance your clinical review productivity for a wide range of care areas, with 3D and 4D real time interaction.
- **“Smart layout” protocol**, which adapts the screen layout to the type of loaded exams.
- **High resolution, real-time rendering modes:** MIP/Min IP/Average; standard volume rendering and enhanced volume rendering with Volume Illumination¹ on the fly; Navigator view for endoluminal fly-through; lumen and curved views.
- **Compare mode** for reviewing multiple exams and modalities, simplified with **Dynamic Load to load and register²** new series into the current review session.
- Rich set of **2D/3D ROI tools**, including **AutoContour** for 3D semi-automatic contouring of structures of interest (CT, MR, PET), and **ROI color coded** for tissue classification based on voxel values.
- **Advanced tools**, to take advantage of full 3D capabilities: **AutoSelect**, for easy point-and-click segmentation; One or two clicks **Quick Vessel Trace** to analyze all vessels, in curved reformat, lumen, or MPR view; **Advanced 3D processing tools**.
- **Summary Table** extends reporting capabilities, while collecting measurements as they are deposited on the images.

- **Multiple Export capabilities**, in different DICOM and non-DICOM formats, including **Key Image Notes**.
- **3D Suite³:** STL, OBJ, VRML and 3MF file formats to export single or multi mesh objects, for 3D printing. 3D Suite is not intended for clinical use.
- **Save state**, which allows to save and restore to pause and resume the review at a later time or to facilitate collaboration workflows.

Modality Features

- **Dual energy and GSI⁴** GE Healthcare CT images supported with dedicated review protocols.
- **Bone VCAR⁵** combines spine labelling on the fly for all CT review and a dedicated layout to improve spine review.
- **MR review support:** Recognition of PSD name, and specific MR parameters at image loading, which simplifies protocol customization.
- **Dedicated review workflows for PET/CT and PET/MR images**, including **PET SUV measurements**.
- **Q.Freeze⁶** combines the quantitative benefits of **4D PET gated** imaging into an image that provides both frozen patient motion and reduced image noise.
- **Dedicated protocol-driven workflows for Innova 3D XA images⁹**.

Platform features

- **“Smart Compression” technology** for fast and responsive client performance on AW Server.
- **Fast automated Pre-processing⁷**

- Multi-sessions switch on AW workstation with **AutoLaunch⁸**.
- Support of **Chinese, Japanese, Korean and Russian languages**.

System Requirements

- **AW Server 3.2 and above**, with recommended monitor resolution up to dual 2MP or single 3MP.
- **AW 4.7 Workstation and above** (some performances limitations can occur with Z800)
- **CentricityTM Universal Viewer**

Intended use / Indications for use

Volume Viewer is a medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography and PET Scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, surgery and treatment planning.

Regulatory compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its feature may not be available in some other countries or regions. Please contact your sale associate.

Rx Only

¹Requires Volume Illumination license

Volume Illumination may not be available in all countries or regions.

Please contact your sale associate.

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBoneTM Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

Volume Viewer


General Features

Unified user interface

- The user interface provides one click access to favorite tools and more space for clinical images.
- The page layout is organized to facilitate navigation through the review steps. Guided protocols walk the unfamiliar user through a clear workflow of creating and manipulating volumetric data, helping to reduce the learning curve by providing a consistent workflow.
- The customizable toolbar can be docked at the top, bottom, left, or right, by a simple drag & drop. This preferred location is saved for future use for each user.
- An extensive set of advanced tools are grouped into categories on the toolbar to facilitate quick retrieval when needed.
- Each category and tool are easily configurable by a simple drag & drop in the Toolbar so that they are only one click away for your future reviews
- To reduce long mouse trips, up to 6 favorite tools can be accessible from each image via a right mouse menu.
- Customization of the toolbar and tool behaviors with a user's preferences is provided in a single intuitive Preferences menu.
- Keyboard shortcuts and quick access tools for on-viewport controls

The following mouse modes are available to help you manipulate the images:

- Interactive VR adjustment lets you adjust VR opacity interactively with the mouse as Window Level/Window Width
- Direct Paging (continuous or contiguous)
- Free hand rotation on 3D and oblique views
- Percentage fusion, to easily change the transparency of objects on fused images. There are several fusion modes available for your use.

- Mag glass displays zoomed-up imagery around the cursor.

Protocol management and loading

- Customizable protocol-driven workflows for a wide range of clinical applications, including support for multiple review steps. You are guided through the review steps by a Navigation bar at the top of the User Interface.
- Review protocols created by user with Volume Viewer can be saved as favorite protocols and accessed easily in routine.
- Progressive Load enhances image loading performance. This feature has two modes:
 - Images are loaded in sequential mode (e.g., Reformat)
 - Images are loaded in interlaced mode (e.g., 3D/VR)In either case, image review can begin as soon as the first image is displayed.
- Navigation through series and exams can be performed without exiting the patient list.

Smart layout

- Volume Viewer General Review protocol presents an adaptive layout for single and multi-modality PACS-like reading of data. This protocol distributes and displays all loaded series of a single exam in an intelligent layout to maximize use of available monitor space. With multiple exams, the protocol automatically detects and distributes "like" series for optimal comparison.

GE imagination at work

Advanced visualization capabilities

Volume Viewer offers various capabilities to display advanced rendering modes:

- 3D Volume Rendering / Volume Illumination¹ images
- MPR views with different slice thickness. The following rendering modes are available for all Thick Slab: MIP, Min IP, Volume Rendering, Average.
- Navigator views which display interactive endoluminal views exportable as movies.
- Curvilinear reformatting allowing display of curved, lumen and cross section views for various structures (vessels, spine, etc).
- Merge multiple volume rendered models into a single view.

Volume Rendering

Volume Illumination


Predefined cut planes

Volume Viewer lets the user define cut planes to isolate specific structures in the VR model, in 16 configurations, including Left / Right / Inferior / Superior / Anterior / Posterior / Front. You can then display the voxels values on each plane.

Multiple Volume Rendering and Volume Illumination¹ objects

- Volume Viewer allows you to merge up to 8 Volume Rendered models from the same series into a single 3D view.
- You can adjust independently the threshold, colors, and transparency of each Volume Rendering or Volume Illumination

Auto-center fly through with smart cursor

- In fly-through studies (airways, colon, angiography) navigation along the centerline of the structure of interest can be recorded step by step and can be exported as a movie.
- In addition to the regular 180° or less Navigator view, Volume Viewer enables a Fish-Eye View with any wide camera angle value from 180° to 360°. This provides a view of structures both in front of and behind the users' virtual location on the same image.

Lumen view

- Lumen View provides an unfolded 3D view around a user defined centerline. The lumen view can be interactively adjusted (rotation around the centerline, width, field of view).

3D and 4D capabilities

- Volume Viewer allows 3D cursor synchronization in any orientation on the fly.

GE imagination at work

- Real time reformat in oblique planes is available for your exams, as well as simple and double oblique interactive modes.
- Quick access to cross reference display gives you the ability to display the current slice location or all slices, or other planes.
- Real time interaction in 4D mode, by using Cine tool on all planes and 3D views for CT, MR, and PET multiphase data.

Compare mode


- With Volume Viewer several volumes can be loaded either from the same exam or from different exams. This is useful for exam comparison (Compare Mode) or for Multiphase examinations (of the liver for example). The Compare Mode is also available for PET/CT examinations and takes advantage of a dual screen configuration.
- Multi-phase CT and MR data are listed as separate phases in the Series Selection screen, allowing selection of all phases or a subset of phase data for review. The Series Selection Panel displays the multi-phase data split into arterial and portal phases, for example.
- Dynamic Load, compatible with both Volume Viewer and Fusion protocols, lets you drag & drop 3D volumes from CT, MR, PET, and 3D XA modalities into a desired viewport. Together with the Integrated Registration² option, Dynamic Load allows new volumes to be registered and loaded on the fly. Save State series can be restored as a separate session. For dual monitor configurations, a Save State series can be displayed on the right monitor together with a current session on the left monitor. For single monitor configurations, a user may simply switch between the two sessions.
- The zoom and pan functions are propagated to all images displayed in the same orientation.

2D/3D ROI for quantitative measurements


- Volume Viewer contains a set of basic 2D/3D tools: distance, angle, report cursor, arrow tool, annotation, free hand ROI tool, that can be placed on fused images as well.
- Easy deposition and labelling of measurements.
- All 3D ROIs can be customized, and color coded to display statistics computed on different intensity range and help assessing tissue classification based on voxel values.
- All measurements are considered bookmarks and can be accessed by navigating through the list of measurements.

AutoContour tool

- AutoContour provides consistent tools for 3D contouring ROI, with the same workflow used for CT, MR, and PET images.
- For CT and MR images, the contouring tool will suggest contour boundaries after defining 1 point or a diameter on a structure of interest, which can be edited prior to validation.

- For PET images, the contouring tool has been adapted for standard SUV calculations.

Advanced 3D tools

Segmentation tools

- AutoSelect segmentation tool helps adding or removing any continuous structures of interest by simple point and click.
- You can also access advanced segmentation tools such as Scalpel, Threshold, Dilate, Erode, Subtraction, Intersection, Addition, Filter Floaters.

Quick Vessel Trace (CT and MR images)

- You can perform vessel analysis with one or two points clicks from any CT or MR Angiography images.
- Quick Vessel Trace automatically extracts the vessel after user clicks and launches vessel tracking for fast review in curved reformat, cross-sectional, best L-section, lumen, and MRP view.


GE imagination at work

- Tracking can be extended proximally or distally for a full view of the vessel.

Summary table

- The Summary Table collects measurements and organizes them according to finding and exam date.
- The Summary Table allows interactive navigation through findings in the image viewports and can be docked or expanded.
- The Summary Table can be copied into the client clipboard and then pasted in email, reports, office software etc. (on AWS only)

Export capabilities

Volume Viewer contains multiple standard options for exporting the results of a review session:

- Save images to the database, as new DICOM series.
- Save the significant images as Key Image Note objects in the database. Key Image Notes and End Review allow you to flag images of interest as Key Image Notes (IHE profile) and push them to the archiving system when you exit the application.
- Save State object are used to save the work in a new series containing all the post processing data (3D Model, displays, measurements, annotations, etc), for future review.
- Save STL, OBJ, VRML or 3MF file from 3D Volume Rendering / Volume Illumination¹ viewport for 3D printing purposes through 3D Suite³

Additional export tools are available to create specific sequences of images to be exported:

- The Batch tool allows creating a sequence of reformatted images or a sequence of rotating 3D views.
- The Movie tool creates a comprehensive movie including different rotations, zooms, and pan of the image, which can be exported as DICOM series or mpeg file.
- Quick Export: Exports in a single click a full batch of contiguous images at the displayed thickness for 2D images, or a batch of rotations of a 3D View.
- Cardiac Review and Export: Processing and reviewing cardiac exams for CT, MR and PET with manual oblique


reformatted protocols can be exported into a multi-phase Cine movie that allows the referring physician to review the exam in a dynamic mode.

- Capture and send images from the viewport to your computer clipboard on the fly (Ctrl+C/Ctrl+V on AW Server only).

Modality Features


CT image quality

- CT filters to enhance or smooth CT images to help improve image quality on the fly.
- The lung filter helps enhance contours of images reconstructed in standard mode for excellent visualization of lung structures

Standard recon w/o filter Standard recon with Lung filter

- A smoothing filter reduces noise while maintaining high image contrast.

No filter Smooth filter

No filter Smooth 3D+ filter

CT Dual energy protocols

- Three protocols are available to enable review of dual energy images acquired on compatible GE Healthcare scanners and enhancing pixels with specific ratio: DE Calcium ($HU_{80}/HU_{140} \geq 1.25$), DE Uric Acid ($HU_{80}/HU_{140} \leq 1.25$) and DE Custom (custom threshold on HU_{low}/HU_{high}).


GE imagination at work

- In addition, all Volume Viewer standard protocols are compatible with GSI⁴ GE Healthcare CT images.

Bone VCAR⁵


- Bone VCAR is available in Volume Viewer toolbar and can be launched on the fly
- Vertebrae are automatically recognized and labelled (deep learning algorithm)
- A dedicated protocol provides a specific layout of the spine including curved views

Review of MR studies

Volume Viewer offers support for review of MR studies, through dedicated protocols.

- The MR-specific parameters (PSD name, image weighting, scan plane, etc) are recognized at loading and are used to customize the layout display.
- Anatomy-based protocols for standard review (e.g. spine, brain stroke...), with dedicated and guided review steps, to help review standard MR exams.
- Dedicated MR Cardiac Viewer for single and multiple Cine Review. Enables comparison between cardiac series such as Time Course and Myocardial Delayed Enhancement (MDE) images.
- Support of MR Multi Echo, Multiphase, Diffusion series.
- Support of parametric series (e.g., ADC, MTT, and other parametric maps) created in READY View, allowing direct functional measurements.


- 2D and 3D ROI propagation to other phases / series / exams allowing you to easily correlate information from multiple sources.
- Improved MR annotation consistency.
- In combination with Integrated Registration², enables direct access to MR image registration (inter/intra exam), image fusion and Whole Body MR Review protocols
- Ability to automatically bind several MR axial series corresponding to sequential axial locations into a single series. This is useful for consolidating multi-stage acquisitions for displaying Whole Body MR scanning.

SUV measurements on PET images

- Volume Viewer supports Standardized Uptake Values (SUV) for image display and measurements. Several SUV scales are available like SUVbw, SUV lbm, SUVbsa, as well as SUV Peak.
- Window/level presets may be user-defined for PET images based on SUV values.
- The SUV values are available in all the basic 2D/3D ROI tools of Volume Viewer and in the AutoContour tool.

PET/CT and PET/MR dedicated protocols

- Volume Viewer contains a list of predefined protocols allowing the review of PET/CT and PET/MR data. The screen layouts contain fused views between the morphological and functional images, as well as 3DMIP rendering of PET images.
- Factory protocols are customizable for your own review.
- PET 4D protocol allows users to load and display gated and dynamic PET series, as well as summing or reframing these series.

GE imagination at work

PET IQ improvements

- Enhanced 3D visualization algorithm to fully support Q.Clear PET images.
- Q.Freeze⁶ option provides a respiratory motion-corrected Gated PET volume, statistically similar to conventional static PET with significantly reduced or eliminated blurring effects due to patient respiration. It provides also a corrected Gated PET series where counts are summed back to each respiratory phase.

Dedicated features for 3D XA images⁹

The optional Innova 3D XR application produces XA images (512x512 images or 256x256 images depending on your settings). This format opens the Volume Viewer world to Innova 3D XA images:

- Single exam: Innova Navigator, Curved Reformat, Multi Oblique, MPVR 3mm Average, Segment Structure
- Compare exam: Compare Axial, Compare Coronal, Compare Sagittal, Reformat
- Layouts and protocols customization
- The Innova options provide access to Workflow booster tools, easy multi-modality review "on-the-fly" dedicated Innova review protocols with dedicated Review Steps.

Smart compression

- The Smart Compression technology automatically displays images with full fidelity when the image is still, then uses the selected compression level for increased interaction speed during user interaction. This allows for diagnostic reads on full fidelity static images with responsive dynamic display even at low bandwidth. On-image visual indicators notify the user when compression is in effect.

Pre-processing⁷ features on AW Server

- Preprocessing is a Volume Viewer feature that automatically performs routine processing tasks and saves the results so they are ready for you. When new exams are transferred to the database, Volume Viewer recognizes user-defined keywords in the Series Description and launches the appropriate pre-processing. The results are stored in a Save State object with the original study so you can load it when you are ready for review.

AutoLaunch⁸ features on AW workstation

- AutoLaunch is a Volume Viewer feature that automatically preloads exams, so they are ready for your review.

- When you are processing a study and a new exam is transferred to the AW Workstation, this feature automatically launches the new exam with an appropriate Volume Viewer protocol in the background. When you are ready, one click in the AutoLaunch window instantly brings up the exam in the Volume Viewer application, eliminating waiting time and extra steps to load the exam into computer memory for processing.
- AutoLaunch is compatible with CT, MR and PET single volume protocols of Volume Viewer.
- When combined with optional applications⁸, this feature gives access to data already preprocessed, that you can review by just clicking on AutoLaunch window.

Summary of Operation

Volumetric models are loaded by selecting the exam or series. The user can select a protocol category from an anatomical selector or go directly to a Review Layout. In either case, images are loaded progressively in the background; this gives control to the user in just a few seconds after selecting the images. Selecting a Review Layout launches a volumetric display protocol with predefined layout preferences. Review Layouts may be customized and combined to suit the user's workflow. Selecting a protocol category unlocks a variety of visual protocols that include the layout, threshold, rendering mode and filming formats. Some of these protocols direct the user through the process providing capabilities to interactively view and manipulate the model, increasing productivity and consistency for all modalities.

References

¹Requires Volume Illumination license

Volume Illumination may not be available in some other countries or regions.

Please contact your sale associate

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBone™ Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

GE imagination at work

GE Medical Systems SCS
283 rue de la Miniere
78530 Buc France

DOC1460639 rev9

General Electric reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Illustrations are provided for informational purposes, and your configuration may differ.

This information does not constitute legal, financial, coding, or regulatory advice in connection with your use of the product or service. Please consult your professional advisors for any such advice. Operation of GE Healthcare products should neither circumvent nor take precedence over required patient care, including human intervention of healthcare providers. GE Healthcare products and services do not code medical procedures. Accurate coding is the responsibility of the provider or billing professional.

GE, the GE Monogram, imagination at work, Centricity and AutoBone are trademarks of General Electric company.



General Electric Company, by and through its GE Healthcare division.

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
General Electric Company, prin intermediul diviziei sale GE Healthcare.

©2019 General Electric Company

[Scrieți aici]

GE imagination at work

SECRET DE AFACERI

CONFIDENTIAL

VessellQ™ Xpress with AutoBone™

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Overview

VessellQ™ Xpress with Autobone is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

What's new

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

SECRET DE AFAGEPI

CONFIDENTIAL

Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- Fast tracking detects entire length of vessel and displays the centerline based on a single start point.
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- Thrombus detection with measurements and comparison over time.
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.
- Save State lets you save and

restore the state of the processed images at any stage

System Requirements

VesselIQ Xpress with AutoBone option can be installed on GE's AW workstation or AW Server with VolumeShare 5, VolumeShare 7 or later software.

Compatible hardware VolumeShare 5:

- 8400, 8600, Z800 workstation
- AW Server 2.0

Compatible hardware VolumeShare 7:

- Z800, Z820 (and later)
- AW Server 3.1 and above (recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048))
- Centricity™ Universal Viewer

Indications for Use

VesselIQ Xpress is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. VesselIQ Xpress is a post processing application option for the Advantage Workstation (AW) platform, CT Scanner or PACS stations, which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and

directional vessel tortuosity visualization. VesselIQ Xpress automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel, analyzing calcified and non-calcified plaque to determine the densities of plaque within a coronary artery, measure areas of abnormalities within a vessel.

VesselIQ Xpress is a software post processing option for the Advantage Workstation (AW) platform, which can be used in the analysis of the 3D Angiography data. It provides a number of display, measurement, and batch filming/archive features and will aid physicians in studying user selected vessels for stenosis analysis, pre/post stent planning, and directional vessel tortuosity visualization.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

Regulatory Compliance

This product complies with the following requirements:

- European Council Directive 93/42/EEC concerning medical devices.

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, imagination at work, AutoBone, VesselIQ are trademarks of General Electric Company.
All other third party trademarks are the property of their respective owners
Dec1471990.

2


SECRET DE AFACERI!

CONFIDENTIAL

Features Detail

Anatomy-based protocols

- The program offers you predefined protocols for fast, easy review of all vascular structures:
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

AutoBone™

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail

- Automatic calcium removal for vessel lumen visualization

Auto Aorta

- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.


- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.

Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.
- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
 - Minimum diameter
 - Maximum diameter
 - Mean diameter
 - C...

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, Imagination at work, AutoBone, VesselIQ are trademarks of General Electric Company.
All other third party trademarks are the property of their respective owners
Doc1471890.


Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.

Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.


Advanced Vessel Analysis Tools

- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for

quick review in curved, oblique longitudinal or cross-sectional views.


- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis. Available measurements include:
 - Distance
 - Cross-sectional area
 - Single or dual reference point comparison
 - Relative percent stenosis / aneurysm
 - Color Identification
 - 3D Angle
 - Thrombus Extraction
 - Volume
 - Mean, min, max diameter

- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.

- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

GE imagination at work

© 2015 General Electric Company.
 All rights reserved. Data subject to change.
 GE, GE Monogram, imagination at work, AutoBone, VesselIQ are trademarks of General Electric Company.
 All other third party trademarks are the property of their respective owners
 Doc1471990.

SECRET DE AFACER

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one or two point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.


Summary Table & Report

- Measurements are automatically captured in summary table for easy review.
- Simply click on any measurement in the summary table and you will be taken to that measurement location.

- Prepare report with a single click or

- generate an advanced report with customized options
- Send the report to database, filmer or PDF.

Multiple Study Comparison

- With VesselIQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click.

Save state captures:

- 3D images
- Vessel tracking
- Measurements
- Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VesselIQ Xpress provides user-friendly tools an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

GE imagination at work

© 2015 General Electric Company. All rights reserved. Data subject to change. GE, GE Monogram, imagination at work, AutoBone, VesselIQ are trademarks of General Electric Company. All other third party trademarks are the property of their respective owners. Doc1471990.

System Requirements

- AW Server version 3.1 or later
- AWS supports single or dual monitor configurations.
- AW Standalone version 4.7 or later

Indications for Use

VessellQ Xpress is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

VessellQ Xpress is a post processing application option for the Advantage Workstation (AW) platform, CT Scanner or PACS stations, which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and directional vessel tortuosity visualization.

VessellQ Xpress automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel,

analyzing calcified and non-calcified plaque to determine the densities of plaque within a coronary artery, measure areas of abnormalities within a vessel.

VessellQ Xpress is a software post processing option for the Advantage Workstation (AW) platform, which can be used in the analysis of the 3D Angiography data. It provides a number of display, measurement, and batch filming/archive features and will aid physicians in studying user selected vessels for stenosis analysis, pre/post stent planning, and directional vessel tortuosity visualization.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

Regulatory Compliance

This product complies with the following requirements:

- European Council Directive 93/42/EEC concerning medical devices.

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, Imagination at work, AutoBone, VessellQ are trademarks of General Electric Company.
All other third party trademarks are the property of their respective owners.
Doc1471990.

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
CONFIDENTIAL

Lung VCAR

For lung nodule analysis - automatic visualization, measurement, reporting and follow-up

CT Scanning often is utilized as the exam of choice for intricate CT exams to visualize and analyze complex lung pathology. The detection of pulmonary nodules and assessment of their evolution with CT are of major importance in chest imaging.

What's new

- Synchronized 2D, Digital Contrast Agent (DCA) and segmentation analysis.
- Automatic nodule visualization.
- Automatic nodule analysis (volume, doubling time, % growth).
- Automatic follow-up
- Customizable and interactive report

Overview

Lung VCAR brings efficient CT pulmonary nodule assessment and diagnosis. The innovative Digital Contrast Agent (DCA) feature automatically visualizes lung nodules to help you confirm the presence or absence of suspicious lesions from 2 to 12 mm in size. Lung VCAR allows automated follow-up for lesion matching by the registration of two or more datasets, automatic lesion classification, and a customizable reporting tool.

Visit us:

www.gehealthcare.com/awi/applications/lung-vcar/

SECRET DE AFACERI

CONFIDENTIAL

Features

Review

- Ability to synchronize multiple images for nodule comparison.
- Ability to review single or multiple exams and compare axial, sagittal, oblique, coronal, and volume-rendered images.
- Automatically propagates previous exam bookmarks to current exam.
- Automatically segments both right and left lung to reduce visual distraction.
- Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

Analysis

- Performs automatic segmentation of all nodule types.
- Provides automatic nodule analysis, including % growth, doubling time
- Automatic Nodule Contour to verify pixels within the volume

Reporting

- Features a customizable, interactive patient reporting tool

System requirements

- AW Workstation
- AW Server

Indications for Use

Lung VCAR/ AdvantageALA is an image analysis software package for AW systems, which allows the user to study suspicious lesions within the lungs using CT helical- and axial-acquired images.

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GE imagination at work

© 2012 General Electric Company
All rights reserved. Data subject to change
GE and GE Monogram are trademarks of General Electric Company
* Trademark of General Electric Company

2

SECRET DE AFACERI
CONFIDENTIAL

Lung VCAR

Features

Scanning

Lung VCAR sensitivity and segmentation performance is a function of slice thickness and dose. Lung VCAR was designed to be compatible with slice thicknesses ranging from 0.625 -1.25 mm with a dose range of >40.

Review

Correlated Workflow


The Lung VCAR user interface is designed to organize routine chest radiology reading and make it more efficient. Viewports and displayed anatomy are synchronized and correlated, facilitating reading by making review and problem solving more transparent.

Lung VCAR gives you the ability to:

- Synchronize multiple images for nodule comparison.
- Review single or multiple exams and compare axial, sagittal, coronal, oblique, and volume rendered images.
- Automatically propagate and display previous exam bookmarks on the current exam when two exams are loaded.


Lung VCAR provides you with tools to help improve your workflow and analysis precision:

Lung Segmentation - Lung VCAR automatically segments both the right and left lungs. Segmentation reduces visual distractions associated with irrelevant anatomy during lung nodule evaluation. Total lung volume is captured in the patient report.


Automatic Follow-up synchronizes and displays current and prior patient exams and their corresponding bookmarks for quick and easy comparison review.

Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

DCA Color Mapping highlights nodules in your choice of red or, for a colorblind-friendly display, orange.


DCA Object Size (0-10mm) gives you full control over visualization sensitivity to a desired nodule size and the number of false positives displayed.

DCA Shape Choice (Sphere / Cylindrical) lets you select from two different filters to help visualize specific shapes and assist in problem solving:

- Cylindrical for vessels
- Spheres for nodules

DCA Toggle enables you to toggle the DCA-highlighted nodules on and off.

X-Reference / Correlation Bar provides you with a quick reference for localizing a nodule's global location when reading 2D axial images. The x-reference bar is synchronized with both review and analysis layouts for immediate screen refresh when deposited at any location or on any anatomy.

Review Controller Bookmark Management

Lung VCAR provides you with custom-designed (select dot or small box bookmark format) bookmark management capability linked to review, paging and Digital Contrast Agent for a more organized nodule tracking and analysis workflow.

GE imagination at work


Lung VCAR -1

SECRET DE AFACERI


CONFIDENTIAL

Review Controller Interface lets you dynamically and easily view, page and deposit bookmarks within an exam. It provides:

- On/Off toggle capability for DCA-highlighted nodules
- Smooth, contiguous paging through large data sets
- Thick slab MIP paging
- Bookmark paging (next / prior and add / delete)

Automatic Bookmark Propagation lets you automatically propagate bookmarks from previous to current exam or current to previous exam.

Deferred bookmarking lets you queue bookmarked nodules for later segmentation (operator initiated segmentation).


Immediate bookmarking automatically segments the nodule and displays the analysis review layouts immediately after you place a bookmark.

Layouts

Layout Preferences let you create and save a layout to match your reading style. Layout elements such as viewport, image type in each viewport, slice thickness, window level/width, magnification, rendering mode, can be easily customized.


Analysis

Analysis mode allows you to maximize use of a second monitor by displaying a segmented nodule of interest to an orientation useful for quick problem solving and analysis.

Automatic segmentation of all nodule types provides automated segmentation and sizing of all nodules types. This feature can:

- Display 3D volume rendered nodule view.
- On follow-up exams, display doubling time and percentage of nodule growth.
- Compare nodules from multiple exams.

Shutter view provides you with a quick Volume Rendered framed view of a point of interest in relationship to surrounding structures for fast visualization and problem solving.

Manual or automatic nodule characterization can automatically characterize and label nodule type and apply the necessary segmentation algorithm for accurate isolation. You can at any time manually override this feature and apply a different segmentation algorithm and labeling for any nodule. **Automatic Nodule Analysis** provides the following nodule information:

- % growth
- Doubling time
- Volumes (solid component, non-solid component, total volume)
- Maximum distance in 3D (long axis)
- 2D measurements
- Hounsfield units (Min, Max, Avg)
- Scan interval

GE imagination at work

Lung VCAR -2

Paging enhancements provide a special scrolling tool that lets you page through large data sets while the image displays seamlessly.

Temporal statistics display versatility lets you display statistics for a single exam, or two or more exams for an efficient review of data so you can make a fast, informed decision.

Nodule ID	Date	Size	Status
1	12/11/11	1.2	Stable
2	12/11/11	1.5	Stable
3	12/11/11	1.8	Stable
4	12/11/11	2.1	Stable
5	12/11/11	2.4	Stable
6	12/11/11	2.7	Stable
7	12/11/11	3.0	Stable
8	12/11/11	3.3	Stable
9	12/11/11	3.6	Stable
10	12/11/11	3.9	Stable

Automatic nodule contour capability lets you verify pixels within the volume. It displays segmentation nodule contour on the 2D image for quick segmentation validation.


Reporting

Patient Reporting

Lung VCAR provides customizable and interactive report capabilities giving you more reporting structure and flexibility. You can specify and include:

- Number of images per sheet
- Facility logo
- Physician signature
- Patient demographics
- Comments / impressions

Reports can be exported to CD, HTTP, or you can choose a printed paper output.

- Patient and physician reports are customizable.
- Layout and entered information is retained for use with subsequent reports.

- A patient questionnaire/survey that can be printed out and given to patients to gain a better understanding of patient history and indications of factors that may impact therapy management is available.
- Elements that you can enter into a report include: no findings (default); volume measurements from each exam; screen capture of the volume rendered nodule; axial images.
- You can generate reports as PDF files, store them as DICOM objects, or print them to a standard postscript printer.
- You can indicate negative findings with just a single click.
- You can archive and restore full nodules histories by saving as DICOM objects.

Interactive Reporting

- Temporal Statistics defined by Automated Nodule Analysis are displayed in a follow-up table that gives you a quick visual comparison between previous and current exams.
- You can select and propagate desired images to the patient report.
- You can select and propagate Screen Captures, complete with global annotation to the patient report.

Summary

With automated nodule visualization and localization, increased detection sensitivity, accurate and repeatable segmentations analysis, Lung VCAR could streamline your reading workflow, enhance your efficiency, and increase your diagnostic confidence.

Instructions for use

Lung VCAR/AdvantageALA is intended to provide an optimized non-invasive application to measure abnormalities in the lung (for example, nodules, lesions, etc.) from a set of computed tomography (CT) images.

The software is designed to support the physician in confirming the presence or absence of physician identified lung lesions (e.g. nodules). The software allows measurement of volume over time using a consistent standardized measurement protocol, thus providing an estimation of the volume doubling time. The Lung VCAR/AdvantageALA software allows analysis and displays statistics for nodule characterization all the different nodule types.

Lung VCAR Digital Contrast Agent (DCA) module is an automated highlight feature for the visual identification of possible lesions. Digital Contrast

GE imagination at work

AdvantageSim MD-3

Agent (DCA) is a 3D filter that produces images that highlight spherical (S) or cylindrical (C) anatomical regions, such as nodules, cysts, scars, and vessels. Images are made available to the physician to aid in characterization of suspicious nodules and thus, the patient management care decision process. Lung VCAR/AdvantageALA provides the physician with additional information, meant to complement diagnosis based on classical techniques. Lung VCAR/ AdvantageALA is an image analysis software package for AW systems, which allows the user to study suspicious lesions within the lungs using CT helical – and axial-acquired images.

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
CONFIDENTIAL

Colon VCAR EC

A comprehensive reading workflow solution for colonic lesion detection.


Colorectal cancer is a major health concern for both men and women. Recent studies conclude that early detection of colonic polyps can reduce the mortality rate of colon cancer. You need a protocol that lets you quickly and accurately detect and report results of CT colon scans; a program that optimizes your workflow and improves your reading productivity.

Overview

Colon VCAR EC gives you a comprehensive reading workflow solution for detecting colonic lesions with flexibility and efficient performance. With it you can visualize anatomy that would otherwise be obscured by tagged fluid. Plus digital contrast agent filter highlights potentially suspicious regions.

What's new

- The program allows you to read and problem solve in 2D, 3D, or 360 degree dissection views.
- Prone/Supine polyp linking provides quick navigation between datasets.
- Digital Contrast Agent (DCA) provides automatic visualization of shapes characteristic of polyps.
- Electronic Cleansing (EC) subtracts tagged stool and fluid, making it easier to identify lesions.

Visit us:

www.gehealthcare.com/aw/applications/


SECRET DE AFACERI

CONFIDENTIAL

Features

- Prone and supine synchronized reviewing for problem solving.
- You synchronize prone and supine navigation based on polyp bookmark location and distance from center line.
- Provides you with quick, easy to use polyp and center line measurement tools for complete reporting.
- You can leverage any combination of synchronized 2D, 3D, and 360 views.
- Dissection views make primary reading and problem solving more efficient.

- You can vary navigation speed.
- Lets you track entire colon from rectum to cecum for thorough, precise exams.
- You are provided with unrestricted viewing angles of all colon surfaces.
- You can apply color to bookmarked lesions for better visualization.
- The program allows you to assess luminal, mural and extracolonic anatomy and pelvic structures.

- A selection of application-specific tools ensures ease of use and thorough exams.
- Your workflow benefits from fast, interactive electronic patient reporting.

System Requirements

- AW Workstation
- AW Server

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability.

It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

This product has not been approved for sale in the United States by the United States Food and Drug Administration (FDA).

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

Features Detail

Electronic Cleansing


Electronic Cleansing, accessible in both 2D and 3D review, subtracts tagged stool and fluid from the generated data making it easier to identify lesions that may be visible in only one view.

Tagged material is identified based on Hounsfield units. Processing tasks such as Auto Segmentation, 3D flythrough, and Dissection also take advantage of tagged material subtraction.

Auto Segmentation shows the tagged material while 3D flythrough enables quick anatomy visualization without the tagged material. With this method, the tagged material is removed yet shown in a different color at an adjustable opacity level.

Similarly, the tagged material can be removed from the dissection view as seen in the images below.

User Interface

The Colon VCAR EC user interface has been optimized to facilitate fast, accurate examination of the colon's full extent. It provides seamless integration between automatic and manual interactions for a level of flexibility designed to reduce reading times.

Simultaneous, Indexed Viewing

The program can display both prone and supine

views, which can be manually synchronized for simultaneous and indexed viewing. During review, you can "Lock to Target" with a single button click to quickly problem solve during the primary read. And you are able to bookmark locations of suspected polyps.

360-degree Dissection

This powerful tool provides a 360-degree view of the entire internal large bowel simultaneously displayed in supine and prone positions. This ability speeds viewing and analysis.

This feature lets you leverage all available views: 2D, 3D, and 360-degree dissection during the primary read. For more traditional reading methods, a 2D, 3D, or Dissection view can be laid out and reserved for primary reading with reference to other views for problem solving.


Joystick Navigation

Colon VCAR EC lets you adjust the speed and direction of a flythrough simply by using the mouse. Flythroughs can also be set on autopilot for hands-off reviewing. Additionally, a single mouse click lets you lock onto a region of interest and fly around the target anomaly.

Computer Aided Reading


Colon VCAR EC's Contrast Agent (DCA), a 3D filter, highlights spherical anatomy such as polyps in color. The resulting images facilitate the characterization of potential polyps.

Virtual Biopsy

Problem solving while 3D reading is facilitated by the Virtual biopsy tool, which allows you to view areas in the navigator in color-coded or black/white scales. Different densities, determined by Hounsfield units, can be displayed using the color scale.

Auto-Segmentation


Colon VCAR EC quickly performs the time-intensive task of colon segmentation for both prone and supine data sets. All you do is select the protocol.

The program automatically defines a luminal centerline in both prone and supine positions. This can be applied to a high percentage of cases, depending on the level of insufflation and collapse.

Semi-automatic and manual segmentation modes are also available. Semi-automatic segmentation provides complementary, and often necessary, assistance in cases when the colon is collapsed and operator intervention is required. This mode is typically used to connect colon segments around a collapse to form a single centerline. Manual segmentation is typically used to extend the centerline to the end of the cecum for completeness. Manual editing is also available to refine automatic tracking if needed.

Dissection (centerline) Contouring Verification

Using a color mapping technique, Dissection Contouring Verification outlines colon anatomy included in the centerline trace. You can review verification in any orthogonal plane.

Automatic High Resolution Flythrough

Real time 512-matrix resolution automatic flythrough enhances navigational image quality, and provides operational flexibility, functionality, and simplicity. Anatomy remains sharp with minimal pixilation. Steering remains centered as you navigate through the colon. You can specify your navigation speed, and 3D navigation is correlated with 2D and 360-degree Dissection views for rapid synchronized referencing in primary or secondary reading.

512 matrix resolution

Save State

This feature gives you the ability to save any centerline edits, bookmarks, and colored polyps for later reloading. It is also convenient for quick recall during consultation or exam review.

Prone/Supine Polyp Linking


After bookmarking a polyp in prone and supine datasets, the program provides you with the capability to match and label them for easy tracking. This gives you the ability to read registered prone/supine datasets for a one pass read.

Synchronized Prone/Supine 3D Dissection View Navigation

Once prone/supine polyp linkage is made, the program co-registers and links the datasets and polyps, letting you visualize opposing anatomy for a synchronized review.

Small Bowel Extraction

This feature quickly segments the small bowel, giving you an unobstructed view of the large intestine.

Analysis Mode

The program's Analysis Mode provides a variety of tools to help streamline your workflow:

- Lock Mode lets you target and possible lesion and stay locked on it at all times.
- Synchronized display can be customized to show two oblique images, axial, sagittal, coronal, volume rendered and/or 3D views within the colon in a Lock mode. Leveraging multiple synchronized planes and views in a problem-solving mode may enhance your reading productivity.

GE imagination at work

Colon VCAR EC -2

Color Centerline Tracking

With Colon VCAR EC you can display a real-time 3D model of the colon with a defined centerline to quickly orient the location of the area you are viewing to the overall colon. This synchronization is always linked and available for quick reverence.

Patient Report*

Customizable, interactive reporting capability gives greater structure and flexibility in generating your reports. You can choose how many images to display per sheet, place your facility's logo, apply physician signatures, list patient demographics, and detail your comments and impressions. Export your report to CD, HTTP, or print it to paper.

The Quick Report panel gives you instant access to report input fields and provides access to automatic bookmark annotation on the colon map.

Measurement Tools


With Colon VCAR EC, you can measure polyp distance from distal rectum to polyp location in millimeters, perform 2D measurements, and obtain statistics of ROI in Hounsfield units.

Archive

You can define the number of images per polyp you want to archive to either a networked location or to magnetic optical disk.


Customizable Review and Analysis Layouts

Create layouts that match your reading style and save them as protocols for future recall. Viewports can be shaped into a variety of square or rectangular configurations.

Five of the most popular protocol layouts are pre-built for quick loading. Choose from several dissection configurations:


- 360-dissection gives you convenience and productivity benefits.
- 120-dissection lets you read around distortions if necessary.

- Use the 2D read configuration if the exam has poor insufflation and multiple collapses. 2D can also be used for extracolonic review.

- Use the 3-series layout to load a prone, supine, and decubitus series in parallel to obtain a correlated view.

- You can also customize new layouts to suit your own reading style.
- If you work on a one-monitor Advantage Workstation, a single monitor display is also available.

Movie Generator

You can generate movie clips of the area of interest simply by defining the start and end locations with a mouse click. A cross-sectional movie will be generated. You can also generate lock-to-target movies. Movies can be saved as mpeg, avi, HTTP and saved to CD.

GE imagination at work

Colon VCAR EC -3

SECRET DE AFACERI
CONFIDENTIAL

Entire Colon Movie

Generate a fly-through movie of the entire colon with one button click. Export the result as an mpeg or as screen captures to a PACS workstation for review or consultation.

Additional Features

- Automatic center line tracking
- SmartCursor™ for easier navigation
- Synchronized reformatted views
- Multiple correlations of anatomy with CT colonography
- Exploration of suspicious areas in primary 2D reformatted sections
- Primary 3D reformatted sections. Obliques, Navigator, Quick Virtual Dissection, and axial views are all correlated to the same anatomy of interest.
- Black/white or color display

Summary

Colon VCAR EC helps make colon evaluation a less traumatic procedure for patients. Preparation for the exam is less severe; the procedure is less time intensive and less expensive than optical colonoscopy. Moreover, the tools and features this program puts into your hands will streamline your workflow and potentially enable greater diagnostic confidence.

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability. It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.


*Report feature is only available on AW Workstation

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

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
CONFIDENTIAL

CT Perfusion 4D

Comprehensive perfusion analysis available for neurology, oncology and cardiology. Easy-to-use, protocol-driven workflow for all organs.

Perfusion 4D has been expanded beyond stroke and oncology to now include dynamic myocardial perfusion. For any organ where blood perfusion is of interest Perfusion 4D is a complete package to quantitatively provide the necessary information to guide treatment decisions.

The protocol-driven design leads the user step-by-step through the process, reducing keystrokes and improving repeatability. Get the information you want quickly and reliably.

Overview

GE's Perfusion 4D is a fast, easy-to-use automated software for analyzing CT Perfusion images related to stroke, tumor angiogenesis and dynamic myocardial perfusion. Its simple user interface and automated perfusion post-processing make it easy for you to diagnose quickly and accurately - and make treatment decisions more confidently.

GE leverages its deconvolution expertise, with its innovative delay correction algorithm. Perfusion 4D takes perfusion assessment to the next level, with the addition of CT Dynamic Myocardial Perfusion. Providing the classic color overlays with quantitative measurements to evaluate myocardial blood perfusion.

What's new

- Dynamic myocardial perfusion maps
- Improved IQ of functional maps in presence of noise.
- Improved processing times
- Improved ROI mirroring for brain protocols
- Streamlined workflow for tissue classification.
- Permits injection rates of 4cc/sec.
- Incorporates the GE delay-corrected deconvolution method.
- Dynamic registration

Visit us:
https://www3.gehealthcare.com/en/products/categories/advanced_visualization/applications/ct_perfusion_4d_multi-organs

SECRET DE AFACERI
CONFIDENTIAL

Features

- Visualize all information in true volumetric form with the ability to employ all the volumetric-based image analysis tools in the AW Volume Viewer.
 - Whole heart perfusion with semi-automatic guided workflow to aid in assessment of myocardial ischemia
 - Smart Maps using the intelligent 4D noise suppression algorithm improves functional map image quality in the presence of noise.
 - Streamlined tissue classification and visualization incorporates thresholding of selected functional maps.
 - Dynamic non-rigid registration protocols for dynamic acquisitions within the body.
 - Productivity enhancements include:
 - Ability to select automated vessel detection or manual vessel detection for Brain Perfusion
 - Faster processing times for functional map computations
 - Volumetric visualization of functional maps
 - Access from anywhere using AW Server.
 - Time stamps on functional maps and Regions of interest (ROIs).
 - Interactive large vessels exclusion.
 - Many basic features, including the ability to:
 - Process double phase acquisition data.
 - Create 4D Regions of interest (ROIs).
- Display Regions of interest (ROI) statistical information.
 - Display averaged time intensity information of the Regions of interest (ROI).
 - Simultaneously review any number of functional maps.
 - Customize protocols
 - Save parameters.
- Includes these CT post-processing protocols:
 - Brain stroke
 - Brain tumor
 - Dynamic Myocardial Perfusion
 - Body tumor
 - Liver perfusion
 - Pancreas perfusion
 - Prostate perfusion
 - Kidney perfusion
 - Soft tissue perfusion
 - Spleen perfusion
 - Bone perfusion
 - Liver Dynamic Registration
 - Cardiac Dynamic Registration
 - Functional maps:
 - Regional cerebral blood volume
 - Regional cerebral blood flow
 - Regional mean transit time
 - Contrast arrival delay
 - Transit time to peak of impulse residue function
 - Capillary permeability surface area for lengthened acquisition protocols
 - Mean Slope of Increase
 - Base image
 - Average image

System Requirements

Advantage Workstation Volume Share

4.7 or higher, AW Server 3.2 or higher, Centricity™ Universal Viewer

- Z820, Z440 (and later) with 24GB or higher of RAM for the AW
- Recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048) for the AWS

Intended Use

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. It supports the analysis of CT Perfusion images, obtained by cine imaging (in the head and body) after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability). The results are displayed in a user-friendly graphic format as parametric images. This software will aid in the assessment of the extent and type of perfusion, blood volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

GE imagination at work

© 2017 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
DOC1996217.

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Visit us:
http://www3.gehealthcare.com/en/products/ctsegofsw/advan_ced_visuizatiilor/applications/ct_perfusion_4d_multi-organe

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GE imagination at work

© 2012 General Electric Company
All rights reserved. Details subject to change.
GE and GE Monogram are trademarks of General Electric Company.
© 2012 General Electric Company.
* Trademark of General Electric Company as rights reserved. See subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

SECRET DE AFACERI

CONFIDENTIAL

SmartScore 4.0

Advanced imaging software that detects, quantifies and scores cardiac calcium plaque burden. Instantly.

Cardiovascular disease remains one of the most common health issues in the world today. As with many conditions, early detection and patient risk assessment are vital to preventing or minimizing long-term negative effects. But many cases of this disease aren't diagnosed until the patient presents with symptoms. Conventional procedures to assess risk can be expensive, time consuming, and uncomfortable. An accessible, patient-friendly risk assessment method could help clinicians devise a regimen for their patients that might lessen the chance for serious cardiac events.


Overview

SmartScore 4.0 is designed to identify the presence of regional and global coronary artery calcification from a CT scan, then measure and score the results. Scores can be calculated using a standard Agatston/Janowitz (AJ) method. When correlated with a patient's personal information, the score can yield an estimation of a patient's risk for coronary artery disease.

What's new

- Non-invasive alternative to conventional assessment procedures.
- Score can be correlated to age group cohort to determine patient's risk per population
- Provides information on coronary artery wall calcium plaque buildup.
- Automatically detects calcium and highlights it in green.
- Streamlines workflow by networking patient demographics from CT scanner directly to SmartScore program.
- Free Hand Trace lets you outline specific ROIs.

Visit us:

www.gehealthcare.com/awi/applications/smart-score-4/

Features

- AJ 130 scoring method uses conventional Agatston/Janowitz technique with a threshold of 130 HU which is adjusted to the appropriate image slice thickness.
- Volume scoring calculates volumes (mm³) of calcified plaque above the 130 HU threshold.
- Mass scoring calculates mass (mg) of calcified plaque above the 130 HU threshold.
- Individual and aggregate scores are computed for each artery type.
- Provides two methods of calcium scoring.
- Customizable parameters can be predefined based on your site's preferences.
- User interface streamlines your workflow and reduces reading time.
- Report options allow you to custom tailor your reports and distribute them in a variety of formats.

System Requirements

- EKG monitor with recording device and x-ray translucent lead cable.
- Gantry hardware upgrade kit for those scanner systems already in operation.

Minimum platform release:

- Software for Advantage Workstation 4.2P or higher.
- AW Server 2.

Recommended Options

Postscript Printers:

- Codonics: NP-1660M
- Kodak: 3600 DMI
- Codonics 1660M, 1660MD or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Kodak Dmi3600
- Quantum GL2101HD with film/thick paper
- Quantum GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

Intended Use

SmartScore is a non-invasive software option that can be used to evaluate calcified plaques in the coronary arteries, which may be a risk factor for coronary artery disease. SmartScore may be used to monitor the progression/regression of calcium in coronary arteries over time, which may aid in the prognosis of cardiac disease.

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
CONFIDENTIAL

CardIQ Function Xpress

Real time cardiac review the instant you're ready to read.

Cardiac disease is one of the leading health concerns world wide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and function supported by quantifiable data.


Overview

Cardiac IQ Function Xpress post-processing software helps you evaluate cardiac function and diagnose cardiovascular disease with a high degree of confidence. Providing accurate and reproducible quantification of left and right ventricular volumes, ejection fractions, and myocardial mass, Cardiac IQ Function Xpress is optimized to assess cardiac function using multiphase, multi-slice cardiac CT images. The CardIQ Function Xpress option provides an easy-to-use and time-effective way for you to perform cardiovascular functional analysis.

What's new

- Automatically pre-processes and loads exams.
- Automatically detects all chambers in all phases.
- Automatically performs ejection fraction and chamber volume analysis with 91% reliability on LV and RV automatic ejection fraction.
- Provides single click myocardial analysis with bull's eye images.
- Performs simultaneous case review in real time.
- Lets you load multiphase data with no load time.
- Lets you review multiple cases in real time.
- Lets you generate reports with imported analytics.

Visit us:

www.gehealthcare.com/aw/applications/cardiq-function-xpress/

SECRET DE AFACERI

CONFIDENTIAL

Features

- Apply the CardIQ Function Xpress application to standard axial CT images acquired on a GE scanner using the cardiac CT SnapShot imaging acquisition option.
- Automatically selects each heart chamber for individual volume analysis.
- Performs behind-the-scene function data process for real time ejection fraction review, volume analysis, and myocardial analysis.
- Enables you to extract, render, and display 3D volumetric models of the endocardium for ejection fraction calculation.
- Automatic detection of the endocardial and epicardial walls for wall motion, wall thickness, wall thickening and myocardial mass analysis. Lets you perform volume analysis of all chambers.
- You can automatically calculate left atrium volume while excluding the pulmonary vein.
- With a single click you can activate visual wall motion with short axis images in basal, mid, and distal orientations along with a two-chamber long axis view.
- You can perform myocardial analysis with wall motion, wall thickness, and mass calculations.
- The application's flexible reporting tool lets you include graphical representations.
- Display on screen a table of key functional parameters for instant visualization.

Intended Use

CardIQ Function Xpress is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

Image Requirements

CardIQ Function Xpress accepts standard CT image sets acquired on qualified GE CT scanners using the appropriate cardiac imaging software. These images must meet the same image requirements as those for the basic Volume Share 2 application.

System Requirements

- AW Volume Share Workstation with Volume Share 2 running 4 8400 hardware or higher (running 64 bit).
- Auto launch and preprocessing are available only on 8400 workstations with 16 GB RAM.
- 2 monitor configuration.
- Color landscape monitor.

Recommended Printers

- Codonics NP-1660M
- Kodak 3600 DMi
- Codonics 1660M, 1660MD, or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Quantum GL2101HD with film/thick paper
- Quantum GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

Regulatory Compliance

This product complies with the European CE Marking regulation for Medical Devices Directive: Directive 93/42/EEC, dated 14 June 1993.


GE imagination at work

© 2012 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.


Overview

CardIQ Function Xpress allows you to accurately and reproducibly quantify myocardial mass, and left and right ventricular volumes and ejection fractions. The software automatically detects left ventricular endocardial and epicardial contours so that you can more easily assess left ventricular (LV) and right ventricular (RV) functional parameters. The program is optimized to perform assessment of cardiac function using GE multiphase, multi-slice cardiac CT images.

The CardIQ Function Xpress menu page consists of multiple image processing tools for single or multiphase cardiac studies. These protocols can be accessed either from the AW Volume Viewer application or directly from the application menu page on the browser.

CardIQ Function Xpress benefits to cardiac patients include:

- Minimally invasive alternative to help diagnose heart function.
- Data is acquired during a coronary CTA exam, so no additional dose or contrast is required.
- Offers a diagnostic option to patient who may not be able to undergo other non-invasive exams.

Feature Details

CardIQ Function Xpress uses single or multiphase image data sets to measure, graphically plot, and tabulate a number of ventricular function and wall motion parameters.

Pre-processing

- The program automatically processes multiphase data for functional analysis.
- You can select configurations to process only the chambers you require.
- For more streamlined reading, ejection fraction, volume analysis, and myocardial analysis are completed before you view the case.

Auto Launch


- Load multiple cases into auto launch so they are ready to review when you are.
- Switch between exams without having to quit the application and reload.
- Exams ready for review display in green.

Automatic Chamber Detection

- The program automatically loads single or multiphase data sets and selects the heart's chambers for volume analysis.
- Analyze individual chambers by automatically contouring any of the heart's four chambers.
- The system automatically detects the left and right ventricles and the left atrium and contours the endocardium for ejection fraction calculation.
- Choose single or multi-phase chamber detection.

Automatic Ejection Fraction

- The system features a protocol that, when chosen, automatically detects the end diastolic and end systolic phases and calculates an accurate ejection fraction of the left ventricle, right ventricle, and left atrium.

- Ejection fraction measurements are available for


the right and left ventricles.

- Stroke volume, end diastolic and end systolic volumes are automatically calculated and shown in a graphical representation for assessment of heart failure.
- You can easily send ejection fractions and other quantification results to your report.
- With dual screens, you can display quantification assessments on the left monitor and a qualitative overview of the full chamber from apex to base on the right monitor.
- Customize layout designs to accommodate your viewing preferences:
 - Base to apex
 - Basal, mid, and distal along with a two-chamber long axis.
- You can change the selection of end diastolic and end systolic phase.
- While contour detection is automatic, scalpel and auto select tools are available with which you can edit the contours.

Display Table

- A display table shows volume and provides a graphical update of end diastolic and end systolic locations.

- Click on available tabs to select different chambers for analysis.
- Left and right ventricle measurements include:
 - Stroke volume
 - Ejection fraction
 - Cardiac output
 - Myocardial mass
 - End systolic volume
 - End diastolic volume
 - Pulmonary vascular resistance
 - Systemic vascular resistance
 - Body surface area
 - Heart rate

- Left atrium measurements include:
 - Min and max volume
 - Ejection fraction
 - Left atrium index
 - Body surface
 - Heart rate
- Right atrium measurements include:
 - Min and max volume
 - Heart rate
- You can film, save, and export the volume graphs.

Myocardial Analysis

- CardIQ Function Xpress automatically detects the endocardial and epicardial walls to analyze wall motion, wall thickness, wall thickening, and myocardial mass analysis.
- You can edit the contours with a freehand trace and click to spline tools.
- The program displays multiple bull's eye graphs to graphically represent these cardiac characteristics:
 - Wall motion
 - Wall thickness
 - Myocardial mass
 - Cardiac output
- Choose from four bull's eye graph overlays to show 17, 20, and 49 segment model, or a coronary artery territory model.
- Manual scoring is available for wall motion and thickening in 17 and 20 segment.
- Cord display is available for visual assessment of wall abnormalities.

Reporting

- The application provides you with automatic report generation capabilities.
- You can import hospital or clinic logo and address to make your reports site specific.
- You can capture 3D volume rendering endocardium images with ejection fraction calculations displayed for inclusion in your report.
- You can export the volume graph of all heart chambers to your report.
- The program automatically fills these report fields:
 - Ejection fraction
 - Stroke volume
 - End diastolic volume
 - End systolic volume
 - Pulmonary vascular resistance
 - Systemic vascular resistance
 - Myocardial mass
 - Cardiac output
- You can create predefined conclusion templates for various types of findings.
 - Report templates are fully customizable, so you

GE imagination at work

CardIQ Function Xpress -2

can include logos, diagrams, signatures, patient history, and messages.

- Generate reports in PDF, HTML, or plain text
- You can combine CardIQ Function Xpress reports with CardIQ Xpress 2.0 Elite reports.

Simplified Workflow

CardIQ Function Xpress simplifies your workflow with:

- Automatic ejection fraction protocol
- Streamlined multiphase beating images for easy visualization of cardiac wall motion.
- Single click myocardial analysis with bull's eye images.
- Simplified reporting.

Summary

CardIQ Function Xpress provides you with the tools you need to accurately analyze, evaluate, and report on cardiac anatomy and function. Its ease of use will streamline your workflow and its accuracy will give you greater diagnostic confidence.

Indication for Use

CardIQ Function Xpress is used in conjunction with CT cardiac images to automatically calculate and display various left ventricular and right ventricular functional parameters such as ejection fraction, end systolic and end diastolic volumes, stroke volumes, wall motion, wall thickening, cardiac output, myocardial mass, and systemic and pulmonary vascular resistance. Volume measurement of each chamber of the heart is also available. With CardIQ Function Xpress atrium volumes may be used to determine volume assessment of atrial disease to including, but not limited to, atrial fibrillation. CardIQ Function Xpress is a CT, non-invasive image analysis software package, which aids in the assessment of cardiac function and in determination of cardiovascular disease diagnosis and management. CardIQ Function Xpress is available on the AW Workstation, CT Scanner, PAC, or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CardIQ Function Xpress is an extension of Volume Viewer 4 and is dedicated to cardiovascular imaging. It provides an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of CT Angiographic images.

GE imagination at work

© 2012 General Electric Company
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

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
CONFIDENTIAL

CardIQ Xpress 2.0 Reveal

Bringing integration and automation to your CT Cardiac post processing needs.

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and tools for thorough evaluation into the extent of disease.

Overview

CardIQ Xpress 2.0 Reveal is an integrated post processing imaging analysis application dedicated to cardiovascular imaging on GE Advantage Workstation (AW) and AW Server (AWS). The CardIQ Xpress 2.0 Reveal software option can be used to display, reformat and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from a single or multiple cardiac phase image data sets. Cardiac motion is a very real challenge that can occur at any heart rate. CardIQ Xpress 2.0 Reveal is designed to work with SnapShot Freeze* images to automatically process and display images generated with reduced motion blur artifact.

Highlights

- Automatically segment coronary tree across phases.
- Automatically tracks and labels coronary arteries.
- Improved centerline editing tools for faster edits.
- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines.
- Plaque ID tool assists in visualizing and quantifying plaque burden
- Relative perfusion highlights and quantifies hypo-dense areas of myocardium

Features

- Pre-processing automatically recognizes cardiac datasets and performs all segmentations as data arrives on the system saving processing time.
- Auto Coronary analysis automatically segments the coronary tree, tracks the centerline and labels the coronary arteries.
- For readers who prefer a 2D read, three pre-defined orientations provide quick access to planes for best visualization of the coronaries.
- Quick AVA allows access to vessel tracking at any time.
- Relative perfusion color codes and quantifies percent of hypo-dense areas of myocardium with four selectable color maps and fusion overlays with the coronary tree.
- Plaque ID provides volume measurements for four distinct Hounsfield ranges to aid with identification and visualization of coronary plaque in axial, curved or MPVR views.
- Automatically display SnapShot Freeze processed images for reduced motion blur¹.
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes.
- Portfolio of tools such as stenosis measurements, IVUS views and preset volume rendering models assist in communication of specific findings back to referring physicians.
- Measure ES and ED for ejection fraction & volume with automatic extraction of the left ventricle.

- Create multiphase beating heart images
- Select oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Display 4D aortic valve and mitral valve views with a single click

Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart, heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis. The user can see the different valve

planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and color codes like IVUS images. No new or additional diagnostic information is added.

CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

System Requirements

CardIQ Xpress 2.0 Reveal option can be installed on GE's AW workstation or AW Server with VolumeShare 5, VolumeShare 7 or later software.

Compatible hardware VolumeShare 5:

- 8400, 8600 workstation 16GB RAM
- 2800 24GB RAM
- AW Server 2.0

Compatible hardware VolumeShare 7:

- 2800, 2820 (and later) 24GB of RAM
- AW Server 3.1 and above (recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048))
- Centricity™ Universal Viewer

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European Council Directive 93/42/EEC concerning medical devices.

¹ SnapShot Freeze option is not compatible with Centricity Universal Viewer.

GE imagination at work

SECRET DE AFACERI
CONFIDENTIAL

Feature Details

Pre-Processing

- Automated processing of single or multiphase data
- Customizable menu to choose which protocols are pre-processed
- Pre-processed protocols include: Auto Coronary Analysis, Labeling of vessels, VR heart, Angiographic View, and Auto EF

Auto Launch

- Multiple exams can be loaded into the auto launch to ready them for review (AW only)
- Color-code display of ready-to-read exam listing (AW only)
- The ability to switch between exams without having to quit out of the application and reload


Automated Coronary Vessel Analysis (CVA)

- CVA allows the users to track, extract, visualize and measure coronary arteries from either single or multiple cardiac phase images
- Auto launch case selection and display of automatically tracked and labeled coronary vessels
- Simple, corrections to vessel branch tracking

- Curved, oblique, longitudinal and cross-sectional reformatted views are automatically generated in the software for quick review
- Current-state tracking points within vessel analysis may be saved for future review and/or manipulation
- Color-coded plaque analysis with volume measurements for evaluating the change in plaque size over time
- Generate and save rotational movies from curved planar reformation, best section and cross-sectional (lumen) views
- Quantitative or qualitative coronary vessel assessment on user-selected vessel segments.
- Vessel measurements including:
 - Distance and volume
 - Cross-sectional area and mean diameter

- Single or Dual reference point comparison
- Relative percent stenosis
- Single-screen filming capabilities with multi-views within the screen to show the entire picture of the vessel
- Automatically display of the coronary vessel tree
- Using pre-set protocols, 2D or 3D coronary vessel tree models can be generate and displayed in an automatic fashion to give users a qualitative overview of coronary vessel structures

- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines
- Improved centerline editing tools for faster edits
- The vessel tree models can be in the form of either 3D volume rendered or Maximum Intensity Projection (MIP)

GE imagination at work

SECRET DE AFACERI
CONFIDENTIAL

2D Review

- Three pre-defined orientation protocols for easy review of the coronary vessels and chambers
- Dual reformat review allows automatic review of the coronaries from an axial image while linking to oblique views in longitudinal and cross-sectional planes
- Reformat review predefines workflow steps to automatically present thick, multi-planer views of cardiac anatomy
- One-click access to quick vessel analysis

Quick/Dynamic AVA

- Single or dual click vessel analysis from any protocol; 3D, reformat or vessel tree
- Real-time tracking of the center line with instant visualization of a vessel for quick inspection
- Ability to extend tracking proximally or distally for full view of the vessel
- Ability to select a point above and below potential lesions for automated vessel analysis

Vessel Data Base

- A vessel name data base for quick access to vessel labels
- Anatomically relevant listings
- Ability to add new vessel names into the data base
- Automated processing of single or multiphase data

Relative Perfusion

- Hypo-dense areas of myocardium can be high-lighted with color to aid in the visualization of ischemic heart disease
- Pre-defined layouts for density enhancement
- A hybrid display view to show the relationship of the vessel tree with the perfusion defect
- Quantification of a hypo-perfused area as it relates to myocardial defects
- Four selectable color maps to help display the hypo-dense areas
- Exportable statistics via one touch

IVUS Views

- Interactive volume rendered images to better define the calcium, lumen and non-calcified plaque related to lumen narrowing or wall abnormalities

- Can be applied to any 2D MIP image to include; cross-section, MPR, best L-section images
- Display cross-sections perpendicular to the vessel centerline to create an IVUS-like image. No new diagnostic information is added

PlaqID


- Customizable color mapping to HU ranges for easy identification of plaque in axial, curved reformat or MPVR images
- Four distinct colors to aid in differentiating vessel lumen, non-calcified plaque and calcified plaque
- Easy modification of the color and transparency of plaque
- Volume, area to track size of plaque over time
- Smooth transition of color from one density to another
- Plaque volume and area automatically captured in summary table.

GE imagination at work

SECRET DE AFACERI

CONFIDENTIAL

GE imagination at work

SECRET DE AFACERI
CONFIDENTIAL

- Load multiple phases from systole to diastole to create a beating heart movie
- Save movies as a DICOM image set or export in JPEG/MPEG format

SnapShot Freeze (SSF)

- Automatically display motion blur-reduced SnapShot Freeze processed images
- Reprocess SnapShot Freeze images for artifacts caused by motion blurring after manual vessel editing
- Motion blurring correction requires the use images produced by a CT using the optional SnapShot Freeze feature

Protocols

- CardIQ Xpress 2.0 Reveal is supplied with a set of pre-defined protocols that are easily adaptable to customized protocols
- The protocols have comprehensive instructions and tools which minimize the need to memorize the procedures or to refer continuously to user documentation.

Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining

treatment paths from a set of Computed Tomography (CT) Angiographic images CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart, heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis.

The user can see the different valve planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and

color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

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GE imagination at work

© 2015 General Electric Company
All rights reserved. Data subject to change
GE, the GE Monogram, Centricity și SnapShot Freeze sunt trademarks of General Electric Company
™ Trademark of General Electric Company
DOC1039477

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GE imagination at work

© 2015 General Electric Company
All rights reserved. Data subject to change
GE, the GE Monogram, Centricity și Snapshot Freeze sunt trademarks of General Electric Company
* Trademark of General Electric Company
DOC1036477

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Copyright © 2016 General Electric Company. General Electric Company își rezervă dreptul de a face modificări în specificații și ale caracteristicilor indicate aici, sau să întrerupă fabricarea produsului desceis în orice moment fără vreo notificare sau obligație.
All rights reserved. Data subject to the
GE, the GE Monogram, Centricity și SnapShot Freeze sunt trademarks of General Electric Corp
® Trademark of General Electric Corp
© 2015 General Electric Corp