MITRA Experience Peace of Mind BLOOD BAG

INSTRUCTION FOR USE

INTRODUCTION: Blood Bag System manufactured by Mitra Industries (P) LTD. are as per the international standards. These are used for collection, preservation, processing, transport, separation and administration of whole blood and its components. Since the entire blood bag system is in a closed loop, all the above can be achieved in sterile condition.

PRESERVATION: Blood Bag System may contain anticoagulant and/or preservative solution depending on the application envisaged as indicated on the label. The Blood Bag system are available in different volume capacity of 100 ml, 250 ml, 350 ml, 450ml and 500ml. in different configuration for various purposes. Each single or multiple sterile blood bag system packed in a transparent over pouch is then packed in aluminium pouches as per the specified unit packing. These aluminium pouches are further packed in corrugated cartons as per the packing specifications.

STORAGE INSTRUCTIONS: Blood Bag System should be stored at temperature between 4 °C to 30 °C. 

PRECAUTIONS:

- Read the instructions for use, Aluminum pouch label and blood bag label carefully before use and follow all instructions enunciated there in.
- Do not use if there is a visible sign of damage or deterioration to the bag or solution.
- Do not use unless the anticoagulant is clear.
- Do not collect blood after 10 days of opening the aluminium pack.
- Do not reuse. ② Any such attempt would lead to deterioration in quality of blood or its components.
- Do not touch the blood taking needle.
- Safe guard against needle stick injury by using Needle Stick Protection Device (NPD), provided with the bag system.
- Keep away from children.
- To be used by suitably qualified and trained skilled personnel only.
- Do not use blood bags after the expiry date indicated on the label.

ADDITIONAL EQUIPMENT: Weighing scales with mixing or Haemo scales, a stripper-cutter-sealer, aluminium rings are required during blood collection. For separation of blood components, a blood bag centrifuge, a plasma expessor and other related equipment may additionally be required.

PROCEDURE:

- Open the aluminium pouch and the transparent over pouch by tearing along the notch and remove the blood bag just before use.
- Identify Blood Bag System unit using appropriate donor identification system. Confirm that all numbered tubing of each Blood bag unit contains own identical segment number.
- Place bag on weighing scale kept at about 60cm. below donor's arm. Adjust the indicator of balance at zero.
- Apply blood pressure cuff in the upper arm and disinfect site of phlebotomy.
- Clamp the donor tubing and make a loose knot in the tubing at about 10 cm. from the bag end of tubing.
- Inflate blood pressure cuff to 60 mm Hg.
- Remove the Temper Evident Needle Protective Cap as described below and accomplish phlebotomy.
 - Hold the Temper Evident Needle Protective Cap with the finger as close to the seal breaking point as comfortable. Keep the hub stationary, twist the Temper Evident Needle Protective Cap in clockwise or anti-clockwise direction to break the seal at the pinch of position.
 - After the seal has been broken remove the Temper Evident Needle Protective Cap of the needle in a parallel manner, being careful not to drag the Temper Evident Needle Protective Cap across the needle point.
- On venepuncture, hold the needle in position with a piece of tape.
- Remove the clamp to permit flow of blood.
- If pre donation sampling bag is included in the system, close white clamp, provided in the donor line & direct of initial blood into PDS bag by Means of Closure, provided above the PDS bag. Close white clamp of PDS bag when required volume of initial blood collected. Open white clamp, provided in the donor tube to continue collection of blood into the mother bag.
- Position the PDS bag with the tube holder assembly down wards. Ensure the air in the PDS bag is at the top & will not enter the blood collection tube.
- Insert blood collection tube firmly into the tube holder, remove from holder when filled up to its designated capacity. Repeat with other tubes for multiple sample collection. The PDS bag may be remove after the donor samples are collected. The tube above PDS bag must be sealed before removing the pouch.
- During blood collection the blood should be mixed well with anticoagulant solution using weighing scale/ haemo scale mixing mechanism. In case instrumental mixing is not available then mixing has to be accomplished manually by gently inverting the blood bag at several intervals. Never squeeze the bag since this may force air into the donor vein. In the first five minutes of collection, proper mixing of blood with anticoagulant plays a critical role and hence it should be strictly followed.

215.00 mm

145.00 mm

14. Ensure constant observation of donor throughout donation.
15. Collect appropriate amount of blood indicated on the primary bag label with the help of blood volume- weight table below.

Blood collection volume	100 ml.	250 ml.	350 ml.	450 ml.	500 ml.
Net weight	106 gm.	264 gm.	370 gm.	475 gm.	528 gm.

16. Clamp the collection tube of the bag artery forceps or plastic clamp after collection of appropriate amount of blood.
17. Deflate the cuff. Place sterile swab on venepuncture site, apply little pressure and with draw the needle from the vein.
18. Ask the donor apply pressure over the swab on venepuncture site by the fingers of the other hand to plug the blood flow.
19. Tighten knot firmly after collection .clamp between knot and needle and cut tube distal to the knot to separate the needle. Where the needle stick protection device is provided, pull the donor collection tube and slide the needle stick protection device over the hub with a gentle 'pull & click' activity which locks the needle ,in to needle stick protection device thereby 'guarding' sharp point and preventing risk of transmission of the blood borne infections through needle stick injuries.
20. Strip blood from collection tube into the bag, mix and allow collection tube to refill.
21. Collect pilot sample of anticoagulated blood for serological typing or microbiological screening in pre-numbered (and having identical number) pilot tubes either by knotting or clamping or RF sealing.
Begin sealing at needle end of the collection tube and work towards the bag. The last segment of the collection tube that is mostly used for cross matching should not be detached from the main bag.
22. Store blood filled units between 2-6 °C immediately after collections. In case platelet rich plasma (PRP) or platelet concentrate (PC) has to be separated, the bag should kept at 22 °c till platelets are separated taking not more than 4 hours for the purpose. Platelets can be stored up to 5 days in TOTM bags at 20-24 °C.
23. Infuse the blood within 21 days if collected in Blood Bag Systems containing CPD and within 35 days if collected in Blood Bag containing CPDA.
24. RBC preservative solution (SAG-M) can extend RBC storage upto 42days.
For further processing and separation of Blood components using our multiple blood bag system, standard component processing techniques should be followed.

CAUTIONS: ⚠

- (a) The use of aseptic technique is strongly recommended.
- (b) Mix Blood and anticoagulant at several intervals during collection.
- (c) Blood bag filled with Blood should be stored between 2-6 °C.
- (d) Quality of anticoagulant solution may be deteriorate if stored at temperature beyond the specified limits.
- (e) Cross-match before transfusion.
- (f) Do not add medication to the Blood except for red cell preservative solutions. If provided in the multiple blood bag system.
- (g) Mix blood thoroughly before use and do not warm administration .
- (h) Use of sterile transfusion set, with filter, is a must.
- (i) Do not dispense without prescription.
- (j) Report any reactions to Mitra Industries (P) Ltd.
- (k) **Medical device containing phthalate.**  **to be use under supervision of registered medical practitioner when treating children or pregnant or nursing women.**

DISPOSAL: Dispose of the needle safety using a needle cutter in a puncture proof container. The rejected /used Blood bag shall be dispose of as per national regulations.

LIMITATIONS: Blood bags are susceptible to damage inspite of utmost care in the design, selection of components, manufacture and rigorous quality testing prior to sale. Improper handling during transportation, storage, etc. may damage the Blood Bags. Complications due to various reasons other than the quality of the Blood Bag may be arise during Blood Transfusion.

Manufactured in India by:
MITRA INDUSTRIES (P) LTD.
 14/4, Delhi-Mathura Road, Faridabad - 121003, Haryana, India.
 Customer Care No.: +91-129-2271068, 2271069
 Toll Free No.: 1800 1020 307 E-mail: customersupport@mipl.in
 Website: www.mitraindustries.com

EC REP European Representative
Advena Ltd.
 Tower Business Centre, 2nd Flr., Tower Street,
 Swatar, BKR 4013, Malta

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- (j) Raportați orice reacție la Mitra Industries (P) Ltd.