A total of 488 potentially cross-reacting specimens were tested, the results are presented in Table 2.

Table 2

Potentially cross- reacting specimens	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2	
		Negative	Positive
Rheumatoid factor	10	10	0
Cytomegalovirus	5	5	0
Epstein Barr Virus	13	11	2
Malaria	85	85	0
Syphilis	52	51	1
Herpes	5	5	0
HBSAg	124	124	0
Anti-HCV	115	114	1
Anti-HBc	15	15	0
Anti-HTLV ¹ / ₂	15	15	0
Anti-HEV	5	5	0
Specimens from STI clinics	10	10	0
Trypanosomiasis	24	24	0
E Coli infected	10	10	0
Total	488	484	4

Sensitivity

A total of 750 HIV-1 and HIV-2 antibody positive specimens were tested by KHB® Diagnostic Kit for HIV (1 + 2) Antibody (Colloidal Gold) V2 (Table 3).

Table 3

Population	No. of specimens tested	KHB® Diagnostic Kit for HIV (1 + 2) Antibody (Colloidal Gold) V2		Sensitivity
		Positive	Negative	
HIV-1 Ab positive specimens	360	360	0	100%
non- B subtypes HIV-1 Ab positive specimens	40	40	0	100%
HIV-2 Ab positive specimens	100	100	0	100%
HIV Ab positive specimens	250	250	0	100%
Total	750	750	0	100%

Performance characteristics in comparison with an agreed reference standard

	Initial(95% CI)	Final(95% CI)
Sensitivity%	100% (99.2%-100%)	100% (99.2%-100%)
Specificity%	100% (99.4% -100%)	100% (99.4% -100%)
Invalid rate%	0%	
Inter- reader variability%	0.09%	

	Reference testing results			
Diagnostic Kit for HIV (1+2) Antibody (Colloida Gold) v2 Final results		HIV+	HIV-	Total
	HIV (1&2)+	460	0	460
	HIV-1+	0	0	0
	HIV-2+	0	0	0
	HIV-	0	657	657
	Total	460	657	1117

- Sensitivity (95% CI): 460/460=100% (99.2% -100%)
- Specificity (95% CI): 657/657=100% (99.4%-100%)
- Indeterminate results: 0/1117=0%
- False positive results: 0/657=0% (see Annex3)
- False positive results: 0/460=0% (see Annex3)



Seroconversion panels

The mean seroconversion index was 0.25 specimens compared to the benchmark assay, a HIV-1/2 antibody detection EIA assay. Thus Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) v2 detected HIV-1/2 antibodies in 8 different seroconversion panels, on average, 0.25 specimen earlier than the benchmark assay.

CE certification ensures good quality. Products are sold at home and aboard and its preformance is widely recognized.

KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2


Features and Benefits

- Superior Raw Material**
Our HIV 1+2 rapid test who utilized superior gp 160 antigens as the major raw material, which is expressed in mammalian cells in our labs. Unlike the antigen expressed in E. Coli, this antigen has great solubility in Water, high sensitivity, and very low cross-reactivity. With this superior ntigen, the sensitivity and specificity of the kit is ensured.
- “Sandwich” Testing Principle**
Based on the double antigen “Sandwich” principle, this kit can easily detect the HIV antibody (both IgG and IgM) within 2-3 Weeks after infection.
- Easy to Use**
Small Sample Size: 40μL
Reads the result within 15-25 minutes
No equipment required
- Suitable for whole blood, serum, and plasma samples

A Brief of Overview

KIT COMPONENTS

Component (product code)	50 tests R-401-50-C-2	50 tests R-401-50-C-3	50 tests KH-R-02	50 tests KH-R-03	50 tests KH-R-02-CE	1 tests R-401-1-C-CE	25 tests R-401-25-C-1-CE	Accessories kit A-GOLD-01
HIV (1+ 2) antibody test cassette	50	50	50	50	50	1	1	N/A
Sample diluent	4ml x 1 bottle	4ml x 2 bottle	4ml x 1 bottle	4ml x 2 bottle	4ml x 1 bottle	0.5ml x 1 bottle	4ml x 1 bottle	N/A
Safety lancet	50	50	N/A	N/A	N/A	1	25	50
Alcohol pad	50	50	N/A	N/A	N/A	1	25	50
Disposable transfer pipette	50	50	N/A	N/A	N/A	1	25	50
Instruction for use	1 pcs	1 pcs	1 pcs	1 pcs	1 pcs	1 pcs	1 pcs	1 pcs


Principle of the Procedure

The Diagnostic Kit for HIV (1+ 2) Antibody (Colloidal Gold) V2 from KHB adopts the solid phase colloidal gold immunochromatographic technology for the qualitative detection of antibodies against HIV-1/HIV-2. The gold HIV-1 gp 160 conjugate and the gold HIV-2 gp36 conjugate are coated to the conjugate pad in advance. The test line (recombinant gp36 and gp41 antigens) and the control line (monoclonal antibody against gp160) are pre-coated on the surface of Nitrocellulose (NC) membrane. When the specimen is added to the sample pad, it migrates through the conjugate pad, reconstitutes and mixes with the colloidal gold-antigen complex and moves on till the end of the strip. A reddish-purple test line will be visible in the strip if there are enough antibodies against HIV-1/HIV-2 in the specimen. If the specific antibodies are absent, or present at a very low level, no test line appears.

Test Procedure

Equilibrate all specimens and the devices to room temperature (18- 28℃) before testing.
Sample diluent from different batches can not be mixed.

- Take out a test cassette from the foil pouch, and place it on a horizontal surface.
- Use the specimen of either serum/plasma or whole blood: Add 40μl of specimen to the specimen area firstly (Figure 1), then slowly add 1 drop of sample diluent vertically to the same area (Figure 2). Please avoid direct contact of the diluent bottle with the specimen area.
- Incubate the cassette at room temperature (18- 28℃) and read the result after 15 minutes but not more than 25 minutes. Reading the test before 15 minutes or after 25 minutes may cause false result.

Results

Positive result. Clear visible reddish-purple bands appear both at the C-line and the T-line of the cassette.

Negative result. A reddish-purple band appears only at the C-line of the cassette.

Invalid result. A reddish-purple band appears only at the T-line of the cassette. Repeat the test and contact KHB® supplier if the C- line remains invisible.

Invalid result. A reddish reddish-purple band appears at neither the C-line nor the T-line of the cassette. Repeat the test of the sample. If the result remains the same, then the result is invalid and contact KHB® supplier.

Precaution

The test is for in vitro diagnostic use only. The test cassette is for single use only. The precautions are included, but not limited to the following :

[Safety Precautions]

Wear gloves during the entire testing process.

- Do not pipette by mouth.
- Do not eat or smoke while handling specimens.
- Clean and disinfect all the areas that may be contaminated by spills of specimens or reagents with appropriate disinfectant.
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials as if they were infectious wastes, in a biohazard container.

[Handling Precautions]

- Do not mix and interchange different specimens.
- The sample diluent is specific for each lot of reagent. Do not use diluents from different lots.
- Do not perform the test under environment which leads to rapid evaporation (e.g. high temperature (>45℃) and low humidity, fast flow).
- Use of hemolytic specimens, rheumatoid factors-contained specimens and lipidemic, icteric specimens may cause invalid results or false results.
- Pooled specimens or specimens other than specified (e.g. saliva, urine) are strictly forbidden to be tested.
- Ensure that the specified volume of specimen is added by either precision pipette or specimen transfer pipette.

SHELF LIFE: 24 Months
STORAGE: 4-30℃

PERFORMANCE CHARACTERISITCS

Specificity

A total of 3400 specimens including serum/plasma and whole blood were tested by KHB® Diagnostic Kit for HIV (1+ 2) Antibody (Colloidal Gold) V2 (Table 1).

Table 1

Population	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+ 2) Antibody (Colloidal Gold) V2		Sensitivity
		Positive	Negative	
Blood donor specimens	2500	2500	0	100%
Whole blood donation specimens	500	500	0	100%
Hospitalized patient specimens	200	200	0	100%
Pregnant women specimens	200	200	0	100%
Total	3400	3400	0	100%

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Product Service

EC Certificate

EC Design-Examination Certificate

Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex IV (4)
(List A)

No. V7 052591 0011 Rev. 00

Manufacturer:

**Shanghai Kehua Bio-engineering
Co., Ltd.**

1189 North Qinzhou Road
200233 Shanghai
PEOPLE'S REPUBLIC OF CHINA

EC-Representative:

DIAnearing Diagnostics Engineering & Research GmbH
Friedrichstrasse 26, 69221 Heidelberg-Dossenheim, GERMANY

Product:

**Reagents for determination of infection markers
for HIV-1 / HIV-2: Anti-HIV-1 / HIV-2 Rapid Tests**

The Certification Body of TÜV SÜD Product Service GmbH declares that a design examination has been carried out on the respective devices in accordance with IVDD Annex IV (4). The design of the devices conforms to the requirements of this Directive. See also notes overleaf.

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Product Service

EC Certificate

EC Design-Examination Certificate

Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex IV (4)
(List A)

No. V7 052591 0011 Rev. 00

Model(s):

Diagnostic Kit for HIV (1+2) antibody
(Colloidal Gold) V2

Facility(ies):

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233 Shanghai, PEOPLE'S
REPUBLIC OF CHINA

Parameters:

Brand Name:

Diagnostic Kit for HIV (1+2) antibody (Colloidal Gold) V2

Cat.no.:

KH-R-02-CE, R-401-1-C-CE, R-401-25-C-1-CE

**DIAGNOSTIC KIT FOR HIV
(1+2) ANTIBODY
(COLLOIDAL GOLD)
V2**

Instructions for use

REF

R-401-50-C-2 (50 Tests)

***KHB®* Shanghai Kehua Bio-engineering Co., Ltd.**

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DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2

REF R-401-50-C-2 (50 Tests)

PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE

Intended Use

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB is an *in vitro*, visually read test for the qualitative determination of antibodies against HIV-1 and HIV-2 in human serum, plasma, venous and capillary whole blood. This product is intended for use as an aid to detect antibodies against HIV (1+2) from individuals suspected to be infected by HIV.

It is intended to be used by laboratory professional and health care workers at point-of-care.

Summary and Explanation of the Test

AIDS (Acquired Immunodeficiency Syndrome) is an immunosuppressive disease characterized by depletion of the helper T-cell, which leaves the infected individual susceptible to opportunistic infections and some malignancies. Viruses that have been detected to cause AIDS are HIV-1, first isolated in 1983⁽¹⁾ and HIV-2, a distinctive but related type isolated in 1985. Presence of the virus elicits the production of the antibodies specific to either HIV-1 or HIV-2, and this test is designed to detect the antibodies to indicate the presence of HIV infection.

Principle of the Procedure

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB adopts the solid phase colloidal gold immunochromatographic technology for the qualitative detection of antibodies against HIV-1/HIV-2. The gold HIV-1 gp160⁽²⁾ conjugate and the gold HIV-2 gp36⁽³⁾ conjugate are coated to the conjugate pad in advance. The test line (recombinant gp36 and gp41 antigens) and the control line (monoclonal antibody against gp160) are pre-coated on the surface of Nitrocellulose (NC) membrane. When the specimen is added to the sample pad, it migrates through the conjugate pad, reconstitutes and mixes with the colloidal gold-antigen complex and moves on till the end of the strip. A reddish-purple test line will be visible in the strip if there are enough antibodies against HIV-1/HIV-2 in the specimen. If the specific antibodies are absent, or present at a very low level, no test line appears^(4, 5).

Materials Provided

- | | |
|---|------------------|
| • HIV (1+2) antibody test cassette with desiccant | 1 × 50 cassettes |
| • Sample diluent | 4ml ×1 bottle |
| • Safety lancet | 1 × 50 pieces |
| • Alcohol pad (70% Isopropyl alcohol) | 1 × 50 pieces |
| • Disposable transfer pipette | 1 × 50 pieces |
| • Instruction for use | 1 piece |

Cassette

Sample Diluent

Materials Required But Not Provided

- Specimen transfer device (graduated)
- Tips (Pipette)
- Protective gloves
- Specimen and test waste container
- Timer

- Collection tube (with or without anticoagulant, see Specimen Collection and Storage Requirements)

Warnings and Precautions

The test is for *in vitro* diagnostic use only. The test cassette is for single use only. The precautions are included, but not limited to the following:

[Safety Precautions]

- Wear gloves during the entire testing process.
- Do not use the accessories if seal is broken or package is damaged.
- All the accessories are for single use. Do not reuse the accessories.
- Do not pipette by mouth.
- Do not eat or smoke while handling specimens.
- Clean and disinfect all the areas that may be contaminated by spills of specimens or reagents with appropriate disinfectant.
- Decontaminate and dispose of all specimens, reagents, accessories and other potentially contaminated materials as if they were infectious wastes, in a biohazard container.

[Handling Precautions]

- The accessories are only used for specimen collection, do not store the specimen in the pipette.
- The assay is validated for use with 40uL of specimen, obtained when the specimen reaches the second graduation of the pipette. Use of other volumes may not produce the correct result.
- Do not mix and interchange different specimens.
- The sample diluent is specific for each lot of reagent. Do not use diluents from different lots.
- Do not perform the test under environment which leads to rapid evaporation (e.g. high temperature (> 45°C) and low humidity, fast air flow).
- Use of hemolytic specimens, rheumatoid factors-contained specimens and lipidemic, icteric specimens may cause invalid results or false results.
- Pooled specimens or specimens other than specified (e.g. saliva, urine) are strictly forbidden to be tested.
- Ensure that the specified volume of specimen is added by either precision pipette or specimen transfer pipette provided within the test kit.

[Warning]

- For *in vitro* diagnosis of the infection of HIV only, this assay will not provide a diagnosis for any other disease.
- The cassette, individually enclosed in a foil pouch and packed 50 tests per kit box, is only for "single use".
- Do not use the test cassette if the pouch is damaged or the seal is broken. Contact the supplier when it happens as this is an important indicator for quality and performance of the IVD.
- The sample diluent in this kit contains Sodium Azide (<0.1%) as the preservative, be cautious while using with the kit. In case of contact with sample diluent to eyes and/or skin, flush affected area with copious amount of water. Do not dispose the sample diluent directly in drains.
- Care should be taken to avoid the contact between the sample area of the cassette and the diluent bottle to avoid the contamination of the diluent.

Storage Requirements

- The test cassettes and the sample diluent must be stored at 4-30°C until expiration date.

- Do not use expired devices.
- Test cassette should be used immediately upon opening the pouch.
- Sample diluent should be used within 3 months of opening.

Specimen types

Human serum, plasma, and whole blood specimens should be collected through venipuncture (or whole blood through fingerstick).

Plasma specimens and venous whole blood specimens can be prepared with EDTA, heparin or sodium citrate as anticoagulant.

Specimen Collection and Storage Requirements

[Whole Blood]

- Collection by venipuncture

Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture.

Venous whole blood specimens can be stored up to 8 hours at 18-24°C or up to 7 days at 2-8°C if not tested at the same day of specimen collection, and should not be frozen.

Whole blood specimens should be brought to room temperature (18-28°C) prior to use.

Haemolytic or turbid whole blood specimens should not be used.

- Collection by fingerstick⁽⁶⁾

Figure a


Figure b


Figure c


Figure d

1. Choose the ring finger of the patient.
2. Open the alcohol pad and clean the fingertip (Figure a). Allow the fingertip dry.
3. Remove the cap of the safety lancet. (Figure b)
4. Place the lancet on the fingertip, press the lancet against the finger and puncture the skin. (Figure c)
5. Discard the lancet into a bio-hazard waste container immediately.
6. Squeeze the fingertip to form a bulb and use the disposable transfer pipette to collect the drop of the blood. (Figure d)
7. Allow the blood to reach the second graduation and the specimen volume will be 40ul. (Figure e) Avoid formation of bubbles during the collection of specimen.
8. Transfer the specimen to the sample well immediately after collection. Avoid formation of bubbles.
9. Discard the disposable transfer pipette into a bio-hazard waste container immediately.

Figure e

[Plasma or Serum]

- Plasma

Draw the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture and lightly mix, then centrifuge blood to get plasma supernatant.

- Serum

Draw the whole blood into the collection tube by venipuncture, leave to settle for 2 hours

for blood coagulation and then centrifuge blood to get serum supernatant.

Serum and plasma specimens may be stored at 2-8°C for up to 7 days from time of draw, store at -18°C or below for long time storage. Several freeze-thaw cycles should be avoided.

Plasma or Serum specimens should be brought to room temperature prior to use.


Plasma or serum specimens with visible bio-contamination should not be used.

Test Procedure


Equilibrate all specimens and the devices to room temperature (18-28°C) before testing.

Sample diluent from different batches cannot be mixed.

1. Take out a test cassette from the foil pouch, and place it on a horizontal surface.
2. Use the specimen of either serum/plasma or whole blood: Add 40µl of specimen to the sample area firstly (Figure 1), then slowly add 1 drop of sample diluent vertically to the same area (Figure 2). Please avoid direct contact of the diluent bottle with the sample area.
3. Incubate the cassette at room temperature (18-28°C) and read the result after 15 minutes but not more than 25 minutes. Reading the test before 15 minutes or after 25 minutes may cause false result.


Interpretation of Results


Positive Result (Figure 3)

Clear visible reddish-purple bands appear both at the control line (C-line) and the test line (T-line) of the cassette. The test result is valid even the band of the control line (C-line) appears much lighter or darker than the band of test line (T-line).


Negative Result (Figure 4)

A reddish-purple band appears only at the control line (C-line) of the cassette (Figure 4).

Invalid Result 1 (Figure 5)

Reddish-purple band appears at neither the control line nor the test line of the cassette (Figure 5). Repeat the test of the specimen. If the result remains the same, then the result is invalid and contact KHB supplier.

**Invalid Result 2 (Figure 6)**

A reddish-purple band appears only at the test line (T-line) of the cassette (Figure 6). Repeat the test and contact KHB supplier if the C-line remains invisible.

Limitations of the Procedure

- The kit is designed to detect antibodies against HIV-1 and HIV-2 in human serum, plasma, and whole blood. Specimens other than the specified ones may not supply accurate results.
- The intensity of the T-line does not necessarily correlate to the titer of antibody in the specimen.
- The presence of the control line (C-line) only means the flow of the conjugate.
- When specimens contain high concentration of antibody to HIV-1 are tested on the device, the control line could be absent. In this case, please perform further analysis according to the instructions of Figure 6 in Interpretation of Results.
- As this product is intended to detect antibodies against HIV (1+2) from individuals, clinical diagnosis on AIDS should not be made only based on the test results.
- A non-reactive result should not exclude the possibility of infection caused by HIV-1 or HIV-2. A false negative result can occur in the following circumstances:
 - Negative result may not have detected recently acquired HIV infection.
 - Low levels of antibody (e.g., early seroconversion samples) are below the detection limit of the test.
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration, in exceptional cases this may lead to observation of false negative results.
 - Specimens are not properly stored.
 - False negative result due to very high concentration of antibody(ies) to HIV (which is known as 'hook effect').

For these reasons, care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- The product is not validated on specimens from infants, children, patient on anti-retroviral treatment.
- The product is not specifically designed for detection of antibody to HIV-1 group O.**

Performance Characteristics

- Specificity**

A total of 3400 specimens including serum/plasma and whole blood were tested by KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 1)

Table 1

Population	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Specificity
		Negative	Positive	
Blood donor	2500	2500	0	100%
Whole blood donation	500	500	0	100%
Hospitalized patient specimens	200	200	0	100%
Pregnant Women	200	200	0	100%
Total	3400	3400		100%

A total of 488 potentially cross-reacting specimens were tested, the results are presented in Table 2.

Table 2

Potentially cross-reacting specimens	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2	
		Negative	Positive
Rheumatoid factor	10	10	0
Cytomegalovirus	5	5	0
Epstein Barr Virus	13	11	2
Malaria	85	85	0
Syphilis	52	51	1
Herpes	5	5	0
HBsAg	124	124	0
Anti-HCV	115	114	1
Anti-HBc	15	15	0
Anti-HTLV1/2	15	15	0
Anti-HEV	5	5	0
Specimens from STI clinics	10	10	0
Trypanosomiasis	24	24	0
E Coli infected	10	10	0
Total	488	484	4

- Sensitivity

A total of 750 HIV-1 and HIV-2 antibody positive specimens were tested by KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 3)

Table 3

Population	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Sensitivity
		Positive	Negative	
HIV-1 Ab positive specimens	360	360	0	100%
non-B subtypes HIV-1 Ab positive specimens	40	40	0	100%
HIV-2 Ab positive specimens	100	100	0	100%
HIV Ab positive specimens	250	250	0	100%
Total	750	750	0	100%

A total of 40 HIV-1 Non-B Subtypes positive specimens (in Table 3) were tested by KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2, the results are presented in Table 4.

Table 4

HIV Subtype	No. of specimens tested	Positive on KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2
A	1	1
A1	1	1

A2	1	1
C	3	3
CRF01_AE	3	3
CRF02_AG	3	3
CRF06_cpx	3	3
CRF36_cpx	3	3
D	3	3
F1	2	2
F2	1	1
G	3	3
H	3	3
J	3	3
K	3	3
Group O	3	3
Recombinant Form (recombination of subtype B, D, CRF19 and CRF21_A2D)	1	1
Total	40	40

A total of 100 paired HIV-1 antibody positive whole blood and plasma specimens were tested by KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 5)

Table 5

Population	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Sensitivity
		Positive	Negative	
HIV-1 Ab positive whole blood	100	100	0	100%
HIV-1 Ab positive plasma*	100	100	0	100%

*The 100 HIV-1 Ab positive plasma are included in the HIV-1 Ab positive plasma in Table 3.

Seroconversion panels

The sensitivity was evaluated on 30 commercially available seroconversion panels. Overall the Diagnostic kit for HIV (1+2) Antibody (colloidal gold) V2 detect 38 out of the 165 seroconversion panel members, which is 13 specimens more than the least sensitive CE marked antibody test and 15 specimens less than the most sensitive CE marked antibody test.

References

1. Gallo, R. C., S. Z. Salahuddin, M. Popovic, G. M. Shearer, M. Kaplan, B. F. Haynes, T. J. Palker, R. Redfield, J. Oleske, B. Safai, and et al. 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* **224**:500-503
2. Viscidi, R., E. Ellerbeck, L. Garrison, K. Midthun, M. L. Clements, B. Clayman, B. Fernie, and G. Smith. 1990. Characterization of serum antibody responses to recombinant HIV-1 gp160 vaccine by enzyme immunoassay. NIAID AIDS Vaccine Clinical Trials Network. *AIDS research and human retroviruses* **6**:1251-1256.
3. Chen, Y. H., Y. Xiao, W. Wu, Q. Wang, G. Luo, and M. P. Dierich. 2000. HIV-2 transmembrane protein gp36 binds to the putative cellular receptor proteins P45 and P62. *Immunobiology* **201**:317-322..

4. **Brandt, J., L. O. Andersson, and J. Porath.** 1975. Covalent attachment of proteins to polysaccharide carriers by means of benzoquinone. *Biochimica et biophysica acta* **386**:196-202.
5. **Osikowicz, G., M. Beggs, P. Brookhart, D. Caplan, S. Ching, P. Eck, J. Gordon, R. Richerson, S. Sampedro, D. Stimpson, and et al.** 1990. One-step chromatographic immunoassay for qualitative determination of choriogonadotropin in urine. *Clinical chemistry* **36**:1586.
6. **Spielberg, F., C. M. Kabeya, R. W. Ryder, N. K. Kifuani, J. Harris, T. R. Bender, W. L. Heyward, and T. C. Quinn.** 1989. Field testing and comparative evaluation of rapid, visually read screening assays for antibody to human immunodeficiency virus. *Lancet* **1**:580-584.

Key to symbols used

	Caution		Temperature limitation (4-30°C)
	Batch code		Manufacturer
	In vitro Diagnostic use Medical device		Consult instructions for use
	Catalogue number or order number		Do not use if package is damaged
	Use by		Sufficient for "n" use
	Do not reuse		

Shanghai Kehua Bio-engineering Co., Ltd.

1189 North Qinzhou Road, 200233, Shanghai, P.R.China

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Certificate

No. Q5 052591 0010 Rev. 02

Holder of Certificate: **Shanghai Kehua Bio-engineering Co., Ltd.**
1189 North Qinzhou Road
200233 Shanghai
PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate: **Design and Development,
Production and Distribution of
In Vitro Diagnostic Reagents using ELISA,
PCR, Clinical Chemistry,
Chemiluminescence and Rapid Test Technologies**

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:Q5 052591 0010 Rev. 02

Report No.: SH2131901

Valid from: 2021-05-20
Valid until: 2024-05-19

Date, 2021-05-05

Christoph Dicks
Head of Certification/Notified Body

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Design and Development, Production and Distribution of In Vitro
Diagnostic Reagents using ELISA, PCR, Clinical Chemistry,
Chemiluminescence and Rapid Test Technologies

Shanghai Kehua Bio-engineering Co., Ltd.
701, Guiping Road, 200233 Shanghai, PEOPLE'S REPUBLIC
OF CHINA

Production of In Vitro Diagnostic Reagents using PCR and
Chemiluminescence Technologies

Shanghai Kehua Bio-engineering Co., Ltd.
F2, Building 2, No.26 West Zhangjing Road, Dongjing Town,
Songjiang District, 201619 Shanghai, PEOPLE'S REPUBLIC
OF CHINA

Warehouse

Shanghai Kehua Bio-engineering Co., Ltd.
Room 106, Room 301, Room 302, Room 303, Room 304, Room
402, Building 2, 99 Yinxi Road, Songjiang District, 201615
Shanghai, PEOPLE'S REPUBLIC OF CHINA

Warehouse