iChemVELOCITY Strips

REF 800-7204

TABLE OF CONTENTS

English	2
Français (FR)	8
Deutsch (DE)	15
Italiano (IT)	22
Español (ES)	28
Português Portugal (PT-PT)	34
Dansk (DA)	41
Svenska (SV)	47
Norsk (NO)	53
Ελληνικά (EL)	59
日本語 (JP)	66
中文 (ZH-CN)	72
Magyar (HU)	78
Polski (PL)	84
Čeština (CZ)	90
Slovenčina (SK)	96
한국어 (KO)	102
Türkçe (TR)	108
Русский (RU)	114
Hrvatski (HR)	121
Български (BG)	127
Română (RO)	134
Србија (SR)	140
Українська (UK)	146
Português Brasil (PT-BR)	152
Tiếng Việt (VI)	159
REFERENCES	165

iChemVELOCITY Strips

REF 800-7204

INTENDED USE

The iChemVELOCITY Strips are manufactured for use only with the iChemVELOCITY and DxU Urine Chemistry Systems and are intended for use exclusively by healthcare professionals. These qualitative or semi-quantitative measurements can be used along with other diagnostic tools in the evaluation of renal, urinary, and metabolic disorders.

PRINCIPLE

Bilirubin: This test is based on the coupling of bilirubin with diazonium salt in an acidic medium. A pinkish tan color proportional to bilirubin concentration is generated.

Urobilinogen: This test is based on the coupling reaction of urobilinogen with a stable diazonium salt in buffer. A pink to red color proportional to the urobilinogen concentration is generated.

Ketone: This test is based on Legal's method^{1,2} in which the test pad contains sodium nitroprusside and glycine in an alkaline medium. A violet color proportional to methylketone is generated.

Ascorbic Acid: This test is based on Tillman's reaction³ in which the presence of ascorbic acid leads to the decolorization of the test pad from gray-blue to orange.

Glucose: This two-step enzymatic reaction uses glucose oxidase, peroxidase⁴ and a chromogen. Glucose oxidase catalyzes the formation of gluconic acid and hydrogen peroxide via the oxidation of glucose. Peroxidase then catalyzes the reaction of hydrogen peroxide with a chromogen via the oxidation of chromogen to colors ranging between green to gray-blue.⁵ Other sugar compounds are not detected.

Protein: This test is based on the "protein error" of pH indicators⁶ on the green color developed from the presence of protein. This dye-binding test is particularly strong with albumin.⁷

Blood: This pseudo-enzymatic test contains organic peroxide and a chromogen. The peroxidase effect of hemoglobin and myoglobin causes a color change to green.⁸

pH: This test contains a mixed indicator which assures a marked change in color between pH 5 and pH 9. Colors range from orange through yellow and green to cyan.

Nitrite: This test is based on modified Griess reaction in which nitrite in the urine reacts with amide to form a diazonium compound.⁹ The subsequent coupling reaction yields a pink color in the presence of nitrite. Some Gram positive and non nitrite-forming bacteria are not detected in this test.

Leukocytes: This enzymatic test pad contains an indoxyl ester and a diazonium salt. Granulocyte esterases react with indoxyl ester and diazonium salt to generate a violet color.

Specific Gravity: Specific gravity is physically determined by refractometry in the Urine Chemistry System and is not analyzed by strip chemistry.

Color: Color is measured directly from scattered light in the sample using proprietary algorithms embedded in the Urine Chemistry System.

WARNING AND PRECAUTIONS

- Each test strip is for single use only. Do not test urine samples with a developed (used) test strip.
- The iChemVELOCITY Strips have been determined to be nonhazardous under the guidelines issued by OSHA 29 CFR 1910. 1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens and used test strips should be exercised accordingly.
- Consult local or state authorities with regard to recommended method of disposal.
- Do not touch test strip pads.
- Discard any discolored strips as they may have deteriorated.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckmancoulter.com/techdocs

REACTIVE INGREDIENTS*

Bilirubin:

2,4 dichlorobenzene diazonium salt 2.9%

Urobilinogen:

Fluorodiazonium tetrafluoroborate 0.6%

Ketone:

Sodium nitroprusside 2.0%

Glycine 70.5%

Ascorbic Acid:

2,6 dichloro-phenolindophenol 0.7%

Glucose:

Glucose oxidase 1.8%

Peroxidase 0.9%

Tetramethylbenzidine dihydrochloride 3.5%

Protein:

Tetra-bromophenol blue 0.2%

Blood:

Cumene hydroperoxide 14.5%

Tetramethylbenzidine dihydrochloride 2.9%

pH:

Bromthymol blue 18.8%

Cresol red 1.2%

Methyl red 1.3%

Nitrite:

Sulfanilamide 3.7%

N-(naphthyl)-ethylenediammonium dihydrochloride 0.3%

Leukocytes:

Indoxylcarbonic acid ester 0.4%

Diazonium salt 0.2%

*Concentrations obtained above are based on reagent composition (w/w) at time of manufacture and may vary within manufacturing tolerances.

STORAGE AND STABILITY

- Store at 2° C to 30° C (36° F to 86° F) under dry conditions.
- iChemVELOCITY Strips are stable for 12 months from the date of manufacture.
- Open expiration: Unused strips that remain in the original capped container are stable until the expiration date printed on the label.
- Do not freeze.
- Protect the strips against light and moisture.
- Recap the container tightly after removal of the strips.
- Do not use test strips after the expiration date.

MATERIALS PROVIDED

- 1 x 100 iChemVELOCITY Urine Chemistry Strips

MATERIALS REQUIRED, BUT NOT PROVIDED

- IRISpec CA/CB/CC Controls

- iChemVELOCITY or DxU Urine Chemistry Analyzer
- Personal Protective Equipment

SPECIMEN COLLECTION AND PREPARATION

Collect fresh urine in a clean, dry and clearly labeled container. Run the test as soon as possible after the collection. If the specimen cannot be tested within one hour after collection, refrigerate the specimen immediately and allow it to return to room temperature prior to testing. A fresh first-morning urine specimen is the most concentrated and is best for the detection of protein, nitrite and bilirubin. It is important to test using a fresh urine specimen for the determination of bilirubin and urobilinogen as these analytes are highly unstable at room temperature and under light exposure. Unpreserved urine specimens at room temperature may undergo pH change and in turn interfere with protein determination.

PROCEDURE

Follow the directions in the system's Instructions for Use. The system will automatically read each analyte and the results are displayed or printed as soon as they are available.

QUALITY CONTROL

Clinical performance of the iChemVELOCITY Strips should be confirmed regularly using IRISpec Urine Chemistry Controls CA/CB/CC. Each individual laboratory should establish its own guidelines for acceptable clinical performance. Follow Local, State and Federal guidelines concerning quality controls.

REPORTABLE RESULTS

Analyte	Report Values
Bilirubin	NEG., 0.5, 1.0, 2.0, 3.0, 4.0, 6.0, ≥8.0 mg/dL
	NEG., 8.5, 17, 34, 50, 70, 100, ≥140 µmol/L
	NEG., +, ++, +++
Urobilinogen	NORM., 2.0, 3.0, 4.0, 6.0, ≥8.0 mg/dL
	NORM., 34, 50, 70, 100, ≥140 µmol/L
	NORM., +, ++, +++
Ketones	NEG., 5, 10, 20, 40, 60, 80, 100, 150, >150 mg/dL
	NEG., 0.5, 1, 2, 4, 6, 8, 10, 15, >15 mmol/L
	NEG., Trace, +, ++, +++, +++++
Ascorbic Acid	NEG., 20, 40 mg/dL
	NEG., +, ++
Glucose	NEG., 50, 100, 200, 500, ≥1000 mg/dL
	NEG., 2.8, 5.6, 11, 28, ≥56 mmol/L
	NEG., Trace, +, ++, +++, +++++
Protein	NEG., 10, 50, 100, 300, ≥600 mg/dL
	NEG., 0.1, 0.5, 1, 3, ≥6 g/L
	NEG., Trace, +, ++, +++, +++++
Blood	NEG., 0.03, 0.1, 0.5, ≥1.0 mg/dL
	NEG., 0.3, 1, 5, ≥10 mg/L
	NEG., Trace, +, ++, +++
pH	5, 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9
Nitrite	NEG., 0.1, 0.2 mg/dL
	NEG., +, ++
Leukocytes	NEG., 25, 75, 250, 500 WBC/µL
	NEG., Trace, +, ++, +++
SG [*]	1.000 – 1.060

*Read by refractometer. Values provided for user convenience.

LIMITATION OF THE PROCEDURE

Note: Diagnostic or therapeutic decisions should not be based on any single result or method.

Bilirubin:

- **Expected Values:** Expected Values: In normal urine, no detectable level of bilirubin should be obtained. Positive results require further investigation.

- **Sensitivity:** 0.5 mg/dL bilirubin
- **Performance Characteristics:** The test is specifically developed for bilirubin. Biliverdin does not react with this test pad.
- **Limitations:** Some urine specimens may contain impurities such as food dyes and therapeutic pigments to produce a yellowish or reddish discoloration of the test pad that may lead to the interference. Elevated concentrations of ascorbic acid and nitrite may inhibit the reaction. Bilirubin is light sensitive and prolonged exposure of urine specimens to light may result in diminished or false negative values. Elevated urobilinogen concentrations may slightly enhance the response of this test pad.

Urobilinogen:

- **Expected Values:** In normal urine, urobilinogen is usually present at concentration up to 1 mg/dL. A result of 2 mg/dL represents the transition from normal to abnormal state and the specimen should be investigated further for possible liver disease and hemolytic disorders.
- **Sensitivity:** 2 mg/dL urobilinogen
- **Performance Characteristics:** The diazonium-based test is more specific for urobilinogen than Ehrlich's reagent-based test. Test strips cannot determine the absence of urobilinogen, which may be significant in biliary obstruction¹⁰.
- **Limitations:** This test is inhibited by elevated concentrations of formaldehyde. Food dyes and medications that have an intrinsic red color in acidic medium such as red beets, azo dyes, phenazopyridine and *p*-aminobenzoic acid may produce false positive results. Prolonged exposure to light may lead to diminished or false negative values.

Ketone:

- **Expected Values:** Detectable amounts of ketone do not appear in the urine of normal specimen. Positive ketone values may result from the following conditions: starvation, dietary imbalance, diabetes mellitus, eclampsia, insulin dosage monitoring, vomiting and other metabolic disorders.
- **Sensitivity:** 5 mg/dL ketone
- **Performance Characteristics:** The test does not measure β -hydroxybutyric acid and is only slightly sensitive to acetone.
- **Limitations:** Elevated concentrations of phenylpyruvic acid may interfere with the test pad and produce a variety of colors. Phthaleins and anthraquinone derivatives exhibit a red color in alkaline medium and this may mask the response. Large amounts of levodopa and medications containing sulfhydryl groups may produce atypical color reactions.¹¹

Ascorbic Acid:

- **Expected Values:** Ascorbic acid is found in various food supplies and dietary supplements. Concentrations of ascorbic acid greater than 20 mg/dL can be expected to cause strong interference with glucose, blood and nitrite.
- **Sensitivity:** 20 mg/dL ascorbic acid
- **Performance Characteristics:** The oxidized form, dehydroascorbic acid, does not react with this test pad.
- **Limitations:** No interferences are reported.

Glucose:

- **Expected Values:** A small amount of glucose (up to 20 mg/dL) may be present in the normal urine. The detectable sensitivity of this test has been adjusted to exclude the minute amounts of glucose. Therefore, any positive response should be further evaluated.
- **Sensitivity:** 20 mg/dL glucose
- **Performance Characteristics:** The test pad does not react with other reducing sugars such as fructose and galactose.
- **Limitations:** Ascorbic acid concentrations of up to 30 mg/dL did not interfere with glucose assay test results (no false negative results). Acetoacetic Acid concentrations of up to 150 mg/dL did not interfere with glucose assay test results (no false negative results). High specific gravity, acidic pH values and gentisic acid may inhibit color formation. Cleaning agents such as hypochlorite and peroxide may lead to false positive results.

Protein:

- **Expected Values:** Normal urine contains very little protein: usually less than 10 mg/dL is excreted. Positive protein results are considered pathological and should be investigated.
- **Sensitivity:** 10 mg/dL protein
- **Performance Characteristics:** The test pad detects primarily albumin. A negative result does not exclude the presence of other protein molecules such as Bence Jones proteins, globulin and mucoproteins.
- **Limitations:** Food dyes such as red beets and therapeutic pigments such as methylene blue and pyridium may mask the coloration of the test pad. Interference may occur with high specific gravity. Interference may also occur with disinfectants, wetting agents and blood substitutes (quaternary ammonium compounds, polyvinylpyrrolidone, chlorohexidine). There is no interference from pH.

Blood:

- **Expected Values:** Normal urine contains no detectable hemoglobin or intact red blood cell. Any positive results should be further evaluated.
- **Sensitivity:** 0.015 mg/dL blood
- **Limitations:** Reducing agents such as ascorbic acid, uric acid, glutathione and gentisic acid may cause false negative results. Preservatives (formalin) and cleaning agents such as hypochlorite may result in false positives. High concentrations of nitrite and specific gravity can delay the reaction. There is no interference from pH for this test.

pH:

- **Expected Values:** The normal pH of urine can vary between pH 5.0 and pH 9.0.
- **Performance Characteristics:** pH values are determined to within 0.5 unit over the range from 5.0 to 9.0.
- **Limitations:** No interferences are known for pH test pad.

Nitrite:

- **Expected Values:** Normal urine contains no detectable nitrite. However, a negative result does not rule out a urinary tract infection.
- **Sensitivity:** 0.06 mg/dL nitrite
- **Performance Characteristics:** This test is specific for nitrite. Results may depend on the ability of bacteria to reduce nitrate to nitrite, the number of bacteria and the retention time of the urine in the bladder.
- **Limitations:** Food dyes and therapeutic pigments such as red beets and pyridium may cause false positive responses. A negative response in the presence of bacteriuria may be caused by the following: non-nitrite producing microorganisms, low nitrate diet, antibiotic therapy, strong diuresis, high levels of ascorbic acid or insufficient urinary retention time in the bladder. There is no interference from pH or specific gravity for this test.

Leukocytes:

- **Expected Values:** Normal urine specimen should not produce a positive result.
- **Sensitivity:** 25 WBC/ μ L
- **Performance Characteristics:** Leukocyte test detects the presence of esterase in the granulocytic white blood cells. The test result is most frequently accompanied by the presence of bacteria that may or may not produce a nitrite positive reaction.
- **Limitations:** False positive results may occur in the presence of preservatives such as formaldehyde and formalin. High concentrations of protein, glucose, cephalexin and gentamicin may diminish the color response. Test results may be positive in the absence of observable cells if the granulocytes have lysed. The test can be negative in the presence of visible leukocytes if they have not lysed and/or are not granulocytes. There is no interference from pH or specific gravity for this test.

ANALYTICAL SENSITIVITY

Analyte	Analytical Sensitivity
Bilirubin	0.5 mg/dL
Urobilinogen	2 mg/dL
Ketones	5 mg/dL

Analyte	Analytical Sensitivity
Ascorbic Acid	20 mg/dL
Glucose	20 mg/dL
Protein	10 mg/dL
Blood	0.015 mg/dL
PH	NA
Nitrite	0.06 mg/dL
Leukocytes	25 WBC/ μ L

TECHNICAL SUPPORT

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at beckmancoulter.com/techdocs.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

REVISION HISTORY

Revision AA, 11/2018

Initial Release

Revision AB, 04/2019

Correct storage temperature error

Revision AC, 09/2020

Correct product name. Replace "Same as Manual." with "No interferences are known for pH test pad." in pH Limitations. Add glucose to substances that may diminish color response and remove glucose from non interfering substances in Leukocyte Limitations. Change SDS and Symbol Glossary web site address to beckmancoulter.com/techdocs.

Revision AD, 08/2021

Updated Proprietary Name

Updated copyright year

iChemVELOCITY Strips (Bandelettes iChemVELOCITY)

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Phiên bản AD, 08/2021

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REFERENCES

1. Legal, E. A.: New Acetone Reaction and Its Applicability for the Examination of Urine. *Chem. Centr.* 14: 652 (1883)
2. Chertack, M. and Sherrick, J.: Evaluation of Nitroprusside Dip Test for Ketone Bodies. *JAMA* 167: 1621 (1958)
3. Roe, J. H.: Chemical Determination of Ascorbic, Dehydroascorbic and Diketogulonic Acids. *Methods of Biochemical Analysis*, Vol 1:115 (1954) ed. By D. Glick, Interscience Publishers, New York
4. Comer, J.: Semiquantitative Specific Test Paper for Glucose in Urine. *Anal. Chem.* 28:1748 (1956)
5. Appel, W., Nurck, C. and Merkle, U.: A Rapid Test for Urinary Glucose with an Ascorbic Acid Zone. *Medical Laboratory* 6: 29 (1979)
6. Sorenson, S.: The Measurement of Hydrogen Ion Concentration and Its Importance for Enzymatic Process. *Biochem. Z.* 21: 131 (1909)
7. Vonderschmitt, D. and Scholer, A.: Teststreifen fur Screening- Untersuchungen zum semiquantitativen Nachweis von Proteinurien. *J. Clin. Chem. Biochem.* 19: 997 (1981)
8. Leonards, J.: Simple Test for Hematuria Compared with Established Test. *JAMA* 179: 807 (1982)
9. Weltmann, O.: Method for the Simple Detection of Urinary Tract Infections. *Wien. Med. Wschr.* 72: 618 (1922)
10. Strasinger, S. K. and De Lorenzo, M. S.: *Urinalysis and Body Fluids*. 4th edition, F. A. Davis Company, Philadelphia (2001)
11. Simerville, J.A., Maxted, W.C. and Pahira, J.J.: Urinalysis: A Comprehensive Review. *A., Fam. Physician* 71: 1153 (2005)

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Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland +(353) (0) 65 683 1100