Registrieren | Artikulieren
Registration | Articulation
Enregistrement de l'occlusion
Montage sur articulateur
Registrazione | Messa in
articolatore
Registro | Articulación

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User Manual

Mode d'emploi

Istruzioni d'uso

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




AMANNGIRRBACH

Deutsch	04 - 05
English	06 - 07
Française	08 - 09
Italiano	10 - 11
Español	12 - 13

AMANNGIRRBACH

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ARTEX®CR

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

ARTEX®CR

SUPPORT IN INCLINED POSITION

The Artex®CR may be positioned on the table top - front raised

OPEN CENTRIC

Lift **(2)** pivoting lever to stop. The centering axes (2 semi-axes) are moved in, thus releasing the centric.

The centric **(3)** is open. Fixation of the condyle ball is suspended. The top part is released.

CLOSE CENTRIC

Move lever **(4)** downward to stop. The two semi-axes are moved out until softly stopped against the condylar centre, fixing the Artex top section in central position.

Centric closed **(5)**. Visible fixation on condylar centre by the semi-axis.

CONDYLAR BOX

Horizontal condyle path inclination **(6)**, Bennett angle, ISS, protrusion and retrusion are individually adjustable on the CR condylar unit. **One glance is enough.** All values are checked and in view.

CONDYLE PATH INCLINATION (HCI)

Release HCI **(7)** set screw using allen key. Adjust required angle by rotating condylar box **(8)**, followed by locking the HCI adjustment with an allen key.

Adjustment of HCI from -20° to +60°.

Read angle of inclination from front edge of rear section (congruent with scale).

Fig.: HCI = +60°.

SETTING OF BENNETT ANGLE

Open centric. Release locking screw of Bennett angle using allen key **(9)**. The individually set Bennett angle (-5° to +30°) may be read from the scale.

Fig.: Bennett angle = 0°.

Rotate condylar unit **(10)** to required scale reading by using the rotary knob, then lock with allen key.

Fig.: Bennett angle = +30°.

IMMEDIATE SIDESHIFT (ISS)

Open centric. Release locking screw of Bennett angle using allen key **(11)**. ISS individually and infinitely adjustable to 1,5 mm per side.

Push condylar unit **(12)** in mesial direction by applying side pressure to the rotary knob until the scale agrees with the required mm-reading. Then lock screw for rotary knob using allen key.

Fig.: ISS = 1 mm

PROTRUSION

Whilst the centricity is open, the Artex top section is protrusively guided over the adjustment claw **(13)** by rotating the blue knurled screw in clockwise direction.

Protrusion **(14)** infinitely adjustable up to 6 mm.

Fig.: Protrusion = 6 mm.

(15) Example of individually programmed CR condylar unit.

Bennett angle = 30°,
HCI = 60°,
Protrusion = 6 mm,
ISS = 1 mm

RETRUSION

(16) Open centric, release locking screw of retrusion insert using allen key.

Guide retrusion insert **(17)** to required individual value, then relock using allen key.

Fig.: Retrusion = 2 mm.

Retrusively guide Artex top section **(18)** (based on a Bennett angle = 0°) via the adjustment claw by rotating the blue knurled screw in counter-clockwise direction.

Fig.: Retrusion = 2 mm

SEQUENCE OF MOVEMENT

Condyle ball in centered position and movement (view into condylar box).

Attention:

Sequences of movement in Arcon articulators must be performed in permanent contact of condyle balls (bottom part of articulator) with the condylar housing (top part of articulator). Only thus they can be reproduced.

Start position:

(19) Bennett angle adjusted to 30°. Condyle ball in centered position = initial position for sequence of movement.

Sequence of movement:

(20) Bennett angle at 30°. Condyle ball guided subject to permanent contact with the condylar box on all three dimensions. This is an imitation of the adjustment claw of the protrusion pointer. The condyle ball may never lose its contact with the fossa, neither in cranial nor mesial direction.

DISTRACTION

Loose screw from frontal. **(21)** Install distraction to the desired value (0-3 mm) and fix it with the screw.

COMPRESSION

(22) Lift frame top by the required compression reading = distraction. Position models into articulator. Lower frame top for compression and retighten screw.

Fig.: Distraction set to 3 mm.

ARBITRARY PIN/CENTRE OF AXIS

Maxillary cranium/axis-based transfer of the maxilla can be effected directly on the articulator by adapting the facebow **(23)**. Both the arbitrary pin (which reproduces the position of the ear opening/bore of the porus buttons) and the axis mean point can be used as a reference.

! WARNINGS AND SAFETY ASPECTS:

Anatomical anterior guidance

Angle of protrusion track = 40°; on lateral movements the angle varies between 30 and 40°, according to the predetermined Bennett angle

General information

Articulators are precision instruments. Therefore, it is important to handle them with care. Please read the following notes:

- Do not place the articulator in a boil-out unit, a pressure pot or an ultrasonic bath.
- For cleaning and care, use only cleaning agents that are suitable for light metals. Regarding the use of the cleaning agent, please follow the instructions of the according manufacturer. Do not use cleaning agents with a pH <5 or >8 since those might damage the anodized surfaces.

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Amann Girrbach AG
Amann Girrbach GmbH