# **OM-MA**

**For use on IMMULITE® 2000 systems**

**SIEMENS**

# IMMULITE® 2000 OM-MA

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of CA125 antigen in serum, as an aid in monitoring the response to therapy for patients with epithelial ovarian cancer, and in detecting residual ovarian cancer in patients who have undergone first-line therapy and would be considered for diagnostic second-look procedures.

Catalog Number: **L2KOP2** (200 tests)

Test Code: **OV** Color: **Aqua**

**Caution:** In the United States, federal law restricts this device to sale by or on the order of a physician.

CA125 antigen levels in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Results reported by the laboratory to the physician must include the identity of the assay used to measure CA125 antigen levels. Values obtained with different assays cannot be used interchangeably. Before changing assays, the laboratory must confirm baseline values for patients being serially monitored.

## Summary and Explanation

The CA125 determinant was originally identified by a monoclonal antibody, selected for reactivity with a cell line from a patient with serous papillary cystadenocarcinoma of the ovary. This antibody was found to react with cell lines derived from epithelial ovarian carcinomas but not with nonmalignant tissues including normal adult and fetal ovary.<sup>1,2</sup>

Although the precise nature of the CA125 determinant remains unclear, there is agreement that the molecule is a high molecular weight (1000 kDa) glycoprotein with a smaller quantity of carbohydrate than mucins.<sup>2</sup> There is some evidence that

more than one form of the CA125 molecule exists.<sup>3,4</sup>

Epithelial neoplasms of the ovary originate from the single layer of cells covering the ovary. These epithelial cells have a high proliferative capacity, repairing the ovarian surface after ovulation. Suppression of ovulation by oral contraceptives, pregnancy and lactation may therefore reduce the risk of ovarian cancer.<sup>5</sup> Ovarian malignancy has been associated with a variety of peptide growth factors, oncogenes and tumor suppressor genes.<sup>5</sup> According to a recent study comprising approximately 10% of the US population, more than 48% of ovarian cancer can be found in women 65 years of age and older, rather than in younger women. The incidence increases with age, reaching a peak of 54 in every 100,000 women in the age group 75 to 79 years.<sup>6</sup>

Measurement of CA125 before and after cytoreductive surgery for ovarian cancer has been shown to predict the likelihood of a patient being left with residual disease.<sup>11</sup>

IMMULITE 2000 OM-MA uses a murine monoclonal antibody for the capture and a rabbit polyclonal antibody for detection of the CA125 antigen. The monoclonal antibody was established by immunization with human mucin prepared from a pool of patients with epithelial ovarian cancer. The monoclonal antibody, which forms the basis of specificity for the kit, recognizes a repetitive protein determinant expressed in the protein core of the CA125 antigen.

This antibody has specificity for an epitope that overlaps with, or is very near to, that bound by the M11 monoclonal antibody. The M11 monoclonal is incorporated in many commercially available CA125 immunoassays.<sup>15</sup> The polyclonal antibody is affinity-purified against CA125 antigen, resulting in a reagent reacting with multiple epitopes on this antigen.

## Principle of the Procedure

IMMULITE 2000 OM-MA is a solid-phase, two-site chemiluminescent immunometric assay.

**Incubation Cycles:** 2 × 30 minutes

## Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 OM-MA has not been tested with all possible variations of tube types.

**Volume Required:** 50  $\mu$ L serum

**Storage:** 1 day at 2–8°C or 2 months at –20°C.<sup>14</sup>

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

Results reported by the laboratory to the physician must include the identity of the assay used to measure CA125 antigen levels. Values obtained with different assays cannot be used interchangeably.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### OM-MA Bead Pack (L2OP12)

With barcode. 200 beads, coated with murine monoclonal anti-CA125 antibody. Stable at 2–8°C until expiration date.  
**L2KOP2:** 1 pack

### OM-MA Reagent Wedge (L2OPA2)

With barcodes. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-CA125 antibody in buffer, with preservative; 6.5 mL of a buffer, with preservative. Stable at 2–8°C until expiration date.  
**L2KOP2:** 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### OM-MA Adjustors (LOPL, LOPH)

Two vials (Low and High), 3 mL each, of CA125 in a nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KOP2:** 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### OM-MA Sample Diluent (L2OPZ)

For the on-board dilution of high samples. 25 mL of concentrated (ready-to-use) processed, CA125-free nonhuman protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test

tube, so that the barcodes can be read by the on-board reader.

**L2OPZ:** 3 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water, test tubes, controls

## Assay Procedure

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

### Recommended Adjustment Interval:

4 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of CA125.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Interpretation of Results

An IMMULITE 2000 OM-MA result of 21 U/mL or greater is considered elevated and a result less than 21 U/mL is considered normal.

Users should be aware that a result greater than or equal to 21 U/mL may be found in a small percentage of healthy individuals and in patients with nonmalignant conditions, such as

pericarditis, cirrhosis, severe hepatic necrosis, endometriosis (Stages II-IV), first trimester pregnancy, and ovarian cysts, or in patients with non-ovarian malignancies, such as uterine carcinoma, hepatoma, pancreatic adenocarcinoma, and lung cancers.

A result below 21 U/mL does not necessarily indicate the absence of residual or recurrent ovarian cancer because some patients with histopathologic evidence of ovarian carcinoma may have CA125 measurements below 21 U/mL.

The unit (U/mL) adopted by the IMMULITE 2000 OM-MA is arbitrary. Since there are no existing standard reference units, users should **not** make quantitative comparisons between different methods of measuring CA125. Results reported by the laboratory to the physician must include the identity of the assay used to measure CA125 antigen levels. Values obtained with different assays cannot be used interchangeably.

## Expected Values

Based on its relationship to IMMULITE OM-MA (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Sixty-four serum samples from female adults (age range: 17–86 years) in apparent good health were included in the two clinical sites. They were tested by the IMMULITE OM-MA assay. The results ranged from 1.9 to 16.3 U/mL, with a median of 4.8 U/mL. The graph shows the distribution of the CA125 values. (See "Female Normals" Graph.)

The same studies also included adult female patients with ovarian cancer, nonmalignant (gastrointestinal, genitourinary and other nonmalignant diseases) and malignant conditions (breast, gastrointestinal, genitourinary and other cancers). The distributions of these patients, as well as male and female healthy individuals included in the studies, are shown below for the ranges of CA125 measurement indicated.

**Healthy Subjects:**

Healthy Subjects	< 21 U/mL	≥ 21 – < 50 U/mL	≥ 50 – < 100 U/mL	≥ 100 U/mL	n
Female	64	0	0	0	64
<50 years	50	0	0	0	50
≥50 years	14	0	0	0	14
Male	28	0	0	0	28

**Patients:**

Patients	< 21 U/mL	≥ 21 – < 50 U/mL	≥ 50 – < 100 U/mL	≥ 100 U/mL	n
Ovarian Cancer	36	10	7	43	96

**Nonmalignant Conditions**

	< 21 U/mL	≥ 21 – < 50 U/mL	≥ 50 – < 100 U/mL	≥ 100 U/mL	n
GI	5	0	1	0	6
GU	7	0	0	1	8
Others	40	4	0	0	44

**Malignant Conditions**

	< 21 U/mL	≥ 21 – < 50 U/mL	≥ 50 – < 100 U/mL	≥ 100 U/mL	n
Breast	11	2	1	0	14
GI	3	2	1	0	6
GU	4	1	2	4	11
Others	9	3	0	6	18

Elevated levels of CA125 may be associated with endometriosis,<sup>12</sup> the first trimester of pregnancy,<sup>13</sup> menstruation, and cancer of the pancreas, stomach, colon or rectum.<sup>2</sup>

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

**Limitation**

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this

assay should always be used in combination with the clinical examination, patient medical history, and other findings.

**Performance Data**

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 500 U/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 1 U/mL

**High-dose Hook Effect:**

None up to 80,000 U/mL

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three CA125 solutions (1025, 1950 and 4080 U/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for CA125. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:**

Biotin Test Level (ng/mL)	OM-MA (CA125) Concentration (U/mL)	
	15.2	189
	Bias (%)	
5.0	7	0
9.0	14	8
19.0	22	19
32.0	20	26
100	40	26
1500	36	32

Specimens that contain biotin at a concentration of 5 ng/mL demonstrate a less than or equal to 10% change in

results. Biotin concentrations greater than this may lead to falsely elevated results for patient samples.

The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.<sup>17</sup> Subjects who take up to 300 mg of biotin per day may have plasma biotin levels as high as 1160 ng/mL.<sup>18</sup>

**Hemolysis:** Presence of hemoglobin in concentrations up to 192 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

**Method Comparison:** The assay was compared to IMMULITE OM-MA (mono/poly, LKOP) on 161 samples. (Concentration range: 4.16 to 364.7 U/mL. See graph.) By linear regression:

(IML 2000) = 0.97 (IML) + 0.30 U/mL  
r = 0.985

Means:  
106 U/mL (IMMULITE 2000)  
109 U/mL (IMMULITE)

## References

1) Bast RC Jr, et al. Reactivity of a monoclonal antibody with ovarian carcinoma. *J Clin Invest* 1981;68:1331-7. 2) Jacobs I, et al. The CA125 tumour-associated antigen: a review of the literature. *Hum Reprod* 1989;4(1):1-12. 3) Hosono MN, et al. Different antigenic nature in apparently healthy women with high serum CA125 levels compared with typical patients with ovarian cancer. *Cancer* 1992;70(12):2851-6. 4) Kaneko Y, et al. Study for the nature of interfering substances responsible for nonspecific reactive phenomena occurred occasionally in EIA determination of CA125. *Rinsho-Byori* 1992;40(9):948-52. 5) Bast RC Jr, et al. Cell growth regulation in epithelial ovarian cancer. *Cancer* 1993;71(4):1597-601. 6) Yancik R. Ovarian cancer. *Cancer* 1993;71(2):517-23. 7) Nagell JR van. Ovarian cancer screening. *Cancer* 1993;71(4):1523-8. 8) Kramer BS et al. A National Cancer Institute sponsored screening

trial for prostatic, lung, colorectal and ovarian cancers. *Cancer* 1993;71(2):589-93. 9) Bast RC Jr, et al. Perspectives on the future of cancer markers. *Clin Chem* 1993;39(11):2444-51. 10) Jacobs I, et al. Multimodal approach to screening for ovarian cancer. *Lancet* 1988;i:268-71. 11) Brand E, et al. The decline of CA125 level after surgery reflects the size of residual ovarian cancer. *Obstet Gynecol* 1993;81(1):29-32. 12) Pittaway DE. CA-125 in women with endometriosis. *Obstet Gynecol Clin North Am* 1989;16(1):237-52. 13) Niloff JM, et al. CA125 antigen levels in obstetric and gynecologic patients. *Obstet Gynecol* 1984;64(5):703-7. 14) Bonfrer JMG, et al. Technical evaluation of three second generation CA125 assays. *Eur J Clin Chem Clin Biochem* 1994;32:201-7. 15) Nustad K, Bast RC, O'Brien TJ, et al. Specificity and affinity of 26 monoclonal antibodies against CA125 antigen: first report from the ISOBM TD-1 Workshop. *Tumor Biol* 1996;17:196-219. 16) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4. Wayne, PA: NCCLS, 1998. 17) Grimsey P, et al. 2017. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *Int. J. Pharmacokinet.* 2(4), 247-256. 18) Piketty ML, et al. 2017. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 55(6):817-825.

## Technical Assistance

In the United States, contact  
Siemens Healthcare Diagnostics  
Technical Services department.  
Tel: 877.229.3711. Outside the United  
States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare  
Diagnostics Products Ltd. is certified to ISO  
13485.

## Tables and Graphs

### Precision (U/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	5.6	0.22	3.9%	0.61	11%
2	23	1.1	4.8%	1.2	5.2%
3	41	1.8	4.4%	1.9	4.6%
4	92	3.7	4.0%	4.7	5.1%
5	111	4.5	4.1%	5.0	4.5%
6	184	8.2	4.5%	9.2	5.0%
7	292	12	4.1%	15	5.1%
8	394	20	5.1%	21	5.3%

### Linearity (U/mL)

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	99	—	—
	4 in 8	53	50	106%
	2 in 8	27	25	108%
	1 in 8	13	12	108%
2	8 in 8	103	—	—
	4 in 8	57	52	110%
	2 in 8	30	26	115%
	1 in 8	16	13	123%
3	8 in 8	118	—	—
	4 in 8	65	59	110%
	2 in 8	35	30	117%
	1 in 8	17	15	113%
4	8 in 8	137	—	—
	4 in 8	68	69	99%
	2 in 8	35	34	103%
	1 in 8	16	17	94%

### Recovery (U/mL)


	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	11	—	—
	A	52	62	84%
	B	97	108	90%
	C	172	215	80%
2	—	18	—	—
	A	66	68	97%
	B	114	115	99%
	C	206	221	93%
3	—	45	—	—
	A	90	94	96%
	B	136	140	97%
	C	230	247	93%
4	—	67	—	—
	A	116	115	101%
	B	157	161	98%
	C	245	268	91%

### Specificity


Compound <sup>1</sup>	Amount Added <sup>2</sup>	Apparent U/mL <sup>3</sup>	% Cross-reactivity <sup>4</sup>
AFP	10,000 IU/mL	ND	ND
CA15-3	1753 U/mL	7.1	0.41%
CA19-9	4000 U/mL	ND	ND
CEA	10,000 ng/mL	4.8	0.05%
Cisplatin	100 µg/mL	ND	ND
Cyclophosphamide	1000 µg/mL	ND	ND
Doxorubicin hydrochloride	100 µg/mL	ND	ND
5-Fluorouracil	1000 µg/mL	ND	ND
Mitomycin C	100 µg/mL	ND	ND
Vincristine	1000 ng/mL	ND	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2017-2020 Siemens Healthcare Diagnostics. All rights reserved.