OTOCOMPACT STEEL

0476

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English

These instructions apply to all the versions of the OTOCOMPACT STEEL unit and all the related accessories. For this reason, not all the paragraphs of this manual will refer to the equipment you have bought.

PROFESSIONAL version MEDI-CARE SOLUTIONS will provide the authorised service centre with all the information necessary to install and repair the equipment.

No part of this manual is to be reproduced, stored in a retrieval system or transmitted in any form or in any means, i.e. electronic, mechanical, photocopying, translation or otherwise, without the prior written permission of MEDI-CARE SOLUTIONS.

OTOCOMPACT STEEL is a medical device which use is allowed only at specialized staff in ENT field. MEDI-CARE SOLUTIONS refuses any responsibility for damages and / or accidents caused by using the machine by non-authorized people.

The medical device OTOCOMPACT STEEL must be used matching (electrical connection) just other medical devices **and with exclusive use destination the patient area**

The expected useful life of the device is 10 years

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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English

1. INTRODUCTION

- 1.1 Intended and environment of use
 - 1.1.1 General characteristics and identification of the parts
- 1.2 Cautions
- 1.3 Symbology
- 1.4 Possible side effect, contraindications

2. INSTALLATION

3. USE OF THE INSTRUMENTS

- 3.1 Front commands
- 3.2 Insufflator
- 3.3 Syringe
- 3.4 Suction
 - 3.4.1 Suction Accessories
 - 3.4.2 Antimicrobial Filter
- 3.5 Medicaments atomizer
- 3.6 Menu'

4. USE OF THE ACCESSORIES

- 4.1 Germicide compartment
- 4.2 Larynx-mirror warmer
- 4.3 Camera support
- 4.4 Handpiece Led light source
- 4.5 Single or double LED light SOURCE
- 4.6 Articulated arm
- 4.7 Chair Wireless
- 4.8 Fixed holder for endoscope

5. MAINTENANCE

- 5.1 Filling the water tank
- 5.2 Removing the secretion tank
- 5.3 Cleaning of the unit
 - 5.3.1 Oil Pump Maintenance
- 5.4 Chart products for cleaning and disinfections
 - 5.4.1 Secretions tank disinfections
- 5.5 Disposable secretion tank
- 5.6 Empty of the condense
- 5.7 Discharge or scrap
- 5.8 Table products for cleaning or disinfections

6. TECHNICAL FEATURES


- 6.1 Technical data charts – absorbed power
- 6.2 Working periods of the instruments summary – identification of the unit

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1. INTRODUCTION

1.1 Intended of use and environment of use

Unit for use exclusively for the sector otolaryngologist born to perform both cleaning of the nasal cavity and ear through active tools is to assist and support toolbox of ENT specialist who is the only user required by MEDI-CARE SOLUTIONS Ltd.

1.1.1 General characteristic and identification of the parts

Unit for otorhinolaryngology completely autonomous and mobile on castor wheels composed of (see pic. 1):

- A) Instrument tray holder which is able to hold up to a maximum of 4 instruments. It has articulated arm which enable one to position it in the most comfortable way for the professional's activity. On the board is placed the electronic CONSOLLE with LCD display and touch screen.
- B) Aspiration instrument holder
- C) Anti-scratch working top with drop-saving back border embodied in the top; or alternatively desk for instruments trays with plastic cover.
- D) Housing for the auxiliary instruments complete with a series of drawers & spaces of various sizes and dimensions according to one's choice

1.2 CAUTIONS

- 10) To move the unit with the drawers closed and pc-monitor and Microscope within rest position. Unit weight it is about 120kg (based on configuration).
- 11) Do not install the machinery in environments which are exposed to explosions or frost
- 12) The temperature of the environment in which the equipment is installed and used must be between 10 ° C and 30 ° C. The relative humidity environment of use or storage must be understood between 30% and 70%; pressure 500hpa at 1060hpa
- 13) During installing, check that the electrical tension on the sticker on the machine corresponds to that in use in order to avoid damaging.
- 14) The building, repairing, modifying, calibration and all the other operations that require the opening of the unit should be carried out by personnel authorized by MEDI-CARE SOLUTIONS
- 15) The electric circuit in which the unit is connected to, should comply to the IEC 64.8 Sec. 710 regulation (regulation regarding electrical circuits in rooms used for medical purposes).
- 16) The machinery should be used exclusively following the instructions illustrated in this manual.
- 9) MEDI-CARE SOLUTIONS refuses any responsibility for the following: damages or accidents caused by the carelessness of the user, the noncompliance to the safety rules, a non-correct use of the machine, a use of the machine which is different to that of its purpose. MEDI-CARE SOLUTIONS denies all responsibility which has nothing to do with faults in the material and non-correct building.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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10) Before installation and functioning check that there are not loose parts, burnt electric cables, or various anomalies and if they are evident please contact MEDI-CARE SOLUTIONS or authorised technical staff.

11) Mobile phones and radios can influence the good functioning of medical device OTOCOMPACT STEEL.

12) If during the use of unit OTOCOMPACT STEEL happen strange changes of performance or various anomalies we suggest to interrupt all operations, switch off the machine and contact authorised staff MEDI-CARE SOLUTIONS in order to avoid any kind of damage.

13) Do not install this unit in places exposed to the risk of explosions, rich in oxygen or exposed to frost

14) The unit is not protected against splashing water and foreign bodies: **DEGREE IP20**

15) **WARNING!** To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth."

16) **WARNING!** There shall be no modification of this device

17) **WARNING!** The unit in question should not be used adjacent to or stacked with other equipment, if necessary use it close to or stacked with other equipment it is necessary to observe the DM to check whether, in the configuration in which it is used, its operation remains normal

18) The UNIT can be transported and stored with the environment:

Temperature: da -20°C a 60°C;
Humidity: 10% - 90%;
Pressure = by 500hpa at 1060hpa

19) The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment

20) **Warning!** In the event of a serious accident occurring in relation to the device, it is necessary to report to the Manufacturer and the competent authority of your country

• **APPLIED PART OF THE UNIT**

1. Terminal of Auricular wash syringe and medical syringe
2. Suction terminal
3. Insufflation Terminal

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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1.3 Symbols

Classification of the "BF" type machinery. "BF" type machinery (machinery of the type I, II, III or machinery housing an internal electric source with a protection degree against direct and indirect contact and in particular for:

- permitted current dispersion
- reliability of the ground contact)

Class I machinery BF.

This symbol indicates: TO CLARIFICATION CONSULT THE ATTACHED DOCUMENTATION

This symbol indicates: ATTENTION! CONSULT THE ATTACHED DOCUMENTATION BEFORE USE THE UNIT

Symbol that indicate "generic warning".

Earthling contact of the machinery

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1.4 WARNING, POSSIBLE SIDE EFFECTS

Always remember, before you use the MD (medical device) to make a careful medical history and a careful examination of the patient to recognize situations in which it recommended the use of active instruments such as the ear wash syringe, suction, insufflation, that can cause side effects such as redness, bleeding, abrasions, minor burns, etc. .

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CONTRAINDICATIONS

The contraindications are those specific situations in which the use of the tools are preferably not active carried and consequently the device must be used where appropriate with a restriction defined relative to the cost / benefit ratio resulting from the use of the system itself. It is therefore considered that any adverse or residual risk can be defined by MEDI-CARE; is the responsibility of the specialist do the analysis costs / benefits of the use of MD based on their knowledge and experiences.

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2. INSTALLATION

The installation of the equipment must be made respecting the instructions of the following booklet of instructions. It must be performed by authorized personnel MEDI-CARE SOLUTIONS or, given the simplicity of the operations can also be performed by the customer under its responsibility.

The user should always make sure that the electric circuit complies to IEC 60364-7-710 regulation or special rules of your country. The operations are (after unpacking the unit and take IFU):

- Taking unit out from the pallet from basement or solid parts and not from top or instrument tablet,
- Unpacking the instruments and put them on the tablet as picture 1
- Unblock the wheels brake and positioning unit
- Screw, if it's present, the monitor and microscope support pole
- Full the tank but not to the rim
- Insert the plug in a socket ➔ It is advisable to insert the plug in a socket with a bipolar switch compliant to IEC and having the following technical characteristics:

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3. INSTRUMENTS

In pic2 it is showed a general picture of the unit instruments. For the instruments use available on the machine it is necessary first of all to switch on the machine itself by pushing the button (A) placed on the little front of the working top. It will lighten with a green light and the appearing of the screen first menu (B) in the LCD display of the console indicates that the machine is effectively functioning. The switch-off of the machine is effectuated through the same button.

3.1 Console LCD with touch-screen drives

Here following a list of the drives located on the console of the unit (picture 3):

1. Up in the right do appear the switches +, - for the variation of light intensity of the LED hand-piece if taken out (percentage indication);
2. Switch for activating the water heating cycle;
3. Button for switching on germicidal lamp;
4. Button for activating insufflators through pedal on-off;
5. Button "Menu" for entering in regulating section and visualize parameters.
6. Chair section

N.B. THE CORRECT FUNCTIONING OF EACH BUTTON WILL BE ILLUSTRATED IN THE PARAGRAPHS CONCERNING THE INSTRUMENTS.

3.2 Insufflator (see picture 2-E)


For a correct functioning of the instruments we suggest to respect the following using times:

Work = 5 minutes' - Pause = 10minutes

The tool delivers compressed air up to a maximum of 1.8 Bar for insufflation within ENT. It is possible to use the insufflator or directly through its extraction from the seat or through the foot-control. In that case it is enough to push the button 4 on the console LCD (see 3.1); then on the display appears the signal FOOT CONTROL, and consequently press the foot-control.

To return to the direct drive it is enough to re-push the above mentioned button (the signalisation FOOT-CONTROL disappears). During its functioning on the display appears the INSUFFLATOR'S symbol and the approximate value (+/- 5% of precision) in BAR on the exit pressure. To vary the insufflations intensity you must act on the knob placed under the instruments tablet (SEE PIC 2-L). To increase the pressure it is enough to turn the knob to the right.

Combining with it, there is a Politzer olive in autoclavable material (universal measure) directly engaging to the adjusting gate (see picture 3.1)

**Fig.3.1
Pic. 3.1**

The Display shows the value of the pressure in exit of the neck/spout; therefore to obtain the maximum pressure it is necessary to wait for 3-5 seconds and then push the black button on the handpiece!

Notice: use only spare parts MEDI-CARE SOLUTIONS.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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WARNING! In the event that the activation of the instrument, the compressor does not work, press the black button to download the air in the circuit!

Notice: All these parts are cleaned, washed according to its internal procedures with common disinfectants. See Table 5.7

3.3 Ear washings syringe (picture 2 C)

For a correct functioning of the instruments we suggest to respect the following using times: **Work: 5'- Pause:10'**

The tool provides approximately 350 ml / min of water at a pressure of 1.8-2 Bar. The functioning of the instruments is obtained by extracting it from its seat. On the display appears a symbol indicating that the instrument has been extracted; now for letting water flow out, you can press any of the 2 buttons placed on the hand-piece of the washings syringe, or both. The water heating cycle starts after having pushed the button 2 (picture 3) on the console (see chapter 3.1); on the display of the console appears the signalisation "ON 37" in intermittent way; when the water temperature has reached the 37° C the signalisation "ON 37°C" stops lightening and an acoustic signal advises the operator about the end of heating cycle. Now to let water come out you can push one of the two buttons placed in the washings syringe hand piece or both.

NOTE: The first heating cycle can last some minutes (about 10).

WARNING! In the event of a reserve water audible and visual alarm display "RESERVE"; immediately turn off the machine and refill the tank (See 5.1 chapter)

N.B. After a prolonged pause with the machine off, when the temperature is reached, it is always advisable to let come out from the syringe a water flow for about 10 seconds so that to asset the machine.

N.B. In the evidence that the water temperature overcomes the set value, a safety thermostat blocks automatically the washings syringe functioning, until the normal cycle is re-activated. On the display appears the sequence "TERMIC ALARM" to indicate that the thermostat is active. In that case it is necessary to contact an authorised technician of MEDI-CARE SOLUTIONS, which is the only authorised to intervene.

Note: For the maintenance, cleaning please follow the instructions in the Technical Manual here attached and chapter 5.8.

3.4 Suction (picture 2 D)

For a correct functioning of the instruments we suggest to respect the following using times: **Work: 5' - Pause: 10'**

The tool allows you to draw with a maximum vacuum of 0.85 bar +/- 5% and a real flow up to 50 litres / min +/- 5%. The instrument's functioning is obtained by extracting it from its seat. They will display a symbol indicating the successful extraction and vacuum BAR value obtained with a precision of +/- 5%. To vary the intensity of suction it is necessary to operate on the knob placed in front of the unit (picture 2N). The vacuum control on the front of the unit indicates the empty created by the empty pump. Turning the knob to the right the suction power is increased and viceversa.

Note: The display shows the vacuum-control releases the vacuum in function of the resistance that is located in the exit of the silicone neck; therefore to obtain and verifying the maximum depression it is necessary to close the exit.

NOTICE: If the container is full of secret operative system overflow aspiration is stopped, so to drain (see 5.2)

Notice: The useful capacity is stated on the vessel, the unit cannot work on inclined planes

Notice: Use only spare parts by MEDI-CARE SOLUTIONS.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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ATTENTION! In order to avoid the pump block, don't suction large quantities of non-viscous liquids in a short time; aspirate also in the rinsing phase, small quantities of disinfectant liquids with 2/3mm in diameter cannulas.

If required, the unit can be provided with the **automatic discharge** through a system that drives the secret to a final pipe connection to be connected to the discharge net. This pipe (**the connection should be done by the customer**) is sited in the lowest part of the unit and it is suitable for tubes section 6 x 8:

ATTENTION! Do not suction DENATURED ALCOL or sodium ipocloride solutions or pure deionized water (see cap.5)

Notice:Carefully read the manufacturer's manual attached

If required, in addition to the automatic discharge, it can also be possible to add to the unit the automatic rinsing system of the suction pipes (side "PATIENT").

There is a tank of about 750cc where it has to be put the disinfectant solution that should flow through the pipes.

NOTE: In the case of automatic or multipurpose secretion jar MEDI-CARE SOLUTIONS suggests using a solution called PULIJET NEW PLUS or similar; do not use solutions containing sodium hypochlorite. In the case of disposable jar, MEDI-CARE SOLUTIONS suggests to add to one liter of water, about 150cc of solution of sodium hypochlorite or Amuchina and fill the tank of the rinse in.

To turn on the rinsing cycle, please follow these instructions:

1. PUT EVERY INSTRUMENT IN ITS LODGMENT
2. PUSH THE BUTTON ON THE DISPLAY (SEE BELOW)
3. REMOVE THE SUCTION PIPE AND DISCONNECT THE TERMINAL PART
4. INSERT THE CANULA IN ITS CONNECTIONS (see mail below)

The cycle (more or less 40-60 sec) is automatic: at the end it automatically switches off. The customer has to remove the cannula and put it in its support and if necessary to repeat the operation. During the cycle it is not possible to make other operation.

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WARNING! Avoid to drip the solution for rinsing on the basis of the unit

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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Fig 3.4.1
Pic. 3.4.1

Attention! If you require oil suction pump predisposition, it is necessary to pay attention to the unavoidable ordinary maintenance. For this, please, refer to this instruction manual and pump producer instruction.

3.4.1 SUCTION ACCESSORIES

The unit is set with a pulling-out neck headed by a piece of silicone pipe (see pic. 3.4.1); moreover together with the machine is supplied a kit for setting the funnel for secretions collecting (see pic. 3.4.2) which contains a pipe fitting to be coupled to the suction pipe + funnel (see assembling phase).

Notice: All these parts can be cleaned, cleaned or washable according to their internal procedures with common disinfectants (See table 5.7); alternatively, they are autoclavable at 121 °C at 1.1 bar for 10 min (see table 5.8).

Fig 3.4.2
Pic. 3.4.2

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WARNING: Do not keep the rolling shutter always closed because it could cause suction locking.

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N.B. In case it is not used to avoid dispersion, it is advisable to keep the tap closed while the first suction is used.

N.B. During the RINSING cycle, keep the tap open to ensure correct recirculation operation.


Fig 3.4.3
Pic. 3.4.3

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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Work: 5 minutes / Pause: 10 minutes

The tool provides dry compressed air to a maximum pressure of 1.8-2 bar to spray disinfectant powder or liquid or a choice of air blow through the politzer cannula. The functioning of the instrument is obtained by extracting it from its seat. On the LCD display a symbol indicating that you are using the instrument will appear. After the application of the ampoule to the syringe to obtain the vaporization of the medicaments products, operate on the lever of the instrument hand piece (see picture 4). It is possible to adjust the pressure and the flow of air through the tap under the tablet. In counter clockwise to increase the pressure and vice versa.

WARNING! To avoid keeping the black tube too much under pressure, we suggest not to keep the syringe extracted too long without using it.

NOTE: USE ONLY SPARE PARTS supplied by MEDI-CARE SOLUTIONS.

Note: For maintenance, cleaning and sterilization, follow the table 5.8 attached to this document.

3.6 Menu




On the LCD of the console there is the MENU switch (button 5, picture 3) that allows to visualised and set some parameters of the machine working; in fact, entering the STAT section you can visualise the following parameters:

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- Visualise functioning time of the machine
- Visualise functioning time of the UV lamp
- Reception of the chair through RF
- Select the language (Italian, English, French)
- It possible to: Pushing the button **"change name"** you can set the doctor name, by means of the pseudo-keyboard.
- It possible to setting the **DATE**; set the date and time and then push the button **"SAVE"**
- It possible to set the **RF channel** for chair communication

NOTE: To go back to the previous page of the menu, just press again the MENU for any of the sub-pages.


**Fig. 5
Pic. 5**

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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For set the time of Mirror warmer press SETUP button and regulated the time with +/- push button.

WARNING! Don't touch the Poffs and Voffs buttons.

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Supporto Manipolo Led – Handpiece
Led light source Support

Innesto Manipolo Led – Handpiece
Led light connection

4. USE OF ACCESSORIES

4.1 Germicidal Space (picture 2 I)

The germicidal lamp can be switched on and switched off through the button 3 "UV" (see picture 3 section 3.1) placed on the console. A safety switch turn automatically off the lamp in case of accidental opening of the germicidal space, and turn it on when the same is closed. Please note that pushing the button "MENU" and then the button "STATS" (picture 5) on the display appears the time (in hours) of use of the UV lamp; after 500 hours of use on the display appears automatically the sequence "CHANGE UV" that indicates the substitution of the lamp in order to preserve its efficiency.

WARNING! We underline once more that the germicidal space does not constitute element of substitution of the sterilizing process but it is only of help and support.

Notice: Do not use lamps with different features from these imposed by MEDI-CARE SOLUTIONS.

4.2 Larynx mirror warmer (Picture 2 G)

The functioning of the larynx-mirror warmer is obtained by pushing the button placed in the front column of the unit (see photo 2 G). Please note that by pushing 3 times the switch nr. 6 "MENU" on the console it is possible, through the switches +/-, set the lit up time; the range goes from a minimum of 10 sec. To a maximum of 60 sec. The pushing of the button during the functioning causes the interruption of the air flux before the set time.

WARNING! Do not put tools in the grate escape.

WARNING! Do not touch the grate after use.

Notice: The power of element is 110W

4.3 Camera wire.

The UNIT OTOCOMPACT STEEL (and Professional version) can support a camera wire (see enclosed manual). See Pic 2-N

4.4 Handpiece light source LED

On request it is possible to predispose in the column module a white LED light source with temperature of 5500 °K and 180 Lumen with option of single or double connector for eventual connection to endoscopes or laryngoscope. Insert the hand-piece connector in the socket under the module; the ignition is automatic, you only have to take it out from its place (see picture 6) then, with the keys +/- on the touch screen you can regulate the intensity (from 30 to 100%)


Fig.6/Pic. 6

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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4.5 Single or double LED light source

Otocompact steel is equipped with (a YELLOW LED or white) single or double light source, with f.o. cable.

The single light source is automatically switched on when the operator takes out the cable from its site. The regulation of intensity can be done through the knob sited in front of it (see pic. 3.5.1) In the case of double Light source to turn and to adjust the light intensity use the regulator from the front panel see. Pic 7-c

Notice: The LED NEUTRAL or WHITE (4100-5500 °K) their duration to about 50,000 hours.

NOTICE: For replaced LED, call **MEDI-CARE SOLUTIONS** or **MEDI-CARE SOLUTIONS** authorized distributor.


Fig. 7-b
Pic. 7-b

Fonte di luce LED/ LED LIGHT SOURCE

NOTE: One of the sources can be activated directly through the sensor in the tablet

WARNING! DO NOT LOOK THE LED LAMP FOR TOO TIME; MAY BE HARMFUL TO THE EYES.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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4.6 Articulated arm

On request the unit can be set with an articulated arm with rotating bushing for monitor or tray holder with supportable a max weight **of 7 Kg.**

4.7 Wireless Control of the Chair

When the unit is combined to a Wireless activated Chair its display (touch screen) will show on the bottom right side the symbol of the chair; by pushing on the symbol the Chair control will be activated.

All touch screen control buttons have a dual functionality; the instant touch of the button will activate an automatic movement that will stop at the end of the actuator drive. If the same button is pressed for more than one second, the correspondent movement will be manually controlled, and it will stop when the pressure on the button ends.

It is possible to set 3 memorized working positions M1, M2, M3 by pushing the STORE button when the desired position is reached by the Chair. The Display will also show the M0 position which is pre-set and it allows the Chair to automatically reach the Patient Exit Position (called 0 position). Is present a **"Synchro"** button; it can give a contemporaneous and out of phase movement of the backrest with the leg rest. If for one minute you do not make any movement, the system returns to normal (independent controls). To exit the Chair Control section of the Touch Screen and return to the main section it is necessary to push the "Menu" button displayed on the bottom right side.

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4.8 Cleaned and Dirty Optics Housing

On the Unit is standard planned an integrated support for the dirty and cleaned optics. There are 2 closed supports for the resting optics and/or for inserting some disinfectant solutions.

Warning! Do not use aggressive solutions and do not maintain them on the support more than 24 hours, wash the supports with hot water after use.

Warning! To preserve the supports we suggest to use solutions for the sterilization at cold or generators of free in order to conserve the supports.

Attention! For the disinfectants, replenishing is needed, in order to avoid the dripping on the compartment, take the closed support out from its seat, and then proceed with the replenishing.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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Warning! The Optic supports are movable to allows the disinfection control

NOTICE: You can set a timer to indicate the exposure time of the optics inside the supports.

4.9 Microscope pole support

In the unit, if required, there is a support for a microscope (Only model of Kaps or Zumax).

Notice: Installation of the microscope must be carried out by specialized personnel and by participating in the manufacturer's user manual.

WARNING! For use, maintenance and cleaning of microscope see attached instruction manual.

WARNING! It is recommended not to positioning the arm of microscope to the rear of the unit.

Warning! Do not open or take off all drawers when microscope is installed.

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5. MAINTENANCE

Warning! The following proceeding have to be carried out when the machine is not plugged into the mains and with specific PPE aids to avoid biological contamination.

5.1 Filling in the water tank

In case the water in the tank is finished, an acoustic signal and the alarm "RESERVE" appear on the touch screen. It is necessary to fill the tank up, so please proceed in the following way:

- 1- Switch off the device
- 2- Open the door of the auxiliary instruments compartment
- 3- Take the tank out, activate the 3-way valve (IF PRESENT) in order discharge the air, unscrew the cap and add water with a funnel, without reaching the maximum level.
- 4- Close the cap, close the 3 way valve
- 5- Close the door and switch on the device

Fig. 5.1
Pic. 5.1

Notice: It is recommended to use only distilled H2O mixed with disinfectants solutions to avoid problems of limestone deposits in the tubes and proliferation of bacterial charges.

WARNING! Whit automatic discharge system do not use pure distilled water as the level probes may not work and lead to the system blocking due to overflow.

WARNING! Do not leave water in the tank for more than 12 hours; in this case, to avoid bacterial contamination, it is recommended to empty the tank and dry it with a sterile cloth or alternatively, add a commercially available preservative solution to the water

5.2 Cleaning of the secretion tank

At the end of a working day, it is necessary to clean the suction canister. In order to carry out this operation, proceed in the following way (see picture 6):

- 1- Unplug the unit
- 2- Open the door where the tank is housed
- 3- Unscrew the tank by turning it anticlockwise
- 4- Empty the tank and clean it thoroughly
- 5- Re-screw the tank to the cover, close the door and plug the unit again.

Fig. 6
Pic. 6

WARNING! Do not reverse the PATIENT and VACUUM tubes!

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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

Fig. 2
Pic. 2


5.3 General maintenance and spare parts

To guarantee the correct functioning of the instrument MEDI-CARE SOLUTIONS suggests making a periodical revision. MEDI-CARE SOLUTIONS does not take charge of any responsibility for damages, hitches or inefficiencies due to the substitution of non-original spare parts.

5.3.1 OIL PUMP MAINTENANCE

1. Oil Level Check

It is suggested to check the oil level of the pump's engine every 2-3 working months, according to the hours of activity sustained by the Pump:


2. Oil Refill

After approximately 500h of suction use, and in case of dark oil due to its deterioration, proceed to its replacement by following point 5.4 of instructions manual Rif. LB.4 9601057: see page 9-5.2 see pict. 1).

3. Filter Replacement

The pump's filter is able to reduce residual fumes for 2000 working hours, after such period of time it is necessary to replace it with a new filter; see the following pictures 2,3:

B. To unscrew the old filter and to replace it with a new one.
Attention: To tighten gradually without excessive force.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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5.4 Cleaning the unit

The choice of materials used for building the unit has been carefully taken in order to correspond to the cleaning and disinfecting needs. Given the big differences between chemical products used in a medical study, it is possible to ruin the tops. The major damages are caused by the amount of time that these chemical products are left on the tops. It is important to immediately dry the interested parts with a damp cloth.

To clean the painted or metal parts, use a specific product for these materials which doesn't damage their composition.

IMPORTANT! Do not use products with alcohol, ammonia or petrol.

5.4.1 Cleaning and maintenance of secretion tank

The tank which contains the aspirated liquids can be cleaned with water and detergents which are compatible with the materials (see chart) used and/or sterilized in an autoclave at 135°C. Separate the main parts of the tank before the sterilization, insert the glass upside down in the autoclave and do not place anything on top of it whilst disinfecting.

DO NOT USE SOLVENTS OR ALCOHOL FOR CLEANING AND DISINFECTING. THE USE OF THESE PRODUCTS CAN DAMAGE THE PLASTIC PIECES OF THE SYSTEM.

Always substitute the damaged or faulty parts with original spare parts supplied by MEDI-CARE SOLUTIONS.

In case of cleaning in autoclave please stick strictly to the temporary chart shown in following table 5.5

5.5 Vase secreted disposable

The unit can be setted for single-use kit of different type. **Please keep attending during the replacement of the ampoules in order to avoid to invert the suction tubes VACUUM and PATIENT, because it could damage the motor. If the intake of liquids or solids into the engine, immediately turn-off the suction and call the service center for maintenance.**

5.6 Empty of the condense

Near the panel there is a filter for collecting the condense created by the compressor; in case it is filled it is necessary to empty it; place a glass under the filter, turn the knob and pull downside the knob, empty the condense and close again the knob!

Warning! If the tap is not closed the compressor leaks air through it!

5.7 Discharge or scrap

The machine does not contains dangerous elements and can therefore follow normal procedures of discharging and scraping, in harmony with rules of Country where the unit has worked.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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5.8 Instrument sterilization chart

* Instrument and detachable parts sterilization chart		
INSTRUMENT	TEMP	KIND OF STERILIZATION
Ear Wash Syringe handpiece***	134 +/- 1°C	In autoclave at 2.13 Bar >= 4Min /500 cycles
Ear Wash Syringe terminal ***	134 +/- 1°C	In autoclave at 2.13 Bar >= 4Min /500 cycles
Medicament Atomizer terminal ***	134°C	In autoclave at 3 Bar >= 5Min /500 cycles
Suction Cannula	134°C	In autoclave at 2.1 Bar = 10Min/100 cycles
Suction tube VHR (if present)	121°C	In autoclave at 1.1 Bar = 10Min/500 cycles
Canister multipurpose ***	134°C	In autoclave at 2.1 Bar = 15Min/30 cycles
** Cleaning and disinfecting products chart (other parts)		
Disinfectants		Pharmaceutic products and cosmetics
+ bactol 5% + clorammina - DDT + delegol 5% = dimammina 5% = iodine dye - fenic acid + lysoform 2% - TB lysoform + bactol + merfen 2% + octozone + perydrol + resorcina 1% = sagrotane 5% + alcohol to be purely burnt + mercury chloride + trosilina G + oxiginated water = zefirol + sodium ipocloride (amuchina)		+ blood plasma + sun tan lotion + hydroplex = iodine dye + lanolyn = mantolo 90% in alcohol -nail varnish - solvent + mouth wash + periston R + vaselyn + vicks rub
CHART SIGNS		
+ compatible - non compatible = compatible but with caution		

* For additional material autoclavable see in single paragraph of this manual

** For more products, please contact MEDI-CARE SOLUTIONS.

*** See Instruction of use of manufacturer

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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6. TECHNICAL FEATURES

6.1 Technical data chart

Electrical tension	230VAC +/-10%
Frequency	50/60 Hz
Power absorbed by the unit (according to the configuration)	1200 VA
Weight	120 Kg
Maximum number of instruments	4 pcs
Maximum air pressure	2.8 Atm
Air capacity	26 l/m
Maximum water pressure	2.2 Atm
Water capacity	0.4 l/m
Dimensions (in mm):	
Height of the working bench	805 mm
Width	860 mm
Depth	745 mm

Absorbed power chart

DESCRIPT.	ABSORB.	TENSION	SERVICE
Germicide Lamp	11 W	230 V ac	Continuos
Aspirator	120 W	230 V ac	Continuos
Compressor	120 W	230 V ac	Continuos
Pump	12 W	230 V ac	Continuos
Electrovalve	24 W	24 V dc	Continuos
Resistance	800W	230 V ac	Discontinuos
Led lights	3W X 3	5 V dc	Continuos

6.2 Summary of the working periods of the instruments

INSTRUMENT	WORK	OFF
Syringe	5'	10'
Insufflator	5'	10'
Suction	5'	10'
Atomizer	5'	10'

Building Datas

1	2	3	4	5	6	22	23
7	8	9	10	11	12	21	20
N° matricola							

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

6.3 SCHEDA TECNICA/TECHNICAL FEATURES

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NOTE: MEDI-CARE SOLUTIONS on request will provide electrical circuit diagrams, the parts list, descriptions, calibration instructions or other information to assist the AUTHORIZED SERVICE PERSONNEL in the repair of parts of the MD that are considered repairable.

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
In the back part of the unit is located the label of the fuses – its values are listed in the chart 6.2 of the following MI.
Inside the machine, the protective mass is indicated by the following stuck symbol:

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7.2.2 Norme Applicate/ Standard Rules applicated:

- CEI EN 60601-1-2 (2014)
- CEI EN 60601-1-2/EC (2014)
- IEC 60601-1 third edition;
- CEI EN 60601-1 (2007);
- CEI EN 60601-1/EC (2010)
- CEI EN 60601-1/A11 (2012)

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

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www.euroclinic.it ; info@euroclinic.it.

7.3 DISPOSAL

According to the legislative decree 14/03/2014, n. 49, "Implementation of the European Directive 2012/19/UE about decrease of dangerous substances on electric and electronic equipment" and moreover it refers to garbage disposal, the symbol of crossed out garbage bin (on the back of the equipment) indicates that the product at the end of its working life must be collected separately from other garbage according to the following law 31 may 2016 n° 121 and 8 march 2010, n.65 or follow your country laws.

The separate collection of rubbish of this equipment (in disuse) must be arranged and coordinated by the manufacturer. When the user wants to throw away the equipment, he has to contact the Manufacturer Company that will indicate all procedures to follow for the disused equipment separate collection of rubbish.

Therefore, the appropriate separate collection of rubbish, aimed to actuate afterwards disused equipment recycling, treatment and to ambient garbage disposal, contributes to avoid possible bad effects on the environment and on the health. In addition, these procedures raise the reemployment and/or the recycling of the materials of which the equipment is composed.

The unauthorized garbage disposal, on the part of the holder, requires administrative law sanctions application.

MEDI-CARE SOLUTIONS S.R.L. remains available for further information about this matter:

Registered Office: MEDI-CARE SOLUTIONS Ltd

Via Della Zecca 1, 40121 Bologna (BO)

Facility: Via Pietro Nenni 3, 40026 Imola (BO) - ITALY

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The electro-medical device requires special precautions regarding EMC (Electromagnetic Compatibility) and must be installed in accordance with the EMC information contained in the accompanying documents.

Portable and mobile radio-communication equipment can affect the operation of EUT (Equipment under test).

The use of accessories, transducers and cables other than those specified, with the exception of those sold by the manufacturer as spare parts for internal components, may cause an increase in emissions and a decrease in immunity.

The device or the system must not be used near other devices and, if it is necessary to use it near other appliances, the electro-medical device must be observed to check normal functioning in the configuration in which it is used.

The electro-medical device has been tested and founded in compliance with the emission and immunity limits of electro medical equipment according to IEC60601-1-2: 2014. These limits are Specifically designed to provide adequate protection against dangerous interferences in a typical medical room. However, there is no guarantee that interferences will not occur in a particular installation. If the electro-medical device, interacting with another one, causes or receives detectable interferences, the user is invited to limit them by adopting one or more of the following measures:

1. Reorient or relocate the receiving device;
2. Increase the distance between the appliances;
3. Connect the equipment to an outlet on a circuit different from the device (s) causing the interference;
4. Contact the manufacturer or local technician for assistance.

11	Ing. Francesco Marrone	04-10-2021	DG Fabio Casella	04-10-2021
Release	Redatto	Data	Approvato	Data

Guida e dichiarazione del fabbricante – emissioni elettromagnetiche/ *Guidance and manufacturer's declaration - electromagnetic emissions*

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Guidance and manufacturer's declaration - electromagnetic emissions

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