
**TECRES®**  
ADVANCING HIGH TECHNOLOGY

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**INSTRUCTIONS FOR USE**  
*For the attention of medical personnel*

- CEMEX GENTA ID GREEN is a radiopaque bone cement containing Gentamicin that can be used for fixing prosthetic joint implants to bone tissue.
- CEMEX GENTA ID GREEN is a low polymerization temperature cement as it uses a powder to liquid ratio that is higher than the traditional 2:1 ratio. CEMEX GENTA ID GREEN releases Gentamicin over time. Colorants lend the bone cement a green colour which allows it to be easily identified during surgical procedures.
- CEMEX GENTA ID GREEN is a STERILE and SINGLE-USE medical device.
- CEMEX GENTA ID GREEN is an ideal bone cement for both manual and syringe application.

The package contains 1 40g pouch of sterile powder and 1 vial containing 16.7g of sterile liquid. The liquid component is sterilized by filtration and the powder is sterilized via gamma- ray radiation.

|                          | <b>Cemex Genta ID Green</b> |
|--------------------------|-----------------------------|
| Liquid component         | 16.7 g vial                 |
| METHYLMETHACRYLATE       | 98.20% w/w                  |
| N,N-DIMETHYL-p-toluidine | 1.80% w/w                   |
| HYDROQUINONE             | 75 ppm                      |
| Powder component         | 40 g pouch                  |
| POLYMETHYL METHACRYLATE  | 82.68% w/w                  |
| BARIUM SULPHATE          | 10.00% w/w                  |
| BENZOYL PEROXIDE         | 3.00% w/w                   |
| GENTAMICIN SULPHATE      | 4.22% w/w*                  |
| COLOURED PIGMENTS        | 0.10% w/w                   |

\*Equivalent to 1.0 g (1.0 M.I.U.) equalling 2.5% of Gentamicin base in a 40 g unit.

### **INDICATIONS**

The CEMEX GENTA ID GREEN bone cements are indicated for fixing prosthetic joints to bone tissue.

These cements are particularly indicated where there is the risk or presence of infections caused by organisms sensitive to Gentamicin.

### **CONTRAINDICATIONS**

The use of CEMEX GENTA ID GREEN must be carefully considered in the presence of myasthenia gravis or hypersensitivity to Gentamicin, the colouring agents (E102 and E133) or any of the other components of the bone cement. When the loss of musculature or neuromuscular compromise of the unhealthy limb would make the surgical procedure unjustifiable.

### **PRECAUTIONS**

Read the Instructions for Use carefully.

For the safe use of the CEMEX GENTA ID GREEN bone cement, the surgeon should be thoroughly familiar with the product's properties, handling and preparation characteristics, its limitations for use and the correct application technique.

It is therefore recommended that the surgical staff perform practical tests before the operation using the same instruments and environmental conditions. The surgeon must ensure that the prosthesis to be implanted is suitable for use with bone cement.

### PRECAUTIONS FOR PREPARATION OF THE CEMENT

- Make sure that the package is intact and that the components have not undergone any changes such as a yellow or brown coloured powder or syrup-like consistency of the liquid.
- Temperature has a great influence on the preparation characteristics of any bone cement. Temperatures above 23°C for the product, preparation tools or in the operating room lead to an acceleration of the different preparation stages. Similarly, a reduction in temperature causes a deceleration in these times. **Before using CEMEX GENTA ID GREEN, it is strongly recommended to ensure that it has been stored at a temperature of 23°C ± 1° C over the past 24 hours;**
- Prolonged exposure to high humidity (>70%) can cause the viscosity to increase with a subsequent acceleration in the preparation and application times of the cement.
- Ensure that the tools used to prepare the cement are specifically suited for that task.
- In order to avoid introducing glass fragments into the mixture, do not open the vial directly above the mixing vessel.
- Do not mix the cement under a flow of air, as this will cause the liquid to evaporate quickly and therefore modify the cement's performance.

### PRECAUTIONS RELATING TO THE APPLICATION OF THE CEMENT

Data obtained from clinical investigations indicates the necessity for strict compliance with strictly aseptic surgical techniques. It is important to remember that any deep wound infection from the surgery constitutes a serious risk for the success of the prosthetic implant. This infection may be latent and not manifest itself even several years after the operation.

The use of bone cement requires extremely close collaboration between the surgeon and the anaesthetist. The anaesthetist must be told when the bone cement is implanted during the operation. In some cases events defined as "bone implantation syndrome" (BCIS) may occur which are characterized by a number of clinical features that include hypoxia, hypotension, cardiac arrhythmias, increased pulmonary vascular resistance (PVR), and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology. These phenomena are commonly associated with, but is not restricted to, cemented hip arthroplasty which usually occurs at one of the five stages in the surgical procedure: femoral reaming, acetabular or femoral cement implantation, insertion of the prosthesis or joint reduction (Donaldson et al., 2009, Br J Anaesth).

The blood pressure of patients should be monitored carefully during and immediately following the application of the bone cement. In addition, overpressurisation of the bone cement should be avoided during the insertion of the bone cement and implant in order to minimise the occurrence of pulmonary embolism.

To minimize the risk of placing cells, bone marrow fat, bone fragments or other materials into the blood stream, a thorough irrigation of the bone cavities with Ringer's or saline solution is recommended before inserting the cement.

During cementation it is important to maintain the positioning of the prosthetic component by manually applying pressure until the polymerization process has been completed. This must be done in order to achieve proper fixation.

### PRECAUTIONS FOR THE USER

Extended breathing of the vapours from the liquid component can cause a soporific effect; excessive exposure to the concentrated vapours of the monomer may cause irritation of the respiratory tract and the eyes. Avoid contact of the monomer with skin or mucous membranes. Contact dermatitis has been reported in susceptible subjects. Wearing a second pair of gloves and strict adherence to the mixing instructions may reduce the possibility of hypersensitivity reactions. The liquid component of the bone cement is a powerful lipid solvent, therefore, avoid direct contact with the human body. Rubber or latex gloves do not always provide effective protection against the monomer. More suitable gloves made of different materials are however available on the market. Check that these materials are suitable for contact

with bone cement (see the relative technical data sheets).

Do not allow rubber or elastomeric tools to come into direct contact with the liquid component. The vapours from the liquid component may have an adverse reaction with soft contact lenses.

The liquid component is flammable and highly volatile and the operating room must therefore be properly ventilated. The liquid component and/or its fumes must not be exposed to an open flame or incandescent material.

After the mixing of the two components, the consistency of the bone cement changes after just a few minutes, rapidly increasing its viscosity until becoming a hardened mass that securely anchors the prosthesis to the bone. The reaching of this state can be easily noticed by the increase in temperature of the bone cement.

The cement will spontaneously cool after a few minutes, indicating the completion of the reaction and also the moment in which the prosthesis can be released.

### **INTERACTIONS**

The use of CEMEX GENTA ID GREEN should be carefully assessed if administered concurrently with ototoxic or nephrotoxic drugs.

### **SPECIAL WARNINGS**

To use CEMEX GENTA ID GREEN bone cements in a safe and effective manner, the surgeon must be properly trained on the use of bone cements and must also be familiar with the handling characteristics, the limitations of use and the correct application of the cements.

Tests from clinical studies clearly indicate the necessity for strict compliance to correct aseptic surgical techniques.

Correctly cemented prostheses are stable and durable. Loosening and fracture of either the cement or the prosthesis, or both, can occur due to disease, trauma, inadequate cementing technique, mechanical failure of the materials or latent infection. All patients should therefore undergo regular and long term check-ups after the operation.

Do not add foreign substances to the cement's components.

### **Warning CEMEX GENTA ID GREEN bone cements reach temperatures above physiological temperatures during the exothermic polymerization reaction.**

In the event of an infection following the operation, patients should immediately consult their physician in order to reduce the risk of infection to the implant.

The use of CEMEX GENTA ID GREEN as the first choice in the fixation of a prosthetic implant should be carefully assessed since this could increase the risk of developing Gentamicin resistant bacteria.

The use of CEMEX GENTA ID GREEN should be carefully assessed for patients with coagulation disorders and for patients with severe cardiopulmonary failure.

The use of CEMEX GENTA ID GREEN should be carefully assessed for patients with pre-existing renal impairment.

**CAUTION:** Do not resterilize and/or reuse. This is a single-use device intended for use on one patient. Do not divide the product into 2 or more portions for use at different times. This would be considered a reuse of the product which could lead to an error in the ratio between the powder and liquid components as well as a loss of sterility. Resterilization must not be performed as this would create a risk of infection to the patient; resterilization could also alter the morphology, the efficacy of the antibiotic and the mechanical properties of the device, causing the device to malfunction with serious risk to the patient's health.

The residual material should be considered surgical waste and must therefore be discarded at the end of the surgical procedure.

### **USE DURING PREGNANCY AND LACTATION AND IN CHILDREN**

There are currently no tests on the safety of using bone cement during pregnancy or while breast-feeding. Bone cement should not be used during

the first trimester of pregnancy, and during the rest of the pregnancy period it should only be used for life-threatening illnesses.

The CEMEX GENTA ID GREEN bone cements are indicated for use in children only in the case of limb preservation where no other procedure is likely to succeed.

## DATA FOR USE

### PREPARATION

Bone cements are temperature sensitive. Any **increase in temperature**(of either ambient, the cement components or the mixing equipment) above the recommended temperature of 23°C will **reduce the preparation times** of the cement. A **reduction in temperature increases these times**.

Open the outer packaging and place the pouch with the powder and the liquid vial on a sterile shelf in the operating room. Open the vial and transfer all of the liquid into a suitable mixing vessel. Open the pouch and pour all of the powder over the liquid.

To minimize the amount of air bubbles, mix the cement by moving the spatula from the outer rim of the vessel towards the centre. The liquid must moisten all the powder; use the spatula to carefully work any lumps of unmoistened powder into the mass of moistened dough. The amount of mixed cement required for clinical use is determined by the surgeon in each individual case. **Caution: Do not arbitrarily change the ratios of the liquid and solid components.**

The mixing time is between 1 and 1.5 minutes, but the actual time is influenced by the temperature, humidity and the mixing technique and is thus determined based on the surgeon's experience.

Syringe application: after mixing, insert the cement into a suitable sterile device. Based on his experience, the temperature and humidity during storage, in the operating room and the injection accessories, the surgeon will decide when the cement is ready to be extruded.

Manual application: after mixing, continue kneading the dough until it no longer sticks to the gloves. The dough is now ready to be applied. The temperature and humidity during storage of the product, in the operating room and of the mixing tools as well as the surgeon's hands can cause differences in the preparation and application times of the cement, which are determined based on the surgeon's experience.

### APPLICATION

The cement must be inserted into the bone cavity in this stage. Firmly press the material into the bone cavity. If a syringe is being used, begin extruding the cement in a uniform manner by starting in the distal area of the bone cavity and moving towards the proximal area. This application method minimizes the entrapment of air bubbles.

### INSERTION OF THE PROSTHESIS

After the cement has been inserted, the prosthesis is put into position and held in place without movement until the bone cement has polymerized. Remove any excess bone cement before it has completely hardened.

The cement's polymerization time depends on the type of cement, handling, temperature and humidity during storage and in the operating room.


**Caution! The temperature and moisture of the bone cavity accelerates the polymerization of the cement and the application of the prosthesis must therefore be completed as quickly as possible.**

### EFFECTS OF TEMPERATURE ON THE PREPARATION AND APPLICATION TIMES OF CEMEX GENTA ID GREEN BONE CEMENTS

The preparation and application of the cement is greatly influenced by the temperature during storage and in the operating room.

The effect that temperature has on the setting time of cements was evaluated in a laboratory test. Refer to the setting time vs. temperature graph that is reported. (Data obtained at controlled ambient and storage conditions subject to standard deviation). Besides temperature and humidity, other factors can

influence the cement's setting time: mixing technique (speed, use of a mixer), degree of mixing, the using of all the solid and liquid components, foreign substances in the cement (blood, saline solution, etc.), the preheating of the prosthesis.

### RECOMMENDATIONS

To improve the use of CEMEX GENTA ID GREEN bone cement:

- Use cements at a temperature of about 23°C and 60% relative humidity.
- Remove debris and thoroughly rinse the implantation site with physiological saline solution;
- Prevent as much liquid as possible from coming between the bone tissue and cement by drying the bone surface using gauze and/or suction tubes before and during cementation;

### ADVERSE EVENTS

The blood pressure often drops temporarily immediately after implanting the bone cement and the prosthesis.

The following serious and frequent adverse events may occur with the use of bone cement, but are not directly due to acrylic bone cement.

The surgeon must be aware of these events and be ready to treat them if they appear.

Serious:

- Myocardial infarction
- Cerebrovascular accident
- Cardiac arrest
- Sudden death
- Pulmonary embolism

Most frequent:

- Thrombophlebitis
  - Haemorrhage, hematoma
  - Superficial or deep wound infection
  - Trochanteric bursitis
  - Short-term cardiac conduction irregularities
- Other potential adverse events:
- Heterotopic new bone formation
  - Trochanteric separation
  - Loosening or displacement of the prosthesis
  - Breakage of the cement

### PHARMCEUTICAL PRECAUTIONS

Do not store at a temperature above 25°C and a relative humidity higher than 70%. Keep away from sunlight.

Sterility is only guaranteed if the containers are intact.  
Do not attempt to resterilize any of the components.  
Do not use the product if the powder has a yellowish or brownish colour or if the liquid is syrupy.  
These two conditions indicate that the product has not been stored correctly.

**DISPOSAL**

The disposal of the device or its components is carried out in accordance with local waste regulations. Allow the surplus bone cement to set prior to disposal.

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TECRES S.p.A.