



Weck®

## Ligation Solutions

Metal and polymer ligation systems to help you ligate with security and confidence

# Weck® Ligation Solutions

## Metal and polymer ligation systems to help you ligate with security and confidence

Beginning in Edward Weck's Manhattan storefront in 1884, Weck brand products have been there for surgeons — step-by-step — since the earliest days of modern surgery.

In 1963, the Weck brand redefined the closure market with the introduction of the Hemoclip® Traditional Metal Ligating Clip. Over the last fifty years, the Weck brand has continued to refine the science of closure to produce a line of metal and polymer ligation solutions demanded by surgeons around the world.

Whether you prefer the classic feel of metal clips or the modern flexibility of polymer, Teleflex offers ligation products designed to suit your distinct preferences and procedural needs:

- **Weck Hem-o-lok® Polymer Locking Ligation Systems** provide unparalleled security, based on an ex-vivo study, against leading competitors in three common clip failure modes\*, the ability to ligate larger vessels through a smaller port, and cool ligation. Hem-o-lok Clips are ideal for both open and laparoscopic cases.
- **Horizon™ and Hemoclip Metal Ligation Systems** satisfy your need for traditional metal ligating clips for use in cardiovascular and general surgical applications. Available in a range of sizes, surgeons have relied on Weck Metal Clips for generations.



\*Data on file (2013 internal study), Teleflex Incorporated, Report #D001591. Testing conducted on porcine carotids, sample size = 33, p ≤ 0.05. Clinical performance cannot be extrapolated from the data. Testing pressures range beyond physiological pressures. The Weck® Hem-o-lok® medium/large clip was compared to Ethicon® LIGAMAX™ 5, Ethicon LIGACLIP® ERCA 10mm, and Covidien Endo Clip™ III. The Weck Hem-o-lok large clip was compared to Ethicon LIGACLIP ERCA Large Clips and Covidien Endo Clip™ L.

## Table of Contents

4-7

**WECK® HEM-O-LOK® POLYMER LOCKING LIGATION SYSTEM**  
Secure locking polymer technology

8-11

**WECK® HORIZON™ METAL LIGATION SYSTEM**  
Precise easy loading titanium clips

12-15

**WECK® HEMOCLIP® METAL LIGATION SYSTEMS**  
Classic metal clip ligation

16-21

**RESOURCES**  
Clip Loading Methods, Warranty and Repair, Applier Care and Maintenance, Applier Management, Product Use Information

# Weck® Hem-o-lok® Polymer Locking Ligation System

## Secure-locking polymer technology

### Confidence you can feel

Confidence. That's what the Weck® Hem-o-lok® Polymer Locking Ligation System from Teleflex provides. Confidence created by a ligation system designed for clip security — firmly seated in the applier, securely closed around the patient vessel.

### Security by design

Based on an ex-vivo study, the Hem-o-lok System delivers unparalleled ligation security compared to leading competitors in three common clip failure modes,\* with the Weck Hem-o-lok Clip's integrated ridges and proprietary distal locking mechanism.

### Polymer flexibility

The Hem-o-lok Clip's non-absorbable polymer composition and flexible hinge allow it to pass through a smaller port and provide a larger distal opening than comparable metal clips. As a result, Hem-o-lok Clips can ligate up to a 10 mm structure through a 5 mm port, and up to a 16 mm structure through a 10 mm port.


### Cool ligation

Since the Hem-o-lok System is a cool ligation system, there's no chance of thermal spread to vital structures, a common concern of many surgeons when using energy-based ligation solutions.

### Size ranges to meet your needs

The Hem-o-lok System is available in a range of clip sizes: Medium, Medium-Large, Large and Extra-Large. Our broad range of applier lengths (including 45 cm endoscopic appliers for extended reach) allows surgeons to deliver clips with control and confidence.

## WECK® HEM-O-LOK® POLYMER LIGATION CLIP FEATURES



\*Data on file (2013 internal study), Teleflex Incorporated, Report #D001591. Testing conducted on porcine carotids, sample size = 33, p ≤ 0.05. Clinical performance cannot be extrapolated from the data. Testing pressures range beyond physiological pressures. The Weck® Hem-o-lok® medium/large clip was compared to Ethicon® LIGAMAX™ 5, Ethicon LIGACLIP® ERCA 10mm, and Covidien Endo Clip™ III. The Weck Hem-o-lok large clip was compared to Ethicon LIGACLIP ERCA Large Clips and Covidien Endo Clip™ L

Hem-o-lok Ligating Clips are not intended for use as a fallopian contraceptive tubal occlusion device.

Hem-o-lok Ligating Clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.

## MANUAL LOAD APPLIERS

### HEM-O-LOK OPEN APPLIERS

- Open appliers available in Medium, Medium-Large, Large, Extra-Large
- Available in 8" and 11" lengths and curved or right-angle jaws
- Opens at box lock for easy, thorough cleaning



544180



544179




544240

### HEM-O-LOK ENDOSCOPIC APPLIERS

- Endoscopic appliers available in Medium-Large, Large, Extra-Large
- Available with 32 cm or extended-length 45 cm shafts
- Available in both one-piece and take-apart designs



544965



544995T

## AUTOMATIC APPLIER

### WECK AUTO ENDO5® CLIP APPLIER


- 5 mm automatic clip applier with 15 Medium-Large polymer clips
- Medium-Large clips designed for 3-10 mm size vessels
- Clip designed for security – offers greater confidence in secure occlusion



AE05ML

### WECK AUTO ENDO10™ CLIP APPLIER

- 10 mm automatic clip applier with 12 Large polymer clips
- Large clips designed for 5-13 mm size vessels
- Long 35 mm shaft for enhanced reach




AE10LG

# Weck® Hem-o-lok® Polymer Locking Ligation System

Secure-locking polymer technology

Vessel/Tissue Bundle Size Range







HEM-O-LOK POLYMER LIGATING CLIPS

- 544220 Medium
- 544230 Medium-Large
- 544240 Large
- 544250 Extra-Large



## ORDERING INFORMATION

### Weck Hem-o-lok Polymer Locking Ligation System

CLIPS	OPEN MANUAL LOAD APPLIERS/REMOVERS	ENDOSCOPIC MANUAL LOAD APPLIERS/REMOVERS	AUTOMATIC APPLIERS
 M 544220 Hem-o-lok Medium Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile	544113 8" Curved 544114 11" Right-Angle (70°) 544115 11" Curved 544123 11" Remover (for M, ML, L)		
 ML 544230 Hem-o-lok Medium-Large Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile	544170 8" Curved 544171 11" Curved 544172 11" Right-Angle (70°) 544123 11" Remover (for M, ML, L)	<b>Endo5® ML, 5 mm diameter</b> 544965 32.5 cm applier 544965L 45 cm applier 544965T 34 cm applier, take-apart 544121 32 cm remover (for ML and L) 544121L 45 cm remover (for ML and L) 544121T 32 cm remover, take-apart (for ML and L)	<b>AE05ML 35 cm Applier</b> 15 ML clips per applier, 3 appliers per box, sterile, disposable
 L 544240 Hem-o-lok Large Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile	544179 11" Right-Angle (70°) 544180 8" Curved 544181 11" Curved 544182 10" Right-Angle (70°), narrow shaft 544123 11" Remover (for M, ML, L) 544124 11" Remover (for L, XL)	<b>Endo10™ L, 10 mm diameter</b> 544995 32 cm applier, 10 mm 544995L 45 cm applier, 10 mm 544995T 34 cm applier, 10 mm, take-apart 544121 32 cm remover, 5 mm (for ML and L) 544121L 45 cm remover, 5 mm (for ML and L) 544121T 32 cm remover, 5 mm, take-apart (for ML and L) 544130 32 cm remover, 10 mm (for L and XL) 544130L 45 cm remover, 10 mm (for L and XL) 544130T 33 cm remover, 10 mm take-apart (for L and XL)	<b>AE10LG 35 cm Applier</b> 12 L clips per applier, 3 appliers per box, sterile, disposable
 XL 544250 Hem-o-lok Extra-Large Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile	544191 11" Curved 544192 11" Right-Angle (70°) 544124 11" Remover (for L, XL)	<b>Endo10™ XL, 10 mm diameter</b> 544990 32 cm applier 544990L 45 cm applier 544990T 34 cm applier, take-apart 544130 32 cm remover (for L and XL) 544130L 45 cm remover (for L and XL) 544130T 33 cm remover, take-apart (for L and XL)	
<b>STORAGE TRAY</b>		544500 capacity: 3 endoscopic 32 cm appliers P37319 capacity: 2-3 endoscopic 45 cm appliers	

Hem-o-lok Ligating Clips are not intended for use as a fallopian contraceptive tubal occlusion device.  
 Hem-o-lok Ligating Clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.



# Weck® Horizon™ Metal Ligation System

## Precise easy-loading titanium clips

The Horizon System is the premium Weck titanium clip ligation platform, providing surgeons a strong, secure clip coupled with an efficient, convenient easy-load cartridge system.

### The Horizon Clip

The pre-formed chevron shape of the Horizon Clip engulfs tissue with precise tip-to-tip closure. The heart-shaped wire is designed to give each clip a firm grip on vessels, while the triangulated cross-section of the clip leg maximizes surface-to-surface contact between clip and jaw, practically eliminating clip fallout.

### The Horizon Cartridge

The design of the Horizon Cartridge creates a nearly frictionless one-hand load. Horizon Cartridges are color-coded to match applicators and adhesive backing allows convenient placement of the cartridge to any sterile area.

Clip images represent approximate, not exact, sizes.



MICROCLIP™



SMALL



MEDIUM



MEDIUM-LARGE



LARGE



### OPEN APPLIERS

- Available with curved (20°), angled (55°) or right-angled (70°) jaws
- Opens at box-lock for easy, thorough cleaning
- Available in a wide range of applicator lengths including 14"



### ENDOSCOPIC APPLIERS

- Medium, Medium-Large, and Large
- Extended length available (45 cm shaft)
- Large, lightweight rotation knob allows one-finger, 360° rotation of applicator shaft
- Flush ports for thorough cleaning



### CLIP OPTIONS

- The Horizon Clip Cartridge provides easy and precise loading that can enhance surgical procedure efficiency
- Both 6- and 24-clip cartridges are available



002200



002205

### SPECIALTY APPLIER

- The Horizon 11" Multi-Angle Access Open Surgical Applier for small (wide slot) clip cartridges incorporates a 45° angled jaw (or curved jaw) for an optimal angle of access to the surgical target
- The Horizon Endosaphenous Vein Applier incorporates a 5 mm shaft diameter with a 7 mm angled jaw for ligation of vein branches



137117



238200

Weck Ligating Clips are not intended for use as a contraceptive tubal occlusion device.

Weck Ligating Clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.

# Weck® Horizon™ Metal Ligation System

Precise easy-loading titanium clips



## ORDERING INFORMATION

### Weck Horizon Metal Ligation System

#### TITANIUM LIGATING CLIPS

ITEM CODE	CLIP SIZE	CLIPS/CARTRIDGE	CARTRIDGES/BOX	CLIPS/BOX
005200	MicroClip™	6	30	180
001200	Small	6	30	180
001204	Small	24	25	600
001201	Small-Wide	6	30	180
001205	Small-Wide	24	25	600
002200	Medium	6	30	180
002204	Medium	24	25	600
003200	Medium-Large	6	20	120
003204	Medium-Large	24	15	360
004200	Large	6	20	120
004204	Large	24	15	360

## ORDERING INFORMATION

### Weck Horizon Metal Ligation System

#### OPEN SURGICAL APPLIERS

ITEM CODE	APPLIER LENGTH	CLIP SIZE	JAW TYPE*
537061	6" (15 cm)	MicroClip	Curved
537081	8" (20 cm)	MicroClip	Curved
137061	6" (15 cm)	Small	Curved
137081	8" (20 cm)	Small	Curved
137085	8" (20 cm)	Small	Angled
137111	11" (28 cm)	Small	Curved
137062	6" (15 cm)	Small-Wide	Curved
137082	8" (20 cm)	Small-Wide	Curved
137086	8" (20 cm)	Small-Wide	Angled
137112	11" (28 cm)	Small-Wide	Curved
237061	6" (15 cm)	Medium	Curved
237081	8" (20 cm)	Medium	Curved
237085	8" (20 cm)	Medium	Angled
237111	11" (28 cm)	Medium	Curved
237115	11" (28 cm)	Medium	Angled
237117	11" (28 cm)	Medium	Right-Angled
237141	14" (35.5 cm)	Medium	Curved
237145	14" (35.5 cm)	Medium	Angled
337081	8" (20 cm)	Medium-Large	Curved
337085	8" (20 cm)	Medium-Large	Angled
337111	11" (28 cm)	Medium-Large	Curved
337115	11" (28 cm)	Medium-Large	Angled
337117	11" (28 cm)	Medium-Large	Right-Angled
337141	14" (35.5 cm)	Medium-Large	Curved
337145	14" (35.5 cm)	Medium-Large	Angled
437081	8" (20 cm)	Large	Curved
437085	8" (20 cm)	Large	Angled
437111	11" (28 cm)	Large	Curved
437115	11" (28 cm)	Large	Angled
437117	11" (28 cm)	Large	Right-Angled
437141	14" (35.5 cm)	Large	Curved
437145	14" (35.5 cm)	Large	Angled

#### ENDOSCOPIC LIGATING APPLIERS

ITEM CODE	APPLIER LENGTH	DIAMETER	CLIP SIZE	JAW TYPE
238110	12" (31 cm)	10 mm	Medium	Endoscopic
238170	18" (45 cm)	10 mm	Medium	Endoscopic
238200	13" (32 cm)	5 mm Shaft, 7 mm Jaw	Medium	Endosaphenous
338110	13" (32 cm)	10 mm	Medium-Large	Endoscopic
338170	18" (45 cm)	10 mm	Medium-Large	Endoscopic
438110	14" (36 cm)	12 mm	Large	Endoscopic
438170	18" (45 cm)	12 mm	Large	Endoscopic

#### ENDOSCOPIC LIGATING APPLIERS – TAKE-APART

ITEM CODE	APPLIER LENGTH	DIAMETER	CLIP SIZE	JAW TYPE
238110T	11" (28 cm)	10 mm	Medium	Endoscopic
338110T	11" (28 cm)	10 mm	Medium-Large	Endoscopic

#### SPECIALTY OPEN SURGICAL LIGATING APPLIERS

ITEM CODE	APPLIER LENGTH	CLIP SIZE	JAW TYPE
137117	11" (28 cm)	Small-Wide	45° Angled
137118	11" (28 cm)	Small-Wide	Curved

#### METAL CLIP REMOVERS (REMOVES ALL WECK METAL CLIPS)

ITEM CODE	APPLIER LENGTH	STYLE
523120	6.5" (16.5 cm)	Straight
523121	10.5" (26.7 cm)	Bent Shaft, 30° Angle

Weck Ligating Clips are not intended for use as a contraceptive tubal occlusion device.

Weck Ligating Clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.

\*Curved 20°, Angled 55°, Right-Angle 70°

# Weck® Hemoclip® Metal Ligation Systems

Classic metal clip ligation relied on by surgeons for years

## HEMOCLIP TRADITIONAL METAL LIGATING CLIPS

- Friction fit loading system
- Available in titanium and tantalum
- Cartridge base to assist loading

## HEMOCLIP PLUS SYSTEM

The Hemoclip Plus System is an evolution of the classic metal Hemoclip. Hemoclip Plus Metal Ligating Clips feature:

- An easy-load cartridge system
- Available in titanium
- Clip loads from the cartridge with minimal pressure
- Adhesive or non-adhesive clip cartridges

## CLIP OPTIONS

- Chevron-shaped clip for precise tip to tip closure
- Heart-shaped wire for tenacious grip on vessels
- Lateral transverse grooves for increased vessel contact



523130

## OPEN APPLIERS

- Available in a variety of lengths and tip angles
- Opens at box-lock for easy cleaning



Curved tip



Right angle tip



Weck ligating clips are not intended for use as a contraceptive tubal occlusion device.  
Weck ligating clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.

# Weck® Hemoclip® Metal Ligation Systems

Classic metal clip ligation relied on by surgeons for years

## ORDERING INFORMATION

### Weck Hemoclip Traditional Clips and Appliers

#### HEMOCLIP TRADITIONAL-TITANIUM LIGATING CLIPS

ITEM CODE	CLIP SIZE	CLIPS/ CARTRIDGE	CARTRIDGES/ BOX	CLIPS/BOX
523735	Small	25	12	300
523835	Small	10	24	240
523700	Medium	25	10	250
523800	Medium	10	20	200
523760	Medium/Large	10	20	200
523860	Medium/Large	10	16	160
523770	Large	15	10	150
523870	Large	10	12	120

#### HEMOCLIP TRADITIONAL-TANTALUM LIGATING CLIPS

ITEM CODE	CLIP SIZE	CLIPS/ CARTRIDGE	CARTRIDGES/ BOX	CLIPS/BOX
523135	Small	25	12	300
523335	Small	10	24	240
523100	Medium	25	10	250
523300	Medium	10	20	200
523160	Medium/Large	10	20	200
523360	Medium/Large	10	16	160
523170	Large	15	10	150
523370	Large	10	12	120

■ Small ■ Medium ■ Medium/Large ■ Large

#### HEMOCLIP TRADITIONAL LIGATING CLIP APPLIERS

ITEM CODE	CLIP SIZE	APPLIERS LENGTH	JAW TYPE
523140	Small	6" (15 cm)	Curved
523150	Small	8" (20 cm)	Curved
523151	Small	11" (28 cm)	Curved
523105	Medium	8" (20 cm)	Right Angled
523106	Medium	11" (28 cm)	Right Angled
523109	Medium	6" (15 cm)	Curved
523110	Medium	8" (20 cm)	Curved
523111	Medium	11" (28 cm)	Curved
523165	Medium/Large	8" (20 cm)	Curved
523166	Medium/Large	11" (28 cm)	Curved
523167	Medium/Large	8" (20 cm)	Right Angled
523168	Medium/Large	11" (28 cm)	Right Angled
523175	Large	8" (20 cm)	Curved
523177	Large	8" (20 cm)	Right Angled
523178	Large	11" (28 cm)	Right Angled
523180	Large	11" (28 cm)	Curved

#### HEMOCLIP TRADITIONAL LIGATING CLIP PRECISION APPLIERS

ITEM CODE	CLIP SIZE	APPLIERS LENGTH	JAW TYPE
523540	Small	6" (15 cm)	Curved
523550	Small	8" (20 cm)	Curved
523551	Small	11" (28 cm)	Curved

#### HEMOCLIP LIGATING CLIP ACCESSORIES

ITEM CODE	CLIP SIZE	APPLIERS LENGTH	JAW TYPE
523120	Removing Forceps	6" (15 cm)	Straight
523121	Removing Forceps	10" (25.4 cm)	Angled Shanks
523130	Cartridge Base	7" (17.78 cm)	

## ORDERING INFORMATION

### WECK HEMOCLIP PLUS CLIPS AND APPLIERS

#### HEMOCLIP TRADITIONAL-TITANIUM LIGATING CLIPS

ITEM CODE	CLIP SIZE	CLIPS/ CARTRIDGE	CARTRIDGES/ BOX	CLIPS/BOX
533735	Small w/tape	25	12	300
533737	Small	25	12	300
533835	Small w/tape	10	18	180
533837	Small	10	18	180
534735	Small Strongpoint w/tape	25	12	300
534737	Small Strongpoint	25	12	300
534835	Small Strongpoint w/tape	10	18	180
534837	Small Strongpoint	10	18	180
533700	Medium w/tape	25	10	250
533702	Medium	25	10	250
533800	Medium w/tape	10	18	180
533802	Medium	10	18	180
533860	Medium-Large w/tape	10	12	120
533862	Medium-Large	10	12	120
533870	Large w/tape	10	12	120
533872	Large	10	12	120

#### HEMOCLIP PLUS LIGATING CLIP APPLIERS

ITEM CODE	CLIP SIZE	APPLIERS LENGTH	JAW TYPE
533140	Small	6" (15 cm)	Curved
533150	Small	8" (20 cm)	Curved
533151	Small	11" (28 cm)	Curved
533152	Small	8" (20 cm)	Right Angled
534140	Small Strongpoint	6" (15 cm)	Curved
534150	Small Strongpoint	8" (20 cm)	Curved
533105	Medium	8" (20 cm)	Right Angled
533106	Medium	11" (28 cm)	Right Angled
533108	Medium	11" (28 cm)	Angled
533109	Medium	6" (15 cm)	Curved
533110	Medium	8" (20 cm)	Curved
533111	Medium	11" (28 cm)	Curved
533164	Medium-Large	11" (28 cm)	Angled
533165	Medium-Large	8" (20 cm)	Curved
533166	Medium-Large	11" (28 cm)	Curved
533167	Medium-Large	8" (20 cm)	Right Angled
533168	Medium-Large	11" (28 cm)	Right Angled
114527	Medium-Large	13" (32 cm)	Endoscopic
533175	Large	8" (20 cm)	Curved
533177	Large	8" (20 cm)	Right Angled
533178	Large	11" (28 cm)	Right Angled
533180	Large	11" (28 cm)	Curved
533181	Large	11" (28 cm)	Angled

#### HEMOCLIP LIGATING CLIP ACCESSORIES

ITEM CODE	CLIP SIZE	APPLIERS LENGTH	JAW TYPE
523120	Removing Forceps	6" (15 cm)	Straight
523121	Removing Forceps	10" (25.4 cm)	Angled Shanks
523130	Cartridge Base	7" (17.78 cm)	

WECK ligating clips are not intended for use as a contraceptive tubal occlusion device.  
WECK ligating clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.




# Weck® Clip Loading Methods

## Manual applier technique

### RING-HANDLED APPLIER GUIDELINES


To prevent jaw applier damage or dropped clips, the user should avoid any grasp that can compress the two shanks of the applier during loading. Applier is held at an angle appropriate to sliding the jaws into the cartridge.

#### Preferred Loading Methods




**Method #1 Shank and Spring**

One shank and the spring are held between the thumb and index finger, which is placed between shafts to spread jaws and stabilize loading.



**Method #2 Pencil Grip**


Box-lock is held between the thumb and side of index finger, with tips of index finger and other fingers supporting box-lock area for stability. The index finger is behind the box-lock, not beside it.




**Method #3 Angled Appliers**

Thumb should be placed on top of the box-lock while other fingers support one side of the shaft.

#### What to Avoid




Grasping the applier at both rings tends to force applier tips together and “fights” the loading action. This may cause the applier to twist or rock, allowing jaw damage or misloaded clips.



Holding both shanks of applier, even near the box-lock, may also tend to force jaws together during loading. Keep the index finger behind the box-lock or shank.

### RING HANDLE AND ENDOSCOPIC APPLIER CLIP LOADING

#### Proper Jaw Angle



Applier jaws need to be positioned perpendicular to the cartridge body for slide-in of the jaw into the cartridge seat. It is necessary to set the shanks of the applier (the main body of the applier) at an acute angle to the cartridge.

# Warranty/Repair

Send all appliers in need of repair or refurbishment to:

Teleflex® Applier Repairs  
4620-A Industry Lane  
Durham, NC 27713

For applier repairs outside the U.S., contact your local distributor.

### Open Surgical Appliers

All Weck Ring Handled Reusable Open Surgical Appliers are covered by a warranty. Teleflex will repair, or replace, free of charge, these appliers as long as the appliers are used as intended for application of Weck Clips and have not been repaired by unauthorized personnel.

### Endoscopic Surgical Appliers

Weck Manual Load Endoscopic Appliers are covered by a three year warranty on defects in materials and workmanship and a three year repair service on tip alignment and refurbishment as long as the appliers are used as intended for application of Weck Clips and have not been repaired by unauthorized personnel.

Teleflex | Toll Free: 866.246.6990 | [www.teleflex.com](http://www.teleflex.com)

# Care and Maintenance

## Weck® Open and Endoscopic Clip Appliers

### CLEANING

Follow hospital, Association of the Advancement of Medical Instrumentation (AAMI), Association of Perioperative Registered Nurses (AORN), or regional standards for cleaning your appliers.

#### Open Appliers

Opening the box lock:

- Lightly squeeze shanks together with one hand
- Lift the spring off the flag
- Open the applier



#### Endoscopic Appliers — One Piece

Endoscopic appliers are provided with a cleaning port, which permits access to the interior channels and cavities. Uncap the cleaning port and flush the applier with warm distilled or filtered water for approximately two minutes or until the visible gross debris is removed from the device.



#### Lubrication

Lubrication is essential every time instruments are processed. Do not use mineral oil, petroleum, or silicone-based products. Use a non-silicone water-soluble lubricant prior to sterilization.

### Sterilization

STERILIZATION METHOD	INSTRUMENT CONFIGURATION	TEMPERATURE	EXPOSURE TIME (MINIMUM)	DRYING TIME (MINIMUM)
US PARAMETERS				
Pre-vacuum	Wrapped	270°F (132° C)	4 minutes	10 minutes
Pre-vacuum	Wrapped	275°F (135° C)	3 minutes	10 minutes
Gravity Displacement	Wrapped	275°F (135° C)	10 minutes	10 minutes
NON-US PARAMETERS				
Pre-vacuum	Wrapped	273°F (134° C)	3 minutes	10 minutes
Pre-vacuum	Wrapped	279°F (137° C)	3 minutes	10 minutes
Gravity Displacement	Wrapped	270°F (132° C)	10 minutes	10 minutes


### INSPECTION

Please examine endoscopic and open appliers for potential damage before use in surgery. Pay particular attention to the jaws. Damaged or misaligned jaws may not allow clips to close acceptably for occlusion of intended structure.

#### Open and Endoscopic Appliers

##### Alignment:

Always check the alignment of the applier jaws before use. When closed, jaw tips should be directly aligned and not offset.





#### Open Appliers

Check the following components of your open clip appliers:



##### Spring:

Examine to ensure spring is not loose, bent, or rusted



##### Flag Condition:

Examine to ensure flag is not missing or damaged



## Applier Management

### Two convenient options to manage your open appliers

An essential part of accurate ligation is a finely crafted applier, designed specifically for Weck® Clips and maintained to the highest standards. We are proud to provide, service, maintain and replace these appliers for you — ensuring that each clip is applied by a device designed specifically for that clip.

#### TWO WAYS TO OBTAIN HORIZON™, HEM-O-LOK®, HEMOCLIP® PLUS AND TRADITIONAL APPLIERS:

##### Option 1 — Applier Purchase

All Weck Appliers are available for purchase. Purchased open appliers will be maintained, repaired or replaced for as long as you own them, at no charge, through Teleflex® as long as the appliers are used as intended for application of Weck Clips and have not been repaired by unauthorized personnel. Contact your local Surgical Sales Representative for more details.

##### Option 2 — Applier Loaner Program

Hospitals purchasing Weck Clips can take advantage of the Teleflex Applier Loaner Program. Weck Open Metal Clip Appliers are loaned to the hospital for use with Weck Clips. Teleflex will maintain, repair or replace all appliers needed to be sufficient for the volume of Weck Clips being used. These appliers will be provided to the hospital at no charge for exclusive use with Weck Clips. If the hospital discontinues its use of Weck Clips, then the appliers are to be returned to Teleflex (Signed Loaner Agreement required).



## Product Use Information

Weck Metal and Polymer Clips are intended for use in procedures involving ligation of vessels or tissue structures. Surgeons should apply the appropriate size clip for the size of the vessel or tissue structure to be ligated.

Weck Ligating Clips are contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomy. In procedures other than laparoscopic donor nephrectomy, Teleflex recommends ligation of the renal artery with more than one clip on the patient side with a minimum distal renal artery cuff of 2–3 mm beyond the distal clip. Application of a second clip on all other vessels other than the renal artery should be dictated by the surgeon's judgment. Security of the closure should be confirmed after ligation.

Always check the alignment of the applier jaws before use. When closed, jaw tips should be directly aligned and not offset. Alignment of the jaw is critical for safe application of the clip. If this is not done, patient injury may occur. Follow the Instructions for Use for proper maintenance, care and cleaning necessary to ensure proper functionality.

Do not attempt to close the jaws on a vessel or anatomic structure without a clip properly loaded into the jaws. Closure of empty jaws on a vessel or anatomic structure may result in patient injury. Before applying a clip, verify the structural size and condition of the vessel or structure and use the proper size clip. Ligating clip systems differ in closure characteristics according to clip design and other variables. It is the responsibility of the user to select structures for the application of clips and confirm secure grip of the clips after placement, and after the use of other surgical devices in the immediate area of the application.

Your Weck Applier has been designed and calibrated for use with specific Hem-o-lok, Horizon, Hemoclip Plus or Hemoclip Traditional Ligating Clips as indicated. If applicable, applier color coding matches the color of the ligating clip cartridge for which it is to be used. Weck does not assume responsibility for unsatisfactory results caused by the use of any equipment or clips not specifically identified by Weck as an integral part of that specific system.

For clip and associated applier information, please refer to the Instructions for Use included with the applier.





For more information visit [teleflex.com/ligation](https://teleflex.com/ligation)



Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit [teleflex.com](http://teleflex.com).

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüschi® and Weck® – trusted brands united by a common sense of purpose.

### Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

### Regional Offices

**United States:** Phone +1 919 544 8000, Toll Free +1 866 246 6990, [cs@teleflex.com](mailto:cs@teleflex.com), 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

**Latin America:** Phone +1 919 433 4999, [la.cs@teleflex.com](mailto:la.cs@teleflex.com), 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

**International:** Phone +353 (0)9 06 46 08 00, [orders.intl@teleflex.com](mailto:orders.intl@teleflex.com), Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

**Australia/New Zealand** 1300 360 226

**Austria** +43 (0)1 402 47 72

**Belgium** +32 (0)2 333 24 60

**Canada** +1 (0) 905 943 9000

**China (Shanghai)** +86 (0)21 6163 0965

**China (Beijing)** +86 (0)10 6418 5699

**Czech Republic** +420 (0)495 759 111

**France** +33 (0)5 62 18 79 40

**Germany** +49 (0)7151 406 0

**Greece** +30 210 67 77 717

**India** +91 (0)44 2836 5040

**Italy** +39 0362 58 911

**Japan** +81 (0)3 6632 3600

**Korea** +82 2 536 7550

**Mexico** +52 55 5002 3500

**Netherlands** +31 (0)88 00 215 00

**Portugal** +351 22 541 90 85

**Singapore** +65 6439 3000

**Slovak Republic** +421 (0)3377 254 28

**South Africa** +27 (0)11 807 4887

**Spain** +34 918 300 451

**Switzerland** +41 (0)31 818 40 90

**United Kingdom** +44 (0)1494 53 27 61

For more information, please visit [teleflex.com](http://teleflex.com).

Rx Only. Caution: U.S. federal law restricts this device to sale by or on the order of a physician. Hem-o-lok Ligating Clips are contraindicated for use as a fallopian contraceptive tubal occlusion device and contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.

Teleflex, the Teleflex logo, Arrow, Deknatel, Endo5, Endo10, Hemoclip, Hem-o-lok, Horizon, Hudson RCI, LMA, MicroClip, Pilling, Rüschi and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks or registered trademarks appearing herein are the property of their respective owners.

Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2016). Subject to technical changes without further notice.

© 2016 Teleflex Incorporated. All rights reserved. MC-002140

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

**Teleflex Medical**  
**IDA Business and Technology Park**  
Dublin Road  
Athlone  
Co. Westmeath  
Ireland

Holds Certificate Number:

FM 544574

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

The design and manufacture of non-active digestive tract devices; non-active gynaecological devices, non-active regional anaesthesia devices, non-active respiratory devices, non-active surgical devices, non-active urology devices and active surgical devices.

For and on behalf of BSI:

Graeme Tunbridge, Senior Vice President Medical Devices

Original Registration Date: 2009-03-09

Latest Revision Date: 2023-01-26

Effective Date: 2023-02-12

Expiry Date: 2026-02-11

Page: 1 of 1



...making excellence a habit.™

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

**No.****CE 540595**

## Issued To:

**Teleflex Medical  
IDA Business and Technology Park  
Dublin Road  
Athlone  
Co. Westmeath  
Ireland**

In respect of:

**The design and manufacture of non active digestive tract devices; non active gynecological devices; non active regional anaesthesia devices; non active respiratory devices; non active surgical devices; non active urology devices.**

on the basis of our examination of the quality assurance system under the requirements of Council Directive 93/42/EEC, Annex II excluding section 4. The quality assurance system meets the requirements of the directive. For the placing on the market of class III products an Annex II section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):



Gary E Slack, Senior Vice President Medical Devices

First Issued: **2009-01-13**Date: **2020-06-09**Expiry Date: **2024-05-26**

...making excellence a habit.™

Page 1 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

# EC Certificate - Full Quality Assurance System

## Supplementary Information to CE 540595

Issued To:

**Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

Number	Device Name	Intended purpose per IFU
<b>Class III</b>		
---	EpiStar CSE - Spinal-Epidural Anaesthesia Kits	See CE 544836
---	Spinostar Spinal Needles	See CE 560441
<b>Class IIb</b>		
10735	Sterile Percutaneous Nephrostomy Catheter	Puncture and dilation of percutaneous approaches into the upper urinary tract.
35404	Sterile Tracheostomy Tube	Cannulation of tracheostomised patients through an existing tracheostoma.
14099	Sterile Tracheostomy Tube	Cannulation of tracheostomised patients, in whom the stoma was created by percutaneous dilative tracheostomy.
58005	Sterile Ureter Stent	Routine drainage of the renal pelvis via the ureter or a ureter-skin stoma to an external collection site.
34924	Sterile Suprapubic Cystotomy Set	Routine suprapubic drainage of the bladder
31074	Sterile Ureterocutaneostomy Catheter	Routine drainage of urine through a ureterocutaneous stoma site.

First Issued: **2009-01-13**

Date: **2020-06-09**

Expiry Date: **2024-05-26**

...making excellence a habit.™

Page 2 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.



# EC Certificate - Full Quality Assurance System

## Supplementary Information to CE 540595

Issued To:

**Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

Number	Device Name	Intended purpose per IFU
<b>Class IIa</b>		
MD 0106	Sterile Transurethral Catheter	---
MD 0101	Sterile Tracheostomy Retainer Set	---
MD 0106	Sterile Rectal Tube	---
MD 0101 MD 1102	Sterile Breathing Circuit	---
MD 0101 MD 1102	Non-sterile Breathing Circuit	
MD 0101	Sterile Cricothyrotomy Set	---
MD 0102	Sterile Epidural Set	---
MD 0101	Sterile EZ Blocker Kit	---
MD 0106	Sterile Guidewire	---
MD 0106	Sterile Kidney Stone Extractor	---
MD 0101	Sterile Tracheal Tube	---
MD 0101	Non-sterile Tracheal Tube	---
MD 0101	Sterile Laryngeal Mask	---
MD 0101	Non-sterile Laryngeal Mask	---

First Issued: **2009-01-13**

Date: **2020-06-09**

Expiry Date: **2024-05-26**

...making excellence a habit.™

Page 3 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

# EC Certificate - Full Quality Assurance System

## Supplementary Information to CE 540595

Issued To:

**Teleflex Medical  
IDA Business and Technology Park  
Dublin Road  
Athlone  
Co. Westmeath  
Ireland**

Number	Device Name	Intended purpose per IFU
<b>Class IIa</b>		
MD 0106	Sterile Laparoscopy Bag	---
MD 0101	Sterile Bronchial Tube	---
MD 0101	Sterile Suprapubic Cystotomy Set	---
MD 0303	Sterile Drainage Tube	---
MD 0101	Sterile Tracheostomy Tube, Inner cannula	---
MD 0106	Sterile Ureter Catheter	---
MD 0102	Sterile Needle Introducer	---
MD 0101	Sterile Percutaneous Nephrostomy Catheter	---
MD 0106	Non-sterile Temperature Sensor	---
MD 0101	Sterile Breathing Bag	---
MD 0101 MD 0106	Sterile Irrigation System for Ureterocutaneostomy	---

First Issued: **2009-01-13**

Date: **2020-06-09**

Expiry Date: **2024-05-26**

...making excellence a habit.™

Page 4 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
Date: **2020-06-09**  
Issued To: **Teleflex Medical  
IDA Business and Technology Park  
Dublin Road  
Athlone  
Co. Westmeath  
Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
Arrow International CR, a.s. Jamska 2359/47 Zdar Nad Sazavou 59101 Czech Republic	<b>Design Manufacture</b>
Arrow International CR, a.s. Prazska 209 Hradec Kralove 50004 Czech Republic	<b>Design</b>
Arrow Medical Ltd Hatton Garden Industrial Estate Kington Hereford HR5 3RB United Kingdom	<b>Manufacture</b>

...making excellence a habit.™

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
BBF Sterilisationsservice GmbH Willy-Rüsch-Straße 10/1 71394 Kernen Germany	<b>Radiation (Gamma Sterilization)</b>
Chelle Medical Limited Le Rocher P.O Box 221 Victoria Mahe Seychelles	<b>Manufacture</b>
Chemiczna Spółdzielnia Pracy Technochemia ul. Fabryczna 3 05-600 Grójec Poland	<b>ETO Sterilization</b>
Contract Medical International, spol. sr.o. Vážní 848 500 03 Hradec Králové Czech Republic	<b>Manufacture</b>

...making excellence a habit.™



# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
Daqing Medical Device (Tianjin) Co., Ltd 10A & 11A Tianzhi Industrial Center No.12 Hong Yuan Road Xiqing Economic Development Area 300385 Tianjin People's Republic of China	<b>Manufacture</b>
Degania Silicone Limited Kibbutz 1513000 Degania Bet Israel	<b>Manufacture</b>
Forefront (Xiamen) Medical Devices Co., Ltd No. 28 Haijing East Road & No. 61 Haijing South Road Xiamen Area of China (Fujian) Pilot Free Trade Zone 361026 Xiamen, Fujian People's Republic of China	<b>Manufacture</b>

...making excellence a habit.™

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
Forefront Medical Technology (Pte) Ltd 35 Joo Koon Circle Singapore 629110 Singapore	<b>Manufacture</b>
M.E.M., Inc. 8 Bishop Lane Madison Connecticut 06443 USA	<b>Manufacture</b>
Medicoplast International GmbH Heusweilerstrasse 100 DE-66557 Ilingen Germany	<b>ETO Sterilization</b>
Parker Hannifin CSS Merrillville 1201 East 86th Place Merrillville IN, 46410 United States	<b>Manufacture</b>

...making excellence a habit.™

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
Date: **2020-06-09**  
Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
Plaxtron Industrial (M) Sdn. Bhd. Plot 28, Kawasan Perusahaan Jelapang II Zon Perdagangan Bebas, Ipoh Perak 30020 Malaysia	<b>Manufacture</b>
Professional Contract Sterilization Inc. 40 Myles Standish Blvd Taunton Massachusetts 02780-1026 USA	<b>ETO Sterilization</b>
safemed medical devices s.r.o Trabantská 292 19015 Praha 9 Czech Republic	<b>Manufacture</b>

...making excellence a habit.™

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
sfm medical devices GmbH Brückenstraße 5 63607 Wächtersbach Germany	<b>ETO Sterilization Manufacture</b>
SINA-SterilGamma Sdn. Bhd. LOT 88077, Jalan Perigi Nenas 7/1 Taman Perindustrian Pulau Indah 42907 Pelabuhan Klang, Selangor Malaysia	<b>ETO Sterilization</b>
SP Medical A/S Møllevej 1 4653 Karise Denmark	<b>Design Manufacture</b>
SP Medical Sp. z o.o. Ul. Ceramiczna 2K 98-220 Zduńska Wola Poland	<b>Manufacture</b>

...making excellence a habit.™



# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
STERIS AST CZ s.r.o. Prumyslová Zona Kosikov 80 Velka Bites 59501 Czech Republic	<b>ETO Sterilization</b>
Synergy Sterilisation (M) Sdn Bhd. Plot 203 Kuala Ketil Industrial Estate Kuala Ketil Kedah 09300 Malaysia	<b>ETO Sterilization</b>
Synergy Sterilisation Kulim (M) Sdn. Bhd Lot 71, Kulim Industrial Estate Kulim Kedah 09000 Malaysia	<b>ETO Sterilization</b>

...making excellence a habit.™

# EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

## List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

<b>Subcontractor:</b>	<b>Service(s) supplied</b>
Teleflex Medical Asia Pte. Ltd. 21 Merchant Road #04-01 Royal Merukh S.E.A 058267 Singapore	<b>Design</b> <b>Manufacture</b>
Teleflex Medical Sdn. Bhd. Lot PT 2577, Jalan Perusahaan 4 34600 Kamunting Perak Malaysia	<b>Design</b> <b>ETO Sterilization</b> <b>Manufacture</b>
The Laryngeal Mask Company (Malaysia) Sdn. Bhd. Lot 19 & 1920 Industrial Zone Phase 1 Kulim Hi-Tech Park Kulim 09000 Malaysia	<b>Manufacture</b>

...making excellence a habit.™

# EC Certificate - Full Quality Assurance System Certificate History

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road  
 Athlone  
 Co. Westmeath  
 Ireland**

Date	Reference Number	Action
13 January 2009	7245725	First issue.
17 March 2009	7325719	Company address amended. Extension to scope. Addition of Willy Rüsç, Germany as subcontractor for design and manufacture.
25 August 2009	7399879	Addition of 'epidural catheter Epistar and Epistar CSE' to scope. Addition of SFM as significant subcontractor for manufacture. Addition of 'design' to services supplied by Teleflex Medical Malaysia, Arrow International CR, a.s. and Arrow International Inc., Czech Republic.
11 November 2009	7455515	Addition of CeMed GmbH for manufacturing to the list of significant subcontractors.
20 April 2010	7497906	Laryngeal Mask added to scope. Addition of Tianjin Medis Medical Device Co. Ltd as significant subcontractor for manufacture.
08 September 2010	7558508	Scope reworded in accordance with generic device groups. Certificate renewal.
23 May 2012	7778467	Correction of significant subcontractor address and addition of new scope activities for subcontractors.

# EC Certificate - Full Quality Assurance System Certificate History

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

Date	Reference Number	Action
04 February 2013	7932588	The addition of a significant subcontractor SP Medical A/S.
14 May 2014	8134266	Addition of peripheral angioplasty balloon catheters to product family, covered by scope expression 'non-active surgical devices'. Addition of significant subcontractors Hotspur Technologies, Inc and Teleflex Medical Asia Pte Ltd.
09 March 2015	8293488	Addition of 8 crucial suppliers.
28 August 2015	8406490	Certificate renewal. Removal of Hotspur Technologies, Inc. from list of significant subcontractors.
05 August 2016	8571081	Addition of Contract Medical International, spol. sr.o. to the list of significant subcontractors. Addition of EZ Blocker non-active respiratory device.
09 January 2017	8665617	Change to the address of subcontractor (Forefront).
16 July 2018	8939923	Addition of Daqing Medical Device (Tianjin) Co., Ltd to the list of significant subcontractors.



# EC Certificate - Full Quality Assurance System Certificate History

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

Date	Reference Number	Action
4 March 2019	7779566	Traceable to NB 0086.
Current	3124666	<p>Certificate renewal.</p> <p>Addition of supplementary product information table.</p> <p>Removal of Control of Sterilization from Service(s) supplied for Arrow International CR, a.s. (Zdar), Arrow International CR, a.s. (Hradec Kralove), Contract Medical International spol. sr.o., SP Medical A/S, sfm medical devices GmbH, Teleflex Medical Asia Pte. Ltd. and Teleflex Medical Sdn. Bhd.</p> <p>Removal of Crucial Supplier from Service(s) supplied for Arrow Medical Ltd, Chelle Medical Limited, Forefront (Xiamen) Medical Devices Co., Ltd, Forefront Medical Technology (Pte) Ltd, Parker Hannifin CSS Merrillville, Plaxtron Industrial (M) Sdn. Bhd. and The Laryngeal Mask Company (Malaysia) Sdn. Bhd.</p> <p>Addition of Manufacture to Service(s) supplied for Arrow Medical Ltd, Chelle Medical Limited, Forefront (Xiamen) Medical Devices Co., Ltd, Forefront Medical Technology (Pte) Ltd, M.E.M., Inc., Parker Hannifin CSS Merrillville, Plaxtron Industrial (M) Sdn. Bhd., and The Laryngeal Mask Company (Malaysia) Sdn. Bhd.</p> <p>Removal of Manufacture from Service(s) supplied for Arrow International CR, a.s. (Hradec Kralove)</p> <p>Addition of Degania Silicone Limited, safemed medical devices s.r.o and SP Medical Sp. z.o.o. as subcontractors for Manufacture.</p>

...making excellence a habit.™

Page 3 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.  
 This certificate was issued electronically and is bound by the conditions of the contract.

# EC Certificate - Full Quality Assurance System Certificate History

Certificate No: **CE 540595**  
 Date: **2020-06-09**  
 Issued To: **Teleflex Medical**  
**IDA Business and Technology Park**  
**Dublin Road**  
**Athlone**  
**Co. Westmeath**  
**Ireland**

Date	Reference Number	Action
		Addition of STERIS AST CZ s.r.o., Synergy Sterilisation (M) Sdn Bhd., Synergy Sterilisation Kulim (M) Sdn Bhd., Chemiczna Spółdzielnia, Medicoplast International GmbH, Professional Contract Sterilization Inc., SINA-SterilGamma Sdn Bhd and Teleflex Medical Sdn. Bhd. as subcontractors for ETO Sterilization. Addition of BBF Sterilisationservice GmbH as subcontractor for Gamma Sterilization. Removal of CeMed GmbH, Tianjin Medis Medical and Willy Rüschi GmbH Administrative correction of details for Arrow Medical Ltd, Chelle Medical Limited, Contract Medical International spol. sr.o., Daqing Medical Device (Tianjin) Co., Ltd, Forefront (Xiamen) Medical Devices Co., Ltd and SP Medical A/S. Change of address for Teleflex Medical Asia Pte. Ltd. Name change from Süddeutsche Feinmechanik GmbH (SFM) to sfm medical devices GmbH Name change from Parker Medical Systems Division - Merrillville to Parker Hannifin CSS Merrillville

...making excellence a habit.™

Page 4 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

**ARTICLE 120 SELF-DECLARATION**  
 EU MDD DECLARATION OF CONFORMITY GLOBAL ADDENDUM  
 FOR EXTENDED TRANSITION TO MDR

<b>Legal Manufacturer (LM)</b>	<u>Name:</u> Teleflex Medical Inc.  <u>Address:</u> 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA
<b>Authorized Representative</b>	<u>Name:</u> Teleflex Medical  <u>Address:</u> IDA Business and Technology Park Dublin Road Athlone Co. Westmeath Ireland
<b>Incoming Notified Body (MDR NB)</b>	<u>Name:</u> BSI Netherlands  <u>Identification Number:</u> CE 2797
<b>Notified Body That Issued the Certificate Under MDD (MDD NB)</b>	<u>Name:</u> SGS Belgium NV  <u>Identification Number:</u> CE 1639

Teleflex Medical Inc. declares that the product(s) listed in Appendix A meet the provisions of Regulation (EU) 2023/607 of the European Parliament and of the Council of 15 March 2023 amending Regulations (EU) 2017/745 (MDR) as regards the transitional provisions for certain medical devices.

This declaration is made on the basis that the product(s) listed in Appendix A are currently compliant to:

- Council Directive 93/42/EEC dated 14 June 1993, as amended by 2007/47/EC (MDD), the provisions of Article 120(3c), and are intended to or have been confirmed to meet the provisions of Regulation (EU) 2017/745 (MDR)
- Council Directive 93/42/EEC dated 14 June 1993, as amended by 2007/47/EC (MDD), the provisions of Article 120(3c), Article 120(2) points (a) or (b) of MDR and are intended to meet the provisions of Regulation (EU) 2017/745 (MDR).
- Council Directive 93/42/EEC dated 14 June 1993, as amended by 2007/47/EC (MDD), and Article 120(3c) of MDR

Furthermore, Teleflex Medical Inc. declares:


- A formal application has been lodged in accordance with Section 4.3, first subparagraph, of Annex VII MDR before MDD certificate expiration or before 26 May 2024. The LM and MDR NB have signed a written agreement in accordance with Section 4.3, second subparagraph, of Annex VII MDR before 26 September 2024.

A Quality Management System (QMS) will be implemented in compliance with Article 10(9) of Regulation (EU) 2017/745 (MDR) prior to 26 May 2024. Documentation regarding this QMS has been included as part of the application for conformity assessment made to the MDR NB. The MDR NB will complete an assessment of the QMS as part of its conformity assessment activities.

Post-market surveillance, market surveillance, vigilance, and economic operator requirements will be conducted pursuant to Article 120(3e) of the Regulation (EU) 2017/745 (MDR).

The extended transitional period for the product(s) listed in Appendix A shall end on the dates stated in Appendix A, as set forth in Article 120 (3) of the MDR.

### Approvals

<b>Name and Title of Approver:</b>	Kim Campbell, Sr. Regulatory Affairs Manager
<b>Signature of Approver:</b>	
<b>Date Approved:</b>	11-Jul-2023
<b>Place of Issue:</b>	Teleflex Medical Inc. 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

### Change History for Declaration of Conformity Addendum:

Revision	Date	Change Order Number or N/A	Reason for Revision
00	See Agile	DCO-068469	Initial Release through DCO-068469



**APPENDIX A**

<b>Product Name</b>	<b>Product Code(s)</b>	<b>MDD EC Certificate(s) No</b>	<b>MDD EC Certificate(s) Expiry Date</b>	<b>End date for Transition Period</b>	<b>MDR Product Class</b>	<b>MDR Class Rule</b>	<b>MDR Conformity Assessment Route(s)</b>	<b>Scheduled MDR Submission Date</b>
Metal Manual Ligation System (Metal Ligating Clips)	001200	BS EN ISO 13485:2016 – US97/10878  European Directive (Design Examination) 93/42/EEC – US19/819943634	14-JUL-2023	31-DEC-2028	Class III	Rule 8	Annex IX	30-Aug-2021
	001204							
	001201							
	001205							
	002200							
	002204							
	003200							
	003204							
	004200							
	004204							
	005200							
	523100							
	523135							
	523160							
	523170							
	523300							
	523335							
	523360							
	523370							
	523700							
	523735							
	523760							
	523770							
	523800							
	523835							
	523860							
523870								
533700								
533702								
533735								
533737								
533800								
533802								
533835								
533837								
533860								
533862								
533870								

	533872 534735 534737 534835 534837							
--	--	--	--	--	--	--	--	--

EC Design Examination Certificate: Certificate US19/819943634

# Teleflex Medical

3015 Carrington Mill Blvd., Morrisville, NC, 27560, United States

Device Identification:

**Sterile and non-sterile Hemoclip Traditional, Hemoclip Plus,  
Horizon and Vesocclude Metal Ligation Clips**

Intended Purpose of Device:

**Sterile and non-sterile ligating clips intended for use in procedures  
that require the ligation of vessels or anatomic structures,  
including cardiac and neurological procedures**

has been assessed and certified as meeting the requirements of

## Directive 93/42/EEC

on Medical Devices, Annex II, Section 4

It is certified that the manufacturer's design dossier (and product, where applicable) for the above device has been examined and, based on the evidence submitted, it is considered that the device conforms to the relevant Essential Requirements of EC Directive 93/42/EEC.

This certificate is issued in conjunction with a certificate covering the full quality assurance system to Annex II, which must be subject to satisfactory surveillance audits.

This certificate is valid from 01 February 2020 until 25 October 2023  
Issue 1. Certified since 25 October 2018  
and first certified by SGS Belgium NV since 01 February 2020

Certification is based on report number WW/PCI 606960 dated 20 August 2018

Addenda to that report have been issued on the following dates:

Addendum Date  
06 January 2020

Reason for Addendum  
Addition of Vesocclude Metal Ligation Clips

Authorised by

**SGS Belgium NV, Notified Body 1639**

SGS House Noorderaan 87 2030 Antwerp Belgium  
t +32 (0)3 545-48-48 f +32 (0)3 545-48-49 www.sgs.com

LPMD5007 - Certificate CE1639 Annex II-4\_EN rev. 02

Page 1 of 1





Certificate US97/10878.02

The management system of

# Teleflex Medical

375 Forbes Blvd, Mansfield, MA, 02048-1805, United States

has been assessed and certified as meeting the requirements of

## ISO 13485:2016 EN ISO 13485:2016

For the following activities

**The scope of registration appears on page 2 of this certificate.**

This certificate is valid from 13 July 2021 until 14 July 2024  
and remains valid subject to satisfactory surveillance audits.  
Recertification audit due a minimum of 60 days before the expiration date.  
Issue 2. Certified since 26 September 2000

Multiple certificates have been issued for this scope  
The main certificate is numbered US97/10878.00

The validity of this certificate depends on the validity of the main certificate.



Authorised by



0005

SGS United Kingdom Ltd  
Rossmore Business Park, Ellesmere Port, Cheshire CH65 3EN UK  
t +44 (0)151 350-6666 f +44 (0)151 350-6600 [www.sgs.com](http://www.sgs.com)

21HC 13485 2016 0421 M2

Page 1 of 2



This document is issued by the Company subject to its General Conditions of Certification Services accessible at [www.sgs.com/terms\\_and\\_conditions.htm](http://www.sgs.com/terms_and_conditions.htm). Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at <http://www.sgs.com/en/certified-clients-and-products/certified-client-directory>. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.



# Teleflex Medical

## ISO 13485:2016

## EN ISO 13485:2016



Issue 2


Detailed scope

**Design, development, manufacture and distribution of reusable medical and surgical instruments for general and specialty use; sterile and non-sterile disposable surgical, urology, anaesthesia and respiratory medical devices, sterile disposable electrosurgical medical devices. Design of Non-Sterile Nasal and Oral Mucosal Devices. Design and development of sterile single use absorbable and non-absorbable sutures, pledgets and suture guides and manufacturing of non-sterile absorbable and non-absorbable suture material.. Manufacturing of sterile single use absorbable and non-absorbable sutures.**

**Distribution of sterile single use absorbable and non-absorbable sutures and non-sterile suture material. Distribution of medical devices for endoscopy; fiber optic illuminators; sterile single use instruments for cardiovascular and general surgical procedures.**



0005





Certificate US97/10878.00

The management system of

# Teleflex Medical

3015 Carrington Mill Blvd., Morrisville, NC, 27560, United States

has been assessed and certified as meeting the requirements of

## ISO 13485:2016 EN ISO 13485:2016

For the following activities

**The scope of registration appears on page 2 of this certificate.**

This certificate is valid from 15 July 2021 until 14 July 2024 and remains valid subject to satisfactory surveillance audits. Recertification audit due a minimum of 60 days before the expiration date. Issue 22. Certified since 26 September 2000

Multiple certificates have been issued for this scope  
The main certificate is numbered US97/10878.00  
This is a multi-site certification.  
Additional site details are listed on the subsequent page.

Authorised by



0005

SGS United Kingdom Ltd  
Rossmore Business Park Ellesmere Port Cheshire CH65 3EN UK  
t +44 (0)151 350-6666 f +44 (0)151 350-6600 [www.sgs.com](http://www.sgs.com)

21HC 13485 2016 0421 M2

Page 1 of 2



# Teleflex Medical

## ISO 13485:2016

## EN ISO 13485:2016



Issue 22

Detailed scope

**Design, development, manufacture and distribution of reusable medical and surgical instruments for general and specialty use; sterile and non-sterile disposable surgical, urology, anaesthesia and respiratory medical devices, sterile disposable electrosurgical medical devices. Design of Non-Sterile Nasal and Oral Mucosal Devices.**

**Design and development of sterile single use absorbable and non-absorbable sutures, pledgets and suture guides and manufacturing of non-sterile absorbable and non-absorbable suture material..**

**Manufacturing of sterile single use absorbable and non-absorbable sutures.**

**Distribution of sterile single use absorbable and non-absorbable sutures and non-sterile suture material. Distribution of medical devices for endoscopy; fiber optic illuminators; sterile single use instruments for cardiovascular and general surgical procedures.**

Additional facilities

**375 Forbes Blvd, Mansfield, MA, 02048-1805, United States**



0005

