STORZ

KARL STORZ—ENDOSKOPE

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20 1614 20-1 Kaltlicht-Fontäne Power LED 175 SCB

INSTRUCTION MANUAL

20 1614 20-1 Cold Light Fountain Power LED 175 SCB

MANUAL DE INSTRUCCIONES

20 1614 20-1 Fuente de luz fría Power LED 175 SCB

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the cold light fountain power LED 175. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology may occur for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG

Utilization and reproduction by third parties require express authorization by KARL STORZ SE & Co. KG.

All rights reserved.

RxOnly Federal (USA) law restricts this device to sale by or on the order of a physician.

It is recommended to check the suitability of the product for the intended procedure prior to use.

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2. Controls, displays, connectors, and their uses

- ① Power switch, lights up when in operation
- ② Light outlet point (not exchangeable)
- ③ Button to increase brightness *
- ④ Button to reduce brightness *
- ⑤ Current light intensity in percent
- ⑥ Standby button **
- ⑦ SCB indicator display ***
- ⑧ to ⑬ - see next page (symbols)
- ⑭ Potential equalization connector
- ⑮ Fuse holder
- ⑯ Power cord socket
- ⑰ Two SCB ports (SCB®, see next page, (symbols))

* Keeping the buttons depressed engages continuous adjustment of the light intensity (autorepeat function).

** If you press the Standby button ⑥ once, the brightness is reset to the lowest value. The integrated LED illuminates. When you press this button again, the brightness is reset to the last valid value. The integrated LED extinguishes.

***If there are several power LED 175 light sources available in the same SCB system at the same time (see next page), the system assigns each power LED 175 light source its own number (1-9), provided that doing so is supported by the infrastructure. IMAGE1 S™ can only support one light source. Up to 9 light sources can be supported in the OR1™ system. This number is displayed at position ⑦ and makes it clear to the user on which unit he is making the changes when using the remote control via the SCB.

Example: The user changes the light intensity of LED light source no. 2 from '25 %' to '50 %' in the camera head's 'SCB menu'. This changes the light intensity of the power LED 175 for which a '2' is displayed at position ⑦.

If only one power LED 175 light source is connected to the same SCB system, no number is displayed. Similarly, no number appears if there is no SCB system present.

If the light source is being remote controlled by means of a CCU, an 'A' is shown.

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	Follow instructions for use (see ⑧ on page IV)
	The KARL STORZ-SCB interface (KARL STORZ Communication Bus) permits remote control of equipment functions, as well as remote display of equipment parameters.
	Brightness (see ⑨ on page IV)
	Lamp / Light source (see ⑩ on page IV)
	Standby button
	Light outlet (see ⑪ / ⑫ on page IV)
	Type CF applied part according to IEC 60601-1 (Explanation: Limitation of leakage currents; CF = Cardiac Floating. See ⑬ on page IV).
	Potential equalization connector (⑭ on page IV)
	Alternating current
	Electronic information product pollution control (China RoHS). The symbol is found on the identification plate on the underside of the unit.
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE). The symbol is found on the identification plate on the underside of the unit.
	Manufacturer

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Symbols on label and packaging of the device

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. This can be downloaded from www.karlstorz.com.

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4. General information

4.1 Device description

The power LED 175 cold light fountain is designed to supply light for endoscopic diagnostic and surgical procedures.

The device offers a wide range of functions. These include:

- Efficiency: High light intensity at with low power consumption
- Comfortable working environment: Intelligent cooling concept for minimal noise emission
- Durable: No costs or effort required to change the lamp if used properly
- Integration: The SCB function enables the light source to be controlled remotely and the system status information to be retrieved, while providing simple configuration during initialization (see ⑰ on page IV, or ⑰ on page V).

4.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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5. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

5.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

NOTE: Any damage to the device resulting from incorrect operation is not covered by the manufacturer's warranty.

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WARNING: Read this instruction manual thoroughly before putting the device into operation. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

WARNING: The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep device out of reach of patients.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this (cf. Clause 16 of the 3rd Ed. of IEC 60601-1).

WARNING: Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

WARNING: At the end of its useful operating life, dispose of the unit as electronic scrap.

WARNING: Avoid allowing fluids to penetrate the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

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WARNING: Ensure sufficient air circulation for the device.

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WARNING: DO NOT use the device in the presence of flammable anesthetics.

WARNING: Servicing may only be performed by specialists authorized by KARL STORZ.

WARNING: Test this device prior to each procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

WARNING: The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

WARNING: Do not use the device if the input keys are not working properly.

WARNING: Always supervise the start up sequence of the device in order to detect possible key or display errors.

WARNING: The device must not be used in oxygenated environments.

WARNING: The device may only be operated within specified environmental conditions.

WARNING: Always disconnect the device from the mains before cleaning.

WARNING: Always pull out the power plug before any cleaning/disinfection procedures! Make sure the light source system has cooled down.

WARNING: The device must be set up in such a way that the power cord can be simply pulled out at any time.

WARNING: The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.

WARNING: The device may be subject to sudden failures. Keep a spare device at hand or alternate procedures in mind.

WARNING: Applied parts of devices used in combination with this device must be of type BF or CF.

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WARNING: Patient leakage currents may be additive when energized endoscopes are used with energized endotherapy devices. This is particularly important if a CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

WARNING: The compatibility of this device with any accessories and/or energized endotherapy devices should be checked in the individual instructions for use.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Always use the lowest possible settings for your intervention.

WARNING: The light source is transmitting high energy light leading to high surface temperatures on the light cable and on the tip of the endoscope. Tissue damage is possible.

WARNING: Danger of glare! Never look into the light outlet when the light system is switched on.

WARNING: Insufficient ventilation of the unit can cause an internal build-up of heat which leads to a safety shut-off of the light source.

CAUTION: If a continuous tone is emitted when the unit is switched on, an error has occurred during start-up. Switch the unit off and start it up again. If the error occurs again, there is a general error. Send the unit in for repair.

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WARNING: Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see Clause 16 of the 3rd Ed. of IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

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Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

5.2 Normal use

It is only allowed to use the power LED 175 light source in combination with approved accessories from KARL STORZ or with equipment that is approved as safety-related uncritical.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

WARNING: No modification of this equipment is allowed.

5.2.1 Indication

The device is intended to create light for endoscopic diagnostic observation and surgical procedures.

5.2.2 Contraindication

The medical device is contraindicated for ophthalmologic interventions. Mandatory protection filters for these kinds of interventions are not implemented for the medical device.

For other interventions the responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

5.3 Intended use

Cold light sources are designed to supply light for endoscopic diagnostic and surgical procedures.

5.4 User qualification

Only qualified doctors and medical assistant staffs are allowed to use the power LED 175 light source after training. The training has to be conducted by skilled persons only. The area agency or the responsible KARL STORZ sales representative is available for trainings, for advices about further training alternatives.

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5.5 Patient profile

The use of the power LED 175 is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

5.6 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Have been thoroughly trained at least once in the use of the device
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means

5.7 Intended conditions of use

5.7.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual. It is the user's responsibility to make sure the device is safe and operates properly before using it. During treatment using the power LED 175 cold light fountain, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

WARNING: Any treatment may only be performed if visual observation of the action of the device is ensured.

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5.7.2 Other intended conditions

Frequency of use: one or more times a day
 Length of use: from a few minutes to several hours a day.
 Place of installation: Positioning on a level, vibration-free surface
 Mobility: can be moved if positioned on a cart
 Combination: can be used on the patient at the same time as other equipment required for the operations
 Control: can be controlled via the KARL STORZ-SCB.

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5.8 User position

For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$ at a distance of approx. 30 – 70 cm to the front panel.

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
5.9 Safety precautions at the site of installation

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.
 It is not intended for use in hazardous zones.

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This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants. The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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! WARNING: The device should be placed on a flat surface or a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

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5.10 Safety features

The cold light fountain power LED 175 is equipped with the following safety features:

5.10.1 Self-test

A self-test is performed each time the device is switched on. If any of the tests is negative, a corresponding error message is displayed.

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Functions and features

Following installation and switching on of the unit, the light source emits cold light from its light port. The intensity of the light can be regulated by pressing the + and - buttons ③ and ④. The current intensity is displayed by the green LED. When the unit is switched on, the lowest setting is always active.

- ① **NOTE:** The light intensity can also be controlled by remote control or automatically, provided the technical requirements for this are satisfied (see following paragraph).
- ① **NOTE:** Various light sources (up to 9) can be integrated in one SCB installation. The units can then only be differentiated on the remote station (monitor) by their SCB address. For further information: please see the explanation for ① on page V.
- ① **NOTE:** Not all systems support this (IMAGE1 S™ can only support one light source).

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6. Installation and operating instructions

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6.1 Unpacking the equipment

Carefully unpack the cold light fountain power LED 175 and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

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6.2 Basic equipment


KARL STORZ Set No.: 20161401-1

- 1 power LED 175 cold light fountain 20161420-1
- 1 400 A power cord, length 300 cm
- 1 Instruction manual 96206197 D
- 1 Protective wallet for instruction manual
- 1 SCB connection cable 100 cm 200901 70

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6.3 Installation

NOTE: The cold light fountain power LED 175 including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

! WARNING: The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

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Set the device on a horizontal, flat surface (Fig. 1).

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6.3.1 Connecting the ground line

The device is equipped with a connector ⑭ for attaching a ground line.

The device's ground line should be installed by a qualified electrician (Fig. 2).

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6.3.2 Connecting the power cord

Connect power cord and insert the power plug fully into the power cord receptacle ⑯ (Fig. 3).

WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

WARNING: Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

NOTE: Ensure that access to the power plug is ensured at the site of installation.

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6.3.3 KARL STORZ-SCB

NOTE: The term 'SCB' is explained on page VI. The bottom of page V also provides an example of use of the SCB function (see explanation for ⑦).

Connect the SCB cable 20 0901 70 to an (equivalent) SCB socket on the back of the device (Fig. 4). Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ-SCB control NEO System Instruction Manual).

NOTE: The SCB 'power LED 175' system can be operated with any camera control unit (CCU) which has an SCB (e.g., IMAGE1 S™).

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i NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Therefore, to unplug the SCB cable, pull out the connector – do not pull directly on the cable.

Commissioning the SCB system

The units of the SCB system can be turned on in any order. When the camera control unit (e.g., IMAGE1 S™) is turned on, the following SCB message appears: New Device

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6.3.4 Connecting light cable

Plug the light cable into the socket ② as far as it will go (Fig. 5). Only handle the light cable using the handle piece. Never pull on the cable.

i NOTE: The light outlet ② is provided with an antiglare flap, which does not allow any direct outlet of light. Light is only emitted when the light cable is inserted.

! WARNING: Risk of burning!
The light source transmits high energy light which can lead to high surface temperatures. Tissue damage is possible.

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6.3.5 Connecting the light cable to the endoscope or videoscope

Connect the light cable to the endoscope by a quarter turn of the knurled collar on the screw base (Fig. 6 or 7).

i NOTE: The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.

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General applications

The power LED 175 is suitable for use in a wide range of applications. The correct light intensity depends on the particular application, the instrument being used, and the connected light cable. If the camera control unit used is not SCB-compatible, the intensity should always be kept to a minimum.

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**6.4 Commissioning the
power LED 175 light source**


Press the power switch ① (Fig. 8).

ⓘ NOTE: If no error is noted, the unit's self test is ended with a short tone. The light source is then at the lowest intensity setting (approx. 5 %). However, the intensity display ⑤ shows the last setting used before the unit was switched off. That means: the segment which was the highest illuminated segment prior to the unit being switched remains permanently lit.

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6.4.1 Choosing the last setting

The unit starts in standby mode (at the lowest setting). If the standby button is then pressed, the last valid setting is activated (Fig. 9).

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6.4.2 Choosing the lowest setting

The unit starts in standby mode and is already automatically at the lowest setting. Therefore, no further action is necessary.

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6.4.3 Choosing the individual setting

The unit starts in standby mode and is automatically at the lowest setting. To select the desired setting, press the + button (repeatedly, if necessary).

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6.4.4 Setting the brightness manually


The brightness can be set with the buttons ③ and ④ (Fig. 10).
To increase the brightness press the button ③.
To decrease the brightness press the button ④.
The current brightness setting is shown by the LED display ⑤.

! **NOTE:** If the buttons ③ and ④ are held down, the autorepeat function is activated. The value then continues to change in the selected direction until the minimum or maximum setting is reached.

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Standby button ⑥

Pressing this button once resets the brightness to the lowest value. The integrated LED illuminates.
Pressing this button again resets the brightness to the last valid value. The integrated LED extinguishes.
If the plus button is pressed when the device is in standby mode, the light intensity returns to the lowest value (5 %).

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6.4.5 Decommissioning

For decommissioning, switch off the unit at the power switch.

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**6.5 Safety functions of the
equipment**

- The unit includes a range of safety functions. During power-up, the unit performs a self-test. This enables errors or impermissible operating conditions to be detected, right from the start. In the case of errors, the unit enters safety mode immediately. No light is emitted.
- If impermissible operating conditions first become apparent during operation, this will initially be signaled. If the situation deteriorates, the light will switch off automatically. This gives the user the possibility to react to the development either by resolving the cause of the error or by concluding his treatment in a controlled manner.

The following error displays may be shown:**a) Constant tone after turning on and dark intensity display:**

There are two possible causes of this:

- **A hardware problem was discovered.**
The light source is not working properly. This can also be assumed if the error is signaled constantly even when the unit is cold. The unit must be sent in for repair.
- **The unit is too hot:**
It is irrelevant whether the unit has become too hot during operation or solely as a result of the external circumstances (storage). As long as it is too hot, it will not start. The unit must be switched off and allowed to cool.

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**b) Warning sound (two short tones) every
15 seconds; simultaneously, the fields of the
display which show the current light intensity
flash**

- The unit has become too hot during operation. If the signal sounds, the unit can still be used for a while.

CAUTION: Remove the cause of the overheating (e.g., poor convection, effect of heat or light, etc.) or reduce the light intensity (output). If this is not possible, conclude your treatment.

- If the temperature continues to increase, the cold light turns off when the unit reaches the critical temperature. The optical/acoustic signal remains on. It turns off once the temperature drops back down to an admissible value. However, the cold light remains turned off (restart block).
- If the critical temperature is not exceeded or if the unit temperature drops as a result of the measures taken, the unit continues to operate. As soon as the temperature reaches admissible values, the signal (optical/acoustic) is deleted.

Notes on errors a) and b)

i NOTE: If the overheating fuse has responded and the unit is subsequently turned off and immediately turned back on, error a) is signaled (constant alarm). The unit must firstly be allowed to cool down to permissible temperatures.

- As soon as the unit has cooled down, turn the unit back on. As long as there are no other problems (e.g., blocked fan), the unit starts again normally.

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7. Maintenance

7.1 Fuse replacement

! **WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Switch off the device and remove the power plug from the electrical outlet (Fig. 12).

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
Remove the fuse holder (15) with a slim screwdriver or other suitable tool (Fig. 13).

! **WARNING:** Only use fuses of the correct rating.

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! **WARNING:** The fuse must fulfill the requirements for a component of high integrity characteristics.

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Refit the fuse holder (15) (Fig. 15).
Reconnect the power cord.
Test the device for proper operation.

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7.2 Reprocessing

7.2.1 General warnings for medical devices

WARNING: Risk of infection: Incorrectly prepared medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: Country-specific laws and regulations must be observed.

NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

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CAUTION: Make absolutely sure that no liquid can get into the device housing.

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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
7.3 **Maintenance and safety check**

WARNING: Always unplug the unit and disconnect the endoscope before carrying out any maintenance or cleaning.

7.3.1 **Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

7.3.2 **Safety checks/repeat inspections as per IEC 62353**

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/ repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 6.4).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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7.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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7.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

WARNING: At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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Servicehotline: +49 7461/708 980
E-Mail: technicalsupport@karlstorz.com

7.6 Repair program

For videoscopes and devices, individual repair is necessary. Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

For repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany, you can refer repairs directly to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 NEUHAUSEN

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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7.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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7.8 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty. The same is valid for any cases of improper use (i.e. operation outside of the specified operating data).

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Issue
Complete failure of the device.

Possible cause

1. Power supply failure.
2. Defective power fuse.
3. Power cord connector is not properly connected to device socket.

Recommended remedy

1. Check that there is electricity to the wall outlet.
2. Change fuses as described in the instruction manual. Use only fuses of the correct rating.
3. Push power cord connector firmly into device socket.

Issue
No light is emitted.

Possible cause

1. Electronics are defective.
2. Overheating of the power LED 175 (LED has switched off due to overheating of the unit caused by covered air vent).

Recommended remedy

1. Notify Service.
2. Remove objects from in front of the air vents. Switch off the power LED 175, let it cool down (approx. 10–15 minutes) and switch it back on.
Ensure that it is sufficiently ventilated.

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Issue
No light output, power switch is lit (on).

Possible cause

1. Strong impacts or vibrations during operation.
2. Power supply or LED is defective.

Recommended remedy

1. Turn the unit off and then on again.
2. Have the unit repaired.

Issue
Too little light.

Possible cause

1. Ends of light cable and/or endoscope are soiled.
2. Light cable or endoscope is defective.
3. LED service life has been exceeded.

Recommended remedy

1. Clean ends of light cables, and light entrances and exits on the endoscope.
2. Replace light cable or endoscope.
3. Have the unit repaired.

Issue
No light output, power switch is lit (on). LED-brightness indicator is flashing.

Possible cause

1. Defective LED.
2. Defective fan.
3. Defective temperature sensor.
4. Overheating of the unit.

Recommended remedy

- 1.-3. Have the unit repaired.
4. Turn the unit off and then on again.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

According to IEC 60601-1-2:

Observe the information on electromagnetic compatibility in chapter 10 (p. 33-49).

Directive compliance (for 20161420-1)
According to the Medical Device Directive (MDD):

This medical device belongs to Class I
 This medical device bears the CE mark in accordance with MDD 93/42/EEC.

Conformidad con la norma (para 20161420-1)
Según CEI 60601-1, CEI 60601-1-6, ANSI/AAMI ES60601-1, CEI 60601-2-18, UL 60601-1, CAN/CSA-C22.2 No. 60601-1:

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8.3 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Item	Cat. no.
SCB connecting cable, length 100 cm	20 0901 70
All KARL STORZ light cables are suitable for use	(see catalog; corresponding specialities)

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10. Electromagnetic Compatibility (EMC) Information

The cold light fountain power LED 175 20 1614 20-1 is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Accessories and cables used for EMC Compliance				
Type	Shielded	Length [m]	Ferrite	Used for
PA (Potential equalization)	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

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**Electromagnetic
Compatibility
(EMC) Information**

WARNING: The use of an accessory or cable with the power LED 175 other than those specified in this manual may result in increased emissions or decreased immunity of the power LED 175. When using an accessory or cable with power LED 175 other than those specified in this manual, it becomes the responsibility of the user of the power LED 175 to determine compliance with IEC 60601-1-2 when using this item.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the power LED 175, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests			
Guidance and manufacturer's declaration – electromagnetic immunity			
The power LED 175 is intended for use in the electromagnetic environment specified below. The user of power LED 175 should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the power LED 175 require continued operation during power mains interruptions, it is recommended that the power LED 175 be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	100 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the power LED 175 further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	10 V _{rms} on 150 kHz to 80 MHz 1 kHz 80 % AM modulation 6 V _{rms} in ISM bands	

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Table 3 – Test levels for Radiated and conducted Immunity Tests
Guidance and manufacturer's declaration – electromagnetic immunity

The power LED 175 is intended for use in the electromagnetic environment specified below. The user of power LED 175 should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	10 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the power LED 175, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the power LED 175 is used exceeds the applicable RF compliance level above, the power LED 175 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the power LED 175.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group Guidelines and manufacturer’s declaration – electromagnetic emissions		
The power LED 175 is intended for use in the electromagnetic environment specified below. The customer or user of the power LED 175 should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The power LED 175 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The power LED 175 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5			
Recommended separation distances between portable and mobile RF communications equipment and the power LED 175			
The power LED 175 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the power LED 175 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the power LED 175 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.			
Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzrcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada, Del. Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstućo st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraße Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & CO. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +84 28 3823 8000
Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com