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Instruction for use
A solid-phase enzyme immunoassay kit
for the qualitative detection of IgM antibodies
to *Borelia burgdorferi sensu lato*
in human serum or plasma

Borelia burgdorferi IgM EIA

Catalogue number **REF** **K118G**

For 96 determinations

In vitro diagnostic medical device


XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: info@xema.com.ua
www.xema.com.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de

ASSAY PROCEDURE

During performing several independent series of tests, Positive and Negative Control Serum should be used **each time**.

K118M

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY AND CALCULATION OF RESULTS	8
12. INTERPRETATION OF THE RESULTS	8
13. PERFORMANCE CHARACTERISTICS	9
14. LIMITATIONS	9
15. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the qualitative detection of IgM antibodies
to *Borelia burgdorferi sensu lato*
in human serum or plasma

Borelia burgdorferi IgM EIA

1. INTENDED USE

ELISA reagent kit Borelia burgdorferi IgM EIA is a solid-phase enzyme immunoassay for the qualitative detection of IgM antibodies to *Borelia burgdorferi sensu lato* in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Borrelia burgdorferi sensu lato - is a group of borreliosis or Lyme disease pathogens, a common infection, the main host and vector of which is the ixodid tick. The disease is transmitted only through a tick bite.

In the early stages of borreliosis, fatigue, chills and headaches may be observed, and later more serious symptoms may occur, such as joint pain, meningitis, numbness in the extremities, facial nerve paralysis, memory disorders, and eye and heart damage. After the spirochete penetrates the skin, a creeping erythema occurs, and after several days or weeks, it reaches many organs by haematogenous or lymphatic means. In general, the incubation period is from 3 to 45 days.

Early diagnosis of the disease is based on clinical and epidemiological data. The diagnosis is confirmed by laboratory, usually by serological methods - the detection of specific antibodies to *Borrelia burgdorferi* in the blood.

IgM antibodies appear in the blood first, a few days after infection, but can be detected by laboratory tests in 2-3 weeks. After about 6 weeks, the concentration of antibodies reaches a maximum and then gradually decreases. IgG antibodies begin to be detected 4-6 weeks after infection and the maximum amount of IgG antibodies is synthesised 2-3 months after the onset of early symptoms of the disease. Then their number gradually decreases, but they remain in the body for several years.

3. TEST PRINCIPLE

The detection of IgM antibodies to *Borrelia burgdorferi* is based on the indirect enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized recombinant *Borrelia burgdorferi* antigen. The analysis procedure includes three stages of incubation:

- during the first stage specific to *Borrelia burgdorferi* antibodies from the specimen are bound onto the microwell surface;
- during the second stage horseradish peroxidase-conjugated specific monoclonal anti-IgM antibodies bind to the antigen-antibody complexes, fixed in the formed at the previous stage complexes;
- during the third stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density (OD) in the microwell is directly related to the concentration of the measured IgM antibodies to *Borrelia burgdorferi* in test specimen.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P118MZ	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with recombinant antigen of <i>Borrelia burgdorferi</i> , ready to use
CN118MZ	CONTROL -	Negative Control Serum K-	0.5 mL	1	Solution based on human serum, free of IgM antibodies to <i>Borrelia burgdorferi</i> , with preservative, ready to use (yellow liquid)
CP118MZ	CONTROL +	Positive Control Serum K+	0.2 mL	1	Solution based on human serum, containing of IgM antibodies to <i>Borrelia burgdorferi</i> , with preservative, ready to use (red liquid)
T118MZ	CONJ HRP	Conjugate Solution	12 mL	1	Solution of monoclonal antibodies to IgM conjugated to the horseradish peroxidase, ready to use (red liquid)
SP118MZ	DIL SPE	EIA Buffer	12 mL	1	Buffer solution with detergent and preservative, ready to use (purple liquid)
R055Z	SUBS TMB	Substrate Solution	12 mL	1	Tetramethylbenzidine (TMB) substrate solution, ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	30 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	12 mL	1	5.0% solution of sulphuric acid, ready to use (colourless liquid)

The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs).

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450/620-680 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 1^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The positive and negative control serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, avoid no more than three cycles of thawing-freezing samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The *Borrelia burgdorferi* IgM EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The *Borrelia burgdorferi* IgM EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Positive and Negative Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- the remaining strips should be immediately resealed in the bag along with the silica gel, closed with the zip-lock, and stored at +2...+8°C within 2 months;
- EIA Buffer, Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Positive and Negative Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. The remaining strips should be immediately resealed in the bag along with the silica gel and closed with the zip-lock to prevent moisture from affecting the plate's strips.

9.3. Washing solution preparation

Add the contents of the 30 mL washing solution concentrate vial to 750 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the washing solution concentrate, mL	2.5	5	7.5	10	12.5	15	17.5	20	22.5	25	27.5	30
Volume of water, mL	62.5	125	187.5	250	312.5	375	437.5	500	562.5	625	687.5	750

10. ASSAY PROCEDURE

- 10.1. Put the desired number of strips into the frame based on the number of test samples and 4 wells for Positive and Negative Control Serum (1 well for Positive Control (CP) and 3 wells for Negative Control Serum (CN)).
- 10.2. Dispense **90 µL of EIA Buffer** to all wells.
- 10.3. Dispense **10 µL of Positive and Negative Control Serum as well as 10 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Positive and Negative Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Positive and Negative Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CP	SAMP5	SAMP13	SAMP21								
B	CN	SAMP6	SAMP14	SAMP22								
C	CN	SAMP7	SAMP15	SAMP23								
D	CN	SAMP8	SAMP16									
E	SAMP1	SAMP9	SAMP17									
F	SAMP2	SAMP10	SAMP18									
G	SAMP3	SAMP11	SAMP19									
H	SAMP4	SAMP12	SAMP20									

- 10.4. Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **30 minutes at +37°C**.
- 10.5. At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well **5 times** using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5 µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6. Add **100 µL of Conjugate Solution** to all wells.
- 10.7. Cover strips with a plate sealing tape and incubate for **30 minutes at +37°C**.
- 10.8. At the end of the incubation period, aspirate and wash each well **5 times** as described in 10.5.
- 10.9. Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.10. Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.11. Read the optical density (OD) of the wells at 450 nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the stop solution. Set photometer blank on air.

11. TEST VALIDITY AND CALCULATION OF RESULTS

11.1. The test results are valid only if Positive and Negative Control Serum are within the specified ranges and if all other test parameters are also within the given assay specifications, namely:

- OD of CONTROL- < 0.15;
- OD of CONTROL+ > 1.5;
- $OD(CN) \times 0,5 < OD(CN) < OD(CN) \times 2$.

11.2. Calculate the mean OD value of the Negative Control Serum:

$$\text{meanOD(CN)} = (\text{OD1(CN)} + \text{OD2(CN)} + \text{OD3(CN)})/3$$

11.3. Calculate the Cut Off value by adding to the mean OD value of the Negative Control Serum the coefficient 0.3.

$$\text{Cut off} = \text{meanOD(CN)} + 0.3$$

11.4. Calculate Positivity Index (PI) for each sample by dividing the OD of the sample by Cut off value:

$$\text{PI} = \text{ODsample}/\text{Cut off}$$

12. INTERPRETATION OF THE RESULTS

- If PI value > 1.1 the result is **POSITIVE**,
- If PI value is between 0.9 and 1.1 the result is **EQUIVOCAL**,
- If PI value < 0.9 the result is **NEGATIVE**.

If equivocal results are obtained, it is recommended to conduct a reexamination of the sample in several replicates. If the result is equivocal again, a new sample should be obtained within 5-7 days and retested. If the result remains equivocal, the sample should be considered negative.

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1. Precision of Measurement

Reproducibility (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples, with different levels of IgM antibodies to the *Borrelia burgdorferi sensu lato* antigen, during 1 day in 43 replicates on one series of ELISA kit.

Nº serum	mean OD	mean PI	CV PI, %
1	0.31	1.02	8.29
2	1.08	3.59	6.65

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation (CV) for 2 samples of each serum for 4 days in 8-replicate determinations.

Nº serum	mean OD	mean PI	CV PI, %
1	0.27	0.9	8.1
2	1.1	3.66	6.8

13.1.2. Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL, hemoglobin in a concentration of up to 10 mg/mL and triglycerides in a concentration of up to 10 mg/mL.

13.2. Diagnostic performance characteristics

The clinical sensitivity and specificity of the assay were evaluated using a serum panel with 8 positive and 8 negative clinical serum samples and were 100%. The relative sensitivity and specificity of the assay were investigated in a sample of 96 donor sera characterised for the content of IgM antibodies to *Borrelia burgdorferi* antigen in commercial Kits, and the results were 99.7% and 97.5%, respectively.

14. LIMITATIONS

A positive result is evidence of the presence of IgM antibodies to *Borrelia burgdorferi sensu lato* antigen. The diagnosis cannot be based on the results of an IgM antibody test to *Borrelia burgdorferi* alone and requires confirmation, including an assessment of the patient's clinical presentation and history, the detection of IgG antibodies to *Borrelia burgdorferi* and conducting an immunoblot test.

A negative result indicates the absence of IgM antibodies to *Borrelia burgdorferi sensu lato* or antibody levels below the limit of sensitivity of the kit.

The results of serum tests in patients with immunosuppression and immunological disorders should be interpreted with caution.

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LOT _____

DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua