BIOCHEMISTRY CONTROL URINE I, II

BioSystems

ORINA CONTROL DE BIOQUÍMICA I, II

NIVEL I	NIVEL II
COD 18054 1 x 5 mL	COD 18066 1 x 5 mL

USO PREVISTO

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INTENDED USE

The Biochemistry Control Urine is intended for quality control in the laboratory only and is supplied with suggested intervals of acceptable values. For *in vitro* diagnostic professional use in the clinical laboratory only.

COMPOSITION

Biochemistry Control Urine. For 5 mL. Lyophilized human urine containing component concentrations suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests.

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

PREPARATION AND USE

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 5 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. Analyze the reconstituted material in the same way as the patient samples.

STORAGE AND STABILITY

Store at 2-8°C.

Biochemistry Control Urine is stable until the expiration date given in the label.

The components of the reconstituted material are stable for at least 7 days at 2-8°C and 30 days at -20°C (when frozen once).

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

SPECIFIED VALUES

The assigned concentration values for components and their traceability are shown in the enclosed value sheets. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

The suggested intervals of acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only. Each laboratory should establish its own precision parameters.

NOTICE TO USER

The Summary of Safety and Performance is available from EUDAMED database: <https://ec.europa.eu/tools/eudamed>.

Language versions for IFU and SDS are available in the BioSystems website <https://einfo.bio>.

M18054i-05

First version updated according to IVDR 2017/746 requirements.

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