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N
TearFlo™ Sterile Tear Flow Test Strips (Schirmer Test)
Instructions for use
 Intended use: Intended for use during a tear secretion test.

5. **Schirmer test:** the oldest method of diagnosis of the dry eye syndrome (DES). **2. Product description** The TearFlo™ strip is a rounded, blunt end that should be placed under the lower eyelid without touching the cornea. The printed scale (0-35 mm), graduated in 1mm increments (± 0.03), has the accuracy ± 1 mm. Each strip is individually packaged in a sealed and coated medical package and in film to ensure product sterility. **3. Contraindications:** Do not use in patients with damaged cornea, conjunctivitis, conjunctival inflammation, corneal or conjunctival ulceration, and a case of hypersensitivity to any of the product components. **5. Users:** The product is designed exclusively for professional users. **6. Medical device class** An Im class medical device bearing the CE mark and having an EU declaration of conformity. **7. Use** Do not apply any topical medication or perform any manipulation on the eye prior to performing the examination. The examination is performed by placing the strip under the lower eyelid. The examiner holds the examination hair with the head resting against the forehead. Be careful not to touch the rounded end of the strip when opening the package. With the patient looking up, carefully insert the rounded end under the lower part of the eye lid (about 10mm from the temporal edge or corner of the eye) and then bend the strip down slightly. Wait 5 minutes, remove the strip, and record the length of the moistened area, measuring from the folded end. **Interpretation of results:** ≥ 15 mm - normal result; 10-15 mm - indicates an initial or mild dry eye; < 10 mm - advanced DES. **8. Disposal:** Because of the possibility of contamination, dispose of only with medical waste. The product is not recyclable.

R
TearFlow™ Sterile Tear Flow Test Strips – bandelettes de diagnostic
(Test de Schirmer)
Notice d'utilisation

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Madhu Instruments Pvt. Ltd

F-90/3D, Okhla Industrial Area,
Phase -1, New Delhi – 110 020, India

incl. site: A 260, Okhla Industrial Area,
Phase -1, New Delhi, 110 020, India

Development, Manufacture and Supply of Diagnostic & Surgical devices used in ophthalmology like Ophthalmic Diagnostic Strips, Ophthalmic Devices for Tissue Manipulation, Surgical Blades & Knives, Ophthalmic Cannula, Ophthalmic Microsurgery Instruments, Silicon Device (Silicon Band, Silicon Sponge, & Silicon Tire), Diagnostic & Surgical Lenses, Ophthalmic Devices for Teaching & Training Devices, Ophthalmic Surgical Instruments, Cataract Protection.

The validity of the **qualityaustria** certificate will be maintained by annual surveillance audits and one renewal audit after three years.

This **qualityaustria** certificate confirms the application and further development of an effective

QUALITY MANAGEMENT SYSTEM

complying with the requirements of standard

EN ISO 13485:2016

Medical devices - Quality management systems -
Requirements for regulatory purposes

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Konrad Scheiber
General Manager

Dr. Mag. Anni Koubek
Specialist representative

Quality Austria - Trainings,
Zertifizierungs und
Begutachtungs GmbH is
accredited according to
the Austrian Accreditation
Act by the BMWFVW
(Federal Ministry of
Science, Research and
Economy)

Quality Austria is
accredited as an
organisation for
environmental verification
by the BMLFUW (Federal
Ministry of Agriculture,
Forestry, Environment and
Water Management)

Quality Austria is
authorized by the VDA
(Association of the
Automotive Industry)

For accreditation
registration details please
refer to the applicable
decisions or recognition
documents

Quality Austria is the
Austrian member of IQNet
(International Certification
Network)

Dok. Nr. FO_24_028

31f01c68-7e5b-4e42-
82cd-2bc899b8354

The current validity of the certificate is documented exclusively on the Internet under
<http://www.qualityaustria.com/en/cert> EAC: 19.2

EC CERTIFICATE

Full Quality Assurance System

Certificate No.: 11305-2017-CE-IND-NA-PS Rev 4.0

Initial certification date:
16 October 2017

Valid Until: 27 May 2024

This is to certify that the quality system of:

Madhu Instruments Pvt. Ltd.

A - 260, Okhla Industrial Area, Phase-1, New Delhi – 110 020, INDIA

For design, production and final product inspection/testing of:

STERILE DISPOSABLE MEDICAL DEVICES FOR OPHTHALMOLOGY PROCEDURES

Has been assessed with respect to:

The conformity assessment procedure described in Annex II excluding section 4 of Council Directive 93/42/EEC on Medical Devices, as amended

and found to comply

Further details of the product(s) and conditions for certification are given overleaf.

Place and date:
Høvik, 18 May 2021

For the issuing office:
Notified Body 2460
DNV Product Assurance AS

Alessandra Rinna
Assessor

Notice: The Certificate is subject to terms and conditions as set out in the Certification Agreement. Failure to comply may render this Certificate invalid.

NOTIFIED BODY 2460: DNV Product Assurance AS, Veritasveien 3, 1363 Høvik, Norway, Tel +47 67 57 88 00, www.dnv.com

ICP-4-5-i1-MDD-f2, rev.0

Jurisdiction

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as “Forskrift om Medisinsk Utstyr” by the Norwegian Ministry of Health and Care Services.

Certificate history:

Revision	Description	Issue Date
0.0	Supersedes DNVGL (NB 0434) Certificate no: 8911-2016-CE-IND-NA following transfer to notified body functions to DNV GL Nemko Presafe AS (NB 2460)	16 Oct 2017
1.0	EU Representative Change of Address	18 Jul 2018
2.0	EU Representative Change	10 Feb 2020
3.0	Recertification	13 May 2021
4.0	Editorial changes	18 May 2021

Products covered by this Certificate:

Product Description	Product Name	Class
Ophthalmic Diagnostic Strips	MIPL/A1 Schirmer Tear Test Strips Brand: Tear Touch MIPL/A6 Schirmer Tear Test Strips with Blue Mark Brand: Tear Touch Blue	Im /Is
Ophthalmic Devices for Tissue Manipulation	MIPL/B1 Cellulose Sponge Spears Brand : High Soak MIPL/B2 PVA Sponge Spears Brand : Soft Soak	Is

Ophthalmic Devices for Tissue Manipulation	MIPL/D1 Iris Retractors Brand: Iris Care MIPL/D2 Capsular Tension Ring Size 12-10mm, 13-11mm, 14-12mm MIPL/D3 Capsular Tension Ring with Scleral Fixation Arm Size 12-10mm, 13-11mm, 14-12mm (for Left, Right, Both) MIPL/D4 Capsular Tension Ring Injector MIPL/D8 Capsule Hook Brand : Capsule Care MIPL/D10 Pupil Dilator 3 to 8 loop Brand : Gupta Ring MIPL/D11 Capsule Support Segment Brand : Gupta Segment	Ila
Surgical Blades and Knives	MIPL/C1 Corneal Trephine - 1 to 20mm Brand: Nano Edge MIPL/C2 Suction Trephine - 6mm to 20mm Brand: Nano Edge MIPL/C4 Trephine Punch (vacuum & non-Vacuum) 6mm to 20mm Brand: Nano Edge MIPL/C3 Micro Incision Blade – - Lance Tip (15° to 45°) Brand: Nano Edge - Slit (1.2 to 4.0mm) Brand: Nano Edge - MVR (19G, 20G, 23G & 25G) Brand: Nano Edge - Crescent (1 to 2.5 mm) Brand: Nano Edge - Enlarger (5.0 to 5.5mm) Brand: Nano Edge	Ila

Ophthalmic Cannula	MIPL/E2 Silicon Tip Cannula - Regular/Brush Tip (20G, 23G, 25G,27G) MIPL/E4 Anterior Chamber Maintainer 19G, 20G,23G MIPL /E5 Lacrimal Cannula (DCR) Regular / Olive Tip (4.5, 7.5cm, 11cm, 15cm, 17.5cm) MIPL/E6 Anaesthesia Cannula 19G,21G, 23G, 25G, 27G MIPL/E7 Cystotome Cannula 25G, 26G 27G, 30G MIPL /E8 Air injection & Irrigation Cannula 16G ,19G ,20G,21G,23G ,25G ,26G ,27G, 28G,30G MIPL /E9 Hydrodissection / Hydrodelineation Cannula 23G, 25G, 27G, 30G MIPL /E10 Cortex Aspirating Cannula 21G, 23G, 25G, 26G MIPL /E11 Irrigating & Aspirating Cannula 23/23,22/23, 23/22 MIPL /E12 Capsule Polishers 21G,23G,25G,27G MIPL /E15 Bimanual Cannula 21G,22G,23G MIPL /E16 Lasik Cannula 16G, 23G, 25G, 27G Brand: Maxiflo (All cannulas)	Ila
Cataract Protection	MIPL/O3 Eye Shield Brand : Clear Shield	Is

The complete list of devices is filed with the Notified Body

Sites covered by this certificate

Site Name	Address
Madhu Instruments Pvt. Ltd. (Unit II)	A - 260, Okhla Industrial Area, Phase-1, New

Certificate No.: 11305-2017-CE-IND-NA-PS Rev 4.0
Place and date: Høvik, 18 May 2021

	Delhi – 110 020, INDIA
--	------------------------

EU Representative

Obelis s.a. Bd , General Wahis 53, B-1030 Brussels , Belgium

Terms and conditions

The certificate is subject to the following terms and conditions:

- Any producer (see 2001/95/EC for a precise definition) is liable for damage caused by a defect in his product(s), in accordance with directive 85/374/EEC, as amended, concerning liability of defective products.
- The certificate is only valid for the products and/or manufacturing premises listed above.
- The Manufacturer shall fulfil the obligations arising out of the quality system as approved and uphold it so that it remains adequate and efficient.
- The Manufacturer shall inform the Notified Body of any intended updating of the quality system and the Notified Body will assess the changes and decide if the certificate remains valid.
- Periodical audits will be held, in order to verify that the Manufacturer maintains and applies the quality system. the Notified Body reserves the right, on a spot basis or based on suspicion, to pay unannounced visits.

The following may render this Certificate invalid:

- Changes in the quality system affecting production.
- Periodical audits not held within the allowed time window.

Conformity declaration and marking of product

When meeting with the terms and conditions above, the producer may draw up an EC declaration of conformity and legally affix the CE mark followed by the Notified Body identification number.

End of Certificate

APPENDIX TO EC CERTIFICATE

Appendix to Certificate no.:
11305-2017-CE-IND-NA-PS Rev. 4.0

Valid Until:
27 May 2024

This is an Appendix issued to EC Certificate issued for manufacturer:
Madhu Instruments Pvt. Ltd.

originally issued in compliance with:
the Council Directive 93/42/EEC on Medical Devices, as amended

Based on assessment audit performed, the following changes to the certification has been accepted as compliance with Council Directive 93/42/EEC on Medical Devices has been established.

A device already covered by the certification is placed on the market under an additional name or Identifier.

Products covered by this Certificate (replaces information on certificate):		
Product Description	Product Name	Class
Ophthalmic Diagnostic Strips	MIPL/A1 Schirmer Tear Test Strips Brand: Tear Touch, Tear Flo MIPL/A6 Schirmer Tear Test Strips with Blue Mark Brand: Tear Touch Blue	Im / Is
Ophthalmic Devices for Tissue Manipulation	MIPL/B1 Cellulose Sponge Spears Brand: High Soak MIPL/B2 PVA Sponge Spears Brand : Soft Soak	Is
Ophthalmic Devices for Tissue Manipulation	MIPL/D1 Iris Retractors Brand: Iris Care MIPL/D2 Capsular Tension Ring Size 12-10mm, 13-11mm, 14-12mm Brand: CapSafe	Ila

Place and date:
Høvik, 02 February 2022

For the issuing office:
DNV Product Assurance AS - Notified Body
2460
Veritasveien 3, 1363 Høvik, Norway

Hazem Tinawi
Technical Reviewer

	MIPL/D3 Capsular Tension Ring with Scleral Fixation Arm Size 12-10mm, 13-11mm, 14-12mm (for Left, Right, Both) Brand: CapSafe MIPL/D4 Capsular Tension Ring Injector Brand: CapSafe MIPL/D8 Capsule Hook Brand : Capsule Care MIPL/D10 Pupil Dilator 3 to 8 loop Brand: Gupta Ring MIPL/D11 Capsule Support Segment Brand: Gupta Segment, Canabrava Gupta Segment	
Surgical Blades and Knives	MIPL/C1 Corneal Trephine - 1 to 20mm Brand: Nano Edge MIPL/C2 Suction Trephine - 6mm to 20mm Brand: Nano Edge MIPL/C4 Trephine Punch (Vacuum & non Vacuum) 6mm to 20mm Brand: Nano Edge MIPL/C3 Micro Incision Blade – a. Lance Tip (15° to 45°) Brand: Nano Edge b. Slit (1.2 to 4.0mm) Brand: Nano Edge c. MVR (19G, 20G, 23G & 25G) Brand: Nano Edge d. Crescent (1 to 2.5 mm) Brand: Nano Edge e. Enlarger (5.0 to 5.5mm) Brand: Nano Edge	Ila
Ophthalmic Cannula	MIPL/E2 Silicon Tip Cannula - Regular/Brush Tip (20G, 23G, 25G, 27G) MIPL/E4 Anterior Chamber Maintainer 19G, 20G, 23G MIPL/E6 Anaesthesia Cannula 19G, 21G,	Ila

	23G, 25G, 27G MIPL/E7 Cystotome Cannula 25G, 26G, 27G, 30G MIPL/E8 Air injection & Irrigation Cannula 16G, 19G, 20G, 21G, 23G, 25G, 26G, 27G, 28G, 30G MIPL/E9 Hydrodissection / Hydrodelineation Cannula 23G, 25G, 27G, 30G MIPL/E10 Cortex Aspirating Cannula 21G, 23G, 25G, 26G MIPL/E11 Irrigating & Aspirating Cannula 23/23, 22/23, 23/22 MIPL/E12 Capsule Polishers 21G, 23G, 25G, 27G MIPL/E15 Bimanual Cannula 21G, 22G, 23G MIPL/E16 Lasik Cannula 16G, 23G, 25G, 27G Brand: Maxiflo (All above cannulas) MIPL/E5 Lacrimal Cannula (DCR) Regular / Olive Tip (4.5, 7.5cm, 11cm, 15cm, 17.5cm) Brand: Lacri Care	
Cataract Protection	MIPL/O3 Eye Shield Brand: Clear Shield	Is

Appendix History -		
Revision	Description	Issued Date
0	Brand Addition in Bold & Editorial changes	02 February 2022

TECHNICAL CONSTRUCTION FILE
Ref. No.: MIPL/TCF/STTS

Issue No.: 03
Issue Date: 10.06.2020
Revision No.: 03
Revision Date: 02.02.2022

DECLARATION OF CONFORMITY

Section No.- 19

DECLARATION OF CONFORMITY

According to annex –II excluding section IV of the MDD 93/42/EEC & as amended directive 2007/47/EC concerning medical devices we:

Madhu Instruments Pvt. Ltd.

Registered. Address :- F-90/3D, Okhla Industrial Area Phase-I, New Delhi-110020 India.

Manufacturing Facility: - A-260, Okhla Industrial Area Phase-I, New Delhi-110020 India.

Declare under our sole responsibility that the product:

S. No.	Product Group	Commercial Name	Model / Variants	Brand	Rule & Class	Sterile/ Non sterile	Lot No	Mfg. date	Exp. Date
1.	Ophthalmic Diagnostic Strips	Schirmer Tear Test Strips	N/A	Tear Flo	Rule 5, Class Im /Is	Sterile	N/A	N/A	N/A

Meets the provisions of the MDD 93/42/EEC & by amended directive 2007/47/EC concerning medical devices which apply to them:

The MD is in Class Im/ Is, according to annex-IX of the Medical Device Directive and certified as per Rule 5. The conformity procedure referred to Annex II-excluding section 4 has been applied in order to affix the CE marking on the device.


Following standards were used to prove the products conformity with the essential requirements of the Directive:

[EN ISO 14971:2012], [EN ISO 10993-1 : 2009/AC :2010], [EN ISO 10993-5:2009], [EN ISO 10993-7: 2008/AC 2009], [ISO 10993-10:2010], [EN ISO 10993-12 :2012], [EN ISO 10993- 18 :2009], [EN ISO 11607-1:2009], [EN ISO 11607-2:2006], [EN ISO 13485:2016], [ISO 14644-1:2015], [ISO 14644-2 : 2015], [ISO 14644-3: 2005], [ISO 14644-4 : 2001], [ISO 14644-5 : 2004], [ISO 14644-7 : 2004], [ISO 14644-8:2013], [ISO 14644-9:2012], ISO 14644-10:2013], [EN ISO 11737-1 :2006/AC :2009], [EN ISO 11737-2 :2009], [EN ISO 11140-1 :2009], [ISO 11138-1 : 2017] , [EN ISO 11138 (Part 2&3) :2009], [EN 1041:2008], [EN ISO 15223-1:2021], [ISO 15223-2:2010], [EN ISO 11135-1:2007 & ISO 11135:2014], [EN 62366: 2008], [EN ISO 15004-1 :2009], [I.P. 2018], [USP 42]

The certification company Quality Austria has issued the following certificate: EN ISO 13485:2016, No. 00154/0 issued on the 25th March 2014 valid upto 31 March 2024.

Signatory established within the EU who has been empowered to enter into commitments on our behalf:

Specimen Copy

Prepared By  Executive QA	Checked & Issued By  Head QA	Approved By  Director
--	---	--

TECHNICAL CONSTRUCTION FILE
Ref. No.: MIPL/TCF/STTS

Issue No.: 03
Issue Date: 10.06.2020
Revision No.: 03
Revision Date: 02.02.2022

DECLARATION OF CONFORMITY

Section No.- 19

NOTIFIED BODY:

DNV Product Assurance AS,
Veritasveien 3,
1363 Høvik, Norway,
Tel +47 67 57 88 00,
www.dnv.com

NOTIFIED BODY No.: 2460

Obelis s.a.

Bd. Général Wahis 53
B-1030 Brussels, Belgium
Phone: 32.2.732.59.54
Fax: 32.2.732.60.03
E-mail: mail@obelis.net

CE Certificate No: 11305-2017-CE-IND-NA-PS Rev 4.0

Date of Issue: 18 May 2021

Valid till: 27 May 2024

For Madhu Instruments Pvt. Ltd.

Name: **Dinesh Goel** (Employee ID code-099)

Position: Manager, Quality Assurance

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Dla Madhu Instruments Pvt. Ltd.