







Metrohm Autolab B.V.  
Kanaalweg 29/G  
3526 KM, Utrecht  
The Netherlands  
+31302893154  
autolab@metrohm.com  
<https://www.metrohm.com/en/products/electrochemistry/>

# **NOVA**

2.1.4

## **User Manual**

Metrohm Autolab Teachware  
Metrohm Autolab B.V.  
3526 KM, Utrecht

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above or at [autolab@metrohm.com](mailto:autolab@metrohm.com).

This documentation is protected by copyright. All rights reserved.

# Table of contents

<b>1</b>	<b>NOVA installation</b>	<b>1</b>
1.1	Software compatibility .....	1
1.2	Hardware compatibility .....	2
1.3	Software installation .....	2
1.4	External devices .....	5
1.4.1	Metrohm Devices support .....	6
1.4.2	Metrohm Devices installation .....	7
1.4.3	Spectrophotometer support .....	7
1.4.4	Spectrophotometer installation .....	8
1.4.5	Autolab RHD Microcell HC support .....	8
1.4.6	Autolab RHD Microcell HC installation .....	9
1.5	Powering the instrument .....	9
1.6	Autolab hardware installation .....	10
1.7	Software license .....	12
1.8	Intended use .....	13
1.9	Options .....	13
<b>2</b>	<b>Conventions</b>	<b>16</b>
2.1	Scientific conventions .....	16
2.2	Software conventions .....	16
2.3	Numbering conventions .....	17
2.4	Warning label conventions .....	18
2.5	NOVA information, warnings and errors .....	18
2.6	NOVA menus and controls .....	20
<b>3</b>	<b>Release notes</b>	<b>23</b>
3.1	<b>Version 2.1.4 release</b> .....	<b>23</b>
3.1.1	Chrono charge discharge galvanostatic .....	24
3.1.2	Editable overlay plot properties .....	25
3.1.3	Disable event logging .....	26
3.1.4	Import table in the repeat loop during a running procedure ...	27
3.2	<b>Version 2.1.3 release</b> .....	<b>28</b>
3.2.1	Procedure information .....	28
3.2.2	Export to RelaxIS .....	29
3.2.3	OCP command .....	30
3.2.4	Import table in Repeat command .....	31
3.2.5	Light source shutter control .....	31
3.2.6	Demo database file .....	32



3.2.7	Colored data file command tiles .....	33
<b>3.3</b>	<b>Version 2.1.2 release .....</b>	<b>34</b>
<b>3.4</b>	<b>Version 2.1.1 release .....</b>	<b>34</b>
3.4.1	Signal names, identity and locations .....	34
3.4.2	Current range logging .....	37
3.4.3	Event logging .....	38
3.4.4	Export options for Spectrophotometer control panel .....	40
3.4.5	Export options for Spectrophotometer control panel .....	41
3.4.6	Spectroelectrochemistry procedure .....	42
<b>3.5</b>	<b>Version 2.1 release .....</b>	<b>44</b>
3.5.1	Search function .....	44
3.5.2	Check cell .....	46
3.5.3	Current interrupt .....	46
3.5.4	Spectrophotometer manual control .....	47
3.5.5	Spectroelectrochemical measurements .....	48
3.5.6	Repeat number in Repeat command .....	50
3.5.7	Custom command name .....	51
3.5.8	Zoom function .....	52
3.5.9	Electrochemical Frequency Modulation .....	53
3.5.10	Corrosion rate analysis .....	54
3.5.11	New Plots frame controls .....	55
3.5.12	Device drivers installation .....	57
<b>3.6</b>	<b>Version 2.0.2 release .....</b>	<b>58</b>
3.6.1	Managed schedules .....	58
3.6.2	New color picker .....	60
3.6.3	Data handling command shortcut button .....	61
3.6.4	Library filters .....	62
3.6.5	Extended Sampler information .....	63
3.6.6	Zoom function for data analysis commands .....	64
3.6.7	Custom name for Build signal command .....	65
<b>3.7</b>	<b>Version 2.0.1 release .....</b>	<b>66</b>
3.7.1	Procedure and data tags .....	67
3.7.2	New plot options .....	68
3.7.3	Import data .....	70
3.7.4	Number of recent items .....	70
3.7.5	Plot preview .....	71
3.7.6	Print plot .....	72
3.7.7	Region insensitivity .....	74
3.7.8	Electrochemical interface toggle .....	74
3.7.9	ECl10M measurements .....	75
3.7.10	Library column display .....	77
3.7.11	Data grid column display .....	77
3.7.12	Estimated duration .....	78
3.7.13	Interpolate command .....	79
3.7.14	Hydrodynamic analysis .....	79
<b>3.8</b>	<b>Version 2.0 release .....</b>	<b>80</b>
3.8.1	Dynamic data buffers .....	80

3.8.2	Value of Alpha .....	81
3.8.3	Autolab RHD Microcell HC support .....	81
3.8.4	PGSTAT204 and M204 combination with Booster10A support .....	82
3.8.5	ECI10M module support .....	82
3.8.6	AC voltammetry .....	83
<b>4</b>	<b>Dashboard</b> .....	<b>84</b>
4.1	<b>Actions</b> .....	<b>85</b>
4.2	<b>Recent items</b> .....	<b>86</b>
4.3	<b>What's going on</b> .....	<b>88</b>
4.4	<b>Instruments panel</b> .....	<b>90</b>
<b>5</b>	<b>Instruments panel</b> .....	<b>92</b>
5.1	<b>Change the default instrument</b> .....	<b>95</b>
5.2	<b>Autolab control panel</b> .....	<b>96</b>
5.2.1	Instrument information panel .....	97
5.2.2	Tools panel .....	98
5.2.3	Autolab display panel .....	127
5.3	<b>Autolab RHD Microcell HC control panel</b> .....	<b>134</b>
5.3.1	Autolab RHD Microcell HC hardware setup .....	134
5.3.2	Autolab RHD Microcell HC manual control panel .....	135
5.4	<b>Autolab Spectrophotometer control panel</b> .....	<b>136</b>
5.4.1	Autolab Spectrophotometer hardware setup .....	137
5.4.2	Autolab Spectrophotometer manual control panel .....	140
5.5	<b>Metrohm devices control panel</b> .....	<b>149</b>
5.5.1	Metrohm Dosino control panel .....	149
5.5.2	Metrohm Sample Processor control panel .....	155
5.5.3	Metrohm Stirrer control panel .....	164
5.5.4	Metrohm Remote box control panel .....	169
<b>6</b>	<b>Library</b> .....	<b>173</b>
6.1	<b>Default procedures</b> .....	<b>175</b>
6.2	<b>Add location</b> .....	<b>177</b>
6.3	<b>Default save Location</b> .....	<b>179</b>
6.4	<b>Moving files to a new location</b> .....	<b>180</b>
6.5	<b>Remove location</b> .....	<b>180</b>
6.6	<b>Load from Library</b> .....	<b>182</b>
6.7	<b>Edit name and remarks</b> .....	<b>183</b>
6.8	<b>Rating and tagging</b> .....	<b>184</b>
6.9	<b>Preview plot</b> .....	<b>186</b>
6.10	<b>Column visibility</b> .....	<b>187</b>



6.11	<b>Filtering the Library</b>	188
6.12	<b>Sorting the Library</b>	191
6.13	<b>Rearranging Library columns order</b>	192
6.14	<b>Locating files</b>	193
6.15	<b>Delete files from Library</b>	194
6.16	<b>The data repository</b>	195
6.17	<b>Merge data</b>	197
6.18	<b>Search function</b>	200
<b>7</b>	<b>NOVA commands</b>	<b>203</b>
7.1	<b>Control commands</b>	<b>204</b>
7.1.1	Message	204
7.1.2	Send email	205
7.1.3	Repeat	206
7.1.4	Increment	224
7.1.5	Play sound	227
7.1.6	Build text	228
7.1.7	.NET	229
7.1.8	Procedure Information	233
7.2	<b>Measurement - general</b>	<b>235</b>
7.2.1	Autolab control	236
7.2.2	Apply	239
7.2.3	Cell	240
7.2.4	Wait	240
7.2.5	OCP	246
7.2.6	Set pH measurement temperature	247
7.2.7	Reset EQCM delta frequency	248
7.2.8	Autolab R(R)DE control	250
7.2.9	MDE control	251
7.2.10	Synchronization	254
7.3	<b>Measurement - cyclic and linear sweep voltammetry commands</b>	<b>257</b>
7.3.1	CV staircase	257
7.3.2	CV linear scan	260
7.3.3	LSV staircase	262
7.4	<b>Measurement - voltammetric analysis commands</b>	<b>264</b>
7.4.1	Sampled DC voltammetry	265
7.4.2	Normal pulse voltammetry	267
7.4.3	Differential pulse voltammetry	270
7.4.4	Differential normal pulse voltammetry	273
7.4.5	Square wave voltammetry	276
7.4.6	PSA (Potentiometric stripping analysis)	279
7.4.7	AC voltammetry	283



<b>7.5</b>	<b>Measurement - chrono methods commands</b>	<b>286</b>
7.5.1	Record signals	286
7.5.2	Chrono methods	290
<b>7.6</b>	<b>Measurement - impedance commands</b>	<b>302</b>
7.6.1	FRA measurement	303
7.6.2	FRA single frequency	305
7.6.3	Additional properties	307
7.6.4	Electrochemical Frequency Modulation	327
<b>7.7</b>	<b>Data handling commands</b>	<b>332</b>
7.7.1	Windower	332
7.7.2	Build signal	337
7.7.3	Calculate signal	345
7.7.4	Get item	358
7.7.5	Import data	358
7.7.6	Export data	362
7.7.7	Generate index	364
7.7.8	Shrink data	365
<b>7.8</b>	<b>Analysis - general commands</b>	<b>367</b>
7.8.1	Smooth	368
7.8.2	Peak search	372
7.8.3	Regression	373
7.8.4	Derivative	376
7.8.5	Integrate	378
7.8.6	Interpolate	379
7.8.7	FFT analysis	379
7.8.8	Convolution	381
7.8.9	Calculate charge	385
7.8.10	Hydrodynamic analysis	386
7.8.11	ECN spectral noise analysis	387
7.8.12	iR drop correction	391
7.8.13	Baseline correction	392
7.8.14	Corrosion rate analysis	398
<b>7.9</b>	<b>Analysis - impedance</b>	<b>403</b>
7.9.1	Electrochemical circle fit	403
7.9.2	Fit and simulation	405
7.9.3	Kronig-Kramers test	449
7.9.4	Include all FRA data	452
7.9.5	Potential scan FRA data	453
<b>7.10</b>	<b>Metrohm devices commands</b>	<b>454</b>
7.10.1	Dosino	455
7.10.2	Sample Processor	460
7.10.3	Stirrer	468
7.10.4	Remote	469
<b>7.11</b>	<b>External devices commands</b>	<b>472</b>
7.11.1	Spectroscopy	472
7.11.2	External device control	479
7.11.3	RHD control	483



<b>8</b>	<b>Default procedures</b>	<b>486</b>
8.1	<b>Cyclic voltammetry</b>	<b>487</b>
8.1.1	Cyclic voltammetry potentiostatic	487
8.1.2	Cyclic voltammetry galvanostatic	490
8.1.3	Cyclic voltammetry potentiostatic current integration	493
8.1.4	Cyclic voltammetry potentiostatic linear scan	496
8.1.5	Cyclic voltammetry potentiostatic linear scan high speed	499
8.2	<b>Linear sweep voltammetry</b>	<b>503</b>
8.2.1	Linear sweep voltammetry potentiostatic	503
8.2.2	Linear sweep voltammetry galvanostatic	506
8.2.3	Linear polarization	509
8.2.4	Hydrodynamic linear sweep	512
8.2.5	Hydrodynamic linear sweep with RRDE	518
8.2.6	Spectroelectrochemical linear sweep	525
8.3	<b>Voltammetric analysis</b>	<b>529</b>
8.3.1	Sampled DC polarography	530
8.3.2	Normal pulse voltammetry	535
8.3.3	Differential pulse voltammetry	539
8.3.4	Differential normal pulse voltammetry	543
8.3.5	Square wave voltammetry	547
8.3.6	AC voltammetry	551
8.4	<b>Chrono methods</b>	<b>555</b>
8.4.1	Chrono amperometry ( $\Delta t > 1$ ms)	556
8.4.2	Chrono coulometry ( $\Delta t > 1$ ms)	559
8.4.3	Chrono potentiometry ( $\Delta t > 1$ ms)	561
8.4.4	Chrono amperometry fast	564
8.4.5	Chrono coulometry fast	566
8.4.6	Chrono potentiometry fast	570
8.4.7	Chrono amperometry high speed	573
8.4.8	Chrono potentiometry high speed	576
8.4.9	Chrono charge discharge potentiostatic	580
8.4.10	Chrono charge discharge galvanostatic	583
8.5	<b>Potentiometric stripping analysis</b>	<b>587</b>
8.5.1	Potentiometric stripping analysis	587
8.5.2	Potentiometric stripping analysis constant current	588
8.6	<b>Impedance spectroscopy</b>	<b>590</b>
8.6.1	FRA impedance potentiostatic	591
8.6.2	FRA impedance galvanostatic	594
8.6.3	FRA potential scan	597
8.6.4	FRA current scan	601
8.6.5	FRA time scan potentiostatic	605
8.6.6	FRA time scan galvanostatic	608
8.6.7	Electrochemical Frequency Modulation	612
<b>9</b>	<b>Additional measurement command properties</b>	<b>615</b>
9.1	<b>Sampler</b>	<b>616</b>

<b>9.2</b>	<b>Automatic current ranging</b>	<b>618</b>
<b>9.3</b>	<b>Cutoffs</b>	<b>620</b>
9.3.1	Cutoff configuration	621
9.3.2	Combining cutoffs	623
<b>9.4</b>	<b>Counters</b>	<b>624</b>
9.4.1	Counter configuration	625
9.4.2	Counter action - Pulse	627
9.4.3	Counter action - Autolab control	629
9.4.4	Counter action - Shutter control	630
9.4.5	Counter action - Get spectrum	631
9.4.6	Combining counters	632
<b>9.5</b>	<b>Plots</b>	<b>633</b>
9.5.1	Default plots	634
9.5.2	Custom plots	634
9.5.3	Plot options	636
<b>9.6</b>	<b>Automatic integration time</b>	<b>641</b>
<b>9.7</b>	<b>Value of Alpha</b>	<b>642</b>
<b>10</b>	<b>Procedure editor</b>	<b>644</b>
<b>10.1</b>	<b>Creating a new procedure</b>	<b>645</b>
<b>10.2</b>	<b>Global options and global sampler</b>	<b>647</b>
<b>10.3</b>	<b>End status Autolab</b>	<b>651</b>
<b>10.4</b>	<b>Procedure tracks</b>	<b>652</b>
<b>10.5</b>	<b>Procedure wrapping</b>	<b>653</b>
<b>10.6</b>	<b>Procedure zooming</b>	<b>654</b>
<b>10.7</b>	<b>Command groups</b>	<b>655</b>
10.7.1	Grouping commands	655
10.7.2	Ungrouping commands	656
10.7.3	Renaming groups	657
<b>10.8</b>	<b>Enabling and disabling commands</b>	<b>658</b>
10.8.1	Disabling commands	658
10.8.2	Enabling commands	659
<b>10.9</b>	<b>Adding and removing commands</b>	<b>660</b>
10.9.1	Adding commands	660
10.9.2	Removing commands	667
<b>10.10</b>	<b>Moving commands</b>	<b>668</b>
10.10.1	Moving commands using the drag and drop method	669
10.10.2	Using the drag and drop method to move commands to a command group or a sub-track	670
<b>10.11</b>	<b>Moving multiple commands</b>	<b>672</b>
<b>10.12</b>	<b>Stacking commands</b>	<b>674</b>
10.12.1	Creating command stacks	675



10.12.2	Remove commands from stacks .....	677
<b>10.13</b>	<b>Links .....</b>	<b>678</b>
10.13.1	Viewing links .....	679
10.13.2	Creating links .....	681
10.13.3	Editing links .....	689
<b>10.14</b>	<b>My commands .....</b>	<b>692</b>
10.14.1	Saving a My command .....	693
10.14.2	Editing My commands .....	695
<b>11</b>	<b>Running measurements .....</b>	<b>698</b>
<b>11.1</b>	<b>Starting procedure .....</b>	<b>698</b>
<b>11.2</b>	<b>Procedure validation .....</b>	<b>701</b>
<b>11.3</b>	<b>Procedure cloning .....</b>	<b>702</b>
<b>11.4</b>	<b>Plots frame .....</b>	<b>704</b>
11.4.1	Displaying multiple plots .....	706
<b>11.5</b>	<b>Real time modifications .....</b>	<b>708</b>
11.5.1	Real-time properties modification .....	708
11.5.2	Procedure control .....	711
11.5.3	Reverse scan direction .....	712
11.5.4	Display the Manual control panel .....	713
11.5.5	Enable and disable plots .....	715
11.5.6	Q+ and Q- determination .....	718
<b>11.6</b>	<b>End of measurement .....</b>	<b>719</b>
11.6.1	Procedure time stamp .....	719
11.6.2	Post validation .....	720
<b>11.7</b>	<b>Specify plot preview .....</b>	<b>721</b>
<b>11.8</b>	<b>Detailed plot view .....</b>	<b>722</b>
11.8.1	Plot properties .....	723
11.8.2	Toggle the 3D view .....	725
11.8.3	Toggle the step through data mode .....	725
11.8.4	Add an analysis command .....	727
11.8.5	Zooming options .....	728
11.8.6	Print plot .....	729
11.8.7	Export plot to image file .....	731
11.8.8	Relocate plots .....	733
<b>11.9</b>	<b>Viewing the data grid .....</b>	<b>737</b>
11.9.1	Current range logged in the data grid .....	739
11.9.2	Events logged in the data grid .....	739
11.9.3	Formatting the data grid .....	740
11.9.4	Sorting the data grid .....	741
11.9.5	Changing the order of the columns in the data grid .....	742
11.9.6	Exporting the data from the data grid .....	743
<b>11.10</b>	<b>Convert data to procedure .....</b>	<b>745</b>

<b>12 Data analysis</b>	<b>747</b>
<b>12.1 Smooth analysis</b>	<b>748</b>
12.1.1 SG mode	750
12.1.2 FFT mode	752
<b>12.2 Peak search</b>	<b>754</b>
12.2.1 Automatic search mode	757
12.2.2 Manual peak search	757
12.2.3 Manual adjustments	768
12.2.4 Results	770
<b>12.3 Regression analysis</b>	<b>771</b>
<b>12.4 Integrate</b>	<b>775</b>
<b>12.5 Interpolate</b>	<b>779</b>
<b>12.6 Hydrodynamic analysis</b>	<b>784</b>
<b>12.7 Baseline correction</b>	<b>788</b>
12.7.1 Zooming in/out	793
12.7.2 Fine tuning the baseline correction	794
<b>12.8 Corrosion rate analysis</b>	<b>796</b>
12.8.1 Tafel Analysis	799
12.8.2 Polarization Resistance	803
<b>12.9 Electrochemical circle fit</b>	<b>805</b>
12.9.1 Zooming in/out	809
12.9.2 Fine tuning the baseline correction	810
12.9.3 Copy as equivalent circuit	812
<b>12.10 Fit and simulation</b>	<b>813</b>
12.10.1 Direct fitting or simulation	814
12.10.2 Fitting or simulation using the dedicated editor	816
12.10.3 Viewing the result	823
<b>13 Data handling</b>	<b>827</b>
<b>13.1 Get item</b>	<b>827</b>
<b>13.2 Shrink data</b>	<b>830</b>
<b>14 Data overlays</b>	<b>834</b>
<b>14.1 Create an overlay</b>	<b>834</b>
<b>14.2 Adding data to an overlay</b>	<b>836</b>
<b>14.3 Changing overlay plot properties</b>	<b>838</b>
<b>14.4 Changing overlay plot settings</b>	<b>840</b>
<b>14.5 Hiding and showing plots</b>	<b>842</b>
<b>14.6 Remove data from overlay</b>	<b>845</b>
<b>14.7 Additional Overlay controls</b>	<b>847</b>



<b>15 Procedure scheduler</b>	<b>850</b>
<b>15.1 Remove instrument from schedule</b>	<b>852</b>
<b>15.2 Creating a procedure schedule</b>	<b>853</b>
15.2.1 Open procedures	853
15.2.2 Recent procedures	855
15.2.3 Search Library	856
15.2.4 Remove procedure	857
<b>15.3 Using synchronization points</b>	<b>858</b>
<b>15.4 Naming and saving the schedule</b>	<b>861</b>
<b>15.5 Running the schedule</b>	<b>864</b>
15.5.1 Starting the complete procedure schedule	865
15.5.2 Starting the schedule sequentially	866
15.5.3 Procedure schedule control	867
<b>15.6 Inspecting procedures or data</b>	<b>869</b>
<b>15.7 Schedule zooming</b>	<b>871</b>
<b>16 Hardware description</b>	<b>873</b>
<b>16.1 General considerations on the use of the Autolab potentiostat/galvanostat systems</b>	<b>874</b>
16.1.1 Electrode connections	874
16.1.2 Operating principles of the Autolab PGSTAT	877
16.1.3 Environmental conditions	897
16.1.4 Noise considerations	897
16.1.5 Cleaning and inspection	899
<b>16.2 Instrument description</b>	<b>900</b>
16.2.1 Autolab N Series (AUT8) instruments	900
16.2.2 Autolab F Series (AUT8) instrument	912
16.2.3 Autolab MBA N Series (AUT8) instruments	925
16.2.4 Autolab Compact Series (AUT4/AUT5) instruments	927
16.2.5 Multi Autolab Series (MAC8/MAC9) instruments	941
16.2.6 Autolab 7 Series (AUT7) instruments	951
16.2.7 $\mu$ Autolab Series instruments	963
<b>16.3 Module description</b>	<b>972</b>
16.3.1 Common modules	973
16.3.2 Optional modules	998
<b>17 Diagnostics</b>	<b>1181</b>
<b>17.1 Connecting the instrument</b>	<b>1181</b>
<b>17.2 Running the Diagnostics</b>	<b>1183</b>
<b>17.3 Integrator calibration</b>	<b>1186</b>
<b>17.4 Diagnostics options</b>	<b>1186</b>
<b>17.5 Firmware update</b>	<b>1187</b>



<b>18 Warranty and conformity</b>	<b>1189</b>
<b>18.1 Warranty</b> .....	<b>1189</b>
<b>18.2 Spare part availability</b> .....	<b>1190</b>
<b>18.3 Safety information</b> .....	<b>1190</b>
<b>18.4 EU Declaration of conformity</b> .....	<b>1192</b>
18.4.1 EU Declaration of Conformity .....	1192
18.4.2 EU Declaration of Conformity .....	1195
18.4.3 Declaration of Conformity .....	1197
18.4.4 Declaration of Conformity .....	1199
<b>18.5 Environmental protection</b> .....	<b>1201</b>
<b>Index</b>	<b>1202</b>





# 1 NOVA installation

This chapter describes how to install the NOVA software on the host computer and to how connect Autolab and external devices to the host computer.

The NOVA installation package is supplied on CD-ROM or USB support provided with the Autolab instrument. It can also be downloaded from the Metrohm Autolab webpage.



## NOTICE

Leave the Autolab disconnected from the computer when installing NOVA for the first time.

## 1.1 Software compatibility

NOVA requires Windows 7 or later as operating systems in order to run properly. NOVA can be installed on 32 bit and 64 bit versions of Windows.



## NOTICE

Previous versions of Windows are not supported.

The minimum and recommended specifications are reported in *Table 1* and *Table 2*, respectively.

*Table 1 Overview of the minimum specifications for NOVA*

CPU	1 GHz or faster 32-bit (x86) or 64-bit (x64) processor
RAM	2 GB RAM
HD	20 GB available hard disc space
GPU	DirectX 9.0c compliant display adapter with 64 MB RAM

*Table 2 Overview of the recommended specifications for NOVA*

CPU	Intel Core i5 or equivalent AMD processor
RAM	8 GB RAM






Figure 1 The installation wizard

Click the **INSTALL** button to start the installation. The files will be copied on the computer. If needed, the installation folder can be changed using the installation wizard.

When prompted to do so, please click the **Install** button provided in the Metrohm Autolab Driver installation window (see figure 2, page 4).




Figure 2 Install the Metrohm Autolab device drivers

**NOTICE**

Make sure that the Always trust software from Metrohm Autolab B.V. check box is ticked.

When prompted to do so, please click the **Install** button provided in the Avantes Driver installation window (see figure 3, page 4).




Figure 3 Install the Avantes device drivers

**NOTICE**

Make sure that the Always trust software from Avantes BV check box is ticked.

The installer will indicate when the installation is completed, as shown in Figure 4.




Figure 4 The installation is complete

## 1.4 External devices

The following additional external devices can be connected to the host computer:

- **Metrohm liquid handling devices:** these devices can be used to handling liquid samples and to automate the handling thereof.
- **Spectrophotometers:** the Autolab or the supported Avantes spectrophotometers can be used to perform spectroelectrochemical measurements in combination with the Autolab potentiostat/galvanostat.
- **Autolab RHD Microcell HC:** this device can be used to perform temperature-controlled measurements.



## 1.4.2 Metrohm Devices installation

Connecting a **USB** controlled Metrohm device to the host computer will trigger the installation of the instrument (see figure 5, page 7).




Figure 5 The Metrohm Device Driver installer

The installation will complete automatically.



### NOTICE

In order to control the supported Metrohm Sample Processors, an additional Windows component must be present on the computer. The controls for the Metrohm Sample processors use the **Microsoft msxml6.0.dll** library for the configuration files (XML file format). This component may not be preinstalled on every Microsoft operating system. Please ensure the availability of this dll on the operating. If this package is missing, please download the installation package from the Microsoft website.

## 1.4.3 Spectrophotometer support

NOVA provides support for Autolab spectrophotometers.

The following Autolab spectrophotometer are supported through a **USB** connection to the host computer:

- Autolab Spectrophotometer UA
- Autolab Spectrophotometer UB

Additionally, compatible Avantes spectrophotometers are also supported when connected to the host computer through a **USB** connection. The following devices are supported:

- AvaSpec ULS2048-USB2





### 1.4.6 Autolab RHD Microcell HC installation

No driver is required to control the Autolab RHD Microcell HC. When this type of device is connected to the host computer, it is immediately recognized by NOVA and listed in the **Instruments** panel.



#### NOTICE

The Autolab RHD Microcell HC is only detected by NOVA if a stage is connected to the controller with a cell mounted on the stage.



#### NOTICE

The Autolab RHD Microcell HC is connected to the host computer through a serial port. If no serial port is present on the computer, a USB to Serial port adapter can be installed. The drivers required for this adapter are not included in the installation package of NOVA and need to be installed separately.

## 1.5 Powering the instrument

In order to use the instrument, it must be connected to the mains using the mains connection socket, located on the back plane of the instrument. Before connecting the instrument to the mains make sure that the mains output voltage matches the value indicated on the main voltage indicator, located above the connector (*see figure 7, page 9*).




Figure 7 The required mains voltage is indicated above the connector






Figure 9 The Driver Manager application

The Driver Manager can be used at any time to select the driver to use to control the Autolab. Two drivers are available:

- **Nova only (recommended setup):** this is the latest driver for the Autolab, allowing up to 127 instruments to be connected to the host computer. This driver is compatible with 32 bit and 64 bit versions of Windows.
- **Legacy driver:** this is an older driver version which can be used in combination with the GPES or FRA software. No further developments are planned for this driver. The maximum number of devices connected to the host computer is 8. Data transfer may be slower than with the NOVA only driver. This driver is **only** compatible with 32 bit versions of Windows.



## NOTICE

The **Nova only** driver will not work with previous versions of NOVA (version 1.10 and older). In order to use previous versions of NOVA, it is necessary to start the Driver Manager application provided with the previous version and select one of the available drivers provided with this previous version (please refer to the **Getting Started** manual of the previous version of NOVA for more information).

Click the installation button for the required driver to change the device driver and follow the instructions on screen. The selected driver will be installed for all connected Autolab instruments. New instruments connected to the host computer will be configured using the selected driver.



## NOTICE

The **Driver Manager** application can be used to change the device driver at any time.

## 1.7 Software license

The Autolab NOVA software, and all its components, provided in conjunction with the Metrohm Autolab potentiostat/galvanostat instruments is copyrighted and owned by Metrohm Autolab.

The software is provided as a **Free Licensed Closed-Source** product with limited warranty. The software can be installed on any computer without specific authorization from Metrohm Autolab.

Metrohm Autolab retains the copyright to the software. You may neither modify nor remove references to confidentiality, proprietary notices or copyright notices. Modifications of the software in part or as a whole is not permitted.

Metrohm Autolab warrants that the software, when operated properly, is suitable for the specified use with the electrochemical instrumentation from Metrohm Autolab or compatible external instrumentation.

Metrohm Autolab is exempt from further warranty or liability. Metrohm Autolab is neither liable for third-party damages or consequential damage not for loss of data, loss of profits or operating interruptions, etc.

## 1.8 Intended use

All Metrohm Autolab products are designed for electrochemical research and development within the normal environment of a laboratory. The instrumentation shall therefore only be used for this purpose and within the specified environmental conditions. All other uses fall out of the scope of the instrumentation and may lead to voiding of any warranty.

## 1.9 Options

The application options can be defined by selecting the *Options* from the **Edit** menu. A window will be displayed, showing two different sections (see figure 10, page 13).




Figure 10 The application Options window

The following properties are available in the **General** section (see figure 10, page 13):

- **Auto save measured data:** specifies if measured data should be saved automatically at the end of each measurement, using the provided  toggle. This option is on by default.
- **Embedded executable IF030:** specifies the path to the embedded application for instrument fitted with the **IF030** controller, using the provided  button. This is a system property, do **not** change this unless instructed by Metrohm Autolab.
- **Embedded executable IF040:** specifies the path to the embedded application for instrument fitted with the **IF040** controller, using the provided  button. This is a system property, do **not** change this unless instructed by Metrohm Autolab.



- **Logging:** specifies if error logging should be used and which level of logging should be used, if applicable, using the provided drop-down list. This is a system property, do **not** change this unless instructed by Metrohm Autolab.
- **Check for updates at:** specifies the URL for version checks of NOVA. This is a system property, do **not** change this unless instructed by Metrohm Autolab.
- **Number of recent items:** defines the number of recent items shown in the **Recent items** panel of the dashboard. The default value is 5. Please refer to *Chapter 4.2*, for more information on the recent items.
- **Defined tags:** provides a list of tags used in NOVA. This list is empty by default and is automatically populated by user-defined tags through the tagging feature of NOVA. If needed, tags can be removed or added to this list directly. More information on the use of tags can be found in *Chapter 6.8*.
- **Event logging:** specifies if event logging should occur during the measurement, using the provided  toggle. This option is on by default. For more information about **Event logging** please refer to *Chapter 11.9.2*



## CAUTION

Modifying the system properties shown in the **General** section can interfere with the operation of the instrument. Do **not** change these properties unless instructed by Metrohm Autolab.

The **Plots** section displays the default plot options used in NOVA for all plots (see *figure 11, page 15*).




Figure 11 The plot options

In this section all the default options for plots can be specified. Clicking the **Reset** button will reset all the options to the factory default values.



## NOTICE

Please refer to *Chapter 9.5.3* for more information on the Plot options.



## 2 Conventions

Throughout NOVA and all Metrohm Autolab products, the conventions detailed in this chapter are used.

### 2.1 Scientific conventions

The following scientific conventions are used:

- All units are specified in the International System of Units (SI), unless otherwise specified. See *B. N. Taylor, A. Thomson, The International System of Units (SI), NIST Special Publication 330, 2008 Edition* for more information.
- Electrochemical values like potential and current are indicated according to the International Union of Pure and Applied Chemistry (IUPAC) convention. Positive currents and (over)potentials are associated with oxidation processes. Negative currents and (over)potentials are associated with reduction processes. See *A. D. McNaught, A. Wilkinson, IUPAC, Compendium of Chemical Terminology: IUPAC Recommendations, Blackwell Science: Oxford, England; Malden, MA, USA, 1997* for more information.

### 2.2 Software conventions

The following standard interaction conventions are used in NOVA:

- A right-handed mouse where the left button is used for selecting items and the right button may open context-related menus is assumed.
- Quickly pressing and releasing the mouse button is called 'Clicking'. A click of the left mouse button on a menu option, a button, an input item on the screen, will result in an action.
- Quickly pressing and releasing the right mouse button is called 'Right-clicking'. A click of the right mouse button on a suitable location on the screen opens a context-sensitive menu, if applicable.
- By clicking and holding down the left mouse button you can 'Drag' items from one window and 'Drop' it in another by releasing the button. This action will be called 'Drag and Drop' and it is the key mechanism for creating a procedure.
- Quickly pressing and releasing the mouse button twice is called 'Double-clicking'. A double-click of the left mouse button is used to perform particular actions, and mainly is applied through standard usage in window actions.

The following selection methods are used in NOVA:



- To select any item on screen, click the item.
- To select consecutive items on screen, click the first item, press and hold down **[SHIFT]**, and then click the last item.
- To select nonconsecutive items on screen, press and hold down **[CTRL]**, and then click each item.

## 2.3 Numbering conventions

All **numeric** values are defined in NOVA according to the local culture defined for the Windows operating system. Depending on these settings, the decimal separator symbol can either be **.** or **,**.

Improper use of the local culture settings defined in Windows may lead to wrong values. For example, typing 0,3 in NOVA on a computer which uses the **.** as decimal separator will be validated as 3.



### NOTICE

It is recommended to consult the local culture settings defined in Windows before using the NOVA software.

Scientific (exponential) numbering is done using the **e** or **E** symbol. A value of **1e2** or **1E2** is converted to **100**.

The following prefixes are using in NOVA for engineering notation:

- **T**, for *Tera* (1000000000000).
- **G**, for *Giga* (1000000000).
- **M**, for *Mega* (1000000).
- **k**, for *Kilo* (1000).
- **m**, for *Milli* (0.001).
- **μ**, for *Micro* (0.000001).
- **n**, for *Nano* (0.000000001).
- **p**, for *Pico* (0.000000000001).






Figure 12 Information is highlighted in blue

- **Warning:** any item highlighted in yellow indicates that an issue has been identified and that user intervention is recommended in order to resolve the issue (see figure 13, page 19). Whenever possible, the cause and a possible solution will be offered. It is possible to ignore the warning and continue working with the software however this may lead to invalid data.




Figure 13 Warnings are highlighted in yellow

- **Error:** any item highlighted in red indicates that a problem has been identified and that user intervention is required in order to resolve the error (see figure 14, page 20). Whenever possible, the cause and a possible solution will be offered. It is **not** possible to ignore the error. No measurements are possible until the error is resolved.



<b>Redo 'Action name'</b>	Redoes the specified action. Keyboard shortcut: <b>[CTRL] + [Y]</b>
<b>Cut</b>	Cuts the selected item(s) to the clipboard. Keyboard shortcut: <b>[CTRL] + [X]</b>
<b>Copy</b>	Copies the selected item(s) to the clipboard. Keyboard shortcut: <b>[CTRL] + [C]</b>
<b>Paste</b>	Pastes the items in the clipboard at the specified location. Keyboard shortcut: <b>[CTRL] + [V]</b>
<b>Select All</b>	Selects all visible items. Keyboard shortcut: <b>[CTRL] + [A]</b>
<b>Options</b>	Specifies the default options used in the application.
<hr/>	
<b>View</b>	
<b>Zoom in</b>	Zooms in on a plot. Keyboard shortcut: <b>[CTRL] + [=]</b>
<b>Zoom out</b>	Zooms out on a plot. Keyboard shortcut: <b>[CTRL] + [-]</b>
<b>Fit all</b>	Adjusts the plot area to the best possible scale. Keyboard shortcut: <b>[F4]</b>
<b>Manual control</b>	Displays the Manual control panel for the default instrument. Keyboard shortcut: <b>[F10]</b>
<hr/>	
<b>Measurement</b>	
<b>Run</b>	Starts the procedure defined in the selected tab on the <b>default</b> instrument. Keyboard shortcut: <b>[F5]</b>
<b>Run on ►</b>	Starts the procedure defined in the selected tab on the <b>specified</b> instrument.
<b>Instrument #1</b>	
<b>Instrument #2</b>	
<b>...</b>	



<b>Pause</b>	Pauses the running command in the selected tab.
<b>Skip</b>	Skips the running command in the selected tab.
<b>Stop</b>	Stops the measurement running in the selected tab.
<hr/>	
<b>Help</b>	
<b>User manual</b>	Displays the NOVA User Manual. Shortcut key: <b>[F1]</b>
<b>About</b>	Displays the About dialog.

## 3 Release notes

This chapter describes the release notes of the current and previous versions of NOVA. The release notes are provided in reverse chronology. The following version have been released:

- **Version 2.1.4:** minor update of NOVA 2.1.
- **Version 2.1.3:** minor update of NOVA 2.1. *Version 2.1.3 release (see chapter 3.2, page 28)*. This version was released on March 9<sup>th</sup>, 2018.
- **Version 2.1.2:** minor update of NOVA 2.1. *Version 2.1.2 release (see chapter 3.3, page 34)*. This version was released on May 5<sup>th</sup>, 2017.
- **Version 2.1.1:** minor update of NOVA 2.1 *Version 2.1.1 release (see chapter 3.4, page 34)*. This version was released on March 24<sup>th</sup>, 2017.
- **Version 2.1:** the current major release of NOVA *Version 2.1 release (see chapter 3.5, page 44)*. This version was released on November 15<sup>th</sup>, 2016.
- **Version 2.0.2:** minor update of NOVA 2.0 *Version 2.0.2 release (see chapter 3.6, page 58)*. This version was released on July, 6<sup>th</sup>, 2016.
- **Version 2.0.1:** minor update of NOVA 2.0 *Version 2.0.1 release (see chapter 3.7, page 66)*. This version was released on April, 1<sup>st</sup>, 2016.
- **Version 2.0:** the original major release of NOVA 2 *Version 2.0 release (see chapter 3.8, page 80)*. This version was release on October, 7<sup>th</sup>, 2015.

### 3.1 Version 2.1.4 release

Version 2.1.4 adds the following new functionality:

1. A new **Default Procedure** "Chrono charge discharge galvanostatic" was added, along with a new **Demo Database** file that demonstrates this procedure. *Chrono charge discharge galvanostatic (see chapter 3.1.1, page 24)*.
2. Possibility to edit plot properties within the **Overlay** tab. *Editable overlay plot properties (see chapter 3.1.2, page 25)*.
3. Option to disable event logging. *Disable event logging (see chapter 3.1.3, page 26)*
4. Possibility to import a table in the Repeat for multiple values command while the procedure is running. *Import table in the repeat loop during a running procedure (see chapter 3.1.4, page 27)*








Figure 16 Demo data file 22 shows galvanostatic charge discharge of a cylinder cell battery

### 3.1.2 Editable overlay plot properties

The possibility to edit plot properties has been added to the **Overlay** tab. This new functionality allows the plotted data series to be configured, including the point or line styles and colors. Each data series in the overlay plot can now be given a title that is propagated automatically to the legend. The overlay plot can be given a title as well. For more detailed information, please proceed to *Changing overlay plot properties* (see chapter 14.3, page 838).




Figure 17 Overlay tab with a new plot properties panel

### 3.1.3 Disable event logging

The possibility to disable event logging has been added to the **Application options**. When event logging is disabled, events such as cutoffs, user intervention, and the current range in which data points were collected will not be logged. Disabling event logging reduces slightly the size of the measurement data .nox file and reduces the RAM used when creating the .nox data file. Event logging is enabled by default. It is not recommended to disable event logging unless the measurement is expected to create a very large data file and RAM availability during data file creation is expected to be a problem. For more detailed information, please proceed to *Options* (see chapter 1.9, page 13).




Figure 18 The application options window



## NOTICE

Modifying the options shown in the **General** section of the application options can interfere with the operation of the instrument. Do **not** change these properties unless instructed by Metrohm Autolab.

### 3.1.4 Import table in the repeat loop during a running procedure

The possibility to import a table while the procedure is running has been added to the **Repeat** command. This new functionality allows the repeat values to be imported at the moment the Repeat command is encountered in the procedure sequence. Therefore, changes can be made to the values in the table during the measurement sequence, until the moment the command is executed. Use of this feature is recommended for advanced NOVA users only because validation cannot be performed on the values within the table when it is imported during run time. For more detailed information, please proceed to: *Import a table using the Import table option (see chapter 7.1.3.2.4, page 214)*.

Figure 19 The repeat for multiple values command properties with a toggle to repeat a table at runtime





## NOTICE

Use of the **Procedure information** command tile is optional.



## NOTICE

More information on the **Procedure information** command is available in *Chapter 7.1.8*.

### 3.2.2 Export to RelaxIS

NOVA 2.1.3 offers the possibility to export electrochemical impedance spectroscopy data in the RelaxIS file format. This functionality has been added to the existing **Export data** command *Chapter 7.7.6*.



## NOTICE

RelaxIS is a comprehensive software solution for the analysis of impedance spectra. This software is a product of **rhd instruments**. For more information, please visit <https://www.rhd-instruments.de/index.php/en/relaxis-en.html>.

The RelaxIS file format option is located in the File format drop down menu of the **Export data** command, as shown in (*see figure 21, page 29*).





Figure 21 The export data command with RelaxIS file format

When the RelaxIS file format is selected, the list of link-able properties is pre-configured to those supported by the RelaxIS file type, as shown in (*see figure 22, page 30*). Data from a **FRA measurement** can be linked



### 3.2.4 Import table in Repeat command

The **Repeat** for multiple values command now provides the possibility to import a table. The table can be imported from a CSV file, as shown in (see figure 23, page 31).



The screenshot shows a control panel titled 'Values' with a table and two configuration sections. The table has three columns: 'Scan Rate', 'Concentration', and 'Temperature'. Below the table are two sections: 'Add range' and 'Import Table'.

	Scan Rate	Concentration	Temperature
1	0,05	0,1	25
2	0,01	0,1	25
3	0,005	0,1	25
4	0,001	0,1	25

**Add range**

Begin value:

End value:

Number of values:

Distribution:

**Import Table**

Column delimiter:

Decimal separator:

Rows to skip:

Figure 23 Import a table in the Repeat for multiple values command



## NOTICE

More information on how to import a table in the Repeat command can be found in .

### 3.2.5 Light source shutter control

NOVA 2.1.3 adds the possibility to open and close the shutter of an Auto-lab light source from the Spectrophotometer control panel. A toggle is provided that allows the light source shutter to be opened and closed using a TTL pulse from the default connected Autolab instrument.




Figure 24 The light source shutter open toggle is located in the Spectrophotometer control panel



## NOTICE

More information on the Spectrophotometer manual control panel of NOVA 2.1.3 can be found in *Chapter 5.4.2*.

### 3.2.6 Demo database file

A new file has been added to the **Demo database**. Demo data file 21 demonstrates linear sweep voltammetry with in-situ UV/Vis spectroscopy, as shown in *Demo database file* (see chapter 3.2.6, page 32).






Figure 25 Demo data file 21 shows LSV of ferrocyanide oxidation with in-situ UV/Vis spectroscopy

### 3.2.7 Colored data file command tiles

The appearance of the command tiles within a data file has been altered to provide a visual contrast to the command tiles within a procedure file. This style change aims to provide a visual cue that a data file is open in the **Procedure editor**. The tiles in the data file are shaded green. In addition, the word **DATA** now appears at the top left corner of the **Procedure editor**, as shown in (see figure 26, page 33).





Figure 26 Command tiles within a data file are shaded green



For example, applying a **Smooth** command on  $i$  vs  $E$  data (WE(1).Current vs Potential applied), as shown in *Figure 27*, now produces two new signals called *Smoothed WE(1).Current* and *Potential applied*, as shown in *Figure 28*.



*Figure 27* Adding a Smooth command to the  $i$  vs  $E$  plot

The calculated signals are automatically plotted (see *figure 28*, page 36).




Figure 28 Smoothed WE(1).Current vs Potential applied plot is created. Additionally, the calculated data is copied to the parent measurement command (CV staircase), as shown in Figure 29.





Figure 29 The calculated signals are also available in the parent command

The calculated *Smoothed WE(1).Current* signal is known by the NOVA procedure as a valid current signal, obtained by applying a **Smooth** command on the WE(1).Current signal of the **CV staircase** command.

### 3.4.2 Current range logging

NOVA 2.1.1 now logs the active current range for each data point recorded in all measurement commands. This information is stored in the data file and is reported in the data grid (see figure 30, page 37).



Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Index	Current range
0,288086	10,4533	2,86011E-7	0,287781	118	100 nA
0,290527	10,4777	2,88239E-7	0,290253	119	100 nA
0,292969	10,5021	2,97638E-7	0,292664	120	100 nA
0,29541	10,5265	2,98279E-7	0,295074	121	100 nA
0,297852	10,5509	2,96173E-7	0,297607	122	100 nA
0,300293	10,5753	2,99164E-7	0,299957	123	100 nA
0,302734	10,5997	3,02643E-7	0,302338	124	100 nA
0,305176	10,6241	3,08868E-7	0,305298	125	1 µA
0,307617	10,6485	3,0954E-7	0,30777	126	1 µA
0,310059	10,6729	3,07281E-7	0,310394	127	1 µA
0,3125	10,6973	3,1076E-7	0,312897	128	1 µA
0,314941	10,7217	3,21289E-7	0,315277	129	1 µA
0,317383	10,7461	3,21442E-7	0,317657	130	1 µA
0,319824	10,7705	3,18298E-7	0,32016	131	1 µA
0,322266	10,7949	3,19824E-7	0,322601	132	1 µA

Figure 30 The active current range is now reported in the data grid



## NOTICE

This new feature only applies to measurements carried out with NOVA 2.1.1 or later. Measurements carried out with earlier versions of NOVA may not display the active current range properly.



## NOTICE

More information on current range logging is available in *Chapter 11.9.1*.



Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range	Overload	Cutoffs
0.476516	12.0811	4.61322E-7	0.479523	1	196	1.51173E-6	0	100 nA	Current overload	
0.460957	12.1055	4.62459E-7	0.461964	1	197	1.51173E-6	0	100 nA	Current overload	
0.463398	12.1299	4.65382E-7	0.464436	1	198	1.51173E-6	0	100 nA	Current overload	
0.46584	12.1543	4.69227E-7	0.466816	1	199	1.51173E-6	0	100 nA	Current overload	
0.468261	12.1787	4.91069E-7	0.469258	1	200	1.51173E-6	0	100 nA	Current overload	
0.490723	12.2031	4.92767E-7	0.49176	1	201	1.51173E-6	0	100 nA	Current overload	
0.493164	12.2275	4.93378E-7	0.49408	1	202	1.51173E-6	0	100 nA	Current overload	
0.495605	12.2519	4.9762E-7	0.496552	1	203	1.51173E-6	0	100 nA	Current overload	
0.498047	12.2763	4.99329E-7	0.499115	1	204	1.51173E-6	0	100 nA	Current overload	
0.500488	12.3007	5.03754E-7	0.501495	1	205	1.51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.50293	12.3251	5.05707E-7	0.503937	1	206	1.51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.505371	12.3495	5.05615E-7	0.506256	1	207	1.51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.507813	12.3739	5.10681E-7	0.508759	1	208	1.51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.510254	12.3983	5.12939E-7	0.511047	1	209	1.51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A

Figure 32 The details of the cutoff condition are reported in the data grid

The same applies to user intervention, where the action performed by the user is reported in the grid (see figure 33, page 39).

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range	User events
0.710449	16.7506	7.02209E-7	0.71106	1	529	1.72456E-5	-1.44562E-5	1 µA	
0.708008	16.763	7.06177E-7	0.708679	1	530	1.72456E-5	-1.44562E-5	1 µA	
0.705566	16.7752	6.9458E-7	0.706207	1	531	1.72456E-5	-1.44562E-5	1 µA	
0.703125	16.7874	6.99463E-7	0.703766	1	532	1.72456E-5	-1.44562E-5	1 µA	
0.700684	16.7997	7.01599E-7	0.701324	1	533	1.72456E-5	-1.44562E-5	1 µA	
0.698242	16.8119	6.95801E-7	0.698833	1	534	1.72456E-5	-1.44562E-5	1 µA	
0.695801	16.8241	6.9519E-7	0.696411	1	535	1.72456E-5	-1.44562E-5	1 µA	
0.693359	16.8363	6.93359E-7	0.694061	1	536	1.72456E-5	-1.44562E-5	1 µA	
0.690918	16.8485	6.88035E-7	0.691528	1	537	1.72456E-5	-1.44562E-5	1 µA	Reverse scan direction toggled
0.688477	16.8607	6.91833E-7	0.689056	1	538	1.72456E-5	-1.44562E-5	1 µA	
0.686035	16.8729	6.84204E-7	0.686707	1	539	1.72456E-5	-1.44562E-5	1 µA	
0.683594	16.8851	6.88782E-7	0.684204	1	540	1.72456E-5	-1.44562E-5	1 µA	
0.681152	16.8973	6.81763E-7	0.681793	1	541	1.72456E-5	-1.44562E-5	1 µA	
0.678711	16.9095	6.72302E-7	0.679382	1	542	1.72456E-5	-1.44562E-5	1 µA	

Figure 33 User events are reported in the data grid



## NOTICE

This new feature only applies to measurements carried out with NOVA 2.1.1 or later. Measurements carried out with earlier versions of NOVA may not display the recorded events properly.



## NOTICE

More information on event logging is available in *Chapter 11.9.2*.

Furthermore, NOVA now provides indications whenever the measurement conditions can be improved, by highlighting the affected command in blue in the procedure editor (see figure 34, page 40).





Figure 34 Commands are highlighted in blue when the measurement conditions can be improved

### 3.4.4 Export options for Spectrophotometer control panel

NOVA 2.1.1 adds the possibility to export data measured in the **Spectrophotometer** control panel. The data can be exported to ASCII or Excel using the  button located in the top right corner of the control panel (see figure 35, page 40).





Figure 35 The measured data can be exported to ASCII or Excel

It is also possible to export the chart displayed in the **Spectrophotometer** control panel using the  button, as shown in Figure 36.







Figure 36 The chart can also be exported



## NOTICE

More information on the manual control of the Autolab and Avantes spectrophotometers can be found in *Chapter 5.4*.

### 3.4.5 Export options for Spectrophotometer control panel

NOVA 2.1.1 adds the possibility to toggle the *Step through data* option on or off in the Spectrophotometer control panel using the  button in the top right corner of the control panel (see *figure 37, page 42*).




Figure 37 The Step through data option can be used in the Spectrophotometer control panel

When the *Step through data* mode is on, an additional indicator is added to the plot, showing the X and Y coordinates of the point indicated by the arrow. The indicator can be relocated anywhere in the plot area.



## NOTICE

More information on the manual control of the Autolab and Avantes spectrophotometers can be found in *Chapter 5.4*.

### 3.4.6 Spectroelectrochemistry procedure

A procedure for spectroelectrochemical measurements in combination with support spectrophotometers has been added to the **Default procedures**. This procedure is based on a synchronized measurement, using a **LSV staircase** command. The procedure is shown in *Figure 38*.




Figure 38 The default procedure for spectroelectrochemical measurements

This procedure includes a counter that is used to trigger a spectroscopy measurement every points. The **Spectroscopy** command stacked on the **LSV staircase** command uses the data from the two preceding **Spectroscopy** commands to calculate the absorbance and transmittance automatically.



## CAUTION

This procedure requires an Autolab spectrophotometer or a supported Avantes spectrophotometer.



## NOTICE

More information on this procedure can be found in *Chapter 8*.



## 3.5 Version 2.1 release

Version 2.1 adds the following functionality:

1. Application-wide search function for Procedures, Data and Schedules *Search function (see chapter 3.5.1, page 44)*.
2. Check cell tool *Check cell (see chapter 3.5.2, page 46)*.
3. Cell off after current interrupt *Current interrupt (see chapter 3.5.3, page 46)*.
4. Manual control for Autolab and Avantes spectrophotometers *Spectrophotometer manual control (see chapter 3.5.4, page 47)*.
5. New command and command options for spectroelectrochemical applications *Spectroelectrochemical measurements (see chapter 3.5.5, page 48)*.
6. Repeat number added to **Repeat** command *Repeat number in Repeat command (see chapter 3.5.6, page 50)*.
7. Custom names for commands *Custom command name (see chapter 3.5.7, page 51)*.
8. Zoom function for procedure editor and schedule editor *Zoom function (see chapter 3.5.8, page 52)*.
9. New Electrochemical Frequency Modulation (EFM) measurement command available *Electrochemical Frequency Modulation (see chapter 3.5.9, page 53)*.
10. Corrosion rate analysis command expanded with Linear polarization analysis mode *Corrosion rate analysis (see chapter 3.5.10, page 54)*.
11. Improved **Plot** frame controls *New Plots frame controls (see chapter 3.5.11, page 55)*.
12. All device drivers are now included in the NOVA installer *Device drivers installation (see chapter 3.5.12, page 57)*

### 3.5.1 Search function

NOVA now provides the possibility to search for Procedures, Data or Schedules. A dedicated input field is located in the top right corner of the application (see figure 39, page 45).




Figure 39 A search box is provided in the top right corner

The search function can be used to specify a string, with or without wild-cards. When triggered, the search function will look for all Procedures, Data and Schedule items in all the *Locations* specified in the **Library**, except the Default procedures.

The results of the search will be reported in a dedicated tab, grouped by item type. The table controls used to display the results are the same as those used by the **Library**. The results can therefore be sorted or filtered as required.



## NOTICE

The search function will look for all items that match the specified search string in the Name or Remarks.



## NOTICE

More information on the search function can be found in *Chapter 6.18*.



### 3.5.2 Check cell

The Check cell tool is now available from the instrument control panel. This tool can be used to check the electrode connections and the noise level by performing five consecutive current or potential measurements and determining the average value and standard deviation of each measurement (see figure 40, page 46).



Figure 40 The Check cell tool can now be used to check the noise level

The tool can therefore be used to assess the instrument noise pickup and optimize the measurement conditions.



## NOTICE

More information on the **Check cell** tool can be found in *Chapter 5.2.2.4*.

### 3.5.3 Current interrupt

The current interrupt tool has been modified and now allows the possibility to switch the cell off at the end of the measurement (see figure 41, page 47).




Figure 41 The current interrupt tool now provides the possibility to set the cell end state

The *Cell state after measurement* toggle, located in the **Settings** panel, can be used to specify the state of the cell at the end of the measurement. This toggle is off by default.



## NOTICE

More information on the current interrupt tool can be found in *Chapter 5.2.2.2*.

### 3.5.4 Spectrophotometer manual control

NOVA now provides a complete manual control interface for Autolab and Avantes spectrometers. This interface can be used to setup the hardware configuration of the connected spectrophotometer and manually control the spectrophotometer (see figure 42, page 48).




Figure 42 Autolab and Avantes spectrophotometers can be manually controlled

Using this interface it is possible to acquire spectra using the specified properties. It is also possible to save measured spectra as dark and reference (blank) spectra and convert the measured data to absorbance, transmittance or reflectance.



## NOTICE

More information on the manual control of the Autolab and Avantes spectrophotometers can be found in *Chapter 5.4*.


### 3.5.5 Spectroelectrochemical measurements

New measurement command and command options have been added to NOVA in order to facilitate spectroelectrochemical measurements. The **Avantes** command is now replaced with the **Spectroscopy** command, which can be used to control Autolab (and Avantes) spectrophotometers.

It is no longer necessary to initialize and close this type of device in a procedure and the new **Spectroscopy** command now supports a software acquisition mode, which can be used at any time without triggers (*see figure 43, page 49*).







**Properties** ➔

**Software trigger**

Command name

Action

Device name

Start wavelength  nm

Stop wavelength  nm

Integration time  ms

Number of averages

Enable light source shutter control

DIO connector

Shutter open

More

*Figure 43 Software and hardware control is now possible to Autolab and Avantes spectrophotometers*

New measurement options are available for all measurement commands that support them. These options can be used to control the light source shutter position or to acquire a spectrum on the spectrophotometer connected to the DIO port (see figure 44, page 50).






Figure 45 The Repetition number is now available in the Repeat command



## NOTICE

More information on the Repeat command can be found in *Chapter 7.1.3*.

### 3.5.7 Custom command name

For improved readability in the procedure editor, it is now possible to specify a name for all commands in a procedure. Providing a custom name will overrule the default name of the command (*see figure 46, page 52*).




Figure 46 Custom names can now be given to all commands



## NOTICE

When a custom name is provided for a command, the content of this command is no longer updated during a measurement, if applicable.

### 3.5.8 Zoom function

The procedure editor and the schedule editor now offer the possibility to zoom in or out at any time to increase or decrease the size of the items shown on screen. The controls for this new zoom function are located in the top right corner of the editor frame (see figure 47, page 52).




Figure 47 Zoom controls are now available

Using this function will either scale the size of the items and the text up or down (between 200 % and 50 % of the original size), as shown in Figure 48.







Figure 48 Zooming in on the procedure editor

The following zooming controls are available:

- **Zoom out:** decreases the scaling of the items and text shown on screen. The  button or **[CTRL] + [-]** keyboard shortcut can be used to do this.
- **Zoom to 100%:** resets the scaling of the items and text shown on screen to the default size. The  button or **[F4]** keyboard shortcut can be used to do this.
- **Zoom in:** increases the scaling of the items and text shown on screen. The  button or **[CTRL] + [=]** keyboard shortcut can be used to do this.



## NOTICE

More information on the zoom controls of the procedure editor and the schedule editor can be found in *Chapter 10.6* and *Chapter 15.7*, respectively.

### 3.5.9 Electrochemical Frequency Modulation

This version of NOVA provides support for Electrochemical Frequency Modulation (EFM) measurements. These measurements are based on the application of a small amplitude voltage perturbation and recording of the electrochemical response of the cell. Using the measured data, corrosion rate information can be determined.

EFM measurements use a special two component sinewave modulation. During this type of measurements, the response from the cell at the applied frequency, higher harmonics of these frequencies and intermodulated frequencies are recorded. *Figure 49* shows a typical measurement.




Figure 49 Example of an EFM measurement



## CAUTION

Electrochemical Frequency Modulation measurements require a **FRA32M** module.



## NOTICE

More information on Electrochemical Frequency Modulation command can be found in *Chapter 7.6.4*.

### 3.5.10 Corrosion rate analysis

The **Corrosion rate analysis** command has been complemented with a new mode: **Polarization Resistance**. This analysis method is based on the **ASTM G59** standard and it uses the *Stern-Geary* equation to determine the corrosion current and the corrosion rate (see figure 50, page 55).




Figure 50 The Linear polarization method has been added to the Corrosion rate analysis command

Provided that the analysis is carried out in a low overpotential range with respect to the corrosion potential, the Linear polarization analysis method can provide a direct estimation of the corrosion current and corrosion rate, using user-defined Tafel slopes.



## NOTICE

More information on the **Corrosion rate analysis** command can be found in *Chapter 7.8.14*.

### 3.5.11 New Plots frame controls

The **Plots** frame now provides new controls that can be used to disable plots (see figure 51, page 56).




Figure 51 Disabling plots in the Plots frame



## NOTICE

Disabling plots can be done at any time.

It is also now possible to relocate the plot order or overlay plots by dragging the plots in the frame (see figure 52, page 57).






Figure 52 Rearranging the plot order



## NOTICE

It is not possible to relocate plot during a measurement.



## NOTICE

More information on the disabling of plots and the relocation of plots can be found in *Chapter 11.5.5* and *Chapter 11.8.8*, respectively.

### 3.5.12 Device drivers installation

The installation package of NOVA 2.1 now installs all required device drivers during the installation process, as described in *Chapter 1.3*.

The following drivers are installed:

- **Autolab device drivers:** required for using the Autolab potentiostat/galvanostat.
- **Metrohm device driver:** required for using any supported Metrohm liquid handling instrument.
- **Spectrophotometer device driver:** required for using any supported Autolab (or Avantes) spectrophotometer.






Figure 53 Schedules are now managed through the Library

The most recent **Schedules** are now also listed in the **Recent items** panel on the **Dashboard**, as shown in Figure 54.




Figure 54 The most recent Schedules are now listed in the Recent items panel



## NOTICE

More information on the **Library** can be found in *Chapter 6*.







Figure 56 The Custom tab provides additional controls for specifying the color

On the Custom tab, colors can be specified using RGB values or by changing the hue of the selected color or by selecting any available color in the provided RGB color matrix.

### 3.6.3 Data handling command shortcut button

A shortcut button has been added to this version of NOVA allowing data handling commands to be added to a procedure or data. The  shortcut button, located in the top right corner of the procedure editor, works in the same way as the data analysis shortcut button already available in the NOVA (see figure 57, page 62).




Figure 57 The data handling shortcut button can be used to add data handling commands easily



## NOTICE

The data handling commands shown in the popout menu depend on the selected command.



## NOTICE

More information on the use of the data handling shortcut button can be found in *Chapter 13*.

### 3.6.4 Library filters

To facility data management, the **Library** now provides filtering options that can be used to force the **Library** to display items that fit within the specified filter conditions. *Figure 58* shows an example using two filter conditions, one on the instrument serial number and one on the rating.

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 04 - Hydrodynamic linear sweep	Fe <sup>2+</sup> /Fe <sup>3+</sup> - NaOH 0.2 M	AUT71848	31-8-2015 13:55:57	14-6-2016 12:47:46	★★★★★	☐ ★ ☆ ☆ ☆ ☆
Demo 05 - Fe(II) - Fe (III) on pCP	Fe <sup>2+</sup> /Fe <sup>3+</sup> Reversibility Test - LSV with i	AUT71848	31-8-2015 14:40:14	14-6-2016 12:47:47	★★★★★	☐ ★ ★ ☆ ☆ ☆
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanostatic	AUT71848	31-8-2015 11:27:11	14-6-2016 12:47:46	★★★★★	☐ ★ ★ ★ ☆ ☆

Figure 58 The Library now provides filtering options for better data handling



## NOTICE

More information on the **Library** filters can be found in *Chapter 6.11*.

### 3.6.5 Extended Sampler information

To provide more information on how and when the **Sampler** records the signals during any electrochemical measurement, the **Sampler** editor has been extended with a table that provides more details on the sampling conditions. In the **Sampler** editor, shown in *Figure 59*, an additional **More** button is now available.

Signal	Sample	Average	d/dt
WE(1).Current	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
WE(1).Potential	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Power	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Resistance	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Charge	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
External(1).External 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrator(1).Charge	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrator(1).Integrated Current	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Time	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Sample alternating

More

Figure 59 The Sampler editor now provides more information








Figure 61 Zooming in and out is now possible for relevant analysis commands



## NOTICE

More information on this new option can be found in *Chapter 12.7.1* and in *Chapter 12.9.1* for the **Baseline correction** and **Electrochemical circle fit** commands, respectively.

### 3.6.7 Custom name for Build signal command

It is now possible to define a custom name to the signals generated by the **Build signal** command, as shown in *Figure 62*.




Figure 62 The Build signal command now offers the possibility to provide custom names



### 3.7.1 Procedure and data tags

It is now possible to assign tags to procedures and data files. This is a convenient tool for bookkeeping purposes. The controls for tagging data or procedures are provided in **Tags** sub-panel of the **Properties** panel of the procedure editor (see figure 63, page 67).

The screenshot shows the 'Properties' panel for a procedure named 'Demo 20 - Iron screw in seawater'. It features several input fields and buttons. The 'Procedure name' field contains the text 'Demo 20 - Iron screw i...'. The 'Remarks' field contains 'Linear sweep', 'voltammetry', and 'potentiostatic'. The 'Estimated duration' field is set to '0' with a unit 's'. Below these fields are two buttons: 'Edit End status' and 'More'. At the bottom, there is a 'Tags' section with a 'Rating' field showing five stars and an 'Add' button.

Figure 63 The Tags sub-panel provides controls for tagging data and procedure files

Two types of tags can be assigned to data or procedures:

- **Rating:** a rating based on a stars system can be assigned to each data or procedure file. By default, no stars are assigned to a file, but it is possible to change this at any time.
- **Tags:** text tags can be added to each data or procedure file. By default, no tags are assigned to a file, but it is possible to change this at any time.



## NOTICE

The rating and tags are updated when the file is saved.

It is also possible to provide a rating and tags directly from the Library (see figure 64, page 68).

Demo Database							+ -	
Name ▲	Remarks	Instrument	Rating	Tags	Measurement date	Last modified		
Demo 01 - Copper	CuSO4 0.01 M	AUT71848	☆☆☆☆☆	Add	31-8-2015 11:07:50	11-2-2016 14:3		
Demo 02 - Lead de	Pb(ClO4)2 0.0		★★★★☆	Add	4-2-2009 11:04:15	11-2-2016 14:29		
Demo 03 - Bipoten	RRDE measurt	MAC80064#3	☆☆☆☆☆	Add	15-7-2013 13:45:21	11-2-2016 14:30		
Demo 04 - Hydrod	Fe2+/Fe3+, N	AUT71848	☆☆☆☆☆	Add	31-8-2015 13:53:57	11-2-2016 14:3		
Demo 05 - Fe(II) - F	Fe2+/Fe3+ Re	AUT71848	☆☆☆☆☆	Add	31-8-2015 14:40:14	11-2-2016 14:3		
Demo 06 - Galvanc	Lead depositic	AUT71848	☆☆☆☆☆	Add	31-8-2015 11:27:11	11-2-2016 14:3		
Demo 07 - Chrono	Example of fa:	AUT71848	☆☆☆☆☆	Add	1-9-2015 13:20:24	11-2-2016 14:3		
Demo 08 - Superca	Supercapacito	AUT71848	☆☆☆☆☆	Add	1-9-2015 13:29:23	11-2-2016 14:3		
Demo 09 - Superca	Supercapacito	AUT50229	☆☆☆☆☆	Add	1-9-2015 13:50:29	11-2-2016 14:3		
Demo 10 - Differen	Differential pu	AUT50477	☆☆☆☆☆	Add	18-8-2015 15:11:45	11-2-2016 14:3		

Figure 64 Tags can also be defined in the Library



## NOTICE

More information on the rating and tagging of procedures and data can be found in *Chapter 6.8*.

### 3.7.2 New plot options

The plot options have been expanded in order to allow for additional control of the plotting of data (see *figure 65, page 69*).




Figure 65 New plot options have been added to the Axes sub-panel

The following options have been added to the **Axes** sub-panel:

- **Fixed scale (on/off):** defines if the axis should be automatically scaled or if a fixed scale should be used. When this property is switched on, it is possible to define a minimum and maximum value for the axis.
- **Custom ticks (on/off):** defines if custom major and minor ticks should be used for the axis. When this property is switched on, it is possible to define the distribution of major and minor ticks.






Figure 67 The number of recent items can be edited in the NOVA options

The default number of items is 5 and can be edited at any time.

### 3.7.5 Plot preview

NOVA now offers the possibility to assign one of the plots of a data file as a preview plot to display in the Library. This provides a quick preview of the data contained in each data file. The preview plot is shown in a tooltip (see figure 68, page 72).




Figure 68 A plot preview is now added to each data item in the Library



## NOTICE


Data measured with previous versions of NOVA will create a preview plot when changes to the file are saved in the current version.



## NOTICE

More information on the plot previews can be found in *Chapter 6.9*.

### 3.7.6 Print plot

It is now possible to print plots using the provided  button (see figure 69, page 73).






Figure 69 Plots can now be printed

A print preview dialog will be displayed, allowing finetuning of the print output (see figure 70, page 73).




Figure 70 A print preview dialog is shown






Figure 72 A tooltip shows the active electrochemical interface in bold



## NOTICE

More information on the optional ECI10M module can be found in *Chapter 16.3.2.8*.

### 3.7.9 ECI10M measurements

In order to more easily identify measurements carried out with the **ECI10M** as the active electrochemical interface, the **(ECI10M)** suffix will be shown in the procedure editor, next to the serial number of the active instrument, below the procedure title (*see figure 73, page 76*).




Figure 73 The ECI10M suffix is shown in the procedure editor

For measurements that have been carried out with the **ECI10M**, the same suffix will be added to the instrument serial number in the **Library**(see figure 74, page 76).




Figure 74 The ECI10M suffix is shown in the Library



## NOTICE

This suffix is not shown for measurements carried out with previous version of NOVA.

### 3.7.10 Library column display

The display settings of the grids used in the **Library** are now non-volatile. This change affects the following settings of the Library:

- **Column order:** the order in which the columns appear.
- **Column visibility:** the visibility of the available columns.
- **Sorting options:** the sorting options used in the Library.



#### NOTICE

The display settings used in the Library are stored on the local computer and can be defined for each type of Library location.



#### NOTICE

More information on the display settings used in the Library can be found in *Chapter 6.10* and *Chapter 6.13*.

### 3.7.11 Data grid column display

The display settings of the grids used in the **Data grid** are now non-volatile. This change affects the following settings of the data grid:

- **Column order:** the order in which the columns appear.
- **Column formatting:** the data formatting used in each column.
- **Sorting options:** the sorting options used in the data grid.



#### NOTICE

These settings are stored for each command in the data file.



#### NOTICE

More information on the display settings used by the data grid can be found in *Chapter 11.9*.






Figure 76 When commands are stacked, the Estimated duration takes underlying commands into account



## NOTICE

The Estimated duration is determined based on the interval time and the estimated number of points as well as the duration of underlying commands, if applicable.

### 3.7.13 Interpolate command

The new Interpolate command, available in the Analysis - general group of commands, is now available. This command can be used to determine Y or X value based on a user-defined X or Y value, by linear interpolation.



## NOTICE

More information on the Interpolate command can be found in *Chapter 7.8.6*.

### 3.7.14 Hydrodynamic analysis

The Hydrodynamic  $i \propto \sqrt{\omega}$  command has been renamed to **Hydrodynamic analysis** and it has been extended with the Koutecký-Levich analysis method .





### 3.8.2 Value of Alpha

For the **CV staircase** command and the **LSV staircase** command, the *Alpha value* advanced property is now available (see figure 77, page 81).




Figure 77 The Alpha value property is available for the CV staircase and LSV staircase command



## NOTICE

More information on the use of the value of Alpha can be found in Chapter 9.7.

### 3.8.3 Autolab RHD Microcell HC support

NOVA 2.0 introduces support for the Autolab RHD Microcell HC support (see figure 78, page 81).



Figure 78 The Autolab RHD Microcell system





Figure 80 NOVA introduces the support of the ECI10M module



## NOTICE

More information on the ECI10M module can be found in *Chapter 16.3.2.8*.

### 3.8.6 AC voltammetry

This new version of NOVA provides a new command for AC voltammetry measurements *AC voltammetry* (see *chapter 7.4.7, page 283*). A default procedure for this electrochemical method is also available *AC voltammetry* (see *chapter 8.3.6, page 551*).

## 4 Dashboard

The **Dashboard** is the home screen of NOVA. Whenever NOVA starts, the **Dashboard** is always shown to the user (see figure 81, page 84).




Figure 81 The Dashboard

At any time, when NOVA is used, it is possible to show the display by clicking the home tab (🏠). The Dashboard provides four different panels:

- **Actions:** this panel provides a list a shortcut buttons to trigger a common task in NOVA *Actions* (see chapter 4.1, page 85).
- **Recent items:** this panel provides a list of the last procedures, data and schedule items *Recent items* (see chapter 4.2, page 86).
- **What's going on:** this panel provides messages to the user about ongoing or finished events in NOVA *What's going on* (see chapter 4.3, page 88).
- **Instruments:** this panel provides a list of connected instruments *Instruments panel* (see chapter 4.4, page 90).

## 4.1 Actions

The **Actions** panel provides a series of buttons that can be used to quickly trigger a common action or control of the NOVA software (see figure 82, page 85).




Figure 82 The Actions panel provides shortcut to the most common actions in NOVA

The following shortcut buttons are provided:

- **Open library:** this button open the **Library**. More information on the **Library** can be found in *Chapter 6*.
- **New procedure:** creates a new blank procedure. More information on the procedure editor can be found in *Chapter 10.1*.
- **Import procedure:** imports a procedure from a *.nox* file in the **Library**. More information on the **Library** can be found in *Chapter 6*.
- **Import data:** imports NOVA data from a *.nox* file in the **Library**. More information on the **Library** can be found in *Chapter 6*.
- **Import command:** imports a command from a *.noi* file in the **My Commands** group of command. More information on the **My Commands** can be found in *Chapter 10.14*.
- **New schedule:** creates a new procedure schedule. More information on the **Procedure scheduler** can be found in *Chapter 15*.
- **Import schedule:** imports a procedure schedule from a *.nos* file in the **Library**. More information on the **Procedure scheduler** can be found in *Chapter 15*.






Figure 84 The Recent items panel is automatically updated when data or procedures are saved

It is possible to remove items from the **Recent items** panel by right-clicking an item and selecting the *Remove from recent items* option from the context menu (see figure 85, page 88).








Figure 86 The What's going on panel is used to provide messages to the user

At any time, it is possible to clear the **What's going on** panel using the **X** button (see figure 87, page 89).




Figure 87 Clearing the panel content







Figure 90 Clicking the refresh button will update the content of the Instruments panel

The content of the Instrument panel is updated when the  button is clicked (see figure 91, page 91).










Figure 91 The Instruments panel is refreshed




 <p>VIRT00001</p>	<p>This symbol is used to identify all <b>Autolab PGSTAT204</b> instruments. These instruments have a serial number starting with <b>AUT5</b>.</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all <b>Autolab PGSTAT101</b> instruments. These instruments have a serial number starting with <b>AUT4</b>.</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all <b>Multi Autolab Series</b> instruments (M101 and M204). These instruments have a serial number starting with <b>MAC8</b> (for the M101 Multi Autolab systems) and <b>MAC9</b> (for the M204 Multi Autolab systems).</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all <b>µAutolab type II</b> and <b>µAutolab type III</b> instruments. The µAutolab type II instruments are identified by a serial number starting with <b>µ2AUT7</b> and the µAutolab type III instruments are identified by a serial number starting with <b>µ3AUT7</b>.</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all <b>Autolab 7 Series</b> instruments (PGSTAT302, PGSTAT30, PGSTAT12, PGSTAT100) as well as the older <b>Autolab 9 Series</b> instruments (PGSTAT30, PGSTAT20, PGSTAT10 and PGSTAT100). These instruments are identified by a serial number starting with <b>AUT7</b> or <b>USB7</b>.</p>

The following tiles are used to identify the connected **Autolab RHD Microcell HC** controllers:





 <p>12</p>	<p>This symbol is used to identify all <b>Autolab RHD Microcell HC controllers</b> connected to the computer through a RS232 connection. These instruments are identified by their serial number (or device name).</p>
---	--

The following tiles are used to identify the connected **Autolab or Avantes** spectrophotometers:



 ASM80002	This symbol is used to identify all <b>Autolab or Avantes spectrophotometers</b> connected to the computer through a USB connection. These instruments are identified by their serial number (or device name).
---	--

The following tiles are used to identify the connected **Metrohm** devices:

 2429	This symbol is used to identify all <b>Metrohm 800 Dosino</b> devices connected to a USB controlled Metrohm device. These instruments are identified by their serial number (or device name).
 2358	This symbol is used to identify all <b>Metrohm 801 Magnetic Stirrers</b> or <b>804 Titration Stands</b> with a stirrer connected to it (either <b>Metrohm 802 Rod Stirrer</b> or <b>Metrohm 741 Magnetic Stirrer</b> ). These instruments are identified by their serial number (or device name).
 4589	This symbol is used to identify all <b>Metrohm 814, 815 or 858 Sample Processor</b> devices connected by USB to the host computer. These instruments are identified by their serial number (or device name).
 41774	This symbol is used to identify all <b>Metrohm 6.2148.010 Remote Box</b> devices connected to a USB controlled Metrohm Device. These instruments are identified by their serial number (or device name).

In *Figure 92*, two instruments are connected (a Multi Autolab system with Serial Number MAC91234 and an Autolab N Series instrument with serial number AUT81234). A virtual instrument, with serial number VIRT00001 is also connected. This instrument is identified as a PGSTAT204.

The serial number of the N series instrument is shown in bold underlined font (**AUT81234**) indicating that this is the *default* instrument.



## NOTICE

When a measurement is started, it will always be executed on the *default* instrument, unless otherwise specified.

The following actions can be performed in the **Instruments** panel:

- Change the *default* **Autolab** instrument.
- Open the instrument control panel.

## 5.1 Change the default instrument

The default instrument, displayed in the bold underline in the **Instruments** panel, is the instrument used in any measurement, unless otherwise specified. This is also the instrument used for procedure validation purposes.

To change the default instrument, right-click any instrument tile in the **Instruments** panel and select the *Make [Instrument serial number] the default instrument* from the context menu (see figure 93, page 95).




Figure 93 Defining the default instrument



## NOTICE

Only one instrument can be set as default instrument.





- **Instrument information panel:** this panel displays information about the instrument.
- **Tools panel:** this panel provides quick access to the hardware setup and a number of direct measurement tools like current interrupt and positive feedback.
- **Autolab display panel:** this panel provides a number of manual controls of the instrument.




## NOTICE

The available channels in a multi channel Autolab are listed in the **Channels** sub-panel. Each channel is identified by a letter or a number. More information is provided in *Chapter 16.2.5*.

### 5.2.1 Instrument information panel

The **Instrument information** panel shown in the instrument control panel provides information on the selected instrument (*see figure 96, page 97*).



*Figure 96 The Instrument information panel*

This information is updated in real time and is provided for information only. The following items are listed:

- **State:** indicates the state of the instrument (idle or measuring).
- **Instrument type:** indicates the type of instrument.
- **Modules:** shows the extension modules of the instrument.
- **Embedded processor:** shows the type of embedded processor installed in the instrument (IF030 or IF040).
- **Embedded software:** shows the embedded application name and version number.






## NOTICE

The i-Interrupt, Positive feedback, pH calibration and Reset integrator drift tools are only shown on if the instrument provides the functionality used by these tools.

### 5.2.2.1 Hardware setup

The  button can be used to edit the **Hardware setup**. The hardware setup screen shows three panels (*see figure 98, page 99*).




Figure 98 The hardware setup

The following panels are provided:

- **Autolab module panel:** used to specify the type of Autolab and additional properties of this instrument.
- **Additional modules panel:** provides a list of compatible extension modules which can be installed in the instrument or connected to the instrument.
- **Properties panel:** provides additional parameters for the extension modules.

#### 5.2.2.1.1 Autolab module panel

The **Autolab module** panel can be used to specify the following properties (*see figure 99, page 100*):

- **Main module:** specifies the type of Autolab using the provided drop-down list.
- **Power supply frequency:** specifies if the mains frequency is 50 or 60 Hz.







Figure 100 Optional module can be selected in the optional modules panel



## NOTICE

Only the FRA2, BA, FI20 - Integrator, Booster10A, Booster20A, External Devices, IME303, IME663 and MUX module have additional properties to display in the **Properties** panel.

### 5.2.2.2 Current interrupt

The  button can be used to perform a current interrupt (i-Interrupt) measurement. This tool can be used to determine the uncompensated resistance,  $R_U$ .



## NOTICE

This tool is not available for  $\mu$ Autolab type II and type III instrument as well as the Autolab PGSTAT10.

During a current interrupt measurement, a constant potential is applied on the cell before the current interrupt circuit is triggered. This circuit interrupts the current flow in the cell and measures the potential decay. From the measured potential decay, the uncompensated resistance ( $R_u$ ) value is determined, using a linear and an exponential regression.

Two values of the uncompensated resistance,  $R_u$ , are determined automatically at the end of the measurement:

- **Ru linear:** this value is obtained from a linear regression performed on the initial segment of the voltage decay.
- **Ru exponential:** this value is obtained from an exponential regression performed on the initial segment of the voltage decay.

Proper determination of this value requires an accurate measurement of the current. The measurements must therefore be carried out at a potential value where the current is high enough to be measured properly and the current range must be adjusted in accordance.



## NOTICE

For accurate measurements, the current should be at least in the order of 1 mA.

When the i-Interrupt tool is used, the control screen for this tool will be displayed (*see figure 101, page 103*). The control screen provides two panels and one plot area.




Figure 101 The *i-Interrupt* tool

The **Settings** panel shows the properties used in the current interrupt measurement (see figure 102, page 103).




Figure 102 The *i-Interrupt* Settings panel

The following properties and controls are available:

- **High speed:** a  toggle that can be used to switch the high speed ADC module (ADC10M or ADC750) off or on (default off).








Figure 103 The measured and fitted data

The measured data points are shown as a point plot. The linear regression is shown using a green line and the exponential regression is shown as a red line. The start and end value of the two regression methods are shown using vertical lines with matching colors.

It is possible to hide or show the measured data or the regression data by checking or unchecking the check boxes shown in the legend (see figure 104, page 105).




Figure 104 Using the check boxes to show or hide the data measured during the current interrupt





## NOTICE

The copied value can be pasted in a suitable property field in NOVA.

### 5.2.2.3 Positive feedback

The Positive feedback button can be used to perform a positive feedback measurement. This tool can be used to determine the uncompensated resistance,  $R_u$ .



## NOTICE

This tool is not available for  $\mu$ Autolab type II and type III instrument as well as the Autolab PGSTAT10.

During a positive feedback measurement, a potential pulse is applied on the cell and the potential is recorded. The  $iR$  compensation value can be adjusted upwards manually until its value is close to the actual value of the uncompensated resistance,  $R_u$ . When the compensated resistance reaches a value close to the actual value of  $R_u$ , potentiostatic loop will start to ring. When the compensated resistance exceeds the  $R_u$  value, the potentiostatic loop is no longer stable and the instrument will oscillate.

When the positive feedback tool is used, the control screen for this tool will be displayed (see figure 108, page 107). The control screen provides two panels and one plot area.




Figure 108 The positive feedback tool

The **Settings** panel shows the properties used in the current interrupt measurement (see figure 109, page 108).






Figure 110 The measured data

The measured data shows the potential profile applied on the cell. Since the positive feedback tool uses an iterative approach, it is possible to adjust the value of the iR compensation value property and repeat the measurement (see figure 111, page 109).




Figure 111 The measured data and previous results shown in overlay

The Show previous results  toggle provided in the **Results** panel can be used to show or hide the data from the previous measurement (see figure 112, page 109).



## NOTICE

The tool only stores the current measured data and the data from the previous measurement.

### Results

Show previous results

Figure 112 The previous data can be enabled and disabled



When the check cell tool is used, the control screen for this tool will be displayed *Figure 114*. The control screen provides two panels and one plot area.




Figure 114 The check cell tool

The **Settings** panel shows the properties using the cell check measurement (see *figure 115, page 111*).

Figure 115 The check cell Settings panel

The following properties and controls are available:

- **Current range:** the current range in which the cell check is performed.
- **Bandwidth:** a drop-down control that can be used to specify the bandwidth of the instrument (high stability, high speed or ultra-high speed).
- **Potential/Current:** a numeric field that can be used to specify the applied potential (in potentiostatic mode) or the applied current (in galvanostatic mode).

