Prospektüs

Prospectus

İPEK - SILK

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Instructions for Use

SILK

NONABSORBABLE, SURGICAL SILK SUTURE, USP/ EP/TS5505

DESCRIPTION

Silk suture is a nonabsorbable, sterile, surgical suture composed of an organic protein called fibrin. This protein is derived from the domesticated species Bombyx mori (B. mori) of the family Bombycidae. Silk sutures are processed to remove the natural waxes and gums. Above USP 7/0 it is braided, USP 7/0 and 8/0 available braided or twisted. Silk suture is dyed black and coated with paraffin wax or silicone. Black dye is Hematein HCK (C.I. color index no. 75290). Silk suture is also available in its natural color, ivory. Virgin silk suture is available in which the sericin gum is not removed and serves to hold the filaments together. Virgin silk coloured blue by Methylene Blue (C.I.# 52015) or black by Hematein.

Silk suture meets requirements established by the United States Pharmacopoeia (U.S.P.) for Nonabsorbable Surgical Suture Monogram and the requirements of Braided Silk Suture Monogram established by European Pharmacopoeia (E.P.).

INDICATIONS

Silk suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures.

ACTIONS

Silk suture elicits an acute inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. While silk sutures are not absorbed, progressive degradation of the proteinaceous silk fiber in vivo may result in gradual loss of all of the suture's tensile strength over time.

CONTRAINDICATIONS

The use of this suture is contraindicated in patients with known sensitivities or allergies to silk. Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, silk should not be used where permanent retention of tensile strength is required.

WARNINGS

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Silk suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Do not re-sterilize. Discard opened packages and unused sutures.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

Skin sutures which remain in place for periods of longer than seven days may cause localized topical irritation and the extended portion of suture may be snipped off after five to seven days.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Silk sutures, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of Silk suture include, allergic response in patients that are known to be sensitive to silk, initial inflammatory tissue reaction and transitory local irritation at the wound site and infection.

HOW SUPPLIED

Silk sutures are available in U.S.P sizes 10/0 to 5 (metric sizes 0.2 - 7.0) in a variety of lengths and affixed to needles or non-needled. Silk sutures are available in one, two, and three dozen boxes. Silk suture is supplied sterile. Silk suture is for single use only.

STORAGE

Store below 25°C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

Do not reuse

Do not resterilize

Do not use if package is damaged

Manufacturer

YYYY-MM Date of Manufacture, Year

YYYY-MM Expiry Date, Year-Month

STERILE EO Sterile EO: Ethylene oxide

LOT Batch Number

Dyed, Non-absorbable, Braided, Coated

Undyed, Non-absorbable, Braided, Coated

Undyed, Non-absorbable, Twisted, Coated

Dyed, Non-absorbable, Twisted, Coated

REF Catalogue Number

Store below 25°C

Keep away from sunlight

Protect from humidity

Recyclable pack

Attention, See instruction for use

CE 1014

IFU-SB-rev-04-05-05-2017

Issue date: 11.09.2012

"EASSI (The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

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