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Installation of an Eye protection filter into a microscope

An instruction how to install the eye protection filter offered by ARC laser.

Parts : 1 Filter A.R.C. Laser BG 03017 / BG 03018 /
BG 06063 / BG06064

Instruction how to install the filter

WARNING:

Do not use the filter as an eye protection device without having the filter properly installed into a microscope as described as follows.

The main body of the eye protection filter is made of aluminum. At both ends, a standard dovetail fitting is available for various kind of microscopes. Take care when removing the filter from the box. Do not touch the glasses inside the main body. Fingerprints may change reflectivity of the filter glass. In case of an unintentional touch of the glass, use Methanol (ref: Lab. grade) and lens tissue to wipe the fingerprint or any other dirt of the glass.

WARNING:


Do not use any soap or detergent, which contains oil or fat.

At the time of installation of the filter, the laser has to be switched off. No attachment or the laser should be available.

Remove the main body with the filter glass from the box and place it close to the microscope. Remove the Binoculars from the microscope, by turning the upper silver screw counterclockwise on the main body of the microscope. Different microscopes may have different screws to hold the binoculars. Remove the binoculars from the main body and place the filter into the main body. Tighten the screw of the main body to hold the filter firmly in position.

Take the binoculars, place the binoculars into the filter main body, and use the hex screw on the filter body to tighten the binoculars firmly to the filter body. At this time, the filter is firmly placed into the microscope.

The laser can be switched on and use in connection with this filter.

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
Removing the protection filter at a laser slit lamp.

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4) Open screw #2 to remove the protection filter (5)

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6) Add the ocular again and close screw #2

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Open the Joint of the ocular to equip your microscope with a protection filter.

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2) Open the screw #1 to remove the ocular. Insert the protection filter according to the image (3)

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4) Close screw #1 to fix the filter. Open screw #2 to give space for the ocular for easy fit in (5)

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