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02 | 03 LENSTAR 900

Lenstar 900

The all-in-one optical biometer & IOL planning platform

From automated and precise measurement acquisition to IOL calculation and surgical planning, Lenstar 900 provides a seamless workflow that boosts your practice's efficiency and optimizes your patients' refractive outcomes.

With a single scan, Lenstar 900 captures axial dimensions, keratometry, and pupillometry, while simultaneously generating a cornea map* to assess the regularity and symmetry of astigmatism. The dual zone keratometry is equivalent to manual keratometry, providing precise K values, axis, and astigmatism measurements^{1,2}.

Equipped with state-of-the-art formulas such as Hill-RBF, Olsen, Barrett, and the EyeSuite IOL toric planner, Lenstar 900 ensures accurate IOL predictions for any patient.

**Availability may vary by region. Please contact your local Haag-Streit Distributor for details.*

Precise measuring data for the Hill-RBF Method

The Lenstar 900's precise measurement technique in combination with the unique Hill-RBF Method allows accurate prediction of the IOL power for all types of eye.

Perfect K values – best toric results

The Lenstar 900 features dual zone keratometry or T-Cone topography for precise astigmatism and axis measurement^{1,2}. Its toric calculator, featuring Barrett³, Hill-RBF⁴, and Olsen⁵, predicts toric IOL power, considering the posterior cornea for best refractive outcomes.

For post-refractive cases – quick & reliable

Barrett True-K, Shammass No-History, and Masket IOL calculation methodologies can be used for post-refractive patients even without any clinical history available^{6,7,8}.


04 | 05 LENSTAR 900 MEASUREMENTS


Optical biometry **Fast & precise measurement of the anterior & posterior segments of the eye**

The Lenstar 900 delivers fast and precise optical biometry for a variety of eye types, including regular, aphakic, pseudophakic, and silicone oil-filled eyes.

Its optional Automated Positioning System (APS) utilizes dynamic eye-tracking to facilitate automated and accurate measurement acquisition.

If necessary, each measurement can be evaluated and adjusted to ensure complete biometric accuracy.

Axial dimensions

Using Optical Low Coherence Reflectometry (OLCR), Lenstar 900 precisely measures all key axial dimensions such as central corneal thickness (CCT), anterior chamber depth (ACD), lens thickness (LT), and axial length (AL) with high precision and reproducibility.

The Dense Cataract Measurement (DCM) mode uses advanced digital signal processing to achieve high success rates in AL measurements, even in patients with dense cataracts.

White-to-white & pupillometry

Based on high-resolution color photography of the eye, Lenstar 900 provides reliable white-to-white measurements to be used in advanced IOL calculation formulas. In addition, Lenstar 900 includes measurement of the pupil diameter, an indicator of the patient's suitability for premium IOLs, as well as for laser refractive procedures.

Keratometry

Lenstar 900's unique dual zone keratometry provides axis and astigmatism measurements equivalent to the "gold-standard" manual keratometry^{1,2}. The closely spaced 32 measurement point pattern improves precision, both delivering more data and minimizing the need for software data interpolation.

Cornea map*

Simultaneously with the keratometry measurement, the Lenstar 900 generates an axial curvature map of the central 3 mm zone of the anterior cornea using 48 measurement points. The cornea map indicates the astigmatism's regularity and symmetry and therefore provides decision support in choosing the appropriate IOL type for toric candidates.

T-Cone topography

For even greater confidence in toric procedures, the optional T-Cone toric platform offers a mountable true 11-ring Placido topographer, delivering detailed axial, tangential, and elevation maps of the central 6 mm optical zone. Integrated with the EyeSuite IOL toric planner, the T-Cone toric platform ensures optimal intervention planning, improving surgical outcomes and accuracy.


**Availability may vary by region. Please contact your local Haag-Streit Distributor for details.*

06 | 07 IOL CALCULATION


EyeSuite IOL The comprehensive IOL calculation platform

EyeSuite IOL is designed to meet the needs of all patients, including those with toric or post-refractive requirements. Its intuitive dashboard ensures accurate IOL power calculations and streamlines lens selection. The platform enables side-by-side comparisons of IOL calculation methods and anticipated outcomes for monofocal, toric, or multifocal lenses across manufacturers.

The IOL editors streamline management and customization of the IOL database, allowing surgeons to adjust IOL constants based on their outcomes and improve calculation accuracy. The database can also be expanded with IOLs from new manufacturers or external libraries like IOL Con, ensuring maximum flexibility and adaptability to individual preferences.

90°

45°

0°

315°

270°

Hill-RBF 3.0

Hill-RBF 3.0 uses the predictive power of artificial intelligence, built on high-quality data collected by leading cataract surgeons worldwide. By leveraging pattern recognition and big data, it delivers highly accurate IOL predictions for any anatomical configuration⁴. Its unique reliability check identifies unusual eye configurations and labels them as outliers, providing confidence beyond simple power values.

Barrett Universal II & Olsen

Standard IOL calculations work well for average eyes but often fall short in cases like short or long eyes or those with unique anatomical features, as they tend to misinterpret lens positioning. To enable excellent refractive results for any patient anatomy, EyeSuite IOL features latest generation IOL calculation formulas, i.e., Barrett Universal II³ and Olsen⁵. For improved IOL prediction, they include anterior chamber depth (ACD) and lens thickness (LT) in their calculation besides the standard axial length measurements and K-readings.

Post-refractive solutions

IOL calculation for post-refractive patients presents unique challenges. To support surgeons managing patients with prior RK, LASIK, RK or PRK, EyeSuite IOL offers a comprehensive set of cutting-edge formulae for post refractive patients. For patients without surgical history, the integrated Barrett True-K and Shammas No-History methods are employed. If the change in refraction caused by the refractive procedure is known, then the Barrett True-K with history, Masket and modified Masket formulae may also be used^{6,7,8}.

Toric solutions

For planning toric IOLs, correct calculation of the cylinder power and orientation are key. To achieve this, the optional EyeSuite IOL toric planner features the Barrett Toric Calculator, the Hill-RBF, Koch/Abulafia toric planner⁹, and Olsen Toric⁵ for accurate IOL prediction.

08 SURGICAL PLANNING

EyeSuite IOL toric planner Enhance efficiency & refractive outcomes

EyeSuite IOL toric planner is an intuitive platform designed to simplify toric intervention planning. It displays the implantation axis, the incision location and overlays of critical data such as IOL orientation or curvature maps. The incision optimization tool helps surgeons identify the ideal incision location to minimize residual astigmatism.

To enhance intraoperative orientation, surgeons can manually mark anatomical landmarks on the iris and sclera images. They either serve as a baseline point for the intraoperative orientation or as a fallback strategy if external marking is not successful.

All critical planning details—including IOL data, calculations, and custom markings—are consolidated into a tailored printout, providing clear, actionable guidance for the operating room.


Connectivity is key Seamless integration for optimal workflow

With EyeSuite software, the Lenstar 900 seamlessly integrates into your practice network, providing real-time access to all data.

This integration not only saves time but also reduces the risk of transcription errors, ensuring efficiency and accuracy throughout your workflow.

Standardized interfaces, such as GDT or DICOM, enable effortless connection to most electronic medical record (EMR) systems. Additionally, EyeSuite's open data interface and Lenstar 900's separate computer automatically populate data fields with imported values from refraction devices or your EMR.

“The Lenstar 900 is a remarkably easy to use all-in-one IOL power calculation tool that delivers exceptionally accurate axial length, anterior chamber depth & lens thickness by optical biometry. At the same time, its dual zone autokeratometry feature is precise & uniformly consistent. The Lenstar 900 is an excellent choice for surgeons migrating towards torics & other premium channel IOLs where highly accurate outcomes are critical for success.”

WARREN E. HILL, MD, FACS

Technical specifications

Lenstar LS 900

Measured variables & modes

Corneal thickness ^{CT}

Measurement range	300–800 µm
Display resolution	1 µm

Anterior chamber depth ^{ACD}

Measurement range	1.5–6.5 mm
Display resolution	0.01 mm

Lens thickness ^{LT}

Measurement range	0.5–6.5 mm
Display resolution	0.01 mm

Axial length ^{AL}

Measurement range	14–32 mm
Display resolution	0.01 mm

White-to-white distance ^{WTW}

Measurement range	7–16 mm
Display resolution	0.01 mm

Keratometry ^K

Measurement range	
for radius	5–10.5 mm
Display resolution	0.01 mm
Measurement range	0–180°
for axis angle	
Display resolution	1°

Pupillometry ^{PD}

Measurement range	2–13 mm
Display resolution	0.01 mm

Measurement modes

- 'Normal' eye
- Aphakic eye
- Pseudophakic eye
- Silicone-filled eye
- Combination of the above.

Laser safety

Class 1 laser product

Onboard IOL calculation formulae

Hill-RBF Method, Hill-RBF/Abulafia-Koch Toric Calculator, Olsen, Barrett Universal II, Barrett True-K, Barrett True-K Toric, Barrett Toric, Haigis, Hoffer Q, Holladay 1, SRK/T, SRK II, Masket, Modified Masket, Shammas No-History.

IOL calculation data interfaces

- Holladay IOL Consultant Professional Edition (Holladay 2 formula & Holladay toric calculator)¹⁰
- PhacoOptics (Olsen formula)¹¹
- Okulix (Ray-Tracing by Prof. Preussner)¹².

Electronic medical record system interfaces

- DICOM (SCU)
- EyeSuite Script Language
- GDT
- EyeSuite command line interface.

The above-mentioned measurement ranges are based on the standard settings of the device for automatic measurement and analysis.

Lenstar LS 900 optical biometer intended use

The Lenstar LS 900 is a non-invasive, non-contact OLCR (optical low-coherence reflectometry) biometer used for obtaining ocular measurements and performing calculations to assist in the determination of the appropriate power and type of IOL (intraocular lens) for implantation after removal of the natural crystalline lens.

The Lenstar LS 900 measures:

- Axial length
- Corneal thickness
- Anterior chamber depth
- Aqueous depth
- Lens thickness
- Corneal curvature
- Radii for flat and steep meridian
- Axis of the flat meridian
- White-to-white distance
- Pupil diameter.

1 Hill W, Osher R, Cooke D, Solomon K, Sandoval H, Salas-Cervantes R, Potvin R. Simulation of toric intraocular lens results: manual keratometry versus dual zone automated keratometry from an integrated biometer. *J Cataract Refract Surg.* 2011 Dec; 37(12): 2181-7. **2** Gundersen KG, Potvin R. Prospective study of toric IOL outcomes based on the Lenstar LS 900 dual zone automated keratometer. *BMC Ophthalmol.* 2012 Jul 16; 12:21. **3** Barret G. Flight of the arrow: improving outcomes with toric intraocular lenses, XXXII Congress of the ESCRS 2014 London, Video Prize Winner, Category Innovation, available on ESCRS on demand, <http://escrs.conference2web.com/content/23476/> accessed March 16, 2015. **4** W.E. Hill; IOL Power Selection by Pattern Recognition; ASCRS EyeWorld Corporate Education; ASCRS 2016. **5** Olsen T. Improving IOL power Calculation by measurement of the lens thickness with the Lenstar LS 900 presented at the ESCRS in Paris 2010. **6** Wang L, Hill WE, Koch DD. Evaluation of intraocular lens power prediction methods using the American Society of Cataract and Refractive Surgeons Post-keratorefractive Intraocular Lens power Calculator. *J Cataract Refract Surg.* 2010 Sep; 36(9): 1466-73. **7** McCarthy M, Gavanski GM, Paton KE, Holland SP. Intraocular lens power calculations after myopic laser refractive surgery: a comparison of methods in 173 eyes. *Ophthalmology.* 2011 May; 118(5): 940-4. **8** Wang L, Tang M, Huang D, Weikert MP, Koch DD. Comparison of Newer Intraocular Lens Power Calculation Methods for Eyes after Corneal Refractive Surgery. *Ophthalmology.* 2015 Dec;122(12):2443-9. **9** Abulafia A, Barrett GD, Kleinmann G, Ofir S, Levy A, Marcovich AL, Michaeli A, Koch DD, Wang L, Assia EI. Prediction of refractive outcomes with toric intraocular lens implantation. *J Cataract Refract Surg.* 2015 May;41(5):936-44. **10** <https://www.hicsoap.com/> **11** <https://www.phacooptics.net/> **12** <http://okulix.de>

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6. Edition / 2013 – 09

Introduction

Thank you for choosing a HAAG-STREIT appliance. Provided you comply carefully with the regulations in these instructions for use, we can guarantee the reliable and unproblematic use of our product.

WARNING!

Read the instructions for use carefully before commissioning this device. They contain important information concerning the safety of the user and patient.

Purpose of use

The Lenstar LS 900 is a non-invasive, non-contact OLCR (optical low-coherence reflectometry) biometer for measuring various ocular parameters necessary for calculating and determining suitable IOLs (intraocular lenses) for patients, for implantation after removal of the natural crystalline lens. The Lenstar LS 900 measures:

- Axial eye length
- Corneal thickness
- Anterior chamber depth
- Aqueous depth
- Lens thickness
- Corneal curvature
- Radius of corneal curvature meridian from flat to steep
- Axes of the meridians
- White-to-white distance
- Pupil diameter

Contents

1. Safety	4
1.1 Areas of application of the device.....	4
1.2 Patient population.....	4
1.3 Ambient conditions.....	4
1.4 Shipment and unpacking.....	4
1.5 Installation warnings.....	4
1.6 Operation, environment.....	4
1.6.1 Plausibility of the measurements.....	5
1.6.2 IOL calculation.....	6
1.6.3 References.....	6
1.6.4 IOL constants.....	6
1.6.5 IOL constants derived from data from an optical biometer.....	7
1.6.6 IOL constants derived using data from an immersion ultrasound biometer.....	7
1.6.7 IOL constants derived from data from a contact ultrasound biometer.....	7
1.7 Optical radiation.....	7
1.8 Disinfection.....	7
1.9 Warranty and product liability.....	7
1.10 Symbols.....	7
2. Introduction	8
2.1 Basic construction.....	8
2.2 Examination components (LS 900).....	8
2.3 Control component (PC).....	9
2.4 Instrument table (option).....	9
3. Appliance assembly / installation.....	9
3.1 Safe system configuration according to EN 60601-1.....	9
3.2 Instrument table (option) assembly.....	10
3.2.1 Mounting the roller rails/slide plate.....	10
3.2.2 Mounting the sticker for automatic position recognition.....	11
3.3 Transporting the appliance.....	11
3.3.1 Device connection, cables under the device.....	11
3.3.2 Device connection, cables at side (option).....	12
3.4 Computer connection.....	12
3.4.1 Installing a short forehead rest when using the optional T-Cone.....	12

4. Operation	12	7. Software / Help menu / Error messages.....	17
4.1 Position of patient during measurement	12	8. Maintenance.....	17
4.2 Measuring with the optional T-Cone	13	8.1 Function check / zero adjustment.....	18
4.2.1 Fitting the T-Cone.....	13	8.2 Cleaning	18
4.2.2 Measuring with the optional T-Cone	13	8.3 Maintenance of T-Cone (option).....	18
4.2.3 Removing the T-Cone.....	14	A. Appendix.....	18
4.3 Fixation	14	A.1 Accessories / spare parts	18
4.4 Measured variables	14	A1.1 LENSTAR LS 900.....	18
4.4.1 A-Scan	14	A1.2 T-Cone (option).....	18
4.4.2 Keratometry.....	14	B. Legal regulations.....	18
4.4.3 White-to-white distance.....	14	C. Classification	19
4.4.4 Pupillometry & visual axis.....	14	C.1 LENSTAR LS 900 biometer.....	19
5. Commissioning	15	C.2 T-Cone (option).....	19
5.1 Switching on the appliance.....	15	D. Disposal	19
5.2 Switching off the appliance.....	15	E. Standards.....	19
6. Technical data.....	15	F. EMC supplement	20
6.1 Power supply (primary side)	15	F.1 General.....	20
6.1.1 Device (secondary side).....	15	F.2 Emitted interference (standard table 1).....	20
6.2 Eye length measurement (A-Scan) and central fixation	15	F.3 Immunity (standard table 2).....	21
6.2.1 Keratometry	15	F.4 Immunity on non-life support devices (standard table 4).....	22
6.2.2 Illumination	15	F.4 Safe distances on non-life support devices (standard table 6).....	23
6.2.3 Positioning aid (as of serial number 2000).....	15		
6.2.4 Technical data of T-Cone (option).....	15		
6.3 Reproducibility	15		
6.4 Measuring accuracy	16		
6.4.1 Reproducibility table B.....	16		
6.5 Measured variables	16		
6.5.1 Corneal thickness (CT).....	16		
6.5.2 Anterior chamber depth (ACD).....	16		
6.5.3 Lens thickness (LT)	16		
6.5.4 Axial length (AL)	16		
6.5.5 Keratometry (R).....	16		
6.5.6 White-to-white distance (WTW).....	16		
6.5.7 Pupillometry.....	16		
6.5.8 Study design.....	17		

1. Safety

FORBIDDEN!

Failure to comply with the instructions can result in material damage and pose a danger to users and patients.

WARNING!

These warnings must absolutely be complied with to guarantee safe operation of the device and to avoid any danger to users and to patients.

NOTE!

Important information: please read carefully.

1.1 Areas of application of the device

The device is intended for use in doctor's practices, hospitals and optometrists' and opticians' premises.

1.2 Patient population

The patient must be capable of sitting up straight and keeping his head still. He must be physically and mentally able to cooperate well and mentally capable of following the examination. Patients must be at least 6 years old.

1.3 Ambient conditions

Transport:	Temperature	from -40°C	to +70°C
	Air pressure	from 500 hPa	to 1060 hPa
	Relative humidity	from 10%	to 95%
Storage:	Temperature	from -10°C	to +55°C
	Air pressure	from 700 hPa	to 1060 hPa
	Relative humidity	from 10%	to 95%
Use:	Temperature	from -10°C	to +35°C
	Air pressure	from 800 hPa	to 1060 hPa
	Relative humidity	from 30%	to 75%

1.4 Shipment and unpacking

- Before unpacking the appliance, check whether the packaging shows traces of improper handling or damage. If this is the case, notify the transport company that has delivered the goods to you.
- Unpack the equipment together with a representative of the transport company. Make a report of any damaged parts. This report must be signed by you and by the representative of the transport company.
- Leave the device in the packaging for a few hours before unpacking it (condensation).

- Check the appliance for damage after it is unpacked.
- Return defective appliances in the appropriate packaging.
- Store packaging material carefully, so that it can be used for possible returns or when moving.
- Check that the contents of the packaging correspond to the contents specified on the leaflet included in the packaging.

1.5 Installation warnings

WARNING!

- Only use a HAAG-STREIT approved external medical power supply (EN 60601-1, UL 60601-1, CAN/CSA C22.2 No 601.1-M90).
- The plug, cable and socket must be functioning perfectly.
- Before carrying out maintenance or cleaning work, the device must always be disconnected from the mains by unplugging the power supply or plug.
- Computers and further ancillary devices (printers, etc.) must comply with the EN 60601-1 standard, or be connected through galvanic isolation to external networks (isolating transformer, galvanic Ethernet isolator, etc.), see also "Safe system configuration in accordance with EN 60601-1".
- For connection to PC, only use the supplied USB cable (2 m).
- The power supply must be positioned in such a way that proper heat dispersion is guaranteed.

1.6 Operation, environment

FORBIDDEN!

- It is expressly forbidden to open the device!
- Never use the device in potentially explosive environments where volatile solvents (alcohol, benzine, etc.) and combustible anaesthetics are in use.

WARNING!

- The doctor or the operator is obliged to inform the patient about the safety instructions concerning him and to ensure that these instructions are complied with.
- Only staff trained and experienced in the generation of measuring data and calculations and the manual entry, editing and deletion of data may examine patients, operate the device and interpret results.
- All users must be appropriately trained and familiarised with the contents of the instructions for use, especially with regard to the safety instructions contained therein.
- Measurements can be carried out with dilated or undilated pupils. Dilatation of the pupil only influences pupillometry.

- We recommend checking the calibration of the T-Cone when fitting it and repeating it if necessary (see section "Function check").

NOTE!

- This appliance may only be used for the purpose described in these instructions for use.
- Installation only by trained specialists.
- The PC on which the EyeSuite software is installed may not contain any other software which could restrict the correct operation of EyeSuite.
- Please switch the computer off if it is not used for a long time.
- The appliance may not be transported, stored or operated outside the specified ambient conditions (see section "Ambient conditions").
- The appliance is used in a room in the medical area with attenuated light.
- If exposed to the effects of an external force (e.g., by being accidentally knocked or dropped), this type of sensitive measuring device must be promptly checked according to section "Function check" and, if necessary, returned to the factory for repair.
- If the device is moved/transported, a functional check must be performed in accordance with section "Function check".
- Perform a zero adjustment and function check when prompted by the software.
- The selected measuring mode must be checked before every measurement.
- Keep these instructions for use in a place where they are accessible to those working with the device at all times.
- Warranty claims can only be made if the instructions in these instructions for use are complied with.
- The manufacturer of the appliance is not liable for loss or damage due to unauthorised handling of the same. All warranty claims arising in this case are null and void.
- Always remove the dust cover before switching the appliance on. The light sources can be destroyed by overheating. Likewise, make sure that the appliance is switched off before attaching the dust cover.
- Repairs may only be conducted by suitably trained and authorised specialist personnel. Incorrect repairs can pose considerable risks for operating staff and patients.
- Only original spare parts and original accessories may be used for repairs.
- The software must be installed by trained personnel.

- The optional T-Cone may only be used with the Lenstar.
- The optional T-Cone may only be used with Lenstar devices with a serial number of ≥ 2000 or Lenstar devices converted to white light illumination.
- The optional T-Cone must be checked for damage before use.
- The T-Cone must not be exposed to strong direct sunlight.

1.6.1 Plausibility of the measurements

WARNING!

- Users must check measurement readings for plausibility. This includes the checking of the A scan and the cursors, which automatically adjust to the signal, the keratometry values, the white to white distance and the pupillometry, whenever one of the measurements displays an unusually high standard deviation. The operator must also take into account the type (e.g., posterior subcapsular cataract) and density of the cataract when evaluating plausibility.
- Prior to the measurement, the user must verify that the patient is not wearing contact lenses. Wearing contact lenses will result in erroneous measuring results.
- It may not be possible, under certain circumstances, to carry out measurements on persons with fixation problems.
- In cases of thick cataracts and uncertain measurement of the axial length, ultrasound biometry should be performed as a control examination.
- Dense lenticular opacities may make it impossible to measure the axial eye length and lens thickness.
- Pronounced opacities of the central cornea can likewise make it impossible to measure corneal thickness, anterior chamber depth, lens thickness or axial eye length.
- Blood in the vitreous may make it impossible to measure the axial eye length.
- Keratometry may be erroneous in eyes that have undergone keratorefractive surgery because such eyes may significantly deviate from spherical surfaces.
- The user should make a visual check, when carrying out the measurement procedure, to ensure that all light spots are present.
- If the appliance repeatedly generates error messages, stop using it and contact the customer service.
- You are recommended always to examine both of the patient's eyes five times. The user should subject the measurement readings to extra scrutiny if there is a notable difference between the right and left eye.

- The following are classed as notable differences:
 - More than 1 D with respect to central corneal refractive power
=> 0.18 mm difference with respect to the corneal curvature radius
 - More than 0.3 mm with respect to axial eye length
 - More than 1 D with respect to emmetropic IOL refractive power
- The user must check the A-scan when measuring anterior chamber depth in pseudophakic mode. If only one IOL signal is visible, it is not clear whether this signal relates to the front or back of the IOL. Uncertainty in this case can lead to the displayed reading for anterior chamber depth being inaccurate by the thickness of the IOL (approx. ± 1 mm).
- The measured corneal thickness is not intended as a basis for the correction of intraocular pressure readings for diagnosing glaucoma.
- An excessively tilted or decentered IOL may make it impossible to measure the anterior chamber depth, lens thickness and aqueous depth.
- Measurement readings obtained from patients with a non-intact cornea (e.g., due to a corneal transplant, corneal opacity or corneal scarring, etc.) may possibly be inaccurate (this applies to keratometry in particular), and the user should check the data for plausibility.
- Ambient light has a bearing on pupil-diameter measurement readings. The user is responsible for ensuring the correct level of ambient light when carrying out pupillometry. The LS 900 cannot monitor ambient light, so do not use pupillometry as the decisive factor when considering keratorefractive surgery.
- Keratometry may not be accurate in eyes with keratoconus, and so the user must verify its plausibility.
- The user must verify that the eye assignment (OD, OS) is correct for the measured eye.
- The white-to-white distance reading is merely an indirect measurement of the inner lateral dimensions of the anterior ocular section. It therefore provides only approximate indications of the actual inner lateral dimensions of the anterior ocular section and of the size of the implant used.
- The measured results of patients with asteroid hyalosis may not be accurate (in particular the axial length measurement) and must therefore be verified by the user in terms of their plausibility.
- If the optional T-Cone is used, the user should verify the even illumination of the cone.

1.6.2 IOL calculation

The measurements taken with the LS 900 are a central element of every IOL calculation. A further important parameter in calculating the lens to be implanted is the IOL constant. When using the Lenstar LS 900, only IOL constants optimised for optical biometers should be used. Please contact your IOL manufacturer for information on optimised IOL constants for optical biometry. An alternative source of information for IOL constants optimised for optical biometry is the website of the "User Group for Laser Interference Biometry" (ULIB) at the University of Wuerzburg, Germany. Even though the constants published there have been optimised for a different optical biometer, published data [1, 2, 3] indicate that these IOL constants can also be used for calculating IOL power with the Lenstar.

To further improve the clinical result, HAAG-STREIT recommends that every surgeon creates personally optimised IOL constants based on pre-operative measurement data generated with the Lenstar and reliable postoperative data (e.g., 3 months after the operation).

1.6.3 References

- [1] Buckhurst P J, Wolffsohn J S, Shah S, Naroo S A, Davies L N, Berrow E J, "A new optical low coherence reflectometry device for ocular biometry in cataract patients", British Journal of Ophthalmology 2009;93:949-953
- [2] Holzer M P, Mamusa M, Auffarth G U, „Accuracy of a new partial coherence interferometry analyser for biometric measurements”, British Journal of Ophthalmology 2009;93: 807-810
- [3] Rohrer K, Frueh B E, Wälti R, Clemetson I A, Tappeiner C, Goldblum D, "Comparison and Evaluation of Ocular Biometry Using a New Noncontact Optical Low-Coherence Reflectometer", Ophthalmology 2009, 116:2087-2092

1.6.4 IOL constants

HAAG-STREIT recommends using personalised IOL constants to achieve maximum prediction accuracy of the IOL calculation. Using personalised IOL constants minimises the effects of individual surgical techniques, individual measurement and surgical equipment and individual physiological differences in the patient cohort under treatment on the IOL calculation.

1.6.5 IOL constants derived from data from an optical biometer

If there are optimised or personalised IOL constants available from a different optical biometer, these may also be used with the Lenstar [1, 2, 3]. Optimised constants for optical biometers can also be found on the website of the "User Group for Laser Interference Biometry" (ULIB) at the University of Würzburg, Germany.

1.6.6 IOL constants derived using data from an immersion ultrasound biometer

If there are optimised IOL constants available from an immersion ultrasound biometer, they may be used as a starting point for further optimisation with IOL calculation with the Lenstar. Differences in IOL calculation may still occur as the keratometry data are collected in different ways. The "User Group for Laser Interference Biometry" (ULIB) at the University of Würzburg, Germany has published on its website a manual explaining how to correct the effect of keratometry on the IOL constants. Still, the IOL constants acquired in this way should only be used as a starting point for further optimisation/personalisation.

1.6.7 IOL constants derived from data from a contact ultrasound biometer

If there are optimised IOL constants available based on data from a contact ultrasound biometer, they must be stringently recalculated for use with the Lenstar. The "User Group for Laser Interference Biometry" (ULIB) at the University of Würzburg, Germany has published on its website a manual explaining how to correct the effect of ultrasound biometry and keratometry on the IOL constants. Constants acquired in such a way should only be used as a starting point for further optimisation/personalisation.

1.7 Optical radiation

In accordance with EN 60825-1, the limiting values for class 1 lasers are respected when the device is used in the defined manner.

WARNING!

The light from this device may be dangerous. The risk of eye damage increases with the irradiation period. An exposure time with this device at maximum intensity of over 100 individual measurements per patient eye with dilated pupil per day exceeds the risk guideline value.

NOTE!

The device complies with the limit values for risk group 1 in accordance with EN 62471 insofar as no more than 100 individual measurements are performed per day and per patient's eye with dilated pupil. If this limit value of 100 individual measurements is exceeded, the patient's eye may be damaged by the white illumination. (Radiance 1,800 Wm⁻²sr⁻¹)

1.8 Disinfection

NOTE!

- The device does not need to be disinfected.
- If the optional T-Cone is used, its tip must be cleaned with max. 70% alcohol after each patient.

For more information on cleaning, please refer to the «Maintenance» section.

1.9 Warranty and product liability

- The product should be treated as described in the "Safety" chapter.
- Improper handling can damage the product. This would void all guarantee claims.
- Continued use of a product damaged by incorrect handling may lead to personal injury. In this case, the manufacturer accepts no liability.
- Repairs and modifications to the product must only be performed by HAAG-STREIT service technicians or authorised personnel.
- The optional T-Cone must be stored in its original dust cover and protected from direct sunlight.

1.10 Symbols

Read the instructions for use carefully

General warning: Read the accompanying documentation

Disposal instructions
See "Disposal" section.

CSA certification mark
with approval for USA

European certificate
of conformity

Manufacturer

Year of production

Direct current

 Alternating current

 HS reference number

IP20 Housing protection

 Logo of manufacturer:
HAAG-STREIT AG

 Serial number

2. Introduction

2.1 Basic construction

The system is divided into two parts: one concerned with examination (LS 900) and the other with control (Notebook, PC). The examination part communicates via a USB connection with the external PC. The LS 900 is operated using the "EyeSuite" software installed on the PC. Integral, automatic error recognition for measurements guarantees reliable examination results.


2.2 Examination components (LS 900)

Overview

1. Front ring
2. Housing
3. Service cover
4. Cable cover
5. Control lever
6. Type plate
7. Side identification sticker


Head rest (option)

8. Head rest
9. Headband
10. Mark for optimum eye height
11. Chin rest
12. Chin rest height adjustment
13. Hand grips for patient

T-Cone (option)


14. Dust cover of storage packaging
15. T-Cone topography add-on
16. "Top" – shows which side of the T-Cone is up.
17. Locking/release clip
18. Type plate
19. Base plate of storage packaging
20. Short forehead band with four Phillips screws

Device state

The device state indicator allows device monitoring without PC software:


21. State indicator	Dark	Off
	Orange	Standby
	Green	On
	Blue	Light source on
	RED	ERROR

Control lever


The joystick is used to position the device in relation to the patient's eye.

22. Trigger

Connections

23. USB-device connection
24. DC-device connection
25. USB 2.0 cable
26. DC cable
27. Cable tension compensator

2.3 Control component (PC)

A commercial PC is used as the control component for the biometer.


WARNING!

The software must be installed by trained personnel in accordance with the separate installation instructions. For further information, please contact your HAAG-STREIT representative.

2.4 Instrument table (option)

An adjustable instrument table (option) allows the height of the device to be set comfortably to the height of the individual patient.

- 28. Table top
- 29. Left-hand drawer (for external medical power supply) / Switch box SB01
- 30. Right-hand drawer (empty)
- 31. Elevator column (mechanical with spring)
- 32. Stand base with castors

3. Appliance assembly / installation

3.1 Safe system configuration according to EN 60601-1

FORBIDDEN!

- It is not permitted to make modifications to the measuring device.

WARNING!

- The device must only be installed on a stable, horizontal surface.
- All externally connected devices must comply with the standards relevant to safety.

NOTE!


- After installation, a zero adjustment and function check must be performed in accordance with the section "Function check".


NOTE!

- The device should be set up in such a way that the plug is always easily accessible and the device can be disconnected from the power supply easily.

- 33. Power supply (LS 900)
- 34. Instrument table
- 35. Isolation Transformer
- 36. Laptop / PC
- 37. Power supply (PC)
- 38. Printer

Galvanic
isolation 4 kV

3.2 Instrument table (option) assembly


The instrument table is delivered in a separate package.

To assemble the table, please consult the instructions for use enclosed with it.

3.2.1 Mounting the roller rails/slide plate


The roller rails/slide plate must be positioned in such a way that, when the device is centered, the patient's eye will be at a distance of 68 mm from the front of the device.

- 39. Rail
- 40. Side identification sticker
- 41. Slide plate
- 42. Headrest columns
- 43. Aperture for cables

3.2.2 Mounting the sticker for automatic position recognition

- Take off rail cover and move the device to the side, if already in position.
- Clean surface of table.
- Remove protective film (44) from back of the sticker. Start lifting it off carefully from the corner opposite the white area.
- Position the sticker with the white side against the left roller rail (44) and the edge of the table (46). Press the black/white area down firmly and stroke out any air bubbles.
- Carefully tear off the rest of the sticker (47) (Positioning Tool) along the perforations
- Check whether the position recognition function in the software works correctly.

3.3 Transporting the appliance


Transport the appliance over longer distances in its original packaging.

WARNING!

In its original packaging, the device can withstand the transport conditions as per EN ISO 15004-1.


For short distances, grasp the device by the base with both hands and lift (48), then support the cross slides with your second hand (49).

3.3.1 Device connection, cables under the device


The device can be invisibly wired by passing cables through the cross slide and using the cable channel underneath. To do this, you need either the HAAG-STREIT instrument table HSM 901 (option) or an aperture in the table-top used.

- Remove cable cover (50).
- Expose USB and supply cables (table aperture, (51)).
- Insert the supply cable through the opening beneath the cross slide until the cable emerges from the opening at the top at the front (52).
- Insert the USB cable through the same opening (52) on the cross slide.
- Insert both cables in the appropriate box and fix with the tension compensator (53).
- Replace the cable cover (54).
- Replace the rail cover.

3.3.2 Device connection, cables at side (option)

- Remove cable cover (55).
- Expose USB and supply cables (left or right aperture on rail cover).
- With a knife, cut carefully into the snap-off point provided on the left or right of the cable cover (56).
- Carefully use fingers to snap off the cut part (57).
- Insert both cables in the appropriate box and fix with the tension compensator.
- Join both cables together with a spiral wrap and fix in place using the cable channel on the rail cover (58).
- Replace cable cover
- Replace the rail cover.

3.4 Computer connection

WARNING!

For connection to PC, only use the supplied USB cable (2 m)

- Connect the electric power supply cable. Integral mains components work with the voltages specified under section A.1.1 "Electrical data". It is not necessary to select the voltage on the appliance.
- If an instrument table HSM 901 (option) has been supplied, the power supply of the LS 900 can be connected to the Switchbox SB01 (left drawer). Use the instructions for use enclosed with the switchbox and the instrument table.


3.4 Installing a short forehead rest when using the optional T-Cone

NOTE!

If the Lenstar LS 900 biometer with the T-Cone topography add-on is used on an instrument table with a HAAG-STREIT head rest (*HS art. no. 7200123*), the short forehead band delivered with the T-Cone must be used on the head rest for better measurability. If the Lenstar is operated without the T-Cone, the short forehead band does not need to be replaced with the long one.


- Detach the forehead band by removing the four Phillips screws (60) using a size 1 Phillips screwdriver.
- Remove the forehead band (59).
- Insert the short forehead band (*HS art. no. 1021653*) into the head rest and align the holes with the holes in the head rest.
- Attach the short forehead band using the four Phillips screws supplied.

4. Operation

4.1 Position of patient during measurement

The positioning of the device is done manually by the user. The patient must be positioned in such a way that the distance from the measuring head to the eye is approx. 68 mm. A steady head position is promoted by resting the patient's head well against the chin rest and forehead band, and by the patient holding on to the handles supplied. This can positively influence adjustment time and measuring accuracy. The patient should sit up as straight as possible.


**NOTE!**

In order to get the best possible results, the patient should be requested to keep the eye as wide open as possible during the measurement and to focus on the measuring beam. Blinking is permitted, but should be kept to a minimum.


4.2 Measuring with the optional T-Cone

4.2.1 Fitting the T-Cone


- Remove the dust cover from the storage packaging.
- Hold the T-Cone with the locking clip (17) pressed down next to the metal ring and fit it on the front ring of the Lenstar as shown on the right.
- Keep the locking clip depressed so that the T-Cone lies flush with the front ring. Make sure that "TOP" on the metal ring of the T-Cone is at the top. The T-Cone is held in place on the front ring of the Lenstar LS 900 biometer by strong magnets.

- Ensure that the red safety mark is no longer visible on the locking clip. If the red safety mark is still visible, remove the T-Cone and fit it again.

The T-Cone is fitted correctly. The red safety mark is not visible.

The T-Cone is not fitted correctly. The red safety mark is visible. Please remove the T-Cone and fit it again.

**WARNING!**

The T-Cone can only be used with the Lenstar with which it was calibrated. To use the T-Cone on another Lenstar or to use another T-Cone on a Lenstar with an already saved T-Cone calibration, the calibration must be performed again.

**WARNING!**

The measuring distance from the T-Cone to the eye (apex) is approx. 6 mm. Depending on the anatomy of the patient, the tip of the T-Cone may touch the eyelid or bridge of the nose. To avoid injuries, care must always be taken when moving the Lenstar with fitted T-Cone towards the patient.

**NOTE!**

When you use the T-Cone for the first time with the Lenstar LS 900 biometer, the software will prompt you to calibrate the T-Cone with the Lenstar LS 900 biometer. To do so, follow the software wizard and read the software operating instructions (F1 key).

**NOTE!**

We recommend performing a test measurement every time the T-Cone is fitted and removed. The test measurement can be invoked in the software's biometry menu. To do so, follow the software wizard and read the software operating instructions (F1 key).

4.2.2 Measuring with the optional T-Cone

- Fit the T-Cone as described in 4.2.1.
- Prior to every measurement procedure, the tip of the T-Cone should be cleaned with a lint-free cloth soaked in 70% alcohol. To do so, use a moistened cotton bud or lint-free cloth. The cloth or cotton bud must not be so wet that it drips. Ensure that the T-Cone is dry after the cleaning. Information on cleaning the T-Cone can be found under "Cleaning".
- Retract the Lenstar / T-Cone completely before positioning the patient in the head rest. Always start the measuring procedure with the Lenstar / T-Cone in the position furthest away from the patient.
- Explain to the patient that he should focus on the red, flashing light (measuring beam) in the centre of the T-Cone. The second eye (the one not being measured) can be covered with the eye patch on the optional head rest.
- Start the measuring procedure by pressing the button on the joystick and following the instructions on your PC screen. Detailed information on the measuring procedure can be found in the operating instructions for the software (F1 key).

5. Commissioning

Energy to the Lenstar LS 900 biometer is supplied by the medical power supply FRIWO 3288-DT12/12 HED (*HS no:1020392*).

5.1 Switching on the appliance

- Connect the power supply up to the mains.
- Switch on the PC.
- Start the software on the PC.

5.2 Switching off the appliance

- Exit the software on the PC.
- Switch off the PC.
- Disconnect the Lenstar from the mains if you do not intend to use it for an extended period of time.

6. Technical data

Type designation	LS 900
Dimensions (W x D x H):	310 x 260 x 420 mm
Weight:	6.2 kg

6.1 Power supply (primary side)

Power supply:	100 – 240 V / 50 – 60 Hz 280 – 140 mA
	<i>HS no:1020392</i> , FRIWO 3288-DT12/12 HED

6.1.1 Device (secondary side)

Voltage:	12 V ===
Current:	1 A

6.2 Eye length measurement (A-Scan) and central fixation

Light source:	Superluminescent diode
Wavelength:	820 nm
Power on patient's eye:	< 0.6 mW

6.2.1 Keratometry

Light source:	LED
Wavelength:	950 nm

6.2.2 Illumination

	Serial number of the system	
	Up to 1999*	As of 2000
Light source:	LED	LED
Light colour:	Green	White

* Devices with a serial number < 1999 can be subsequently converted to white light illumination

6.2.3 Positioning aid (as of serial number 2000)

Light source:	LED
Wavelength:	940 nm

6.2.4 Technical data of T-Cone (option)

Type designation:	T-Cone
Diameter:	ø 84 mm
Length:	63.5 mm
Weight:	0.2 kg
Placido rings:	11
Covered optical zone:	≤ 6 mm

6.3 Reproducibility

- The T-Cone complies with the requirements of the ISO 19980:2012 standard for the central and medial measuring zones; the T-Cone does not cover the peripheral measuring zone.
- International standard for ophthalmological instruments.
- Corneal topograph, type B.
- The T-Cone complies with the requirements of the ANSI Z80.23-2012 standard for the central and medial measuring zones; the T-Cone does not cover the peripheral measuring zone.
- US national standard for ophthalmological instruments
Corneal topography system, type B.

6.4 Measuring accuracy

Exact test area (axial curvature difference in mm / elevation difference in μm)

	Central zone (Diameter ≤ 3 mm)			Medial zone (Diameter ≤ 6 mm)		
	Average value Curvature diff.	2 st. dev. Curvature diff.	2 st. dev. Elevation diff.	Average value Curvature diff.	2 st. dev. Curvature diff.	2 st. dev. Elevation diff.
Torus (mm), $R_1=7.987$, $R_2=7.584$	0.006	0.065	0.66	0.004	0.044	0.33
Ellipsoid (mm), $R=7.79$, $k=-0.255$	0.001	0.083	0.78	0.020	0.038	1.51
Sphere 1 (mm), $R=6.448$,	0.013	0.037	0.61	0.012	0.025	1.12
Sphere 2 (mm), $R=7.804$,	-0.008	0.034	0.33	-0.008	0.026	0.75
Sphere 3 (mm), $R=8.844$,	-0.017	0.056	0.46	-0.024	0.038	1.61
Sphere 4 (mm), $R=10.501$,	-0.040	0.082	0.48	-0.069	0.044	0.58

Sphere 1: $\varnothing 11.5$ mm

Torus, asphericity, sphere 2,3,4: $\varnothing 14$ mm

All surfaces: Precision $\pm < 1$ μm

Test surfaces were centred within ± 0.1 mm and their symmetry axes within $\pm 0.5^\circ$ with regard to their measuring axis.

6.4.1 Reproducibility table B

Human cornea (axial curvature difference, D / elevation difference, μm)

Zone 1: Centre (Diameter ≤ 3 mm)			Zone 2: Average value (Diameter ≤ 6 mm)		
Average value	1 st. dev.	2 st. dev.	Average value	1 st. dev.	2 st. dev.
-0.001 / -	0.38 / 0.66	0.76 / 1.31	-0.020 / -	0.36 / 3.22	0.72 / 6.45

n = 42 eyes

6.5 Measured variables

Measurement ranges are based on the 'Phakic' measurement mode.

6.5.1 Corneal thickness (CT)

Measurement range:	300 – 800 μm
Display resolution:	1 μm
In vivo repeatability (1.SD):	± 2.3 μm

6.5.2 Anterior chamber depth (ACD)

Measurement range:	1.5 – 6.5 mm
Display resolution:	0.01 mm
In vivo repeatability (1.SD):	± 0.04 mm

6.5.3 Lens thickness (LT)

Measurement range:	0.5 – 6.5 mm
Display resolution:	0.01 mm
In vivo repeatability (1.SD):	± 0.08 mm

6.5.4 Axial length (AL)

Measurement range:	14 – 32 mm
Display resolution:	0.01 mm
In vivo repeatability (1.SD):	± 0.035 mm

6.5.5 Keratometry (R)

Measurement range radius:	5 – 10.5 mm
Display resolution:	0.01 mm
In vivo repeatability (1.SD):	± 0.03 mm
Measuring area of axial angle:	$0-180^\circ$
Display resolution:	1°
In vivo repeatability (1.SD):	$\pm 11^\circ$

6.5.6 White-to-white distance (WTW)

Measurement range:	7 – 16 mm
Display resolution:	0.01 mm
In vivo repeatability (1.SD):	± 0.04 mm

6.5.7 Pupillometry

Measurement range:	2 – 13 mm
Display resolution:	0.01 mm

The above mentioned measurement ranges correspond to the default setting for the automatic analysis. The in vivo reproducibility was evaluated in a clinical study of cataract patients (see Tables 1 "All eyes" and 2 "Special eyes")

6.5.8 Study design

- The clinical trial for deriving in vivo reproducibility was approved by the local ethics committee. The study was planned and conducted as a prospective, non-randomised comparative study.
- Two phases were scheduled in the approved study protocol. In the first phase, measurements of axial length (AL), central corneal thickness (CT), anterior chamber depth (ACD), central lens thickness (LT), median corneal radius (R) and the axis position of the flat meridian (axis) were taken.
- The white-to-white distance (WTW) was measured in the second phase.
- A total of 144 eyes in 80 subjects were included in study stage 1 and a total of 40 eyes in 20 subjects in study stage 2.
- Subjects with different conditions of the anterior and posterior segment of the eye (cataract in different stages, pseudophakia with different IOLs, aphakia, silicone-oil fill) and subjects with healthy eyes were included in this clinical study. Data were analyzed for all eyes (see Tab. 1) and for a subgroup of eyes with special health conditions (see Tab. 2), with a complete set of 5 repeat measurements on both eyes of each subject.
- The special eyes group included eyes with one or more of the following conditions: pseudophakia, aphakia and silicone-oil fill.

Tab. 1: All Eyes

[unit]	n	Mean _{grand}	SD _{repeat}	CV
AL [mm]	45 / 90	23.973	0.035	0.00145
CT [µm]	53 / 106	557.1	2.3	0.00407
ACD [mm]	34 / 68	3.19	0.04	0.01220
LT [mm]	27 / 54	4.56	0.08	0.01784
R [mm]	34 / 68	7.67	0.03	0.00396
Axis [°]	27 / 54	72	11	0.14191
WTW [mm]	9 / 18	12.27	0.04	0.00337

Tab. 2: Special Eyes

[unit]	n	Mean _{grand}	SD _{repeat}	CV
AL [mm]	10 / 20	24.087	0.056	0.00234
CT [µm]	11 / 22	564.4	2.8	0.00496
ACD [mm]	5 / 10	7.75	0.03	0.00333
Axis [°]	3 / 6	80	13	0.16092

Abbreviations

n	Number of subjects / number of eyes		
Mean_{grand}	Overall mean of results on all eyes		
SD_{repeat}	Repeatability standard deviation		
CV	Coefficient of variation		
SD	Standard deviation	LT	Lens thickness
AL	Axial length	R	Corneal radius of curvature
CT	Corneal thickness	Axis	Axis of flat meridian
ACD	Anterior chamber depth	WTW	White-to-white distance

7. Software / Help menu / Error messages

The software's help section contains instructions and help for performing an examination and descriptions of the error messages. The help can be opened via the F1 key or in the [?] – [Help] menu.

WARNING!


The software must be installed by trained personnel in accordance with separate installation instructions.

8. Maintenance

The LS 900 is practically maintenance-free and requires only minimum care to work for as long as possible to your complete satisfaction. However, we recommend instructing a service engineer to inspect the biometer periodically. HAAG-STREIT or your local agent will be happy to provide further information.

8.1 Function check / zero adjustment

The inspection gauge (65) supplied with the device will assist the user in checking the proper functioning of the unit. The type plate with the serial number can be found on the rear of the gauge. This check takes place for the first time when the device is commissioned. Subsequent test intervals (1 week) will be specified by the software (message displayed). For exact procedure see software instructions for use. Should the test be outside the tolerance limits, a software message appears and the unit must be taken out of use. Notify your HAAG-STREIT customer support center.

8.2 Cleaning

Regular dusting of the device with a soft cloth is sufficient. More stubborn dirt can be removed using a soft, lint-free cloth dampened with water or alcohol at maximum 70%.

WARNING!

Avoid making the device wet and use only the means listed above. Under no circumstances use solvent or any abrasives.

WARNING!

The T-Cone must never be immersed in cleaning fluid or disinfectant. No abrasive or otherwise aggressive agents must be used for cleaning/disinfection. Suitable media include water and up to 70% alcohol. Ensure that the T-Cone is dry after every cleaning/disinfection.

A dust cover is included in the accessories of the LS 900. Cover the appliance when the room is being cleaned or if it is not used for longer periods of time. Always remove the dust cover before switching on the power.

WARNING!

The appliance must not be switched on when covered! (Heat build-up, fire hazard).

8.3 Maintenance of T-Cone (option)

The T-Cone requires practically no maintenance apart from the cleaning of the tip between patients (see Measuring with the T-Cone).

A. Appendix

A.1 Accessories / spare parts

A1.1 LENSTAR LS 900

Component	HS art. no.	Note
Electric power supply lead CH	1001319	Max. length 2.50 m
USA power supply lead	1001316	Max. length 3.10 m
Dust cover (small)	1001395	
Checking gauge (replacement)	1021124	
Head rest (LS 900)	7200123	
Eye patch	1400113	
Instrument table HSM 901 Workstation, manual	7220085	
Instrument table HSM 901 Workstation, electric lifting column 230 V	7220149	

A1.2 T-Cone (option)

Component	HS art. no.	Note
T-Cone topography add-on	1021534	
Dust cover of storage packaging	1021665	
Base plate of storage packaging	1021666	
Short forehead band with Phillips screws for forehead band	1021653 1005072	4 pcs

B. Legal regulations

- HAAG-STREIT maintains a quality management system in accordance with EN ISO 13485. The device has been developed and designed taking into consideration all the standards listed in section E 'Standards'.
- This is a Class IIa device in accordance with Appendix IX of Directive 93/42/EEC. By affixing the CE mark we confirm that our device complies with the applicable standards and directives.
- You can request a copy of the declaration of conformity for the appliance from HAAG-STREIT at any time.

C. Classification

C.1 LENSTAR LS 900 biometer

EN 60601-1:	Continuous operation
EN 60825-1:	Laser class I
EN 60529:	Enclosure protection IP20
EN 62471:	Risk group 1
EN ISO 15004-2	Group 2
93/42EEC (medical devices)	Class IIa
FDA:	Class II

C.2 T-Cone (option)

ISO 19980:2005	Corneal topograph, type B
ANSI Z80.23-2008	Corneal topograph, type B
93/42EEC (medical devices)	Class I
FDA:	Class I

D. Disposal

- Electrical and electronic devices must be disposed of separately from household waste! This appliance was made available for sale after the 13th August 2005.
- For correct disposal, please contact your HAAG-STREIT representative. This will guarantee that no hazardous substances enter the environment and that valuable raw materials are recycled.

E. Standards

EN 60601-1	EN ISO 15004-1
EN 60601-1-2	EN ISO 15004-2
EN 62471:2008	EN ISO 10343
EN 60825-1:2007	EN ISO 22665:2013
EN 60529	ISO 19980:2005
ANSI Z80.23:2008	

F. EMC supplement

F.1 General

The Lenstar LS 900 complies with the standard EN 60601-1-2. The instrument is constructed such that the generation and emission of electromagnetic interference is limited so as not to disrupt the intended use of other devices, while the instrument itself possesses an appropriate immunity to electromagnetic interference.

WARNING!

- Electrical medical devices and systems are subject to special EMC measures and must be installed in accordance with the EMC instructions contained in this accompanying document.
- Portable and mobile HF communication systems may interfere with electrical medical devices.

WARNING!

The use of cables or equipment other than those listed may lead to a higher emission or to reduced interference immunity of the Lenstar LS 900 biometer system.

F.2 Emitted interference (standard table 1)

Guidance and manufacturer's declaration – electromagnetic emissions

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Emission of harmonics according to EN 61000-3-2	Class A	

F.3 Immunity (standard table 2)

Guidance and manufacturer's declaration – electromagnetic immunity

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Immunity test standard	EN 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst EN 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	± 1 kV for symmetrical voltages ± 2 kV for asymmetrical voltages	± 1 kV for symmetrical voltages ± 2 kV for asymmetrical voltages	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines EN 61000-4-11	< 5% U_T (> 95% drop in U_T) for ½ cycle < 40% U_T (> 60% drop in U_T) for 5 cycles < 70% U_T (> 30% drop in U_T) for 25 cycles < 5% U_T (> 95% drop in U_T) for 5 s	< 5% U_T (> 95% drop in U_T) for ½ cycle < 40% U_T (> 60% drop in U_T) for 5 cycles < 70% U_T (> 30% drop in U_T) for 25 cycles < 5% U_T (> 95% drop in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this product requires continued function even in the event of interruptions in the energy supply, this product should be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field EN 61000-4-8	3 A/m	0.3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T = the AC mains voltage prior to application of the test level.

F.4 Immunity on non-life support devices (standard table 4)

Guidance and manufacturer's declaration – electromagnetic immunity

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Electromagnetic environment – guidance

Portable and mobile RF communications equipments should be used no closer to any part of this product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Immunity test standard	EN 60601 test level	Compliance level	Recommended distance:
Conducted RF EN 61000-4-6	3 V _{rms} 150 kHz – 80 MHz	3 V _{rms}	$D = 1.2 \sqrt{P}$
Radiated RF EN 61000-4-3	3 V/m 80 MHz – 2.5 GHz	3 V/m	$D = 1.2 \sqrt{P}$ 80 MHz – 800 MHz
		80 MHz – 800 MHz	$D = 2.3 \sqrt{P}$ 800 MHz – 2.5 GHz

Where **P** is the maximum output power rating of the transmitter in watts (**W**) according to the transmitter manufacturer and **D** is the recommended separation distance in metres (**m**). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: At 80 MHz and 800 MHz the higher frequency applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this product is used exceeds the applicable RF compliance level above, this product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this product.
- Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

F.4 Safe distances on non-life support devices (standard table 6)

Recommended safe distances between portable and mobile HF communication devices and this device.

This product is designed to be operated in an electromagnetic environment in which radiated HF interference is controlled. The customer or user of this product can help to prevent electromagnetic interference by maintaining minimum distances between portable and mobile HF communication systems (transmitters) and this product, as recommended below in accordance with the maximum output of the communication system.

Safe distance according to transmission frequency (m)

Nominal output of the transmitter (W)	150 kHz – 80 MHz	80 MHz – 800 MHz	800 MHz – 2.5 MHz
	$D = 1.2 \sqrt{P}$	$D = 1.2 \sqrt{P}$	$D = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters with a nominal output not listed in the table above, the distance D can be calculated in meters (**m**) using the equation for the respective column, in which P is the nominal output of the transmitter in watts (**W**) according to the specifications of the transmitter manufacturer.

NOTE 1: To calculate the recommended safe distance of transmitters in the frequency range of 80 MHz to 2.5 GHz an additional factor of $10/3$ was used to reduce the probability of a mobile/portable communication device causing interference if inadvertently brought into the patient area.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic wave propagation is influenced by absorption and reflection of buildings, objects and people.

Should you have any further questions, please contact your HAAG-STREIT dealer at:


<http://www.haag-streit.com/contact/contact-your-distributor.html>

HS HAAG-STREIT
DIAGNOSTICS

PRODUCTS CERTIFIED FOR BOTH THE U.S.
AND CANADIAN MARKETS, TO THE APPLI-
CABLE U.S. AND CANADIAN STANDARDS

Certified Quality Management System

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