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REF 07P7120

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BIBLIOGRAFIE

- Gotto AM. Lipoprotein metabolism and the etiology of hyperlipidemia. *Hosp Pract* 1988;23(Suppl 1):4.
- Crouse JR, Parks JS, Schey HM, et al. Studies of low density lipoprotein molecular weight in human beings with coronary artery disease. *J Lipid Res* 1985;26(5):566.
- Badimon JJ, Badimon L, Fuster V. Regression of atherosclerotic lesions by high-density lipoprotein plasma fraction in the cholesterol-fed rabbit. *J Clin Invest* 1990;85:1234-1241.
- Castelli WP, Doyle JT, Gordon T, et al. HDL cholesterol and other lipids in coronary heart disease. *Circulation* 1977;55(5):767.
- Barr DP, Russ EM, Eder HA. Protein-lipid relationships in human plasma. *Am J Med* 1951;11(4):480.
- Gordon T, Castelli WP, Hjortland MC, et al. High density lipoprotein as a protective factor against coronary heart disease. *Am J Med* 1977;62(5):707.
- Williams P, Robinson D, Baily A. High density lipoprotein and coronary risk factor. *Lancet* 1979;1:72.
- Kannel WB, Castelli WP, Gordon T. Cholesterol in the prediction of atherosclerotic disease: new perspectives based on the Framingham study. *Ann Intern Med* 1979;90:85.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Tietz NW, editor. *Clinical Guide to Laboratory Tests*. Philadelphia, PA: WB Saunders; 1986:256.

- National Cholesterol Education Program. *Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III)*. National Institutes of Health. National Heart, Lung, and Blood Institute. NIH Publication No. 01-3670. May 2001.
- Calam RR, Bessman JD, Ernst DJ, et al. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. (H18-A3). Wayne, PA: Clinical and Laboratory Standards Institute, 2004.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 3rd ed. Washington, DC: AACC Press; 1990:3-104-3-106.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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BIBLIOGRAFIE

- Gotto AM. Lipoprotein metabolism and the etiology of hyperlipidemia. *Hosp Pract* 1988;23(Suppl 1):4-13.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III)—Final Report. National Institutes of Health. National Heart, Lung, and Blood Institute. NIH Publication No. 02-5215. September 2002; I-1-II-22.
- Badimon JJ, Badimon L, Fuster V. Regression of atherosclerotic lesions by high density lipoprotein plasma fraction in the cholesterol-fed rabbit. *J Clin Invest* 1990;85(4):1234-41
- Castelli WP, Doyle JT, Gordon T, et al. HDL Cholesterol and other lipids in coronary heart disease. The cooperative lipoprotein phenotyping study. *Circulation* 1977;55(5):767-72.
- Gordon T, Castelli WP, Hjortland MC, et al. High density lipoprotein as a protective factor against coronary heart disease. *Am J Med* 1977;62(5):707.
- Williams P, Robinson D, Bailey A. High density lipoprotein and coronary risk factors in normal men. *Lancet* 1979;1(8107):72-5.
- Kannel WB, Castelli WP, Gordon T. Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study. *Ann Intern Med* 1979;90(1):85-91.
- Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on detection, evaluation, and treatment of high blood cholesterol in adults (Adult Treatment Panel III). *JAMA* 2001;285:2486-97.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- National Cholesterol Education Program. *Recommendations on Lipoprotein Measurement, from the Working Group on Lipoprotein Measurement*. National Institutes of Health. National Heart, Lung, and Blood Institute. NIH Publication No. 95-3044. September 1995.

14. Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:22–23.
15. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
16. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
17. Camps J, Simó JM, Guaita S, et al. Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogenous Methods for HDL-Cholesterol. *Clin Chem* 1999;45(5):685–88.
18. Roberts WL, Leary ET, Lambert TL, et al. Falsely low direct HDL-cholesterol results in a patient with dysbetalipoproteinemia. *Clin Chem* 2000;46:560–2.
19. Lackner KJ, Schmitz G. Beta-VLDL of patients with type III hyperlipoproteinemia interferes with homogenous determination of HDL-cholesterol based on polyethylene glycol-modified enzymes. *Clin Chem* 1998;44:2546–8.
20. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
21. Cholesterol Reference Method Laboratory Network. *HDL Cholesterol Certification Protocol for Manufacturers*. November 2002. Accessed July 11, 2005 from: <http://www.cdc.gov/labstandards/pdf/crmln/MFRHDLNov2002final.pdf>.
22. Krouwer JS, Astles JR, Cooper WG, et al. *Estimation of Total Analytical Error for Clinical Laboratory Methods; Approved Guideline (EP21-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 2003.
23. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
24. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
25. Young DS, *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000:3-399–3-414.
26. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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BIBLIOGRAFIE

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Randox Laboratories Ltd
Ardmore, 55 Diamond Road,
Crumlin, Co Antrim, BT29 4QY, UK

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FOR USE WITH

Alinity c

Read Highlighted Changes: Revised March 2022.

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Iron2

INTENDED USE

The Iron2 assay is used for the direct colorimetric determination of iron without deproteinization in human serum or plasma on the Alinity c system.

The Iron2 assay is to be used as an aid in the diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis (a disease associated with widespread deposit in the tissues of two iron-containing pigments, hemosiderin and hemofuscin, and characterized by pigmentation of the skin), and chronic renal disease.

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Iron exists in biological fluids as a component of hemoglobin and myoglobin and is bound in serum and plasma to transferrin, which acts as a carrier protein. Increased iron concentrations are seen in hemolytic anemias, hemochromatosis, and acute liver disease. Decreased iron concentrations are seen in iron deficiency and anemia of chronic disease, such as in chronic renal disease.¹ Major causes of iron deficiency include gastrointestinal and menstrual bleeding. For the assessment of the body's iron status, the measurement of transferrin and ferritin can provide more accurate information.²

PRINCIPLES OF THE PROCEDURE

The Iron2 assay is an automated clinical chemistry assay. At an acidic pH, iron is released from transferrin to which it is bound, and then quantitatively reduced to a ferrous state. The iron forms with ferene-S (3-(2-pyridyl)-5,6-bis-[2-(5-furylsulfonic acid)]-1,2,4-triazine), a stable colored complex of which the color intensity is proportional to the amount of iron in the sample.

Particular reaction conditions and a specific masking agent almost entirely eliminate the interference from copper.

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Warnings and Precautions

- IVD**
- For *In Vitro* Diagnostic Use
- Rx ONLY**

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials and all consumables contaminated with potentially infectious materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate regional, national, and institutional biosafety practices should be used for materials that contain, are suspected of containing, or are contaminated with infectious agents.³⁻⁶

The following warnings and precautions apply to: R1	
WARNING	Contains guanidine hydrochloride, acetic acid, thiourea and methylisothiazolones.
H302	Harmful if swallowed.
H332	Harmful if inhaled.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H317	May cause an allergic skin reaction.
H351	Suspected of causing cancer.
H361	Suspected of damaging fertility or the unborn child.
H402*	Harmful to aquatic life.
H412	Harmful to aquatic life with long lasting effects.
Prevention	
P202	Do not handle until all safety precautions have been read and understood.
P261	Avoid breathing mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P271	Use only outdoors or in a well-ventilated area.
P273	Avoid release to the environment.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response	
P301+P330+P312	IF SWALLOWED: Rinse mouth. Call a POISON CENTER or doctor / physician if you feel unwell.
P302+P352	IF ON SKIN: Wash with plenty of water.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P308+P313	IF exposed or concerned: Get medical advice / attention.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P337+P313	If eye irritation persists: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

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The following warnings and precautions apply to: R2	
	
WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
H402*	Harmful to aquatic life.
H412	Harmful to aquatic life with long lasting effects.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

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Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.

- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

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Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position.
Onboard	System Temperature	30 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

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INSTRUMENT PROCEDURE

The Iron2 assay file must be installed on the Alinity c system prior to performing the assay.

The Alinity ci-series system software version 3.2.0 or higher must be installed on the Alinity c system prior to performing the assay.

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■ SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

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Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum
	Serum separator
Plasma	Lithium heparin
	Lithium heparin separator
	Sodium heparin

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Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - pooled specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
 - specimens with fungal growth
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

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NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

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- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

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Specimen Storage

Specimen Type	Temperature	Maximum Storage Time
Serum/Plasma	Room temperature (20 to 25°C)	10 hours ⁷
	2 to 8°C	7 days ⁸
	-20°C	12 months ⁹

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For additional information on sample handling and processing, refer to CLSI GP44-A4.¹⁰ The storage information provided here is based on references.

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Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

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- Iron2 assay file
- 04V6201 Consolidated Chemistry Calibrator
- Controls containing iron
- Saline (0.85% to 0.90% NaCl) for specimen dilution

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- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- Minimum sample cup volume is calculated by the system and printed on the Order List report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
 - Sample volume for single test: 20.0 µL.

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Sample Dilution Procedures

Sample dilutions have not been evaluated for the Iron2 assay. Samples with an iron value exceeding 1143 µg/dL (204.6 µmol/L) are flagged with the code "> 1143 µg/dL" (> 204.6 µmol/L). The standard dilution factor for the Iron2 assay is 1:1.45.

For details on configuring automated dilutions, refer to the Alinity ci-series Operations Manual, Section 2.

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Calibration is stable for approximately 15 days (360 hours) but is required with each change in reagent lot. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

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- At least 2 levels of controls (low and high) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

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Quality Control Guidance

Refer to "Basic QC Practices" by James O. Westgard, Ph.D. for guidance on laboratory quality control practices.¹¹

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RESULTS

Calculation

The Iron2 assay utilizes the Linear data reduction method to generate a calibration and results.

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Reportable Interval

Based on representative data for the limit of quantitation (LoQ) and the limit of detection (LoD), the ranges over which results can be reported are provided below according to the definitions from CLSI EP34, 1st ed.¹²

	µg/dL	µmol/L
Analytical Measuring Interval (AMI) ^a	7 - 1143	1.3 - 204.6
Reportable Interval ^b	4 - 1143	0.7 - 204.6

^a AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in µg/dL (µmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

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NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the reportable interval of 4 µg/dL (0.7 µmol/L). To flag values using the lower limit of the analytical measuring interval of 7 µg/dL (1.3 µmol/L), the operator must edit the Low Linearity value, adjusted by the standard dilution factor.

For detailed information on editing the result settings of assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- Falsely elevated iron results may be observed at the low end of the analytical measuring interval in samples with triglyceride concentrations above 200 mg/dL.
- Falsely elevated iron results may be observed at the low end of the analytical measuring interval in samples with unconjugated bilirubin concentrations above 25 mg/dL.
- Iron dextran treatment can result in elevated total iron results.
- Use of the Iron2 assay for patients undergoing treatment with deferoxamine or other iron chelating compounds is not recommended.
- Transiently elevated iron levels can be observed post ingestion of supplements/vitamins that contain iron.¹³
- Rifampicin levels above 5 mg/L may produce artificially low results with the Iron2 assay.
- Substances that demonstrated interference with the Iron2 assay are listed in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- Potential interference has not been evaluated for substances other than those described in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.

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Reference Range

Age	Range (µg/dL)	Range (µmol/L)
Pediatric¹⁴		
0 to < 14 years	16 - 128	2.9 - 22.9*
14 to < 19 years (Female)	20 - 162	3.6 - 29.0*
14 to < 19 years (Male)	31 - 168	5.5 - 30.1*
Adult¹⁵		
Female	50 - 170	9.0 - 30.4
Male	65 - 175	11.6 - 31.3

* Alternate result units were calculated by Abbott and are not included in the citation provided.

Abbott has not evaluated reference ranges in the pediatric population.

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Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A3.¹⁶ Testing was conducted using 3 lots of the Iron2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 1 instrument. Two controls and 3 human serum panels were tested in a minimum of 2 replicates twice per day on 20 days on 3 reagent lot/calibrator lot combinations, where a unique reagent lot and a unique calibrator lot are paired. The performance from a representative combination is shown in the following table.

Sample	n	Mean (µg/dL)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	245	0.7	0.3	1.4 (1.4 - 4.0)	0.6 (0.6 - 1.7)
Control Level 2	80	69	0.5	0.7	0.8 (0.7 - 1.2)	1.1 (1.0 - 1.8)
Panel A	80	14	0.7	4.6	0.8 (0.6 - 0.9)	5.6 (3.8 - 6.1)
Panel B	80	375	0.7	0.2	1.6 (1.6 - 4.2)	0.4 (0.4 - 1.2)
Panel C	80	975	2.4	0.2	2.8 (2.8 - 8.8)	0.3 (0.3 - 0.9)

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Reproducibility

A study was performed based on guidance from CLSI EP05-A3.¹⁶ Testing was conducted using 1 lot of the Iron2 reagents, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 3 human serum panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

Sample	n	Mean (µg/dL)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	245	0.7	0.3	0.9	0.4	1.0	0.4
Control Level 2	90	69	0.4	0.5	0.4	0.6	0.9	1.3
Panel A	90	14	0.4	2.8	0.4	2.9	0.6	4.3
Panel B	90	374	0.9	0.2	1.0	0.3	1.8	0.5
Panel C	90	970	2.4	0.2	2.7	0.3	7.2	0.7

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Accuracy

A study was performed to estimate the bias of the Iron2 assay relative to standard reference material NIST SRM 3126. Testing was conducted using 3 concentrations of the standard across 3 lots of the Iron2 reagents, 2 lots of the Consolidated Chemistry Calibrator, and 1 instrument. The bias ranged from -3.7% to 2.4% across concentrations of the standard, calibrator, and reagent lots.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.¹⁷ Testing was conducted using 3 lots of the Iron2 reagents on each of 2 instruments over a minimum of 3 days. The limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized below. These representative data support the lower limit of the analytical measuring interval.

	µg/dL	µmol/L
LoB ^a	1	0.2
LoD ^b	4	0.7
LoQ ^c	7	1.3

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^b The LoD presented in the table is in alignment with the Iron2 assay on the ARCHITECT c System. The observed LoD on the Alinity c system was 3 µg/dL (0.5 µmol/L) and represents the lowest concentration at which the analyte can be detected with 95% probability based on n ≥ 60 replicates of low-analyte level samples.

^c The LoQ presented in the table is in alignment with the LoQ for the Iron2 on the ARCHITECT c System. The observed LoQ on the Alinity c system was 5 µg/dL (0.9 µmol/L). This LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from n ≥ 60 replicates of low-analyte level samples.

Linearity

A study was performed based on guidance from CLSI EP06-A.18
This assay is linear across the analytical measuring interval of 7 to 1143 µg/dL (1.3 to 204.6 µmol/L).

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A study was performed based on guidance from CLSI EP07, 3rd ed.¹⁹ Each substance was tested at 2 levels of the analyte (approximately 50 µg/dL and 180 µg/dL).

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A study was performed based on guidance from CLSI EP07, 3rd ed.¹⁹ Each substance was tested at 2 levels of the analyte (approximately 50 µg/dL and 180 µg/dL).

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No Significant Interference (Interference within ± 10%)

Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Mebendazole	2 mg/L	6.78 µmol/L
Methyldopa	25 mg/L	118 µmol/L
Metronidazole	130 mg/L	759 µmol/L
Phenylbutazone	330 mg/L	1069 µmol/L
Rifampicin	5 mg/L	6.10 µmol/L
Sodium heparin	4 U/mL	N/A
Stanozolol	60 mg/L	183 µmol/L
Theophylline (1,3-dimethylxanthine)	60 mg/L	333 µmol/L

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Interferences from medication or endogenous substances may affect results.²⁰

Method Comparison

A study was performed based on guidance from CLSI EP09-A3,²¹ using the Passing-Bablok regression method.

Iron2 on the Alinity c system vs Iron on the ARCHITECT c System						
	n	Units	Correlation Coefficient	Intercept	Slope	Concentration Range
Serum	113	µg/dL (µmol/L)	1.00	-1 (-0.2)	1.05	16 - 879 (2.9 - 157.2)

Iron2 on the Alinity c system vs Iron2 on the ARCHITECT c System						
	n	Units	Correlation Coefficient	Intercept	Slope	Concentration Range
Serum	115	µg/dL (µmol/L)	1.00	1 (0.1)	1.00	16 - 933 (2.9 - 166.9)

BIBLIOGRAPHY

- Chen TK, Knicely DH, Grams ME. Chronic kidney disease diagnosis and management: a review. *JAMA* 2019;322(13):1294-1304
- Jacobs DS, DeMott WR, Oxley DK. *Laboratory Test Handbook*, 5th ed. Hudson, OH: Lexi-Comp; 2001:203-204.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 6th ed. Washington, DC: US Government Printing Office; June 2020.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Henriksen LO, Faber NR, Moller MF, et al. Stability of 35 biochemical and immunological routine tests after 10 hours storage and transport of human whole blood at 21°C. *Scand J Clin Lab Invest* 2014;74(7):603-610.
- Klat RL, Stein CDS, Santos MD, et al. Effects of freeze-thaw cycles and assessment of short-term storage stability on serum iron, ferritin, and transferrin. *Clin Lab* 2021;67(2).
- Chen H, Sternberg MR, Schleicher RL, et al. Long-term stability of 18 nutritional biomarkers stored at -20 °C and 5 °C for up to 12 months. *J App Lab Med* 2018;3(1):100-108.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking*. 1st ed. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Burtis CA, Bruns DE, editors. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, MO: Saunders Elsevier; 2015.
- Colantonio DA, Kyriakopoulou L, Chan MK, et al. Closing the gaps in pediatric laboratory reference intervals: a CALIPER database of 40 biochemical markers in a healthy and multiethnic population of children. *Clin Chem* 2012;58:854-868.
- Wu AHB, editor. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Elsevier Saunders; 2006.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
REF	List Number
SN	Serial number
Other Symbols	
DISTRIBUTED IN THE USA BY	Distributed in the USA by
FOR USE WITH	Identifies products to be used together
INFORMATION FOR USA ONLY	Information needed for United States of America only
PRODUCT OF IRELAND	Product of Ireland
R1	Reagent 1
R2	Reagent 2
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

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BIBLIOGRAFIE

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
5. Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 4th ed. St. Louis, MO: Elsevier Saunders; 2006:1912.
6. Taylor EC, Sethi B. Stability of 27 biochemistry analytes in storage at a range of temperatures after centrifugation. *Br J Biomed Sci* 2011;68:147-157.
7. Wu W, Yang D, Tiselius HG, et al. Collection and storage of urine specimens for measurement of urolithiasis risk factors. *Urology* 2015;85:299-303.
8. Remer T, Montenegro-Bethancourt G, Shi L. Long-term urine biobanking: storage stability of clinical chemical parameters under moderate freezing conditions without use of preservatives. *Clin Biochem* 2014;47:307-311.

9. Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
10. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
11. Wu AHB. *Tietz Clinical Guide to Laboratory Tests*, 4th ed. Philadelphia, PA: WB Saunders; 2006:706-708.
12. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
13. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
14. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
15. Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.
16. Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
17. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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BIBLIOGRAFIE

- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:807-810, 812-813.
- Jacobs DS, DeMott WR, Oxley DK, et al. *Laboratory Test Handbook*. Hudson, OH: Lexi-Comp; 2001:144-147, 258-260, 275-279.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Kaplan LA, Pesce AJ, Kazmierczak SC. *Clinical Chemistry – Theory, Analysis, and Correlation*, 3rd ed. St. Louis, MO: Mosby, 1996:461–463.
- McPherson RA, Bidkorpeh EK, Castellani WJ, et al. *Analysis of Body Fluids in Clinical Chemistry; Approved Guideline (C49-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 2007.
- Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:22,40,44,52.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Pesce AJ, Kaplan LA, editors. *Methods in Clinical Chemistry*. St. Louis, MO: CV Mosby, 1987:75.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:2139,2164,2168.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000:3-176, 3-644, 3-714.

18. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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For laboratory professional use only.

NAME

Phosphorus2 (also referred to as Phos2)

INTENDED USE

The Phosphorus2 assay is used for the quantitation of phosphorus in human serum, plasma, or urine on the Alinity c system.

Measurements of phosphorus (inorganic) are used in the diagnosis and treatment of various disorders, including parathyroid gland and kidney diseases, and vitamin D imbalance.

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Phosphorus is a major component of bone mineral, phospholipids in cell membranes, and nucleic acids. Phosphorus acts as a major pH buffer in serum and urine.¹ Phosphorus is essential for normal muscle contractility, neurologic function, electrolyte transport, and oxygen carrying by hemoglobin (2,3 diphosphoglycerate).²

Serum phosphate is filtered at the glomerulus and is reabsorbed primarily in the proximal tubule. Parathyroid hormone (PTH) increases renal phosphate excretion by inhibiting the sodium-phosphate cotransporter in the proximal tubule, whereas vitamin D enhances intestinal phosphate absorption.¹

Hypophosphatemia may be caused by decreased intake, impaired intestinal absorption, redistribution into cells or bones, and renal losses (as in Fanconi syndrome).¹ Clinical manifestation of severe hypophosphatemia include encephalopathy, dilated cardiomyopathy, generalized muscle weakness leading to respiratory failure, destruction of muscles (rhabdomyolysis), and hemolysis.¹ Hemolysis of red blood cells may be seen with serum phosphorus levels < 0.5 mg/dL.³

Hyperphosphatemia is caused by excessive phosphate intake, increased intestinal absorption, redistribution from intracellular stores, or impaired renal excretion. Acute levels can increase the risk for precipitation of calcium phosphate in the kidney and soft tissues.¹ In chronic hyperphosphatemia (renal insufficiency) with a phosphate level > 6.5 mg/dL, mortality is higher due to the increased risk for the development of coronary and other vascular calcification, leading to increased systolic blood pressure and left ventricular hypertrophy.¹ A rapid increase in serum phosphorus can result in hypocalcemia, which can cause tetany, hypotension, seizures, and cardiac arrhythmias.¹

A urine phosphorus-creatinine ratio and urine fractional excretion of phosphorus may provide useful information about therapies altering intestinal absorption or urine phosphorus handling.⁴

PRINCIPLES OF THE PROCEDURE

The Phosphorus2 assay is an automated clinical chemistry assay. Inorganic phosphate reacts with ammonium molybdate to form a heteropolyacid complex.⁵ The absorbance at 340 nm is directly proportional to the inorganic phosphorus level in the sample.

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Kit Contents

Phosphorus2 Reagent Kit 04U03

NOTE: Some kit sizes may not be available. Please contact your local distributor.

Volumes (mL) listed in the following table indicate the volume per cartridge.

REF	04U0320	04U0330
Tests per cartridge	280	700
Number of cartridges per kit	4	4
Tests per kit	1120	2800
R1	20.7 mL	44.4 mL
R2	20.0 mL	44.0 mL

R2 Active ingredient: ammonium molybdate tetrahydrate 2.90 g/L.
Preservative: ProClin 300.

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Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials and all consumables contaminated with potentially infectious materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate regional, national, and institutional biosafety practices should be used for materials that contain, are suspected of containing, or are contaminated with infectious agents.⁶⁻⁹

The following warnings and precautions apply to: R1	
	
DANGER	Contains sulfuric acid.
H314	Causes severe skin burns and eye damage.
H290	May be corrosive to metals.
Prevention	
P234	Keep only in original container.
P260	Do not breathe mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.

Response	
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.
P310	Immediately call a POISON CENTER or doctor / physician.
P390	Absorb spillage to prevent material damage.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: **R2**

WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
H402*	Harmful to aquatic life.
H412	Harmful to aquatic life with long lasting effects.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P273	Avoid release to the environment.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EC 1272/2008 (CLP) has been implemented.

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Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 2 hours before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

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Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position.
Onboard	System Temperature	28 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

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Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

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INSTRUMENT PROCEDURE

The Phosphorus2 assay file must be installed on the Alinity c system prior to performing the assay.

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Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types, collection tube types, and anticoagulants have not been verified with this assay.

Specimen Type	Collection Vessel	Special Conditions
Serum	Serum tubes	
	Serum separator tubes	
Plasma	Dipotassium EDTA tubes	
	Lithium heparin tubes	
	Lithium heparin separator tubes	
	Sodium heparin tubes	

Specimen Type	Collection Vessel	Special Conditions
Urine (random specimens or timed specimens collected over intervals shorter than 24 hours) ¹⁰	Clean glass or plastic container	
Urine (24 hour) ¹⁰	Clean glass or plastic container	Urine specimens should be collected in 6 mol/L HCl, 20 to 30 mL, to avoid precipitation of phosphate complexes. ¹⁰

- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.

The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - pooled specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
 - specimens with fungal growth
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation. For accurate results, plasma specimens should be free of platelets and other particulate matter. Ensure centrifugation is adequate to remove platelets.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low-speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low-speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen Type	Temperature	Collection Vessel	Maximum Storage Time	Special Instructions	
Serum/ Plasma	Room temperature (20 to 25°C)	Serum tubes	24 hours		
		Serum separator tubes			
		Lithium heparin plasma tubes Lithium heparin separator tubes Sodium heparin tubes Dipotassium EDTA tubes			
	2 to 8°C	Serum tubes	3 days		
		Serum separator tubes			
		Lithium heparin plasma tubes Lithium heparin separator tubes Sodium heparin tubes Dipotassium EDTA tubes			
Urine	Room temperature (20 to 25°C)	Clean glass or plastic container	4 days	Acidify to pH < 5	
		Clean glass or plastic container		7 days	Acidify to pH < 5
		Clean glass or plastic container		1 month	Acidify to pH < 5 Avoid multiple freeze/thaw cycles.

It is the responsibility of the individual laboratory to determine specific specimen stability criteria for their laboratory per their laboratory workflow.

For additional information on sample handling and processing, refer to CLSI GP44-A4.¹² The storage information provided here is based on references or data maintained by the manufacturer.

Each laboratory may establish a range around -20°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Stored specimens must be inspected for particulates. If present, mix with a low-speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

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- Phosphorus2 assay file
- 04V6201 Consolidated Chemistry Calibrator
- Controls containing phosphorus
- Saline (0.85% to 0.90% NaCl) for specimen dilution

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- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- Minimum sample cup volume is calculated by the system and printed on the Order List report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
 - Sample volume for single test: 3.6 μ L (serum/plasma); 16.0 μ L (urine).

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the Consolidated Chemistry Calibrator package insert **[REF]** 04V6201 and/or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Serum/Plasma

The standard dilution factor, applied automatically by the system software to all results for the Phosphorus2 assay, is 1:1.73. Samples with a phosphorus value exceeding 25.3 mg/dL (8.17 mmol/L) are flagged with the code "> 25.3 mg/dL" (> 8.17 mmol/L). Automated or manual sample dilutions have not been evaluated for the Phosphorus2 assay.

Urine

The standard dilution factor, applied automatically by the system software to all results for the Phosphorus2 assay, is 1:17.33. Samples with a phosphorus value exceeding 186.2 mg/dL (60.14 mmol/L) are flagged with the code "> 186.2 mg/dL" (> 60.14 mmol/L). Automated or manual sample dilutions have not been evaluated for the Phosphorus2 assay.

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Calibration is stable for approximately 14 days (336 hours) but is required with each change in reagent cartridge. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

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- At least 2 levels of controls (low and high) are to be run every 24 hours.
- Run both levels of control with each cartridge change.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

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Quality Control Guidance

Refer to "Basic QC Practices" by James O. Westgard, Ph.D. for guidance on laboratory quality control practices.¹³

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RESULTS

Calculation

The Phosphorus2 (Phos2) assay utilizes the Linear data reduction method to generate a calibration and results for both the serum/plasma and urine applications.

Urine sample quantification (Phos2-U) is achieved using the calibration generated with the Phosphorus2 (Phos2) assay parameter file.

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Reportable Interval

Based on representative data for the limit of quantitation (LoQ) and the limit of detection (LoD), the ranges over which results can be reported are provided below according to the definitions from CLSI EP34, 1st ed.¹⁴

Serum

	mg/dL	mmol/L
Analytical Measuring Interval (AMI) ^a	0.5 - 25.3	0.16 - 8.17
Reportable Interval ^b	0.2 - 25.3	0.06 - 8.17

^a AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in mg/dL (mmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

^b The reportable interval extends from the LoD to the upper limit of the AMI.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the AMI. Samples with a phosphorus value below 0.5 mg/dL (0.16 mmol/L) are reported as "< 0.5 mg/dL" (< 0.16 mmol/L).

Urine

	mg/dL	mmol/L
Analytical Measuring Interval (AMI) ^a	2.4 - 186.2	0.78 - 60.14
Reportable Interval ^b	1.0 - 186.2	0.32 - 60.14

^a AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in mg/dL (mmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

^b The reportable interval extends from the LoD to the upper limit of the AMI.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the AMI. Samples with a phosphorus value below 2.4 mg/dL (0.78 mmol/L) are reported as "< 2.4 mg/dL" ("< 0.78 mmol/L").

LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- Substances that demonstrated interference with the Phosphorus2 assay are listed in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- Potential interference has not been evaluated for substances other than those described in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- In very rare cases, gammopathy may cause unreliable results.¹⁵

EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

Reference Range (Serum)

Category	Range (mg/dL)	Range (mmol/L)
Children ³	4.0-7.0	1.29-2.26
Adults* ³	2.5-4.5	0.81-1.45

* Adult plasma reference intervals are approximately 0.2 mg/dL (0.06 mmol/L) lower than for serum.¹⁶

NOTE: The system will flag results based only on the serum reference range in the assay file.

Reference Range (Urine)¹⁷

Category	Range	Range
24 Hour (adults)	0.4-1.3 g/day	12.9-42.0 mmol/day
Random (adult males)	5-189 mg/dL	1.6-61 mmol/L
Random (adult females)	7-148 mg/dL	2.3-48 mmol/L

24-Hour Urinary Excretion

To convert results from mg/dL to g/day (24-hour urinary excretion)

$$24\text{-hour excretion} = [(V \times c) \div 100\,000] \text{ g/day}$$

Where:

V = 24-hour urine volume (mL)

c = analyte concentration (mg/dL)

To convert results from mmol/L to mmol/day (24-hour urinary excretion)

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Representative performance data are provided in this section.

Results obtained in individual laboratories may vary.

The Alinity c system and the ARCHITECT c System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity c system.

Precision

Within-Laboratory Precision

Serum

A study was performed based on guidance from CLSI EP05-A3.¹⁸ Testing was conducted using 3 lots of the Phosphorus2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 1 instrument. Two controls and 3 human serum panels were tested in a minimum of 2 replicates, twice per day on 20 days on 3 reagent lot/calibrator lot combinations, where a unique reagent lot and a unique calibrator lot are paired. The performance from a representative combination is shown in the following table.

Sample	n	Mean (mg/dL)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	3.4	0.03	0.8	0.03 (0.03-0.05)	1.0 (1.0-1.4)
Control Level 2	80	7.4	0.03	0.3	0.06 (0.06-0.06)	0.8 (0.8-0.8)
Panel A	80	1.1	0.02	2.0	0.03 (0.02-0.05)	2.5 (1.8-4.2)
Panel B	80	8.9	0.04	0.5	0.06 (0.06-0.06)	0.7 (0.7-0.7)
Panel C	80	22.8	0.13	0.6	0.16 (0.14-0.19)	0.7 (0.6-0.8)

^a Includes within-run, between-run, and between-day variability.

^b Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

Sample	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	1.10	0.006	0.5	0.009 (0.009-0.009)	0.8 (0.8-0.8)
Control Level 2	80	2.38	0.008	0.3	0.019 (0.015-0.019)	0.8 (0.6-0.8)
Panel A	80	0.36	0.005	1.5	0.008 (0.007-0.008)	2.4 (2.1-2.4)
Panel B	80	2.88	0.009	0.3	0.021 (0.017-0.021)	0.7 (0.6-0.7)
Panel C	80	7.37	0.044	0.6	0.053 (0.046-0.061)	0.7 (0.6-0.8)

^a Includes within-run, between-run, and between-day variability.

^b Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

Urine

A study was performed based on guidance from CLSI EP05-A3.¹⁸ Testing was conducted using 3 lots of the Phosphorus2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 1 instrument. Two controls and 3 urine panels were tested in a minimum of 2 replicates, twice per day on 20 days on 3 reagent lot/calibrator lot combinations, where a unique reagent lot and a unique calibrator lot are paired. The performance from a representative combination is shown in the following table.

Sample	n	Mean (mg/dL)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	24.8	0.12	0.5	0.24 (0.24-0.26)	0.9 (0.9-1.0)
Control Level 2	80	48.1	0.26	0.5	0.31 (0.26-0.41)	0.7 (0.5-0.8)

Sample	n	Mean (mg/dL)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Panel A	80	5.2	0.08	1.5	0.21 (0.18-0.21)	4.1 (3.4-4.1)
Panel B	80	92.2	0.63	0.7	1.09 (1.09-2.15)	1.2 (1.2-2.3)
Panel C	80	161.3	0.93	0.6	1.11 (1.11-1.78)	0.7 (0.7-1.1)

^a Includes within-run, between-run, and between-day variability.

^b Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

Sample	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	8.02	0.037	0.5	0.076 (0.076-0.082)	0.9 (0.9-1.0)
Control Level 2	80	15.52	0.083	0.5	0.099 (0.085-0.132)	0.6 (0.5-0.8)
Panel A	80	1.66	0.026	1.6	0.066 (0.057-0.066)	4.0 (3.4-4.0)
Panel B	80	29.78	0.202	0.7	0.351 (0.351-0.695)	1.2 (1.2-2.3)
Panel C	80	52.10	0.301	0.6	0.360 (0.360-0.574)	0.7 (0.7-1.1)

^a Includes within-run, between-run, and between-day variability.

^b Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot combinations.

Reproducibility

Serum

A study was performed based on guidance from CLSI EP05-A3.¹⁸ Testing was conducted using 1 lot of the Phosphorus2 reagents, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 3 human serum panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

Sample	n	Mean (mg/dL)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	3.4	0.02	0.5	0.02	0.7	0.02	0.7
Control Level 2	90	7.4	0.03	0.4	0.04	0.5	0.04	0.5
Panel A	90	1.1	0.02	1.4	0.02	1.6	0.02	1.7
Panel B	90	8.9	0.03	0.4	0.04	0.4	0.05	0.5
Panel C	90	22.8	0.06	0.3	0.15	0.6	0.15	0.7

^a Includes repeatability (within-run), between-run, and between-day variability.

^b Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

Sample	n	Mean (mmol/L)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	1.11	0.006	0.5	0.008	0.7	0.008	0.7
Control Level 2	90	2.39	0.010	0.4	0.012	0.5	0.013	0.5

Sample	n	Mean (mmol/L)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Panel A	90	0.35	0.006	1.6	0.006	1.7	0.006	1.8
Panel B	90	2.88	0.010	0.4	0.011	0.4	0.016	0.5
Panel C	90	7.36	0.020	0.3	0.047	0.6	0.048	0.7

^a Includes repeatability (within-run), between-run, and between-day variability.

^b Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

Urine

A study was performed based on guidance from CLSI EP05-A3.¹⁸ Testing was conducted using 1 lot of the Phosphorus2 reagents, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 3 urine panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

Sample	n	Mean (mg/dL)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	25.6	0.25	1.0	0.26	1.0	0.39	1.5
Control Level 2	90	48.7	0.35	0.7	0.37	0.8	0.59	1.2
Panel A	90	4.8	0.16	3.2	0.18	3.8	0.25	5.2
Panel B	90	91.4	0.49	0.5	0.53	0.6	1.38	1.5
Panel C	90	160.0	0.74	0.5	0.91	0.6	2.51	1.6

^a Includes repeatability (within-run), between-run, and between-day variability.

^b Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

Sample	n	Mean (mmol/L)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	8.26	0.080	1.0	0.085	1.0	0.125	1.5
Control Level 2	90	15.72	0.113	0.7	0.121	0.8	0.189	1.2
Panel A	90	1.56	0.050	3.2	0.059	3.8	0.080	5.1
Panel B	90	29.52	0.158	0.5	0.174	0.6	0.446	1.5
Panel C	90	51.67	0.237	0.5	0.294	0.6	0.809	1.6

^a Includes repeatability (within-run), between-run, and between-day variability.

^b Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

Accuracy

A study was performed to estimate the bias of the Phosphorus2 assay relative to standard reference material (NIST SRM 2186I Standard). Testing was conducted using 3 concentrations of the standard across 1 lot of the Phosphorus2 reagents, 1 lot of the Consolidated Chemistry Calibrator, and 1 instrument. The bias ranged from -0.2% to 0.0% for serum and -5.1% to -0.1% for urine.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.¹⁹ Testing was conducted using 3 lots of the Phosphorus2 reagents on each of 2 instruments over a minimum of 3 days. The maximum observed limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized below.

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^c The LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from n ≥ 60 replicates of low-analyte level samples.

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^c The LoQ presented in the table is in alignment with the LoQ for the Phosphorus2 assay on the ARCHITECT c System. The observed LoQ on the Alinity c system was 2.2 mg/dL (0.71 mmol/L). This LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from n ≥ 60 replicates of low-analyte level samples.

Linearity

A study was performed based on guidance from CLSI EP06-A.²⁰

This assay is linear across the analytical measuring interval of 0.5 to 25.3 mg/dL (0.16 to 8.17 mmol/L) for serum, and 2.4 to 186.2 mg/dL (0.78 to 60.14 mmol/L) for urine.

Analytical Specificity

Interference

These studies were performed on the ARCHITECT c System.

Serum

Potentially Interfering Endogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.²¹ Each substance was tested at 2 levels of the analyte (approximately 3.5 mg/dL and 5.0 mg/dL).

No significant interference (interference within ± 10%) was observed at the following concentrations.

Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Bilirubin - conjugated	30 mg/dL	356 µmol/L
Bilirubin - unconjugated	60 mg/dL	1026 µmol/L
Hemoglobin	125 mg/dL	1.25 g/L
Total protein	15 g/dL	150 g/L
Triglycerides	750 mg/dL	8.47 mmol/L

Interference beyond ± 10% (based on 95% Confidence Interval [CI]) was observed at the concentrations shown below for the following substances.

Potentially Interfering Substance	Interference beyond ± 10% (based on 95% Confidence Interval [CI])				% Interference (95% CI)
	Interferent Level		Analyte Level		
	Default Units	Alternate Units	Default Units	Alternate Units	
Bilirubin - conjugated	40 mg/dL	474 µmol/L	3.5 mg/dL	1.13 mmol/L	-12% (-13%, -11%)
Hemoglobin	250 mg/dL	2.50 g/L	3.5 mg/dL	1.13 mmol/L	14% (14%, 14%)
Triglycerides	1500 mg/dL	16.9 mmol/L	3.5 mg/dL	1.13 mmol/L	15% (14%, 16%)
Triglycerides	1500 mg/dL	16.9 mmol/L	5.0 mg/dL	1.62 mmol/L	10% (9%, 11%)

Potentially Interfering Exogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.²¹ Each substance was tested at 2 levels of the analyte (approximately 3.5 mg/dL and 5.0 mg/dL).

No significant interference (interference within ± 10%) was observed at the following concentrations.

Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Acetaminophen	160 mg/L	1059 µmol/L
Acetylcysteine	150 mg/L	920 µmol/L
Acetylsalicylic acid	30 mg/L	167 µmol/L
Aminosalicylic acid (p-aminosalicylic acid)	47 mg/dL	3074 µmol/L
Amiodarone	5 mg/dL	77.5 µmol/L
Ampicillin-Na	80 mg/L	215 µmol/L
Ascorbic acid	60 mg/L	341 µmol/L
Biotin	4250 ng/mL	17.4 µmol/L
Ca-dobesilate	60 mg/L	143 µmol/L
Cefotaxime	53 mg/dL	1166 µmol/L
Cefoxitin	6600 mg/L	15.4 mmol/L
Cyclosporine	2 mg/L	1.66 µmol/L
Desacetylcefotaxime	6 mg/dL	145 µmol/L
Doxycycline	20 mg/L	45.0 µmol/L
Ibuprofen	220 mg/L	1067 µmol/L
Levodopa	8 mg/L	40.6 µmol/L
Mannitol	4 g/L	22.0 mmol/L
Methicillin (sodium)	8 mg/dL	198 µmol/L
Methotrexate	140 mg/dL	3080 µmol/L
Methyl dopa	25 mg/L	118 µmol/L
Metronidazole	130 mg/L	759 µmol/L
Naproxen	36 mg/dL	1562 µmol/L
Nitroglycerin	0.015 mg/L	0.0660 µmol/L
Omeprazole	0.9 mg/dL	26.1 µmol/L
Phenylbutazone	330 mg/L	1069 µmol/L
Primidone	6 mg/dL	275 µmol/L
Promethazine	0.03 mg/dL	1.06 µmol/L
Ranitidine	2 mg/dL	63.6 µmol/L
Rifampicin	50 mg/L	61.0 µmol/L
Sodium heparin	4 U/mL	NA
Theophylline (1,3-dimethylxanthine)	60 mg/L	333 µmol/L

NA = Not applicable

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI]) was observed at the concentrations shown below for the following substances.

Potentially Interfering Substance	Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI])				
	Interferent Level		Analyte Level		% Interference (95% CI)
	Default Units	Alternate Units	Default Units	Alternate Units	
Mannitol	5 g/L	27.5 mmol/L	3.5 mg/dL	1.13 mmol/L	-11% (-11%, -10%)
Mannitol	10 g/L	54.9 mmol/L	5.0 mg/dL	1.62 mmol/L	-10% (-11%, -10%)

Urine

Potentially Interfering Endogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.²¹ Each substance was tested at 2 levels of the analyte (approximately 30 mg/dL and 90 mg/dL).

No significant interference (interference within $\pm 10\%$) was observed at the following concentrations.

Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Acetic acid (8.5N)	6.25 mL/dL	531 mmol/L
Ascorbate	220 mg/dL	12.6 mmol/L
Boric acid	250 mg/dL	40.5 mmol/L
Glucose	1000 mg/dL	55.5 mmol/L
Hydrochloric acid (6N)	2.5 mL/dL	150 mmol/L
Nitric acid (6N)	5.0 mL/dL	300 mmol/L
Protein	50 mg/dL	0.500 g/L
Sodium carbonate	1.25 g/dL	118 mmol/L
Sodium fluoride	400 mg/dL	95.2 mmol/L
Sodium oxalate	60 mg/dL	4476 μ mol/L

Potentially Interfering Exogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.²¹ Each substance was tested at 2 levels of the analyte (approximately 30 mg/dL and 90 mg/dL).

No significant interference (interference within $\pm 10\%$) was observed at the following concentrations.

Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Acetaminophen	16 mg/dL	1059 μ mol/L
Acetylcysteine	15 mg/dL	920 μ mol/L
Biotin	4250 ng/mL	17.4 μ mol/L
Ibuprofen	22 mg/dL	1067 μ mol/L
Mannitol	18 g/L	98.8 mmol/L

Interferences from medication or endogenous substances may affect results.²²

Method Comparison

A study was performed based on guidance from CLSI EP09c, 3rd ed.²³ using the Passing-Bablok regression method.

Phosphorus2 on the Alinity c system vs. Phosphorus on the ARCHITECT c System						
	n	Units	Correlation Coefficient		Concentration Range	
			Intercept	Slope	Intercept	Slope
Serum	108	mg/dL (mmol/L)	1.00	0.0	1.05	1.0-22.8 (0.34-7.36)
Urine	110	mg/dL (mmol/L)	1.00	-0.6 (-0.21)	1.05	6.2-167.9 (2.00-54.22)

Phosphorus2 on the Alinity c system vs. Phosphorus2 on the ARCHITECT c System						
	n	Units	Correlation Coefficient		Concentration Range	
			Intercept	Slope	Intercept	Slope
Serum	108	mg/dL (mmol/L)	1.00	0.0	1.00	1.1-23.7 (0.36-7.67)
Urine	111	mg/dL (mmol/L)	1.00	-0.3 (-0.08)	1.01	4.4-174.1 (1.42-56.23)

BIBLIOGRAPHY

- Goldman L, Schafer AI, editors. *Goldman-Cecil Medicine*. 25th ed. Elsevier/Saunders; 2016.
- McPherson RA, Pincus MR, editors. *Henry's Clinical Diagnosis and Management by Laboratory Methods*. 23rd ed. St. Louis, Missouri: Elsevier; 2017.
- Burtis CA, Bruns DE, editors. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, MO: Saunders Elsevier; 2015.
- Dominguez JR, Kestenbaum B, Chonchol M, et al. Relationships between serum and urine phosphorus with all-cause and cardiovascular mortality: the osteoporotic fractures in men (MrOS) study. *Am J Kidney Dis* 2013;61(4):555-563.
- Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry*. 5th ed. WB Saunders; 2001:chap 38.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 6th ed. Washington, DC: US Government Printing Office; June 2020.
- World Health Organization. *Laboratory Biosafety Manual*. 4th ed. Geneva: World Health Organization; 2020.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Rifai N, Horvath AR, Wittwer CT editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. Saunders Elsevier; 2018:1444.
- Kachhawa K, Kachhawa P, Varma M, et al. Study of the stability of various biochemical analytes in samples stored at different predefined storage conditions at an accredited laboratory of India. *J Lab Physicians* 2017;9(1):11-15.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking*. 1st ed. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Valentina Molinaris, Mario G. Bianchetti, Gregorio P. Milani, et al. Interferences in the measurement of circulating phosphate: a literature review. *Clin Chem Lab Med* 2020; 58(12): 1971–1977.
- Rifai N, Horvath AR, Wittwer CT, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. St. Louis, MO: Elsevier; 2018.
- Wu AHB, editor. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Elsevier Saunders; 2006.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3rd ed. CLSI Guideline EP09c. Wayne, PA: CLSI; 2018.

Key to Symbols

ISO 15223 Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
REF	List Number
SN	Serial number

Other Symbols

DISTRIBUTED IN THE USA BY	Distributed in the USA by
FOR USE WITH	Identifies products to be used together
INFORMATION FOR USA ONLY	Information needed for United States of America only
PRODUCT OF IRELAND	Product of Ireland
R1	Reagent 1
R2	Reagent 2
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

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For customers in the European Union: if, in the course of using this device, you have reason to believe that a serious incident has occurred, report it to the manufacturer and to your national authority.

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Revised March 2022.

REF 08P4320

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

WARNING: The Alinity c Hemoglobin A1c assay has significant interference with the fetal hemoglobin (HbF). Hemoglobin A1c results are invalid for patients with abnormal amounts of HbF, including those with known Hereditary Persistence of Fetal Hemoglobin. For more information regarding the specific concentrations of HbF that were found to interfere with the Alinity c Hemoglobin A1c assay, refer to the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity section of this package insert.

NAME

Alinity c Hemoglobin A1c Reagent Kit (also referred to as HbA1c)

INTENDED USE

The Alinity c Hemoglobin A1c assay is used in clinical laboratories for the quantitative *in vitro* measurement of percent hemoglobin A1c (NGSP) or HbA1c fraction mmol/mol (IFCC) in human whole blood and hemolysate on the Alinity c analyzer.

Hemoglobin A1c measurements are used as an aid in the diagnosis of diabetes mellitus, as an aid to identify patients who may be at risk for developing diabetes mellitus, and for the monitoring of long-term blood glucose control in individuals with diabetes mellitus.

SUMMARY AND EXPLANATION OF THE TEST

HbA1c is the fraction of hemoglobin A that is first reversibly, then irreversibly glycosylated at one or both *N*-terminal valines of the β -chain.¹ The longer red blood cells are in circulation and the higher the ambient glucose levels, the higher the concentration of HbA1c. HbA1c reflects the average blood glucose level during the preceding 2 to 3 months. The HbA1c assay is useful as an aid in the:

- diagnosis of diabetes mellitus,
- identification of patients at risk for developing diabetes, and
- monitoring of patients with diabetes mellitus.²⁻⁶

For monitoring diabetic patients, it is recommended that glycemic goals are individualized following current professional society recommendations.⁷ As recommended by the American Diabetes Association (ADA), patients in the range of 5.7 to 6.4 %HbA1c (39 to 46 mmol/mol) would be in the category of increased risk for diabetes and results \geq 6.5% (48 mmol/mol) may aid in the diagnosis of diabetes.⁷

Several studies, including the Diabetes Control and Complications Trial (DCCT), have shown that long-term control of diabetes can prevent complications such as cardiovascular disease, retinopathy, nephropathy, and neuropathy. Measurement of HbA1c can be invaluable in the monitoring of glycemic control of diabetic patients.⁸⁻¹⁰

This method is certified to the National Glycohemoglobin Standardization Program (NGSP), standardized to International Federation of Clinical Chemistry and Laboratory Medicine (IFCC), and traceable to DCCT.

PRINCIPLES OF THE PROCEDURE

The Alinity c Hemoglobin A1c assay consists of two separate concentration measurements: glycosylated hemoglobin (HbA1c) and total hemoglobin (THb). The two concentrations are used to determine the percent HbA1c (NGSP units) or the hemoglobin fraction in mmol/mol (IFCC units).

The individual concentration values of HbA1c and THb generated by the Alinity c Hemoglobin A1c assay are used only for calculating the percent hemoglobin A1c or HbA1c fraction, and must not be used individually for diagnostic purposes.

The anticoagulated whole blood specimen is lysed automatically on the system for the Whole Blood application or may be lysed manually using the Hemoglobin A1c Diluent (**A1cDIL**) for the Hemolysate application.

Glycosylated Hemoglobin (HbA1c)

The Alinity c Hemoglobin A1c assay utilizes an enzymatic method that specifically measures *N*-terminal fructosyl dipeptides of the β -chain of HbA1c.

- In the pretreatment process, the erythrocytes are lysed and the hemoglobin is transformed to methemoglobin by reaction with sodium nitrite.
- With the addition of Reagent 1 (**R1**) to the sample, the glycosylated *N*-terminal dipeptide (fructosyl-VH) of the β -chain of hemoglobin is cleaved by the action of protease. The hemoglobin is transformed to stable methemoglobin azide by the action of sodium azide and the concentration of the hemoglobin is determined by measuring absorbance.
- Addition of Reagent 2 (**R2**) starts a reaction and fructosyl peptide oxidase (FPOX) is allowed to react with fructosyl-VH. The HbA1c concentration is measured by determining the resultant hydrogen peroxide.

Total Hemoglobin (THb)

The hemoglobin is oxidized to stable methemoglobin azide by the action of sodium nitrite and sodium azide and the concentration of the hemoglobin is determined by measuring absorbance (sample + **R1**).

Hemoglobin A1c Calculations¹¹

The final result is expressed as %HbA1c (NGSP) or mmol/mol HbA1c (IFCC) and is automatically calculated by the system from the HbA1c/THb ratio as follows:

mmol/mol HbA1c IFCC:

$$\text{HbA1c (mmol/mol)} = (\text{HbA1c/THb}) \times 1000$$

%HbA1c DCCT/NGSP:

$$\text{HbA1c (\%)} = \text{IFCC} \times 0.09148 + 2.152$$

Methodology: Enzymatic

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity c Hemoglobin A1c Reagent Kit 08P43

NOTE: This product is composed of 3 components, which are packaged as a 2 cartridge reagent set. Both cartridges are required to perform the assay.

Volumes (mL) listed in the table below indicate the volume per cartridge set.

REF	08P4320
Tests per cartridge set	260
Number of cartridge sets per kit	5
Tests per kit	1300
R1	46.0 mL
R2	18.3 mL
A1cDIL	67.6 mL
R1 Active ingredients: 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino) phenothiazine sodium salt (0.000817%), protease (bacterial) (< 1 MU/dL). Inactive ingredient: sodium azide (< 0.1%) as a stabilizer. Preservatives: sodium azide (< 0.1%), ProClin 300 (0.05%).	
R2 Active ingredients: peroxidase (horseradish) (5 to 15 kU/dL), fructosyl-peptide-oxidase (<i>E. coli</i> , recombinant) (300 to 900 U/dL). Preservative: ofloxacin (0.001%).	
A1cDIL Active ingredient: sodium nitrite (> 0.05 to < 0.3%). Preservative: ProClin 300 (0.01%).	

Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹²⁻¹⁵

The following warnings and precautions apply to: R1	
	
DANGER	Contains morpholinoethanesulfonic acid, monohydrate*, N,N-dimethylformamide, methylisothiazolones and sodium azide.
H317	May cause an allergic skin reaction.
H360	May damage fertility or the unborn child.
H316*	Causes mild skin irritation.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P201	Obtain special instructions before use.
P261	Avoid breathing mist / vapors / spray.
P280	Wear protective gloves / protective clothing / eye protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P308+P313	IF exposed or concerned: Get medical advice / attention.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

The following warnings and precautions apply to: R2	
H316*	Causes mild skin irritation.
P332+P313*	If skin irritation occurs: Get medical advice / attention.

* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

The following warnings and precautions apply to: A1cDIL	
	
WARNING	Contains sodium nitrite*, maleic acid and methylisothiazolones.
H317	May cause an allergic skin reaction.
H402*	Harmful to aquatic life.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P273*	Avoid release to the environment.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

Follow local chemical disposal regulations based on your location along with recommendations and content in the Safety Data Sheet to determine the safe disposal of this product.

For the most current hazard information, see the product Safety Data Sheet.

Safety Data Sheets are available at www.corelaboratory.abbott or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Reagents are shipped refrigerated or on cold packs.
- Upon receipt, place reagent cartridges in an upright position for 8 hours before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 8 hours before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position.
Onboard	System Temperature	50 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

A1cDIL can be replaced independently of **R1** and **R2**.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range.

Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity c Hemoglobin A1c assay files must be installed on the Alinity c analyzer prior to performing the assay.

For the Whole Blood application, the following assay files are required:

- THbWB
- HbA1cWB
- A1cWB (for results in mmol/mol (IFCC) units) or %A1cWB (for results in %HbA1c (NGSP) units)

For the Hemolysate application, the following assay files are required:

- THbH
- HbA1cH
- A1cH (for results in mmol/mol (IFCC) units) or %A1cH (for results in %HbA1c (NGSP) units)

The Alinity ci-series software version 2.0 or higher must be installed prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen type used for this assay is whole blood only.

The specimen type listed below was verified for use with this assay.

Other specimen types and anticoagulants have not been verified with this assay.

Specimen Type	Collection Vessel	Special Conditions
Whole Blood	Collection tubes Acceptable anticoagulants are: Dipotassium EDTA Lithium heparin Sodium heparin Sodium fluoride/disodium EDTA Tripotassium EDTA	Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated.

- Use only whole blood specimens collected by standard venipuncture techniques into plastic tubes.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Analyze fresh specimens if possible.
- Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. Refer to the Alinity ci-series Operations Manual, Section 10.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- **Do not centrifuge samples.**
- Follow the tube manufacturer's collection instructions for specimen collection tubes.
- For testing whole blood samples less than 600 µL, use 12 x 75 mm polypropylene conical bottom tubes.
- Visually inspect the specimens. If fibrin clots or particulate matter is observed, remove with a clean applicator stick.
- Mix all specimens thoroughly by low speed vortexing or gently inverting 10 times prior to loading onto the Alinity c analyzer.

Prepare frozen whole blood specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.

Do not use Alinity c Sample Cups for whole blood samples. Refer to the Assay Procedure section of this package insert for further information.

Specimen Storage

Analyze fresh specimens if possible.

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Whole Blood	Room temperature	8 hours	
	2 to 8°C	7 days	
Hemolysate	Room temperature	4 hours	Do not freeze hemolyzed specimens.
	2 to 8°C	24 hours	

Whole Blood: If testing will be delayed more than 7 days, store at -70°C or colder.

CAUTION: Whole blood specimens that require freezing must be stored at -70°C or colder.

Avoid more than one freeze/thaw cycle.

Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

08P43 Alinity c Hemoglobin A1c Reagent Kit

Materials Required but not Provided

- Alinity c Hemoglobin A1c assay files
- 08P4301 Alinity c Hemoglobin A1c Calibrators

Whole Blood application:

- Commercially available whole blood controls
- 12 x 75 mm polypropylene conical bottom tubes

Hemolysate application:

- 08P4310 Alinity c Hemoglobin A1c Controls or other commercially available controls
- Calibrated adjustable pipette capable of measuring 222 μL
- Calibrated micropipette capable of measuring 10 μL
- Vortex (optional)

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

Refer to the following table.

NOTE: No more than one replicate can be sampled from a sample cup or tube.

For information on ordering patient samples and controls, refer to the Alinity ci-series Operations Manual, Section 5.

	Sample Type	
	Whole Blood	Hemolysate
1. Order the application:	Whole Blood	Hemolysate
2. Prepare samples according to the instructions from this section of the package insert:	Assay Procedure, Whole Blood Preparation	Assay Procedure, Hemolysate Preparation
3. Load, then run:	Whole Blood samples	Prepared, Hemolyzed samples
4. Report results in percent (NGSP) or mmol/mol (IFCC) hemoglobin A1c.		

Whole Blood Preparation

- WARNING: Do not use Alinity c sample cups.**
- WARNING: Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. Refer to Section 10 of the Alinity ci-series Operations Manual.**
- Select the appropriate sample vessel using the following table:

If sample is:	Then:
< 200 μL	Do not use the Whole Blood application. Follow the instructions in the Assay Procedure, Hemolysate Preparation section of this package insert.
200 μL to 600 μL	Use 12 x 75 mm polypropylene conical bottom tubes only.
> 600 μL	Use suitable tubes or polypropylene conical bottom tubes.
> 78 mm in height in a large tube	Pipette 600 μL of sample into a suitable tube or polypropylene conical bottom tube.

Hemolysate Preparation

- The minimum sample volume requirement is 150 μL for Alinity c sample cups.
 - NOTE:** This volume supersedes the calculated minimum sample volume displayed on the User Interface.

Prepare the hemolysate samples as follows:

- Using a calibrated pipette, dispense 222 μL **A1cDIL** into a tube or sample cup.
- Using a calibrated micropipette, withdraw 10 μL of the well-mixed whole blood patient specimen.
- Wipe excess blood from the exterior of the pipette to ensure accurate transfer of the sample.
- Insert the pipette into the tube or sample cup containing the **A1cDIL** allowing the tip of the pipette to just make contact with the surface of the **A1cDIL** and dispense the 10 μL sample (1:23.2 dilution).
- Withdraw and dispense twice to rinse the pipette, always keeping the tip of the pipette in contact with the fluid in the tube.
- Mix hemolysate thoroughly by low speed vortexing or by gently inverting 10 times. Avoid foaming.
- Allow the hemolysate to stand for a minimum of 1 minute at room temperature prior to testing.
- If the hemolysate is prepared in a tube, transfer to a sample cup and place the cup on the instrument.

NOTE: The number of tests per kit is based on the 222 μL of **A1cDIL** and 10 μL of specimen volumes stated in steps 1 and 2 above. However, alternate volumes may be used for the 1:23.2 dilution, such as 555 μL of **A1cDIL** and 25 μL of specimen.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Do not use Alinity c Sample Cups for whole blood samples.
- Minimum sample volume requirements:
 - Sample volume for single test: 200 μL (Whole Blood), 150 μL (Hemolysate).

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the Alinity c Hemoglobin A1c Calibrators package insert and the Alinity c Hemoglobin A1c Controls package insert or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples cannot be diluted for the Alinity c Hemoglobin A1c assay.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration is stable for approximately 50 days (1200 hours), but is required with each change in reagent lot. Verify calibration with 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Both the Whole Blood and Hemolysate applications use the Alinity c Hemoglobin A1c Calibrators which are supplied separately.

Alinity c Hemoglobin A1c Calibrators are traceable to NGSP and IFCC reference methods.

Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls are to be run every 24 hours. The Alinity c Hemoglobin A1c assay uses:
 - Commercially available whole blood controls for the Whole Blood application.
 - Alinity c Hemoglobin A1c Controls or other commercially available controls for the Hemolysate application.
Refer to the Alinity c Hemoglobin A1c Controls Value Sheet for NGSP and IFCC ranges of the Alinity c Hemoglobin A1c Controls.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to “Basic QC Practices” by James O Westgard, Ph.D. for guidance on laboratory quality control practices.¹⁶

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

- The individual concentrations of the Whole Blood and Hemolysate applications are measured by the analyzer.
NOTE: HbA1c or THb concentrations must not be used individually for clinical purposes.
- Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert.

Conventional Units (NGSP)

The percent HbA1c (%HbA1c) is automatically calculated by the system per the calculation provided in the Principles of the Procedure section.

SI Units (IFCC)

The hemoglobin A1c fraction (mmol/mol HbA1c) is automatically calculated by the system per the calculation provided in the Principles of the Procedure section.

Calculation

The Alinity c Hemoglobin A1c assay utilizes the Linear data reduction method to generate a calibration and results.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in %HbA1c (mmol/mol) which meets the limits of acceptable performance for linearity, imprecision, and bias.

The measuring interval of the Alinity c Hemoglobin A1c assay is 4.0 to 14.0 %HbA1c (20.22 to 129.51 mmol/mol HbA1c).

LIMITATIONS OF THE PROCEDURE

- Do not use Alinity c Sample Cups for whole blood samples. Refer to the Assay Procedure section of this package insert for further information.
- This assay must be performed by qualified laboratory personnel, under appropriate laboratory conditions, solely for the intended use of the assay.
- Do not centrifuge samples.
- Do not freeze specimens that have been hemolyzed with the **A1cDIL**.
- Whole blood specimens that require freezing must be stored at -70°C or colder.
- Do not overfill specimen collection tubes. Whole blood samples greater than 78 mm in height from the bottom of tube will result in an instrument error and results will not be generated. Refer to Section 10 of the Alinity ci-series Operations Manual.
- Use specimen collection tubes, or for sample volumes < 600 µL use the 12 x 75 mm polypropylene conical bottom tubes as recommended in the PROCEDURE, Materials Required but not Provided section of this package insert.
- **WARNING:** The Alinity c Hemoglobin A1c assay should not be used to diagnose diabetes during pregnancy. Hemoglobin A1c reflects the average blood glucose levels over the preceding 3 months (i.e., the average life span of a red blood cell) and therefore may be falsely low during pregnancy or any other condition associated with recent onset of hyperglycemia and/or decreased red blood cell survival.¹⁷⁻²⁰
- Blood transfusions may impact the HbA1c concentration in the patient sample.
- The Alinity c Hemoglobin A1c assay should not be used to diagnose or monitor diabetes in patients with the following conditions:¹⁷⁻²⁰
 - hemoglobinopathies except as demonstrated to produce acceptable performance (e.g., sickle cell trait - refer to the SPECIFIC PERFORMANCE CHARACTERISTICS section of this package insert)
 - abnormal red blood cell turnover (e.g., anemias from hemolysis and iron deficiency)
 - malignancies, and severe chronic hepatic and renal disease
- In cases of rapidly evolving Type 1 diabetes, the increase of HbA1c values might be delayed compared to the acute increase in glucose concentrations. In these conditions, diabetes mellitus must be diagnosed based on plasma glucose concentrations and/or the typical clinical symptoms.
- This test should not replace glucose testing for patients with Type 1 diabetes, pediatric patients, or pregnant women.
- The Alinity c Hemoglobin A1c assay is susceptible to interference effects from conjugated bilirubin at > 15.0 mg/dL (180 µmol/L) and unconjugated bilirubin at > 10.0 mg/dL (171 µmol/L).
- The observed bias for samples containing HbC, HbD, HbE, HbS and HbA2 may be impacted by the method used to determine the reference Hemoglobin A1c concentration.

- The Alinity c Hemoglobin A1c assay is susceptible to interference effects from HbF at > 5%. Glycated HbF is not detected by the Alinity c Hemoglobin A1c assay as it does not contain the β-chain that characterizes HbA1c. However, HbF is measured in the total hemoglobin constituent assay and as a consequence, specimens containing high amounts of HbF (> 5%) may result in lower than expected mmol/mol HbA1c values (IFCC) and %HbA1c values (NGSP).
- Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section of this package insert for specimen limitations.

EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics. For monitoring diabetic patients, it is recommended that glycemic goals are individualized following current professional society recommendations.⁷ The American Diabetes Association (ADA) recommendations⁷ are summarized in the following table.

HbA1c Value	Glycemic Goal
< 8 %HbA1c (64 mmol/mol)	Less stringent
< 7 %HbA1c (53 mmol/mol)	General (non-pregnant adults)
< 6.5 %HbA1c (48 mmol/mol)	More stringent

HbA1c values above 6.5 %HbA1c (48 mmol/mol) are an indication of hyperglycemia during the preceding 2 to 3 months or longer. According to the recommendations of the ADA, HbA1c values above 6.5 %HbA1c (48 mmol/mol) are suitable for the diagnosis of diabetes mellitus. Patients with HbA1c values in the range of 5.7 to 6.4 %HbA1c (39 to 46 mmol/mol) may be at a risk of developing diabetes.^{3, 21}

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity c analyzer and the ARCHITECT c System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity c analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.²² Testing was conducted using 3 lots of the Alinity c Hemoglobin A1c Reagent Kit, 3 lots of the Alinity c Hemoglobin A1c Calibrators, and 1 lot of commercially available controls and 3 instruments. Three controls and 3 human whole blood panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

Whole Blood (%HbA1c)

Sample	Instrument	n	Mean (%HbA1c (NGSP))	Within-Run (Repeatability)		Within-Laboratory (Total) ^a		Precision with Additional Component of Between-Lot ^b	
				SD	%CV	SD	%CV	SD	%CV
Control 1	1	240	5.1	0.01	0.1	0.05	1.0	0.10	1.9
	2	240	5.1	0.04	0.7	0.06	1.2	0.10	1.9
	3	240	5.1	0.02	0.5	0.06	1.1	0.11	2.1
Control 2	1	240	9.2	0.03	0.3	0.07	0.7	0.12	1.3
	2	240	9.3	0.03	0.4	0.07	0.7	0.13	1.5
	3	240	9.3	0.03	0.4	0.06	0.7	0.11	1.2
Control 3	1	240	13.5	0.03	0.2	0.06	0.5	0.15	1.1
	2	240	13.6	0.04	0.3	0.08	0.6	0.16	1.2
	3	240	13.6	0.04	0.3	0.06	0.5	0.12	0.9
Panel near 4.0 %HbA1c	1	240	4.4	0.03	0.7	0.06	1.3	0.07	1.7
	2	240	4.4	0.03	0.7	0.05	1.2	0.07	1.5
	3	240	4.4	0.03	0.7	0.06	1.3	0.08	1.8
Panel near 6.5 %HbA1c	1	240	6.2	0.02	0.3	0.03	0.5	0.07	1.2
	2	240	6.2	0.02	0.4	0.04	0.6	0.07	1.2
	3	240	6.2	0.03	0.4	0.04	0.6	0.08	1.3

Sample	Instrument	n	Mean (%HbA1c (NGSP))	Within-Run (Repeatability)		Within-Laboratory (Total) ^a		Precision with Additional Component of Between-Lot ^b	
				SD	%CV	SD	%CV	SD	%CV
Panel near 8.0 %HbA1c	1	240	7.9	0.03	0.3	0.05	0.6	0.09	1.2
	2	240	7.9	0.04	0.5	0.06	0.7	0.09	1.2
	3	240	7.9	0.04	0.5	0.05	0.7	0.11	1.3

^a Within-Laboratory variability contains within-run, between-run, and between-day variance components.

^b Contains within-run, within-day, between-day, and between-lot variance components.

Whole Blood (mmol/mol)

Sample	Instrument	n	Mean (mmol/mol (IFCC))	Within-Run (Repeatability)		Within-Laboratory (Total) ^a		Precision with Additional Component of Between-Lot ^b	
				SD	%CV	SD	%CV	SD	%CV
Control 1	1	240	32.15	0.127	0.4	0.534	1.7	1.009	3.1
	2	240	32.40	0.195	0.6	0.546	1.7	0.959	3.0
	3	240	32.38	0.169	0.5	0.537	1.7	1.089	3.4
Control 2	1	240	77.17	0.153	0.2	0.627	0.8	1.332	1.7
	2	240	77.64	0.241	0.3	0.639	0.8	1.344	1.7
	3	240	77.80	0.182	0.2	0.582	0.7	1.204	1.5
Control 3	1	240	124.49	0.197	0.2	0.640	0.5	1.641	1.3
	2	240	125.22	0.373	0.3	0.734	0.6	1.693	1.4
	3	240	125.52	0.297	0.2	0.605	0.5	1.220	1.0
Panel near 20.0 mmol/mol	1	240	24.81	0.195	0.8	0.457	1.8	0.745	3.0
	2	240	24.90	0.146	0.6	0.428	1.7	0.682	2.7
	3	240	24.93	0.107	0.4	0.364	1.5	0.781	3.1
Panel near 47.5 mmol/mol	1	240	43.82	0.152	0.3	0.400	0.9	0.811	1.9
	2	240	43.93	0.214	0.5	0.422	1.0	0.804	1.8
	3	240	43.97	0.192	0.4	0.382	0.9	0.891	2.0
Panel near 63.9 mmol/mol	1	240	62.94	0.145	0.2	0.396	0.6	0.886	1.4
	2	240	63.07	0.238	0.4	0.456	0.7	0.905	1.4
	3	240	63.13	0.175	0.3	0.385	0.6	0.959	1.5

^a Within-Laboratory variability contains within-run, between-run, and between-day variance components.

^b Contains within-run, within-day, between-day, and between-lot variance components.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.²³ Testing was conducted using 3 lots of the Alinity c Hemoglobin A1c Reagent Kit on each of 2 instruments over a minimum of 3 days. The maximum observed Limit of Blank (LoB) and Limit of Detection (LoD) are summarized below.

Whole Blood

	Glycated Hemoglobin (HbA1c)	Total Hemoglobin (THb)
	μmol/L	μmol/L
LoB ^a	43.8095	0.0000
LoD ^b	48.1568	6.7745

^a The LoB represents the 95th percentile from n ≥ 60 replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on n ≥ 60 replicates of low-analyte level samples.

Hemolysate

	Glycated Hemoglobin (HbA1c)	Total Hemoglobin (THb)
	µmol/L	µmol/L
LoB ^a	39.0741	12.2147
LoD ^b	41.0713	16.3614

^a The LoB represents the 95th percentile from n ≥ 60 replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on n ≥ 60 replicates of low-analyte level samples.

Linearity

A study was performed based on guidance from CLSI EP06-A.²⁴

This assay is linear across the measuring interval of 4.0 to 14.0 %HbA1c (20.22 to 129.51 mmol/mol HbA1c).

Analytical Specificity

Hemoglobin Derivatives

This study was performed on the ARCHITECT c System.

A specificity study was performed based on guidance from CLSI EP07-A2.²⁵ Specificity was assessed by comparing test samples containing the potential interferents listed below to reference samples.

No interference was observed for:

- Acetylated Hemoglobin with ≥ 50 mg/dL of acetylsalicylic acid (ASA) (aspirin)
- Carbamylated Hemoglobin with ≥ 10 mmol/L of Cyanate
- Labile Hemoglobin with ≥ 1000 mg/dL of Glucose

IFCC

The Alinity c Hemoglobin A1c assay had a difference within ± 7% for samples with concentrations ≥ 38.78 mmol/mol HbA1c.

NGSP

The Alinity c Hemoglobin A1c assay had a difference within ± 5% for samples with concentrations ≥ 5.7 %HbA1c.

Hemoglobin Variants

A study was performed based on guidance from CLSI documents EP07-A2.²⁵

Specificity was assessed by comparing samples containing abnormal hemoglobins to a comparator assay. Heterozygous hemoglobin variants (HbS, HbC, HbD, HbE, HbA2) do not interfere with the Alinity c Hemoglobin A1c assay.

IFCC

The Alinity c Hemoglobin A1c assay had a difference within ± 7%.

NGSP

The Alinity c Hemoglobin A1c assay had a difference within ± 5%.

Heterozygous Hemoglobin Variant	Relative % Difference from Reference Concentration			
	~ 6.0 %HbA1c (5.5 to 6.5 %HbA1c)		~9.0 %HbA1c (8.0 to 10.0 %HbA1c) ^a	
	Relative % Difference	Range ^b	Relative % Difference	Range ^b
HbC	-4.2	-8.9, -1.6	-2.5	-6.2, 0.0
HbD	1.9	-1.7, 4.8	2.3	1.1, 3.5
HbE	1.4	-3.2, 6.6	-0.9	-4.1, 2.5
HbS	-0.3	-5.5, 1.8	-1.9	-4.9, 1.1
HbA2	-0.8	-4.9, 3.6	-0.6	-3.7, 5.4
HbF	Difference exceeds -5% when the amount of HbF in the sample exceeds ~5%. ^c			

See the Heterozygous Hemoglobin Variant HbF interference in the following table.

^a The HbA2 results at ~9.0 %HbA1c consisted of samples between 7.1 to 14.0 %HbA1c.

^b The range is defined as the minimum and maximum relative % difference at each concentration level (~6.0 and ~9.0 %HbA1c).

^c Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

NOTE: The presence of multiple variants in a sample may impact the % difference.

Heterozygous Hemoglobin Variant HbF interference is summarized in the following table.

A negative % difference with HbF is proportional in magnitude to the % HbF present in the sample. For example, when the amount of HbF in the sample was 20.4%, the % difference was -21.3% on the Alinity c analyzer.

Heterozygous Hemoglobin Variant	% HbF	Individual % Interference from Comparator
	3.2	-3.1
	4.6	-4.9
HbF	6.2	-4.6
	8.6	-7.4
	11	-9.8

Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information regarding the specific concentrations of HbF that were found to interfere with the Alinity c Hemoglobin A1c assay.

Interference

This study was performed on the ARCHITECT c8000 and c4000 Systems.

Potentially Interfering Substances

A study was performed based on guidance from CLSI EP07-A2.²⁵ Interference effects were assessed by comparing test samples containing the potential interferents listed below to reference samples.

IFCC

The Alinity c Hemoglobin A1c assay had a difference within ± 7% for samples with concentrations ≥ 38.78 mmol/mol HbA1c.

NGSP

The Alinity c Hemoglobin A1c assay had a difference within ± 5% for samples with concentrations ≥ 5.7 %HbA1c. The data in %HbA1c (NGSP) showing the highest level of interference observed on the ARCHITECT c8000 and c4000 Systems is presented in the table.

Potentially Interfering Substance	Interferent Concentration		% Interference ^a	
	Conventional Units	SI Units	6.0 to 7.0 %HbA1c	≥ 8.0 %HbA1c
Ascorbic Acid	3.0 mg/dL	0.15 mmol/L	0.0	0.0
Bilirubin (Conjugated) ^b	15.0 mg/dL	180 µmol/L	-3.2	-3.3
Bilirubin (Unconjugated) ^c	10.0 mg/dL	171 µmol/L	-3.0	-2.7
Glucose	1000 mg/dL	55.5 mmol/L	0.0	0.0
Rheumatoid Factor	200 IU/mL	200 IU/mL	0.0	0.0
Triglycerides	3000 mg/dL	33.9 mmol/L	-1.6	-4.5
Total Protein	22 g/dL ^d	220 g/L	0.0	-1.1
Urea	667 mg/dL	111.06 mmol/L	0.0	0.0
Vitamin E	8.6 mg/dL	200 µmol/L	1.6	0.0

^a

$$\% \text{ Interference} = \frac{\text{Test Result} - \text{Control Result}}{\text{Control Result}} \times 100$$

^b Samples containing conjugated bilirubin at > 15.0 mg/dL (180 µmol/L) demonstrated interference. Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

^c Samples containing unconjugated bilirubin at > 10.0 mg/dL (171 µmol/L) demonstrated interference. Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

^d The total protein concentration of 22 g/dL includes serum protein as well as hemoglobin.

Potentially Interfering Substance	Interferent Concentration		% Interference ^a	
	Conventional Units	SI Units	6.0 to 7.0 %HbA1c	≥ 8.0 %HbA1c
Acarbose	50 mg/dL	0.77 mmol/L	0.0	0.0
Acetaminophen	200 mg/L	1324.5 µmol/L	-0.59	-0.80
N-acetyl-4-benzoquinone imine	20 mg/L	134.2 µmol/L	-1.5	-0.8
N-acetyl-L-cysteine	1600 mg/L	9816 µmol/L	-3.4	-2.7
Acetylsalicylate	50.8 mg/dL	2.82 mmol/L	0.0	0.0
Atorvastatin	0.06 mg/dL	600 µg Eq/L	1.6	0.0
Captopril	0.5 mg/dL	23 µmol/L	-1.5	-1.1
Chlorpropamide	74.7 mg/dL	2.7 mmol/L	0.0	-1.1
Cyanate	50 mg/dL	6.16 mmol/L	0.0	1.1
Dipyrrone	100 mg/L	300.3 µmol/L	0.2	0.1
Furosemide	6.0 mg/dL	181 µmol/L	0.0	1.1
Gemfibrozil	7.5 mg/dL	300 µmol/L	0.0	0.0
Ibuprofen	50 mg/dL	2425 µmol/L	0.0	1.1
Insulin	450 micro units per mL	450 micro units per mL	0.8	0.0
Losartan	5 mg/dL	0.11 mmol/L	0.0	0.0
Metformin	5.1 mg/dL	310 µmol/L	0.0	0.0
Nicotinic Acid	61 mg/dL	4.95 mmol/L	-1.5	-0.5
Propranolol	0.2 mg/dL	7.71 µmol/L	0.0	-0.5
Repaglinide	0.006 mg/dL	132.57 nmol/L	0.8	0.0
4-acetamido antipyrine	40 mg/L	163.3 µmol/L	-0.4	-0.4
4-aminoantipyrine	40 mg/L	197.0 µmol/L	-0.4	-0.2
4-formylamino antipyrine	40 mg/L	173.2 µmol/L	-0.1	-0.6
4-methylamino antipyrine	40 mg/L	184.3 µmol/L	-0.3	-0.4

a

$$\% \text{ Interference} = \frac{\text{Test Result} - \text{Control Result}}{\text{Control Result}} \times 100$$

Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert for further information.

Method Comparison

A study was performed based on guidance from CLSI EP09-A3 using the Passing-Bablok regression method.²⁶

	Units	n	Correlation			Concentration Range
			Coefficient	Intercept	Slope	
Alinity c Hemoglobin A1c vs NGSP Secondary	Whole Blood %HbA1c (NGSP)	134	0.99	-0.10	1.00	4.4 - 14.0
Alinity c Hemoglobin A1c vs ARCHITECT Hemoglobin A1c	Hemolysate %HbA1c (NGSP)	133	1.00	0.10	1.00	4.5 - 13.1
	mmol/mol (IFCC)	130	1.00	0.10	1.01	25.72 - 119.99
	Hemolysate %HbA1c (NGSP)	133	1.00	0.00	1.00	4.7 - 13.1
	mmol/mol (IFCC)	133	1.00	-0.39	1.02	28.22 - 119.44

Predicted Bias

IFCC

The Alinity c Hemoglobin A1c assay is designed to have a predicted bias of ≤ 5% at 42.06, 47.53, and 53.00 mmol/mol HbA1c using Deming regression.

NGSP


The Alinity c Hemoglobin A1c assay is designed to have a predicted bias of ≤ 3% at 6.0, 6.5, and 7.0 %HbA1c using Deming regression. The predicted bias in %HbA1c (NGSP) ranged from -2.0% to -1.5%.

Allowable Total Difference (ATD) Zone

The Alinity c Hemoglobin A1c assay is designed to have > 95% of observations in the ATD zone and the low limit of the two-sided 95% Confidence Interval (CI) > 89.5%. The percentage of observations in the ATD zone was 99.3% (133/134) and the lower limit of the two-sided 95% CI was 95.9%. The ATD zone plots are presented below.

ATD Zone- Alinity c analyzer


Hemoglobin A1c in Whole Blood, mmol/mol (IFCC)

x = NGSP Secondary Reference Laboratory (mmol/mol)
y = Alinity c Hemoglobin A1c (mmol/mol)

ATD Zone- Alinity c analyzer

Hemoglobin A1c in Whole Blood, %HbA1c (NGSP)

x = NGSP Secondary Reference Laboratory (%HbA1c)
y = Alinity c Hemoglobin A1c (%HbA1c)

BIBLIOGRAPHY

- Bry L, Chen PC, Sacks DB. Effects of hemoglobin variants and chemically modified derivatives on assays for glycohemoglobin. *Clin Chem* 2001;47(2):153-163.
- Sacks DB. The diagnosis of diabetes is changing: How implementation of hemoglobin A1c will impact clinical laboratories. *Clin Chem* 2009;55(9):1612-1614.
- American Diabetes Association Workgroup Report: International Expert Committee report on the role of the A1c assay in the diagnosis of diabetes. *Diabetes Care* 2009;32(7):1327-1334.
- Fonseca V, Inzucchi SE, Ferrannini E. Redefining the diagnosis of diabetes using glycated hemoglobin. *Diabetes Care* 2009;32(7):1344-1345.
- American Diabetes Association Executive Summary: Standards of medical care in diabetes - 2010. *Diabetes Care* 2010;33(Suppl 1):S4-S10.
- Sacks DB, Bruns DE, Goldstein DE, et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002;48(3):436-472.
- American Diabetes Association. Position Statement: Standards of Medical Care in Diabetes - 2012. In: *Diabetes Care* 2012;35(Suppl 1):S11-S63.
- Rochman H. Hemoglobin A1c and diabetes mellitus. *Ann Clin Lab Sci* 1980;10(2):111-115.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med* 1993;329(14):977-986.
- Larsen ML, Hørder M, Mogensen EF. Effect of long-term monitoring of glycosylated hemoglobin levels in insulin-dependent diabetes mellitus. *N Engl J Med* 1990;323(15):1021-1025.
- Geistanger A, Arends S, Berding C, et al. Statistical methods for monitoring the relationship between the IFCC reference measurement procedure for hemoglobin A1c and the designated comparison methods in the United States, Japan, and Sweden. *Clin Chem* 2008;54(8):1379-1385.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Goldstein DE, Wiedmeyer HM, England JD, et al. Recent advances in glycosylated hemoglobin measurements. *Crit Rev Clin Lab Sci* 1984;21(3):187-228.
- Peacock I. Glycosylated haemoglobin: measurement and clinical use. *J Clin Pathol* 1984;37:841-851.
- Horton BF, Huisman TH. Studies on the heterogeneity of haemoglobin, VII. Minor haemoglobin components in haematological diseases. *Br J Haematol* 1965;11(3):296-304.
- Lind T, Cheyne GA. Effect of normal pregnancy upon the glycosylated haemoglobins. *Br J Obstet Gynaecol* 1979;86(3):210-213.
- American Diabetes Association. Position Statement: Diagnosis and classification of diabetes mellitus. *Diabetes Care* 2010;33(Suppl 1):S62-S69.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.

- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
REF	List Number
SN	Serial number
Other Symbols	
A1cDIL	Hemoglobin A1c Diluent
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF CANADA	Product of Canada
R1	Reagent 1
R2	Reagent 2
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

Alinity, ARCHITECT and related brand marks are trademarks of Abbott. Other trademarks are the property of their respective owners.

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown
Prince Edward Island
C1E 2B9 Canada

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.corelaboratory.abbott

Revised March 2022.

©2017, 2022 Abbott Laboratories

TECHNOPATH CLINICAL DIAGNOSTICS

MULTICHEM A1c

REF 04V0610

(HB001X)

(HB002X)

Bi-Level

Level: 1

Level: 2

2 x 6 x 1mL

6 x 1mL

6 x 1mL

LOT A2710220

Level 1: A2710221

Level 2: A2710222

 2024-10-31

IVD

INSTRUCTIONS FOR USE

www.technopathclinicaldiagnostics.com

INTENDED USE

Multichem A1c control is intended for use as an assayed quality control material to monitor the precision of laboratory testing procedures for the analyte, HbA1c as listed in the package insert.

SUMMARY AND PRINCIPLE

The use of quality control material is indicated as an objective assessment of the precision of methods and techniques and is an integral part of good laboratory practices. Two levels of control are available to allow performance monitoring within the clinical range.

COMPOSITION

This product is prepared from human whole blood to which stabilizers have been added. The control is used in liquid form for convenience.

PRECAUTIONS

Warning: Biological source material. Treat as potentially infectious.

1. For In Vitro Diagnostic Use.
2. Each blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV, and non-reactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases.
3. This product may also contain other human source material for which there is no approved test. It is recommended that this material be handled in accordance with the OSHA Standard on Bloodborne pathogens¹, Biosafety Level 2² or other appropriate biosafety practices³.
4. Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.
5. Material Safety Data Sheet (MSDS) is available for professional user on request.

PROCEDURE FOR USE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit or reagent being used. Frozen Control should be thawed at 2 to 8°C. Mix by gentle inversion several times prior to use. Remove the cap before use. Replace the cap after use and return to 2 to 8°C storage immediately. Prior to subsequent use, remove from 2 to 8°C storage and mix by gentle inversion several times before use. Do not shake vigorously.

In the event of damage to packaging, contact Abbott Customer Service.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened between -20°C and -80°C. For optimum performance, avoid storage in a frost-free freezer or in a freezer which is opened regularly for other applications such as storage of clinical samples. Once the control material is thawed and opened, it should be stored tightly capped at 2 to 8°C. HbA1c is stable for 30 days under these conditions.

LIMITATIONS

1. Do not use the product past the expiration date.
2. If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
3. This product is not intended for use as a standard or calibrator.
4. If aliquoting this product, please follow the thawing and mixing procedure as outlined in the 'Procedure for Use'. Aliquot for single use and freeze immediately and use within 30 days.

ASSIGNMENT OF VALUES

The values provided in the data sheet were derived from replicate analyses and are specific for a particular lot of product. These values have been generated using third party manufacturers' instrument systems and are specific to one measurement procedure. Technopath make no accuracy claims regarding these values. Tests were performed by the control manufacturer and/or by independent laboratories, for various methods and instrument systems. As a tool to assist in establishing their own mean, laboratories can import the values into their instrument specific system. For more details and to register for access to this file please visit www.technopathcd.com.

Values are provided only as guidelines, each laboratory should establish its own statistical limits. Laboratory means may vary from the values listed during the shelf life of the control. Technopath monitors the values over the shelf life of the control and provides update(s) at www.technopathcd.com or contact your local Abbott customer service representative.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organisation. Laboratory Biosafety Manual. 3rd ed. Geneva: World Health Organisation; 2004.

Customer Service: For questions, including suspected shipping damage, contact your local area Abbott Customer Service representative or find country specific contact information at www.abbottdiagnostics.com.

GLOSSARY

ANALYTES (Value Assigned)

HbA1c

TERMS

Mean

Range

Units

Level

Data Sheet

Bi-Level

ALINITY/ARCHITECT Chemistry Systems			LEVEL 1		LEVEL 2		LEVEL 1		LEVEL 2			
ANALYTES		Arch	Alinity	Mean	Range	Mean	Range	Mean	Range	Mean	Range	
HbA1c			{4P52}	{08P43}	NGSP UNITS		IFCC UNITS					
			(%)	5.50	4.95-6.05	9.95	8.96-10.9	mmol/mol	36.6	30.6-42.7	85.3	74.4-96.2

Footnotes

1. ARCHITECT and Alinity are a trademark of Abbott Laboratories in various jurisdictions.
 2. Some of the assays listed may not be available in all countries.
 3. Value assignment is performed using the ARCHITECT list numbers which have been shown to be equivalent to the Alinity list numbers and are considered family members.
- Values are not provided.

 Techno-path Manufacturing Ltd, Fort Henry Business Park, Ballina, Co. Tipperary, Ireland.
Ph: +353 (0) 61 525700 Web: www.technopathclinicaldiagnostics.com

Distributed by
Abbott Laboratories Inc, Abbott Park, IL 60064 USA and ABBOTT 65205 Wiesbaden, Germany

CONTROL	REF			LOT	IVD			
Control	Catalog Number	Manufacturer	Use By (yyyy-MM-DD)	Lot Number	In Vitro Diagnostic Medical Device	Consult Instructions for Use	Biological Risks	Temperature Limitation
SN	INFORMATION FOR USA ONLY			DISTRIBUTED BY		PRODUCT OF IRELAND		GTIN
Serial Number	Information needed for United States of America only			Distributed By		Product of Ireland		Global Trade Item Number

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BIBLIOGRAFIE

- McPherson RA, Pincus MR, editors. *Henry's Clinical Diagnosis and Management by Laboratory Methods*. 23rd ed. St. Louis, Missouri: Elsevier; 2017.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 6th ed. Washington, DC: US Government Printing Office; June 2020.
- World Health Organization. *Laboratory Biosafety Manual*. 4th ed. Geneva: World Health Organization; 2020.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Rifai N, Horvath AR, Wittwer CT, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. Elsevier; 2018:1439.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking*. 1st ed. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Elfatih A, Anderson NR, Fahie-Wilson MN, et al. Pseudo-pseudohypercalcaemia, apparent primary hyperparathyroidism and Waldenström's macroglobulinaemia. *J Clin Pathol* 2007;60:436-437. doi: 10.1136/jcp.2006.036756
- Adeli K, Ceriotti F, Nieuwesteeg M. Reference information for the clinical laboratory. In: Rifai N, Horvath AR, Wittwer CT, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. Elsevier; 2018:1745-1818.
- Bourdel-Marchasson I., Laksir H., Puget E. Interpreting routine biochemistry in those aged over 65 years: a time for change. *Maturitas* 2010;66(1):39-45. doi.org/10.1016/j.maturitas.2010.02.004
- Wu AHB, editor. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Elsevier Saunders; 2006.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures*. 2nd ed. CLSI Guideline EP06. Wayne, PA: CLSI; 2020.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3rd ed. CLSI Guideline EP09c. Wayne, PA: CLSI; 2018.

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REFERINȚE

1. Friedman RB, Young DS. *Effects of Disease on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press, 1989.
2. Thomas L. Immunochemical technics. In: Thomas L, editor. *Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results*, 1st ed. Frankfurt/Main, Germany: TH-Books Verlagsgesellschaft; 1998:1427–41.
3. Simó JM, Joven J, Clivillé X, et al. Automated latex agglutination immunoassay of serum ferritin with a centrifugal analyzer. *Clin Chem* 40;1994:625–9.
4. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne pathogens.
5. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
6. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
7. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
8. Westgard JO, Barry PL. *Cost-effective Quality Control: Managing the Quality and Productivity of Analytical Processes*. Washington, DC: AACC Press, 1986.
9. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000.

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Configure assay parameters — General				
<input checked="" type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: FER	Type: Photometric	Version: †		
Number: 2906	Run controls for onboard reagents by:†† Lot			
<input checked="" type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks		
Reaction mode: End up	Primary Secondary	Read times		
Wavelength: 572 /	None	Main: 32 – 33		
Last required read: 33				
Absorbance range: ___ - ___	Color correction: ___ - ___			
Sample blank type: Self	Blank: 19 – 20			

<input type="radio"/> Reaction definition	<input checked="" type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks	
Reagent: FEROB	Reagent volume: 160	R1	R2
Diluent: Saline	Water volume: ___	50	50
Diluent dispense mode: Type 0	Dispense mode: Type 0	Type 0	Type 2
Dilution name	Sample	Diluted sample	Dilution factor
STANDARD	: 25.0	___	= 1:1.00
Dil 1	: 20.0	15.0	= 1:9.57
___	: ___	___	= ___
			Default dilution
			<input checked="" type="radio"/>
			<input type="radio"/>
			<input type="radio"/>

<input type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input checked="" type="radio"/> Validity checks	
Reaction check: Rate Ratio		A	B
	Read time: 24 – 26	18 – 20	
	Calculation limits: 0.5000 – 9.0000		
	Minimum: 0.0750		
Maximum absorbance variation: ___			

Configure assay parameters — Calibration				
<input type="radio"/> General	<input checked="" type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: FER	Calibration method: Spline			
<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks	
Calibrator set: FER	Blank: Water	Calibrator level: 0.0†	Concentration: ___	
	Cal 1: FER1	††		
Replicates: 2 [Range 1 – 3]	Cal 2: FER2	††		
	Cal 3: FER3	††		
	Cal 4: FER4	††		

<input type="radio"/> Calibrators	<input checked="" type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks		
Calibrator: FER	Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	25.0	___	___	___	___
Cal 1: FER1	25.0	___	___	___	___
Cal 2: FER2	25.0	___	___	___	___
Cal 3: FER3	25.0	___	___	___	___
Cal 4: FER4	25.0	___	___	___	___

<input type="radio"/> Calibrators	<input type="radio"/> Volumes	<input checked="" type="radio"/> Intervals	<input type="radio"/> Validity checks	
Calibration intervals:				
Full interval: 720 (hours)				
Calibration type:				
Adjust type: None				

<input type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input checked="" type="radio"/> Validity checks	
Blank absorbance range: ___ - ___				
Span: Blank - Blank				
Span absorbance range: ___ - ___				
Expected cal factor: 0.00				
Expected cal factor tolerance %: 0				

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Configure assay parameters — Results				
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results	<input type="radio"/> Interpretation
Assay: FER		Assay number: 2906		
Dilution default range:		Result units: ng/mL		
		Low-Linearity: 10.0		
		High-Linearity: 500.0		
Gender and age specific ranges:*				
GENDER	AGE (UNITS)	NORMAL	EXTREME	

Configure result units	
Assay: FER	Version: †
Result units: ng/mL	
Decimal places: 1	[Range 0 – 4]
Correlation factor: 1.0000	
Intercept: 0.0000	

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Configure assay parameters — Results				
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results	<input type="radio"/> Interpretation
Assay: FER		Assay number: 2906		
Dilution default range:		Result units: µg/L		
		Low-Linearity: 10.0		
		High-Linearity: 500.0		
Gender and age specific ranges:*				
GENDER	AGE (UNITS)	NORMAL	EXTREME	

Configure result units	
Assay: FER	Version: †
Result units: µg/L	
Decimal places: 1	[Range 0 – 4]
Correlation factor: 1.0000	
Intercept: 0.0000	

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ARCHITECT

INTENDED USE

The Quantia Ferritin Standard is intended for use in establishing the calibration for the Quantia Ferritin assay by turbidimetry. This standard is intended to be used together with the Quantia Ferritin reagents. Carefully read the instructions with the corresponding reagents before use.

PRINCIPLE

When the standards, which contains ferritin, are mixed with the Quantia Ferritin reagent, a clear agglutination occurs which can be measured by turbidimetry.

CONTENTS

REF 6K49-03 Quantia Ferritin Standard, CAL 1-4 1 x 1 mL

This product is of human origin and is standardized to the WHO International Standard. The concentrations in ng/mL are indicated on the vial labels. Contains sodium azide (< 0.1%).

NOTE: Substituting components from other manufacturers may affect the results of the assay.

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Do not mix fresh reagents with in-use reagents.
-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these standards and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents. Donor units of human sourced material have been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, and anti-HCV.
- The following warnings and precautions apply to Quantia Ferritin Standard:
 - Contains sodium azide.
 - EUH032 Contact with acids liberates very toxic gas.
 - P501 Dispose of contents / container in accordance with local regulations.
- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.
- For a detailed discussion of safety precautions during system operation, refer to the **ARCHITECT System Operations Manual, Section 8.**

STORAGE

Unopened standards are stable until the expiration date shown on the label if stored at 2 to 8°C.

All Quantia products contain preservatives, but are still susceptible to contamination. Handle with normal precautions.

PREPARATION

The standards are ready-to-use.

QUALITY CONTROL

After calibration, the control values should be within the established range. Otherwise, corrective measures should be taken by the user.

REFERENCES

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

CAL 1	Calibrator 1
CAL 1-4	Calibrators 1 through 4
CAL 2	Calibrator 2
CAL 3	Calibrator 3
CAL 4	Calibrator 4
CONTAINS: AZIDE	Contains sodium azide. Contact with acids liberates very toxic gas.
DISTRIBUTED BY	Distributed by
FOR USE WITH	Identifies products to be used together
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Batch code/Lot number
PRODUCT OF SPAIN	Product of Spain
REF	Catalog number/List number
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date

QUANTIA FERRITIN STANDARD

REF 6K49-03

  Biokit, S.A.
Av. Can Montcau 7
08186 Lliçà d'Amunt
Barcelona, Spain

DISTRIBUTED BY

ABBOTT
65205 Wiesbaden, Germany

 **Abbott**

September 2017

QUANTIA FERRITIN/MYOGLOBIN/IgE CONTROL

October 2017

INTENDED USE

The Quantia Ferritin/Myoglobin/IgE Control is intended for use in monitoring the quality control of results obtained with the Quantia (Ferritin, Myoglobin, and IgE) assays by turbidimetry. This control is intended to be used together with the Quantia (Ferritin, Myoglobin, and IgE) reagents. Carefully read the instructions with the corresponding reagents before use.

PRINCIPLE

When the controls, which contain ferritin, myoglobin, and IgE are mixed with the corresponding Quantia reagents, a clear agglutination occurs which can be measured by turbidimetry.

CONTENTS

REF 6K56-02 Quantia Ferritin/Myoglobin/IgE Control:

CONTROL I	low level	(lyophilized)	1 x 3 mL
CONTROL II	high level	(lyophilized)	1 x 3 mL

This product is of human origin. Refer to the control data sheet for assayed values. Contains sodium azide (0.2% lyophilized, <0.01% reconstituted).

NOTE: Substituting components from other manufacturers may affect the results of the assay.

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Do not mix fresh reagents with in-use reagents.
-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these controls and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.
Donor units of human-sourced materials used in the controls have been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, and anti-HCV.
- The following warnings and precautions apply to **CONTROL I** and **CONTROL II**:
Contains sodium azide.
EUH032 Contact with acids liberates very toxic gas.
P501 Dispose of contents / container in accordance with local regulations.
- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.
- For a detailed discussion of safety precautions during system operation, refer to the **ARCHITECT System Operations Manual, Section 8**.

en

Quantia Ferritin/Myoglobin/IgE Control

REF 6K56-02

G10834/R01

C6K5T0

FOR USE WITH

ARCHITECT

STORAGE

Quantia Ferritin/Myoglobin/IgE Control is stable until the expiration date shown on the label if stored at 2 to 8°C.

All Quantia products contain preservatives, but are still susceptible to contamination. Handle with normal precautions.

Once reconstituted, the controls can be used for 7 days if stored at 2 to 8°C.

PREPARATION

The controls should be reconstituted with 3 mL of distilled water. Allow reconstituted material to stand for 5 minutes. Prior to use, gently swirl the vial.

QUALITY CONTROL

The control values should be within the established range. Otherwise, corrective measures should be taken by the user.

REFERENCES

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

CONTAINS: AZIDE	Contains sodium azide. Contact with acids liberates very toxic gas.
CONTROL	Control
CONTROL I	Control I
CONTROL II	Control II
DISTRIBUTED BY	Distributed by
FOR USE WITH	Identifies products to be used together
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Batch code/Lot number
MEAN	Mean Value
PRODUCT OF JAPAN	Product of Japan
RANGE	Range
REAGENT	Reagent
REF	Catalog number/List number
SI UNITS	SI Units
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date

QUANTIA FERRITIN/MYOGLOBIN/IgE CONTROL

REF 6K56-02

  Biokit, S.A.
Av. Can Montcau 7
08186 Lliçà d'Amunt
Barcelona, Spain

DISTRIBUTED BY

ABBOTT
65205 Wiesbaden, Germany

 **Abbott**

October 2017

Multichem® S Plus

Replace up to four competitor chemistry QC products with Multichem S Plus

Features

- Human-based matrix proprietary formulation.
- Mimics performance of patient samples.
- Targeted at key clinical decision points.
- Liquid for ease of use.

Specifications

- 36 month closed vial stability at -20°C to -80°C.
- 10 day open vial stability at 2°C to 8°C.

Instruments

Multichem S Plus, independent clinical chemistry QC, is suitable for use on all major platforms.

Order Information		
Multichem S Plus	Config.	Part Code
Assayed Single Level (Level 1)	15 x 10mL	CH101PLA
Assayed Single Level (Level 2)	15 x 10mL	CH102PLA
Assayed Single Level (Level 3)	15 x 10mL	CH103PLA
Unassayed Single Level (Level 1)	15 x 10mL	CH101CRP
Unassayed Single Level (Level 2)	15 x 10mL	CH102CRP
Unassayed Single Level (Level 3)	15 x 10mL	CH103CRP
Mini kit (contains one vial of each level)	3 x 1 x 10mL	CH100MCRP

Analyte List

Chemistry

- Albumin
- Bilirubin, Direct
- Bilirubin, Total
- Calcium
- Carbon Dioxide (Bicarbonate)
- Chloride
- Creatinine
- Glucose
- Iron
- Lactic acid
- Magnesium
- Phosphorous
- Potassium
- Protein, Total
- Sodium
- Total Iron Binding Capacity*
- Unsaturated Iron Binding Capacity (UIBC)
- Urea
- Uric Acid

Esoterics

- ACE*
- Bile Acids
- Bili, Indirect*
- Caffeine*
- Calcium, Ionized*
- Copper*
- Cortisol
- Ethanol
- Fructosamine*
- NT-proBNP*
- Osmolality*
- Protein Electrophoresis*
- Triiodothyronine*
- Thyroxine, (Total T4)
- Troponin T*
- Zinc*

Immunoproteins

- Alpha-1 Aminoglycoprotein
- Alpha-1 Antitrypsin
- Alpha 2-macroglobulin*
- ASO*
- ADNase B (Anti-Streptococcal DNase B)*
- Antithrombin III*
- Apolipoprotein A1
- Apolipoprotein B
- Beta-2 Microglobulin
- C1 Inhibitor*
- CH50 (Total hemolytic Complement)*
- Cystatin C*
- Complement C3
- Complement C4
- Ceruloplasmin
- C-Reactive Protein
- Ferritin*
- Haptoglobin
- Hemopexin*
- Immunoglobulin A
- Immunoglobulin G
- Immunoglobulin M
- IgE*
- IgG1, Subclass*
- IgG2, Subclass*
- IgG3, Subclass*
- IgG4, Subclass*
- Kappa Light Chain*
- Lambda Light Chain*
- Lipoprotein (a)*
- Prealbumin
- Properdin Factor B*
- Retinol Binding Protein*
- Rheumatoid Factor*
- Transferrin
- sTfR (Soluble Transferrin Receptor)*

Enzymes

- Acid Phosphatase
- Alanine Aminotransferase (ALT)
- Alkaline Phosphatase (ALP)
- Amylase (Pancreatic)
- Amylase (Total)
- Aspartate Aminotransferase (AST)
- Alpha-hydroxybutyrate Dehydrogenase (α-HBDH)*
- Beta-hydroxybutyrate Dehydrogenase*
- Cholinesterase
- Creatine Phosphokinase
- CK-MB*
- Gamma Glutamyltransferase
- Lactate Dehydrogenase (L to P)
- Lipase
- Prostatic Acid Phosphatase*

Lipids

- Cholesterol, HDL
- Cholesterol, LDL
- Cholesterol, Total
- Phospholipids*
- Triglycerides

Therapeutic Drugs

- Acetaminophen
- Amikacin
- Carbamazepine
- Digoxin
- Gentamicin
- Lithium
- Phenobarbital
- Phenytoin
- Salicylic Acid
- Theophylline
- Tobramycin
- Valproic Acid
- Vancomycin

*Please refer to lot specific package inserts for stability and performance claims.

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- 04R1001 Alinity ci-series Calibrator/Control Replacement Caps

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BIBLIOGRAFIE

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.

4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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