

CT Clinical Image Gallery

Aquilion ONE

GENESIS Edition

Transforming CT

Transforming Clinical Confidence



Table of contents

| | |
|---|-------|
| TIA Follow-up | 4-5 |
| Routine Brain | 6-7 |
| SURESubtraction™ Lung | 8-9 |
| Low Dose Cardiac - FIRST | 10 |
| Cardiac - FIRST | 11 |
| Cardiac Pacemaker SEMAR™ | 12 |
| Pediatric Cardiac - FIRST | 13 |
| Liver Metastases | 14-15 |
| Abdominal Aortic Aneurysm - FIRST | 16-17 |
| Motor Vehicle Accident | 18-19 |
| Wrist CT with the ease of X-ray - FIRST | 20-21 |
| Pediatric Brain in ONE Volume | 22-23 |


We gratefully acknowledge the following institutions that contributed to this brochure: Leiden University Medical Center, the Netherlands




TIA Follow-up

This 61-year-old woman was admitted to the hospital following a TIA. A CT scan was performed to rule out stroke. The scan automatically reconstructed with AIDR 3D Enhanced (Adaptive Iterative Dose Reduction) shows excellent distinction between cortical grey matter and white matter.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|-------------|--------|-----|-----|-------------------|-----------------|-------------------|---------------------------|--------------|----------------------|--------|
| Helical | — | 0.5 mm x 40 | Detail | 120 | 240 | 0.75 | 240 | AIDR* 3D Enhanced | 59.1 | 948.5 | 1.99 | 0.0021 |





* Adaptive Iterative Dose Reduction



Routine Brain

This brain scan of a 26-year-old woman permits clear identification of the gyral pattern, basal cistern, and foramen magnum. No abnormalities are seen.



| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|-------------|--------|-----|-----|-------------------|-----------------|------------------|---------------------------|--------------|----------------------|--------|
| Helical | — | 0.5 mm x 40 | Detail | 120 | 240 | 0.75 | 240 | AIDR 3D Enhanced | 59.1 | 948.5 | 1.99 | 0.0021 |



SURESubtraction Lung

This 75-year-old man with a long history of smoking was sent for a CT scan. Multiple tumors are seen in both lungs, which also show evidence of COPD. SURESubtraction Lung was performed to provide additional blood flow information.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|---------------|----------|-------------|----------|-----|------------------|-------------------|-----------------|------------------|---------------------------|--------------|----------------------|-------|
| Ultra-Helical | CE | 0.5 mm x 80 | Standard | 100 | SUREExposure™ 3D | 0.275 | 295 | AIDR 3D Enhanced | 0.8 | 30.9 | 0.4 | 0.014 |
| Ultra-Helical | CE | 0.5 mm x 80 | Standard | 100 | SUREExposure 3D | 0.275 | 295 | AIDR 3D Enhanced | 2.7 | 97.7 | 1.3 | 0.014 |




Low Dose Cardiac



This 68-year-old woman with a heart rate of 88 bpm presented with atypical chest pain. A one beat prospective scan was performed for a low radiation dose of just 0.36 mSv. The coronary arteries are normal.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|--------------|-------|-----------|-----------------|-------------------|-----------------|----------------|---------------------------|--------------|----------------------|-------|
| Volume | CE | 0.5 mm x 260 | - | SUREkV 80 | SUREExposure 3D | 0.275 | 130 | FIRST | 2.0 | 25.6 | 0.36 | 0.014 |



Cardiac



This 68-year-old man with a BMI 28.7 presented with atypical chest pain. This scan was performed with FIRST (Forward projected model based Iterative Reconstruction SoluTioN) ensuring automatic dose reduction for the prescribed level of image quality.


| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|--------------|-------|------------|-----------------|-------------------|-----------------|----------------|---------------------------|--------------|----------------------|-------|
| Volume | CE | 0.5 mm x 280 | — | SUREkV 100 | SUREExposure 3D | 0.275 | 140 | FIRST | 6.3 | 88.8 | 1.24 | 0.014 |



Cardiac Pacemaker SEMAR

This 65-year-old man with a heart rate of 83 bpm and a BMI of 27.4 presented following a pacemaker insertion. Metallic artifacts from the pacemaker leads are significantly reduced utilizing Single Energy Metal Artifact Reduction (SEMAR) reconstruction technology.


| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|--------------|-------|------------|-----------------|-------------------|-----------------|------------------|---------------------------|--------------|----------------------|-------|
| Volume | CE | 0.5 mm x 280 | — | SUREkV 100 | SUREExposure 3D | 0.275 | 140 | AIDR 3D Enhanced | 5.4 | 74.9 | 1.05 | 0.014 |



Pediatric Cardiac

This 3-month-old baby had suspected congenital defects of the heart. With a heart rate of 131 beats per minute, an ECG gated ultra low dose volume scan was performed in one rotation with a total scan time of just 0.275 s.



| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|--------------|-------|-----------|-----------------|-------------------|-----------------|----------------|---------------------------|--------------|----------------------|-------|
| Volume | CE | 0.5 mm x 320 | — | SUREkV 80 | SUREExposure 3D | 0.275 | 80 | FIRST | 0.7 | 5.8 | 0.22 | 0.039 |



Liver Metastases

This 70-year-old man with liver metastases underwent a follow-up CT scan. An arterial scan of the liver and a portal venous scan of the chest and abdomen were performed. Excellent depiction of the patient's disease is clearly demonstrated.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|---------------|----------|-------------|----------|-----|-----------------|-------------------|-----------------|------------------|---------------------------|--------------|----------------------|--------|
| Ultra Helical | CE | 0.5 mm x 80 | Standard | 120 | SUREExposure 3D | 0.5 | 240 | AIDR 3D Enhanced | 6.3 | 179.7 | 2.69 | 0.014 |
| Ultra Helical | CE | 0.5 mm x 80 | Standard | 120 | SUREExposure 3D | 0.5 | 645 | AIDR 3D Enhanced | 6.6 | 456.1 | 6.61 | 0.0145 |





Abdominal Aortic Aneurysm



This 64-year-old woman presented with an abdominal aortic aneurysm. A CT was requested for stent graft planning.



| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|---------------|----------|-------------|----------|------------|-----------------|-------------------|-----------------|----------------|---------------------------|--------------|----------------------|-------|
| Ultra Helical | CE | 0.5 mm x 80 | Standard | SUREkV 100 | SUREExposure 3D | 0.275 | 444 | FIRST | 4.2 | 209 | 3.13 | 0.015 |



Motor Vehicle Accident

This 45-year-old man presented to the emergency department following a motor vehicle accident. A CT scan of the chest, abdomen and pelvis was requested to evaluate the extent of his injuries. In the chest, multiple rib fractures are seen, with associated pneumothorax and subcutaneous emphysema. The abdominal organs, spine and pelvis are normal.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|---------------|----------|-------------|----------|-----|------------------|-------------------|-----------------|------------------|---------------------------|--------------|----------------------|--------|
| Ultra Helical | CE | 0.5 mm x 80 | Standard | 120 | SURE Exposure 3D | 0.275 | 444 | AIDR 3D Enhanced | 6.3 | 475.2 | 6.89 | 0.0145 |



Wrist CT with the ease of X-ray



This 44-year-old woman presented with fracture of the radius following a fall 3 weeks ago. A follow up CT scan was requested after X-ray to provide more details of the fracture. Laser collimation* was used, eliminating the need for scanograms. The patient was sitting in a comfortable position and the scan was performed in 0.275 s.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|--------------|-------|-----|------------------|-------------------|-----------------|----------------|---------------------------|--------------|----------------------|--------|
| Volume | — | 0.5 mm x 200 | — | 120 | SURE Exposure 3D | 0.275 | 100 | FIRST | 3.5 | 34.8 | 0.02 | 0.0008 |



Volumetric scans of extremities and pediatric examinations can be performed as easily as taking a conventional X-ray. Laser collimation allows the field of view and scan range to be set directly from the touch panel control on the gantry.



*Option



Pediatric Brain in ONE Volume

This 9-year-old boy presented for a CT following a fall. The CT scan was performed in a single rotation with a total scan time of only 0.5 s. The low dose Volume CT images show excellent image quality with clear gray/white matter visualization. The cistern and ventricles have a normal size. No intracranial blood or bone fractures were detected.

| Scan Mode | Contrast | Collimation | Pitch | kVp | mAs | Rotation Time (s) | Scan Range (mm) | Dose Reduction | CTDI _{vol} (mGy) | DLP (mGy-cm) | Effective Dose (mSv) | k |
|-----------|----------|--------------|-------|-----|-----|-------------------|-----------------|------------------|---------------------------|--------------|----------------------|--------|
| Volume | — | 0.5 mm x 320 | — | 120 | 215 | 0.5 | 160 | AIDR 3D Enhanced | 27.0 | 431.7 | 1.38 | 0.0032 |



Disclaimer: Any reference to X-ray exposure is intended as a reference guideline only. The guidelines in this document do not substitute for the judgment of a healthcare provider. Each scan requires medical judgment by the healthcare provider about exposing the patient to ionizing radiation. In clinical practice, the use of the AIDR 3D and FIRST (Forward projected model-based Iterative Reconstruction SoluTion) features may reduce CT patient dose depending on the clinical task, patient size, anatomical location and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task.

Due to local regulatory processes, some of the products included in this brochure may not be available in each country. Please contact your sales representative for the most current information.



CANON MEDICAL SYSTEMS CORPORATION

<https://global.medical.canon>

©Canon Medical Systems Corporation 2016-2018. All rights reserved.
Design and specifications are subject to change without notice.
MCACT0288EAA 2018-02 CMSC/D/Printed in Japan

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485.
Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Aquilion ONE, Aquilion ONE GENESIS, ^{SURE}Subtraction, SEMAR, ^{SURE}Exposure and Made for Life are trademarks of Canon Medical Systems Corporation.

Made For life