



# Flow-c Anesthesia Machine

Compact solution for a streamlined workflow

This document is intended to provide information to an international audience outside of the US.

**GETINGE** 







# A smoother workflow for busy ORs

Crowded workspaces. Packed schedules. Different patients. Complex techniques. These are the challenges you face every operation, every day.

That is why we've developed the Flow-c: a compact anesthesia machine where every detail has been designed to ease your daily work.

With innovation and smart design, Flow-c helps you to create an efficient working environment to ensure the highest standards of care.

Experience the Flow.



# Designed with you

– to get work flowing

Every detail of the Flow-c has been thoughtfully designed in collaboration with clinicians to ease your daily work in the fast-paced OR.

The rear of an anesthesia machine often becomes a tangled mess. But not with the Flow-c. Cables and hoses are routed inside panel arms and behind the back covers.

## Easy to use

The intuitive touch screen gives you one point of control for all functions. The system's simplicity of operation saves time and contributes to safety in busy ORs. Tools are right where you want them, so you can work in an ergonomic and comfortable position.

## Flexible mounting


Mounting capabilities are one area where our focus on design innovation stands out. Despite its compact footprint, Flow-c packs in a greater rail length than other machines on the market. We were determined that its compact form should add flexibility, not limit it, and we have used every millimeter of space to the full. Rails are stepless, so you can personalize your Flow-c, adding monitors, tables and other accessories where it best suits you.

## Cable management without clutter

The back of the Flow-c minimizes the clutter of hoses and cables. These are neatly routed and covered by specially designed panels, contributing to improved hygiene and safety.

## Convenient details

Small features often have a big impact on the experience of using a machine every day. After listening carefully to users, we've introduced some details that will make a difference to your workflow. We know your hands are often full, so we've added practical hooks for convenient placement of tubes and more. You'll also find surprisingly generous table and drawer space, as well as USB ports for data transfer and charging.



Configurable quick access presets for key functions such as agent, oxygen and FGF adjustments.






The interface presents key information at a glance, and requires minimal training – so you can focus on the patient.

The Flow-c has a convenient shelf, giving easy access to the USB ports. The work-space is lit by a dimmable LED lamp.





# Next-generation technology at the core

## – for effective and gentle treatment

The Flow family was created by the engineers behind the world-class Servo ventilator platform. The innovative technology at the heart of your Flow-c ensures superior ventilation performance<sup>1</sup>, and outstanding agent efficiency—reducing the use of costly and harmful anesthetic agents.

### The power to care

Ventilation performance is not only about modes. Most importantly, it's about ensuring the power and precision needed to ventilate all patient categories. Power when you need it, yet gentle on the lungs – that's a key strength of our next-generation Flow core technology.

### Servo gas modules

The Servo gas modules enable ICU-quality ventilation. They deliver up to 200 l/min inspiratory flow and are capable of adjusting pressure and flow several times within every breath, according to each patient's needs.

### Volume Reflector rebreathing system

Our patented Volume Reflector is a smart rebreathing system. In combination with Servo gas modules, it enables accurate tidal volumes down to 5 ml, providing better ventilation performance compared to bellows, turbine and piston-operated systems.<sup>1</sup>

The rigid Volume Reflector is never empty, ensuring uninterrupted ventilation, and compensates effectively for any leakage.<sup>2</sup> The Volume Reflector is oxygen driven, minimizing the risk of hypoxic mixtures. It has a small system volume for fast wash-in and wash-out and a rebreathing fraction of 98%.

### Electronic injection vaporizers

Electronic injection technology enables precise delivery of agents, primarily during the inspiratory phase, with minimal waste. The lightweight and maintenance-free vaporizers can be refilled and exchanged while the machine is running, and do not require annual calibration.




# Low-flow anesthesia the safest<sup>5</sup> way

## Active inspired O<sub>2</sub>Guard protects your patients

O<sub>2</sub>Guard is designed to prevent hypoxia.<sup>4</sup> This unique safety mechanism overrules the clinician's settings and increases the flow of fresh gas and oxygen should the O<sub>2</sub> level drop below 21%. Conventional guards will only trigger an alarm. The O<sub>2</sub>Guard is a standard feature on all Flow models. Learn more at [www.getinge.com/o2guard](http://www.getinge.com/o2guard)

»O<sub>2</sub>Guard is the only commercially active inspired hypoxic guard available.«<sup>5</sup>




## Agent concentration in target organ – MAC Brain guides you

Due to pharmacokinetics there is a time delay in agent concentrations between the lungs and the target organ, the brain. The unique MAC Brain tool visualizes the difference to support better dosing and planning of agent delivery.

## Visual support when lowering the flows

The VRI (Volume Reflector Indicator) is a useful visual guide that enables you to optimize the rebreathing fraction and thus save anesthetic agent. The tool makes it simple to set the optimal Fresh Gas Flow (FGF) and volume ratio. Agent consumption can easily be monitored via the interface.





# See the bigger picture

– ownership with less stress  
and easier flow



## **Minimizing your long-term costs**

The purchase price is just one part of the total cost of owning an anesthesia machine over time, so we have designed the Flow-c to minimize cost of ownership. From an intuitive interface that optimizes workflows and minimizes staff training time, to innovations that significantly reduce consumption of anesthetic agents.

## **Getinge Care: protecting your investment**

Optimizing uptime is an excellent way to boost productivity and reduce costs. A Getinge Care service plan ensures your equipment always performs to its full potential, allowing you to focus on saving lives.

## **Smart fleet management reduces risks**

We make it easy to manage a large fleet of Flow machines. They all share the same user-friendly interface, and they have many components in common, such as the Volume Reflector and vaporizers. Getinge Online gives you fleet overview and can be accessed from any device.

We also offer an extensive range of readily available consumables. These are designed for the highest possible level of patient safety and outstanding ease of use.

## **Extensive training programs**

Keeping skills updated improves patient outcomes, reduces risks and boosts productivity. We tailor training to meet your needs, which includes e-learning in addition to hands-on training courses.



#### **Optimized equipment with Getinge Online**

Getinge Online gives you access to real-time information about your Flow-c machines, e.g. agent consumption, running hours, next preventive maintenance. Through the portal a service technician can resolve most issues remotely, saving time and costs.

#### **Connectivity between devices with MSync**

MSync makes it easy to connect the Flow-c to patient monitor, HIS and patient data management system. Patient data is transferred via HL7 (MSync) in real time to support decision-making.



# Flow-c at a glance

— everyday work that simply flows

## 15" tilttable and rotatable touch screen

Flexible and intuitive to use

## USB ports and LED light

Data transfer and charging easily accessible. Three-step dimmable LED lights up the workspace

## Open platform

Customize to match your needs

## Flexible workspaces

Practical for your everyday needs

## Electronic vaporizer

Precise and cost-efficient

## Smart cable and hose management

Less clutter, easier cleaning and improved mobility

## Central 4-wheel brake and cable pushers

Enhances device mobility

## 215 cm multipurpose rails

Extensive mounting options



# Flow-c Anesthesia Machine

## Technical specifications

Gas volumes, flows and leakages associated with the breathing system are stated in the technical specifications and adhere to BTPS reference conditions. (Body temperature, ambient pressure, Saturated).

All gas concentration readings are normally referenced to dry gas conditions, ambient room temperature and atmospheric pressure (ATPD).

The condition for measured inlet gas pressures and flows is STPD (Standard Temperature and Pressure Dry); 20 °C, standard pressure at 101.3 kPa and 0 % relative humidity (dry).

### Weight and dimensions

Base system weight (out of the box weight)	120 kg (249 lbs)
System nominal weight *	134 kg (295 lbs)
* Equipped with control panel, patient cassette, one full vaporizer, one CO <sub>2</sub> absorber.	
System max weight including maximum load	270 kg (595 lbs)
Dimensions of base plate	697 × 863 mm (27.4" × 34.0")
Drawers	1
Vertical rail	4
Wheels	Four wheels (diameter 125 mm/4.9")
Working surface/writing table	380 mm × 480 mm (15.0" × 18.9")
Reading lamp	Adjustable LED light integrated into the shelf tower.

### Display

Type	LED touch screen, complete with 11 membrane switches and one rotary knob
Size	432 mm × 295 mm (17.0" × 11.6")
Placement	Attached to display arm
Viewing area	15"
Waveforms	Up to 6 waveforms, user configurable
Trends	<ul style="list-style-type: none"><li>• Graphic display, 1 to 24 hour resolution</li><li>• Numeric display, 1 to 60 minute resolution</li></ul>



## Essential performance (term defined in IEC 60601-1)


- Oxygen flow under all conditions except the failure of the oxygen supply or generation of a clinical and/or technical alarm.
- Delivery of a non-hypoxic gas mixture to the patient or generation of a clinical and/or technical alarm.
- Non-delivery of excessive concentrations of a volatile anesthetic agent or generation of a clinical alarm.
- Airway pressure monitoring and associated clinical alarms (Ppeak, PEEP).
- Gas measurement accuracy (for Isoflurane, Desflurane, Sevoflurane, CO<sub>2</sub>, N<sub>2</sub>O, O<sub>2</sub>) and generation of gas measurement associated clinical alarms or generation of a technical alarm.
- Delivery of ventilation at the patient connection port within the alarm limits set by the operator or generation of a clinical or technical alarm.

## Environment

Environment	Operating conditions	Non-operating conditions
Ambient temperature	+15 °C to +35 °C (+60 °F to +95 °F) (Desflurane: +15 °C to +30 °C, +60 °F to +85 °F)	-25 °C to +60 °C (-15 °F to +140 °F)
CO <sub>2</sub> absorber	15 °C to 35 °C (60 °F to 95 °F)	0 °C to 35 °C (32 °F to 95 °F)
Storage at higher temperatures can result in reduced efficiency and service life due to moisture loss. When correctly stored, canisters will maintain absorption capacity for a period of two years.		
Battery	+15 °C to +35 °C (+60 °F to +95 °F)	+5 °C to +40 °C (+40 °F to +104 °F)
When the system is disconnected from a mains power supply, a fully charged battery can be stored in the anesthesia system for up to six weeks at temperatures between +5 °C and +40 °C (+40 °F and +105 °F). At temperatures between +50 °C and +60 °C (+125 °F and +140 °F) storage time is one week. If these limits are exceeded, battery performance can no longer be guaranteed.		
Relative humidity (non-condensing)	15% to 95%	<95%
Atmospheric pressure	700 hPa–1060 hPa	470 hPa–1060 hPa

## Maximum weight, number, and position of accessories

- Accessories must be installed according to any installation and safety guidelines given in the accessories installation instruction. Additional local, regional, and/or national guidelines related to occupational safety may apply.
- The following illustrations show a typical configuration. The setup given in the table has been verified by the manufacturer. The manufacturer assumes no responsibility for other configurations.
- The functionality of the system is extended by installing accessory carriers with appropriate accessories using the vertical rails.



Ref.no	Accessory carrier	Accessory	Rail	Max. load
1	Control panel arm	Control panel	A, B	
2	<ul style="list-style-type: none"> <li>GCX arm with VESA interface</li> <li>Monitor arm slide-in plate</li> <li>Monitor arm 2 pin</li> </ul>	Patient monitor	A, B	12.5 kg (27.5 lbs)
3	Height adjustable arm VESA	PDMS system	B	13.6 kg (30.0 lbs)
Remark: Installed in lowest rail position. Requires Extra mains power outlet option.				
4	<ul style="list-style-type: none"> <li>CPU mounting small *</li> <li>CPU mounting large **</li> </ul>	<ul style="list-style-type: none"> <li>CPU</li> <li>CPU</li> </ul>	<ul style="list-style-type: none"> <li>D</li> <li>D</li> </ul>	<ul style="list-style-type: none"> <li>11.3 kg (24.9 lbs)</li> <li>18.2 kg (40.1 lbs)</li> </ul>
* CPU width/max.height/ max. depth: 76-114/356/356 mm (3.0-4.9/14.0/14.0")				
** CPU width/max.height/ max. depth: 114-178/457/457 mm (4.9-7.0/18.0/18.0")				
5	N/A	Drawer with load	N/A	10 kg (22.0 lbs)
Remark: Incl. vaporizer, etc.				
6	<ul style="list-style-type: none"> <li>Horizontal short rail DIN</li> <li>Horizontal short rail Duoflex</li> </ul>	Suction container	C, D	3 kg (6.6 lbs)
7	Equipment arm	<ul style="list-style-type: none"> <li>Downward pole short</li> <li>Upward pole short</li> <li>Quad hook for cable management</li> <li>Parameter box</li> </ul>	C	7.5 kg (16.5 lbs)
8	<ul style="list-style-type: none"> <li>Backup gas rack O<sub>2</sub> Pin index</li> <li>Backup gas rack Air Pin index</li> <li>Backup gas rack N<sub>2</sub>O Pin index</li> </ul>	Backup gas cylinders	N/A	7 kg (15.4 lbs)/cylinder
Accessory			Equipment weight	Max load
Additional writing table			4.5 kg (9.9 lbs)	5 kg (11.0 lbs)
GCX arm with VESA interface			2 kg (4.4 lbs)	12.5 kg (27.5 lbs)
Monitor arm slide-in plate				
Monitor arm 2 pin				
Equipment arm			3 kg (6.6 lbs)	7.5 kg (16.5 lbs)
CPU mounting			2.5 kg (5.5 lbs)	18.2 kg (40.1 lbs)
Height adjustable arm VESA			10 kg (22.0 lbs)	13.6 kg (30.0 lbs)



## Standards – safety and functionality

Safety	IEC 60601-1:2005 + A1:2012 IEC 60601-1-2:2014 IEC 60601-1-8:2006 + A1:2012 ISO 80601-2-13:2011 + A1:2015 IEC 62304:2006 + Cor1:2008 + A1:2015 ISO 5360:2016	
Electromagnetic compatibility	IEC 60601-1-2:2014 Refer to <b>Electromagnetic Compatibility, Flow-c Anesthesia System</b>	
Respiratory gas monitoring	ISO 80601-2-55:2011	
Anesthetic gas delivery	ISO 80601-2-13:2011 + A1:2015	
Usability	IEC 62366-1:2015	
Cleaning	IEC 60601-1:2005 + A1:2012 ISO 80601-2-13:2011 + A1:2015	
Classification according to IEC 60601-1:		
Class I equipment	According to the type of protection against electrical shock	
Type B equipment	According to the degree of protection against electrical shock	
Continuous operation	According to the mode of operation	
Classification according to EU Medical Directive 93/42/EEC:		
The anesthesia system is classified as IIb		
Classification according to IEC 60529:		
Ingress Protection	IP21 Valid when the patient cassette is in place and the patient cassette lid is closed. Make sure any fluid has been wiped from the connections in the vaporizer slots before connecting a vaporizer.	
IP number	First digit – Solids	Second digit – Liquids
IP21	Protected against solid foreign objects of 12.5 mm diameter and greater.	Protected against vertically falling water drops.

## Power supply

### Mains power

Mains power	100–240 V, AC 50–60 Hz (without auxiliary power outlets) 100–120 V, 220–240 V, AC 50–60 Hz (with auxiliary power outlets)
Power consumption	300 VA (auxiliary power outlets not included) 1500 VA (maximum auxiliary configuration)

### Battery

Type	Sealed acid-lead rechargeable
Capacity	18 Ah
Operating time	Approx. 90 minutes
Charging time	Approx. 6 hours

### Auxiliary power outlets

All auxiliary power outlets are connected to an isolation transformer.  
Voltage depends on mains power supply.

Voltage	Type of electrical outlet	Max load total	Max load from each outlet
<b>220–240 V</b>	4 × IEC	• 5 A	• 5 A
	• 4 × CEE 7/3 (EU)	• 5 A	• 2 A
	• 4 × BS 1363 (UK)	• 5 A	• 2 A
<b>100–120 V</b>	• 4 × IEC	• 10 A	• 10 A
	• 4 × CEE 7/3 (EU)	• 10 A	• 4 A
	• 4 × BS 1363 (UK)	• 10 A	• 4 A
	• 4 × NEMA 5-15R (US)	• 10 A	• 4 A



## Gas supply

All gases and anesthetic agents must conform to the European and American Pharmacopeia.

Required minimum concentration of oxygen supplied to the system is 99.5 vol. % of O<sub>2</sub>.

When using Oxygen 93, the required O<sub>2</sub> concentration is 90.0–96.0 vol.%.

### Central gas

Supply pressure:	
• O <sub>2</sub>	• 250–600 kPa (2.5–6.0 bar, 36–87 psi) *
• Air	• 250–600 kPa (2.5–6.0 bar, 36–87 psi) **
• N <sub>2</sub> O	• 250–600 kPa (2.5–6.0 bar, 36–87 psi)

Hospital central gas supply must be able to deliver a flow of at least 60 l/min at a supply pressure of 280 kPa (2.8 bar, 41 psi)

Connection standards	AGA DISS NIST French standard British standard
----------------------	--

Maximum levels	Air	O <sub>2</sub>
	• H <sub>2</sub> O <7 g/m <sup>3</sup>	• H <sub>2</sub> O <20 mg/m <sup>3</sup>
	• Oil <0.5 mg/m <sup>3</sup>	
	• Chlorine must not be detectable	

If the compressed air is generated by a liquid ring compressor there is a potential risk of chlorine in the supplied air.

Maximum inlet gas temperature	<35 °C (<95 °F)
-------------------------------	-----------------

External vacuum source pressure -0.9 to -0.6 bar (-90 to -60 kPa)

### Backup gas

Cylinder connection standards	• Pin index safety system (PISS) • DIN
-------------------------------	---

DIN connections 6, 9, 12, 13

Backup gas rack, excluding valves	
Size	655 mm × 140 mm (25.8" × 5.5") (H × Ø)
Weight	Max 7 kg per cylinder, including gas.

Cylinder configuration	Max. two cylinders and only one of each
	• O <sub>2</sub>
	• Air
	• N <sub>2</sub> O

Cylinder pressure	
• O <sub>2</sub>	Max. 20,000 kPa (200 bar, 2900 psi)
• Air	Max. 20,000 kPa (200 bar, 2900 psi)
• N <sub>2</sub> O	Max. 8,000 kPa (80 bar, 1160 psi)

Pressure measurement Electronically measured cylinder pressure.

Cylinder safety valve opening pressure	
• O <sub>2</sub>	650 kPa (6.5 bar, 94 psi)
• Air	650 kPa (6.5 bar, 94 psi)
• N <sub>2</sub> O	650 kPa (6.5 bar, 94 psi)

\* The auxiliary O<sub>2</sub> device will not meet specified performance if the supply pressure is less than 300 kPa (3 bar, 44 psi).

\*\* The Venturi vacuum ejector pump will not meet specified performance if the supply pressure is less than 300 kPa (3 bar, 44 psi).

## Suction unit

Type	High vacuum/high flow rate
<b>Vacuum ejector pump – Venturi</b>	
Compressed air consumption (Suction unit)	50–90 NI/min* at a supply pressure equivalent to patient suction supply pressure (Air).
Max free flow (suction flow)	28 NI/min
Max. vacuum (suction)**	-0.9 to -0.6 bar (-90 to -60 kPa), at a supply pressure equivalent to patient suction supply pressure (Air).
<b>External vacuum source – Medical vacuum system</b>	
External vacuum source pressure	-0.9 to -0.6 bar (-90 to -60 kPa)

\* Normal liter (NI) – volume of gas given ambient conditions, for example current atmospheric pressure.

\*\* Max. vacuum varies as a function of atmospheric pressure and supply pressure. Highest performance is obtained at sea level when the supply pressure is approx. 4 bar. Performance decreases with increased altitude.

## Anesthetic Gas Scavenging System (AGSS)

Type	Passive system (including a flow indicator) integrated into the system
Scavenging flow	Minimum 25 l/min (STPD), or 10 l/min (STPD) over the set minute volume, whichever is greater.
Outlet connections	<ul style="list-style-type: none"> <li>• 30 mm ISO taper</li> <li>• DISS EVAC</li> <li>• 12.7 mm/½" in hose Barb</li> <li>• 25 mm/1" Barb</li> <li>• AGA EVAC</li> <li>• WAGD-to-Vacuum connector</li> <li>• 22 mm out. diam. connector and 22 mm int. diam. connection tube</li> </ul>

## Fresh gas flow

Gas mix	Air/O <sub>2</sub> O <sub>2</sub> /N <sub>2</sub> O O <sub>2</sub>	Electronic Servo controlled Electronic Servo controlled Electronic Servo controlled
Fresh gas flow range	<ul style="list-style-type: none"><li>• MAN = 0.1–20 l/min</li><li>• AUTO = 0.1–20 l/min (FGF delivery depending on set MV)</li><li>• AFGO = 1.0–20 l/min</li></ul>	
Fresh gas O <sub>2</sub> /Air Flow (numerical/bargraph)	Selectable	
Fresh gas O <sub>2</sub> /N <sub>2</sub> O Flow (numerical/bargraph)	Selectable	
O <sub>2</sub> concentration accuracy in the fresh gas:		
<ul style="list-style-type: none"><li>• Air/O<sub>2</sub> (21%–100%)</li><li>• O<sub>2</sub>/N<sub>2</sub>O (28%–100%)</li></ul>	<ul style="list-style-type: none"><li>• Fresh gas flow 0.3–20 l/min: ±5%</li><li>• Fresh gas flow 0.3–20 l/min: ±5%</li></ul>	<ul style="list-style-type: none"><li>• Fresh gas flow &lt;0.3 l/min: ±20%*</li><li>• Fresh gas flow &lt;0.3 l/min: ±20% **</li></ul>
Setting resolution, O <sub>2</sub>	1%	
Setting resolution, Flow	0.1 l/min	
O <sub>2</sub> Flush	<ul style="list-style-type: none"><li>• Approximately 56 l/min</li><li>• 2 cmH<sub>2</sub>O expiratory resistance when APL is set to SP</li></ul>	

### Auxiliary O<sub>2</sub>

Auxiliary O <sub>2</sub> flow range	0–15 l/min***
-------------------------------------	---------------

\* Specification valid in typical clinical range for minimal flow anesthesia: Respiratory Rate 5–35 breaths/min, Tidal Volume 100–700 ml, Minute Volume <10 l/min, Pressure 5–40 cmH<sub>2</sub>O, Set O<sub>2</sub> 50–100%.

\*\* Specification valid in typical clinical range for minimal flow anesthesia: Respiratory Rate 5–35 breaths/min, Tidal Volume 100–700 ml, Minute Volume <10 l/min, Pressure 5–40 cmH<sub>2</sub>O, Set O<sub>2</sub> 50–100%.

\*\*\* The apparent gas flow will increase if the ambient pressure decreases.



## Breathing system

Type	Circle system with volume reflector
System volume (incl. absorber, without patient tubings and manual breathing bag)	Approx. 2.8 l
The patient circuit configurations are intended to provide the following range of inspired tidal volumes	<ul style="list-style-type: none"> <li>• Adult: 22 mm tubing, Tidal Volumes 100–2000 ml.</li> <li>• Infant: 15 mm tubing, Tidal Volumes 25–350 ml.</li> <li>• Infant: 10–12 mm tubing, Tidal Volumes 5–100 ml *</li> </ul> <p>* VC: 20–100 ml, PC: 5–100 ml</p>
Drive gas	O <sub>2</sub>
CO <sub>2</sub> absorber, volume	<ul style="list-style-type: none"> <li>• Approx. 700 ml (reusable and disposable version 1)</li> <li>• Approx. 800 ml (disposable version 2)</li> </ul>
CO <sub>2</sub> absorber, absorbent material	Sofnolime™
Patient tube connections	22/15 mm ISO cone
Type of material (breathing circuit system)	<ul style="list-style-type: none"> <li>• PPSU (Polyphenylsulphone)</li> <li>• SBC (Styrene-butadiene copolymer)</li> <li>• PP (Polypropylene)</li> </ul>
System compliance (volume of gas lost due to internal compliance – manual mode only)	Approx. 3 ml/cmH <sub>2</sub> O, i.e. 90 ml at a pressure of 30 cmH <sub>2</sub> O
Inspiratory/expiratory flow resistance of the system (the figures here apply to the breathing tubes recommended by the manufacturer)	<ul style="list-style-type: none"> <li>• 10 mm breathing circuits (including Y-piece): <ul style="list-style-type: none"> <li>– &lt;2.8 cmH<sub>2</sub>O at a flow of 2.5 l/min</li> <li>– &lt;5.5 cmH<sub>2</sub>O at a flow of 15 l/min</li> </ul> </li> <li>• 15 mm breathing circuits (including Y-piece): <ul style="list-style-type: none"> <li>– &lt;1.3 cmH<sub>2</sub>O at a flow of 2.5 l/min</li> <li>– &lt;2.1 cmH<sub>2</sub>O at a flow of 15 l/min</li> <li>– &lt;4.4 cmH<sub>2</sub>O at a flow of 30 l/min</li> </ul> </li> <li>• 22 mm breathing circuits (including Y-piece): <ul style="list-style-type: none"> <li>– &lt;1.1 cmH<sub>2</sub>O at a flow of 2.5 l/min</li> <li>– &lt;2.0 cmH<sub>2</sub>O at a flow of 15 l/min</li> <li>– &lt;3.7 cmH<sub>2</sub>O at a flow of 30 l/min</li> <li>– &lt;6.0 cmH<sub>2</sub>O at a flow of 60 l/min</li> </ul> </li> </ul>

### Manual ventilation

Electronic APL valve	Spontaneous breathing (SP) and adjustable pressure up to 80 cmH <sub>2</sub> O
----------------------	--

### AFGO – Additional Fresh Gas Outlet (option)

Type	<ul style="list-style-type: none"> <li>• 22 mm coaxial/15 mm conical outlet connections</li> <li>• Pneumatic powered SW controlled (from control panel)</li> </ul>
------	--

### Emergency backup ventilation

Emergency APL valve	SP–80 cmH <sub>2</sub> O, SP = 2 cmH <sub>2</sub> O
O <sub>2</sub> emergency flow	0–10 l/min

## Breathing circuits and accessories

Note that the table applies to the breathing circuits recommended by the manufacturer.

Compliance	<ul style="list-style-type: none"> <li>• 10 mm breathing circuits: &lt;0.4 ml/cmH<sub>2</sub>O</li> <li>• 15 mm breathing circuits: &lt;0.7 ml/cmH<sub>2</sub>O</li> <li>• 22 mm breathing circuits: &lt;1.8 ml/cmH<sub>2</sub>O</li> </ul>
Internal volume	<ul style="list-style-type: none"> <li>• 10 mm breathing circuits: 0.4 l</li> <li>• 15 mm breathing circuits: 0.7 l</li> <li>• 22 mm breathing circuits: 1.8 l</li> </ul>
Flow resistance in each limb including Y-piece	<ul style="list-style-type: none"> <li>• 10 mm breathing circuits: <ul style="list-style-type: none"> <li>– At 2.5 l/min: &lt;3.0 cmH<sub>2</sub>O/(l/s)</li> <li>– At 15 l/min: &lt;6.0 cmH<sub>2</sub>O/(l/s)</li> <li>– At 30 l/min: &lt;10 cmH<sub>2</sub>O/(l/s) (30 l/min not applicable for intended patient tidal volume range)</li> </ul> </li> <li>• 15 mm breathing circuits: <ul style="list-style-type: none"> <li>– At 2.5 l/min: &lt;1.5 cmH<sub>2</sub>O/(l/s)</li> <li>– At 15 l/min: &lt;1.5 cmH<sub>2</sub>O/(l/s)</li> <li>– At 30 l/min: &lt;2 cmH<sub>2</sub>O/(l/s)</li> </ul> </li> <li>• 22 mm breathing circuits: <ul style="list-style-type: none"> <li>– At 2.5 l/min: &lt;0.5 cmH<sub>2</sub>O/(l/s)</li> <li>– At 15 l/min: &lt;0.7 cmH<sub>2</sub>O/(l/s)</li> <li>– At 30 l/min: &lt;1.0 cmH<sub>2</sub>O/(l/s)</li> </ul> </li> </ul>
Flow resistance for angled adapter	<ul style="list-style-type: none"> <li>• At 2.5 l/min: &lt;0.2 cmH<sub>2</sub>O/(l/s)</li> <li>• At 15 l/min: &lt;0.6 cmH<sub>2</sub>O/(l/s)</li> <li>• At 30 l/min: &lt;1.0 cmH<sub>2</sub>O/(l/s)</li> </ul>
• Flow resistance for 22 mm joint adapter	<ul style="list-style-type: none"> <li>• At 2.5 l/min: &lt;0.2 cmH<sub>2</sub>O/(l/s)</li> <li>• At 15 l/min: &lt;0.2 cmH<sub>2</sub>O/(l/s)</li> <li>• At 30 l/min: &lt;0.2 cmH<sub>2</sub>O/(l/s)</li> </ul>

# Ventilator

Type	Pneumatic powered Servo controlled
Patient range	Neonatal to Adult
Ventilation modes	<ul style="list-style-type: none"> <li>• Manual/Bag</li> <li>• Additional Fresh Gas Outlet (AFGO, option)</li> <li>• Volume Control (VC)</li> <li>• Pressure Control (PC, option)</li> <li>• Pressure Support (PS, option)</li> <li>• Pressure Regulated Volume Control (PRVC, option)</li> <li>• Synchronized Intermittent Mandatory Ventilation (SIMV, option)</li> <li>• Low VT ventilation (option)</li> <li>• High performance ventilation (option)</li> </ul>
Tidal volume range (volume controlled modes)	20–350 ml, $\pm 10\%$ or 10 ml, whichever is greater * 50–1600 ml, $\pm 10\%$ or 10 ml, whichever is greater 50–2000 ml, $\pm 10\%$ or 10 ml, whichever is greater ** Note that the Flow-c delivers tidal volumes down to 5 ml when using Pressure Control
Tidal volume setting range	Infant range: <ul style="list-style-type: none"> <li>• 20–350 ml, resolution 1 ml *</li> <li>• 50–350 ml, resolution 1 ml</li> </ul> Adult range: <ul style="list-style-type: none"> <li>• 100–1600 ml, resolution 10 ml</li> <li>• 100–2000 ml, resolution 10 ml **</li> </ul>
Minute volume setting range	Infant range: 0.3–20 l/min Adult range: 0.5–60 l/min
Inspiratory pressure (pressure controlled modes)	<ul style="list-style-type: none"> <li>• 0–80 cmH<sub>2</sub>O <math>\pm 15\%</math> or <math>\pm 2</math> cmH<sub>2</sub>O, whichever is greater</li> <li>• 0–120 cmH<sub>2</sub>O <math>\pm 15\%</math> or <math>\pm 2</math> cmH<sub>2</sub>O, whichever is greater **</li> </ul>
Inspiratory pressure setting range	Infant range: <ul style="list-style-type: none"> <li>• 0–80 cmH<sub>2</sub>O, resolution 1 cmH<sub>2</sub>O</li> </ul> Adult range: <ul style="list-style-type: none"> <li>• 0–80 cmH<sub>2</sub>O, resolution 1 cmH<sub>2</sub>O</li> <li>• 0–120 cmH<sub>2</sub>O, resolution 1 cmH<sub>2</sub>O **</li> </ul>
Compressible volume compensation	Yes
Inspiratory flow	200 l/min (3.3 l/s)
Breathing frequency	4–100 $\pm 1$ breaths/minute
I:E (VC, PC)	1:8.3–4:1
PEEP	0–50 cmH <sub>2</sub> O
Trigger	Flow/Pressure
Inspiratory pause (VC)	0 to 30% or 0–1.5 s

\* Option Low VT ventilation is required.

\*\* Option High performance ventilation is required.



## Respiratory monitoring

Administered breaths	1–100 ±1 breaths/minute
Loops	Flow–Volume Volume–Pressure
Lung characteristics	Airway resistance (R <sub>dyn</sub> ) Compliance (C <sub>dyn</sub> ) Elastance (E <sub>dyn</sub> )
Inspiratory Minute Volume	0.3–60 l/min
Accuracy Insp. Minute Volume	±15% or ±15 ml multiplied by the breathing frequency, whichever is greater
Expiratory Minute Volume	0.3–60 l/min
Accuracy Exp. Minute Volume	±15% or ±10 ml multiplied by the breathing frequency, whichever is greater
Inspiratory Tidal Volume	5–2000 ml
Accuracy Insp. Tidal Volume	±4 ml (5–20 ml range) * ±15% or 15 ml, whichever is greater (20–2000 ml range) <small>* Accuracy valid for O<sub>2</sub>/Air gas mix, O<sub>2</sub> concentration at 60%, RR at 30 and I:E ≥1:2</small>
Expiratory Tidal Volume	5–2000 ml
Accuracy Exp. Tidal Volume	+7/-4 ml (5–20 ml range) * ±15% or 10 ml, whichever is greater (20–2000 ml range) <small>* Accuracy valid for O<sub>2</sub>/Air gas mix, O<sub>2</sub> concentration at 60%, RR at 30 and I:E ≥1:2</small>
Mean Airway Pressure	0–100 cmH <sub>2</sub> O
Peak Airway Pressure	0–140 cmH <sub>2</sub> O
End Expiratory Airway Pressure	-40–100 cmH <sub>2</sub> O
Airway Pressure	-30–140 cmH <sub>2</sub> O
Airway pressure accuracy (applicable to all pressure measurements)	±5% or ±2 cmH <sub>2</sub> O, whichever is greater

## Alarms

Expiratory Minute Volume: High	0.5–60 l/min
Expiratory Minute Volume: Low	0.01–40 l/min
Excessive leakage	The difference between the maximum and minimum pressures during inspiration is too low
Airway pressure: High	10–120 cmH <sub>2</sub> O
Continuous APL pressure (manual mode only)	Activated when the measured airway pressure exceeds predefined values for more than 15 seconds. Predefined values depend on current APL setting.
High continuous pressure (automatic mode only)	Airway pressure is constant above set PEEP level +15 cmH <sub>2</sub> O more than 15 seconds
Negative airway pressure	Measured airway pressure is below -10 cmH <sub>2</sub> O for more than one second
Regulated Pressure Limited (PRVC mode only)	Permissible pressure limits pre-set tidal volume
PEEP: High	0–55 cmH <sub>2</sub> O
PEEP: Low	0–47 cmH <sub>2</sub> O
Respiratory Rate: High	1–140 B/min and OFF
Respiratory rate: Low	1–140 B/min and OFF
Apnea	5–45 s and OFF
Long apnea (manual mode only)	Infant: No breath detection for up to 60 s Adult: No breath detection for up to 120 s
Check breathing circuit	Activated when inspiratory and expiratory pressures fail to meet preset requirements because of blocked or disconnected tubing
Limited battery capacity	Less than 18 minutes left of battery operation
No battery capacity	Less than 3 minutes left of battery operation
Water trap missing/Replace water trap	The gas analyzer has detected that a water trap replacement is needed
<b>Gas alarms</b>	
FiO <sub>2</sub> : High	23–99% and OFF
FiO <sub>2</sub> : Low	18–99%
EtO <sub>2</sub> : High	13–99% and OFF
EtO <sub>2</sub> : Low	10–99% and OFF
FiCO <sub>2</sub> : High	0.1–10%
EtCO <sub>2</sub> : High	0.1–10%
EtCO <sub>2</sub> : Low	0.1–9.9% and OFF
FiAA: High	<ul style="list-style-type: none"> <li>• 0.1–5.0% and OFF (ISO)</li> <li>• 0.1–8.0% and OFF (SEV)</li> <li>• 0.1–18% and OFF (DES)</li> </ul>
FiAA: Low	<ul style="list-style-type: none"> <li>• 0.1–5.0% and OFF (ISO)</li> <li>• 0.1–8.0% and OFF (SEV)</li> <li>• 0.1–18% and OFF (DES)</li> </ul>
EtAA: High	<ul style="list-style-type: none"> <li>• 0.1–5.0% and OFF (ISO)</li> <li>• 0.1–8.0% and OFF (SEV)</li> <li>• 0.1–18% and OFF (DES)</li> </ul>
EtAA: Low	<ul style="list-style-type: none"> <li>• 0.1–4.0% and OFF (ISO)</li> <li>• 0.1–6.0% and OFF (SEV)</li> <li>• 0.1–12% and OFF (DES)</li> </ul>
Agent mixture: MAC >3	The MAC <sub>40</sub> of the secondary agent is ≥0.6 and the total MAC <sub>40</sub> value is ≥3

## Alarms (continued)

Agent mixture	The second agent is MAC $\geq 0.6$ and the total MAC value is $< 3$
High continuous MAC	Measured MAC exceeds time limit: <ul style="list-style-type: none"> <li>MAC <math>&gt; 2.2</math>; from starting a new case, until 15 minutes after the first vaporizer activation.</li> <li>MAC <math>&gt; 1.8</math> otherwise</li> </ul>
FiN <sub>2</sub> O: High	Inspiratory N <sub>2</sub> O gas supply $> 80\%$
Occlusion in sampling line	Detected occlusion reported from Y-piece gas analyzer

## Vaporizer

Active vaporizers slots	1						
Agents	Isoflurane, Sevoflurane and Desflurane						
Type	Electronic Injector						
Weight (full)	Approx. 3.2 kg (7.1 lbs)						
Dimensions	70 × 215 × 178 mm (2.8" × 8.5" × 7.0")						
Agent capacity	300 ml						
Residual capacity	30 ml (triggering the low level alarm)						
Setting range	<table> <tr> <td>Isoflurane</td><td>0, 0.3–5%, OFF</td></tr> <tr> <td>Sevoflurane</td><td>0, 0.3–8%, OFF</td></tr> <tr> <td>Desflurane</td><td>0, 1.0–18%, OFF</td></tr> </table>	Isoflurane	0, 0.3–5%, OFF	Sevoflurane	0, 0.3–8%, OFF	Desflurane	0, 1.0–18%, OFF
Isoflurane	0, 0.3–5%, OFF						
Sevoflurane	0, 0.3–8%, OFF						
Desflurane	0, 1.0–18%, OFF						
Accuracy	±15% of set value or ±5% of maximum possible user setting (whichever is greater)						
Filling system	<table> <tr> <td>Isoflurane</td><td>• Maquet filling adapter</td></tr> <tr> <td>Sevoflurane</td><td>• Quik Fil®, Maquet filling adapter and SAFE-T-SEAL filling adapter attached to anesthetic agent bottle</td></tr> <tr> <td>Desflurane</td><td>• SAFE-FIL™</td></tr> </table>	Isoflurane	• Maquet filling adapter	Sevoflurane	• Quik Fil®, Maquet filling adapter and SAFE-T-SEAL filling adapter attached to anesthetic agent bottle	Desflurane	• SAFE-FIL™
Isoflurane	• Maquet filling adapter						
Sevoflurane	• Quik Fil®, Maquet filling adapter and SAFE-T-SEAL filling adapter attached to anesthetic agent bottle						
Desflurane	• SAFE-FIL™						
Emptying system	Maquet drain adapter for SAFE-T-SEAL vaporizer						
Vaporizer filling speed	Approx. 4 ml/s						
Tank liquid level	Optical and electronic						



## Gas analyzer

Measuring technology	O <sub>2</sub> Agents, CO <sub>2</sub> , N <sub>2</sub> O	Paramagnetic sensor IR sensor
Warm-up time	ISO standard accuracy Full accuracy	Within 60 s Within 10 minutes
Sampling flow and tolerance	225 ml/min ±10% (Return to circuit), BTPS condition	
Sampling line	Length: 2.0 m 2.5 m 3.5 m 4.5 m	Inner diameter: 1.3 mm 1.5 mm 1.5 mm 1.5 mm

### Measured parameters

Resp. rate	2–100 breaths/minute	
Respiration rate measurement accuracy	<60 breaths/minute >60 breaths/minute	±1 breath/minute Unspecified
Inspiratory and End-Tidal O <sub>2</sub> Concentration	Yes	
Inspiratory and End-Tidal CO <sub>2</sub> Concentration	Yes	
Inspiratory and End-Tidal N <sub>2</sub> O Concentration	Yes	
MAC Y (age dependent)	Yes	
MAC Brain (age dependent)	Yes	
Inspiratory and End-Tidal Agent Concentration	Yes	

Gas measurement accuracy	Gas conc.	Accuracy [%ABS]	Interference	[%ABS]
O <sub>2</sub>	0–25%	±1	N <sub>2</sub> O	0.2
	25–80%	±2	CO <sub>2</sub>	0.2
	80–100%	±3	Any agent	1.0
N <sub>2</sub> O	0–20%	±2	CO <sub>2</sub>	0.1
	0–40%	±3	O <sub>2</sub>	0.1
	40–80%	±5	Any agent	0.1
CO <sub>2</sub>	0–1%	±0.3	N <sub>2</sub> O	0.1
	1–5%	±0.2	O <sub>2</sub>	0.1
	5–7%	±0.3	Any agent	0.3
	7–10%	±0.5		
	>10%	Unspecified		
Isoflurane	0–1%	±0.15	CO <sub>2</sub>	0
	1–5%	±0.2	N <sub>2</sub> O	0.1
	>5%	Unspecified	O <sub>2</sub>	0.1
			2nd agent	0.2
Sevoflurane	0–1%	±0.15	CO <sub>2</sub>	0
	1–5%	±0.2	N <sub>2</sub> O	0.1
	5–8%	±0.4	O <sub>2</sub>	0.1
	>8%	Unspecified	2nd agent	0.2
Desflurane	0–1%	±0.15	CO <sub>2</sub>	0
	1–5%	±0.2	N <sub>2</sub> O	0.1
	5–10%	±0.4	O <sub>2</sub>	0.1
	10–15%	±0.6	2nd agent	0.2
	15–18%	±1.0		
	>18%	Unspecified		
Drift of measurement accuracy	The accuracy includes stability and effects of device drift during operation between calibrations.			

- The respiration rate limit for accurately measured end-tidal values is <60 breaths/minute for I:E = 1:1, <40 breaths/minute for I:E = 1:2 and <30 breaths/minute for I:E = 1:3.
- The accuracy of the gas measurements may be affected if the Ethanol concentration is higher than 0.1%, the Methane concentration is higher than 1% or the Acetone concentration is higher than 1%.
- The partial pressure and the percentage volume of CO<sub>2</sub>, N<sub>2</sub>O, O<sub>2</sub> and anesthetic agent depend on the amount of water vapor in the breathing gas. A partial H<sub>2</sub>O pressure of 11 cmH<sub>2</sub>O is automatically compensated for by the analyzer. Higher H<sub>2</sub>O partial pressures will further dilute the gas sample; at 30 cmH<sub>2</sub>O the general error of all measured gases is -2%.


## External communication

	Quantity	Type	Description
Serial ports	2	RS232	FCI (Flow Communication Interface) protocol
USB	2	USB 1.1	<ul style="list-style-type: none"> <li>• One port for communication</li> <li>• One port for power supply</li> </ul>
Video out	1	VGA	Interface for slave monitor
Ethernet	1	RJ45	Network connection for use with Remote service

## Ordering information

Flow-c Anesthesia Machine and accessories:  
See separate information in "System flowchart, Flow-c", Order no MX-7880.

## Dimensional drawings



### Next-generation Flow core technology

Reduces the need for an ICU ventilator in the OR:

- Servo gas modules deliver up to 200 l/min inspiratory flow and adjusts pressure and flow several times within every breath.
- Innovative Volume Reflector rebreather ensures accurate and precise gas delivery.
- Small system volume (2.7 l) allows fast wash-in and wash-out, saving time and agent consumption.

### Electronic injection vaporizer

- Fast and precise delivery of anesthetic agents.
- Lightweight and possible to fill during use. Holds 300 ml.
- No heating time for Desflurane.

### Low-flow anesthesia

- Fresh Gas Flow (FGF) limit: 0.1 l/min
- VRI (Volume Reflector Indicator) displays the FGF and minute volume ratio to optimize the rebreathing fraction.

### Active hypoxic guard

- O<sub>2</sub>Guard actively intervenes when there is risk of hypoxia, providing added safety at low flows.

### Pause function

- Temporarily stops gas flows and ventilation, giving you time to focus on the patient.

### Space-saving design

- Compact and lightweight (86 x 68 cm, 115 kg).
- Easy to move and adapt to different situations.

### Battery backup

- 90 minutes for added safety in case of power failure.

### Easy cleaning and service

- Just seven parts need to be dismantled for cleaning, saving costs and supporting infection control.
- Preventive maintenance is optimized to reduce complexity with few and easy accessibility parts that are only changed every two years.

### Low cost of ownership

- Modern, easily upgradable platform.
- Maintenance free components, the oxygen and flow sensors are non-consumptive.
- Lower anesthetic agent consumption.
- Reduced training requirements.



### Smart clips, hooks and arms

For convenient placement of suction device, manual breathing bag and other accessories.



### Large drawer

With slots for spare vaporizers and an optional lockable compartment.



# EC DECLARATION OF CONFORMITY FOR MEDICAL DEVICES

acc. to Directive 93/42/EEC on Medical Devices

**Name and Address of the Manufacturer:** Maquet Critical Care AB  
Röntgenvägen 2  
SE-171 54 Solna, Sweden

On our sole responsibility, we hereby declare that the product(s)

**Product Description:** Anesthesia System

**Product Identification:** Flow-c Anesthesia System,  
System version 4.9

**Product-No.:** See Annex I of this document

**Classification (acc. to Annex IX of MDD):** Class IIb

comply with the relevant provisions of the following Directive(s):

## Directive 93/42/EEC on Medical Devices

**Notified Body:** TÜV SÜD Product Service GmbH  
ID No. 0123  
Zertifizierstelle, Ridlerstraße 65  
80339 München, Germany

**Conformity Assessment Procedure:** acc. to Annex II excluding (4) of Directive  
93/42/EEC

## Directive 2011/65/EU on the restriction of the use of certain substances in electrical and electronic equipment

This declaration is valid until the next change of the product or until the expiration date of the Certificate (No. G1 072017 0014, **dated 2020-01-22**, valid to **2024-05-26**) issued by the notified body. Any modification of the medical device not authorized by us will invalidate this declaration.

For the signature of the Manager for Product Compliance on behalf of Maquet Critical Care AB see left corner of page.

# ANNEX I

## TO THE EC DECLARATION OF CONFORMITY

### FOR MEDICAL DEVICES

Product(s) Covered:

Trade Name	Article No.
<b>Base unit</b>	
Flow-c Anesthesia System	6887700
<b>Additional ventilation modes</b>	
PRVC	6688202
PS with backup	6688200
SIMV (VC) and SIMV (PC)	6688197
Pressure Control	6888006
Low VT ventilation	6887892
High performance ventilation	6887889
<b>Additional options</b>	
HLM/CBP	6887895
Insp hold/Exp hold	6887896
Predicted body weight (PBW)	6887893
Agent usage	6887888
Optional module, Additional Fresh Gas Outlet (AFGO)	6681887
AFGO adapter kit, onsite upgrade	6887999
AFGO adapter kit	6887704
O2 Only	6890455
Oxygen 93	6890925

# ANNEX I

## TO THE EC DECLARATION OF CONFORMITY

### FOR MEDICAL DEVICES

<b>AGC (Automated Gas Control)</b>	<b>6881990</b>
<b>Recruitment Maneuvers</b>	<b>6886518</b>
<b>Absorber</b>	
CO2 Absorber Disposable	6677845
CO2 Absorber disposable 12 pcs	6887738
CO2 Absorber Reusable	6692582
Dust Filter for CO2 Absorber Reusable	6880430
<b>Vaporizer</b>	
Vaporizer Sevoflurane, QUIK FIL	6682285
Vaporizer Desflurane, Saf-T-Fil	6682287
Vaporizer Isoflurane, Maquet filling	6682280
Vaporizer Sevoflurane, Maquet filling (Abbott, Baxter)	6682282
Vaporizer Sevoflurane Maquet filling	6886601
Vaporizer Sevoflurane QUIK FIL	6886611
Vaporizer Isoflurane Maquet filling	6886621
Vaporizer Desflurane SAF-FIL	6886631
Vaporizer Sevoflurane SAFE-T-SEAL	6887135
Vaporizer Sevoflurane SAFE-T-SEAL	6887523
MAQUET filling adapter tube for isoflurane bottle	6682290
MAQUET filling adapter tube for sevoflurane bottle (Abbott, Baxter)	6682292
Cover for Vaporizer Slot	6880158



# ANNEX I

## TO THE EC DECLARATION OF CONFORMITY

### FOR MEDICAL DEVICES

Vaporizer Holder	6880155
<b>O2 flowmeter / suction equipment</b>	
External suction vacuum	6887709
Venturi for Suction	6887710
Hydrofobic Bacterial Filter (5 pcs)	6887498
<b>Backup gas supply</b>	
Backup gas rack O2, Pin index	6887416
Backup gas rack Air, Pin index	6887417
Backup gas rack N <sub>2</sub> O, Pin index	6887418
Connector DIN 6	6886811
Connector DIN 9	6886812
Connector DIN 12	6886813
Connector DIN 13	6886814
<b>EVAC outlet Connections</b>	
EVAC Restrictor	6880156
22 mm out. Diam. Connector	6692740
30 mm out. Diam. ISO taper connection	6673980
12.7 mm/½" in house Barb connector	6679627
25 mm/1" in house Barb connector	6679630
AGA EVAC connector	6679645
WAGD-to-Vacuum connector	6692742

# ANNEX I

## TO THE EC DECLARATION OF CONFORMITY

### FOR MEDICAL DEVICES

Connection tube 22 int. diam, (disposable), 50 m	6189646
DISS EVAC connector	6679647
<b>Gas analyzer parts</b>	
Water trap Dryline A, 10 pcs	6522747
Water trap Dryline 1, 10 pcs	6889352
<b>Extra breathing system parts</b>	
Patient cassette complete, extra	6680325
Volume reflector, extra	6887614
<b>Reusable parts</b>	
Tube joint/adaptor 22 mm	6425581
<b>Hardware mounting kits</b>	
Flexible Man bag arm	6887717
Tubing hook small	6887713
Tubing hook large	6887712
Suction tip holder	6887714
Hook and y-piece holder	6880548
Cable holder	6693747
Lockable compartment for drawer	6887718
Additional table	6676962
Hook and hook holder	9368952

## ANNEX I

### TO THE EC DECLARATION OF CONFORMITY

### FOR MEDICAL DEVICES

Clamp, 1 slot	9358193
Horizontal short rail DIN	6887115
Horizontal short rail Duoflex	6887834
Clamp, for Duoflex rail, 10 x30 mm	6038819
Clamp, for DIN rails, 10 x 25 mm	6038801
Monitor arm 2 pin (Philips)	6887066
Panel interface Height adjustable	6886960
Long divider for big drawer	6886664
Channel cover kit, 10 pcs	6887960
Gas bottle support Flow-c/e	6890148
<b>Top Mount</b>	<b>6892715</b>
<b>Horizontal channel slide in plate</b>	<b>6892716</b>
<b>Horizontal channel mount VESA 75</b>	<b>6892717</b>
<b>For upgrade of installed Flow-c</b>	
Upgrade via freight, Flow-c	6888360
Upgrade via E-delivery, Flow-c	6888359
<b>Cleaning adapter</b>	
Cleaning adapter kit	6681910



# Certificate

No. Q5 072017 0013 Rev. 01

**Holder of Certificate:**

GETINGE ✱

**MAQUET CRITICAL CARE AB**

Röntgenvägen 2  
171 54 Solna  
SWEDEN

**Facility(ies):**

MAQUET CRITICAL CARE AB  
Röntgenvägen 2, 171 54 Solna, SWEDEN

See Scope of Certificate

**Certification Mark:**



**Scope of Certificate:**

**Design, development and manufacturing  
of Anaesthesia, Monitoring, Ventilator  
Systems and Perfusion Systems**

**Applied Standard(s):**

EN ISO 13485:2016  
Medical devices - Quality management systems -  
Requirements for regulatory purposes  
(ISO 13485:2016)  
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see:

[www.tuvsud.com/ps-cert?q=cert:Q5 072017 0013 Rev. 01](http://www.tuvsud.com/ps-cert?q=cert:Q5 072017 0013 Rev. 01)

**Report No.:**

713193727

**Valid from:**

2020-12-30

**Valid until:**

2023-12-29

**Date,**

2020-12-28

Christoph Dicks

Head of Certification/Notified Body





Benannt durch/Designated by  
Zentralstelle der Länder  
für Gesundheitsschutz  
bei Arzneimitteln und  
Medizinprodukten  
www.zlg.de  
ZLG-BS-244.10.08



Product Service

# EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

(Devices in Class IIa, IIb or III)

**No. G1 072017 0014 Rev. 01**

**Manufacturer:**

**MAQUET CRITICAL CARE AB**

Röntgenvägen 2

171 54 Solna

SWEDEN

**Facility(ies):**

MAQUET CRITICAL CARE AB

Röntgenvägen 2, 171 54 Solna, SWEDEN

**Product Category(ies): Anaesthesia, Monitoring, Ventilator and  
Perfusion Systems**

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

**Report No.:**

713161324

**Valid from:**

2020-01-22

**Valid until:**

2024-05-26

**Date,**

2020-01-22

Christoph Dicks

Head of Certification/Notified Body



Benannt durch/Designated by  
Zentralstelle der Länder  
für Gesundheitsschutz  
bei Arzneimitteln und  
Medizinprodukten  
www.zlg.de  
BS-MDR-099



Product Service

## EU Quality Management System Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapters I and III  
(Class IIa and Class IIb Devices)

**No. G10 072017 0015 Rev. 00**

**Manufacturer:**

**MAQUET CRITICAL CARE AB**

Röntgenvägen 2  
171 54 Solna  
SWEDEN

The Certification Body of TÜV SÜD Product Service GmbH certifies that the manufacturer has established, documented and implemented a quality management system as described in Article 10 (9) of the Regulation (EU) 2017/745 on medical devices. Details on device categories covered by the quality management system are described on the following page(s). The Report referenced below summarises the result of the assessment and includes reference to relevant CS, harmonized standards and test reports. The conformity assessment has been carried out according to Annex IX Chapter I and III of this regulation with a positive result. The quality management system assessment was accompanied by the assessment of technical documentation for devices selected on a representative basis. The certified quality management system is subject to periodical surveillance by TÜV SÜD Product Service GmbH. The surveillance assessment shall also include an assessment of the technical documentation for the device or devices concerned on the basis of further representative samples.

**Report No.:** 713170827

**Preceding certificate No.:** This certificate is issued for the first time

**Valid from:** 2020-02-17

**Valid until:** 2025-02-16

**Date of initial issuance / Rev.00:** 2020-02-17

Christoph Dicks  
Head of Certification/Notified Body

**Issue date:** 2020-02-17



Benannt durch/Designated by  
Zentralstelle der Länder  
für Gesundheitsschutz  
bei Arzneimitteln und  
Medizinprodukten  
www.zlg.de  
BS-MDR-099



Product Service

## EU Quality Management System Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapters I and III  
(Class IIa and Class IIb Devices)

**No. G10 072017 0015 Rev. 00**

Device(s):	Risk classification	CND Code	Intended Purpose
INSTRUMENTS FOR ANESTHESIA AND PULMONARY VENTILATION SUPPORT	IIb	Z120301	Intended for respiratory support, administration of anesthetic and treatment of neonatal, pediatric and adult patients.

The validity of this certificate  
depends on conditions and/or  
is limited to the following:

Revision History including  
Changes:

Revision / Issue Date / Report  
Rev. 00 / 2020-02-17 / 713170827