

Passive Middle Ear Implant


mXACT PRO

Partial Prosthesis Kit



1


2


3

4

5

English	3
Deutsch	7
Français	12
Italiano.....	17
Español	21
Português	25
Nederlands	29
Dansk.....	34
Svenska	38
Suomi.....	42
Norsk	46
Polski	50
Türkçe	55
Magyar.....	59
Čeština.....	63
Slovenčina.....	67
Hrvatski	71
Slovenščina	75
Srpski jezik	79
Български език.....	83
Română.....	88
Русский.....	93
Ελληνικά.....	98
Lietuvių kalba	102
Latviešu valoda.....	106
Eesti keel	110
Українська	114
Қазақша.....	119
한국어	124
Tiếng Việt	127
中文.....	131
中文.....	134
العربية.....	140

Instructions for use

Passive Middle Ear Implant mXACT PRO Partial Prosthesis Kit

Device description

The MED-EL passive middle ear implant – tympanoplasty prosthesis is an ossicular replacement prosthesis made of titanium. The ossicular chain can be reconstructed to restore the conductive hearing mechanism within the middle ear.


The provided prosthesis is attached to a metal bar. This bar serves as a mount and a length adjustment fixture for the prosthesis. With this device the length of the prosthesis can be adjusted according to the distance between the tympanic membrane and the stapes head determined intraoperatively. The tympanoplasty Sizers are optional templates of ossicular replacement prostheses. They are made of polypropylene. The Sizers are only to be used during a MED-EL tympanoplasty surgery.

Caution: Do not implant or lose the tympanoplasty Sizers in the middle ear.

Contents of the package

Refer to the corresponding figure at the beginning of these instructions for use.

- 1 Tympanoplasty prosthesis: mXACT PRO partial prosthesis
- 2 Length adjustment fixture: The functional length of the partial prosthesis can be steplessly adjusted from 0.75 mm to 3.5 mm.
- 3 Sizers for partial prosthesis (serving as templates for the prosthesis): The package contains 5 Sizers with functional lengths of 1.0 mm, 1.5 mm, 2.0 mm, 2.5 mm and 3.5 mm. The length (in millimetres) is marked on the individual Sizer and indicates the functional length of the prosthesis.

Sizer for partial prosthesis

Tympanoplasty prosthesis

Intended use

The tympanoplasty prosthesis is intended to be used as a partial replacement of the ossicles of the human middle ear. This passive implant restores the mechanical sound transmission from the tympanic membrane to the oval window. The partial prosthesis is intended to be used for replacement of the ossicular chain if at least the stapes and its footplate are still present and functional. The tympanoplasty prosthesis is intended to treat patients of all ages. The tympanoplasty Sizers are intended to be used temporarily during a partial tympanoplasty replacement of the ossicles of the human middle ear, to determine the functional length of the needed prosthesis. The tympanoplasty Sizers for partial ossicular replacement are intended to intraoperatively measure the distance between stapes head and tympanic membrane (to define the functional length of the tympanoplasty prosthesis). The tympanoplasty Sizers are single-use surgical instruments for transient use.

Indications

The tympanoplasty prosthesis is indicated in case of congenital or acquired defects of the ossicular chain due to e.g. chronic otitis media, traumatic injury, malformation, cholesteatoma. The tympanoplasty prosthesis is indicated to treat patients with inadequate conductive hearing from previous middle ear surgery. The tympanoplasty Sizers are indicated to be used optionally for patients undergoing a tympanoplasty surgery.

Contraindications

The prosthesis and the tympanoplasty Sizers should not be used if:

- the patient is known to be intolerant to the materials
- an acute middle ear infection is present

Caution: The physician must fully assess the potential risks and benefits for the patient and his/her realistic expectations with the device prior to the decision to implant. The physician must exercise medical judgement and consider the patient's complete medical history.

Information about use

Determination of the length of the prosthesis

1. The distance between the tympanic membrane and the stapes head can be measured by using the tympanoplasty Sizers. The Sizers have 5 different lengths from 1.0 mm to 3.5 mm with steps of 0.5 mm or 1.0 mm. For covering the headplate of the prosthesis, the measure of a slice of cartilage (recommended: 0.5 mm) must be factored in to finally determine the length of the prosthesis.
2. It is recommended to remove the Sizers from the package with a microsurgical tool (e.g. micro suction tube, forceps or any other delicate surgical instrument).
3. If a Sizer fits exactly between the tympanic membrane and the stapes head, the length can be calculated directly. If a Sizer is too short and the next size is too long, the length of the prosthesis must be calculated to a size in between.
4. The Sizers must be removed from the patient prior to the implantation of the prosthesis. The Sizers are for transient use only and must not be implanted.

Length adjustment of the prosthesis

Refer to the corresponding figures at the beginning of these instructions for use. All illustrations are representative only.

- ① The shaft of the prosthesis can be rotated to adjust the orientation of the coupling structure. The length of the prosthesis must be adjusted by moving the shaft according to the length determination/calculation as described.
- ② To establish a secure connection between the head plate and the shaft, the head plate must be closed (e.g. using delicate tweezers).
- ③ The excess shaft of the prosthesis shall be trimmed close to the head plate using a scalpel (recommended blade #15). The residual shaft must not exceed a length of approximately 0.2 mm (= shaft diameter).
- ④ The fixture must be lifted carefully.
- ⑤ The prosthesis is ready to be removed from the length adjustment fixture within the packaging.

Possible adverse effects

Surgery of the middle ear involves manipulation of the fragile ossicular bones and exposes the inner ear to the risk of surgical trauma. Serious complications may arise either during or after middle ear surgery that may result in irreparable damage to otologic structures causing irreversible partial or total

loss of hearing. Subsequent surgical procedures may be required to correct these conditions, if possible.

Complications that may occur include: sensorineural hearing loss due to trauma during surgery; vertigo; tinnitus; granuloma; perforation of the tympanic membrane; post-surgery displacement of the prosthesis due to the development of scar tissue; post-surgical translocation of the prosthesis due to trauma or incorrect position; extrusion of the prosthesis; and post-surgical infection, including otitis media. It is also possible that the prosthesis could cause a significant loss in residual hearing. The responsible physician should ensure that the patient is fully apprised of the possible adverse effects of the implant procedure.

Additionally, patients are exposed to the normal risk of surgery and anesthesia. Major ear surgery may result in numbness, swelling, or discomfort around the ear, the possibility of facial palsy, disturbance of balance or taste, or neck pain. If these occur, they are usually transient and resolve within a few weeks after surgery.

Warnings and precautions

Pre-operative

- All devices are intended for single use only. Do not reuse or resterilize.
- Do not use if the sterile packaging is damaged, or if the use-by date (expiration date) has passed.
- Open the package only in a sterile environment to avoid contamination.
- Do not alter the devices and only use them as intended.
- Do not use damaged devices.
- Do not implant the prosthesis if the patient is known to be intolerant to the implant material (titanium). Consider patient intolerances for surgical tool selection.
- Do not implant the prosthesis in direct contact with passive middle ear implants other than from MED-EL.

Intra-operative

Ossiculoplasty is a sophisticated surgical intervention and should be performed by a trained surgeon. Inadequate surgical technique can lead to irreparable damage.

MED-EL offers a variety of ossicular prostheses. The correct type and length of the prosthesis to be implanted is determined intra-operatively. Care must be taken when choosing the proper type and length of the prosthesis to be implanted to avoid post-operative complications (e.g. damage to the inner ear if the selected prosthesis is too long; insufficient contact between the prosthesis and the foreseen middle ear structures if the selected prosthesis is not in direct contact with the stapes head or stapes footplate).

Avoid accidental bending of the prosthesis in order to prevent functional damage.

Avoid excessive pressure on auditory structures of the middle ear, especially on the stapes head or stapes footplate.

The middle ear structures must be treated with care, particularly the area of the oval window and the stapes footplate. Avoid traumatization of the inner ear due to luxation of the stapes footplate.

A slice of cartilage (recommended: 0.5 mm) must be placed onto the implant, just between the head of the implant and the tympanic membrane.

Do not implant or lose the tympanoplasty Sizers in the middle ear.

Post-operative

The patient shall be advised regarding proper hygiene of the operated ear. After surgery, both surgeon and patient must be vigilant for any signs of complications (such as development of a granuloma or loss of benefit over time).

The physician must inform the patient that the patient should avoid situations such as sudden, loud noises or severe ambient pressure variations and that the patient should consult his/her physician with questions regarding post-operative activities.

The Implant Card shall be filled out and handed over to the patient.

MRI safety information

Non-clinical testing has demonstrated that the tympanoplasty prosthesis is MR Conditional.

A patient with this device can be safely scanned in an MR system meeting the following conditions: Static magnetic field of up to 7T, maximum spatial field gradient of 79T/m (7 900 G/cm), maximum force product of 356T²/m (356 000 000 G²/cm).

In non-clinical testing, the worst-case image artefact caused by the device extends approximately 4.59 mm if imaged with a gradient echo pulse sequence in a 7T MR system.

Storage and disposal

The product must be stored in its sterile packaging in a dry place and in a temperature range between -29 °C (-20.2 °F) and +60 °C (+140 °F). Devices must not be used after the use-by date specified on the package. Packaging should be disposed of in accordance with local requirements.

Warranty statement

Please contact your local MED-EL representative for information on our warranty provisions.

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MED-EL Elektromedizinische Geräte GmbH