BunnyLens AF

Providing Aberration Free Vision. Aspheric design, also available as EASY PRELOADED IOL

Technical Specifications

Overall diameter.....	11 mm (10D and above); 11.5mm for lower D
Optic diameter.....	6.00 mm
Haptic angulation.....	5°
Optic design.....	Aspheric
Edge design.....	360° Continuous Square Edge
Power range.....	-5.0 to +10.0 (1D increments) +10.5 to +30.0 (0.5D increment) +31.0 to +40.0 (1D increment)

Additional powers within the range can be supplied by special make to order*

Material.....	Hydrophilic Acrylic HEMA/EOEMA copolymer
Filtration.....	UV blocker and Violet Light Filter
Refractive Index.....	1.46 (hydrated @ 35° c)
Y.A.G laser.....	Compatible
A constant (SRK/T)	Optical / Immersion US biometry: 118.5* Contact US biometry: 118.16*
Placement.....	Capsular Bag

CE Approved

* It is recommended that surgeons personalize their A. Constant based on their surgical techniques and equipment, experience and post-operative results. For more information please visit Hanita Lenses web.

Attributes

Advanced Optical Design

The aspheric BunnyLens AF was designed using the most advanced tools, by a professional R&D team of optical and mechanical engineers. The optical profile was calculated using ZEMAX™ software – a simulating tool for the optical design optimization. Calculations were aimed in order to minimize all aberrations, including the spherical aberration of the cornea, and to optimize the MTF (Modulated Transfer Function) of the IOL.

Eye Model


The Optical design of BunnyLens AF was performed using the advanced Arizona Eye model [1]. The parameters and dimensions of the eye model are consistent with average human data. The model was designed to match clinical levels of aberrations, both on and off axis. The retina curvature is designed to split the tangential and sagittal foci off-axis. The result is an accurate simulation of the visual performance of the BunnyLens AF in the Post-operative eye. [1] Field Guide to Visual and Ophthalmic Optics; Jim Schwiegerling; Nov. 2004.

Geometrical Design

- BunnyLens AF ensures excellent stability and centration due to four- point-fixated mechanical design of the haptics.
- 360° double square edge in order to minimize PCO
- Excellent memory – slow gentle release, superior foldability

Material

- The BunnyLens AF is made from hydrophilic acrylic HEMA/ EOEMA copolymer material with a UV and Violet Light Filter, having a proven excellent reputation and many years of clinical experience
- The BunnyLens AF is characterized by excellent biocompatibility and mechanical quality.
- The BunnyLens AF material incorporates a natural yellow, violet filtering, chromophore for better protection of the retina.

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Design, manufacturing, and sales of: soft and hard contact lenses, intraocular lenses, Capsular Tension Rings and foldable Intra Ocular Lens Delivery System

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Hanita Lenses
Kibbutz Hanita 2288500 ISR

est conforme aux exigences des normes internationales
complies with the requirements of the international standards

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
GMED N° 34047-1

Ce certificat est délivré selon les règles de certification GMED / This certificate is issued according to the rules of GMED certification

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ANNEX II excluding section 4 Directive 93/42/EEC concerning medical devices

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For class III devices, a EC design certificate is required

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Soft and hard contact lenses, intraocular lenses, soft contact lenses 54% water content, capsular tension ring

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LNE/G-MED certifies that, on the basis of the results contained in the file referenced T000216-5-S1, the quality system - for design, manufacturing, and final inspection - of medical devices listed here above complies with the requirements of the Directive 93/42/EEC, annex II excluding section 4

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The validity of the certificate is subject to periodic or unexpected verification

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Hanita Lenses
Kibbutz Hanita 2288500 Israel
Conception, fabrication, contrôle final / *Design, manufacturing, final control*

LNE/G-MED

0459

For the General Director
Lionel DREUX
G-MED Certification Director

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