Wip'Anios Excel,

20753N

9049N

Descriere

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WIP'ANIOS Excel

The choice of **safety** into a **suitable format**

Make your choice!

CLEANSING DISINFECTANT WIPES

Wip'Anios Excel Travel 25

- 1 24 packets 25 wipes each
Size : 130 x 210 mm
Texture : 100% viscose
Ref. 2617.466

2 Wip'Anios Excel 60

- 12 packets 60 wipes each
Size : 240 x 300 mm
Texture : 40% viscose 60% polyester
Ref. 2446.492

3 Wip'Anios Excel 50 et 100

- 12 packets 50 wipes each
Size : 180 x 200 mm
Texture : 100% viscose
Ref. 2446.424

4 6 packets 100 wipes each

- Size : 180 x 200 mm
Texture : 100% viscose
Ref. 2446.655

5 Wip'Anios Excel King Size

- 1 bucket containing 1 refill of 350 wipes
Size : 200 x 360 mm
Texture : 70% viscose 30% polyester
Ref. 2446.672

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EN2016-08 - Non contractual pictures. Dangerous- respect the precautions of use. Medical Devices class I/II. Use biocides with caution. Before using, read the label and information relative to the product. Biocidal product.

WIP'ANIOS Excel

WIP'ANIOS EXCEL wipes combine cleaning and disinfection in a single step making it possible to ensure **macroscopic and microscopic hygiene of high surfaces and medical devices**, at any time, in any environment (operating room, dialysis department, environment near the patients, medical office, interiors of medical vehicles, etc.).

✓ Cleansing disinfectant wipes

- Impregnation solution: Surfa'safe Premium.
- Efficiency ensured until the expiry date.
- Complies with EN 16615 phase 2/step 2.

✓ Ready to use

- Uniform and controlled impregnation, no dilution in qualitative water which is not controlled.
- Assured traceability.

✓ Suitable for all types of surfaces

- In particular, haemodialysis generators, including those that may come into contact with foodstuffs (Group 1 - PT2 & PT4).
- Class IIb non-immersible medical devices.
- Wiping endoscope sheaths in preliminary treatment.
- Wide compatibility with a wide range of materials.

✓ Formats suited to different uses

- Available in packets and buckets.
- 4 different sizes.

SAFETY GUARANTEED

MULTI-USE

ENSURED VERSATILITY

For the safety of your patients, operators and equipment.

SAFETY GUARANTEED

FOR THE USER

The impregnation solution is formulated :

- Without CMR substances
- Alcohol free
- Fragrance free

for **optimal use tolerance and comfort.**

FOR EQUIPMENT

The wipes are **impregnated with a Surfa'safe Premium solution**: its neutral pH and lack of amphoteric surfactants ensure maximum compatibility with a wide range of materials (metals, polymers, polycarbonate, aluminium, etc.).

FOR THE ENVIRONMENT

The **100% biodegradable viscose** of WIP'ANIOS EXCEL (50 or 100) and WIP'ANIOS Excel Travel wipes makes it possible to limit the impact on the environment.

VALIDATED ANTIMICROBIAL EFFICACY

- Cleansing disinfectant wipes validated according to EN 16615.
- Bactericide, fungicide and virucide impregnation solution.

ACTIVE AGAINST	STANDARDS	CONTACT TIME
BACTERIAS	EN 13727	2 minutes
	EN16615	2 minutes
	EN 13727 on VRE, A.baumannii βLSE, E.coli OXA 48, E.cloacae OXA 48, K.pneumoniae OXA 48, MRSA	1 minute
	EN 13697 sur ERV, A.baumannii βLSE, E.coli OXA 48, MRSA	2 minutes
YEASTS	EN 13624	2 minutes
	EN16615	2 minutes
VIRUS	PRV, (surrogate of HBV), BVDV (surrogate of HBV), Vaccinia virus, HSV, Rotavirus	1 minute

STANDARD EN 16615

➤ The application standard phase 2/step 2 **EN 16615** (May 2015) makes it possible to claim activity combining the action of the wipe and the solution by measuring the **bactericidal and yeasticidal efficacy** of a wipe pre-impregnated with a disinfectant cleansing solution.

Test condition: Quantitative surface test with mechanical action using wipes (4 zone test) in 1 to 5 minutes (or 60 minutes), at a temperature between +4°C and +30°C, diluted in 30°f hard water for concentrated products and in the presence of interfering substances (CP: 0.3g/L; CS: 3g/L albumin + 3mL/L erythrocytes).

A significant reduction must be established to claim the "cleansing disinfectant wipe" status: 5 log reduction for bacterias - 4 log reduction for yeast.

MAKE YOUR CHOICE !

Wip'Anios Excel Travel 25

Size : 130 x 210 mm
Texture : 100% viscose

Wip'Anios Excel 50 ou 100

Size : 180 x 200 mm
Texture : 100% viscose

Wip'Anios Excel 60

Size : 240 x 300 mm
Texture : 40% viscose 60% polyester

Wip'Anios Excel King Size

Size : 200 x 360 mm
Texture : 70% viscose, 30% polyester

➤ Microbiological quality managed :
Radio-sterilized wipes for an optimum safety.

Particularly suitable for haemodialysis generators

PERFECT FOR ...

- **Small medical devices** (pressure sensors, blood sugar testers, tensiometres, reflex hammers, conventional ultrasound probes)
- **Wiping endoscopes** preliminary treatment before disinfection
- **Surfaces in the vicinity of patients** (adaptable, bed, stretcher)

PERFECT FOR ...

- **Medical equipment and devices**
- Surfaces requiring wipes in a larger format (external device structures, examination table, dialysis equipment, etc.)

GMED certifie que le système de management de la qualité développé par
GMED certifies that the quality management system developed by
GMED certifica que el sistema de gestión de la calidad adoptado por

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Design, development, manufacturing, servicing and sales of antimicrobial products, detergents and materials for disinfecting processes of medical devices dedicated to hospitals and clinical sectors.

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Valable jusqu'au / Expiry date / Fecha de expiración June 10th, 2021 (included)

Etabli le / Issued on / Fecha de preparación August 3rd, 2018

On behalf of the President
Béatrice LYS
Technical Director

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ANNEX II excluding section 4 Directive 93/42/EEC concerning medical devices

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For class III devices, a EC design certificate is required

Fabricant / Manufacturer / Fabricante

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Detergent disinfectant or disinfectant for invasive and/or non invasive devices.

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LNE/G-MED certifies that, on the basis of the results contained in the file referenced P177315, the quality system - for design, manufacturing, and final inspection - of medical devices listed here above complies with the requirements of the Directive 93/42/EEC, annex II excluding section 4.

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On behalf of the G-MED Certification Director

G-MED Certification Technical Director

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*Medical devices covered by this certificate are referenced on the manufacturer's list of products
dated April 10th, 2018 (11 pages) authenticated by LNE/G-MED on May 3rd, 2018*

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This certificate covers the following site and activities:

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Detergents disinfectants or disinfectants for invasive and/or non-invasive devices :

- Disinfectants for dialysis circuits
- Disinfectants for manual use and/or automatic machines for invasive devices
- Detergent disinfectant for manual use and/or automatic machines for invasive devices
- Detergent disinfectant or disinfectant sprays for non-invasive devices
- Detergent disinfectant or disinfectant sprays for non-invasive devices
- Detergent disinfectant or disinfectant wipes for invasive devices
- Detergent disinfectant or disinfectant wipes for non-invasive devices

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LNE/G-MED 0459

On behalf of the G-MED Certification Director

Béatrice LYS

G-MED Certification Technical Manager

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DESIGN, PRODUCTION, STORAGE AND DELIVERY OF ANTIMICROBIAL AND DETERGENT PRODUCTS FOR HOSPITALS AND THE INDUSTRY, DENTAL, MEDICAL AND PARAMEDICAL SECTORS AND LOCAL INSTITUTIONS.

DESIGN, PRODUCTION, STORAGE AND SHIPPING OF EQUIPMENT FOR DISINFECTION PROCESSES, HOSPITALS AND THE INDUSTRY, DENTAL, MEDICAL AND PARAMEDICAL SECTORS AND LOCAL INSTITUTIONS. RELATED AFTER-SALES SUPPORT. DESIGN AND PROVISION OF RELATED TRAINING.

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- Disinfectants for dialysis circuits
- Disinfectants for manual use and/or automatic machines for invasive devices
- Detergent disinfectant for manual use and/or automatic machines for invasive devices
- Detergent disinfectant or disinfectant sprays for non-invasive devices
- Detergent disinfectant or disinfectant sprays for invasive devices
- Detergent disinfectant or disinfectant wipes for invasive devices
- Detergent disinfectant or disinfectant wipes for non-invasive devices

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AFNOR est une marque déposée - AFAQ is a registered trademark. CERTIF 1313.2.11-2014.