Medtronic

Fortrex™

HP PTA Balloon Catheter
Балонен катетър за HP PTA
Balónkový katéter pro HP PTA
HP PTA-balloonkateter
HP-PTA-Balloonkatheter
Καθετήρας μπαλονιού HP PTA
Catéter de balón para ATP HP
HP PTA balloonkateeter
HP PTA -pallokatetri
Cathéter à ballonnet pour ATP HP
Balonski kateter za HP PTA
HP PTA-balloonkatéter
Kateter Balon HP PTA
Catetere a palloncino HP per PTA
HP TTA баллондық катетері
HP PTA 풍선 카테터
HP PTA balioninis kateteris
HP PTA balonkatetrs
Катетер со балон за HP PTA
HP PTA-balloonkatheter
HP PTA-ballongkateter
Cewnik balonowy HP do PTA
Cateter de balão para PTA HP
Cateter cu balon pentru ATP HP
Катетер баллонный для чрескожной транслюминальной ангиопластики (ЧТА) HP
Balónkový katéter HP PTA
Balonski kateter HP PTA
HP PTA balon-kateter
HP PTA-ballongkateter
HP PTA Balon Kateteri
Балонний катетер HP для ЧТА
Bóng nong mạch vành HP PTA

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttööhjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөніндегі нұсқаулар • 사용 지침 • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanim Talimatları • Інструкція з використання • Hướng dẫn sử dụng

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
1 Device description

Dilatation balloon catheters are used to exert radial force to dilate narrow vessel segments. The Fortrex HP PTA balloon catheter is an over-the-wire (OTW) 0.89 mm (0.035 in) dual lumen catheter (1). A semi-compliant inflatable balloon (2) is mounted at the distal end. The distal atraumatic tip is tapered (3).

The luer marked THRU on the manifold (4) is the proximal opening for the central lumen (5). The central lumen ends at the distal tip of the catheter. This lumen is used to pass the catheter over a guidewire with a maximum diameter of 0.89 mm (0.035 in). The luer marked BALLOON (6) is used to inflate and deflate the dilatation balloon, using a mixture of contrast medium and saline solution. The balloon has 2 radiopaque marker bands (7) for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating or working section of the balloon (8).

Figure 1. Fortrex HP PTA balloon catheter

The Fortrex balloon catheter is available in multiple balloon sizes. Nominal balloon diameter and length are printed on the manifold (4). Refer to the package label for information about catheter length and sheath compatibility.

2 Intended purpose

The intended purpose of the Fortrex balloon catheter is to restore lumen patency and blood flow through native or synthetic arteriovenous dialysis fistulae, or to facilitate stent post-dilatation in native or synthetic arteriovenous dialysis fistulae.

The Fortrex balloon catheter is not for use in direct contact with the heart or the central circulatory system.

2.1 Intended patient population

The Fortrex balloon catheter is intended for adult patients in need of percutaneous transluminal angioplasty to improve and maintain luminal diameter in native or synthetic arteriovenous dialysis fistulae.

2.2 Indications for use

The Fortrex HP PTA balloon catheter is indicated for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in native or synthetic arteriovenous dialysis fistulae.

2.3 Contraindications

The Fortrex HP PTA balloon catheter is contraindicated for use in coronary arteries or the neurovasculature, or when unable to cross the target lesion with a guidewire.

2.4 Clinical benefits

The clinical benefit of treating obstructive lesions of native or synthetic arteriovenous dialysis fistulae with the Fortrex balloon catheter is restoring lumen patency and blood flow to preserve hemodialysis access sites.

2.5 Intended users

Interventional procedures using the Fortrex balloon catheter should only be performed by physicians who have experience with interventional techniques in the vascular system.

This device should only be used within a clinical setting in a sterile environment.

3 Performance characteristics

The Fortrex HP PTA balloon catheter is a dual-lumen catheter with a balloon that is bonded onto the distal tip. The Fortrex balloon catheter balances push and track to deliver high pressure, targeted treatment to challenging lesions. The balloon design is optimized for shape retention and peak performance at the rated high pressure.

4 Warnings

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- This device should only be used by physicians who have experience with and who understand the clinical and technical aspects of percutaneous transluminal angioplasty.
- To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Never use air or any gas to inflate the balloon.
- Consider precautions to prevent or reduce clotting. Physician experience and discretion determines the appropriate anticoagulation therapy for each patient.
- Do not manipulate the balloon when it is inflated.
- Use the catheter with caution for procedures involving synthetic vascular grafts.
- Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.
- Use direct fluoroscopic observation while manipulating the catheter in the vascular system.
- Do not retract the balloon catheter unless the balloon is free from the stent and fully deflated under vacuum, as vessel trauma or product damage may occur.

5 Precautions

- Do not use contrast media that are contraindicated for intravascular use with this device.
- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.
- To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Ensure that the guidewire is in place before changing the position of the balloon catheter.
- Balloons should not be inflated in excess of the rated burst pressure.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

6 Potential adverse events

The potential adverse events (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

Abrupt or subacute closure	Fever	Regional ischemia
Allergic reaction to device materials or procedure medications	Hematoma	Sepsis
Amputation	Hemorrhage	Shock
Aneurysm	Hypertension or hypotension	Stenosis or restenosis
Arrhythmia	Iatrogenic arteriovenous fistula	Stroke
Bleeding requiring transfusion	Infection	Transient ischemic attack
Contrast medium reaction or renal failure	Inflammation	Venous thromboembolism
Death	Intraluminal thrombus	Vessel injury (such as dissection, perforation, or rupture)
Embolism	Myocardial infarction	Vessel spasm
Endocarditis	Pseudoaneurysm	

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

7 How supplied

Warning: This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

The catheter is packaged with a protective sheath positioned over the balloon.

The device is supplied sterile and is sterilized using ethylene oxide.

8 Storage

Store the device in a dry area away from sunlight.

Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

9 Device selection

Consider the following details when selecting a catheter:

- Consult the balloon compliance card that is packaged with the device for the diameters of the balloons at given pressures.
- The diameter of the expanded balloon should not exceed the diameter of the vessel immediately distal to or proximal to the stenosis.
Warning: To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Verify that the selected accessories accommodate the balloon catheter based on sizing information on the packaging label. Refer to the product labeling for size information.
Caution: Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.

10 Recommended items

Prepare the following items using sterile technique:

- 10 mL syringe filled with sterile heparinized saline
- Contrast medium (The standard inflation medium is a 1:1 mixture of contrast medium and sterile saline.)
Caution: Do not use contrast media that are contraindicated for intravascular use with this device.
- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

11 Preparation

1. Carefully remove the catheter from the inner package.
2. Remove and discard the protective tubing.

11.1 To displace air:

Note: As supplied, the balloon and the guidewire lumen of the catheter contain air. This air must be displaced to make certain that only liquid fills the balloon while the catheter is in the bloodstream.

1. Flush out the lumen for the guidewire with heparinized saline solution through the wire luer marked THRU.
2. Mix equal volumes of contrast medium and normal saline.
3. Fill a 10 mL syringe with approximately 4 mL of the mixed contrast medium solution.
Warning: Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.
Caution: Do not use contrast media that are contraindicated for intravascular use with this device.
4. Evacuate air from the balloon segment:
 - a. Attach the syringe to the luer marked BALLOON. Hold the syringe pointing down with the plunger pointed upward.
 - b. Apply negative pressure and aspirate the lumen for 15 seconds. Slowly release the pressure to neutral so that the contrast medium fills the shaft of the dilatation catheter.
 - c. Disconnect the syringe from the BALLOON luer of the dilatation catheter.
 - d. Remove all air from the syringe barrel. Reconnect the syringe to the BALLOON luer of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
 - e. Slowly release the device pressure to neutral.
Note: On first use, avoid any positive pressure to prevent the balloon from unfolding and partially inflating outside of the vessel.
5. Repeat steps 4b to 4e as necessary to remove all air from the system.
6. Replace the syringe with a manometer-controlled dilatation system. Do not introduce air into the system.

12 Insertion and dilatation

1. If applicable, use an appropriately sized introducer sheath to introduce the catheter percutaneously. Refer to the product labeling for sizing information.
Caution: Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.

2. After the guidewire is in place through the lesion, place the prepared catheter over the guidewire and advance the catheter tip to the introduction site. Use a guidewire that is of suitable length to maintain control and position of the guidewire at all times.
- Warning:** Use direct fluoroscopic observation while manipulating the catheter in the vascular system.
- Warning:** If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.
- Caution:** To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
3. Position the catheter with the center of the balloon in the stenotic area of the vessel. The radiopaque marker bands indicate the working length of the balloon.
- Caution:** Ensure that the guidewire is in place before changing the position of the balloon catheter.
4. When the radiopaque marker bands are appropriately positioned, inflate the balloon to dilate the target area.
- Warning:** Never use air or any gas to inflate the balloon.
- Warning:** Do not manipulate the balloon when it is inflated.
- Caution:** Balloons should not be inflated in excess of the rated burst pressure.
- Note:** Repeat inflation of the balloon (maximum 10 times) until the desired result is achieved.
5. Using the inflation device, deflate the balloon by aspirating the catheter.
6. Withdraw the catheter while maintaining a vacuum in the balloon.
- Warning:** Do not retract the balloon catheter unless the balloon is free from the stent and fully deflated under vacuum, as vessel trauma or product damage may occur.
- Note:** Multiple balloon inflations and deflations may cause resistance during device withdrawal. If resistance is encountered, use a gentle counterclockwise twisting motion to ease withdrawal through the sheath. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath together, while maintaining vessel access with the guidewire.
7. Use angiography to verify dilatation results.
8. Repeat *Step 1* through *Step 7* as needed to achieve the desired results.

13 Disposal

Caution: Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

14 Disclaimer of warranty

Important: This disclaimer of warranty does not apply in any countries where such a disclaimer is not permitted by law.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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The diagram illustrates the Fortrex HP PTA balloon catheter. It is a long, thin medical device with a handle on the left and a catheter tip on the right. The handle has a control knob (5) and a connector (4). The catheter body has a central lumen (1) and a side lumen (2). The distal end features a balloon (8) with two radiopaque marker bands (7). The catheter is shown in a slightly curved position, with a break in the middle indicated by two parallel diagonal lines. The components are numbered 1 through 8, corresponding to the descriptions in the text.

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Medtronic and Medtronic logo are trademarks of Medtronic. ™™ Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic, Inc.
Medtronic B.V.
Medtronic Limited