ALBU-0600	R 2 x 125 mL + Std 1 x 2 mL
ALBU-0700	R 4 x 250 mL + Std 1 x 2 mL
ALBU-0250	R 12 x 20 mL
☞ ALBU-5600	R 125 mL
☞ ALBU-5700	R 250 mL
☞ ALBU-5220	R 20 mL

FTRO-ALBU-v22 (07/2024) (PIT-ALBU-4-v22)

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BIBLIOGRAFIE

1. Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 66.
2. Dufour D. R., The Liver: Function and Chemical Pathology. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
3. Johnson, A. M., Amino Acids and Proteins. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 286.
4. Doumas, B. T., et al., Albumin standards and the measurement of serum albumin with bromocresol green, Clin Chim Acta (1971), **31**, 87.
5. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
6. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACC Press (1997).
7. Young D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACC Press (1995).

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 64.
3. Dufour, R., *The liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
4. Schumann, G., et al., *Clin Chem Lab Med.* (2002), **40**, 718.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Panteghini, M. & Biais, R., *Enzymes*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Dufour, D.R., *The Pancreas: Function and Chemical Pathology*, Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 651 and appendix.
3. Winn-Deen, E.S., et al., Clin. Chem., (1988), **34**, 2005
4. Guder W.G, et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples*, (2002). WHO/DIL/LAB/99.1 Rev.2.
5. Schumann, G., et al., Clin. Chem. Lab Med., (2006), **44**, 1146.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACCC Press, (1997).
8. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACCC Press, (1995).

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BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFIA/BIBLIOGRAFIA

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 64.
3. Dufour, R., *The liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
4. Schumann, G., *et al., Clin Chem Lab Med.*, (2002), 40, 718.
5. Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples.* (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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- Special Programming instructions :
see § INSTALLATION AND USE

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English - EN

INTENDED USE

AST/GOT 4+1 SL is an *in vitro* diagnostic reagent intended for the quantitative determination of aspartate aminotransferase (AST) in human serum and plasma samples on analyzers or semi-automatic analyzers. This *in vitro* diagnostic device is for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Aspartate aminotransferase (AST) is a transaminase also known as glutamate oxaloacetate transaminase (GOT). AST catalyzes the transfer of the aminogroup of L-aspartate to α-ketoglutarate to give L-glutamate. AST is present in both mitochondria and cytoplasm and is widely distributed in the body, with high concentrations in the heart, liver, skeletal muscles and kidneys. Damage to any of these tissues may cause raised levels. In case of acute forms of hepatitis, especially viral hepatitis, the enzyme level is markedly elevated. In case of myocardial infarction, AST activity increases and reaches a peak after 18-24 hours. The activity falls back to normal after 4-5 days, provided no new infarct has occurred.
The following pathological states are examples of disorders also resulting in an increase of enzyme activity: other liver injury (e.g. chronic hepatitis, liver carcinoma, intake of alcohol, administration of drugs), muscular dystrophy, other damage to cardiac or skeletal muscles and acute pancreatitis.

AST measurement is indicated to screen patients at risk of liver diseases and to help diagnose or monitor liver diseases.

LIMITATION OF USE



The quantitative assay of aspartate aminotransferase (AST) alone cannot be used to diagnose a disease or a specific pathology.

The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁴⁾

IFCC method without pyridoxal phosphate (P-5'-P).

Kinetic.

MDH = Malate dehydrogenase.

COMPOSITION

Reagent 1 : R1

Tris buffer, pH 7.80 (30°C)

L-Aspartate 330 mmol/L

LDH ≥ 2 000 U/L

MDH ≥ 1 000 U/L

Sodium azide < 0.1% (w/w)

Reagent 2 : R2

α-Ketoglutarate 78 mmol/L

NADH 1,1 mmol/L

Sodium azide < 0.1% (w/w)

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2

- CONT-0060 ELITROL I

- CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L).

- Analyzers or semi-automatic analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.

- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Do not interchange reagent vials from different kits.

STABILITY

Store at 2-8 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES

Specimen ^(2,8)

- Serum.

- Plasma (lithium heparin).

- Using any other specimen type should be validated by the laboratory.

Warnings and precautions

- Samples must be free from hemolysis. ^(1,3)

- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

Storage and stability ^(2,5)

- 24h at room temperature.

- 7 days at 2-8 °C.

- 3 month at -20 °C.

REFERENCE VALUES ⁽⁴⁾

<i>Serum/plasma</i>	U/L	μkat/L
Men	≤ 35	≤ 0.58
Women	≤ 31	≤ 0.52

Reference values for infants are higher than for adults.

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

INSTALLATION AND USE

For use on Selectra Pro analyzers :

- Consult operator manual.

- Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer. Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

PROCEDURE

Manual Procedure

Wavelength : 340 nm

Optical path : 1 cm

Sample/reagent ratio : 1:20

Temperature: 37 °C

Read against distilled water.

Working reagent (4 volumes of R1 + 1 volume of R2)	1000 μL
Sample	50 μL

Mix and after 1 minute incubation, read absorbance at 1 minute intervals during 3 minutes. Calculate the change of absorbances per minute (ΔA/min)

Automatic Procedure

- These reagents may be used in several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

- High AST values may induce falsely low results due to the depletion of the substrate (total consumption of NADH before reading of the result). For Selectra Analysers, the application contains a specific alarm to warn the users.

CALCULATION

Activity (U/L) = ΔA /min x 3333

Conversion factor : U/L x 0.0167 = μkat/L

CALIBRATION

ELICAL 2 is traceable to IFCC reference method.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,

- at least once per day,

- after every calibration,

- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range

10.0 to 450.0 U/L (0.17 to 7.50 μkat/L).

Samples having greater concentrations should be diluted 1:10 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 4500.0 U/L (75.00 μkat/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

LoD = 2.5 U/L (0.04 μkat/L)

LoQ = 5.0 U/L (0.08 μkat/L)

- Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		Mean		Within-run	Total
	n	U/L	μkat/L	CV (%)	
Level 1	80	34.1	0.57	1.7	3.4
Level 2	80	67.9	1.13	0.8	1.9
Level 3	80	353.6	5.89	0.4	2.0

- Correlation

A comparative study has been performed between AST/GOT 4+1 SL reagent on a Selectra ProM analyzer and a similar commercially available system on 114 human serum samples. The sample concentrations ranged from 10.0 to 529.8 U/L (0.17 - 8.83 μkat/L). The results are as follows : Correlation coefficient : (r) = 0.999 Linear regression: y = 0.927 x - 0.3 U/L (0.01 μkat/L)

- Limitations/Interferences

- AST can be underestimated in case of severe vitamin B6 deficiency. ⁽³⁾

- Hemolyzed samples should not be used since significant hemolysis may increase AST concentration because of high levels of AST in erythrocytes. ⁽³⁾

- Studies have been performed to determine the level of interference from different compounds. The following AST levels were tested: 35.0 U/L and 350.0 U/L.

No significant interference is defined by a recovery ±10% of the initial value.

Unconjugated Bilirubin : No significant interference up to 30.0 mg/dL (513 μmol/L).

Conjugated Bilirubin : No significant interference up to 29.5 mg/dL (504 μmol/L).

Triglycerides : No significant interference up to 2400 mg/dL (27.1 mmol/L).

Pyruvate : No significant interference up to 3.0 mg/dL (341 μmol/L).

Ascorbic acid : No significant interference up to 20.0 mg/dL.

Acetylsalicylic acid : No significant interference up to 20.0 mg/dL.

Acetaminophen : No significant interference up to 30.0 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results. ⁽⁶⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ^(7,8)

- On board stability/Calibration frequency
On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.
☛The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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BIBLIOGRAFIE

1. Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R. Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
2. Dufour, D.R. *The liver: Function and Chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
3. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
4. Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 172.
5. Howard, C.R., *Breastmilk & Breastfeeding Jaundice*, Pediatric Clinical Advisor: Instant Diagnosis and Treatment, 2nd Ed., Garfunkel, L.C., Kaczorowski, J.M. & Christy, C., (Mosby Inc. eds.), (2007), 82.
6. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACC Press (1997).
8. Young D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACC Press (1995).

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BIBLIOGRAFIE

1. Itani, O., Tsang, R.C., *Bone disease, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 614 and appendix.
2. Endres, D.B., Rude, R.K., *Disorders of Bone. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 711.
3. Foley, K. F., Boccuzzi, L., *Urine Calcium: Laboratory Measurement and Clinical Utility, Labmedicine* (2010), **41**, 683.
4. Bauer, P. J., *Anal. Biochem.*, (1981), **110**, 61.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 202.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
2. Burnett, J.R., *Coronary Artery Disease: Lipid metabolism*. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds), (2010), 691 and appendix.
3. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), JAMA, (2001), **285**, 2486.
4. Ahmadraji, T. & Killard, A., Analytical Methods, (2013), **5**, 3612.
5. Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 564.
6. Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
8. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
9. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
2. Burnett, J.R., *Coronary Artery Disease: Lipid metabolism*. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 691 and appendix.
3. Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
4. Allain, C.C., *et al.*, Clin. Chem., (1974), **20**, 470.
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
8. Young D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Wu, A. H. B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 306
2. Sanhai, W.R. et al., Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis. Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 677 and appendix.
3. Apple, F. S. & Jaffe, A. S., Cardiovascular Disease, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 614.
4. Boot, E., de Haan, L., Psychopharmacology (Berl), (2000), 150, 347.
5. Schumann, G., et al., Clin. Chem. Lab. Med., (2002), 40, 635.
6. Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), 59, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
8. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995)

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BIBLIOGRAFIE

1. Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
2. Burnett, J.R., *Coronary Artery Disease: Lipid metabolism*. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds), (2010), 691 and appendix.
3. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), JAMA, (2001), **285**, 2486.
4. Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 684.
5. Ahmadraji, T. & Killard, A., Analytical Methods, (2013), **5**, 3612.
6. Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
8. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
9. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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Immuno-inhibition/IFCC method - Kinetic.

CK-MB SL/CK-MB reagent contains an antibody inhibiting specifically CK-M subunits (i.e. 100% of CK-MM and 50% of CK-MB). The remaining activity, corresponding to CK-B fraction activity, is quantified according to IFCC reference method for total CK measurement. CK-MB activity is then obtained by multiplying by 2 the remaining activity.

CK-B

Creatine Phosphate + ADP $\xrightarrow{\text{CK-B}}$ Creatine + ATP

ATP + D-Glucose $\xrightarrow{\text{Hexokinase}}$ D-Glucose-6-Phosphate + ADP

G-6-P + NADP⁺ $\xrightarrow{\text{G-6-PDH}}$ D-Gluconate-6-Phosphate + NADPH + H⁺

G-6-P : D-Glucose-6-Phosphate

G-6-PDH : Glucose-6-Phosphate Dehydrogenase

The rate of increase of NADPH concentration measured at 340 nm is directly proportional to the enzymatic activity.

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- Stabilitatea la bord

Refer to § PREPARATION AND STABILITY OF WORKING REAGENT.

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BIBLIOGRAFIE

1. Wu, A. H. B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 306
2. Sanhai, W.R. *et al.*, Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 677 and appendix.
3. Apple, F. S. & Jaffe, A. S., Cardiovascular Disease, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 614.
4. Schumann, G., *et al.*, Clin Chem Lab Med., (2002), **40**, 635.
5. Neumeier, D., *et al.*, Clin.Chim.Acta., (1976), **73**, 445.
6. Klein, G., *et al.*, Clin. Chem., (2001), **47**, Suppl. A30.
7. Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
8. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
9. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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English - EN

INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of quantitative tests listed in the value sheet.

These *in vitro* diagnostic devices are for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Concentrations for each analyte to test are lot-specific.

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- VitalScientific reagents listed in the value sheet.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated on the instructions for use of the reagent used in combination with the recommended calibrator.

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.
- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.

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PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.
- Carefully open the vial, avoiding the loss of lyophilizate.
- Add exactly 5 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

Important: To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.
- Do not use after the expiry date stated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.
- Stability of the components:
Between 2 and 8 °C : 5 days
Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light) :
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of direct bilirubin (when stored protected from light) :
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of alkaline phosphatase:
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 4 weeks (when frozen once)

Note:

- Store the vials tightly capped and protected from light after reconstitution.
- Mix content thoroughly prior to use.
- Thawing frozen specimen at 37°C, and vigorous mixing (vortex), should be avoided as this may denature ferritin.

PROCEDURE

To use ELITROL I & II, follow the procedure described in the instructions for the use of the reagent used.

LIMITATIONS

Using ELITROL I & II has been validated with the VitalScientific systems (Analyzers and reagents listed in the value sheet).

Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.

Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values.

For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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BIBLIOGRAFIE

1. Schreiber, W.E., *Iron and Porphyrin Metabolism, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 755 and appendix.
2. Higgins, T., Beutler, E. and Dumas, B.T., *Hemoglobin, Iron and Bilirubin, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
3. Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 634.
4. Higgins, T., *Clin. Chem.*, (1981), 27, 1619.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Sacks, D.B., *Carbohydrates. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus. Clinical Chemistry : Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 729 and appendix
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 444.
4. World Health Organization (WHO), *Diagnosis and Management of Type 2 Diabetes*, (2020).
5. Neeley, W.E., *Clin. Chem.*, (1972), **18**, 509.
6. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
9. Young D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 470.
3. Dufour, R., *The Liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
4. Schumann, G., *et al.*, *Clin Chem Lab Med.*, (2002), **40**, 734
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations*, WHO/DIL/LAB/99.1 Rev.2, (2002).
6. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.* (2004), **59**, 263.
7. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995)

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BIBLIOGRAFIE

1. Sacks, D.B., *Carbohydrates. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus. Clinical Chemistry : Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 729 and appendix
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 444.
4. World Health Organization (WHO), *Diagnosis and Management of Type 2 Diabetes*, (2020).
5. Trinder, P., *Ann. Clin. Biochem.*, (1969), 6, 24.
6. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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FTNA-HBACCa-v8 (04/2024)

CAUTION: Federal Law restricts this device to sale by or on the order of a licensed healthcare practitioner (Rx ONLY)

INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of HbA1c reagent, used for HbA1c quantitative determination on Selectra Pro Series Analyzers according to the values indicated in the value sheet included in the kit.

COMPOSITION

Calibrators are lyophilised hemolysates prepared from human erythrocytes. Stabilizers are added to maintain hemoglobin in the reduced state for the accurate calibration of the hemoglobin A1c procedure.

TRACEABILITY

Values for HbA1c Calibrator Set have been standardized from NGSP/DCCT (%) system. These values have been converted by calculation to traceable values from IFCC reference method (mmol/mol)⁽¹⁾.

REQUIRED MATERIAL BUT NOT PROVIDED

- VitalScientific HbA1c reagent
- Selectra Pro Series Analyzers
- Distilled or deionized water
- General Laboratory equipment.

PRECAUTIONS AND WARNING

- This *in vitro* diagnostic device is for professional use only.
- Take normal precautions and adhere to good laboratory practice.
- All human material should be considered as potentially infectious. All products derived from blood are prepared exclusively from the blood of donors tested individually by FDA - approved methods and found to be negative for HbsAg and antibodies to HCV and HIV1/HIV2. However as no test method can rule out the potential risk of infection with absolute certainty, handle cautiously as potentially infectious.
- Use clean or single use laboratory equipment only to avoid contamination.
- For more information, material Safety Data Sheet (SDS) is available on request for the professional user.

DISPOSAL OF WASTE

Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PREPARATION

Reconstitute each calibrator with 0.5 mL of deionized water. Gently mix for 10 minutes, or until all material has dissolved.

Note : As all samples, these calibrators should be hemolysed before use.

REAGENT DETERIORATION

Reconstituted calibrators which show any evidence of contamination should be discarded.

STABILITY AND STORAGE

Prior to reconstitution : When stored at -2-8 °C and protected from light, the calibrators are stable until the expiry date stated on the label.

After reconstitution : Calibrators are stable for 30 days when stored at -2-8 °C.

- After opening, the vials should be kept correctly and tightly capped to prevent contamination and evaporation.
- Do not freeze.

PROCEDURE

To calibrate with HbA1c Calibrator Set follow the procedure described in the instructions for use of the HbA1c reagent. HbA1c results for samples and controls are determined using the calibration curve obtained with 4 levels of calibrators.

CALIBRATION VALUES

The concentrations indicated in the value sheet enclosed in the kit are lot-specific. The values were determined and defined by VitalScientific on Selectra Analyzers using HbA1c reagent.

For users with Selectra TouchPro software, use the values included in the barcode available on the value sheet.

BIBLIOGRAPHY

1. Weykamp, C. et al., *Clinical Chemistry*. (2008), **54**, 240.
2. Occupational Safety and Health Standards : bloodborne pathogens. (29 CFR Part 1910.1030). Fed.Register. July 1, 2001; 17:260-273.
3. Directive 2000/54/EC. Official Journal of the European Communities N° L262 from October 17, 2000.

SYMBOLS

	Batch code
	Consult instruction for use
	<i>In vitro</i> diagnostic medical device
	Manufacturer
	Temperature limitation
	Use by
	Catalogue number
	Content
	Reconstitute with x mL
	Calibrator 1
	Calibrator 2
	Calibrator 3
	Calibrator 4
	Add exactly 0.5 mL of distilled or deionised water
	Keep at room temperature for 10 minutes before use
	These products are stable x days at x °C
	These products should be immediately and tightly capped to prevent contamination and evaporation
	Modification from previous version

HbA1c CONTROL L+H

FTNA-HBACCn-v9 (04/2024)

References :
HBAC-0049

Kit composition:
Control L : 2 x 0.5 mL Control H : 2 x 0.5 mL

CAUTION: Federal Law restricts this device to sale by or on the order of a licensed healthcare practitioner (Rx ONLY)

INTENDED USE

This *in vitro* diagnostic device is intended for accuracy control of the quantitative HbA1c reagent on Selectra Pro Series Analyzers according to the values indicated in the value sheet included in the kit.

COMPOSITION

Controls are lyophilised hemolysates prepared from human erythrocytes. Stabilizers are added to maintain hemoglobin in the reduced state providing complete control of HbA1c procedure.

TRACEABILITY

Values for HbA1c Control L + H have been standardized from NGSP/DCCT(%) system. These values have been converted by calculation to traceable values from IFCC reference method (mmol/mol)⁽¹⁾.

REQUIRED MATERIAL BUT NOT PROVIDED

- VitalScientific HbA1c reagent
- Selectra Pro Series Analyzers
- Distilled or deionized water
- General Laboratory equipment

PRECAUTIONS AND WARNING

- This *in vitro* diagnostic device is for professional use only.
- Take normal precautions and adhere to good laboratory practice.
- All human material should be considered as potentially infectious. All products derived from blood are prepared exclusively from the blood of donors tested individually by FDA - approved methods and found to be negative for HbsAg and antibodies to HCV and HIV1/HIV2. However as no test method can rule out the potential risk of infection with absolute certainty, handle cautiously as potentially infectious.
- Use clean or single use laboratory equipment only to avoid contamination.
- For more information, Material Safety Data Sheet (SDS) is available on request for the professional user.

DISPOSAL OF WASTE

Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PREPARATION

Reconstitute each control with 0.5 mL deionized water. Gently mix for 10 minutes, or until all material has dissolved.

Note : As all samples, these controls should be hemolysed before use.

REAGENT DETERIORATION

Reconstituted controls which show any evidence of contamination should be discarded.

STABILITY AND STORAGE

Prior to reconstitution : when stored at -2-8 °C and protected from light, the controls are stable until the expiry date stated on the label.

After reconstitution : The controls are stable for 30 days when stored at 2-8 °C.

- After opening, the vials should be kept correctly and tightly capped to prevent contamination and evaporation.
- Do not freeze.

PROCEDURE

HbA1c Control L + H are handled as samples according to the instructions for use of HbA1c reagent. These controls should be assayed together with patient samples, at least once a day and after each calibration. The control frequency should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local, state, and/or federal guidelines.

CONTROL VALUES

The concentrations indicated in the value sheet enclosed in the kit are lot-specific.

These values have been adjusted to be in the normal range (HbA1c Control L) or in the pathological range (HbA1c Control H). The target and range values were determined under strictly standardized conditions on Selectra Pro Series Analyzers using the HbA1c reagent and HbA1c Calibrator Set.

For users with Selectra TouchPro software, use the values included in the barcode available on the value sheet.

BIBLIOGRAPHY

- 1 - Weykamp, C. et al., *Clinical Chemistry*, (2008), **54**, 240.
2. Occupational Safety and Health Standards : bloodborne pathogens. (29 CFR Part 1910.1030). Fed.Register. July 1, 2001; 17:260-273.
3. Directive 2000/54/EC. Official Journal of the European Communities N° L262 from October 17, 2000.

SYMBOLS

	Batch code
	Consult instruction for use
	<i>In vitro</i> diagnostic medical device
	Manufacturer
	Temperature limitation
	Use by
	Catalogue number
	Content
	Reconstitute with x mL
	Control Low
	Control High
	Add exactly 0.5 mL of distilled or deionised water
	Keep at room temperature for 10 minutes before use
	These products are stable x days at x °C
	These products should be immediately and tightly capped to prevent contamination and evaporation
	Modification from previous version

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4 Weykamp, C., *et al.*, *A review of the Challenge in Measuring Hemoglobin A1c*, *J Diabetes Sci Technol*, (2009), 3, 439.

5 Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263

6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).

7. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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
English - EN

❖INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of VitalScientific reagents, used for HDL or LDL cholesterol quantitative determination. This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum.
- The HDL and LDL concentrations of this calibrator are lot-specific.

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- VitalScientific reagent, used for HDL or LDL cholesterol quantitative determination (Ref. CHDL-XXXX or CLDL-XXXX)
- General Laboratory equipment.

TRACEABILITY

HDL LDL CALIBRATOR is traceable to the reference method approved by the CDC (Centers for Disease Control and Prevention).

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.
- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved or CE compliant. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.

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PREPARATION

- Carefully open the vial, avoiding the loss of lyophilized.
- Add exactly 1 mL of distilled or deionised water.
- Carefully close the vial and dissolve the contents completely by occasional gentle stirring avoiding the formation of foam.
- Keep at room temperature for 30 minutes before use.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particles would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

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After reconstitution :

- These products are stable 2 weeks at 2-8 °C or 10 weeks at -20°C.
- These products should be immediately and tightly capped to prevent contamination and evaporation.

❖PROCEDURE

To use HDL LDL CALIBRATOR, follow the procedure described in the instructions for use of the VitalScientific reagent used.

❖LIMITATIONS

Using HDL LDL CALIBRATOR has been validated with the VitalScientific systems (Analyzers and reagents used for HDL or LDL cholesterol quantitative determination (Ref. CHDL-XXXX or CLDL-XXXX)). Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit. For users of Selectra instruments allowing data import for tests, calibrators and controls, use the corresponding barcode available on the value sheet.

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