**Microalbumin Kit
(Immunoturbidimetric Assay Method)****Order Information**

Cat. No.	Package size
105-006173-00	R1: 1×30 mL + R2: 1×7 mL
105-002242-00	R1: 3×36 mL + R2: 1×20 mL
105-002243-00	R1: 3×50 mL + R2: 3×11 mL
105-006174-00	R1: 1×30 mL + R2: 1×7 mL+Cal: 5×1 mL
105-004503-00	R1: 3×36 mL + R2: 1×20 mL+Cal: 5×1 mL
105-004504-00	R1: 3×50 mL + R2: 3×11 mL+Cal: 5×1 mL

Intended Purpose

In vitro test for the quantitative determination of Microalbumin (MALB) concentration in human urine and cerebrospinal fluid (CSF) on Mindray BS series chemistry analyzers. It is intended to be used for diagnosing and monitoring the therapeutic effect of albuminuria, and it could also be used for aiding to diagnose renal diseases, blood brain barrier disorders.

Summary¹⁻⁴

Albumin is the only non-glycosylated plasma protein with 66.3kDa, is the critical component in plasma, and is the most important binding and transport protein. It is synthesized in the liver and has the 19 days half-life. Determination of albumin is clinically relevant in serum, CSF and urine.

In urine, it is known as microalbumin. A healthy kidney infiltrated and excreted unnecessary albumin which sustain in blood as valuable protein. Microalbumin (MALB) is defined as a condition characterized by urinary albumin excretion above 20 mg/L in normal adults. Investigations have established a range of 20 mg/L-200 mg/L as diagnostic for MALB.

The presence of MALB in the urine is a reliable and sensitive evidence to diagnose earliest stage of nephropathy. MALB can be considered as an indicator of deteriorating renal function in diabetic subjects. It also may beneficial to control diabetes treatment program.

In CSF of healthy individuals, albumin is completely originated from the serum, and the change of ALB ratio is a good maker of blood-brain barrier. The albumin CSF/serum ratio is dependent on age, especially during infancy and old age. The albumin CSF/serum ratio is a reliably quantified aid and increased in blood brain barrier disorders, such as viral meningitis, spinal cord tumor. By simultaneously determining IgG in CSF and serum while taking into account the individual albumin ratios, it is possible to differentiate between IgG originating from the blood and CNS-synthesized immunoglobulin.

Assay Principle

Immunoturbidimetric Assay Method

Anti-human albumin antibody + albumin \rightleftharpoons Immunocomplex (agglutination)

When the sample is added to the R1 and R2, turbidity changes caused by antigen-antibody reaction. The turbidity is proportional to albumin quantity in the sample. The actual concentration is then determined by a calibration curve prepared from calibrators of known concentration.

Reagents Components

R1:	Phosphate Buffer	3.030 g/L
	Sodium azide	0.095%
R2:	Phosphate Buffer	2.420 g/L
	Goat anti-human MALB antibody	0.500 g/L
	Sodium azide	0.095%
	Phosphate Buffer	3.030 g/L
Calibrator:	Human ALB	0.025 g/L~0.400 g/L
	Sodium azide	0.095%

MALB Calibrator is a multi-point calibrator, which contains 5 different levels: a, b, c, d, e. The calibrator values are lot-specific with the matched models listed in the value sheet. For CSF MALB test, please scan reagent bar code on the extra label in the box for bar-coded reagent loading.

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8 and protected from light.

On board in use, the reagents are stable for 30 days when refrigerated on the analyzer. Once opened in use, the calibrator is stable for 30 days when refrigerated at 2-8 and protected from light.

Contamination must be avoided.

Do not freeze the reagent.

Specimen collection and preparation**■ Specimen types**

Urine and cerebrospinal fluid (CSF) are suitable for samples.

■ Preparation for Analysis

1. Use the suitable tubes or collection containers and follow the instruction of the manufacturer. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases, avoid effect of the materials of the tubes or other collection containers.

- Specimens should be tested as soon as possible after sample collection and pre-analytical treatment.
- CSF sample should be better to test in one hour.

■ Sample Stability

- Urine⁵: 7 days at 15-25°C
 1 month at 2-8°C
 6 months at (-25)-(-15)°C
 CSF: 3 days at 2-8°C
 6 months at (-25)-(-15)°C

For longer storage periods, samples should be frozen at (-20°C)⁶. Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

Reagent Preparation

R1 and R2 are ready to use.

Calibrator is ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

- General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
- Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.
- Mindray BS series chemistry analyzers and General laboratory equipment.

Assay procedure

Parameters Item	BS-800 chemistry analyzers	
Sample types	Urine	CSF: Dilution First Sample for dilution: 10 µL Diluent (saline): 110 µL
Assay type	Endpoint	
Wavelength	340 nm	
Reaction direction	Increase	
R1	187.5 µL	
Sample or Calibrator	12 µL	
Mix, incubate at 37°C for 4 min, read the absorbance A1, then add:		
R2	30 µL	
Mix thoroughly, incubate at 37°C for 5 min, then read the absorbance A2, Then calculate $\Delta A=(A2-A1)$		

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

Calibration

1. It is recommended to use the Mindray Calibrator (within the kit or MALB Calibrator: 105-002315-00 or other suitable calibrators) and 9 g/L NaCl (saline) for multi-point calibration. The calibrator values are assigned by Mindray standard transfer procedure and routine method. As the values of calibrator are those of urine and CSF, the calibrator for CSF is ready for use. The traceability process is based on ISO 17511⁷, the analyte (measurand) in this calibrator is traceable to the reference material ERM DA470K.

2. Calibration frequency

Calibration is stable for approximately 14 days on BS-800 chemistry analyzers. The calibration stability may vary on different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.

Recalibration may be necessary when the following occur:

- As changed reagent lot.
- As required following quality control procedures or out of control.
- As executes specific maintenance or troubleshooting procedure of chemistry analyzers.

3. The calibrator values are lot-specific with the matched models listed in the value sheet.

Quality control

1. It is recommended to use the Mindray Control (MALB Control: 105-002316-00 or other suitable control) to verify the performance of the measurement procedure; other suitable control material can be used in addition.

2. Single level of control material is recommended to analyze each batch of samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual.

3. Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control doesn't recover within the acceptable tolerances.

Calculation

$\Delta A = (A_2 - A_1)$ sample or calibrator.

$C_{\text{sample}} = (\Delta A_{\text{sample}} / \Delta A_{\text{calibrator}}) \times C_{\text{calibrator}}$.

The BS series chemistry analyzer detects the change of absorbance(ΔA) and calculates the MALB concentration of each sample automatically after calibration.

Conversion factor: $\text{mg/L} \times 0.0152 = \mu\text{mol/L}$

Dilution

If the value of urine sample exceeds 300 mg/L, the sample should be diluted with 9 g/L NaCl solution (saline) (e.g. 1+ 29) and rerun; the result should be multiplied by 30.

If the value of CSF sample exceeds 3600 mg/L, the sample should be diluted with 9 g/L NaCl solution (saline) (e.g. 1+ 3) and rerun; the result should be multiplied by 4.

Expected values²

Sample Type	Units
Urine	$\leq 20 \text{ mg/L}$
Cerebrospinal fluid	$\leq 350 \text{ mg/L}$

The expected value is provided from reference, and Mindray has verified it by 140 urine samples and 246 CSF samples of people from southern China. Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

Performance Characteristics

■ Analytical Sensitivity

The Microalbumin kit has an analytical sensitivity of 4 mg/L (urine) and 48 mg/L (CSF) on BS-2000. Analytical sensitivity is defined as the lowest concentration of analyte that can be distinguished from a sample that contains no analyte. It is calculated as the value lying 3 standard deviations above that of the mean from 20 replicates of an analyte-free sample.

■ Measuring range

The Mindray BS series systems provide the following linearity range:

Sample Type	Units
Urine	4-300 mg/L

A high concentration MALB urine sample (approximately 300 mg/L) is mixed with a low concentration sample ($<4 \text{ mg/L}$) at different ratios, generating a series of dilutions. The MALB concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.990$. The reportable range is 4-9000 mg/L.

Sample Type	Units
CSF	48-3600 mg/L

A high concentration MALB CSF sample (approximately 3600 mg/L) is mixed with a low concentration sample (<48 mg/L) at different ratios, generating a series of dilutions. The MALB concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.990$. The reportable range is 48-14400 mg/L.

■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3⁸, each sample was assayed 2 times per run, 2 runs per day, a total of 20 days. The precision data of urine control and samples on BS-800, CSF control and samples on BS-2000 are summarized below*.

Urine:

Specimen Type (N=80)	Mean (mg/L)	Repeatability		Within-Lab	
		SD (mg/L)	CV %	SD (mg/L)	CV %
Control	20.86	0.27	1.31	0.66	3.18
Urine 1	180.90	1.30	0.72	1.68	0.93

CSF:

Specimen Type (N=80)	Mean (mg/L)	Repeatability		Within-Lab	
		SD (mg/L)	CV %	SD (mg/L)	CV %
Control	225.23	2.53	1.12	5.21	2.32
CSF 1	129.98	2.35	1.81	4.29	3.30
CSF 2	326.31	3.71	1.14	5.77	1.77
CSF 3	454.11	6.60	1.45	9.91	2.18

*Representative data, results in different instruments, laboratories may vary.

■ Analytical Specificity

The samples with different concentration interfering substance were prepared by addition of interferent to human sample pools, and recovers are within $\pm 10\%$ or ± 1.1 mg/L for urine samples, within $\pm 10\%$ or ± 45.0 mg/L for CSF samples, of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The interference studies data of urine samples on BS-800, CSF samples on BS-2000 are summarized below.

Urine:

Interfering Substance	Interferent Concentration (mg/dL)	Analyte Concentration (mg/L)	Relative Deviation (%)*
Ascorbic acid	30	285.45	-0.98
Hemoglobin	300	285.26	-1.38
Bilirubin	66	153.01	-3.98
Glucose	2000	286.92	-0.93
Acetone	350	294.46	+0.57
Creatinine	500	289.75	-0.07
Urea	4200	289.79	-3.77
Urobilinogen	20	285.13	-0.28

CSF:

Interfering Substance	Interferent Concentration (mg/dL)	Analyte Concentration (mg/L)	Relative Deviation (%)*
Hemoglobin	1000	290.74	+5.26
Conjugated bilirubin	40	290.39	+2.29

*Representative data, results in different instruments, laboratories may vary.

■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP09-A3⁹. The Mindray System (Mindray BS-800/Mindray MALB Reagent) (y) was compared with comparison system (HITACHI 7180/Roche MALB Reagent) (x) using the same urine specimens. The statistical data obtained by linear regression are shown in the table below *

Regression Fit	Correlation Coefficient (r)	Sample (N)	Concentration Range (mg/L)
$y=1.045x-0.7978$	0.9990	100	5.65-278.76

*Representative data, results in different instruments, laboratories may vary.

The Mindray System (Mindray BS-2000/Mindray MALB Reagent) (y) was compared with comparison system (Siemens BN II/Siemens MALB Reagent) (x) using the same CSF specimens. The statistical data obtained by linear regression are shown in the table below *:

Regression Fit	Correlation Coefficient (r)	Sample (N)	Concentration Range (mg/L)
$y=0.9671x-5.2243$	0.9982	115	39.45-3255.94

*Representative data, results in different instruments, laboratories may vary.

Result interpretation

The results could be affected by drugs, disease, or endogenous substances¹⁰. When the reaction curve is abnormal, it is recommended to retest and check the result.

Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The reagents avoid direct exposure to sunlight and freezing. The results can't be assured when stored at inappropriate condition.
4. If unintentionally opened before used, store the reagents tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix reagents with different lots and bottles.
Do not use the reagents beyond the expiration date and the in-use date.
Do not mix new opened reagents with in-use reagents.
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the reagent accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. The calibrator was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen.
13. All human material should be considered potentially infectious.
14. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.

15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

References

1. Gutman AB, Moore DH, Gutman EB, et al. Fractionation of serum proteins in hyperproteinemia, with special reference to multiple myeloma. *J. Clin. Invest.* 1941, 20:765-783.
2. Thomas L. *Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results*. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 382-400, 1308-1326.
3. Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. Elsevier Saunders 2012;527-530.
4. Reiber, H. Aktuelle Methoden der Liquoranalytik. *Lab med* 1988; 12:101-109.
5. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissler H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:46pp.
6. CLSI. *Urinalysis; Approved Guideline-Third Edition*. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
7. ISO 17511:2003. *In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials*.
8. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
10. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press; 2000:3-28,3-31.

Graphical symbols

In Vitro Diagnostic
medical device

Unique device
identifier

European
Conformity

Consult Instructions
For use

Use-by
date

Authorized representative
in the European

Batch Code

Temperature
limit

Manufacture

Catalogue
number

Biological

Keep away from sunlight
Indicates a medical device that
needs

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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1. Gutman AB, Moore DH, Gutman EB, et al. Fractionation of serum proteins in hyperproteinemia, with special reference to multiple myeloma. J. Clin. Invest, 1941, 20:765-783.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 382-400, 1308-1326.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;527-530.
4. Reiber, H. Aktuelle Methoden der Liquoranalytik. Lab med 1988; 12:101-109.
5. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of

- Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:46pp.
6. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
 7. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
 8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
 9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
 10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-28,3-31.

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Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Referências

- 1.Gutman AB, Moore DH, Gutman EB, et al. Fractionation of serum proteins in hyperproteinemia, with special reference to multiple myeloma. J. Clin. Invest, 1941, 20:765-783.
- 2.Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 382-400, 1308-1326.
- 3.Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;527-530.
- 4.Reiber, H. Aktuelle Methoden der Liquoranalytik. Lab med 1988; 12:101-109.
- 5.Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:46pp.
- 6.CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
- 7.ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
- 8.CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- 9.CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-28,3-31.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Morada: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Referencias

1. Gutman AB, Moore DH, Gutman EB, et al. Fractionation of serum proteins in hyperproteinemia, with special reference to multiple myeloma. *J. Clin. Invest.* 1941, 20:765-783.
2. Thomas L. *Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results*. 1.ª edición, Frankfurt: TH-Books Verlagsgesellschaft 1998; 382-400, 1308-1326.
3. Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5.ª edición, Elsevier Saunders 2012; 527-530.
4. Reiber, H. Aktuelle Methoden der Liquoranalytik. *Lab med* 1988; 12:101-109.
5. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1, revisión 2: pág. 46.
6. CLSI. *Urinalysis; directriz aprobada, tercera edición*. Documento GP16-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
7. ISO 17511:2003. *In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials*.
8. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. documento EP05-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition*. Documento EP09-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

10.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5.ª edición
Washington, DC: AACC Press; 2000: 3-28, 3-31.

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Bibliografia

1. Gutman AB, Moore DH, Gutman EB, et al. Fractionation of serum proteins in hyperproteinemia, with special reference to multiple myeloma. *J. Clin. Invest.* 1941, 20:765-783.
2. Thomas L. *Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results*. 1st ed. Francoforte: TH-Books Verlagsgesellschaft 1998; 382-400, 1308-1326.
3. Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. Elsevier Saunders 2012;527-530.
4. Reiber, H. Aktuelle Methoden der Liquoranalytik. *Lab med* 1988; 12:101-109.
5. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:46pp.
6. CLSI. *Urinalysis; Approved Guideline-Third Edition*. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
7. ISO 17511:2003. *In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials*.
8. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

10.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed.
Washington, DC: AACC Press; 2000:3-28,3-31.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech
Industrial Park, Nanshan, Shenzhen, 518057, Repubblica popolare
cinese

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Referanslar

1. Gutman AB, Moore DH, Gutman EB, et al. Fractionation of serum proteins in hyperproteinemia, with special reference to multiple myeloma. J. Clin. Invest, 1941, 20:765-783.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 382-400, 1308-1326.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;527-530.
4. Reiber, H. Aktuelle Methoden der Liquoranalytik. Lab med 1988; 12:101-109.
5. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:46pp.
6. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
7. ISO 17511:2003. In vitro diagnostic medical devices- Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and Kontrol materials.
8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-28,3-31.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. Çin

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- Reiber, H. Aktuelle Methoden der Liquoranalytik. Lab med 1988; 12:101-109.
- Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:46pp.
- CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
- ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
- CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-28,3-31.

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© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés.
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)