


TOPAZ

Mobile DR System

Operation Manual



DRGEM Corporation
7F E-B/D Gwangmyeong Techno-Park, 60 Haan-ro,
Gwangmyeong-si, Gyeonggi-do, Rep. of Korea, 14322
TEL: +82-2-869-8566, FAX: +82-2-869-8567

D/N: RMD1311-009, Rev.7




(This page intentionally left blank)



REVISION HISTORY

Revision Number	Date	Description
0	NOV 25, 2013	First Edition
1	MAR 17, 2015	GUI modified, Configuration section added
2	NOV. 02, 2015	Changes of section 1.4, 3 and 4.
3	JAN 28, 2016	Corrections of section 1.3.
4	AUG 24, 2016	Add the Calibration Step of Varian 4336W V4 Detector Add the Service mode
5	JUL 28, 2017	Add Engineer mode 4.4.4, Changes of section 4.4.3 And 4.5.3, Changes SW images
6	MAR 23, 2018	Change of Software Image,
7	NOV 29, 2018	Change the Standard(EMC 4 , Safety 3.1)



ADVISORY SYMBOLS

The following advisory symbols are used throughout this manual.
Their application and meaning are described below.

WARNING

Warning symbol is used to indicate a potential hazard for operators and service personnel that can lead to serious injury, death or radiation exposure.

CAUTION

Caution symbol is used to indicate a potential hazard for operators and service personnel that can lead to injury or damage of equipment.

NOTE

Note symbol is used to indicate important information needed for proper use and correct operation of equipment.

NOTE

Keep this Operation Manual with the equipment at all times, and review the important information whenever required.

Copyright DRGEM Corporation. All rights reserved.

This document is the property of DRGEM Corporation and contains confidential and proprietary information owned by DRGEM Corporation. Any unauthorized copying, use or disclosure of it without the prior written permission of DRGEM Corporation is strictly prohibited.

NOTE

Consult Accompanying Documents - As Applicable



TABLE OF CONTENTS

1. INTRODUCTION	11
1.1 FEATURES	11
1.1.1 STANDARD	12
1.1.2 OPTION	12
1.2 SAFETY INFORMATION	13
1.2.1 STATEMENT OF LIABILITY	13
1.2.2 SYMBOL DEFINITIONS	14
1.2.3 SAFETY GUIDELINES	17
1.2.4 X-RAY PROTECTION	21
1.2.5 RADIATION SAFETY	24
1.2.6 MANUFACTURER'S RESPONSIBILITY	27
1.2.7 MONITORING PERSONNEL	28
1.2.8 RADIATION PROTECTION SURVEY	29
1.3 SPECIFICATIONS	30
1.4 PART DESCRIPTION	35
1.5 APPLICABLE STANDARDS	37
1.6 CUSTOMER SUPPORT	40
2. OPERATION PROCEDURE	41
2.1 TURN ON THE MOBILE X-RAY SYSTEM	42
2.2 OPERATION	42
2.3 TURN OFF THE MOBILE X-RAY SYSTEM	45
3 GRAPHIC USER INTERFACE	47
4 IMAGING SOFTWARE	49
4.1 MAIN GUI	49
4.1.1 LOGIN	49
4.1.2 MAIN MENU	50
4.1.3 INDICATORS	51
4.1.3.1 IMAGE SENDING INDICATOR	51
4.1.3.2 X-RAY INDICATOR	51
4.1.3.3 DETECTOR STATUS INDICATOR	52
4.1.3.4 DETECTOR THE CABLE INDICATOR	52
4.1.3.5 DETECTOR WIRELESS STATUS INDICATOR	53
4.1.3.6 DETECTOR BATTERY STATUS INDICATOR	54



4.1.3.7 HDD STATUS AND STATUS BAR	54
4.2 WORKLIST	55
4.2.1 SEARCH (QUERY FROM WORKLIST SERVER)	55
4.2.1.2 MODALITY FILTER	56
4.2.1.3 SCHEDULE (TIME/DATE) FILTER	56
4.2.1.4 SEARCH BUTTON	57
4.2.1.5 STOP BUTTON	57
4.2.2 EMERGENCY	58
4.2.3 NEW	59
4.2.4 OPEN	59
4.2.5 SELECTION CHECK BOX	60
4.2.6 SELECT ALL	61
4.2.7 SORT	62
4.2.8 COLUMN RESIZE	63
4.3 STUDY LIST	64
4.3.1 SEARCH (QUARY FROM DATABASE)	64
4.3.1.1 EDITOR FILTER	65
4.3.1.2 MODALITY FILTER	65
4.3.1.3 SCHEDULE (TIME/DATA) FILTER	65
4.3.1.4 SEARCH BUTTON	66
4.3.1.5 CLEAR BUTTON	66
4.3.2 EMERGENCY	67
4.3.3 NEW	68
4.3.4 MODIFY STUDY	69
4.3.5 DELETE STUDY	70
4.3.6 TRANSFER	70
4.3.7 EXPORT	70
4.3.8 OPEN	71
4.3.9 SELECTION CHECK BOX	71
4.3.10 SELECT ALL	72
4.3.11 MULTI OPEN	73
4.3.12 SORT	74
4.3.13 COLUMN RESIZE	75
4.3.14 STUDY LIST STATISTICS	76
4.4 PROCEDURE	77
4.4.1 STUDY TAB	78
4.4.2 TOOLS	79
4.4.2.1 CURSOR	79





4.4.5.1 INDICATORS	94
4.4.5.2 RESET BUTTON.....	95
4.4.5.3 DR BUTTON.....	95
4.4.5.4 X-RAY CONDITION CONTROL BUTTON	95
4.4.5.5 ms/mAs SELECT BUTTON	95
4.4.5.6 SAVE	96
4.4.5.7 PATIENT BODY SIZE SELECTION BUTTONS.....	96
4.5 CONFIGURATION	97
4.5.1 LOGIN	97
4.5.1 ACCOUNT SETTING.....	97
4.5.1.1 ADD	98
4.5.1.2 REMOVE	98
4.5.2 NETWORK SETTING.....	99
4.5.2.1 WORKLIST SETTING	99
4.5.2.2 PACS SETTING	100
4.5.2.3 MPPS SETTING	101
4.5.2.4 STORAGE COMMITMENT SETTING.....	102
4.5.2.5 DICOM PRINTER SETTING.....	103
4.5.3 GUI SETTING.....	104
4.5.2.6 BASE GUI SETTING	104
4.5.2.7 MARKER SETTING.....	105
4.5.2.8 ROI SETTING	106
4.5.2.9 OVERLAY SETTING	107
4.6 SERVICE MODE	108
4.6.1 LOGIN.....	108
4.6.2 DRIVING SENSITIVITY	109
4.6.3 MAX. DRIVING SPEED	109
4.6.4 MAX. INCH MOVER SPEED.....	109
4.6.5 HANDLEBAR INCHMOVE	109
4.6.6 DRIVING BEEP ALARM.....	110
4.6.6 BRAKE RELEASE	110
4.6.7 PROGRAM VERSION	110
4.6.8 DRIVING SLOPE CALIBRATION.....	111
4.6.9 LOAD-CELL CALIBRATION.....	111
APPENDIX A. EXPOSURE TABLE	113
APPENDIX B. EXPOSURE INDEX	117
APPENDIX C. GENERATOR SETUP	121



APPENDIX D. VAREX DETECTOR CALIBRATION STEP	127
D1. 4336W V4 DETECTOR.....	127
APPENDIX E. CONFIGURATION SOFTWARE	137
E1. XGRADMAXCONFIG LOGIN	138
E2. ACCOUNT SETTING.....	139
E3. NETWORK SETTING	141
E3.1 WORKLIST SETTING.....	142
E3.2 PACS SETTING	145
E3.3 MPPS SETTING.....	147
E3.4 STORAGE COMMITMENT SETTING	150
E3.5 OPTIONS SETTING.....	153



1.3 SPECIFICATIONS


Item		Specification	
		TOPAZ-32D	TOPAZ-40D
X-ray Generator	Nominal Output	32KW	40KW
	mA Range	Max. 400mA	Max. 500mA
	KV Range	40 ~ 125 KV (option: 150KV)	40 ~ 125 KV (option: 150KV)
	mAs Range	0,1 ~ 500mAs	0,1 ~ 500mAs
	Operation Duty Cycle	100kV, 320mA, 100ms / Rest time : 1min	100kV, 400mA, 100ms / Rest time : 1min
	Type	High Frequency	
System Design & Transport	Drive Type	Motor Driven(Rear-wheel drive)	
	Speed of Movement	Max. 5Km/h	
	Movement Brake/ Safety	Deadman type Handle, Front Safety Bumper, Spring Caster	
	Maximum Step Height	50mm	
	Maximum ramp gradient	10 degrees	
	Main power supply	100 ~ 240VAC, 9 – 4 A	
	Cable Length(Mains, retractable)	3m	
	Driving distance after fully charged	30km	
	Inch mover	Max. 50mm/s	
X-Ray Tube Support	Colum rotation range	± 325 degrees	
	Tube axis rotation range	-30 ~ +90 degrees	
	Tube support axis rotation range	± 180 degrees	
	Collimator rotation range	±120 (option ±180) degrees	
	Max./Min. X-ray focal spot height	1,390mm(Option : 1200mm)	
	Max. horizontal extension	560mm	
Option & Accessories	Remote Controller	Ready-Exposure, Lamp On/Off	
	DAP(Dose Area Product) Recording	Thermal Printer	



● X-ray Tube

Tube Model	E7239X	DXT-8M	E7242X	DXT-11M
Manufacturer	CANON	DRGEM	CANON	DRGEM
Focal Spot Size	1.0/2.0mm	1.0/2.0mm	0.6/1.5mm	0.6/1.5mm
Rating(0.1s)	22.5/47kW@60Hz	22.5/47kW@60Hz	18/50kW@60Hz	18/50kW@60Hz
Max. Anode HU	140kHU(100kJ)	140kHU(100kJ)	200kHU(142kJ)	200kHU(142kJ)
Target Angle	16°	16°	14°	14°
Max. kV	125kV	125kV	125kV	125kV
Weight	16kg(35.3lbs)	16kg(35.3lbs)	16kg(35.3lbs)	16kg(35.3lbs)
Inherent Filtration	0.9mmAl/75kV, IEC60522			
Additional Filtration	1.0mmAl			
Half Value Layer	More than 2.9mmAl eq. at 80kVp			
Leakage Radiation	Less than 100mR/hr			

Tube Model	E7299X	E7876X	E7884X
Manufacturer	CANON	CANON	CANON
Focal Spot Size	0.3/1.0mm	0.6/1.2mm	0.6/1.2mm
Rating(0.1s)	3.7/39kW@60Hz	22/54kW@60Hz	22/54kW@60Hz
Max. Anode HU	150kHU(111kJ)	230kHU(163kJ)	300kHU(210kJ)
Target Angle	12°	12°	12°
Max. kV	150kV	150kV	150kV
Weight	16kg(35.3lbs)	16kg(55.1lbs)	16kg(35.3lbs)
Inherent Filtration	1.3mmAl/75kV	1.3mmAl/75kV	0.9mmAl/75kV
Additional Filtration	1.0mmAl		
Half Value Layer	More than 2.9mmAl eq. at 80kVp		
Leakage Radiation	Less than 100mR/hr		



Tube Model	DXT-10M	DXT-12M
Manufacturer	DRGEM	DRGEM
Focal Spot Size	0.6/1.2mm	0.6/1.2mm
Rating(0.1s)	17/48kW@60Hz	22/54kW@60Hz
Max. Anode HU	150kHU(111kJ)	300kHU(210kJ)
Target Angle	12°	12°
Max. kV	125kV	150kV
Weight	16kg(35.3lbs)	16kg(35.3lbs)
Inherent Filtration	1.0mmAl/75kV	1.0mmAl/75kV
Additional Filtration	1.0mmAl	
Half Value Layer	More than 2.9mmAl eq. at 80kVp	
Leakage Radiation	Less than 100mR/hr	


* Total filtration including X-ray tube assembly and collimator will be matched by appropriate additional filters to within the range from 2.9 to 3.3mmAl. eq.

- DAP meter (Option)

DAP Resolution	0.01 μGym^2
Interface	RS485
Active area	115 x 115mm
Dimension	158 x 134.5 x 17mm

- Digital Flat Panel Detector

Manufacturer / Model	VAREX / 4336W V4	
Active Pixel Area / Matrix	14 x 17 inch	
Pixel Pitch	139um	
Limiting Resolution	3.6 lp/mm	
Scintillator	CsI,	DRZ+
Energy Range	40 – 150kVp	
A/D Conversion	16-bits	
Interface	Gigabit Ethernet	
Weight(including Battery)	3.8kg (8.4lbs)	3.6kg (7.9lbs)



● Collimator

Model / Maker	R108/R108F	DXC-RML/DXC-RMH
Manufacturer	RALCO	DRGEM
Control	Manual with 30sec. lamp timer	
Field Shape	Rectangular	
Max. Field Size	More than 43x43cm(17x17inch) at 100cm SID	
Leakage Radiation	Less than 100mR/hr	Less than 40 mR/h
Max. kVp shield	150kV	
Inherent Filtration	Min. 2.0mmAl eq.	
Luminosity	Over 160LUX at 100cm SID (Typ. 250LUX)	Over 160LUX at 1cm SID
Light source	Single LED	LED and Halogen
Standard	Tape measure, Rotating flange	Rotating flange with fixing knob
Option	Laser line, Near port moving shutters(Only R108F), Manual rotating filters(Only R108F)	Line laser+shutter, Measure tape Near port moving shutters, Additional filter Mounting flange mechanical detent Accessory guides spring DAP rail
Electrical Rating	20-30VAC, 30VA, 50/60Hz	Halogen type - 24 V DC/AC - 50~60Hz 160VA LED type - 12~45V DC 35VA / 20~30VAC 35VA - 50~60Hz
Dimension / Weight	223(W) x 246(D) x 140(H) mm / 6.6kg(14.6lb)	196(W) x 250(D) x 171(H) mm / 7.1kg(15.6lb)



- Image Workstation

Item	Specification
CPU	Intel(R) Core i3-6100(3.70GHz,3MB Cache) or higher
Memory	4GB (1x4GB) Non-ECC DDR3L 1600MHz SDRAM or higher
Display	Intel® HD Graphics 530 (Integrated HD graphics)
Storage	128 GB x 1 SSD (Samsung) + 500GB x 1 SATA HDD or higher
Maker	DELL International
Touch Monitor	17 inch Color LCD, IR Infrared Ray, Display resolution: 1280 x 1024 pixels

OPERATING ENVIRONMENT

Ambient temperature range	10 to 35 °C (50 to 95 °F)
Relative humidity	30% to 75%, non-condensing
Atmospheric pressure range	700 hPa to 1060 hPa
Altitude Limit	This product is rated to operate at an altitude ≤3000m

TRANSPORT AND STORAGE ENVIRONMENT

Ambient temperature range	-10 to 70 °C (14 to 158 °F)
Relative humidity	10 to 90%, non-condensing.
Atmospheric pressure range	500 hPa to 1060 hPa

WARNING

Do not operate this system except in accordance with information included in this section, and any additional information provided by the manufacturer and / or competent safety authorities.

