

VISION® VIRAL TRANSPORT MEDIUM

REF VVTM001

PACKAGING

100T/Set

INTENDED USE

VISION® VIRAL TRANSPORT MEDIUM is produced in accordance with the CDC formulation. It provides active and safe transportation of samples with its enriched buffer solution, and antibiotic and antifungal agents help the sample to remain pure.

DESCRIPTION

VISION® VIRAL TRANSPORT MEDIUM collects and transports viruses. The antibiotic and antifungal agents in it ensure that the sample remains pure.

The liquid in this tube preserves the viral pathogens found in clinical samples. It contains bacterial and fungal antibiotics.

It is suitable for nucleic acid extraction.

Positive samples taken with VISION® VIRAL TRANSPORT MEDIUM can be stored for at least two years under appropriate conditions.

KIT CONTENT

For 100 pieces,

Content	Quantity
VISION® VIRAL TRANSPORT MEDIUM	100
Swab	100

**** Each tube contains 1.5 mL of liquid.**

**** Swabs are suitable for vaginal/cervical sampling.**

PROTOCOL

CLSI MM13-A steps are followed for sample storage and preservation.
















The swab with the sample is put into the VISION® VIRAL TRANSPORT MEDIUM tube.

The swab is completely synthetic (dacron/rayon/flocked polyester) and sterile.


SUGGESTIONS

After sampling, tubes should be stored at +4° C.

SYMBOLS

	Manufacturer		In vitro Diagnostics medical device		Consult the operating instructions
	Best before date		Temperature Limit (2-30°C)		Authorized representative in the European Community
	Batch code		Do not re-use		Keep dry
	Catalog number		Do not use if the package is damage		Protect from sunlight
	Manufacturing date		CE marking		For <n> experiments Contains sufficient quantity

CONTACT

 Vision Biyoteknoloji Araştırma Geliştirme Laboratuvar Sistemleri Sanayi ve Ticaret Limited Şirketi

Address: Yukarı Dudullu Mah. Alemdağ Cad. No:742/A, 34775

Ümraniye/ İstanbul/Türkiye

Phone: +90 2167840298

Email: info@visionbiotechnology.com

Web: www.visionbiotechnology.com