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Telescopes and hr Telescopes

Telescoape si telescoape hr

Anleitung

Manual

Manual

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Telescopes and hr Telescopes

Telescoape si telescoape hr

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning and sterilization of the telescopes. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

The recommended testing and preparation methods for telescopes are described in this manual.

KARL STORZ telescopes are compatible with all light sources and fiber optic light cables sold by KARL STORZ.

Warning: KARL STORZ telescopes are provided unsterilized, and must be cleaned, disinfected and/or sterilized prior to the initial use and before each subsequent use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply, design, features and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations and descriptions in this manual.

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




HAMOU Optik

DCI-Anschluß

Adapter for fiber light cable, unscrewable, suitable for KARL STORZ, Olympus/Winter and Ibe

Adapter for fiber light cable, unscrewable, suitable for KARL STORZ, Wolf (9 mm plug-type connection)

Light inlet piece, without adapter suitable for ACM, Eder and Zimmer (7.5 mm plug-type connection)

Color coding ring for direction of view of telescope

green = 0 degrees,

straight forward = 12 degrees/ 45
degrees

black = 30 degrees, forward-

red oblique = 70 degrees, lateral

yellow = 90 degrees, lateral

blue = 120 degrees,

white = 120 degrees,

retrospective

retrospective

Eyeiece lens

Jacket tube

"Autoclave" marking ring

Type designation and direction of view data.

The meanings of the letter coding after the type

designation are as follows:

A = 0/6 degrees,

straight forward = 30 degrees,

B = forward-oblique = 70 degrees,

C = lateral = 90 degrees,

D = lateral = 120 degrees,

E = retrospective

F = 12 degrees/ 45

degrees

W after the first letter = Wide angle

Object lens and distal exit of the fiber optic light

carrier.

C-mount coupling mechanism

Degree of focus for the eyepiece drive

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HOPKINS rod-lens telescopes are reusable
and consist of an eyepiece lens, a connection for
fiber optic light cables with screw-on adapters for
other makes of fiber optic light cables and a jacket
tube made of a non-corrosive material, which en-
closes the rod-lens system and a built-in fiber optic
light carrier.

HOPKINS telescopes are designed to minimize
patient trauma as much as possible. They are used
in hundreds of thousands of cases annually with
great benefit to patients.

However, these **HOPKINS** telescopes are more
delicate than the instruments normally used by
surgeons.

They must be handled with care, be inspected
upon removal from the patient for any sign of dam-
age, and be inspected before the next case to be
certain that they were properly handled between
cases and are functioning properly. Failure to prop-
erly maintain and use these instruments may result
in trauma to delicate tissues, perforation, infection
or breakage of the instruments which can produce
foreign bodies.

By following the instructions below, your
HOPKINS telescopes will provide you with a
long service life.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words Warning, Caution, and Note convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words Warning and Caution stand out more clearly, they are accompanied by a pictogram.

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

Note: A Note indicates special information about operating the product, or clarifies important information.

Warning: Before using the unit, read the following safety instructions carefully, to avoid putting your patients, personnel or yourself at risk.

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Note: Any damage to the unit resulting from incorrect operation is not covered by the guarantee.

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Review of warnings

Warning: Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Warning: The telescope is delivered unsterilized and must be cleaned, disinfected an/or sterilized prior to initial use and subsequent reuse.

Warning: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy, or if either no image or only parts of the image are visible.

Warning: Adjust the luminous intensity of the light source to a minimum to achieve optimum illumination of the operating field. The higher the intensity of the light source is set, the greater the heat emission at the tip of the telescope.

Warning: Never place the end of a fiber optic light cable or telescope on or under a sterile cover cloth. The intensity of the light may lead to burns on the patient and/or the sterile cover cloth.

Warning: Danger of blindness! Never look into the open end of a light cable

Warning: The named cleaning, disinfection and sterilization methods alone do not guarantee disinfection or sterilization. This can only be achieved if the personnel proceeds according to recognized and validated preparation methods.

Warning: Before any disinfection/sterilization, the telescopes must be carefully cleaned and any organic material, blood and cleaning solution must be removed. Only clean surfaces can be reached by the sterilization.

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Warning: Any deviation from our recommended preparation parameters must be validated by the user.

Warning: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and calibrated.

Warning: High-level disinfection for use on the patient is recommended only for telescopes which come into contact with mucous membrane or relatively small skin injuries.

High-level disinfection for use on the patient is not recommended for telescopes used in laparoscopic, arthroscopic, gynecological, reconstructive, or cosmetic surgery. Sterilization is recommended in these cases.

The disinfection mentioned here does not concern disinfection done for the purpose of protecting the personnel.

Caution: For miniature endoscopes with removable eyepieces, the eyepiece must either be closed with a stop plug or fixed firmly onto the endoscope before soaking in disinfectant solution.

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Normal use

Indication: KARL STORZ **HOPKINS**® and **HOPKINS II**® telescopes are used for visualization of the operating field during diagnostic and surgical endoscopic interventions.

HOPKINS telescopes for extracorporeal visualization (VITOM) are used exclusively outside the body at a predefined working distance. The minimum working distance is set using a corresponding spacer bar.

Contraindication: Use of KARL STORZ **HOPKINS**® and **HOPKINS II**® telescopes is contraindicated if such an application would endanger the patient in the opinion of a qualified physician.

Unauthorized conversions or modifications to the instrument are forbidden for reasons of safety.

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User qualification

KARL STORZ **HOPKINS**® telescopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique.

The information given in these instructions only serves to instruct in the correct handling, cleaning and sterilization of the telescopes. It is not intended as an introduction to the endoscopic technique.

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Safety precautions when using the telescopes

Telescopes must be used in unison with the medical rules and procedures of endoscopy recognized for this method.

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Unpacking the instrument

Check for missing items and evidence of shipping damage. File any complaints to KARL STORZ or the supplier immediately.

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Operating the unit

Warning: The telescope is unsterilized and must be cleaned, disinfected and/or sterilized prior to initial use and subsequent reuse.

Warning: The instructions and interface specifications for medical devices used in combination must be observed precisely.

Warning: If endoscopes are operated with powered endoscopic accessories, the patient leakage currents may accumulate. This is particularly important when using CF-type equipment. In such a case the other units being used should also be of the CF type in order to minimize total patient leakage current.

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
Inspection

Upon receipt, immediately inspect the telescope for signs of damage. This should also be done before and after every procedure.

Image waveguide

Check the optical end faces (distal tips and eyepiece) for scratches and surgical or cleaning residues. The surfaces should appear smooth and shiny.

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To test the image quality, the telescope should be turned slowly while looking through it.

If the image is completely or partially defective, either the rod-lens system in the jacket tube has broken or the lens in the eyepiece is defective. In either case, the telescope must be replaced.

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
A cloudy or spotty image may be caused by moisture or disinfectant residue on the optical end faces.

Note: When looking through **HOPKINS III** PDD telescopes, the image appears yellowish. This effect is normal, caused by the use of special filters.

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In such cases, carefully clean the optical end faces with a soft cloth or swab soaked in alcohol.

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Light guide
Hold the distal tip in the direction of a light source and check the number of dark points at the light guide connector. These dark points represent broken optical fibers in the light guide bundle. Individual broken optical fibers do not necessarily mean a noticeable reduction in image quality.

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However, above a certain amount of approx. 20 to 25%, light transmission and therefore image quality is significantly reduced. The telescope should no longer be used now.

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Warning: If damage is apparent, the telescope should not be used. Telescopes must be replaced if the image is cloudy, there is no image or only a partial image can be seen.

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Sheath
Inspect the entire sheath for straightness, dents and scratches. The sheath must not be bent. LUER LOCK stopcocks, the telescope lock and other movable parts should operate smoothly and be free of defects. Check plastic parts for discoloration, porosity and flexibility. If the parts deviate substantially from their new condition, they must be replaced. This especially applies to seals.

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Precautionary measures

As with all optical instruments, the **HOPKINS**® rod-lens telescope is a very delicate piece of endoscopic equipment. To ensure the long-lasting quality of the telescope, please follow the handling instructions described below:

To prevent breakage, telescopes should be supported by firmly grasping the eyepiece end. Never handle the distal end alone. Never bend the stainless steel shaft. This could lead to breaks or cracks in the rod-lens system.

Handle the telescope with care. Hard knocks, particularly at the distal end, may result in damage or cracks in the telescope and allow liquid, steam and other material to penetrate. If so, the damage will appear as a foggy and/or unclear visual area.

When cleaning, disinfecting and sterilizing telescopes, it is recommended that they be handled separately from other instruments. This prevents the **HOPKINS**® telescopes from being damaged by other instruments (see "Accessories").

Protective cases are available in various lengths for **HOPKINS**® telescopes, which are suitable for storage, transport and sterilization (see "Accessories").

Note: If damage is apparent, the telescope should not be used. Telescopes must be replaced if the image is cloudy, there is no image or only a partial image can be seen.

Caution: **HOPKINS**® telescopes without the „Autoclave“ labeling ring must not be steam sterilized. This can cause irreparable damage.

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
How to use accessories

Before application, make certain the telescope has been correctly cleaned, disinfected and/or sterilized.

Test for proper operation of the telescope and accessories before application.

Operate the video and light source equipment according to the specifications in the appropriate operating instructions.

Connection to camera system and light source

The camera heads usually have an integrated instrument coupling.

Squeeze the levers (and) of the instrument coupler together and insert the endoscope's eyepiece into the coupler. The coupler will latch automatically when the eyepiece is inserted.

The endoscope can be locked by pressing lever to lever .

Connect the light cable to the light source output and the endoscope (quarter turn of the thumb-screw on the screw base).

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Warning: Never place the end of a fiber optic light cable or telescope on or under a sterile cover cloth. The high intensity of the light can lead to burns on the patient and/or the sterile cover cloth.

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
Image sharpness and focal distance
Operate the video equipment according to the specifications in the appropriate operating instructions.

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White balance
The white balance resets the camera system's chrominance controls to conform to the color temperature of the light source in use.
When the camera is next used a new white balance only has to be performed if a light source with a different color temperature is used or if a light cable or endoscope with different diameter is used.

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To perform this readjustment, switch on the light source to be used and point the camera and endoscope at a white surface (the tip of a clean cotton swab, or similar).
The white view should fill at least 70% of the area of the screen. If an image size of 70% cannot be achieved, the endoscope should be detached from the camera head, a white surface illuminated with the endoscope, and this recorded with the camera head.

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


Connecting the sheath and telescope
Sheaths, sheath inserts and telescopes are provided with a 0 or line mark.
Open telescope lock by turning counterclockwise (up to the stop, line marks aligned with one another).
Make certain during assembly that the marks are aligned with one another (line mark and 0 mark are lined up). Insert the telescope into the sheath lock, and close the telescope lock by turning clockwise.

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General information

Warning: The listed cleaning, disinfection and sterilization methods alone do not guarantee disinfection or sterility. This can only be achieved if the staff follows recognized and validated preparation methods.

Warning: The instrument is delivered unsterilized and must therefore be thoroughly cleaned, disinfected and/or sterilized before using for the first time as well as before any further use.

Warning: High-level disinfection for use on the patient is recommended only for instruments which come into contact with mucous membrane or relatively small skin injuries. High-level disinfection for use on the patient is not recommended for instruments used in laparoscopic, arthroscopic, gynecological, reconstructive, or cosmetic surgery. Sterilization is recommended in these cases. The disinfection mentioned here does not concern disinfection done for the purpose of protecting the personnel.

Warning: For all cleaning of contaminated endoscopes and instruments, liquid-proof protective gloves must be worn. When using a compressed air gun, safety goggles are necessary. Aerosols should be removed using a suction device.

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
 Caution: Constantly switching between various preparation methods is especially hard on materials and must therefore be avoided. A one-time change of method, due to procurement of new cleaning equipment for example, is safe.

Caution: Telescopes must not be cleaned in an ultrasonic bath.

 Caution: Cleaning telescopes in an instrument washing machine is only possible using suitable procedures.

 Caution: Please consult the Water and steam quality section in order to prevent damage to the instrument set.

 Warning: Remember that sterilization is only then successful if

- the instruments have been freed of all organic materials and cleaning residues. Only clean surfaces can be reached by the sterilization.
- work proceeds according to recognized and validated methods.
- the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.

Observe the guidelines of the Employer's Liability Insurance Association and equal ranking organizations.

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Water quality requirements

Water with at least drinking water quality must be used to make cleaning and disinfectant solution. We would, however, point out that even water of drinking water quality may have too high a concentration of components in it which might damage the instruments. Corrosion, spotting and discolorations may then occur. It is therefore recommended that completely demineralized water be used.

The steam for sterilization must also satisfy the requirements of the standard in order to prevent damage and discolorations.

Please determine the water quality available in your building with the maintenance department, since resulting damage may often only be recognized years later and then no longer be repairable. KARL STORZ will gladly help with any questions you may have in this regard.

Required water quality (maximum limits)

Total water hardness	5 °d (0.9 mmol/l)
Chlorides	100 mg/l
Silicic acid	15 mg/l
Iron	0.05 mg/l
Manganese	0.05 mg/l
Copper	0.05 mg/l
	(mg/l = ppm)

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
Preparations for cleaning

Separate the fiber optic light cable and camera head from the telescope (refer to How to use accessories).

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To prevent contaminants from drying onto the telescopes, these should be wiped off immediately after use with cleaning/disinfectant solution or immersed in a container with cleaning/disinfectant solution.
To prevent damage, never place more than one telescope or telescopes and other instruments on top of one another. The best policy is to immerse them individually.

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Both adapters on the light inlet piece must be removed before cleaning, to allow easy access to the fiber surface.

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Note: It is recommended that a plastic basin (Art. no. 27645, 27646) be used to prevent scratching, particularly on the optical end faces. A plastic basin also prevents electrolytic corrosion, which may occur if different metals are found in one solution.

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For cleaning and disinfecting, the products especially recommended for endoscopes and approved by KARL STORZ are suitable (see Appendix).

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Caution: Never place telescopes in physiological saline solution.

Caution: Only EPA registered hard surface disinfectant solutions capable of high level of disinfection should be used. Solutions with a strongly acidic or alkaline pH are not recommended. Follow solution manufacturer instructions, paying close attention to the mixing ratio and immersion time. Prolonged exposure may result in damage. Observe the microbiological range of action of the chemicals used.

The solutions should be replaced regularly, since the risk of corrosion increases with the standing time and contamination of the solution. If an agent is not permitted for aluminum parts, it should never be used. Contact with a non-permissible agent may lead to irreparable damage.

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
Manual cleaning

The **HOPKINS** telescopes may be cleaned using an endoscope cleanser (Appendix). Carefully remove all blood, protein material and contaminants from both optical ends of the telescope. Wipe the outside of the telescope with a sponge or cloth soaked in cleaning solution.

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Never remove contaminants on the optical end faces and on the jacket tube with hard objects. Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661) (see Removal of deposits on optical end faces).

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Clean the fiber surfaces and optical end faces with a soft cloth, sponge or cotton swab, soaked with 70% alcohol.

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
The instrument must be rinsed thoroughly after every cleaning or disinfection cycle. Completely demineralized water is recommended for rinsing.

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Then dry using a soft cloth or with sterile compressed air.

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Finally, dry the fiber surfaces again, including the optical end faces, with alcohol.

Note: The fiber surfaces in the light inlet piece especially need to be dried again carefully with alcohol. Residues of disinfectant and cleaning solutions in the light inlet piece may cake and harden when the light guide is switched on and considerably impair light transmission.

Check the telescope for a clear view.

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Once the cleaning work has been completed, free the cleaning instruments, especially the brushes, from residue and disinfect them. Never reuse damaged brushes.

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Instructions for telescopes with a working channel and LUER LOCK stopcocks
Sheath lumens and stopcocks require particular attention while cleaning. To maintain a smooth surface, polishing with "DURAGLIT" polishing wool is recommended if required (see Accessories). For the sterilization procedure, the stopcocks must be taken apart in order to guarantee complete circulation of the sterilization medium.

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Cleaning before disassembly
Remove the rubber cap.

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Clean the outer surface of the sheaths with a sponge, brush or a soft cloth.

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Clean the sheath lumen using a suitable flexible brush. If necessary, a water and compressed air gun may be used (see Accessories).

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Clean the nozzle first with the stopcock opened (cock plug parallel to the nozzle) using a suitable flexible brush.

Then rinse out the stopcock. Make certain rinsing continues until clear water runs out the front end.

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Disassembly
After loosening the thumbscrew remove the instrument channel adapter.

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Both adapters on the light inlet piece must be removed before cleaning.

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Loosen the thumbscrew on the stopcock, remove cock plug .

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... nozzle thoroughly with a suitable brush.

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The nozzle and the sheath lumen must be rinsed out well on the inside at every cleaning. Then dry using sterile compressed air. Again, the cleaning gun with various attachments (Art. no. 27660) may be used for this (see Accessories).

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The rubber cap, adapter and light inlet piece may be cleaned separately.

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Assembly

Thinly coat the contact surfaces on the cock plug with special grease (Art. no 27657). This thin layer of grease not only serves to lubricate the metal-to-metal contact surfaces, but also seals the stopcocks.

Warning: Use only special silicone or oil-based lubricants, since sterilants cannot penetrate through every silicone or oil.

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... fix in place with the thumbscrew . Tighten slightly.

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Test the cock plug for proper operation.

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Note: Take stopcock apart again for sterilization.

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Attach the adapters for the instrument channel and tighten using the thumbscrew .

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Screw both adapters onto the light inlet piece.

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Machine cleaning

When preparing **HOPKINS** telescopes in a cleaning and disinfection machine, particular care is necessary in the choice of chemical agents and machine programs. The risk of damage is much higher than for manual cleaning.

HAMOU telescopes with an eyepiece drive mechanism are not suitable for preparation by machine.

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The machine procedure must be selected in agreement with the machine manufacturer and with the manufacturer of the chemical reagents. Only special procedures may be used which have been tested and approved for this purpose. It is preferable to use thermal disinfection. Completely demineralized water, at least in the final rinse cycle, protects the surfaces.

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
Depending upon the water quality, a final rinse after cleaning, with an added neutralization agent, may be necessary. A neutral pH cleanser should be used (e.g. enzymatic cleanser). If an agent is not expressly approved for aluminum parts, it should not be used since this can lead to changes in the surface.

Jolting and bending the telescopes, as occurs during improper storage, absolutely must be avoided.

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Especially the jacket tube of telescopes must be protected from the strong mechanical forces of the water jet. Use of retaining devices which specially protect the telescope is recommended (see Accessories, Art. nos. 39501A, 39501B, 39501C, 39501U). For telescopes with an irrigation channel, connect the channel to the machine via an irrigation connection. After a program has finished, remove the telescopes as soon as possible from the machine. Corrosion damage could otherwise result.

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Inspection

After cleaning, check the telescope for a clear view.

Image waveguide

Check the optical end faces (distal tips and eyepiece) for scratches and surgical or cleaning residues. The surfaces should appear smooth and shiny. To test the image quality, the telescope should be turned slowly while looking through it.

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If the image is completely or partially defective, either the rod-lens system in the jacket tube has broken or the lens in the eyepiece is defective. In this case, the telescope must be replaced.

A cloudy or spotty image may be caused by moisture or disinfectant residue on the optical end faces. In such cases, carefully clean the optical end faces with a soft cloth or swab soaked in alcohol.

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Light guide

Hold the distal tip in the direction of a light source and check the number of dark points at the light guide connector. These dark points represent broken optical fibers in the light guide bundle. Individual broken optical fibers do not necessarily mean a noticeable reduction in image quality. However, above a certain amount of approx. 20 to 25%, light transmission and therefore image quality is significantly reduced. The telescope should no longer be used then.

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Removing deposits on optical end
faces

If the telescopes are autoclaved, stubborn deposits may, in the course of time, develop on the glass surfaces due to foreign matter in the steam and chemical reactions between the foreign matter and the glass surfaces. These deposits can be removed with the special cleaning paste (Art. no. 27661) included with the autoclavable telescopes.

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
Dab some cleaning paste onto a clean, moist cleaning swab and spread it over the glass surfaces. Then polish to remove even stubborn deposits.

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
Finally rinse the glass surfaces with water and dry them using alcohol, or subject them to a complete cleaning cycle.

! Caution: Cleaning with special cleaning paste should only performed if the image as viewed is cloudy (after approx. 10-20 sterilizations) and not as a matter of routine after every cleaning.

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Removing deposits on metallic
surfaces

Deposits on the metal sheaths can be removed with DURAGLIT polishing wool (Art. no. 27658).

Then the part must again go through a complete cleaning cycle.

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Disinfection and sterilization

Warning: High-level disinfection for use on the patient is recommended only for telescopes which come into contact with mucous membrane or relatively small skin injuries. High-level disinfection for use on the patient is not recommended for telescopes used in laparoscopic, arthroscopic, gynecological, reconstructive, or cosmetic surgery. Sterilization is recommended in these cases. The disinfection mentioned here does not concern disinfection done for the purpose of protecting the personnel.

Warning: Before any disinfection/sterilization the telescopes must be carefully cleaned and any organic material, blood and cleaning solution must be removed. Only clean surfaces can be reached by the sterilization.

Warning: Remember that sterilization is only successful if it proceeds according to recognized and validated methods.

Warning: Any deviation from our recommended disinfection parameters must be validated by the user.

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Caution: Make certain all parts are securely positioned during sterilization. To do this, articles from the chapter "Spare parts, recommended accessories" or the catalog chapters "Cleaning and care" may be useful. Packaging should conform to EN 868.

Caution: Please consult the Water and steam quality section in order to prevent damage to the instrument set.

Caution: During sterilization the telescopes must not have any direct contact with metal.

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Quick overview of the sterilization methods

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* STERRAD sterilization is generally not recommended due to limitations regarding the size of the lumina (see the STERRAD instructions in this instruction manual)

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Chemical Disinfection

With the exception of the HAMOU telescopes, all telescopes can be placed in disinfectant solution. Suitable as a cleaning/disinfectant solution are those agents recommended especially for endoscopes, e.g. Cidex, Gigasept or Kohrsolin i.D. (see Appendix).

For sterilization of the HAMOU telescopes see special instructions (see next page).

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
A plastic basin (Art. no. 27645, 27646) should be used to avoid scratching the telescopes and to eliminate electrolytic corrosion which can occur when dissimilar metals are soaked in the same solution.

To rinse off, remove the instruments from the solution using the sieve and rinse in a second basin with sterile water. Air bubbles in the instruments are removed by tipping the sieve slightly.

After removal from the disinfection solution, rinse thoroughly with sterile water. Dry the telescope with a sterile cloth. Apply alcohol to the glass surfaces to complete the drying process.

Alcohol must be applied carefully, paying special attention to the fiber surfaces in the light inlet piece to complete the drying process. Any residues of disinfectant or cleaning agents in the light inlet piece could burn onto it when the fiber optic light cable is connected and this would considerably impair light transmission.

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Disinfection of the
HAMOU I /HAMOU II telescopes

26156 B, 26156 BU, 27156 BU, 28720 BH,
7200 BH / 26001 A, 26157 B, 26157 BT,
27018 BH, 27020 BH, 27157 B, 28157 B,
28301 BF, 28720 BF, 63157 B, 64200 BH,
7200 BF

The HAMOU telescopes are endoscopic telescope and microscope combined. They must be cared for and handled carefully and in a different manner than other KARL STORZ **HOPKINS** telescopes:

HAMOU telescopes must not be steam sterilized. Only HAMOU II telescopes can be gas or plasma sterilized.

The insertable portion of the telescope may be soaked in a disinfectant solution.

The metal examination sheath, operating sheath and accessory operating instruments may be soaked in a disinfectant solution or sterilized in any standard gas or steam autoclave.

Stopcocks should be disassembled for thorough cleaning and lubrication.

The head of the HAMOU I telescopes may be wiped with a soft cloth slightly dampened with a mild soap solution or 70% isopropyl alcohol.

Jacket tube of telescope should be thoroughly cleaned with a mild soap solution, rinsed and placed in a special holder (Cat. No. 26156 S, see accessories)

Do not allow solution above the level of the locking device. Follow disinfectant solution manufacturer's instructions regarding proper dilution and time period for soaking in a disinfectant solution, but do not exceed 60 minutes.

After removing from the disinfectant solution, rinse jacket tube thoroughly with sterile water.

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Steam sterilization for autoclavable
HOPKINS® telescopes
KARL STORZ **HOPKINS**® telescopes that are
marked „AUTOCLAVE“ can be steam sterilized up
to 134°C without reservation.
For steam sterilization (with saturated steam at
134°C) follow the method in DIN 58946 Part 1,
Clause 3.25.2 and EN 285. Fractionated pre-vacu-
um methods are preferred.

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The fractionated fore-vacuum procedure at 134°C
is a component of the test series that KARL
STORZ telescopes are subjected to. Before a
telescope is classified as steam-sterilizable, it must
withstand several hundred sterilizations undam-
aged in a specially designed and freely program-
mable sterilizer.
With endoscopes, steam sterilization always
means greater stress than immersion.

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During steam sterilization at 134°C, a pressure of
3.05 bar is applied to the telescope and all sealing
surfaces. KARL STORZ endoscopes are designed
and tested accordingly. The physical influence is
considerably less for immersion in a sterilization
solution, however the chemical stress can be much
greater and, if unsuitable chemicals are used, it can
lead to early failure.

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It is impossible to say whether an autoclavable
telescope lasts longer than an immersible tele-
scope. The two methods have different effects
on the endoscope and provide different levels of
hygienic safety. With chemicals it is only possible to
disinfect and not sterilize, as with steam.

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Note: If deactivation of Creutzfeldt-Jakob-disease (CJD) is necessary, the sterilization method recommended by WHO may be used (18 minutes at 134°C). Extension of the sterilization time up to one hour at 134°C is possible in special cases.

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Gravitation methods are mild, though they may only be used for telescopes without channels and without stopcocks.

The steam must follow requirements in EN 285, Appendix B, to prevent discolorations and spotting (see Water and steam quality).

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Clean and dry the telescopes. After caring for the instruments, thinly coat the sealing and contact surfaces, such as on the stopcocks or plungers, with special grease (Art. no. 27657). These parts must be disassembled for sterilization and only replaced in the housing after sterilization. Place the telescopes within a suitable sterilization container (see accessories). Cloths, used for packing the telescopes, must not contain any detergent residue, since this also causes discolorations.

Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage. Load the sterilizer according to its instruction manual.

When the sterilization cycle is complete, remove container from autoclave and allow telescopes to cool to room temperature before removing top container.

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-  Caution: Nonautoclavable telescopes must not be steam sterilized.
-  Caution: Steam sterilization at 134°C is gentler than steam sterilization at 121°C since the prolonged exposure time has a negative effect.
-  Caution: During steam sterilization, the telescopes should not come into direct contact with metal instruments, trays etc.
-  Caution: Do not open the sterilization container immediately after the sterilization process has finished. Never attempt to cool telescopes by pouring cool, sterile liquid over them. Forced cooling will cause severe damage to the telescope.
-  Caution: Hot-air sterilization is not permissible.

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Steam sterilization is not permitted for the following instrument categories. Only methods with a maximum temperature of 65°C, such as gas and plasma sterilization, may be used.

- Nonautoclavable telescopes (Cleaning up to a Temperature of 93°C is allowed)
- Miniature Endoscopes (Cleaning up to a temperature of 65 °C is allowed)
- Telescopes with HAMOU eyepiece drive mechanism
- Magnifying lenses with drive mechanism
- Fiberscopes
- Fluid light cable
- Video cameras
- Prismatic light deflectors
- Headbands

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**Validated
steam sterilization methods**

Fractionated vacuum methods in accordance with EN 285, as well as methods in accordance with DIN 58946 Part 1, Clauses 3.25.1.1 and 3.25.2.1 have been validated.

To achieve the desired sterility assurance level (SAL) of 10⁻⁶, KARL STORZ recommends the following EtO and steam sterilization methods.

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The Vacuum Methods

The vacuum methods consists of four basic phases:

conditioning phase, exposure phase, exhaust phase and drying phase.

The conditioning phase removes air once or several times from the chamber by pulling a vacuum and then warms the instruments by injecting steam.

Sterilization occurs during the exposure phase when the chamber reaches a temperature of 121°C or 134°C and pressure of 1 or 2 bar. The exposure phase in a pre-vacuum type of sterilizer is 4.0 minutes.

The sterilization time is between 5 and 20 minutes, depending on the method, temperature and pressure.

The drying phase is accomplished by pulling an additional vacuum and allowing the instruments to dry under vacuum for approximately 20 minutes.

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Pre-vacuum

The pre-vacuum method consists of four phases.

- Conditioning phase
The conditioning phase removes air from the chamber by pulling a vacuum and then warms the instruments by injecting steam.
- Sterilization phase
Sterilization takes place at a temperature of 134°C and a Operating pressure (p_b) of 2 bar over a period of 5-8 minutes.
- Evacuation phase
The steam is released.
- Drying phase
Drying takes place under renewed application of a vacuum over a period of approx. 5-20 minutes.

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Fractionated pre-vacuum procedure

The fractionated pre-vacuum procedure also consists of the four previously mentioned phases. In the conditioning phase, a vacuum is generated up to four times. Steam is then blown in and sterilization takes place at 121°C (250 to 252°F) or 134°C (270 to 272°F) in 20 or 5 minutes respectively (minimum times).

Before the final drying phase, the steam is removed from the chamber during the evacuation phase. Drying takes place under vacuum.

Validated parameters for the fractionated pre-vacuum method:

Temperature:	270 to 272°F; 134°C
Pressure (p_b):	27 psi; 2 bar
Exposure time:	5 minutes for all instruments.

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Gravitation method

Warning: This method is only suitable for smooth and accessible surfaces, not for instruments with channels, gaps or stop cocks.

The gravitation method is also made up of 4 phases:

- Conditioning phase
Steam is pressed into the chamber and heats the instruments, whereby the air in the chamber may be released through an outlet.
- Sterilization phase
Sterilization takes place at a temperature of 121°C and a Operating pressure (p_b) of 1 bar over a period of 30 - 45 minutes. Sterilization may also take place at 134°C and 2 bar, however the time required is 5 - 10 minutes.
- Evacuation phase
Steam is released through the outlet.
- Drying phase
Drying takes place under ambient air pressure by heating through the heating jacket of the autoclave.

Validated parameters for the gravitation method:

Temperature:	270° to 272°F, 121°C
Pressure (p_b):	15 psi/ 1 bar
Required time:	30 - 45 minutes

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Flash sterilization

Warning: This method is only suitable for smooth and accessible surfaces, not for instruments with channels, gaps or stopcocks.

Warning: Flash sterilization is a method which, from the perspective of hygiene, must only be used in exceptional cases.

Warning: This procedure is permitted if it is certain that the telescope will not be subjected to any thermal shock.

Note: If the telescope is cooled down too quickly due to the room temperature being lower when the telescope is removed from the unit directly after autoclaving, the ambient air enclosed in the telescope will cause condensation. This misting will mean the telescope cannot be used until the condensation has been absorbed by the air.

Flash sterilization may take place either in a fore-vacuum or gravitation sterilization unit. flash sterilization in a fore-vacuum unit takes place without a conditioning or drying phase. Flash sterilization in a gravitation sterilization unit takes place without a drying phase.

We have validated flash sterilization with the use of a gravitation sterilization unit under the following conditions:

Temperature:	134 °C
Pressure (p _v):	2 bar
Exposure time:	10 minutes

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Gas sterilization/plasma sterilization

Sterilization using formaldehyde or ethylene oxide and plasma sterilization are gentle and therefore recommended methods of sterilization for the **HOPKINS** telescopes, with the exception of HAMOU I telescopes. The latter telescopes are subject to special sterilization instructions.

Immerse the telescopes in suitable vessels and vacuum seal them in suitable plastic film.

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Note: Observe national laws and regulations for the gas sterilization method.

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For gas sterilization with ethylene oxide, the following airing times must be observed for instruments due to gas absorption of the materials:
When airing in the airing cabinet, the temperature should be maintained between 45°C and 55°C over a period of twelve hours.

	at 20°C	at 42°C	at 55°C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Plastics	5 days	2 days	30 hrs

(minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes).

Maximum acceptable levels (per 21 CFR 812.100; Federal Register, Vol. 43, No. 122, June 23, 1978) of residues following EtO sterilization are as follows:

Ethylene oxide:	25 ppm* / 250 ppm**
Ethylene chlorohydrine:	25 ppm* / 250 ppm**
Ethylene glycol:	250 ppm* / 5000 ppm**

** for instruments in contact with skin and mucous membrane,
* for instruments in contact with blood.

The instruments can be used immediately after sterilization if sterilized using formaldehyde or plasma sterilization. Extra airing time is not necessary.

Plasma sterilization

Plasma sterilization (STERRAD®) is a special method which is particularly suited for thermolabile items. In plasma sterilization, changes to the coloration of the surfaces of aluminum alloys could occur, though this does not influence the function of the instruments. It must be pointed out in this context that anodized aluminum parts are not always recognizable as such.

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Warning: Sterility is only assured if the manufacturer's instructions for the STERRAD® procedure are followed.

Warning: All deviations from the sterilization parameters recommended for the STERRAD® system must be validated by the user.

Caution: HAMOU I telescopes with eyepiece drive mechanism (26156B/BU, 26157BU, 28720BH, 7200BH) are not suitable for the STERRAD® method.

Note: Instruments which KARL STORZ has classified as compatible with the STERRAD® sterilization procedure, have been validated with at least a hundred STERRAD® cycles.

Note: Observe the requirements of the sterilizer manufacturer when packaging. The use of telescope holders 39501A to V is possible without impairing the results of sterilization.

No outgassing time is prescribed for plasma sterilization.

Following Product categories have been tested using the STERRAD® method. They are suitable for sterilization without significant influence on the material or function.

- Rigid telescopes (except telescopes with eyepiece drive mechanism)
- Flexible, semiflexible fiberscopes and miniature endoscopes
- Video cameras¹ (except Endovision XL)
- Ureteral light probe
- Fiber and fluid light cable
- HF cable
- Insulated and noninsulated Take-apart and c instruments
- Insulated surgical instruments
- Non-insulated surgical instruments
- Trocars / sheaths / cannulas / obturators / Working elements¹

¹ KARL STORZ video cameras manufactured after September 5, 1996, are compatible to STERRAD®. Cameras manufactured before this date may be refitted. Please contact the KARL STORZ service department for details regarding refitting.

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Chemical sterilization

Chemical sterilization (STERIS[™]) is a special method which is particularly suited for thermolabile goods such as fiberscopes.

Except for telescopes with HAMOU eyepiece drive mechanism, all KARL STORZ telescopes are suitable for STERIS[™] sterilization. Color changes to the surface of aluminum alloys could occur during sterilization, which do not however affect the function of the instrument. It must be pointed out in this context that anodized aluminum parts are not always recognizable as such.

Warning: Sterility is only guaranteed if the manufacturer's instructions for the STERIS System have been followed. This also applies to use of the sterilization containers.

Warning: Parts with a lumen must be prepared together with the STERIS Quick Connect kit. Parts with a lumen which are open only on one side should never be sterilized with the STERIS System. All parts must be thoroughly cleaned before sterilization with the STERIS PROCESS[™], which also applies to the other sterilization methods.

The following product categories have been validated with the STERIS PROCESS[™]:

- Flexible ureterosopes, pediatric cystoscopes and choledochoscopes manufactured after July 1, 1997
- Flexible cystoscopes for adults, and hysteroscopes
- Rigid and semiflexible telescopes (without working channel) (except HAMOU with eyepiece drive)
- Video cameras (except nonimmersible video cameras)
- Illuminated ureteral probes
- Light cable
- Calcutript connecting cable

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KARL STORZ product categories which have material compatibility, but which have not been validated in reference to sterility:

- Bronchoscopes
- Gastrosopes and veterinary gastrosopes
- Working elements, proctoscopes and resectoscopes
- Trocars / sheaths
- Reusable tube sets
- Insulated and noninsulated Take-apart and Click Line™ instruments
- Insulated and noninsulated forceps, scissors, punches, and suture instruments
- Instruments for coagulation, irrigation and suction

KARL STORZ product groups that are not sterilizable with the STERIS PROCESS™:

- Telescopes with HAMOU eyepiece drive mechanism
- Magnifiers with drive mechanism
- Motors and their handles (e.g. for shavers)
- IMPERATOR handpieces
- Prismatic light deflectors

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Hot-air sterilization


Caution: Hot-air sterilization is not permissible. Since different metals and plastics are used, damage could result e.g. due to differing expansion of the materials used or high temperature.

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Storage

Instruments should always be stored in suitable containers (see Accessories) and must be absolutely dry.

The place of storage must have normal temperature and air humidity values (30 - 50 % r.h.; 20°C).

Note: To dispose, no special measures are necessary. National laws and regulations must be observed.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts only.

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Limitation of liability

KARL STORZ GmbH & Co. shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly operations, system expansions, readjustments, modifications, or repairs to same have been performed by a person or persons duly authorized by KARL STORZ; and that the instrument has been used in accordance with its operating instructions at all times.

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Warranty

For two years as of delivery to the end customer, we agree to replace the goods free of charge if proof can be provided of faulty materials or faulty workmanship. In doing so we cannot accept the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co.
Postfach 230
D-78503 Tuttlingen/Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall release us from any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Directive compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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Repair and replacement program

KARL STORZ maintains a repair and replacement depot which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective instrument. For telescopes, you receive a guarantee of 1 year, and for instruments 6 months. For fiberscopes and equipment, individual repair is necessary. To bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co
Repair Service Dept.

Dr.-Karl-Storz-Strasse 34, 78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfectated before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

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The KARL STORZ Repair, Refurbish, and Exchange Plan for Telescopes

Telescopes rely on a complex, finely balanced system of light transmitting optics for top performance. The KARL STORZ Repair, Refurbish, and Exchange Plan provides unparalleled services to extend the life of your endoscopic surgical system.

Our Repair, Refurbish, and Exchange Plan offers competitive tiered pricing based on the complexity of the repair.

When KARL STORZ refurbishes a telescope, we provide the following services:

- We completely disassemble the telescope, inspect, clean, and polish all parts.
- We replace any damaged parts with new, genuine Karl Storz parts.
- We provide quality control checks at every stage of refurbishment.

A KARL STORZ repaired or refurbished telescope meets the same rigorous standards as our new scopes:

- We utilize exclusive **HOPKINS** rod lens technology.
- We seal our telescopes with proprietary adhesives which are proven to resist moisture and withstand high temperatures associated with frequent sterilization cycles.
- We manufacture our own light transmission components.
- Factory trained technicians complete work in compliance with original specifications.

When KARL STORZ repairs or refurbishes a telescope you can rely that work will be conducted in the same controlled environment in which the **HOPKINS** telescopes were built.

The KARL STORZ Repair, Refurbish, and Exchange Plan covers our entire rigid telescope inventory.

Note: We do not supply our rod lens systems to any outside companies so that no other repair company besides our own is capable of repairing our telescopes with our original parts.

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
Compromised Optics

Telescope lenses require precise seating and alignment. When repairs fail to address this issue the results carry far reaching clinical implications. The following pictures show examples of inadequately repaired optics compared to examples of KARL STORZ factory repaired optics.

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The upper picture (Fig. 1a) clearly shows visible distortions caused by improper repair. In addition, there is a remarkable loss of brightness and sharpness.

The lower picture (Fig. 1b) shows the same view through an identical telescope, repaired by KARL STORZ however.

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A HAMOU telescope in the contact endoscopy setting. The partial loss of focus caused by the improper repair (Fig. 2a) can clearly be seen ...

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...in comparison with a telescope repaired by KARL STORZ (Fig. 2b).

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



Sterilizing and storage
Note: The containers may be used for sterilization with
vapor or gas.

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

Excerpt from the list of chemicals for the preparation of telescopes and endoscopic instrument sets

Manual cleaning and disinfection/ultrasound
The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU* telescopes with ocular focusing, magnifiers with focusing, motors and their handles, IMPERATOR drill handles, prismatic light deflectors. The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to proper dilution and recommended immersion times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Hersteller/ Manufacturer/ Fabricante	Handelsname/ Commercial name/ Nombre comercial	Typ/ Type/ Tipo
Erolab GmbH & Co. OHG	Bakasept Plus Bakasept Pulver Classic Ultrasepht AF	1, 2, 3 1, 2, 3 1, 2, 3
Estlar Pharma GmbH	Ultrasepht Aktiv Ultrasepht Classic Ald	1, 2, 3 1, 2, 3 1, 2
Fassentus AG	Ald plus Sporcid FF	2 1, 2
Hollia Polska Sp. z o.o	Polisept Hollia	1, 2
Indeba Industria E Commerciao LTDA	Letahidol	2
Jose Collado S.A.	Darodor 4000 Liquido Darodor 9000 Darodor Sinialdehyd 2000 Midasan 2000	1 2 1, 2 1, 2
Lysolform Dr. Rossmann GmbH	Almyral Desofom Lysolformin 2000	1, 2 1, 2, 3 1, 2
Mediochem International	MediClo MediZym	2 1, 3
Merz Hygiene GmbH	Edisonite Super Mucadant-IS Mucadant-Zymaktiv	1, 3 2, 3 1, 3
METASYS Medizinische Technik GmbH (formerly LOWZA)	Mucocort-N neu ID 50	1, 2, 3 2
Orochemia Promagant AB	A 20 Wavacide Desoton Extra Desoton forte Parifektan Endo Parifektan Neu	1, 2, 3, 4 2 1, 2 2 2, 3 1, 2
Schumacher, Dr. GmbH	Parifektan TB Gigasapt Gigasapt AF forte Gigasapt FF Gigasapt Instruk AF Gigasapt Med forte Gigazyme	1, 2, 3, 4 2, 3 1, 2 1, 2 1, 2, 3, 4 1, 2 1, 3
Staris	Lysolol V EnyCase 2 neodisher LM 2	2 1, 3 1, 3
Walgart, Dr. GmbH & Co.	neodisher medician neodisher medician forte neodisher medicym neodisher Sapto Clean neodisher Sapto ProClean	1 1, 3 1, 3 1, 4 1, 2
Whitlsey Industries PTY. Ltd	Aldal Plus Maltix	2 1

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* Caution: The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set.

For flexible endoscopes, only immerse the examination sheath, not the housing/handle part.

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Typ

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Type

- 1= Cleaner
- 2= Disinfectant
- 3= Suitable for ultrasound baths
- 4= Neutralizer
- 5= Rinse aid

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine cleaning and disinfection

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU telescopes with ocular focusing, magnifiers with focusing, motors and their handles, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

	Thermosapt Alka Clean	1
	Thermosapt DK	2
Schülke & Mayr GmbH	Thermosapt ED	2
	Thermosapt ER	1
	Thermosapt RRM-zym	1
Staris	Instro-Klean	1, 3
	needsher DacClean	1, A
	needsher FA	1
	needsher FA forte	1, 3
	needsher Odvario	1
	needsher mediklean	1
Walger, Dr. GmbH & Co.	needsher mediklean forte	1, 3
	needsher mediklar	5
	needsher medizym	1
	needsher SeptoClean	1, 2, A
	needsher Septo DA	2
	needsher Septo DW	2

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To be filled out by instrument owner:
Return address/company stamp:

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Firmenstempel/Unterschrift:

WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

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To be completed by the supplier/importer:

Company stamp/signature:

KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postfach 230, 78503 Tuttlingen, Germany
Phone: +49 (0)7461 708-0,
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com

KARL STORZ Endoskope Berlin GmbH
Ohlauer Straße 43
10999 Berlin, Germany
Phone: +49 (0)30 30 69090
Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.
2345 Argentia Road, Suite 100
Mississauga, ON, L5N 8K4, Canada
Phone: +1 905 816-8100, Fax: +1 905 858-0933
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
600 Corporate Pointe
Culver City, CA 90230-7600, USA
Phone: +1 310 338-8100, +1 800 421-0837
Fax: +1 310 410-5527
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 308 A,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, Mexico
Phone: +52 5552 5056 07, Fax: +52 55 5545 0174

KARL STORZ Endoscopia Argentina S.A.
Cerviño 4449 Piso 10°
1425 Buenos Aires C. F., Argentina
Phone: +54 11 4772 4545, Fax: +54 11 4772 4433
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14, 12739 Skärholmen, Sweden
Postaladdress: PoBox8013, 14108KungensKurva,
Sweden
Phone: +46 8 50 56 480, Fax: +46 8 50 56 4848
E-Mail: kundservice@karl-storz.se

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopie Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopy France S. A.
12, rue Georges Guynemer
Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 146/11/18
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Francia – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

KARL STORZ Endoscopia ROMANIA srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

TOV KARL STORZ Ukraine
18b Geroev Stallingrada avenu
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19, -20
Fax: +380 44 42668-14

OOO KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Foenicia Street, Mina Elhosn
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax +961 1 365151

KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kasachstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@itte.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 – 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty. Ltd.
174 Parramatta Rd.
Camperdown, N.S.W. 2050, Australia
Phone: +61 02 8594 9150, Fax: +61 02 8594 9199
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
10-18 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosgpp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germania

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com