Medtronic

Protégé™ RX

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1 Device description


The Protégé RX self-expanding carotid stent system is a self-expanding nitinol stent system intended for permanent implantation. The stent is made of a nickel titanium alloy (nitinol) and comes premounted on a 6 Fr rapid exchange (Rx) delivery system. The stent system is compatible with a 0.36 mm (0.014 in) guidewire. The stent is cut from a nitinol tube in an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent.

The delivery system, as shown in Figure 1, is comprised of an inner shaft and outer sheath, which are locked together with a tuohy-borst valve (1). The inner shaft terminates distally in a flexible catheter tip (2) and originates proximally at the end of the catheter. Two radiopaque markers, one marker distal (3) and one marker proximal (4) to the constrained stent (5), are mounted on the inner shaft. Delivery systems with tapered stents have an additional radiopaque marker that identifies where the stent diameter transition begins.

Tapered stents are mounted onto the delivery catheter with the smallest diameter toward the distal end of the catheter. The outer sheath connects proximally to the Y-connector (6). The space between the inner shaft and outer sheath constrains the self-expanding stent. This space is flushed through the side port (7) of the Y-connector before the procedure.

The catheter is back-loaded through the distal tip onto a 0.36 mm (0.014 in) guidewire, which exits through the Rx port (8). Radiopaque markers (9 and 10), which mark the location of the constrained stent, provide guidance for positioning the stent before deployment. Before deploying the stent, turn the tuohy-borst valve counterclockwise to unlock the outer sheath. The outer sheath is retracted when the Y-connector (6) is pulled toward the proximal grip (11).

Figure 1. Stent delivery system

2 Intended purpose

The intended purpose of the Protégé RX self-expanding stent system is to restore lumen patency and blood flow in the common carotid artery (CCA), internal carotid artery (ICA), or carotid bifurcation. Implanting the Protégé RX stent in the target vessel restores lumen patency and re-establishes blood flow in the artery through the gentle, outward force exerted by the stent.

2.1 Intended patient population

The Protégé RX self-expanding stent system is intended for adult patients in need of permanent implantation to improve and maintain luminal diameter in stenoses in the common or internal carotid artery in patients with high surgical risk.

This device is not intended for pediatric patients.

2.2 Indications for use

The Protégé RX self-expanding carotid stent system is indicated for treatment of stenoses of the common carotid artery (CCA), internal carotid artery (ICA), and carotid bifurcation.

2.3 Contraindications

The Protégé RX self-expanding stent system is contraindicated under the following conditions:

- Patients for whom antiplatelet therapy, anticoagulants, or thrombolytic drugs are contraindicated
- Peripheral vascular disease or anatomy, which preclude safe sheath, guide catheter or stent system insertion
- Patients with known hypersensitivity to nickel titanium
- Total occlusion of the carotid artery
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Presence of intracranial tumor(s), arteriovenous malformations (AVMs), ipsilateral intracranial stenosis that requires treatment, or aneurysm requiring treatment
- Patients with lesions in the ostium of the common carotid artery
- All the customary contraindications for PTA

2.4 Clinical benefits

The clinical benefits of the Protégé RX self-expanding stent are restoration of blood flow and reduction of the risk of stroke.

2.5 Intended users

Protégé RX stent implant procedures should only be performed by physicians who have experience with interventional techniques in the vascular system.

The device should only be used within a clinical setting in a sterile environment. Use of this device is limited to facilities where surgical support is readily available in the event of a serious complication.

2.6 Performance characteristics

The Protégé RX stent is a nitinol self-expanding carotid stent with tantalum radiopaque markers at the proximal and distal ends of the stent. It conforms to tortuous vessel anatomy because of the flexibility and kink resistance enabled by its open-cell design. The design of the stent is tailored to provide sufficient radial strength and to maintain the vessel lumen diameter.

2.7 Expected device lifetime

The Protégé RX stent is designed to be a permanent implant with a minimum expected lifetime of 10 years following implantation.

3 Warnings

- Through contact with the delivery system, patients may have limited exposure (< 24 hours) to cobalt, which is considered carcinogenic, mutagenic, or toxic to reproduction (CMR).
- This device was designed for single use only. Do not reuse, reprocess, or sterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Carotid procedures require the use of an embolic protection device. Refer to Medtronic product labeling for placement instructions of the embolic protection device. The safety and efficacy have only been demonstrated with the combination of the Protégé RX self-expanding stent system and the Medtronic embolic protection devices.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen, or may cause premature deployment or deployment complications. If resistance is felt, carefully withdraw the stent system without deploying the stent.
- If resistance is felt during movement through the sheath, carefully withdraw the stent system without deploying the stent.
- If resistance is felt when pulling back the Y-connector, do not force deployment. Carefully withdraw the stent system without deploying the stent.
- If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath contacts the catheter tip and withdraw the system as one unit.

4 Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the sterile packaging has been unintentionally opened or the device is damaged.
- Do not use the stent if it is partially deployed.
- Always use an introducer sheath during the implant procedure to protect the vessel and the puncture site. Support from an introducer sheath is also necessary to minimize stent lengthening or shortening during stent deployment.
- To manage bradycardia and hypotension, maintain venous access for pharmaceutical intervention or the placement of a temporary pacemaker.
- During dilation, do not expand the balloon to the point that bleeding or dissection could occur.
- To avoid potential entanglement, allow for and maintain adequate distance between the embolic protection device and the stent delivery system or deployed stent.
- The stent is not designed for recapturing at any time and is not designed for repositioning after establishing vessel apposition.
- The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- Use caution when crossing a deployed stent with any adjunct device.
- Do not expand the stent beyond its nominal diameter.
- Dispose of the delivery system and accessories in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.
- Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.

5 Potential adverse events

The potential adverse events (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

Table 1. Potential adverse events

Risks associated with carotid percutaneous procedures	
• Abrupt closure	• Fever
• Allergic reactions to procedural medications, contrast dye, or device materials	• Hematoma
• Amaurosis fugax	• Hyperperfusion syndrome
• Aneurysm	• Hypotension or hypertension
• Angina or coronary ischemia	• Infection or pain at the puncture site
• Arrhythmia	• Intraluminal thrombus
• Arterial occlusion or thrombosis (at puncture site or remote site)	• Ischemia or infarction of tissue or an organ
• Arteriovenous fistula	• Myocardial infarction (MI)
• Bacterial infection or septicemia	• Nerve damage
• Bleeding from anticoagulant or antiplatelet medications	• Pain (in head and neck)
• Bleeding, with or without transfusion	• Renal failure
• Cerebral edema	• Renal insufficiency (new or worsening)
• Cerebral hemorrhage	• Seizure
• Cerebral ischemia or transient ischemic attack (TIA)	• Severe unilateral headache
• Congestive heart failure (CHF)	• Slow or no-flow during procedure
• Death	• Stent or embolic protection device entanglement or damage
• Distal embolism	• Stroke or cerebrovascular accident (CVA)
• Emergent or urgent endarterectomy surgery (CEA)	• Total occlusion of carotid artery
• Femoral pseudoaneurysm	• Vessel dissection, flap, perforation, or rupture

Risks associated with stent placement	
• Detachment of a component of the delivery system	• Stent damage or entanglement
• Embolization (from air, tissue, or thrombus)	• Stent migration
• Failure to deploy stent	• Stent misplacement
• MRI-related stent heating or vessel damage	• Thrombosis or occlusion of the stent
• Restenosis of stented segment	• Vessel spasm or recoil
• Stent collapse or fracture	

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

Note: After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00001066X.

6 Patient counseling information

Physicians should review the instructions for use for applicable information to be shared with the patient. A patient implant card, which contains identifying information about the implanted device, is included in the device package. After device implant, complete the patient implant card and provide it to the patient before they are discharged.

Physicians should communicate the following instructions to their patients:

- Always carry their implant card with them.
- Access additional information about their device on the website that is listed on their patient implant card.
Note: If the patient is unable to access the website, the physician must provide the information from the website to the patient.
- Always inform any healthcare personnel that they have an implanted device before any procedure has begun.
- Contact their physician if they notice any new or changing symptoms.

7 Storage

Store the device in a dry area away from sunlight.

8 How supplied

The device is supplied sterile and is sterilized using ethylene oxide.

9 Device materials

The stent is comprised of nitinol with tantalum markers. Some patients are allergic to, or can become sensitive to, nickel.

Table 2. Materials in surface area of implantable device

Material description	Surface area of implant (%)
Nitinol	≥ 99.8%
Ta	≤ 0.2%

The following substances are known to be carcinogenic, mutagenic, or toxic to reproduction (CMR). The substances to which patients can be exposed and have been identified to be present above the > 0.1% weight by weight threshold are listed in Table 3.

Table 3. CMR substances

Substance	CAS number	Concentration (% weight by weight)
Cobalt (Co)	7440-48-4	> 0.1% w/w

Warning: Through contact with the delivery system, patients may have limited exposure (< 24 hours) to cobalt, which is considered carcinogenic, mutagenic, or toxic to reproduction (CMR).

10 Directions for use

Warning: This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

10.1 Selection and preparation

10.1.1 Recommended items

- 5 mL syringe filled with heparinized saline
- 0.36 mm (0.014 in) exchange guidewire
- Introducer sheath with hemostatic valve
- Embolic protection device

10.1.2 Select stent size

1. Measure the diameter of the reference vessel (proximal and distal to lesion).
2. Measure the length of the target lesion.
3. Select a stent size according to the vessel diameter and a stent length sufficient to extend proximal and distal to the target lesion.

Table 4 provides the stent sizing specifications.

Table 4. Stent sizing specifications

Stent diameter (mm)	Reference vessel diameter (mm)	Stent length (mm)	Catheter working length (cm)
6	4.5–5.5	20, 30, 40, 60	135
7	5.5–6.5	20, 30, 40, 60	135
8	6.5–7.5	20, 30, 40, 60	135
9	7.5–8.5	20, 30, 40, 60	135
10	8.5–9.5	20, 30, 40, 60	135
8–6 (tapered)	(6.5–7.5) to (4.5–5.5)	30, 40	135
10–7 (tapered)	(8.5–9.5) to (5.5–6.5)	30, 40	135

Note: Properly size the tapered stent diameters to the lumen diameters.

10.1.3 Prepare the stent delivery system

1. Carefully inspect the pouch to ensure that there is no damage to the sterile barrier. Open the pouch and remove the tray and contents. Place the tray on a flat surface and carefully pull off the tray lid and remove the stent delivery system.
Caution: Carefully inspect the sterile package and the device before use. Do not use the device if the sterile packaging has been unintentionally opened or the device is damaged.
2. Examine the distal end of the catheter to ensure the stent is contained within the outer sheath.
Caution: Do not use the stent if it is partially deployed.
Note: If there is a gap between the catheter tip and outer sheath, open the tuohy-borst valve and gently pull the inner shaft in a proximal direction until the gap is closed. Secure the tuohy-borst valve after making the adjustment by turning it clockwise.
3. Flush the delivery system with heparinized saline immediately before backloading the guidewire. Do not power inject the saline. Apply constant, steady pressure to the syringe to ensure the delivery system is fully flushed.
 - a. Attach a 5 mL syringe filled with heparinized saline to the Y–connector side port.
 - b. Loosen the tuohy-borst valve (counterclockwise) and inject saline until it seeps from the proximal end of the tuohy-borst valve.
 - c. Close the tuohy-borst valve by turning it clockwise.
 - d. Continue to flush the system by injecting saline until it comes out between the distal end of the outer sheath and the tip.

- e. Cover the distal end of the outer sheath using fingers and continue to inject saline until it exits the Rx port.
- f. Before loading the catheter onto the guidewire, submerge the distal catheter section into saline.

Note: The system is fully flushed when saline exits the distal end of the outer sheath and the Rx port while applying constant steady pressure.

10.2 Deployment procedure

Warning: Carotid procedures require the use of an embolic protection device. Refer to Medtronic product labeling for placement instructions of the embolic protection device. The safety and efficacy have only been demonstrated with the combination of the Protégé RX self-expanding stent system and the Medtronic embolic protection devices.

10.2.1 Insert the introducer sheath and the guidewire

1. Gain access using an introducer sheath with a hemostatic valve. Refer to the product labeling for the required sheath size. The common femoral artery is the access site for carotid treatments.
Caution: Always use an introducer sheath during the implant procedure to protect the vessel and the puncture site. Support from an introducer sheath is also necessary to minimize stent lengthening or shortening during stent deployment.
2. Insert a 0.36 mm (0.014 in) guidewire of appropriate length across the target lesion.
Caution: To manage bradycardia and hypotension, maintain venous access for pharmaceutical intervention or the placement of a temporary pacemaker.

10.2.2 Dilate the lesion

Per physician discretion, predilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire and the embolic protection device.

Caution: During dilation, do not expand the balloon to the point that bleeding or dissection could occur.

10.2.3 Introduce the stent delivery system

Advance the stent delivery system over the guidewire or the embolic protection device until the guidewire exits through the Rx port. Use one hand to hold the proximal end of the guidewire. Use the other hand to advance the delivery system through the hemostatic valve and guide catheter, or through the introducer sheath, to the site of the lesion.

Caution: To avoid potential entanglement, allow for and maintain adequate distance between the embolic protection device and the stent delivery system or deployed stent.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen, or may cause premature deployment or deployment complications. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Warning: If resistance is felt during movement through the sheath, carefully withdraw the stent system without deploying the stent.

10.2.4 Deploy the stent

1. Advance the delivery system until the distal (leading) radiopaque marker on the inner shaft is distal to the target lesion.
Note: For tapered stents, position the middle radiopaque marker on the inner shaft at the location in which the artery begins to widen.
2. Loosen the tuohy-borst valve to unlock the system.
Note: Initiate stent deployment by retracting the outer sheath (pulling the Y–connector) while holding the inner shaft (proximal grip) in a fixed position. During release of the stent, keep the whole length of the flexible deployment system as straight as possible. To ensure the deployment system is stationary and straight, maintain a slight backward tension on the delivery catheter, using the proximal grip. Deployment is completed when the outer sheath passes the proximal stent marker on the inner shaft and the stent is released from the retainer.
Warning: If resistance is felt when pulling back the Y–connector, do not force deployment. Carefully withdraw the stent system without deploying the stent.
Caution: Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.
Caution: The stent is not designed for recapturing at any time and is not designed for repositioning after establishing vessel apposition.
Caution: The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
Note: If more than one stent is needed, place the most distal stent first. If overlap of sequential stents is necessary, keep the amount of overlap to a minimum.

10.2.5 After stent deployment

1. Use fluoroscopic guidance while withdrawing the entire delivery system as one unit, over the guidewire or the embolic protection device, into the catheter access sheath, and out of the body. Remove the delivery system from the guidewire or the embolic protection device.
Warning: If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath contacts the catheter tip and withdraw the system as one unit.
2. Using fluoroscopy, verify that the stent is fully deployed.
 - a. If the stent is not completely expanded at any point along the lesion, balloon dilation (standard rapid exchange PTA technique) can be performed.
 - b. To dilate the stent, select an appropriately sized rapid exchange PTA balloon catheter and dilate the stent. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel.
 - c. Confirm that the stent is fully expanded, then remove the PTA balloon from the patient.
Caution: Use caution when crossing a deployed stent with any adjunct device.
Caution: Do not expand the stent beyond its nominal diameter.
3. Remove the embolic protection device, then the guidewire and introducer sheath from the body.
4. Close the entry wound.
5. Discard the delivery system, guidewire, embolic protection device, and introducer sheath.
Note: Physician experience and discretion determines the appropriate drug regimen for each patient.

11 MR Conditional

Nonclinical testing demonstrated that the stent in single and overlapped conditions is MR Conditional for stents up to 110 mm. A patient with this device can be scanned safely, immediately after placement, under the following conditions:

- Static magnetic field of 1.5 T or 3.0 T
- Maximum spatial gradient magnetic field of 4000 Gauss/cm or less (40 T/m)
- Maximum MR system-reported whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

MR image quality may be compromised if the area of interest is in the exact location or close to the position of the stent.

11.1 Artifact information

In nonclinical testing, the maximum artifact size as seen on the gradient echo pulse sequence at 3.0 T extends less than 5 mm relative to the size and shape of the stent. The lumen of the stent cannot be visualized.

12 Disposal

Caution: Dispose of the delivery system and accessories in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

13 Disclaimer of warranty

Important: This disclaimer of warranty does not apply in any countries where such a disclaimer is not permitted by law.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic and Medtronic logo are trademarks of Medtronic. ™™ Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic, Inc.
Medtronic B.V.
Medtronic Limited