



MINI WELL TORIC

PRODUCT TECHNICAL SHEET

MINI WELL TORIC

1. Trade Name

MINI WELL TORIC is a kit composed by an acrylic aspherical intraocular lens (EC 0123) and a dedicated injection system (EC 0482)

MINI WELL TORIC Series Code: X7560CZ

Individual product code:

MINI WELL TORIC: X7560CZPXXXYA

Where:

P = positive dioptric power

XXX = spherical equivalent between +7.00 to +32.00

Y = cylindrical dioptric power

2. Manufacturer

IOL: SIFI S.p.A.

Via Ercole Patti, 36

95025 – Aci S. Antonio (CT) – Italy

INJECTOR: Medice AG

Dornierstrasse 11

CH-9423 Altenrhein – Switzerland

3. EC Marking

IOL: EC 0123

INJECTOR: EC 0482

4. Compliance

Legislative Decree 46/97 (Medical Device Directive 93/42/CEE)

EN ISO 13485:2016 Medical Devices, Quality Management Systems - Requirements for Regulatory

Purposes UNI EN ISO 9001:2015 Quality System - Requirements

5. European Classification

IIB

6. General Product Description

MINI WELL TORIC is a progressive extended depth of focus intraocular lens intended to be implanted in the posterior chamber of the eye, more precisely in the capsular bag, to replace the natural crystalline lens for surgical correction of aphakia and presbyopia in adult patients, also compensating for corneal astigmatism.

MINI WELL TORIC is pre-loaded, comes in a single-piece design and is intended for mini incision. This surgical implant is made of a hydrophilic-hydrophobic copolymer with a UV-filtering chromophore and a refractive index of 1.46 at 35°C.

The anterior surface of MINI WELL TORIC is based on a patented proprietary optics that provides an extended depth of focus for progressive vision at all distances. The anterior surface of MINI WELL TORIC is aspheric, the posterior surface is toric and is designed to compensate pre-surgical astigmatism. The posterior surface has a double square edge which is intended to reduce the incidence of posterior capsule opacification.

MINI WELL TORIC has a biconvex optical shape for the entire dioptric range.



MINI WELL TORIC

PRODUCT TECHNICAL SHEET

This surgical medical device is a system consisting of an intraocular lens and a disposable injector pursuant to Directive 93/42 EC.

Physical characteristics:

Water content	25%	
Refraction index	at 35°C (in hydrated status)	1.46

Range of spherical dioptric powers:

Spherical Equivalent Min. (D)	Spherical Equivalent Max. (D)	Incremental Step (D)
+7.0	+32.0	0.5

Range of cylindrical dioptric powers on corneal plane:

Power Min. (D)	Power Max. (D)
+0.65	+3.0

Correlation table:

Cylindrical dioptric power (D) on IOL optical plane	Cylindrical dioptric power (D) on corneal plane*
1	0.65
1.5	1
2	1.3
2.5	1.7
3	2
3.5	2.3
4	2.7
4.5	3

*The cylinder values at the corneal plane shown in the table have been calculated based on the average pseudophakic eye.

Technical characteristics:

Diameter of the optical plate (mm) \varnothing_P	Total diameter (mm) \varnothing_T	N° of loops	Vaulting	Equivalent additional power (D)
6.0	10.75	4	5°	+3

Lens model, spherical dioptric power, spherical equivalent power, cylindric power, suggested A constant and main dimensions (diameter of the optical plate, total diameter, vaulting) are reported on the case and the primary and secondary packaging.

7.Composition/Device Material/Technical Characteristics

MINI WELL TORIC intraocular lens model is made of a hydrophilic-hydrophobic copolymer constituted by 2-HEMA (2-hydroxyethyl methacrylate) and EOEMA (2-ethoxyethyl methacrylate) with a chromophore that filters ultraviolet radiation and a refractive index of 1.46 at 35°C.

Component 2-HEMA is a hydrophilic monomer, while component EOEMA is a hydrophobic monomer. The two combined monomers form an acrylic polymer that compared to the individual monomers has a higher refractive index and better mechanical properties.

Furthermore, the polymer is ultrapure, i.e. without any residual acid, and with zero ionicity, thus excluding any possibility that the lens might become opaque due to calcium deposits.



MINI WELL TORIC

PRODUCT TECHNICAL SHEET

8. Indications

MINI WELL TORIC intraocular lens is intended for primary implantation for visual correction of aphakia and pre-existing corneal astigmatism in adults with or without presbyopia after removal of a cataractous lens and aphakia after lensectomy in presbyopic adults with corneal astigmatism who desire increased spectacle independence for distance and near vision. This implant must be positioned within the capsular bag.

9. Instructions for use

1. Appropriate implant selection
 - a. Examine the label of the complete package to verify the model, dioptric power or spherical equivalent, design and date of expiration of the product.
 - b. Suggested constant A: the constant A on the package is intended as a guideline and starting point for calculating the dioptric power of the implant. We suggest that surgeons should develop the constant value based on their clinical experience.
2. Removing the implant from its package
 - a. After opening the case, check that the information on it corresponds to that listed on the label in the lens package.
 - b. Open the lens external package and transfer the lens container in the sterile field in accordance with good practices in sterile field management.
 - c. Open the injector blister and place it on the sterile operating table.
 - d. Carefully open the lens container and remove the loading chamber in which it is lodged.
 - e. Insert the loading chamber holding the IOL into the specific injector slot.
 - f. Apply viscoelastic solution first through the front hole until the tip is filled and then through the rear hole towards the back of the chamber until the solution comes out, in order to obtain a barrier between lens and silicone cushion.
 - g. In addition, apply a drop of viscoelastic solution on the tip of the silicone cushion. It is recommended to use the SIFI viscoelastic solution, in order to allow the correct lens injection procedure.
 - h. Carefully remove the loading chamber block by lifting it from below with thumb finger and holding the upper part with index and middle fingers. Check that the lens is located at the centre of the loading chamber. Close the cartridge sides together until the "click-lock" mechanism engages.
 - i. Push forward the plunger softly and ensure that the silicone cushion correctly enters the loading chamber. Continue to push the plunger bringing the lens toward the tip of the injector. In the event the lens-injector system gets jammed, stop using the device.
3. Lens Injection
 - a. Introduce the injector tip into the tunnel, with the sloping side facing down and inject the intraocular lens into the bag.
 - b. After the injection, pull out the injector from the incision.

Attention: When the loading chamber is introduced into the injector, the plunger should be in the pulled-back position. The injector tip must be handled with care, especially during insertion of the loading charge into the injector. It must not be damaged in order not to jeopardise the injection operation.
4. Positioning of the MINI WELL TORIC lens
 - a. Rotate the lens within the capsular bag and place the indicators in correspondence with the axis of the astigmatism to be corrected. The indicators identify the lower power meridian forming an imaginary line that represents the axis of the positive cylinder; the IOL higher power meridian is located at 90° compared to the axis of the indicators. The marks on the axis of the cylinder of the MINI WELL TORIC lens must appear to be aligned with the corneal meridian to be corrected.

IMPORTANT. Before surgery, it is recommended that the eye to be operated should be marked by using a sclerocorneal marker while the patient is sitting. The marked points identify the reference axis (which may be 90° or 180°, as determined by the surgeon) to determine the refractive axis to be corrected in order to position the lens into the capsular bag. It is also recommended that you should consult the instructions provided by the



MINI WELL TORIC

PRODUCT TECHNICAL SHEET

Toric Calculator available on the website www.minitoriccalculator.com to determine the optimal positioning axis and the correct dioptric power.

Once the lens has been inserted into the capsular bag, carefully align the indicators of the MINI WELL TORIC intraocular lens with the astigmatism axis to be corrected. An incorrect alignment of the axis of the lens compared to the required positioning axis may compromise the corrective power. This incorrect alignment may derive from a keratometry or from an inaccurate marking of the cornea; from the incorrect positioning of the axis of the intraocular lens during surgery; from an unexpected surgically-induced alteration of the cornea; or from a physical rotation of the lens itself after the implant.

10. Sterilization

The intraocular lens package consists of a double-barrier assembly sterilized by steam; the injector and the blister in which it is lodged are sterilized with ethylene oxide.

Steam sterilization is compliant with standard UNI EN ISO 17665-1.

11. Packaging

MINI WELL TORIC intraocular lens is supplied in a hydrated state. The package consists of a double sterile barrier film.

The primary package of this IOL consists of a plastic container sealed by an aluminum foil containing a specific loading chamber where the IOL is inserted. The plastic container is filled with a sterile and apyrogenic saline solution. The plastic container is packaged in turn in a pouch.

The primary package of the injector consists of a Tyvek blister.

The case representing the SKU contains the lens and its injector, the "Product Information" leaflet, the implant card and the labels with the lens ID information.

12. Warnings

1. Patients with one or more of the following conditions are not eligible candidates for implantation of a progressive intraocular lens, as it can worsen the pre-existing condition, interfere with diagnosis or the treatment of a condition or potentially be a risk for the patient's sight:
 - a) zonular laxity
 - b) irregularity and decentration of capsulorhexis
 - c) recurrent inflammation of anterior or posterior segment of unknown etiology (chronic uveitis)
 - d) patients in whom the intraocular lens might interfere with the ability to observe, diagnose or treat posterior segment conditions
 - e) intraoperative challenges which could increase the risk of complications, such as severe bleeding, significant iris damage, uncontrollable elevated intraocular pressure or significant vitreous loss
 - f) lack of capsular support for implantation of the intraocular lens in the capsular bag
 - g) decompensated glaucoma
 - h) dystrophy of corneal endothelium
 - i) proliferative diabetic retinopathy
 - j) microphthalmus
 - k) pediatric patients
 - l) suspected microbial infection
 - m) patients without an adequate support for the lens from the posterior chamber or the zonule
 - n) bilateral congenital cataracts
 - o) history of retinal detachment or susceptibility to this condition
 - p) monocular patients
2. MINI WELL TORIC IOL must be positioned within the capsular bag and not in the ciliary sulcus.
3. MINI WELL TORIC IOL is indicated for the correction of pre-existing corneal astigmatism in the range recommended in the table.
4. The rotation of MINI WELL TORIC IOLs far from their selected axis can reduce astigmatism correction. Misalignment greater than 30° can increase post-operative refractive cylinder. If necessary, lens repositioning should be performed as soon as possible before lens encapsulation.



MINI WELL TORIC

PRODUCT TECHNICAL SHEET

5. Carefully remove all viscoelastic material from the capsular bag.
6. Any residual viscoelastic material can allow the lens to rotate causing misalignment of the lens with the intended axis of placement.
7. Other potential complications following IOL implantation are correlated to cataract surgery.

13. Precautions for use

1. In patients who have undergone implantation of the MINI WELL TORIC intraocular lenses, the auto-refractometer values may not correspond to the actual post-operative refractive value; it is therefore advisable to check the aforementioned values with the subjective examination, through the "fogging" technique.
2. The recent use of contact lenses affects the patient's refraction; therefore, in those wearing contact lenses, surgeons must define corneal stability without contact lenses before determining the power of the intraocular lens.
3. Do not re-sterilise the intraocular lens.
4. Do not reuse the intraocular lens; reuse compromises the sterility of the implant.
5. Do not use the system if the package is damaged or open.
6. Do not use after the expiration date printed on the package.
7. Do not use balanced saline solutions, Ringer lactate and/or hydroxypropyl methylcellulose during the loading procedures of the lens into the injector.
8. Improper handling of the preloaded system as well as improper bending techniques may cause damage to the haptics and/or the optical part of the lens. The surgeon must not attempt to implant lenses with a damaged optical plane or haptics. For a correct handling of the medical device follow the instructions for use being provided.
9. The intraocular lens implant requires adequate surgical skills.
10. Opening the protective packaging requires immediate use of the intraocular lens to prevent contamination or dehydration.
11. MINI WELL TORIC intraocular lenses must be kept at room temperature. It is recommended to implant them after keeping them at least 60 minutes at the temperature of the operating room. In this way it is possible to avoid that the thermal shock, resulting from an abrupt exposure to body temperature, causes a transient loss of transparency of the device.
12. MINI WELL TORIC must be implanted exclusively with the injector Accuject Pro 1.8-1P included in the package (Incision size: Wound-assisted technique: 1.8 mm; Into-the-wound technique: 2.0 mm; Into-the-bag technique: 2.2 mm).

14. Adverse events

The potential adverse events that may occur during or after cataract surgery with implantation of an intraocular lens may include but are not limited to:

General adverse events for intraocular lenses

1. Residual refractive error
2. Increase of intraocular pressure
3. Corneal oedema
4. Inflammation (Endophthalmitis, Iopion, Cystoid macular oedema)
5. Pupillary block
6. Retinal detachment
7. Displacement of the intraocular lens (Tilt and decentralisation)
8. Secondary surgical intervention (including repositioning of the implant, removal and replacement, PCO, or other surgical procedure)
9. Any other adverse event that leads to permanent damage to eyesight or requires surgery or a doctor to avoid this epilogue.

Adverse events for MINI WELL TORIC intraocular lenses

1. Increased visual disorders related to the optical characteristics of the intraocular lens:
 - a) Reduction of contrast sensitivity
 - b) Light artefacts (halos, glare or starbursts)



MINI WELL TORIC

PRODUCT TECHNICAL SHEET

2. Rotation of the positioning axis for the toric intraocular lens.

15. Safety Information

MINI WELL TORIC intraocular lens do not contain latex or phthalates.

16. Compatibility with other products

There are no interactions between the intraocular lens and other devices or therapies.

17. Validity

4 years

18. Storage and Handling Conditions

Store at room temperature.

If the lenses are exposed to temperatures other than those shown in the case, it is recommended to stabilize them at the operating room temperature at least 60 minutes before the implant.

19. Disposal

No special precautions are required for disposal.

END OF DOCUMENT

Rev. 7.0 of 31/03/2022