Herpes I & II IgG

**For the Qualitative Detection of
IgG Antibodies to Herpes simplex
I and II Virus in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Herpes I & II IgG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the qualitative detection of IgG antibodies to herpes simplex virus (HSV) types I and II in human serum, as an aid in determination of serological status to HSV I & II.

Catalog Number: **L2KHVG6** (600 tests)

Test Code: **HVG** Color: **Dark Pink**

Summary and Explanation

Herpes simplex virus (HSV) is an ancient and ubiquitous virus, known to cause acute and recurrent infections in humans. The virus enters the mucous membranes (ocular, genital or oral) and replicates locally.^{4,6} Infection of neonates during passage through the birth canal may result in neurological damage and death.⁵ In a small portion of infected individuals, the virus may enter the sensory root ganglion, resulting in latent, recurrent infections.¹

In the 1960s it was recognized that HSV consisted of two distinct types, HSV-I and HSV-II. HSV-I is considered to be primarily associated with ocular and oral infection, while HSV-II is considered to be a genital infection.⁶ However, HSV types I and II share several common antigens, and the use of specific monoclonal antibodies or restriction endonuclease mapping may be required to type individual strains.^{2,7}

Infections with HSV type I or type II can differ in their clinical manifestations and severity. The immune response of the host can play an important role in controlling the severity of primary or reactivated infections. Those at highest risk are neonates, who contract the infection during passage through the birth canal, and immunocompromised patients.³

While isolation of the virus in tissue culture is recommended for the diagnosis of active infections, serological testing can provide valuable information in the management of at-risk populations, such as pregnant women.

Principle of the Procedure

IMMULITE 2000 Herpes I & II IgG is a solid phase, two-step chemiluminescent enzyme-labeled immunoassay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Herpes I & II IgG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum.

Storage: 3 days at 2–8°C or for 6 months at –20°C.¹¹

Automatic Predilution Factor: 20.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface

antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹²⁻¹⁴

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

The bead is coated with *inactivated* HSV I & II viral antigen. However, caution is advised because of the possible presence of residual virus when working with, or disposing of, the materials supplied.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

The anti-HSV I & II IgG results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Herpes I & II IgG EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Herpes I & II IgG Bead Pack (L2HVG12)

With barcode. 200 beads, coated with HSV type I antigen (strain MacIntyre) and HSV type II antigen (G strain) (both antigens are partially purified, from infected cell lysates). Stable at 2–8°C until expiration date.

L2KHVG6: 3 packs.

Herpes I & II IgG Reagent Wedge (L2HVG42)

With barcode. 11.5 mL buffer solution. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to anti-human IgG antibody in buffer. Stable at 2–8°C until expiration date.

L2KHVG6: 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Herpes I & II IgG Adjustor (LHVGR)

4 mL of human serum with IgG reactive to HSV I & II, in buffer, with preservative. The Adjustor serves as the assay's Cutoff. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHVG6: 2 vials.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Herpes I & II IgG Controls (LHVGC1, LHVGC2)

LHVGC1 (Negative Control): 2.0 mL negative control containing human serum with IgG nonreactive to HSV I and II, with preservative. **LHVGC2 (Positive Control):** 2.0 mL positive control containing human serum with IgG reactive to HSV I and II, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHVG6: 2 sets.

The IMMULITE 2000 software performs automatic on-board dilution of control samples, and the results will be tracked in the QC database. Enter controls as controls.

For the current control ratio ranges, please refer to the Control insert.

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2KHVG6: 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KHVG6: 5 labels.

Kit Components Supplied Separately

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2IGZ2: 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2IGZ2: 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Sample transfer pipets; distilled or deionized water.

Assay Procedure

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

Recommended Adjustment Interval:

4 weeks.

Quality Control Samples: The Herpes I & II IgG Controls supplied with the kit should be used as quality control material to monitor the performance of the assay at the cutoff range. The Positive Control is used to validate the IMMULITE 2000 Herpes I & II IgG assay at a critical level when determining the serological status of a patient.

For the current control ratio ranges, please refer to the Control insert.

Additional controls may be tested in accordance with guidelines or requirements of local, state and/or federal regulations or accrediting organizations.

If the Herpes I & II IgG Controls do not yield their expected results, then readjustment is indicated providing other means of assessing and correcting the unacceptable control results have failed to identify and correct the problem. The procedure for readjustment and for assessing and correcting unacceptable control results is described in detail in the IMMULITE 2000 Systems Operator's Manual.

Patient results should not be reported unless QC results are acceptable. If QC results are unacceptable, the patient samples from that run should be repeated.

Calculation of Cutoff and S/CO Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times \text{P1}}$$

Calculation and reporting of qualitative results (reactive / nonreactive / indeterminate) and s/co ratio results are handled automatically by the IMMULITE 2000.

The result for a sample is reported as "Indeterminate" if the counts per second for that sample fall within $\pm 10\%$ of the cutoff. The result is reported as "reactive" if the sample's counts are *above* the indeterminate range, and "nonreactive" if *below* this range.

Interpretation of Results

- A result of **"Reactive"** (ratio of ≥ 1.1) indicates that IgG antibodies to HSV I & II were detected in the patient sample. (The assay does not differentiate between antibodies to HSV I or HSV II.)
- A result of **"Nonreactive"** (ratio of < 0.9) indicates that IgG antibodies to HSV I & II were not detected in the patient sample.
- Any result of **"Indeterminate"** (ratio between 0.9 and < 1.1) should be retested. Samples which still test as "Indeterminate" should be tested by an alternate method, or a second sample should be taken — if possible — within a reasonable period of time (e.g., one week).

The presence of IgG antibodies to HSV I or II is an indication of previous exposure to the virus. A single specimen can only be used to estimate the serological status of the individual.

The magnitude of the measured results (cps) above the Cutoff is not indicative of the total amount of antibodies detected.

Expected Values

Individuals infected with HSV I or II may not exhibit detectable levels of IgG antibody in the early stages of infection. Levels of IgG antibody to HSV begin to rise one or two weeks after primary infection. Peak levels are reached in 6 to 8 weeks, then gradually decline with time. Antibody levels may decrease to undetectable levels between periods of reactivation; hence, a patient may have a result of "Nonreactive" despite prior exposure to the virus.^{1,3,6} Reactivation of the virus may or may not result in a detectable level of antibodies in the IMMULITE 2000 Herpes I & II IgG assay.

Prevalence may vary due to geographical or population differences. Therefore, each


laboratory should establish its own reference ranges.

Studies with presumed healthy, asymptomatic subjects were conducted with the IMMULITE Herpes I & II IgG assay at two locations in the US. At one location (Site 1) in the northwestern US, the subjects consisted of 18 individuals undergoing a pre-employment screening (10 females and 8 males), with ages ranging from 26 to 42 years (samples were stored at 2–8°C if the assay was performed within 4 days, otherwise samples were frozen at –20°C). At a second location in the southwestern US, the subjects consisted of 100 blood donors (54 females and 46 males), with ages ranging from 17 to 71 years (samples were collected and frozen).

The prevalence of HSV I & II reactivity was found to be 78% (Site 1) and 76% (Site 2) in these two distinct geographical locations in the United States.

These study sites also included 132 (Site 1) and 71 (Site 2) pregnant women, with ages ranging from 16 to 43 years (Site 1) and from 14 to 48 years (Site 2). The prevalence of HSV I & II reactivity in this specific population at these two sites was 66% and 86%, respectively.

These prevalences are similar to those reported in the literature.⁹ However, prevalence may vary due to geographical or population differences. The distribution of patients in the two populations from these two independent sites is presented in the chart below.

Limitations

The results of the test must be taken within the context of the patient's clinical history, symptomology and other laboratory findings.

The presence of IgG antibodies in a single specimen is not sufficient to distinguish between active or past infection.

Primary or acute HSV infection should be confirmed by the culture of HSV from the patient.^{6,7}

Patients suspected of having primary or active infection should be tested for the presence of IgM antibodies to HSV.

For the determination of seroconversion from *nonreactive* to *reactive*, two serum samples should be drawn three to four weeks apart, during the acute and convalescent stages of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

Individuals with acute HSV infection may not exhibit any detectable IgG antibodies at the early stage of infection.

A rise in HSV I & II IgG antibody level may occur in patients with varicella-zoster virus (VZV), Epstein-Barr Virus (EBV) or cytomegalovirus (CMV) infections due to antigenic crossreactivity within the herpesvirus family.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from neonates, cord blood, or pretransplant patients.

No prozone effect has been demonstrated by antibody levels approaching 100,000,000 cps.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of

interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Herpes I & II IgG results are expressed as a signal-to-cutoff ratio. Unless otherwise specified, all results were generated on serum samples collected in tubes without anticoagulants, gel barriers, or clot-promoting additives.

Precision: Samples were assayed in duplicate in 40 runs for a total of 80 replicates. (See "Precision" table.)

Specificity: Studies were conducted in-house to determine crossreactivity or interference which may occur when using the IMMULITE Herpes I & II IgG assay. Serum samples were obtained from reference laboratories and characterized by Indirect Fluorescent Assay (IFA) or Anti-Complement Immunofluorescence (ACIF) as containing IgG antibodies to herpesviruses, including cytomegalovirus (CMV), varicella-zoster virus (VZV), Epstein-Barr virus (EBV), as well as *Mycoplasma pneumoniae*, and parvovirus B19. Ten samples of each type were assayed using the IMMULITE Herpes I & II IgG kit and commercially available HSV I and HSV II IgG assay kits (Kit C). Due to the high incidence of HSV I and HSV II in the general population, crossreactivity studies utilizing *only* nonreactive HSV I and HSV II serum samples were not performed.

The results indicate that approximately 70 to 80% of the samples were reactive for HSV I & II, reflecting the high reactivity rate (up to 80%) in the adult population.¹⁰ It is possible that the additional 20% reactive rate shown in the table below is due to cross-reactions. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 539 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST® vacutainer tubes. All samples were assayed by the IMMULITE 2000 Herpes IgG procedure, with the following results.

(Heparin) = 0.96 (Serum) + 0.64
 $r = 0.990$

(SST) = 0.94 (Plain Tubes) + 0.68
 $r = 0.994$

(EDTA) = 0.95 (Serum) + 0.42
 $r = 0.995$

Means:
18.7 (Serum)
18.5 (Heparin)
18.3 (SST)
18.1 (EDTA)

Method Comparison: The assay was compared to IMMULITE Herpes I & II IgG using samples from 227 blood donors. All samples were freshly collected and frozen for analysis. Four samples fell close to the cutoff. (See "Method Comparison" table.)

Clinical Performance

Based on the correlation between the IMMULITE Herpes I & II IgG and IMMULITE 2000 Herpes I & II IgG assays (see Method Comparison), the customer can reasonably expect comparable clinical performance.

The clinical performance of the IMMULITE Herpes I & II IgG assay was studied in-house using a CMV/HSV proficiency panel obtained from the United States Centers for Disease Control and Prevention (CDC).

The in-house retrospective study was conducted by assaying the CDC serum panel with the IMMULITE Herpes I & II IgG assay and sending the results to the

CDC for unmasking. The results are presented as a means to convey further information on the performance of this assay with a masked, characterized serum panel. This does not imply an endorsement of the assay by the CDC.

The panel consisted of 72% reactive and 28% nonreactive samples. The IMMULITE Herpes I & II IgG assay demonstrated 97% total agreement with the CDC results. Of the results obtained in-house, there was 99% agreement with reactive specimens and 93% agreement with the nonreactive specimens.

Studies on the clinical performance of the IMMULITE Herpes I & II IgG assay and two other commercially available systems (Kit A and Kit B) were also conducted nonconcurrently at two independent sites.

Site 1

In a prospective clinical study conducted at a medical center located in the northwestern United States, a total of 200 serum specimens were tested. Of the 200 specimens, 18 were obtained from presumed healthy, asymptomatic individuals undergoing a pre-employment screening, 132 were obtained from pregnant women and the remainder were obtained from patients with miscellaneous diseases and conditions (AIDS, heart disease, immunocompromised, kidney transplant/dialysis). There were 33 male and 167 female subjects, with ages ranging from newborn to 73 years. All 200 samples were stored at 2–8°C if the assay was performed within 4 days, otherwise samples were frozen at –20°C.

All 200 specimens were evaluated with the IMMULITE Herpes I & II IgG procedure and Kit A, an HSV I & II IgG commercially available microplate enzyme-linked immunosorbent assay (ELISA).

IMMULITE Herpes I & II IgG			
Kit A	Reactive	Indeterminate	Nonreactive
Reactive	139	1	0
Indeterminate	0	0	1
Nonreactive	2	1	56

Agreement: 99.0%

Two specimens were identified as indeterminate by IMMULITE Herpes I & II IgG, and one by Kit A.

The two discrepant cases were repeated with both IMMULITE Herpes I & II IgG and Kit A, with the repeat results remaining the same as those obtained in the initial testing.

Site 2

In the second outside prospective clinical study conducted at a university clinical virology laboratory in the southwestern United States, a total of 202 serum specimens were tested. Of the 202 specimens, 100 were from presumed healthy, asymptomatic individuals (blood donors), 71 were from pregnant women undergoing rubella screening and 31 were from patients undergoing infectious disease serological screening (rubella, CMV, toxoplasma, AIDS). There were 57 male and 145 female subjects, with ages ranging from 2 days to 71 years.

All 202 specimens (previously collected and frozen by the investigators) were evaluated with the IMMULITE Herpes I & II IgG procedure and Kit A, an HSV I & II IgG commercially available microplate enzyme-linked immunosorbent assay (ELISA).

Kit A	IMMULITE Herpes I & II IgG		
	Reactive	Indeterminate	Nonreactive
Reactive	153	0	6
Indeterminate	1	0	0
Nonreactive	2	3	37

Agreement: 96.0%

Four specimens were identified as indeterminate either by IMMULITE or by Kit A.

The eight discrepant cases were investigated by repeat assays with both IMMULITE Herpes I & II IgG and Kit A, and a third method, Kit B (a commercially available microplate ELISA kit). IMMULITE repeat tests confirmed the initial IMMULITE results for 7/8 discrepant cases. Kit A repeat tests confirmed 4/8 initial Kit A results. The third method, Kit B, agreed with 2/8 discrepant IMMULITE cases and agreed with Kit A in 6/8 discrepant cases. The four indeterminate cases were also repeated with the initial method that gave the indeterminate results. The repeat testing resolved the indeterminate status of these four cases, and IMMULITE Herpes I & II IgG was in agreement with Kit A in 3/4 cases.

References

- 1) Adam E. Herpes virus infection. In: Glaser R, Gotlieb-Stematski R, editors. Human herpesvirus infections. New York: Marcel Dekker, 1982: 1–55. 2) Goldstein LC, Corey L, McDougall JK, Tolentino E, Nowinski RC. Monoclonal antibodies to herpes simplex virus: use in antigenic typing and rapid diagnosis. *J Infect Dis* 1983;147(5):829–37. 3) Lycke E, Jeansson S. Herpesviridae: Herpes simplex virus. In: Laboratory diagnosis of infectious diseases: principles and practice. Vol II. Viral, rickettsial, and chlamydial diseases. New York: Springer-Verlag, 1988. 4) Moseley RC, Corey L, Benjamin D, Winter C, Remington ML. Comparison of virus isolation, direct immunofluorescence and indirect immunoperoxidase techniques for detection of genital herpes simplex virus infection. *J Clin Microbiol* 1981;13:913–19. 5) Smith IW, JF Peutherer. Immunological diagnosis of herpes simplex virus. In: Young H, McMillan A, editors. Immunological diagnosis of sexually transmitted disease. New York: Marcel Dekker, 1988: 371–401. 6) Stagno S, Whitely RJ. Current concepts: herpes virus infections of pregnancy. *N Engl J Med* 1985;313:1270–4, 1327–30. 7) Zhao L, Landry ML, Balkovic ES, Hsiung GD. Impact of cell culture sensitivity and virus concentration on rapid detection of herpes simplex virus by cytopathic effects and immunoperoxidase staining. *J Clin Microbiol* 1987;25(8):1401–5. 8) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 9) Baltz ML, Searcy RL. Clinical significance and advanced serologic diagnosis of ToRCH infections. *Am Clin Lab* 1994;March/April:18–23. 10) Arvin AM, Prober CG. Herpes simplex viruses. In: Murray PR, et al, editors. Manual of clinical microbiology, 6th ed. Washington, DC: ASM Press, 1995: 876–83. 11) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358. 12) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8. 13) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 14) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (ratio)

	Within-Run ¹			Total ²	
	Mean ³	SD ⁴	CV ⁵	SD	CV
1	0.24	0.03	12.5%	0.04	16.7%
2	0.25	0.07	28.0%	0.07	28.0%
3	0.85	0.06	7.1%	0.14	16.5%
4	1.01	0.05	5.0%	0.07	6.9%
5	1.24	0.17	13.7%	0.25	20.2%
6	2.03	0.16	7.9%	0.22	10.8%
7	2.26	0.11	4.9%	0.15	6.6%
8	45.7	2.44	5.3%	5.00	10.9%

Method Comparison

IML	IML 2000		
	Reactive ¹	Indeterm ²	Non-reactive ³
Reactive	167	0	0
Indeterm	0	2	0
Nonreactive	0	2	56

Agreement: 99.1%⁴

Specificity

Sample Reactive To: ³	Number Tested ¹	Reactive to HSV I & II ²	
		By IML ⁵	By Kit C ⁶
Cytomegalovirus (CMV)	10	2	3
Varicella-Zoster Virus (VZV)	10	8	9
Epstein-Barr Virus (EBV)	10	7	10
Mycoplasma pneumoniae	10	9	9
Parvovirus B19	10	8	7

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.
© 2009 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB