Medtronic

RapidCross™

RX PTA Catheter • RX катетър за РТА • Katétr typu RX pro PTA • RX PTA-kateter • RX-PTA-Katheter • Καθετήρας ΡΤΑ (Διαδερμικής διαλυτικής αγγειοπλαστικής) RX • Catéter para ATP RX • PTA RX-kateeter • RX-PTA-katetri • Cathéter pour angioplastie transluminale percutanée (ATP) à échange rapide (RX) • Kateter RX za PTA • RX PTA-katéter • Kateter PTA RX • Catetere RX per PTA • RX TTA катетери • RX PTA 카테터 • RX PTA kateteris • RX PTA katetrs • Катетер RX РТА (катетер за перкутана транслуминална ангиопластика) • RX PTA-katheter • RX PTA-kateter • Cewnik do PTA przystosowany do szybkiej wymiany (RX) • Cateter RX para PTA • Cateter RX pentru ATP • Катетер RX для ЧТА • Katéter RX PTA • Kateter RX PTA • RX PTA kateter • PTA-kateter av RX-typ • RX PTA Kateteri • Катетер RX для ЧТА • Ống thông PTA RX

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttööhjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөніндегі нұсқаулар • 사용 지침 • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

! USA

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
English

1 Device description

The RapidCross RX PTA Catheter is a rapid exchange (RX) coaxial lumen catheter (1). The catheter is compatible with 0.014 in (0.36 mm) guidewires (2). A semi-compliant inflatable balloon (3) is mounted at the distal end. The distal atraumatic tip is tapered (4). The distal portion of the catheter has a lubricious hydrophobic coating (5). The luer on the manifold (6) is used to inflate and deflate the dilatation balloon, using a mixture of contrast medium and saline solution. The balloon has 2 radiopaque marker bands (7) for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating or working section of the balloon (8). The 150 mm and 210 mm devices have 2 additional radiopaque marker bands (9), which denote the middle of the balloon body. A guidewire lumen starts at the guidewire port (10) located 35 cm from the catheter tip and extends to the distal tip. The catheters have depth marks for use as a reference during catheter insertion. The devices with a useable length of 90 cm (11) have proximal depth marks printed on the proximal shaft at lengths of 55 cm (12) and 65 cm (13) from the distal tip. The devices with a useable length of 170 cm have depth marks at 90 cm and 100 cm from the distal tip.

Figure 1. RapidCross 0.014 in (0.36 mm) rapid exchange PTA balloon dilatation catheter

The RapidCross RX PTA Catheter is available in multiple balloon sizes. Nominal balloon diameter and length are printed on the strain relief (14). Refer to the package label for information about catheter length and sheath compatibility.

All balloons expand to sizes above the nominal size at pressures greater than the nominal pressure. Consult the balloon compliance card that is packaged with the device for balloon diameters at given pressures.

2 Intended purpose

The intended purpose of the RapidCross RX PTA Catheter is to restore lumen patency and blood flow to peripheral arteries.

2.1 Intended patient population

The RapidCross RX PTA Catheter is intended for adult patients in need of percutaneous transluminal angioplasty to improve and maintain luminal diameter in peripheral arteries.

This device is not intended for pediatric patients.

2.2 Indications for use

The RapidCross RX PTA Catheter is indicated to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, and infrapopliteal arteries.

2.3 Contraindications

The RapidCross RX PTA Catheter is contraindicated for use in the coronary arteries or in the neurovasculature, or when unable to cross the target lesion with a guidewire.

2.4 Clinical benefits

The clinical benefits of restoring lumen patency and blood flow with the RapidCross RX PTA Catheter in the treatment of peripheral arterial occlusions or lesions include: resolution of arterial disease related symptoms, improvement in quality of life, and increased mobility.

2.5 Intended users

Interventional procedures using the RapidCross RX PTA Catheter should only be performed by physicians who have experience with interventional techniques in the vascular system.

This device should only be used within a clinical setting in a sterile environment.

2.6 Performance characteristics

The RapidCross RX PTA Catheter is a rapid exchange catheter with a balloon that is bonded onto the distal tip. The overall design and material selection allow for treatment of occlusions or lesions in the peripheral arteries. The materials provide kink-resistant pushability to easily navigate tortuous anatomy, and track to distal lesions. The balloon is inflated to its predetermined diameter and exerts a radial force to establish patency.

3 Warnings

- Through contact with the device, patients may have limited exposure (< 24 hours) to cobalt, which is considered carcinogenic, mutagenic, or toxic to reproduction (CMR).
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.
- Use the catheter with caution for procedures involving calcified lesions.
- Ensure that the guidewire port remains inside the introducer sheath or guide catheter while advancing the catheter to the treatment site. If the guidewire port is advanced outside of the sheath or guide, use caution to avoid guidewire prolapse. If prolapse cannot be resolved, remove the catheter and the guidewire together to prevent potential damage to the device or vessel walls. If resistance is felt, remove the sheath together with the catheter and guidewire.
- Use direct fluoroscopic observation while manipulating the catheter in the vascular system.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.
- Never use air or any gas to inflate the balloon.
- Do not manipulate the balloon when it is inflated.
- Consider precautions to prevent or reduce clotting. Physician experience and discretion determines the appropriate anticoagulation therapy for each patient.

4 Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the sterile packaging has been unintentionally opened or the device is damaged.
- Do not use contrast media that are contraindicated for intravascular use with this device.
- Do not wipe the catheter surface with dry gauze.
- Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.
- To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Ensure that the guidewire is in place before changing the position of the balloon catheter.
- Balloons should not be inflated in excess of the rated burst pressure.
- Larger and longer models of the catheter may exhibit longer deflation times, particularly on long catheter shafts.
- Do not retract the catheter unless the balloon is fully deflated under vacuum. Prematurely retracting the catheter can cause vessel trauma or product damage.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

5 Potential adverse events

The potential adverse events (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

Table 1. Potential adverse events

Allergic reaction	Fever	Regional ischemia
Aneurysm	Hematoma	Sepsis
Arrhythmia	Hemorrhage	Shock
Arteriovenous fistula	Hypertension or hypotension	Stenosis or restenosis
Bleeding requiring transfusion	Infection	Stroke
Contrast medium reaction or renal failure	Inflammation	Transient ischemic attack
Death	Intraluminal thrombus	Venous thromboembolism

Table 1. Potential adverse events (continued)

Embolism	Myocardial infarction	Vessel injury (such as dissection, perforation, or rupture)
Endocarditis	Pseudoaneurysm	Vessel spasm

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

6 How supplied

The catheter is packaged with a protective sheath positioned over the balloon.

The catheter is packaged with a protective wire that is positioned in the guidewire lumen.

The device is supplied sterile and is sterilized using ethylene oxide.

Warning: This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

7 Storage

Store the device in a dry area away from sunlight.

Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

8 Device materials

Warning: Through contact with the device, patients may have limited exposure (< 24 hours) to cobalt, which is considered carcinogenic, mutagenic, or toxic to reproduction (CMR).

Table 2. CMR substances

Substance	CAS number	Concentration (% weight by weight)
Cobalt (Co)	7440-48-4	> 0.1% w/w

9 Device selection

Consider the following details when selecting a catheter:

- The diameter of the expanded balloon should not exceed the diameter of the artery immediately distal to or proximal to the stenosis.
Warning: To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Verify that the selected accessories accommodate the balloon catheter based on sizing information on the packaging label. Refer to the product labeling for size information.
Caution: Carefully inspect the sterile package and the device before use. Do not use the device if the sterile packaging has been unintentionally opened or the device is damaged.

10 Recommended items

Prepare the following items using sterile technique:

- 10 mL syringe filled with sterile heparinized saline
- Three-way stopcock
- Contrast medium (The standard inflation medium is 1:1 mixture of contrast medium and sterile saline.)
Caution: Do not use contrast media that are contraindicated for intravascular use with this device.
- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

11 Preparation


- Carefully remove the catheter from the inner package.
- Remove and discard the protective sheath from the balloon and the protective wire from the guidewire lumen.
Caution: Do not wipe the catheter surface with dry gauze.

To displace air:

Note: As supplied, the balloon and the guidewire lumen of the catheter contain air. This air must be displaced to make certain that only liquid fills the balloon while the catheter is in the bloodstream.

- Flush out the lumen for the guidewire with heparinized saline solution by carefully inserting the provided flushing needle into the distal tip of the balloon catheter. Flush the lumen. Do not discard the flushing needle until the end of the procedure, in case additional flushing is required.

Figure 2.

- Flushing needle used to flush guidewire lumen

- Mix equal volumes of contrast medium and normal saline.
- Fill a 10 mL syringe with approximately 4 mL of the mixed contrast medium solution.

Warning: Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.

Caution: Do not use contrast media that are contraindicated for intravascular use with this device.

- Evacuate air from the balloon segment:
 - Attach the syringe to the inflation luer. Hold the syringe pointing down with the plunger pointed upward.
 - Apply negative pressure and aspirate the lumen for 15 seconds. Slowly release the pressure to neutral, so that the contrast medium fills the shaft of the dilatation catheter.
 - Disconnect the syringe from the inflation luer of the dilatation catheter.
 - Remove all air from the syringe barrel. Reconnect the syringe to the inflation luer of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
 - Slowly release the device pressure to neutral.
Note: On first use, avoid any positive pressure to avoid unfolding and partially inflating the balloon outside of the vessel.
- Repeat steps 4b to 4e as necessary to remove all air from the system.
- Replace the syringe with a manometer-controlled dilatation system. Do not introduce air into the system.

12 Insertion and dilatation

- If applicable, use an appropriately sized introducer sheath to introduce the catheter percutaneously. Refer to the product labeling for sizing information.
Caution: Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to product labeling for sizing information.
Warning: Use the catheter with caution for procedures involving calcified lesions or synthetic vascular grafts.
- Backload the distal tip of the prepared catheter over a prepositioned 0.014 in (0.36 mm) guidewire which has been placed at the treatment site.
Warning: Ensure that the guidewire port remains inside the introducer sheath or guide catheter while advancing the catheter to the treatment site. If the guidewire port is advanced outside of the sheath or guide, use caution to avoid guidewire prolapse. If prolapse cannot be resolved, remove the catheter and the guidewire together to prevent potential damage to the device or vessel walls. If resistance is still felt, remove the sheath together with the catheter and guidewire.
Warning: Use direct fluoroscopic observation while manipulating the catheter in the vascular system.
Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.
Caution: To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
Note: To improve catheter movement, thoroughly wipe the guidewire with saline soaked gauze before using the guidewire in the catheter.
- Position the catheter with the center of the balloon in the stenotic area of the vessel. The radiopaque marker bands indicate the working length of the balloon. The 150 mm and 210 mm devices have 2 additional radiopaque marker bands, which denote the middle of the balloon body.
Caution: Ensure that the guidewire is in place before changing the position of the balloon catheter.
- When the radiopaque marker bands are appropriately positioned, inflate the balloon to dilate the target area.
Warning: Never use air or any gas to inflate the balloon.
Warning: Do not manipulate the balloon when it is inflated.
Caution: Balloons should not be inflated in excess of the rated burst pressure.
- Using the inflation device, deflate the balloon by aspirating the catheter.
Caution: Larger and longer models of the catheters may exhibit longer deflation times, particularly on long catheter shafts.
- Withdraw the catheter while maintaining a vacuum in the balloon.
Caution: Do not retract the catheter unless the balloon is fully deflated under vacuum. Prematurely retracting the catheter can cause vessel trauma or product damage.
Note: Multiple balloon inflations and deflations may cause resistance during device withdrawal. If resistance is encountered, use a gentle clockwise twisting motion to ease withdrawal through the sheath. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath together, while maintaining vessel access with the guidewire.
- Use angiography to verify dilatation results.

13 Disposal

Caution: Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

14 Disclaimer of warranty

Important: This disclaimer of warranty does not apply in any countries where such a disclaimer is not permitted by law.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic and Medtronic logo are trademarks of Medtronic. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
**Medtronic, Inc.**
Medtronic Limited
Medtronic B.V.