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**Functional Anatomy** ✕

	Type	Description
	Birth Simulation	Observe the key stages of vaginal birth to understand organ movements and fetal positioning during delivery. [ Female - Full Body (Asian) ]
	Pregnancy	Visualize anatomical changes in the mother and fetus during the 31st week of pregnancy. [ Female - Full Body (Asian) ]
	Cardiology	Visualize cardiac physiology using a beating heart and ECG. [ Female - Regional ]
	Kinesiology	Observe the interactions between bones, ligaments, and muscles during movement of the shoulder, hip, and knee joints. [ Male - Full Body (Asian) ]
	Neural Pathways	Visualize the sensory and motor pathways that transmit information between the central nervous system and peripheral structures of the body. [ Female - Full Body (Asian) ]
	Catheterization	Learn the different steps required to insert a vascular catheter. [ Female - Full Body (Asian) ]

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
The image shows a dark-themed dialog box titled "Anatomage Share" with a close button (X) in the top right corner. It contains two input fields: "Email:" and "Password:", each with a clear button (X) on the right. Below the fields is a paragraph of text: "Connect your Anatomage Table online and sign in to your Anatomage Share account. Anatomage Share allows easy storage and sharing of preset and DICOM files between the Anatomage Table, Anatomage VR, and other Anatomage products. To create an account, visit [www.anatomageshare.com](http://www.anatomageshare.com)." Below this is another paragraph: "All content downloaded from Anatomage Share is for informational and educational purposes only and does not constitute medical advice." At the bottom center is an "Okay" button.

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40

## INSTALLATION INSTRUCTIONS

The Anatamage Table application is available from Anatamage. The software is distributed by downloading an installer and requires a license USB to operate. The installer contains both the application files and demonstration content (Navigation, Cadavers, Functional Anatomy, Case Library, Histology, Curriculum, and Prosection). Be sure you have the license USB and your Table meets the minimum system requirements.

1. Insert license USB into PC workstation and run Anatamage Table EDU 10.0 installer. Workstation should be connected to a network with internet connectivity.
2. Open **Anatamage Table EDU 10.0 Upgrade** installer, double click "autorun" and follow on-screen instructions to complete installation.
3. Launch Table application and enter the Authorization Code to activate software license.

## CONTROLS

The following section discusses controls for the Table software. For touchscreen devices, please refer to your specific touchscreen hardware manufacture's calibration instructions and verification procedures before using with Table EDU 10.0.

Touch Commands	
Right click	Touch and hold.
Selecting icons	Tap icon to select. If compatible multi-icons are shown, use a second tap to select desired multi-icon. A double-tap on the same icon will open the icon.
On-screen keyboard	Tap the keyboard icon on the toolbar next to the windows icon. This will open the onscreen keyboard.

Keyboard-Only Commands	
Exiting Full-Screen and viewing application on single display monitor	<p><b>Step 1:</b> Press <b>F11</b> on keyboard or <b>FN + F11</b> on on-screen keyboard.</p> <p><b>Step 2:</b> Press the Windows key and the left/right arrow to snap application window to left/right display monitor.</p>

**WARNING:** Resizing the application window from full-screen to a single monitor will cause the user interface and scan to be rescaled based on the new application window size.

**KEYBOARD AND MOUSE CONTROL**


Number of Touches	Movement	Result	Description
<b>Volume Viewing:</b>			
Left Click	Drag	Rotate	Rendering will rotate about the scanning region's geometric center point.
Shift + Left Click	Drag	Pan	Rendering will pan in the dragged direction.
Ctrl + Left Click	Drag up/down	Zoom in/out	Rendering will become larger or smaller.
Space + Left Click	Drag up/down	Spin clockwise/ counterclockwise	Rendering will spin clockwise or counterclockwise about its geometric center point.
Scroll Wheel	Roll up/down	Adjust Clipping Plane	Can adjust clipping plane by scrolling through volume rendering in parallel with initial cutting plane.
<b>Slice Mode Viewing:</b>			
Shift + Left Click	Drag	Pan	Rendering will pan in the dragged direction.
Ctrl + Left Click	Drag	Zoom in/out	Rendering will become larger or smaller.
Scroll Wheel	Roll up/down	Scroll through slices	Can scroll through cross-sectional slices of selected data.

## INTRODUCTION TO THE ANATOMAGE TABLE APPLICATION

### LAUNCHING TABLE EDU 10.0 APPLICATION


From the desktop, double-tap (double-click) the Table icon to launch the Anatamage Table application. Users will be shown the Application Toolbar below.

<b>Navigation</b>	Opens navigation tool to open Table content related to a specific subject.
<b>Cadavers</b>	Opens user interface for selecting full body Gross Anatomy data and Regional Anatomy scans.
<b>Functional Anatomy</b>	Opens user interface for selecting functional anatomy related content.
<b>Case Library</b>	Opens user interface for selecting educational clinical case data sets including Ct, MRI and segmented dog and cat.
<b>Histology</b>	Opens user interface for selecting histology slides.
<b>Curriculum</b>	Opens user interface for selecting curriculum views. Curriculum views, provided by Anatamage, are single scans with pre-made annotations and view sequences for teaching purposes.
<b>Prosection</b>	Opens user interface for selecting prosection data.
<b>Anatamage Share</b>	Opens sign in window for Anatamage created materials and other user content.
<b>Open File</b>	Opens multi-dicom (DCM) scan files or Invivo (INV) scan files. Opens user interface for establishing PACS connection and downloading files.
<b>Help</b>	Opens a user interface for selecting “how to” documents for aspects of Table.
<b>Minimize</b>	Minimizes the Table application. Available only when case is currently open.
<b>Exit Application</b>	Closes the Table application.
<b>Cancel</b>	Closes the Application toolbar. Available only when case is currently open.

## USER INTERFACE AND LAYOUT

*Red box outlines Rendering Window. Image rendering, surface models, and annotations will appear in this region. Region accepts keyboard, mouse, single-touch, and multi-touch controls.*

Icon	Description
	<p><b>Application Toolbar</b>                      Opens the Application toolbar. (Navigation, Cadavers, Functional Anatomy, Case Library, Histology, Curriculum, Prosection, Anatomage Share, Open File, Help, Minimize, Exit Application, and Cancel).</p>
	<p><b>Learning Assistant</b>                      Tap to view more information about a structure.</p>
	<p><b>Custom Navigation</b>                      Add the current cadaver, case, slide, or simulation to a custom Subject within the Navigation tool.</p>

← 15

**Layout**

Tap the Layout icon to display all Layout options. 3D Layout is selected by default.

**3D Layout**

Tap to view a volume rendering of the data.

**3D-2D Layout**

Tap to view a split-screen view of a volume rendering of the data and a cross-section in the axial, coronal, or sagittal plane.

**2D Flat Color**

Automatically flat color all structures in the cross section.


**2D Plane**

Change the plane of the 2D slice.

**Ultrasound Viewer**


Access the ultrasound viewer tool.


**2D Layout**

Tap to view one cross-section at a time. Users can view cross-sections in the axial, coronal, or sagittal plane.

← 6.3

The user can switch between a single 2D view and two different 2D layouts that contain three 2D slices each. The default slices for these split-screen layouts is Axial, Coronal, and Sagittal. The layouts can be customized to contain slices in Axial, Coronal, Sagittal, and Parallel planes.

← 6.4

**Histology Layout**

Tap to view a split-screen view of a volume rendering of the data and a histology slide. Any histology slides associated with a structure will be shown when the structure is selected from the volume rendering using a single tap.

6.5

6.5

Choose a prosection or CT/MR case. Select the icon for the content then select the slides icon to open the menu to choose a case.

← 6.6


**Flythrough**

Tap to view split-screen view of a volume rendering of the data and a flythrough data set. The user can load in default flythrough data sets when available or create custom flythroughs if desired.

6.7


### **Curved Planar Reformation**

Tap to view split-screen view of a volume rendering of the data and a curved planar reformation. Curved planar reformations can be exported.

### **Spline**

Tap to create a spline on the volume rendering of the data. The user can choose the points of the spline by tapping on a structure.

### **Curved Spline**

Tap to create a curved spline.

### **Straight Spline**

Tap to create a linear spline.

### **Pointer Tool**

Tap to open the Pointer Tool Dialog. Tap to select a particular pointer icon. Tap or drag in rendering window to move pointer.

The size and color of the pointer may be adjusted in the dialog.

### **Screen Capture**

Tap the camera icon to save a screen shot (.jpg, .png, or .bmp) of the application using three different options.

**Screen Capture with Icons**

Tap the camera with icons to save a screenshot (.jpg, .png, or .bmp) of the application that includes the Table icons.

**Screen Capture without Icons**

Tap the camera without icons to save a screenshot (.jpg, .png, or .bmp) of the application that does not include the Table icons.

**Cropped Screen Capture**

Tap the scissor icon and drag across screen to select an area and save a screenshot (.jpg, .png, or .bmp) of the application within the selected area.

**Pen Tool**

Tap to open the Pen Tool Dialog. Using second tap, select a particular Pen Tool.

Draw by dragging in rendering window.

**Predefined Draw Colors**

Tap to select a predefined Pen Tool color. From top left going clockwise: Red, White, Yellow, Blue.

Default width of Draw Stroke for all colors is 2.

**Custom Pen Tool**

The Custom Pen Tool allows the user to customize the color and width of a pen.

**Step 1:** Tap to select a particular pen preset.

**Step 2:** Tap the Custom Pen Tool Settings icon to adjust color and stroke width.

*The Pen Tool will save the latest setting used for each preset.*

**Eraser:** Tap to enable. Drag on display window to remove pen, text, or arrow marks.

**Text Tool:** Tap to activate and then tap on display window to place text. Use on-screen or external keyboard to enter text.


**Arrow Tool:** Tap to activate and then tap on display window to draw arrow tail. Tap again to draw arrow head.

**Undo:** Tap icon to undo last drawing action.

**Redo:** Tap icon to redo last drawing action.

**Minimize:** Tap icon to minimize Draw Tool dialog.

**Clear:** Tap icon to erase all drawings in the Rendering Window.

### Measurement Tool

Tap the icon to show the associated measurement icons.

The font size and color of the measurement value may be adjusted in the Adjust Text window under the Image Control Settings icon (p. 27).

Measurements can be made in both volume and slice mode viewing.


**WARNING:** The default measurement unit will be millimeters (mm). Any measurement that is less than 1 mm will be displayed in micrometers (um).

### Distance Measurement Tool

**Step 1:** Select distance measurement tool.


**Step 2:** Tap on two locations on the volume. Red dots will indicate the selected spots, and a line will appear between them with a distance measurement. To adjust a measurement spot, select and drag the landmark red dot.

### Angle Measurement Tool

**Step 1:** Select angle measurement tool.

**Step 2:** Tap on three locations on the volume. Red dots will indicate the selected spots, and an angle will appear between them with an angle measurement. Distance measurements of the angle's sides will also appear. To adjust a measurement spot, select and drag the landmark red dot.

### Area Measurement Tool

**Step 1:** Select area measurement tool.

**Step 2:** Tap on multiple locations on the volume. Red dots will indicate the selected spots, and a line will connect the dots to outline the area.

**Step 3:** Tap on the area measurement tool when finished and an area measurement will appear. To adjust a measurement spot, select and drag the landmark red dot.

16

**Curved Measurement Tool**

**Step 1:** Select curved measurement tool.

**Step 2:** Tap on multiple locations on the volume. Red dots will indicate the selected spots, and a line will connect the dots.

**Step 3:** Tap on the curved measurement tool when finished and a distance measurement will appear. To adjust a measurement spot, select and drag the landmark red dot.

**Delete Measurement**

**Step 1:** Select measurement. Specified measurement will turn red and become bold.

**Step 2:** Select icon to remove specified measurement.

**Clear All Measurements**

Select icon to clear all measurements from the volume.

*WARNING: Identification of anatomical landmarks and structures are limited in part to image resolution and subject to user error. To ensure correct identification of landmarks and other fine measurements, it is recommended that users plug-in and use a USB computer mouse and keyboard for the most accurate possible placement of measurement landmarks (red dots). All measurement landmarks, including those placed using the touchscreen interface, can be adjusted by selecting and dragging the landmark. It is the responsibility of the user to place or adjust the measurement landmark locations as needed for analysis.*


**Pin Tool**

Used for placing 3D pin models in volume viewing. To place a pin, select the pin style of choice and then select a place on the volume. The pin may be moved by first tapping on the pin to display a red box and again to display a yellow sphere. Holding and dragging the red box moves the pin's XYZ location while holding and dragging the yellow sphere rotates the angle of the pin. The lumbar needle and syringe within the pin menu are used the same way.

This icon can be used to import in a custom model that can be placed like the pins, lumbar needle, and syringe. An action


menu will appear which will allow you to manipulate the model. Double tapping on a model will open the models action menu also.

This icon will turn on/off the 3D widget which allows you to precisely move the model.

This icon allows you to map the model to a structure.


This icon will allow you to edit the dimensions (L x W x H) of the model.

### Volume Orientation

Tap the icon to show orientation icons.

### Coronal, Sagittal, Axial Views

Tap to select a particular orientation. From top left going clockwise: Coronal View, Sagittal View, Axial View.

*If an orientation icon is tapped a second time, the view will be flipped.*

**WARNING:** Image orientation is based on scanner/DCM definitions or redefined orientations from Invivo6 software.

### 90° Rotation

Tap once to rotate the image 90° clockwise.

**1:1 Life Size Scaling**

Tap icon to rescale image to life size.

*WARNING: Exact scaling depends on scan size, scan resolution, and hardware specifications.*

**Dissection Tool (for segmented cadaver data)**

The Dissection Tool allows the user to create custom cuts and remove structures, or parts of structures, within the sculpt area. After tapping the tool, the icon becomes highlighted signifying it is enabled and the Dissection Tool dialog appears.

**Curved Dissection**

Tap to activate the curved dissection tool, as indicated by a blue background.

**Step 1:** Hold and drag anywhere on the selected structure. A red line will appear outlining the sculpt area.

**Step 2:** Tap inside or outside the area to select dissection boundaries. The cut will be made only on the selected structure.

**Linear Dissection**

Tap to activate the linear dissection tool, as indicated by a blue background.

**Step 1:** Hold and drag anywhere on the selected structure. A red line will appear.

**Step 2:** Tap on either side of the line to select which side will be dissected. The linear cut will be made only on the selected structure.

**Point-To-Point Dissection Tool (for segmented cadaver and DICOM data)**

Use the point-to-point dissection tool to select two points on the cadaver and a cut diameter, then dissect within that area.

**Freehand Dissection Tool** (for segmented cadaver and DICOM data)

The Freehand Dissection Tool allows the user to create custom cuts on the segmented cadaver or scan volume. After tapping the tool, the icon becomes highlighted signifying it is enabled. To make a custom sculpt:

13

**Step 1:** Hold and drag the custom sculpt anywhere on the volume. A green line will appear outlining the sculpt area.

11

**Step 2:** Tap inside or outside the area to remove the visible volume on either the inside or outside of the sculpt outline, respectively.

13

To enable another freehand sculpt, tap the Freehand Dissection Tool icon again. Once the icon is highlighted, repeat steps 1 and 2 to make another freehand sculpt.

13 12

**Craniotomy Tool** (For DICOM data)

**Step 1:** Select craniotomy tool.

**Step 2:** Hold and draw custom sculpt across any cranial portion of the CT or MRI Scan.

**Step 3:** Release touch and software will automatically close the custom sculpt and perform craniotomy.

**WARNING:** The Craniotomy Tool is meant to be used as a Demonstration Tool only. The amount of volume removed from a scan is dependent on the scalar range used to view DICOM data.

14

**Undo:** Will undo the last structure removal action.

**Redo:** Will redo the last structure removal action.

Note: The undo and redo icons do not work for the Freehand Dissection Tool.


**Reset Dissection**

Tap to remove any dissections performed on the structure, restoring the volume to its original state and removing all operations.

**Clipping Plane Control**

Tap icon to display Clipping Plane Control all Clipping Plane icons as well as the Flip and Reset icons. Custom Clipping Plane is activated by default.

**Custom Clipping Plane**


**Step 1:** Tap scalpel icon with the straight line to activate custom clipping plane. The custom clipping plane is activated by default (indicated by a blue highlight.)

**Step 2:** Using one touch, drag anywhere across the volume to generate the custom clipping plane, defined by a blue line. The line will update to show the current clipping definition.

**Step 3:** Release touch to finish defining plane.

**Step 4:** Select a side of the plane to remove by tapping the volume on that side.

Repeat the above steps to create up to six (6) Custom Clipping Planes. After the sixth plane, the seventh plane will replace the first defined Custom Clipping Plane. The eighth will replace the second and so forth.

**Predefined Clipping Planes**

Select a predefined clipping plane by tapping on the icon. The volume will automatically be clipped in the designated direction. From top left going clockwise: Sagittal Plane, Coronal Plane, Axial Plane, Parallel Plane.

*Parallel defines the Table surface as the clipping plane.*

*Only one predefined clipping plane may be applied at a time.*

### Flip Clipping Plane

Tap icon to flip visibility between sides of the clipping plane. This affects all predefined clipping planes and the most recently defined custom clipping plane.

21

### Reset Clipping Plane

Tap icon to remove all applied clipping planes (custom and predefined).


### Clipping Plane Slider Bar

When a clipping plane, either a predefined or a custom plane, is active, the slider bar can be used to adjust the location of the plane. The slider bar can be used to adjust only the most recently defined clipping plane.

### Presets

Tap to show presets numbered 1-10. Use a second tap on one of the numbered icons to select a particular preset, or open the Presets Menu by tapping on the cogwheel icon. Users can create an unlimited number of presets using the Preset Menu. Presets will automatically switch the image in rendering window based on saved definitions.

*See Table 10.0 "How to Use Presets" for more information on how to create and save presets and use the Quick-View settings.*

### Anatomage Curriculum

A collection of premade images using the Anatomage cadavers. Search through systems and topics.

**Quiz Mode**

Select the quiz icon then select the icon to create a quiz or take a quiz. Set up the desired structure view and add annotations to define quiz questions.

Adjust any quiz settings (*enabled features, randomize questions, timing, scoring, etc.*) by tapping on settings.

Select the quiz type (Pick Me, Flashcards, Highlight, Multiple Choice, Game, or Tournament).

See Table 10.0 "How to Use Quiz Mode" for more information on how to create and take a quiz.

**Explore Tool**

Used for identifying segmentation and structures.

Tap the icon to open up the Explore Tool Dialog. Next, select a particular icon with a second tap. The icon will become highlighted to indicate it is enabled.

**Removal Tool**

Tap anywhere on volume. Structure closest to touch location will be removed from view.


**Split-Screen Removal Tool**

Tap to view removed structures on the right half of the screen. When the isolated structure appears on the right, you can turn on sync to have the two images rotate together. Turning off sync will allow the two images to rotate individually.

**Highlight Tool**

Tap anywhere on volume. Structure closest to touch location will be annotated and all other volume elements will be made transparent. Tap on the isolation icon again to remove all transparent structures. Tap a third time on the icon to add

structures back in again. Tap on name to enable more detailed annotations.

**Undo & Redo:**

After using the Removal or Highlight Tool, tap Undo to undo the last action or Redo to redo the last action.

The font size and color of the explore text may be adjusted in the Adjust Text window under the Image Control Settings icon (p. 27).

**Cardiovascular Tool**

Use these tools to visualize the cadaver's cardiovascular system.

**Blood Flow Tool**

Tap on a vein or artery. Blood flow will be simulated to or from the Heart. A second tap will apply a flat color to the vein or artery. A third tap will add annotations to all branches of the vein or artery.


**Vascular Grow Tool**

Tap on a vein or artery to see it grow.

**Play Icon**

Tap the icon multiple times to apply different views. You can play and pause this animation.

**Particle Flow Tool**

Tap once on the icon to turn on blood flow. A second tap will turn off the feature.

**Heart Beat Tool**

Tap once on the icon to turn on a beating heart. A second tap will turn off the feature.


**Visibility Control**

*Structures*

Tap icon to open Volume Visibility dialog to adjust rendering window images. User can add/remove structures in data sets or adjust volume rendering view presets for any DCM file data sets. See p. 33 for more information on using the Volume Visibility dialog.

*Annotations*

Tap “A” icon next to any entry to turn on all annotations under that entry.

**Brightness/Contrast Slider Bars**

*When viewing segmented cadaver data:*


Drag the upper slider bar right/left to add/remove large systems or structures from the volume rendering. Drag the lower slider bar left/right to add/remove the cardiovascular, nervous, and lymphatic systems from the volume rendering.

*When viewing DICOM data:*

Drag the upper slider bar right/left to increase/decrease the *Brightness* (density) of the volume in the Rendering Window. Shift the lower slider bar right/left to increase/decrease the *Contrast* of the volume in the Rendering Window.

**WARNING:** *Brightness and Contrast settings are dependent upon the volume rendering range defined in the Image Control Settings (p. 26).*

**Action Menu**

Tap structure and then action menu or double tap a structure to view action menu for selected structure. The functions of each tool within the Action Menu are outlined below.

**Structure On/Off**

Tap to toggle a structure on or off. A white check mark indicates that the structure is toggled on. A dashed check mark indicates that the structure is toggled off.

**Transparency**

Tap to adjust surface and volume transparencies (50% surface opacity, 0% volume opacity) for selected structure. A blue background indicates that a transparency is applied.

**Flat Color**

Tap downward arrow to choose a color from the drop-down menu. Tap the color button to highlight selected structure in chosen color.


**Annotations**

Tap to enable all annotations for selected structure. A blue background indicates that annotations are applied.


**Origin and Insertion**

Tap to display origin and insertion points for a selected bone. A blue background indicates that origin and insertion points are turned on. A grayed-out icon indicates no origin/insertion for a structure. Tap on any of the painted surfaces to turn on/off the annotations.

**Bony Landmarks**

Tap to display bony landmarks for a selected bone. A blue background indicates that the bony landmark points are turned on. A grayed-out icon indicates no bony landmarks for a structure. Tap on any of the painted surfaces to turn on/off the annotations.

**Spin**

Enable the ability to spin volume in a circle using two fingers.

**Structure Pronunciation:** The Table will read aloud a selected structure.

Turn this function on or off by tapping the icon to enable or disable. Tap the R/L icon to read only the structure name without identifying which side of the body it is on

**Dual Res Cadavers** (*Cadavers only*)

When enabled, the Table will automatically update to a higher quality image when zoomed in.

**User Modes**

Users can place the Table into different “locked” modes, which disable select features. Choose from premade modes or customize. Lock the setting with a password if desired. Master password – AnatoTest0

**Language Menu:** The user can toggle between languages.

*After the language has been changed, the user must adjust the system locale.*

**Case Information** (*For DICOM data, Histology, and Prosection*)

View scan and case information.

**Annotation Controls** (*for DICOM data and Prosection*)

Opens Annotation interface that allows for custom annotations to be made. These annotations can be added to a group and new annotations can be put in a specific group.

To place the annotation, tap a location on the scan. Location of annotations can be edited. Once the annotations are made, they can be turned on/off by checking the box or tapping the “A” at the bottom of the menu. A list of annotations can be exported and imported.


## NAVIGATION

Navigation is an easy way to open Table content related to a specific subject.

There are ten subjects with content:

- Cardiology
- Sports Medicine
- Interesting Pathology
- Optometry
- Nervous
- Dental
- Oncology
- Surgical Devices
- Developmental anatomy
- Respiratory

Each thumbnail image is outlined in a different color which represents where the Table scan is located within the Table:

## CADAVERS

(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).

Tap **Cadavers** on Application toolbar to open the *Gross Anatomy* and *Regional Anatomy* menus.


- Dialog lists all available Gross Anatomy and Regional Anatomy scans with resolution, size, and gender.
- Tap entry name, then tap **Open** or double tap entry name.
- Tap **Close** to close dialog.


## IMAGE ADJUSTMENT

### Segmented Cadaver Data

Selecting the Volume Visibility Control Icon will open the following dialog.

- Structures are organized into systems and categories.
- Tap Checkbox icon at the bottom of the system list to turn **all structures** On/Off.
- Tap Checkbox icon next to each entry to turn systems, categories, or individual structures On/Off.
- Tap a system or category to show associated subsystems. Selected entry name will be highlighted blue.
- Tap grayed color wheel icon next to entry names to adjust **Opacity**, **No Clip**, and **Color** settings for volumes and to toggle **Origin/Insertion** and **Bony Landmarks** for bones. If adjusted, gray color wheel icon becomes colored.

**Origin and Insertion:** Toggles On/Off Origin and Insertion for selected bone structure.

**Bony Landmarks:** Toggles on/off bony landmarks for selected bone structure.

**No Clip:** If enabled, structures cannot be sliced through.

**Transparency:** Auto adjust surface and volume transparencies (50% surface opacity, 0% volume opacity)

**Opacity:** Adjust slider bars to adjust surface and volume transparencies.

**Randomized Color**

Apply randomized colors to a structure, category, or system.

**Flat Color:** Choose a color from quick access or the drop-down color menu. Tap button to highlight structure in chosen color.

- L and R buttons in lower left-hand corner allow the user to toggle On/Off the left and right structures for select systems and categories.
- The counter-clockwise, curved arrow color wheel button restores the cadaver to its default view by removing any Opacity, No Clip, Color settings, Origin/Insertion and Bony Landmarks.
- The A button in the lower right corner expands the Volume Visibility dialog to include Annotations. When the dialog is expanded: the circled A button toggles On/Off Annotation Visibility, the downward and upward arrow buttons allows the user to import and export custom annotations, and the counter-clockwise, curved arrow allows the user to restore to default annotations. See p. 35 for additional details on Annotations.
- Search bar in lower right corner allows user to search for a particular structure. Tap **X** to clear all search terms.

## ANNOTATIONS

Annotations can be activated in the Volume Visibility Dialog. Annotations are linked with clipping planes. Chosen annotations will be displayed, but lead lines will only appear when associated volume is close to the clipping plane. If no clipping plane is defined, lead lines will always appear.

- Tap the “A” button at the bottom right of the Volume Visibility Dialog to expand the dialog to include Annotations. This includes structure annotation, Origins/Insertions, Bony Landmarks.
  - Tap the arrow icon next to Annotations, Bony Landmarks, or Origins/Insertions to expand that annotation list.
  - Tap the circled “A” button at the lower right of the dialog to show/hide all enabled annotations.
  - Check the box next to each entry to enable that specific annotation.
  - In the Annotation list, tap “New Annotation” to create your own annotation. Type out the desired annotation text and tap anywhere on the volume to select annotation coordinate.
- Tap on the settings tab and tap “Adjust Text” to adjust where annotation texts will appear (*Top, Bottom, Right, Left*) and how annotation text will appear (*Size, Color*).
- Annotations saved with Invivo will appear with *Comment* and *Marker* as System and Category, respectively.
- User can adjust coordinates of annotations by tapping the Ellipsis Button to the right of the annotation and then “Edit Location”. A dialog will appear prompting “Edit Location” and the user can tap anywhere on the structure to change its location. When the Annotation (.csv) file is exported, the new location(s) will be updated.

**WARNING:** Saving the new coordinates for an annotation will overwrite the associated information on the currently loaded .csv annotation spreadsheet. A back up annotation spreadsheet is available on the Table desktop.

### Load Default Annotations

Load in default annotations (.csv file) from the Presets, Annotations-folder on the Table desktop.

### Load Annotation File

Load in custom annotations (.csv file with character set “UTF-8”, separated by tab, and set to “quoted field as text”) created using Invivo6 software (or other software).

### Customizing Annotations (for DICOM Data)

Annotations can be added to any DICOM scan. There are no default annotations for these scans, but custom annotations can be made by tapping on “New Annotations”. Placing these is this same as with Gross Anatomy.


**WARNING:** Open same data set that was used to create annotation file in order to preserve correct coordinate system.

## DCM/INV FILES

*(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).*

Selecting the Volume Visibility Control Icon will open the following dialog.


<p><b>Volume Renderings</b>  <i>(Gray Scale, X-Ray, Transparent Soft Tissue, Transparent Hard Tissue, Transparent Soft + Hard, etc.)</i></p>	<p>Collection of different volume rendering presets (filters.) Each can be adjusted using the <i>Brightness</i> and <i>Contrast</i> slider bars on the main Table user interface. Users can create their own custom volume rendering presets using the Invivo6 software. This setting can be exported as a volume configure file (.vcf). <b>Custom</b> loads in a .vcf file.</p>
	<p><b>WARNING:</b> <i>Ultra High-Quality Rendering (UHQ) requires an NVIDIA graphics card to function. If Table application detects that some DCM files are missing, corrupt, or otherwise determined inaccurate, an error message will appear about possible inaccurate reconstruction. The user may continue with volume reconstruction and should exercise caution when reviewing any data with possible inaccuracies.</i></p>

## DCM/INV FILES

*(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).*

Selecting the Volume Visibility Control Icon will open the following dialog.


<p><b>Volume Renderings</b>  <i>(Gray Scale, X-Ray, Transparent Soft Tissue, Transparent Hard Tissue, Transparent Soft + Hard, etc.)</i></p>	<p>Collection of different volume rendering presets (filters.) Each can be adjusted using the <i>Brightness</i> and <i>Contrast</i> slider bars on the main Table user interface. Users can create their own custom volume rendering presets using the Invivo6 software. This setting can be exported as a volume configure file (.vcf). <b>Custom</b> loads in a .vcf file.</p>
	<p><b>WARNING:</b> <i>Ultra High-Quality Rendering (UHQ) requires an NVIDIA graphics card to function. If Table application detects that some DCM files are missing, corrupt, or otherwise determined inaccurate, an error message will appear about possible inaccurate reconstruction. The user may continue with volume reconstruction and should exercise caution when reviewing any data with possible inaccuracies.</i></p>

**CT Settings** (Only available in slice mode for DCM image sets)

Tap to show all available radiology presets. Using a second tap, select a particular radiology preset. The brightness and contrast can be adjusted using the slider bars to the right.

**Available CT presets:** Brain, Abdomen, Mediastinum, Bone, Lung, Liver

***WARNING:** CT presets are dependent upon original scanner/DCM HU definitions. User is responsible for adjusting rendering range and settings to ensure all structures are visible in the scan.*


***WARNING:** MRI scans will automatically have their volume rendering range adjusted for optimal viewing. The user can specify a specific rendering range by manually entering the minimum and maximum values in the settings menu.*

**MODEL SETTINGS**

Open an INV file with models created from a DICOM file in Invivo or MedicalDesignStudio. Open Visual Property Dialog for currently selected digital surface model.

*Visibility*

- Turn On/Off model visibility
- Adjust the *Opacity* of a particular model

*Rendering*

- Adjust *Mode*: Surface, Wireframe, or Points
- Adjust shading: Smooth or Flat
- Turn On/Off **Back Face Culling**
- Enable **No Clipping** (clipping planes do not affect model)

*Material*

- Adjust surface model appearance by changing color and light settings

*Texture*

- Turn On/Off texture applied to model
- Turn On/Off shading applied to model

Tapping **OK** will close the *Visual Property* window and save the changed settings.

This will show a list of all folders made.

### Quick Add

Tap this icon to automatically add the current view as the next preset.

### Quick-View Presets


Tap this icon to access the quick-view preset menu.

Select the Pre-Made Presets to view pre-made views of the cadaver based on common anatomy topics.

Select this icon to load whatever presets you have in your Preset Menu into slider-form.

Advance the Slider Bar or click the Previous and Next icons to move through the presets.

Tap Play to initiate a playback of all of the presets. Tap Loop to have the playback repeat automatically.

Select the Movie Capture icon to export a video of the presets.


**WARNING:** Preset behavior is best when using the same particular image data set and volume rendering range that was used when initially creating the presets. Current presets in the menu will be overridden by newly imported presets.

## FUNCTIONAL ANATOMY

*(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).*

Tap **Functional Anatomy** on Application Toolbar to open the *Functional Anatomy* Dialog.

- Dialogue lists all available Functional Anatomy scans with type and description.
  - **Birth Simulation:** Shows a birth simulation of a pregnant cadaver.
  - **Pregnancy:** Shows a pregnant cadaver with fetus.
  - **Cardiology:** Shows a beating heart with ECG.
  - **Kinesiology:** Shows various joint movements.
  - **Neural Pathways:** Shows what nerves innervate specific dermatomes, muscles, and organs.
  - **Catheterization:** Steps through the placement of a few heart cauterization procedures.

25.1  
25.2  
25.3

10  
25.4  
25.6  
25.7  
25.5

- **Dental:** Shows dental arch and pano slice view.
- **Ocular Applications:** Shows vision and ocular movements.
- **Cardiology (Full Body):** Shows a beating heart with ECG in full body cadaver.
- **Pathways:** Outlines a particular physiology pathway. (Available in both male and female).
- Tap entry name, then tap **Open** or double tap entry name.


25.8

25.7

**Birth Simulation Stages** (*in Birth Simulation only*)

Select any node on the birth stage map to jump to that point in the birth simulation. You can select the main stage icons, or any node in between.


**Volume Visibility** (*in Pregnancy only*)

Tap on the pink icon to open Volume Visibility dialog for the pregnant cadaver. Tap on the blue icon to open Volume Visibility dialog for the fetus.

**Joint Movements** (*in Kinesiology only*)

Tap this icon to open joint movements dialog. Tap on a joint icon or tap directly on a joint and use arrows to view movements.

**Shoulder**

From top left going clockwise:  
Abduction/Adduction, Flexion/Extension,  
Internal Rotation/External Rotation,  
Protraction/Retraction, Elevation/Depression

**Femur on Pelvis**

From left to right: Abduction/Adduction,  
Internal Rotation/External Rotation,  
Flexion/Extension

**Pelvis on Femur**

From left to right: Abduction/Adduction, Internal Rotation/External Rotation, Flexion/Extension


**Knee**

From left to right: Flexion/Extension, Internal Rotation/External Rotation


**ECG (in Cardiology or Catheterization only)**

Tap this icon in either of the above scans to open the ECG interface to adjust the controls and view the beating heart.

**ECG (in Pregnancy only)**

Tap the pink icon to open ECG interface for the pregnant cadaver. Tap the blue icon to open the ECG interface for the fetus.

**Nerve Flow Visualization (in Neural Pathways only)**

Tap this icon to open the nerve flow visualization tool. Tap on a dermatome, muscle, or organ to reveal the nerve pathway. Tap on the highlighted region in the legend to see more details.

**Catheter (in Catheterization only)**


Tap to choose a procedure and path. Go step by step through the procedure.

## CASE LIBRARY

*(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).*

The Case Library is a collection of scans showing various clinical case examples. Through the use of contrast media and other imaging techniques, real patient anatomical features are highlighted. Each scan comes with scan information, if available.

Select **Case Library** from Application toolbar to open the *Case Library* dialog.

- Scans are sorted by the following icons: Head and Neck, Thorax, Abdomen and Pelvis, Upper Limb, Lower Limb, Full Body, Animal, Archaeology, Embryo, and 4D.
- Tap on a divisional icon to view all associated scans.
- Select scan name to display additional clinical information (*Info, History, Findings, Impression, Other*) along with scan resolution and size.

- Tap the MRI or CT icons to see scans organized by tag.
- To search the Case Library, type text into the bottom left search bar.

- There are also two folder icons for Recently Viewed and Case Folders.
- Case Folders include Recently Added and Favorites.


- To create your own folder, tap the  button and enter in the desired Folder name.

- To remove your own folder, tap the  button.
- To add a case to a folder, tap the  button and choose the desired folder.
- To remove a case to a folder, tap the  button.
- To export a custom folder (.txt file), open the custom folder and tap the  button. A dialog will appear allowing the user to choose the name and location of the .txt file.
- To import a custom folder (.txt file), tap the  button under the folder's category and locate the .txt file for the desired custom folder.


**Comparison**

Displays three related scans linked in movement, clipping plane, brightness/contrast, and volume rendering mode for comparison. Note: Comparison Cases are not compatible with Vertical Viewing Mode.

**4D scans**

Scans showing movement

**To Play 4D scan:**

**Step 1:** Tap the View Sequencer icon .

**Step 2:** Tap  to start/pause movement. Tap  to loop playback.

**Step 3:** Tap  to move forward one frame. Tap  to move backward one frame.

**Step 4:** Tap  to cut to the last frame. Tap  to return to the first frame in the sequence.

**Step 5:** Tap the red X in the upper right to close dialog.


*User can rotate or clip volume while video loop is playing. Video playback will pause when user is rotating volume. Pause playback when not viewing scan to conserve computer memory.*

## HISTOLOGY LIBRARY

*(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).*

The Histology Library is a collection of slides showing various physiological and pathological examples with the use of various microscopic and staining techniques. Each slide comes with information, if available.

Select **Histology** from Application toolbar to open the *Histology Library* dialog.

- Slides are sorted by tissue type (ET, CT, MT, NT) or by organ system (*Cardiovascular, Connective Tissue, Reproductive, etc.*).
- Tap on a divisional icon to view all associated slides.
-  Tap the Interactive Icon to view a selection of annotated slides for each tissue type or system.
-  Tap the Pathology Icon to view a selection of pathologic slides for each tissue type or system.



- Tap on this icon then double tap on up to 4 slides to compare them side by side.
- Select scan name to display additional information, including a larger preview image, and microscopic technique/stain when available.
- To search the Histology Library, type text into the bottom left search bar.

## CURRICULUM

The Curriculum consists of a set of single scans with pre-made annotations and view sequences, provided by Anatomage. Scans are sorted by region (*Thorax, Upper Limb*, etc) and system (*Skeletal, Muscular*, etc). The curriculum is intended for teaching purposes.

*(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).*

## PROSECTION

Prosections consists of regional prosected 3D real cadaver images. These images consist of the external data so they can be rotated to view at different angles. No cuts can be made through these images. Some prosection images have pathology (*diseased liver, partial knee replacement, etc.*) and some are “healthy” prosections (*brain, lower limb, etc.*)

***(This section contains images featuring educational content and should not be used as clinical reference or for diagnostic purposes).***

## OPENING DCM/INV SCANS

**Step 1:** Select **Open File** from the Application toolbar.

**Step 2:** Tap on browse and use the Windows Explorer interface to navigate to directory of INV file or DCM file series.

- |            |  |
|------------|--|
| INV file   | Select file and press <b>Open</b> .  |
| DCM series | Select a single DCM file and press <b>Open</b> . Software will scan through folder and check each DCM file's metadata prior to loading all DCM files in the same series. |

**Step 3:** Table application will automatically construct image volume based on INV or DCM file.

For INV files, any additional content created and saved with the patient data using **Invivo6** software (surface models, models, etc.) will be loaded as well.

**WARNING:** Table application and Invivo6 software will load in DCM files contained within the same folder and of the same imaging series when reconstructing the volume. It is the responsibility of the user to confirm that all slice information is available and in the same folder when loading onto Table6.0 or saving from Invivo6.

**WARNING:** If Table application detects that some DCM files are missing, corrupt, or otherwise determined inaccurate, an error message will appear about possible inaccurate reconstruction. The user may continue with volume reconstruction and should exercise caution when reviewing any data with possible inaccuracies.

**WARNING:** When loading additional content created (surface models, comments, etc.), content is created by another user and is not part of the original patient image data.

## PACS INTEGRATION

The Anatomage Table can import scans directly from PACS server. Tap on “Open File”, tap on “Import PACS” and then tap on “Configure”. Under Server Setup, type in the PACS server AE Title, Host Name and Port. Under Client Setup, type in the AE Title and Port.

Tap Test Connection to confirm that the connection is successful.

Please contact Anatomage with any questions regarding this process.

The screenshot shows a dark-themed dialog box titled "PACS Import Configuration" with a close button (X) in the top right corner. It is divided into three sections: "Server Setup", "Client Setup", and "TLS Setup".  
- **Server Setup:** Contains three input fields for "AE Title:", "Host Name:", and "Port:". Below these are two buttons: "Configure Tags" and "Test Connection".  
- **Client Setup:** Contains two input fields for "AE Title:" and "Port:". Below them is a "Local Storage Directory:" label with a "Browse" button and a text field containing "C:\PacsDL". At the bottom of this section is a checked checkbox labeled "Delete Files After Open".  
- **TLS Setup:** Starts with an unchecked checkbox labeled "Enable". Below it are three input fields for "Peer Certificate:", "Client Certificate:", and "Client Private Key:", each with a "Browse" button to its right. At the bottom of this section is a "Cipher List:" dropdown menu currently showing "DES-CBC3-SHA".  
At the very bottom of the dialog box is a large "OK" button.

**WARNING:** Anatomage is not responsible for verifying the accuracy of patient information from a PACS system. Anatomage also assumes no responsibility in relation to privacy requirements for viewing and sharing patient scan information.

# **PROSPECT MASA**

# Anatamage TABLE CONVERTIBLE

The Anatamage Table Convertible is the most technologically advanced 3D anatomy visualisation and virtual dissection platform for comprehensive medical education, clinical training, diagnostics, and patient care, catering to a diverse range of learners at all levels of education.


## Content

ANATOMAGE BODIES	2 male, 2 female real life-size cadavers for gross anatomy
REGIONAL ANATOMY	40 regional scans up to 50µm high resolution
SEGMENTATION	Highly detailed segmentation of 2,950+ structures in male and female cadavers
FUNCTIONAL ANATOMY	Pregnancy, Neural Pathways, Kinesiology, Ocular Application, Cardiology, Dental
SIMULATIONS	Blood Flow, ECG, Heart Motion, Catheterisation, Arthroscopy, Endoscopic views
HISTOLOGY	1000+ high definition histology slides (280+ segmented) featuring tissue throughout the human body
PROSECTION	66 3D surface scans of segmented prosections - some with multiple dissection steps
ANIMAL ANATOMY	Real Dog, Cat, Frog, Mouse cadavers. 285 Animal CT/MRI scans
CLINICAL CASE LIBRARY	1,500+ real patient CT/MRI scans of normal anatomy, pathology, and unique diseases
CURRICULUM	500+ presets of anatomical structures available for supplementing curriculum

## Adjustable

Table Convertible can be tilted into a vertical position and is ideal for classroom demonstrations.

# Tools & Features


DISSECTION	3D cutting and surgical simulations. Layer-by-layer dissection. Individual structures can be cut, removed, interacted with, and isolated
LEARNING TOOLS	Visibility, Flat color, Screenshot, Annotation, Exploration, Pin Drop, and Pen Tool Fully customizable structure display
LECTURE SUPPLEMENT	Lecture tools available for recording Anatomage Table's activities for students to review after class
LEARNING ASSISTANT	Instantly look up any anatomy structure's associated system, category, origin, insertion, blood supply, innervation, actions, and functions
ASSESSMENTS	Create anatomy quizzes through traditional lab practicals, group assessments or self study through flashcards, multiple choice tests, and anatomy structure tests
RESOURCE SHARING	Upload or download lesson plans, quizzes, learning materials, preset or quizzes
ECOSYSTEM	Through Anatomage Share, our suite of products becomes a unified ecosystem, enabling users to easily upload and share materials

# Clinical Applications

3D VOLUME RENDERING	Render CT/MRI scans in 3D and visualise scans with Ultra High Quality volume rendering tools
VIRTUAL DISSECTION	Virtual scalpel, Craniotomy, Measurement, Screenshot, and Pin Drop Tool
RADIOLOGY WORKSTATION	Reference 2D slice data from CT/MRI and compare it to 3D rendering side-by-side. Capable of connecting to PACS to retrieve data. Open your own DICOMs

# Hardware

TABLE SIZE	2x46" life-sized screen (84" inch surface)
DISPLAY	High-definition monitors with intuitive multi-touch screen interface and tempered glass overlay
INTEGRATION	HDMI, USB, Ethernet ports. Connect to the internet, projectors, and external monitors
DIMENSIONS/WEIGHT	Length: 215/137 cm (84.8"/54"), Height: 90/215 cm (35.6"/84.5"), Width: 85 cm (33.5") Weight: 150kg (330 lbs)

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Ref: Purchase procedure of Anatomy Table (Virtual) with included software

To whom it may concern,

Analyzing the documentation issued for the referenced tender, regarding the provisions raised, based on the technical and other available informative materials, we would like to point out some of the main features in the Anatmage Table use:

- The device software comes by default in English + Latin terms, however it can be setup to be used in different languages, among which: French, Italian, Spanish or German. 
- The Anatmage Table 10 application can be open in several windows contemporary, even 2 windows working together (one on the right screen and one on the left). Due to possible hardware saturation we advise to keep open no more than 3 windows at a time. 
- The Anatmage Table allows for the use free of charge of Anatmage Share. This is a Cloud based system that allows users to upload/download and share their table presets, DICOMS, presentations, video tutorials, ecc. Some of these contents can be downloaded and used in the table itself. The service is free of charge and it includes unlimited space for content uploaded as visible for other users, and up to 1GB for private content. 
- There are more than 1600 DICOM dataset cases on cases obtained from High Definition DICOM datasets, from which more than 100 have 3D reconstructions (segmentation). Info, history, observations, case assumptions information is also available. There are also more than 20 embryo develop series scans. 
- The case library includes also fully segmented dog and cat, and partially segmented mouse and frog. Besides more than 300 animal CT scans (horse, pig, sheep, owl, alligator, turtle, chimp, ecc). 
- It is also available more than 1180 high definition histology slides (500+ segmented, 250+ pathological). On prosection there are available more than 70, 3D scans (actually 75). 


# Anatamage

- The support materials made available by Anatamage are including User Manuals, tutorials on how to use each feature of the table, as well as video tutorials over different functionalities. Inside Anatamage Share website, the user can take advantage of the Training Portal (for free), and enjoy interactive training online with 3 different levels: basic, intermediate and advanced.

Please feel free to contact us if you have any questions regarding this matter.

Sincerely,

Carlos Martinez  
Channel Sales Manager  
Anatamage Italy srl


# Anatmage

## Anatmage Table Unique Features

The Anatmage Table, virtual dissection Table for anatomy and physiology education, is exclusively produced by Anatmage, Inc.

To our knowledge, the Anatmage Table is the only all-in-one device in the world able to show the anatomy of a real human body from frozen cadaver data in 1:1 life-size with the following:

- Data derived from 8 different frozen human specimens
- Including 5 full frozen body specimens, 3 males and 2 females
- Up to 50  $\mu\text{m}$  resolution on frozen bodies
- Segmented histology slices
- Macro and microanatomy correlation at runtime
- 1600+ DICOM datasets
- UHQ (cinematic) renderings (DICOM data)
- Segmented and annotated DICOM datasets
- Step-by-step segmented prosections - surface imaging data
- Kinesiology on frozen bodies
- Eye motion and optometry on frozen bodies
- Nerve pathways on frozen bodies
- Catetherisation simulation on frozen bodies
- Heart beating and blood flow on frozen bodies
- Heart pathologies simulation on frozen bodies
- Pregnancy and birth simulation on frozen bodies
- Ultrasound visualization simulation on frozen bodies

Contents and features are continuously updated to ensure they remain both current and accurate.


More than 4000 tables have been sold worldwide and there is a thriving global community of educators who have already spent time developing their own contents and ideas on how to best incorporate the Table into a wide range of curricula and disciplines.

Do not hesitate to contact us for further information about the declaration reported above.

Kind regards,

  
Davide Giacomo Tommasi  
Regional Manager - Europe


# **EXTRASE MENIURI**


**Fig 1**

**Fig 2**

**Fig 3**

**Fig 4**

**Fig 5**



Fig 6

9

Fig 7

15

Fig 8

15



Fig 9

20

Fig 10

20

Fig 11

Fig 12


Fig 13


Fig 14


Fig 15

29

Fig 16




Fig 17

31

Fig 18


31

30

30

Fig 19


30

31


32

Fig 20

37

Fig 21

35

35

Fig 22



Fig 23

**MATERIALE INSTUIRE “HOW TO...”**

# *Anatomage*

# TABLE

## **Cadavers**

How to Use Quiz Mode

### **Table of Contents**

<a href="#">Create a Quiz</a>	1
<a href="#">Edit a Quiz</a>	7
<a href="#">Take a Quiz</a>	8

## Create a Quiz

**Step 1:** Select the Quiz icon  then select the Create Quiz icon .

**Step 2:** Enter a title for the quiz.

**Step 3:** Select Settings to customize quiz settings. Check any of the options below to enable them. Uncheck them to disable.


Table 10 - How to Use Quiz Mode

Some of these options include:

- Repeat skipped questions: If a student skips a question during the quiz, they can retake this question at the end.
- Same vessel type only: If you are quizzing on the nervous or cardiovascular system, allow for only the vessel type in question to be selected. For example, if the question is the femoral artery, students will only be allowed to click on arteries. Veins and nerves will not register.
- Randomize: The questions are randomly chosen each time the quiz is taken.
- Protect with PIN: Password protect any quiz file with a PIN. This PIN will be required to load in the file, edit the file, and start a quiz.
- Highlight color: Adjust what color structures on the cadaver will be highlighted during quiz questions.
- Student: Use the "+" icon to import in a student list in either.csv or .txt. When the students are taking the quiz, they will be able to choose their name from a drop down menu. This is optional, students will also be able to type their name in when taking a quiz, or leave this blank.
- Timers: set a timer per question or a total quiz time, and choose whether or not this time is displayed for the student.
- Display Results: At the end of the quiz, display the # correct, incorrect, and percentage for the quiz.
- Scoring: In the case of a game or an anatomy competition, for example, assign a point value to questions, and choose whether or not to display the score. Skipped questions can be positive and negative whole number values.
- # of Questions from Pool: Assign a certain number of questions from a total set. For example, your quiz file has a total of 50 questions, but you want each student to take only 20. Set this value to 20. Leave blank if you want all 50.
  - Note: The 20 questions will be randomly selected from the 50.

Once the settings are selected, tap "Apply" to return to the previous menu.

**Step 4:** Select the quiz type that you would like to create.

Pick Me will ask students to locate the indicated structure on the cadaver, and will reveal right and wrong results at the end.

Flashcards will ask students to locate the indicated structure on the cadaver, and will reveal the right and wrong results immediately after the answer is confirmed.


Tap the ellipse icon to expand and view other quiz options.

Highlight will light up the indicated structure on the cadaver. Students can test

themselves, then tap to reveal what the right answer is.

Multiple Choice will light up the indicated structure on the cadaver and prompt students with an automatically generated multiple choice list.


Game will convert the quiz into a game with 2, 4, 6, or 8 players. *(Tip: set a timer in the quiz settings to make games more high stakes!)*

Tournament Mode If you have ever participated in an Anatomage Tournament, this is the mode you can use to mimic those competitions. If you are interested in an Anatomage Tournament or to host one yourself, please contact us!

**Step 5:** Select the arrow icon  to continue to the next page.

**Step 6:** Read the instructions on the dialog. Then select the arrow icon  to continue to the next page.

**Step 7:** Use the Table as normal to set up your quiz view and questions. For example, use the

Visibility menu  to set up the desired view and add annotations to define questions. (See "How to Use the Visibility Menu" document for more information.)




Table 10 - How to Use Quiz Mode

**Step 8:** Select the plus icon  to add the selected questions to the quiz.

**Step 9:** If you need to delete a question or reorder the questions, select the number key next to each image to select the question(s). Select the check icon  to select all questions. Select the minus icon  to delete the question(s). Engage the triple line icon  to adjust the order of questions.

**Step 10:** Select the Save icon  to expand all saving options. Tap the Quiz History icon  to save to quiz history, the Export icon  to save to the computer, or the Folder icon  to save to a folder within the Table software.

**\*Please note, you can also add bony landmarks and origin/insertions to a quiz by turning them on and annotating as you normally would.**


## Creating Multiple Quizzes Back-to-Back in the same Session

**Step 1:** Follow the steps in **Create a Quiz** to create and save your first quiz.

**Step 2:** Select the Quiz icon  then select the Create Quiz icon .


This will automatically populate the quiz you just made in case you would like to edit that quiz.

**Step 3:** To move forward with a new quiz, delete the title of the old quiz and input the new quiz name. This will not delete the old quiz, as you saved that in Step 1.

## Edit a Quiz

**Step 1:** Select the Quiz icon  then select the Create Quiz icon .


**Step 2:** Select the load icon  to expand all loading options. Tap the Quiz History icon  to load from quiz history, the Import icon  to import from the computer, or the Folder icon  to load from a folder within the Table software.

**Step 3:** Follow Steps 2-10 from "Create a Quiz" to edit and re-save the quiz.

## Take a Quiz

**Step 1:** Select the Quiz icon  then select the Take a Quiz icon .

**Step 2:** Select the load icon  to expand all loading options. Tap the Quiz History icon  to load from quiz history, the Import icon  to import from the computer, or the Folder icon  to load from a folder within the Table software.

Alternatively, take a default quiz based on all structures currently in view on the rendering window without loading a previously saved quiz.

**Step 3:** Select Settings to customize quiz settings. Check any of the options below to enable them. Uncheck them to disable. Then tap "Apply".


Table 10 - How to Use Quiz Mode

**Step 4:** Select the quiz type that you would like to take. Select the Pick Me icon  or

Flashcards icon . Tap the ellipse icon  to expand options and select the Highlight

icon , Multiple Choice icon , Game icon  or Tournament Mode 

**\*Refer to page 2 and 3 for a description of each setting and quiz type.**

**Step 5:** Enter a student name or select the plus icon to load a student list text file.

**Step 6:** Tap the Start icon  to start the quiz.

# *Anatomage*

# TABLE

## Case Library

### How to Use the Visibility Menu


#### **Table of Contents**

How to Use the Visibility Menu	2
Description of each Volume Rendering for CT scans	3
Description of each Volume Rendering for MRI scans	4
How to Use UHQ	5

## How to Use the Visibility Menu

Step 1: Open a scan in the Case Library. Scans will load in with the Opaque Hard Tissue

filter applied. Tap on the Volume Visibility menu  to open a list of all available volume renderings.

Step 2: Tap on any of the volume renderings in the list to change the filter on the scan.


Table 10 - Case Library - How to Use the Visibility Menu

20

## Description of each Volume Rendering for CT scans













<p><b><u>Opaque Hard Tissue:</u></b></p> <ul style="list-style-type: none"> <li>- Default view when a case is first loaded</li> <li>- Good for looking at bones</li> <li>- Bring in more soft tissues (eg. muscles) by increasing the brightness</li> </ul> 	<p><b><u>Transparent Hard Tissue:</u></b></p> <ul style="list-style-type: none"> <li>- Good for looking at teeth and other hard structures, such as bone and medical devices</li> <li>- Bring in more soft tissues (eg. muscles) by increasing the brightness</li> </ul> 	<p><b><u>X Ray:</u></b></p> <ul style="list-style-type: none"> <li>- Good for viewing radiopaque structures such as bone and medical devices</li> <li>- Decrease brightness if the image is washed out</li> </ul> 
<p><b><u>Transparent Soft Tissue:</u></b></p> <ul style="list-style-type: none"> <li>- Good for looking at tissues that contain air (eg. respiratory sinuses, lungs, and GI tract)</li> </ul> 	<p><b><u>Transparent Soft and Hard:</u></b></p> <ul style="list-style-type: none"> <li>- A combination of the transparent soft tissue and transparent hard tissue renderings</li> <li>- Good for looking at the skeleton while having the skin (and other soft tissues) as a reference</li> </ul> 	<p><b><u>Opaque Soft and Hard:</u></b></p> <ul style="list-style-type: none"> <li>- A combination of the opaque soft tissue and opaque hard tissue renderings</li> <li>- Good for looking at the skeleton while having the skin (and other soft tissues) as a reference</li> </ul> 
<p><b><u>Opaque Soft Tissue:</u></b></p> <ul style="list-style-type: none"> <li>- Good for viewing soft tissues</li> <li>- Soft tissues (eg. muscle) will appear opaque</li> <li>- Suggestion: Apply a cut or a clipping plane to view the internal organs</li> </ul> 	<p><b><u>GrayScale:</u></b></p> <ul style="list-style-type: none"> <li>- Good for viewing bones</li> <li>- In some scans, it's also good for viewing blood vessels and the bronchial tree</li> </ul> 	

Table 10 - Case Library - How to Use the Visibility Menu

## Description of each Volume Rendering for MRI scans


<p><b><u>Transparent Hard Tissue:</u></b> -Default view when a case is first loaded -Good for looking soft tissue</p> 	<p><b><u>Soft Tissue A - MRI:</u></b> - Good for looking at soft tissue including soft tissue masses, ligaments and tendons</p> 	<p><b><u>Transparent A - MRI:</u></b> - Good for looking at soft tissue - Intensifies grayscale</p> 
<p><b><u>Soft Tissue B - MRI:</u></b> -Good for looking at soft tissue -Intensifies grayscale</p> 	<p><b><u>Transparent B - MRI</u></b> - Good for looking at soft tissue - Intensifies grayscale by using different colors</p> 	<p><b><u>GrayScale:</u></b> - Good for looking at soft tissue</p> 

2)

## How to Use UHQ

Step 1: UHQ renderings stand for ultra high quality. Tap on any of the filters in this list to change to a UHQ rendering.



## Citrine:

- Useful for muscle origin/insertion visualization


- Helpful to visualize blood vessels


### Quartz/Opal (MRIs)

- Helpful to visualize tumors in MRIs

## Aquamarine/Amethyst

- Helpful to visualize skin or lungs contours

**NOTE: Please refer to Case Library: How to Use the Slider Bars to learn how to adjust brightness and contrast for renderings.**

**NOTE: Please refer to the Volume Rendering Guide for more detail on changing the volume renderings.**

# **DIPLOMA, ISO9001, DECLARATIE CE**

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

**Anatamage, Inc.**  
3350 Scott Blvd #29  
Santa Clara  
California  
95054  
USA

Holds Certificate Number: **FM 701540**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

The design, manufacture and support of software used in the archiving, processing, communication of medical images, assisting in diagnosis, and the planning and simulation of orthodontics, surgical elective treatment and dental implants; and the manufacture of a device for the placement of dental implants.

For and on behalf of BSI:

\_\_\_\_\_  
Graeme Tunbridge, Senior Vice President Medical Devices

Original Registration Date: 2011-09-19

Latest Revision Date: 2022-09-15

Effective Date: 2022-06-26

Expiry Date: 2025-06-25

Page: 1 of 1

...making excellence a habit.™

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.  
An electronic certificate can be authenticated [online](#).  
Printed copies can be validated at [www.bsigroup.com/ClientDirectory](http://www.bsigroup.com/ClientDirectory)

Issuing Body: BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
BSI Group The Netherlands B.V. is registered in The Netherlands under number 33264284 | A Member of the BSI Group Holdings B.V.  
Contact Office: 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA.

## Certificate of Attendance

This is to certify that on March 28-29, 2017

**Mr. Petre Melesteu**

from company Neotech srl

Successfully attended the course  
"Anatamage Table Training"  
in Milan, Italy.

Davide Giacomo Tommasi  
Regional Manager - Europe  
Anatamage Inc.

# *Anatomage*<sup>®</sup>

Anatomage  
Diploma for



**Petre Melesteu**  
**Neotech**

Who has attended and received a refresh training at  
the EMEA Distributor Meeting that took place on May  
2024 in Milan.

Anatomage, Inc.  
Regional Manager - Europe

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## EU DECLARATION OF CONFORMITY

We, **ANATOMAGE INC**  
3350 Thomas Road #150  
SANTA CLARA CA 95054 USA

declare under our sole responsibility that our Anatamage Tables (Models MTC, MTD, ABA) to which this declaration relates is in conformity with the **Low Voltage Directive 2014/35/EU** of 26 February 2014 and the **Electromagnetic Compatibility Directive 2014/30/EU** of 15 December 2004

per the provisions of the following standards.

### Safety Standards

- IEC 62368-1:2014Ed.2+C1, IEC 62368-1:2018 Ed.3, IEC 60950-1:2005Ed.2+C2 – Audio/Video, Information And Communication Technology – Equipment - Part 1: Safety Requirements

### Electromagnetic Compatibility Standards

- EN 55032:2015+A11:2020 – Electromagnetic Compatibility of Multimedia Equipment – Emission Requirements
- EN IEC 61000-3-2:2019+A1 – Electromagnetic Compatibility (EMC)- Part 3-2: Limits - Limits for Harmonic Current Emissions (Equipment Input Current ≤ 16 A per Phase)
- CENELEC EN 55035:2017+A11 – Electromagnetic compatibility of multimedia equipment – Immunity requirements
- FCC 47CFR: (Part 15 Subpart B) – Title 47 CFR Part 15 Subpart B: Unintentional Radiators
- IC ICES-003:2020 Ed.7 – Spectrum Management and Telecommunications - Interference-Causing Equipment Standard - Information Technology Equipment

In addition to the above directives, the products covered under the declaration of conformity also complies with the **Restriction of the use of certain hazardous substances in electrical and electronic equipment Directive 2011/65/EU** of 8 June 2011.

per the provisions of the following standard.

- EN IEC 63000:2018 - Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

*Kris Thomson, Director of Table Applications*  
Nov 13, 2023