


# ACUSON Sequoia Ultrasound System

4D Heart<sup>AI</sup>  
3.5 (VB30)





- **Overview of 4D Heart<sup>AI</sup> application**
- Review 4D Heart<sup>AI</sup> configuration and touch screen layout
- Illustrate imaging for 4D Heart<sup>AI</sup>
- Review 4D Heart<sup>AI</sup> workflow
- Evaluate 4D Heart<sup>AI</sup> analysis and report



# 4D Heart<sup>AI</sup> – clinical relevance and benefit

- Clinical guidelines recommend 3D evaluation for cardiac size and function<sup>1</sup>
- 3D strain analysis is faster with good reproducibility compared to 2D<sup>2</sup>
- 3D strain is more accurate with less variability due to entire LV cavity detection<sup>3</sup>
- AI-enabled analysis:
  - Chamber contouring for size and function
  - 4D strain analysis of left heart
- Assessment:
  - All visible chambers in TTE
  - Visible ventricles in TEE
- Available without ECG



Sources | 1. Mitchell, Carol, et al., *Journal of the American Society of Echocardiography* (2019) 32;1  
2. Hasanfatta, Naeem et al., *ACTA Scientific Clinical Case Reports* (2022) 3;6: 45-56  
3. Kosarju, Ateet, et al., *Left Ventricular Ejection Fraction*. StatPearls Publishing (2023) Jan

# Objectives

- Overview of 4D Heart<sup>AI</sup> application
- **Review 4D Heart<sup>AI</sup> configuration and touch screen layout**
- Illustrate imaging for 4D Heart<sup>AI</sup>
- Review 4D Heart<sup>AI</sup> workflow
- Evaluate 4D Heart<sup>AI</sup> analysis and report



# 4D Heart<sup>AI</sup> configuration

**System Configuration**

- ▶ System Settings
- ▶ Workflow Enhancement
- ▶ Imaging Settings
- ▶ Measurement & Report
- ▶ Annotations
- ▼ Features
  - Protocols
  - Stress Echo
  - 4D Heart<sup>AI</sup>**
  - 2D Heart<sup>AI</sup>
  - AI Assist
- ▶ Transducer
- ▶ Connectivity & Network
- ▶ Peripheral Devices
- ▶ System Management



Contours

- Left Ventricle
- Left Atrium
- Right Ventricle
- Right Atrium

Measurement Results

Chambers

- Left Ventricle (max 6)
- Left Atrium (max 4)
- Right Ventricle (max 4)
- Right Atrium (max 4)

Measurements

- EF (%)
- GLS (%)
- Mass (g)
- CO (L/min)
- EDV (mL)
- ESV (mL)
- SV (mL)

LV Epicardium

- Show
- Hide

Default Direction for Speckle Tracking

- Longitudinal
- Circumferential

Default Curve After Speckle Tracking

- Volume
- Segmental Strain
- Global Strain

Default Timing/Phase After Speckle Tracking

- End Systolic Strain
- Peak Systolic Strain
- Time to Peak
- Post Systolic Strain Indexed

Note: Capture sends all available measurements to the report.

Display measurements as BSA indexed (BSA suffix in above table)

Number of Beats

- 1
- 2
- 3
- 4
- 5

Results Format

- Individual Beat
- Average of Beats

Restore Defaults

# System configuration – left ventricle (LV) chamber options

Measurement Results

Chambers	Measurements
<input checked="" type="radio"/> Left Ventricle (max 6)	<input checked="" type="checkbox"/> EF (%)
<input type="radio"/> Left Atrium (max 4)	<input checked="" type="checkbox"/> GLS (%)
<input type="radio"/> Right Ventricle (max 4)	<input checked="" type="checkbox"/> Mass (g)
<input type="radio"/> Right Atrium (max 4)	<input checked="" type="checkbox"/> CO (L/min)
	<input checked="" type="checkbox"/> EDV (mL)
	<input checked="" type="checkbox"/> ESV (mL)
	<input type="checkbox"/> SV (mL)

Note: Capture sends all available measurements to the report.

Display measurements as BSA indexed (BSA suffix in above table)

Measurements available for the left ventricle include:

- Ejection Fraction (EF)
- Global Longitudinal Strain (GLS)
- Mass
- Cardiac Output (CO)
- End Diastolic Volume (EDV)
- End Systolic Volume (ESV)
- Stroke Volume (SV)
  
- Maximum display of 6 measurements for LV

**Note:** All measurements are populated to the report regardless of the options selected for on-screen display in system configuration.

# System configuration – left atrium (LA) chamber options

Measurement Results

Chambers	Measurements
<input type="radio"/> Left Ventricle (max 6)	
<input checked="" type="radio"/> Left Atrium (max 4)	<input checked="" type="checkbox"/> EF (%)
<input type="radio"/> Right Ventricle (max 4)	<input checked="" type="checkbox"/> GLS (%)
<input type="radio"/> Right Atrium (max 4)	<input checked="" type="checkbox"/> EDV (mL)
	<input checked="" type="checkbox"/> ESV (mL)

Note: Capture sends all available measurements to the report.

Display measurements as BSA indexed (BSA suffix in above table)

Measurements available for the left atrium include:

- Ejection Fraction (EF)
- Global Longitudinal Strain (GLS)
- End Diastolic Volume (EDV)
- End Systolic Volume (ESV)

- Maximum display of 4 measurements for LA

# System configuration – right ventricle (RV) chamber options

Measurement Results

Chambers	Measurements
<input type="radio"/> Left Ventricle (max 6)	<input checked="" type="checkbox"/> EF (%)
<input type="radio"/> Left Atrium (max 4)	<input checked="" type="checkbox"/> EDV (ml)
<input checked="" type="radio"/> Right Ventricle (max 4)	<input checked="" type="checkbox"/> ESV (ml)
<input type="radio"/> Right Atrium (max 4)	<input checked="" type="checkbox"/> SV (ml)

Note: Capture sends all available measurements to the report.

Display measurements as BSA indexed (BSA suffix in above table)

Measurements available for the right ventricle include:

- Ejection Fraction (EF)
  - End Diastolic Volume (EDV)
  - End Systolic Volume (ESV)
  - Stroke Volume (SV)
- 
- Maximum display of 4 measurements for RV

# System configuration – right atrium (RA) chamber options

Measurement Results

Chambers	Measurements
<input type="radio"/> Left Ventricle (max 6)	<input checked="" type="checkbox"/> EF (%)
<input type="radio"/> Left Atrium (max 4)	<input checked="" type="checkbox"/> EDV (ml)
<input type="radio"/> Right Ventricle (max 4)	<input checked="" type="checkbox"/> ESV (ml)
<input checked="" type="radio"/> Right Atrium (max 4)	


Note: Capture sends all available measurements to the report.

Display measurements as BSA indexed (BSA suffix in above table)

Measurements available for the right atrium include:

- Ejection Fraction (EF)
  - End Diastolic Volume (EDV)
  - End Systolic Volume (ESV)
- 
- Maximum display of 4 measurements for RA

# 4D Heart<sup>AI</sup> touch screen layout




# Objectives

- Overview of 4D Heart<sup>AI</sup> application
- Review 4D Heart<sup>AI</sup> configuration and touch screen layout
- **Illustrate imaging for 4D Heart<sup>AI</sup>**
- Review 4D Heart<sup>AI</sup> workflow
- Evaluate 4D Heart<sup>AI</sup> analysis and report




# Volume acquisition

- Optimal 2D imaging is essential to a quality volume acquisition
- Prior to acquisition it's important to evaluate all reference planes for structural boundaries
- Alignment is important for assessment of specific anatomy



Optimize the 2D image

1. Press **4D**
  - Further optimize 4D image as needed
2. To acquire a volume, select **Clip** key



**Tip:** Ensure entire chambers are in field of view and have patient perform a breath hold to reduce motion of the heart.

## Acquisition with ECG (physio)

- Use a beat capture – must have at least one full R-R capture
- Ensure ECG tracing is optimal and clean
- Multiple heart cycles can be analyzed and averaged (up to **five**)

## Acquisition without ECG (non-physio)


- Use a time capture
- ED/ES key frames are set by user and proceeds to analysis
- Only **one** cardiac cycle can be analyzed
  - Exclude Beat, Beats Analyzed and Average Result are therefore not available
- Tracking is available for all visible chambers in TTE, and visible ventricles in TEE

# Acquisition for 4D strain analysis – volume rates

- 4D bullseye analysis on a dataset > 15vps (volumes per second) results in a cautionary display within the viewport

Low Volume Rate < 15vps

- To optimize vps
  - Decrease line density
  - Reduce depth
  - Reduce lateral and/or elevational field of view when needed
- Acquire steady beats, without any changes to field of view – have patient hold breath if needed





# Objectives

- Overview of 4D Heart<sup>AI</sup> application
- Review 4D Heart<sup>AI</sup> configuration and touch screen layout
- Illustrate imaging for 4D Heart<sup>AI</sup>
- **Review 4D Heart<sup>AI</sup> workflow**
- Evaluate 4D Heart<sup>AI</sup> analysis and report



# Launching 4D Heart<sup>AI</sup> from cine



- Ability to launch 4D Heart<sup>AI</sup> while imaging
  - During live 4D imaging, select **Freeze**
    - At least one R-R cycle must be in cine
  - Select **4D Heart<sup>AI</sup>** on touch screen



**Note:** Bullseye analysis is *ONLY* available in Review workflow.

# Launching 4D Heart<sup>AI</sup> from review



- From Review, select one 4D dataset
- Select **4D Heart<sup>AI</sup>** on touch screen
- Visible cardiac chambers are analyzed
  - All four chambers in transthoracic 4D imaging
  - Left and right ventricles in transesophageal 4D imaging



# Non-physio workflow example



## Ability to launch application without physio

- Use soft key to define key frames (ED & ES) and tap **Analyze** to proceed with analysis
- Only 1 cardiac cycle is analyzed with non-physio workflow
- Controls for Beat Analyzed, Exclude Beat and Average Result will not be available





# Analysis display (ED-ES View)

- Contours display on key frames ED & ES:
  - Apical 4
  - Apical 2 or PSAX or Coronal
- Viewports include
  - ED Viewport
  - ES Viewport
  - Volume render with contours
  - Results Table
- Freeze button is disabled; use Cine soft key to scroll between frames or beats



# Analysis display (Detailed View)

- Provides multiple views for one frame, to give better representation of the contours
- Layout can differ depending on chamber(s) selected for display
  - LV: Apical 3, Apical 2, Apical 4 and SAX
  - LA: Apical 3, Apical 2, Apical 4 and SAX
  - RV: Apical 4, Coronal, SAX
  - RA: Apical 4, RH A2, SAX
- Freeze button is available and users can play image or use Cine soft key to scroll between frames or beats




# Changing end diastole (ED) and end systole (ES) frames

ED and ES frames are depicted as purple markers on physio line as well as label on viewports

Move ED or ES frames with soft keys (per heart beat analyzed)


- Rotate ED soft key – press to set
- Rotate ES soft key – press to set




# Deleting contour

Ability to delete any chamber on ED or ES key frames

- If performed during ED-ES View, users are automatically taken to Detailed View and prompted to place manual points
- Deletion occurs for the active key frame
- Deleting a contour removes analysis from results table



RV example





- Available when a chamber contour is missing on ED and/or ES key frame (due to deletion or detection failure)
- Once launched, viewport will change to Detailed View
- On-screen guide at bottom left of monitor prompts user for point placement
- Workflow consist of semi manual contouring with 3 points for LV, LA, RA, and 6 points for RV
  - LV: Anterior Septal, Posterior, Apex
  - LA: Anterior Septal, Posterior, Roof
  - RA: Anterior Septal, Freewall, Roof
  - RV: Septal, Freewall, Apex in *Apical view* and Posterior, Infundibulum, Lateral in *SAX view*

# Adding contour – left ventricle (LV)

Views displayed in Detailed View:

- Apical 3
- Apical 2
- Apical 4
- SAX



Place AntSept, Post, Apex points in A3C view


# Adding contour – left atrium (LA)

Views displayed in Detailed View:

- Apical 3
- Apical 2
- Apical 4
- SAX

Manual contour points:

- Anterior Septal
- Posterior
- Roof



Place AntSept, Post, Roof points in A3C view


# Adding contour – right atrium (RA)

Views displayed in Detailed View:

- Apical 4
- RH Apical 2
- SAX

Manual contour points:

- Septal
- Freewall
- Roof



Place Sept, FW, Roof points in A4C view


# Adding contour – right ventricle (RV)

Views displayed in Detailed View:

- Apical 4
- Coronal
- SAX

Manual contour points:

- Septal (A4)
- Freewall (A4)
- Apex (A4)
- Posterior (SAX)
- Infundibulum (SAX)
- Lateral (SAX)



Pointer


Place Sept, FW, Apex points in A4C view and Post, Infund, Lat points in SAX view

# Moving MPR axis marker

- Available at 4D Heart<sup>AI</sup> launch
- Not available in Edit Contour mode
- In ED-ES View, one axis marker is seen for each view
- In Detailed View, axis markers are displayed for each

 Pan tool moves axis marker left/right, up/down

 Z rotation tool allows Z rotation of axis marker





Go To ES



Pointer



Select **Edit Contour** on touch screen

- Model appears as mesh in edit mode
- Perform edits on ED & ES frames only
- Works in both ED-ES View and Detail View
- Move between ED and ES frames for editing
  - Click in ED or ES viewport
  - Select the update key on control panel
- Selected contour region is highlighted
  - When pointer is closer to contour, a smaller region becomes active
  - When pointer is further from contour, a larger region becomes active






# Editing contours – tips and tricks

- Select **Edit Contour** on control panel
- To perform edits, enable display for right heart and left heart chambers separately
  - Axis markers auto align to chambers
- If making LV contour edits, disable LV Epi
  - Easier to delineate LV contour
- Enable **SAX View** by pressing soft key
  - Enables slicing up/down through chamber(s) when rotating soft key
- Perform edits in Detail View to see better representation of contours



# Beats analyzed

- Ability to change number of beats analyzed
- Supports between 1 to 5 beats
- If more than 1 beat is analyzed, results can be averaged
- Default number of beats analyzed at software launch is dependent on feature configuration
- Increasing number of beats analyzed after software launch results in algorithm re-running






Average (3 beats)		HR (bpm)	94
<b>LV</b>			
EF (%)	22.8	EDV (ml)	147.5
GLS (%)	-3.6	ESV (ml)	112.7
CO (l/min)	3.3	Mass (g)	138.4
<b>RV</b>			
EF (%)	28.5	ESV (ml)	62.9
EDV (ml)	88.5	SV (ml)	25.7

# Excluding beat

- Available when number of beats analyzed is more than 1
- Used to remove specific beat(s) results when multiple beats are analyzed
- Hides contours in model for selected beat for exclusion
- Excludes selected beat(s) from results
- On-screen message displays in bottom left hand of monitor

Current beat excluded. Contours hidden.




Go To ES



Pointer


# Bullseye View workflow

- Bullseye View option is only available when 4D Heart<sup>AI</sup> is launched from Review
- Left heart chambers LV & LA have to be selected for display in order to select Bullseye View
- 4D strain tracking is available for LV and LA chambers
- While in 4D Heart<sup>AI</sup> analysis, select **Bullseye View** on the touch screen



# 4D strain cardiac timing

- Flexibility to adjust cardiac timing
- Rotate soft keys to adjust
  - **AVC** (aortic valve closure)
  - **Pre-A** (just before mitral valve opening as a result of atrial contraction)






**Note:** AVC and Pre-A soft keys are available for timing adjustment when the display is frozen.

AVC	Pre-A
← Move → 326ms	← Move → 522ms

# 4D strain curve display options





- Left ventricle (LV)
  - Volume
  - Segmental Strain
  - Global Strain
- Left atrium (LA)
  - Volume
  - Global Strain



# 4D strain timing display options

4D strain timing options for the left ventricle (LV) are only available with Segmental Strain curve display



- Peak Systolic Strain
- End Systolic Strain
- Time to Peak
- Post Systolic Indexed



# 4D strain direction display options

4D strain direction options for the left ventricle (LV) are only available with Segmental Strain curve display

- Longitudinal
- Circumferential





# 4D strain: Bullseye View

## Exclude Beat

- Used to remove specific beat(s) 4D strain results when multiple beats are analyzed
- Available when number of beats analyzed is more than 1
- Hides contours in model for the excluded beat
- Excludes selected beat(s) from results
- On-screen message displays in bottom left hand of monitor

Current beat excluded. Contours hidden.




Current beat excluded. Contours hidden.

# Bullseye View interactions

Users can use the pointer to interact with bullseye segments and curves

- Highlighting a segment on bullseye will highlight corresponding segment on model
  - Segment(s) can be enabled or disabled by clicking in bullseye
  - If segmental strain is active, disabling a segment removes corresponding curve
- Hovering over a point on the curve displays the value at that location



# Objectives

- Overview of 4D Heart<sup>AI</sup> application
- Review 4D Heart<sup>AI</sup> configuration and touch screen layout
- Illustrate imaging for 4D Heart<sup>AI</sup>
- Review 4D Heart<sup>AI</sup> workflow
- **Evaluate 4D Heart<sup>AI</sup> analysis and report**



# Accessing 4D Heart<sup>AI</sup> report worksheet

To access the 4D Heart<sup>AI</sup> report worksheet:

- Select **Report** on the touch screen
- Select **4D Heart<sup>AI</sup>** on left-hand side of report

The screenshot displays the Siemens Healthineers Patient Report interface for the 4D Heart<sup>AI</sup> application. The left-hand navigation menu is visible, with '4D Heart<sup>AI</sup>' selected. The main report area shows the 'Last Single Beat Analysis' for 'LV End Systolic Strain [%] - Beat 1 of 1'. This includes a circular strain map with a color scale from -20 to 20, and a table of LV and LA parameters. Below the strain map is a 'Summary' table with columns for LV, Value, and five numbered columns (1-5), along with a 'Method' dropdown. At the bottom, there are buttons for 'Summary', 'Export to PDF', 'Store Report', 'Transfer Report', and 'Print Report'. A 'Collapse All' button is also present. The bottom right corner indicates 'Edited Value \*'.

LV	Value	1	2	3	4	5	Method
GLS	-14.7 %	-9.0	-11.3	-13.3	-19.9	-19.9	Avg
EF	51.1 %	46.8	46.3	49.1	56.7	56.7	
EDV	169.0 ml	158.6	157.5	167.1	181.0	181.0	

**Note:** Entering the report prompts user to exit the 4D Heart<sup>AI</sup> application.

# Edit Report | Summary

- Worksheet opens at default for edits
- Selecting **Summary** toggles from edit report to report summary

The screenshot shows the Siemens Healthineers Patient Report interface. The 'Summary' tab is selected and highlighted with an orange box. The interface displays a circular strain map and a table of cardiac parameters.

**Last Single Beat Analysis**  
LV End Systolic Strain [%] - Beat 1 of 1

	LV	LA
GLS	-19.9 %	GLS Res 35.0 %
EF	56.7 %	GLS Cond -15.5 %
EDV	181.0 ml	GLS Cont -19.5 %
ESV	78.4 ml	EF 47.6 %
CO	5.1 l/min	EDV 44.7 ml
SV	102.6 ml	ESV 85.3 ml
Mass	186.5 g	PreA 998 ms
AVC	367 ms	

**Summary**

LV	Value	1	2	3	4	5	Method
GLS	-14.7 %	-9.0	-11.3	-13.3	-19.9	-19.9	Avg
EF	51.1 %	46.8	46.3	49.1	56.7	56.7	
EDV	169.0 ml	158.6	157.5	167.1	181.0	181.0	

The screenshot shows the Siemens Healthineers Patient Report interface. The 'Edit Report' tab is selected and highlighted with an orange box. The interface displays a circular strain map and a table of cardiac parameters.

**Last Single Beat Analysis**  
LV End Systolic Strain [%] - Beat 1 of 1


	LV	LA
GLS	-19.9 %	GLS Res 35.0 %
EF	56.7 %	GLS Cond -15.5 %
EDV	181.0 ml	GLS Cont -19.5 %
ESV	78.4 ml	EF 47.6 %
CO	5.1 l/min	EDV 44.7 ml
SV	102.6 ml	ESV 85.3 ml
Mass	186.5 g	PreA 998 ms
AVC	367 ms	

**Summary**

LV	LA	RV	RA	
GLS	-14.7 %	GLS 35.0 %	GLS -15.4 %	GLS 43.0 %
EF	51.1 %	GLS Res 55.8 %	EF 37.7 %	EF 43.6 %
EDV	169.0 ml	GLS Cond -35.5 %	EDV 157.2 ml	EDV 55.2 ml
ESV	82.2 ml	GLS Cont -20.3 %	ESV 97.9 ml	ESV 97.8 ml
CO	4.4 l/min	EF 51.0 %	SV 59.3 ml	HR 50 bpm
SV	86.9 ml	EDV 38.2 ml	HR 50 bpm	

# 4D Heart<sup>AI</sup> edit report layout - bullseye

- The last bullseye captured is displayed in the report
- Bullseye displays the Timing / Phase value selected in system configuration or during the exam from the touch screen
- Segments cannot be removed from the bullseye in the report
- 4D strain and quantification of left heart chambers are listed
- Bullseye results cannot be edited; can only delete per chamber by selecting **X**



# 4D Heart<sup>AI</sup> edit report layout - worksheet

## Worksheet includes last five instances captured



- Average Result enabled - individual beat results sent to report as a separate instance
- Average Result disabled - only current beat result sent to report
- Results cannot be edited, only deleted
- Method can be set to Last, Average, or select the instance 1 thru 5
- Results are listed by chamber (LV, LA, RV, RA)

The screenshot displays the Siemens Healthineers Patient Report interface. The left sidebar shows 'Patient History' with 'Cardiac' expanded, and '4D Heart AI' highlighted. The main area shows a table of cardiac parameters for LV and LA chambers. The table has columns for 'Value' and five instances (1-5), along with a 'Method' dropdown menu. A 'Delete' dialog box is open at the bottom, asking 'Do you want to delete the measurement?' with 'OK' and 'Cancel' buttons.

Chamber	Parameter	Value	1	2	3	4	5	Method
LV	GLS	-14.7 %	-9.0	-11.3	-13.3	-19.9	-19.9	Avg
	EF	51.1 %	46.8	46.3	49.1	56.7	56.7	Last
	EDV	169.0 ml	158.6	157.5	167.1	181.0	181.0	Avg
	ESV	82.2 ml	84.4	84.6	85.0	78.4	78.4	1
	CO	4.4 l/min	3.8	3.6	4.2	5.1	5.1	2
	SV	86.9 ml	74.2	72.9	82.1	102.6	102.6	3
	Mass	160.9 g	139.8	155.6	136.0	186.5	186.5	4
LA	HR	50 bpm	51	50	51	50	50	5
	GLS	35.0 %					35.0	Method
	GLS Res	55.8 %	78.1	40.3	69.8	35.0		
	GLS Cond	-35.5 %	-57.7	-21.5	-47.3	-15.5		
	GLS Cont	-20.3 %	-20.4	-18.8	-22.5	-19.5		
	EF	51.0 %	48.4	55.7	55.5	47.6	47.6	
	EDV	28.2 ml	36.0	37.8	31.0	44.7	44.7	

# Report control options


- **Summary & Edit Report** – toggle between edit report function and report summary
- **Export to PDF** – exports a pdf version of the report to a USB
- **Store Report** – saves an image of the entire report
- **Print Report** – prints the report to the configured printer
- **Collapse All** – collapses all sections contained within the report



# 4D Heart<sup>AI</sup> report summary

The report summary is displayed when selecting **Summary** from the Patient Report

- 4D strain results
- Cardiac chamber results from 4D dataset
  - GLS, EF, and volumes
  - LV cardiac output (CO), and LV mass



# Summary

- **Overview of 4D Heart<sup>AI</sup> application**
- **Review 4D Heart<sup>AI</sup> configuration and touch screen layout**
- **Illustrate imaging for 4D Heart<sup>AI</sup>**
- **Review 4D Heart<sup>AI</sup> workflow**
- **Evaluate 4D Heart<sup>AI</sup> analysis and report**



## AI-driven analysis

- Semi-automated contour placement for heart chambers
- AI-enabled size and function analysis
- Ability to edit when needed

## 4D chamber quantification plus left heart strain

- Analysis for irregular heart beats can be performed
- Ability to edit ED/ES frames, AVC, and Pre-A
- Average up to 5 beats

## Physio and non-physio workflow

- Analysis for irregular heart beats can be performed
- Ability to set ED/ES frames
- Average up to 5 beats
- ***Only vendor to have non-physio feature!***

1. Mitchell, Carol, et al., *Journal of the American Society of Echocardiography* (2019) 32;1
2. Hasanfatta, Naeem et al., *ACTA Scientific Clinical Case Reports* (2022) 3;6: 45-56
3. Kosarju, Ateet, et al., *Left Ventricular Ejection Fraction*. StatPearls Publishing (2023) Jan

# Trademarks and disclaimers

The products/features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

ACUSON Origin is a trademarks of Siemens Medical Solutions USA, Inc. At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey towards expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated five million patients worldwide everyday benefit from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics and molecular medicine as well as digital health and enterprise services.

We're a leading medical technology company with over 120 years of experience and 18,500 patents globally. With over 50,000 employees in more than 70 countries, we'll continue to innovate and shape the future of healthcare.

---

## **Siemens Healthineers Headquarters**

Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
Phone: +49 9191 18-0  
siemens-healthineers.com

## **Manufacturer**

Siemens Medical Solutions USA, Inc.  
Ultrasound  
22010 S.E. 51st Street  
Issaquah, WA 98029, USA  
Phone: 1-888-826-9702  
siemens-healthineers.com/ultrasound

Thank you for your enthusiasm!

# Questions?