



# OKM 801

## INFANT INCUBATOR

## OKM 801 Infant Incubator

- Two control modes: Air mode and Baby (Servo) mode controlled by microprocessor,
- Trend Information can be stored, controlled and displayed by the LCD screen,
- Optional Oxygen concentration level display and servo control system,
- ">37°C" override heat setting function,
- Heated servo humidification setting and display and drawer-shaped integrated water tank with level indicator,
- Optional inbuilt weighing scale & Under-bed X-ray cassette tray,
- Embedded integrated sensor module, with dual thermistors for operational safety
- Adjustable continuous bed tilt
- Air curtains on both front and rear outlets,
- Internal battery for alarms and monitoring, inbuilt drawers, and optional IV pole and monitor tray
- Optional integrated weighing scale, LED examination lamp with spot phototherapy and electrical height adjustment
- Audiable and visible alarms, dual temperature measurement for twins
- Optional wireless and USB connection for data output, and HL7 ready,
- Optional sliding rotating bed with safety lock, double wall, access door covers, slow dumping door mechanism
- Optional slave display, Nellcor SpO2 and NIBP measurements, Patient circuit holder, access ports,
- Easy access to the incubator from four sides
- Fully removable control panel in case of any malfunction
- Rotary knob available for Manual O2 Control


### 01

#### Air and Oxygen Circulation and Air Curtain Features



- The air circulation system through front and rear air outlets including left-right air return, the incubator has automatic accelerated vertical air curtains when the cabin doors are open, thereby minimizing the heat loss inside the canopy and providing a stable micro environment to infants.
- Optional Oxygen supply can be calibrated with a single step and the oxygen concentration can be quickly increased to the set value. Servo controlled or manual flowmeter setting options are available.

## 02 Heated Humidification Function



- The high temperature steam humidification method is used to eliminate common pathogenic bacteria in the water tank, greatly reducing the risk of infection.
- Water level indicator is available.

- The drawer type water tank, placed on the outside of the infant's bed, can be removed as a whole for easy cleaning and disinfection. Autoclavable water chamber option is available.
- With the humidification system, a stable thermal environment with the required humidity level is provided. It works perfectly to prevent transepidermal water loss in newborns. It allows clinicians to adjust the relative heated humidity level inside the canopy.

|                            |               |
|----------------------------|---------------|
| Water tank capacity        | 1500 ml       |
| Humidity display range     | %0 RH~%100 RH |
| Humidity control range     | %30 RH~%95 RH |
| Humidity control precision | % ±5 RH       |

## 03 User-Friendly, Comfortable Treatment Environment

- Pressure-relief waterproof washable mattress is made of polyurethane material, which is soft and breathable, thereby improving the comfort of newborns.



- With the use of biocompatible materials, sensitivity, irritation and toxicity formation on the baby's skin are prevented and the skin integrity is preserved.

**04 Trendelenburg Feature**



- The angle of the bed can be adjusted through the touch screen and adjust the trendelenburg level. Thus, it provides smooth continuous positioning to reduce patient trauma from shaky manual movement. It effectively reduces the reflux and vomiting problems.

The trendelenburg directions and level can be adjusted as desired by pressing the right and left trendelenburg buttons on the main screen. Bed inclination can be reset with one touch button in the middle.

**Manual trendelenburg optional**

- By turning the left Trendelenburg manual adjustment knobs on the sides of the device, the bed can be easily the desired angle.



Manual Trendelenburg



**05 X-Ray Tray**



- The X-ray cassette can be placed under the bassinet, eliminating the need to move the newborns and reducing the risk of infection.

## 06 Bed Mechanism and Canopy Damping System (Optional)



- Sliding out mattress tray with safety mechanism.
- The panel damping system allows the panels to fall slowly and silently without the need to hold them.
- Double protection design is available to prevent the front door from opening accidentally.

## 07 User Interface and External Display Feature (Optional)



- Optional dual-screen display for alternative user control and monitorization,
- External 12.1 inch TFT color touch screen, 360° adjustable angle, clearly visible from far distance, fast switching between multiple languages, convenient for medical personnel to monitor and use.
- Medical personnel can easily observe the treatment of infants in the incubator remotely.
- It has single piece injection molding, smooth appearance, long service life, and robust structure.
- Screen brightness can be adjusted, and power on/off for LED examination and spot phototherapy lamps.
- Incubator air temperature, Dual skin temperature, heater power, baby's weight, phototherapy timer, SpO2, Pulse rate, humidity, oxygen concentration and alarms can be monitored on the inbuilt 8 inch color Touch LCD screen.
- Display function and external data output. It is also HL7 ready.
- In addition, the data of measured parameters are recorded in 2, 8, 24 hours and upto 7 days period.
- Can be viewed graphically and numerically on the Trend menu screen.

**08 Storage Feature**



- The large-capacity drawer is ergonomic and can conveniently take cords and other accessories.

Large drawer maximum 5 kg  
Two small drawers max 5 kg each

**09 Electrical Height Adjustment Feature (Optional)**



- The height of the device can be adjusted through electrically powered foot pedals. It allows easy use for clinical staffs.

- By pressing the pedal in the direction of the up arrow, the height of the device can be brought to the desired height, or by pressing the pedal in the direction of the down arrow, the height of the device can be reduced.

**10 Swivel Bed**



- The innovative design with the Slewing Bearing allows for seamless rotation of the incubator bed, providing easy access to the baby from any angle for medical professionals and caregivers.
- Its smooth and silent motion ensures a soothing environment for the infant, fostering a sense of security and aiding in their recovery journey.

**11 Examination and Spot Phototherapy Lamp ( Optional )**

- The examination lamp illuminates the bed evenly and helps the medical staff to observe and care even in a dark environment. Also spot LED Phototherapy is available as option with timers in the control panel.



## 12 NIBP (Optional)

- Easy, convenient, reliable measurement and monitoring of parameters such as systolic blood pressure, mean arterial pressure (MAP), diastolic blood pressure and pulse,
- Manual, automatic measurement mode (1/2/3/4/5/10/15/30/60/90/120/180/240/480 min and 5 min continuous mode,
- Adjustable upper and lower alarm limits for systolic, diastolic blood pressure and MAP measurements,
- Real-time transmission of module working status, including hardware, software and sensor status,
- A unique algorithm with the ability to resist motion interference and weak signal measurement,
- It has the function of gas leakage detection to detect whether the gas path is leaking
- With double overvoltage protection (hardware and software overvoltage protection)
- Double timeout protection

NIBP Indicator Panel



## 13 Nellcor SpO2 (Optional)


- It offers clinicians the opportunity to safely detect respiratory complications early and intervene promptly. Incorporates the latest Nellcor™ digital signal processing technology for accurate, reliable readings even during low perfusion, motion and other forms of signal interference.

It reacts to the patient's status with technology that displays the patient's oxygenation and pulse more quickly than other technologies.

Provides real-time information on patients' respiratory status, including continuous SpO2 and pulse rate monitoring and trend data

Displays realtime SpO2 and pulse rate measurements, plethysmographic waveforms and pulse amplitude  
Includes SatSeconds alarm management, a clinician-controlled feature that can distinguish between real, clinically significant events and transient events by taking into account both the severity and the duration of any desaturation event.

SpO2 Indicator Panel





## SatSeconds alarm management: How does it work?

The SatSeconds™ function can be activated by selecting a SatSeconds limit, or “clock,” of 10, 25, 50 or 100 SatSeconds. Clinicians who choose to employ the SatSeconds function should select a limit suited to their clinical environment and patient conditions. The SatSeconds function is the product of magnitude and time a patient’s saturation exceeds SpO2 alarm limits.

If the oxygen saturation is outside the limits, the SatSeconds feature calculates: [SATURATION POINTS] × [SECONDS].

For example, in the graph below the alarm threshold is set at 95 and SatSeconds is set to 25.

The SatSeconds clock is clearly visible on the monitor. The clock begins to fill as the monitor starts to track a desaturation event.



### EVENT 1

The patient’s SpO2 drops to 86 percent and the duration of the event is two seconds before the saturation returns above the low alarm limit of 92 percent.

6 percent drop below the Low Alarm Limit  
x 2 second duration

12 SatSeconds

Because the SatSeconds Alarm Limit is set at 25 and the actual number of “SatSeconds” equals 12, there is no audible alarm.



### EVENT 2

15 SatSeconds




### EVENT 3

During this event, the patient’s saturation drops to 82 percent, which is 10 percent below the low alarm limit of 92 percent. Since the patient does not return within 2.5 seconds, there is an audible alarm.

10 percent drop below the Low Alarm Limit  
x 2.5 seconds (maximum time allowed)

25 SatSeconds

At this level of saturation, the event would only be able to last for 2.5 seconds. However, the patient’s saturation did not return within that amount of time. Therefore, an audible alarm is heard just under three seconds into the event.



| <b>Technical Specifications</b>                                   |  |
|---|--|
| <b>Screen Specifications</b>                                      |  |
| Built-in Screen   | 8" color touch TFT LCD screen with IP2X Protection                                     |
| Optional-Incubator Integrated External Screen                     | 12.1" color touch TFT LCD screen with IP6X Protection                                  |
| <b>Physical Specifications</b>                                    |  |
| Height  | 144 cm (Fixed base ) / 131 cm to 156 cm (E-base); Trolley height: 880 to 1100 mm       |
| Width   | 100.5 cm   |
| Depth   | 57.7 cm  |
| Weight  | 110 kg with (Fixed base) / 122 kg with (E-base)  |
| <b>Canopy Specifications</b>                                      |  |
| Inlet Hole (Grommet)  | 10 pieces  |
| Intervention Window   | 6 pieces   |
| incubator covers  | 2 long covers; 1 front and 1 back cover (Can be opened 180 degrees )                   |
| Mattress Sponge Size  | 75 cm x39cm x 3 cm   |
| Bed Tilt  | ±13 ° base specification   |
| <b>Electrical Specifications</b>                                  |  |
| Power source  | AC 230V±%10V 50Hz or AC 110V±%10V 60Hz   |
| Maximum heater output power                                       | 400W/240V or 400W/120V   |
| Auxiliary power output  | AC 230V±10V 50Hz or AC 110V±10V 60Hz, Max. Current 3A                                  |
| <b>Temperature Control Modes</b>                                  |  |
| Heater indicator  | 0~100% (10% adjustment)  |
| Temperature control mode  | Air mode<br>Baby mode (Servo)  |
| Air temperature control range                                     | 20~37°C with +/- 0.2 C accuracy<br>37~39°C (User controlled mode) Override 37°C        |
| Air temperature display range                                     | 10-50 C with 0.1 resolution  |
| Baby temperature control range                                    | 34~38°C<br>37~38°C (User controlled mode)  |
| Skin temperature display range:                                   | 10-50 C with 0.1 resolution  |
| Skin Mode Control Accuracy  | ±0.1°C   |
| Dual Skin Temperature Monitoring                                  | Yes  |
| Temperature rise time* (ambient temperature: +25°C)               | ≤20 min  |
| Temperature variability   | ≤0.3°C   |
| Temperature uniformity (bed placed horizontally)                  | ≤0.8°C   |
| Sensitivity of the skin temperature sensor                        | ±0.2°C   |
| <b>Oxygen Concentration Control (Optional)</b>                    |  |
| Oxygen concentration display/measurement range                    | 0%~100% with +/- 1% resolution   |
| Oxygen concentration display/measurement accuracy                 | ±2% (preset oxygen concentration: < 25%)<br>±3% O2 (preset oxygen concentration: >25%) |
| Oxygen concentration control accuracy                             | ±1% O2 volume concentration  |
| Oxygen concentration control range                                | %21~%70  |
| <b>Scale Specifications (Optional)</b>                            |  |
| Weight range  | 0 kg~ 10 kg  |
| Weight display sensitivity  | 1g   |
| Weight display accuracy   | ±5 g   |
| <b>Humidity Control</b>   |  |
| Working time of the water tank after filling with distilled water | ≥ 24 hours maximum 85% RH and 36°C, in Air Mode  |
| Water tank capacity   | 1500 ml  |
| Humidity display range  | 0%RH~100%RH  |
| Humidity control range  | 30%RH~95%RH in 1% increments   |
| Humidity control precision  | ±5%RH  |
| Resolution  | +/- 1%   |
| Humidity sensor measurement accuracy                              | +/- 1.8%   |
| Target humidity range   | User selectable from 30%RH to 95%RH  |

| <b>SpO2 Specifications (Optional)</b> |   |
|---------------------------------------|---|
| SpO2 saturation measurement range     | %1-100  |
| PR measurement range                  | 20-250 bpm  |
| Perfusion measurement range           | 0.03-%20  |
| SpO2 measurement accuracy             | 70% to 100% ± 2<br>(Low sat) 60% to 80% ± 3<br>(Low Perfusion) 70% to 100% ± 2<br>(With Motion) 70% to 100% ± 3 |
| PR measurement accuracy               | 20 to 250 bpm ± 3<br>(Low Perfusion) 20 to 250 bpm ± 3<br>(With Motion) 20 to 250 bpm ± 5                       |

| <b>NIBP Specifications (Optional)</b> |                |
|---------------------------------------|----------------|
| SBP                                   | 40-130 mmHg    |
| MAP                                   | 20-100 mmHg    |
| DBP                                   | 10-90 mmHg     |
| Cuff Pressure                         |                |
| Range                                 | 0-300mmHg      |
| Accuracy                              | ±2mmHg or ± %1 |
| Resolution                            | 1 mmHg         |

| <b>Other Specifications</b>  |                                   |
|--|-----------------------------------|
| Noise level inside the incubator   | ≤ 45 dB at a constant temperature |
| Noise level in the cabin when operating at 65% oxygen in servo oxygen mode | ≤ 49 dB                           |
| In-cabin airflow rate  | ≤ 10 cm/s                         |
| Carbon dioxide (CO2) level:  | < 0.5%                            |
| Air flow:  | < 30 l/m                          |
| Baby Cot extendable size:  | 510 x 848 x 455 mm                |
| User interface with multiple languages                                     |                                   |

| <b>Physical Conditions</b> |             |
|----------------------------|-------------|
| Operating Range            | 20 °C-30 °C |
| Ambient Air elocity:       | <0.3 m/s    |

| <b>Trend Specifications</b> |  |
|-----------------------------|--|
| Trend duration              | 0 to 2, 4, 8, 24 hours up to 7 days  |
| Parameters                  | Canopy air temperature, first and second skin temperatures, oxygen rate, humidity rate, heater, weight |
| Display Type                | Numeric and Graphic  |

| <b>Standard and Optional Features</b> |  |
|---------------------------------------|--|
| <b>Standard</b>                       | <b>Optional</b>  |
| Servo Humidity Control                | Servo Oxygen Control   |
| Air Temperature Control               | VHA (Vertical Height Adjustment) Stand   |
| Skin Temperature Control              | Baby Scale   |
| 3x Drawer                             | External Touch Screen  |
| 2x Air Filter                         | IV pole (load capacity 10 kg)<br>Height Adjustable Monitor Tray ( 850 mm to 1627 mm, extendable to 2000 mm +/- 5%) |
| 2x Iris Access Window Cover           | Manual Trendelenburg   |
| 6x Access Window Cover (QT Sleeve)    | Nellcor SpO2 measurement module / SatSeconds   |
| 10x Grommets                          | NIBP   |
| APGAR Timer                           |  |
|                                       | Swivel Bed   |
|                                       | Examination and Phototherapy Lamps   |

### Functional Alarms

|                          |                                 |
|--------------------------|---------------------------------|
| Over Temperature         | Sensor Module is not in Place   |
| High Humidity Rate       | Sensor Module Position is Wrong |
| Low Humidity Rate        | Canopy Door is Open             |
| Water Case not In Place  | Low Air Temperature             |
| Water Level Low          | High Air Temperature            |
| Skin Probe Placed Wrong  | Air Temp Sensor Measure Fault   |
| Skin Probe Removed       | High Oxygen Rate                |
| Skin Probe Measure Fault | Low Oxygen Rate                 |
| High Skin Temperature    | Oxygen Sensors not Connected    |
| Low Skin Temperature     | Oxygen Sensors Measure Fault    |

### System Alarms

|                              |                                   |
|------------------------------|-----------------------------------|
| Skin Probe Fault             | Fan Motor Speed Fault             |
| High Air Heater Temperature  | Battery System Fault              |
| Air Temp Sensor Fault        | Air Heater System Fault           |
| Low Battery Voltage          | Sensor Module Communication Fault |
| Servo Oxygen System Fault    | Humidity Heater System Fault      |
| Power Failure                | Mainboard Communication Fault     |
| Scale Circuit Fault          | Battery is Not Exist              |
| Air Circulation System Fault | Sensor Module Fan Fault           |
| Humidity Sensor Fault        | Oxygen Leakage                    |

### SpO2 Alarms

#### Technical Alarms

|                                   |
|-----------------------------------|
| Replace Sensor                    |
| Pulse Search                      |
| Extended Update                   |
| Pulse Timeout                     |
| Sensor Failure                    |
| Interference                      |
| Sensor Disconnected               |
| Pulseoximeter Communication Fault |
| Sensor Off                        |

#### Physiological Alarms

|                 |
|-----------------|
| PI Value Low    |
| PI Value High   |
| Pulse Rate Low  |
| Pulse Rate High |
| SpO2 High       |
| SpO2 Low        |
| SpO2 INOP       |

### NIBP Alarms

#### Technical Alarms

|                          |
|--------------------------|
| NIBP communication fault |
| NIBP system fault        |
| Manual stop              |
| Measurement timeout      |
| Signal is saturated      |
| Signal is weak           |
| Pressure error           |
| Cuff is too loose        |
| Cuff air leakage         |

#### Physiological Alarms

|                            |
|----------------------------|
| Over pressure              |
| Excessive patient movement |
| Exceed measure scope       |
| High MAP pressure          |
| Low MAP pressure           |
| High diastolic pressure    |
| Low diastolic pressure     |
| High systolic pressure     |
| Low systolic pressure      |



# OKUMAN

OKUMAN Medikal Sistemler A.Ş.

Kazım Karabekir Caddesi 95/95 06060 Iskitler  
Ankara / Turkey

info@okuman.com.tr  
www.okuman.com.tr

Phone: +90 312 384 05 20  
Fax : +90 312 384 19 75

© 2023 **OKUMAN** Medikal Sistemler A.Ş.  
All rights reserved. Technical specifications are subject to change  
by **OKUMAN** Medikal Sistemler A.Ş. without prior notice.

CE  
1984

Catalogue No: OKM 801-000-032024-EN