INgezim® INFLUENZA A

Prod Ref: 10.FLU.K3

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials with Positive Control	1	1mL	2	1mL
Vials with Negative Control	1	1mL	2	1mL
Vials with Peroxidase Conjugate, ready to use	1	30mL	2	30mL
Bottles with Diluent (DE01-01), ready to use	1	125mL	1	125mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter).

TECHNICAL BASIS

The kit is based in the technique of blocking immunoenzymatic assay (**Blocking ELISA**) which is briefly described below:

Plates are coated with the antigen. After adding the sample to the well, if it contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the sample does not contain specific antibodies, they will not bind to the antigen.

If we add a specific monoclonal antibody against the viral antigen coated to the plate (conjugated with peroxidase), it will compete with the antibodies in the serum. If the sera samples contain specific antibodies, they will not allow the binding of the conjugate to the antigen; whereas if it does not contain specific antibodies, the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, presence or absence of labelled Mab can be detected by adding the substrate (TMB) which, in presence of the peroxidase, will develop a colorimetric reaction.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

Control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Samples must be sera. To be used in the kit, they must be diluted 1/2 in diluent supplied in the kit. This dilution can be made directly in the well by adding 50 µL of diluent and 50 µL of sample.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water. When ready, this solution remains stable between +2°C and +8°C.

Positive and Negative controls:

Controls must be tested as sample serum. Add 50 µL of controls to the appropriate wells containing 50 µL of diluent.

Attention! Once opened they are stable for one month between +2°C and +8°C. For long term storage, we recommend to make aliquots and store them at -20°C.

TEST PROCEDURE

1. All reagents must be brought to room temperature before use.
2. If dilution of samples is going to be made directly in the well, add 50 µL to the required wells and then, 50 µL of controls and samples to the appropriated wells containing the diluent. Mix well contents to assure a perfect homogenization. Be careful to prevent cross-contamination. It is recommended to assay controls in duplicate. Cover the plate with a lid and **incubate at 37°C for 1 hour.**
3. Shake out contents and wash 4 times as it is described previously.
4. Add 100 µL of conjugate to each well. Cover the plate with a lid and **incubate for 30 minutes at 37°C.**
5. Shake out contents and wash 5 times as it is described previously.
6. Add 100 µL of substrate to each well. The use of multichannel pipette is recommended to accurate time and uniformity. **Incubate at room temperature for 10 min.**
7. Add 100 µL of stopping solution to each well in the same order and at the same speed as the substrate solution.
8. Read the absorbance of each well using an ELISA reader with a filter of 450 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

VALIDATION CRITERIA:

The test is considered valid if:

$$\frac{\text{OD Negative control}}{\text{OD Positive control}} \geq 3$$

RESULTS INTERPRETATION:

There are two different ways for interpretation depending on the preferences:

$$\% \text{ Competition} = 100 - \frac{\text{OD Samples}}{\text{OD Negative control}} \times 100$$

Results expressed as competition percentage:

Avian:

- Samples will be considered positive if the percentage of competition is higher than 45%.
- Samples will be considered negative if the percentage of competition is lower than 40%.
- In case of percentages between 40-45% a new evaluation is recommended with a new extraction after 2 weeks.

Swine and Horse:

- Samples will be considered positive if the percentage of competition is higher than 30%.
- Samples will be considered negative if the percentage of competition is lower than 30%.
- In case of percentages between 25-35% a new evaluation is recommended with a new extraction after 2 weeks.

Results expressed as ratio S/N:

$$S/N = \frac{\text{OD of the sample}}{\text{OD of negative control}}$$

Avian:

- Samples will be considered positive if the S/P < 0.55
- Samples will be considered negative if the S/N > 0.60
- In case of S/N between 0.55 and 0.60, a new evaluation is recommended with a new extraction after 2 weeks.

Swine and Horse:

- Samples will be considered positive if the S/N < 0.70
- Samples will be considered negative if the S/N > 0.70
- In case of S/P between 0.75 and 0.65, a new evaluation is recommended with a new extraction after 2 weeks.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

Distribuido en
Distributed in

por
by