# On•Call<sup>®</sup> Sure

**Blood Glucose Monitoring System** 

NEW **ISO** 15197



- GDH-FAD Enzyme
- 10-70% HCT Range
- 0.6μL Capillary/Venous/ Arterial/Neonatal Sample
- No-coding
- 5s Testing Time
- US 510(k) & CE










### **Accurate and Reliable**

Extensive clinical studies demonstrate the accuracy of the *On Call® Sure* Blood Glucose Monitoring System with fresh capillary, venous, arterial and neonatal blood samples can comply with **ISO15197:2013/EN ISO 15197:2015.** 





Plasma YSI Reference Glucose (mg/dL)



#### Consensus Error Grid Analysis Clinical Trial - Fingertip Capillary Blood, by Technician Acon On Call® Sure Blood Glucose Monitoring System

710011 OII OUI OUI DIOGE CIECOSC I TOINTOI III G 9 STEELII						
System Accuracy Results for Glucose Concentration ≥100mg/dL						
Within ±5%	Within ±10%	Within ±15%				
335 / 582 ( 57.6% )	530 / 582 ( 91.1% )	582 / 582 ( 100.0% )				
System Accuracy Results for Glucose Concentration <100mg/dL						
Within ±5mg/dL	Within ±10mg/dL	Within ±15mg/dL				
156 / 216 ( 72.2% )	214 / 216 ( 99.1% )	216 / 216 ( 100.0% )				
System Accuracy Results for both Glucose Concentration ≥100mg/dL and <100mg/dL						
Within ±15% or ±15mg/dL						
798 / 798 ( 100.0% )						
<u> </u>	-	·				

#### Consensus Error Grid Analysis Clinical Trial - Venous Blood, by Technician Acon On Call® Sure Blood Glucose Monitoring System

7100/7 Cir Car Carc Bioca Clacose Homesting System						
System Accuracy Results for Glucose Concentration ≥100mg/dL						
Within ±5%	Within ±10%	Within ±15%				
253 / 522 ( 48.5% )	448 / 522 ( 85.8% )	512 / 522 ( 98.1% )				
System Accuracy Results for Glucose Concentration <100mg/dL						
Within ±5mg/dL	Within ±10mg/dL	Within ±15mg/dL				
129 / 240 ( 53.8% )	206 / 240 ( 85.8% )	240 / 240 ( 100.0% )				
System Accuracy Results for both Glucose Concentration ≥100mg/dL and <100mg/dL						
Within ±15% or ±15mg/dL						
752 / 762 ( 98.7% )						

#### Consensus Error Grid Analysis Clinical Trial - Neonatal Blood, by Technician Acon On Call® Sure Blood Glucose Monitoring System

System Accuracy Results for Glucose Concentration ≥ 100mg/dL						
Within ±5%	Within ±10%	Within ±15%				
70 / 120 ( 58.3% )	104 / 120 ( 86.7% )	120 / 120 ( 100.0% )				
System Accuracy Results for Glucose Concentration <100mg/dL						
Within ±5mg/dL	Within ±10mg/dL	Within ±15mg/dL				
134 / 234 ( 57.3% )	204 / 234 ( 87.2% )	234 / 234 ( 100.0% )				
System Accuracy Results for both Glucose Concentration ≥100mg/dL and <100mg/dL						
Within ±15% or ±15mg/dL						
354 / 354 ( 100.0% )						







### **Advanced Features for Diabetes Management**









## **Specifications**

Feature	Specification Sp	
Technology	Biosensor/Electrochemical, Glucose dehydrogenase(GDH-FAD)	
Testing Range	10-600mg/dL (0.6-33.3mmol/L)	
Calibration Method	No coding	
Result Calibration	Plasma-equivalent	
Test Time	5 seconds	
Sample Size	0.6µL	
Sample Type	Fresh whole blood(Capillary, Venous, Arterial and Neonatal)	
Test Sites	Fingertip, Forearm, Palm	
Hematocrit Range	10-70%	
Memory Storage	1,000 results with data and time	
Test Averaging	7,14,30,60,90-day averages	
Data Transfer	USB	
Control Solution	3 levels	
Audio Feature	Optional beep for sample detection, error messages	
Automatic Shutoff	2 minutes after last action	
Battery	Two(2) CR 2032 3.0V coin cell batteries	
Battery Life	3,000 measurements	
Operating Conditions	5-45 C and 10-90% relative humidity	
Strip Storage Temperature	2-35℃	
Expiration Date	24 months (18 months after first opening)	
	Automatic detection of insufficient sample	
Other Functions	2nd sample application allowed	
	Sample perturbation detection	

## **Catalog**

Product Name	Catalog No.	Con	tents
On Call® Sure Blood Glucose Monitoring System	G115-10U √+	1 Meter 10 Test Strips 10 Lancets 1 Lancing Dev 1 Carrying Case 1 Warranty Ca	500 Section 1997 - 140 Section 1997
On Call® Sure Blood Glucose Meter	G115-11U √+	1 Meter 1 Lancing Dev 1 Carrying Case 1 Warranty Ca	
	G115-12U √+	1 Meter 1 Control Solu 1 Manual 1 Carrying Ca	
	G115-13U √+	1 Meter 1 Control Solu 1 Carrying Case 1 Warranty Ca	
	G135-10U √+	25 Test Strips (25/vial)	1 Package Insert
		50 Test Strips (25/vial)	1 Package Insert
		100 Test Strips (25/vial)	1 Package Insert
		50 Test Strips (50/vial)	1 Package Insert
On Call® Sure Test Strips		100 Test Strips (50/vial)	1 Package Insert
		200 Test Strips (50/vial)	1 Package Insert
	G135-11U √+	25 Test Strips (Individually Foil Po	, 0
		50 Test Strips (Individually Foil Po	, ,
		100 Test Strips (Individually Foil F	ouch) 1 Package Insert
On Call® Lancets	G124-10A √ +	100 Lancets (25/bag)	
On Call® Lancing Device	G124-11A √ +	1 Lancing Device	1 Package Insert
On Call® Sure Control Solution	G125-12U √+	1 Control Solution 0 1 Control Solu	•
On Call® Diabetes Management Software Kit	G124-13A +	1 USB Data Transfer Cable	1 Installation Disk





**ACON Laboratories, Inc.** 5850 Oberlin Drive, #340, San Diego, CA 92121, USA Tel: 1-858-875-8000 Fax: 1-858-200-0729 E-mail: info@aconlabs.com Please visit our website for details: <u>www.acondiabetescare.com</u>

© 2022 ACON Laboratories, Inc. 1151029404