

Synergy^{Resection™} System

User Manual

Important Product Information



DFU-0213-EO

Revision 4 2025-08

CE 2797

Symbols glossary can be found at arthrex.com/symbolsglossary



Synergy^{Resection}™

The Synergy^{Resection} system *User Manual* provides safety operation information for all components of the Synergy^{Resection} console (AR-8305), including the accessories. All operating personnel must read this *User Manual* thoroughly prior to using this system and follow all safety warnings, cautions, and notes.

05	General Warnings and Safety Notices - Read This First	25	Buttons
05	Important Symbols and Conventions	25	Buttons - Main Screen
08	Symbol Definitions	27	Speed Test Mode
08	Shipping, Unpacking, and Warranty Information	27	Speed Adjustment
		27	Oscillation Mode Adjustment
09	Product Description		Accessory Handpiece Speed Limits and Deltas
09	Product Description and Intended Use	28	Drill Speeds
10	Product Features	28	Console Speed Controls
10	AR-8305 Console - Front View	28	Attachment Motor Control
11	AR-8305 Console - Rear View	28	Definition of Direction
12	AR-8305 Display Messages	28	Forward, Reverse and Oscillation Modes
14	Applied Parts and Accessory Features	29	Cruise Control
14	Standard Footswitch (AR-8310)	30	Toggle Function
15	Single Pedal Footswitch (AR-S8311)	31	Footswitches
		31	Standard Footswitch
16	Technical Specifications	31	Gas Pedal Footswitch
16	Console	31	Multi-Function Footswitches
16	Standard Footswitch	31	Accessory Handpieces
16	Single Pedal Footswitch		
16	Ambient Conditions for Operation	32	Cleaning and Disinfecting
16	Ambient Conditions for Storage (in shipping packaging)	32	Console AR-8305
		32	Footswitches
17	Safety, Electromagnetic Compatibility (EMC), and Regulatory Requirements	33	Accessory Cable
		33	Accessory Handpieces
18	Setup	34	Sterilization
18	How to Set Up the Console	35	
18	AC Power Safety Considerations		Transmissible Spongiform Encephalopathy (TSE) Agents
18	Replacing the Fuses	36	
19	Electromagnetic Compatibility	36	Maintenance
19	Basic Setup Procedure for the AR-8305 Console	36	Life Expectancy
			Periodic Maintenance
20	Setting up a Footswitch	36	
20	Setting up the Shaver Handpieces	36	Technical Support
20	Setting up Other Accessory Handpieces	36	How to Display the Software Version
20	Shutdown Procedure	37	Additional Technical Information
21	Operation and Frequently Used Functions	37	
21	Main Screen		Troubleshooting
22	Settings Menu Main Screen	38	Troubleshooting Interference with Other Devices
22	Language Selection Screen		
23	Oscillation Mode Selection Screen	38	Repair Policy
23	Footswitch Override Screen		
24	Reset Defaults Screen	39	End of Life, Environmental Directives
24	Information Screen	43	Electromagnetic Emissions
			Contact Information

List of Figures

10	Figure 1: Front Panel of Console
11	Figure 2: Rear Panel of Console
14	Figure 3: Standard Footswitch - Description
15	Figure 4: Single Pedal Footswitch – Description
21	Figure 5: Main Screen Display - One Tool Attached
21	Figure 6: Main Screen Display - Two Tools Attached
22	Figure 7: Settings Menu Display - Operation
22	Figure 8: Language Selection Display - Operation
23	Figure 9: Oscillation Mode Selection Screen - Operation
23	Figure 10: Footswitch Override Selection Screen - Operation
24	Figure 11: Reset Defaults Screen - Operation
24	Figure 12: Information Screen - Operation
25	Figure 13: Single Handpiece Screen Oscillation Mode Location
26	Figure 14: Two Handpiece Buttons - Operation

List of Tables

10	Table 1: Front Panel Elements
11	Table 2: Rear Panel Elements
12	Table 3: AR-8305 Display Messages
13	Table 4: *TOOL FAILURE* Error Code Definitions
14	Table 5: Elements of the Standard Footswitch
15	Table 6: Elements of the Single Pedal Footswitch
16	Table 7: AR-8305 Console Specifications
16	Table 8: Standard Footswitch Specifications
16	Table 9: Single Pedal Footswitch Specifications
16	Table 10: AR-8305 Ambient Conditions for Operation
16	Table 11: AR-8305 Ambient Conditions for Storage
17	Table 12: Safety, EMC, and Regulatory Requirements
21	Table 13: Main Screen Display Elements
22	Table 14: Settings Menu Display - Operation
22	Table 15: Language Menu Display - Operation
23	Table 16: Oscillation Mode Display - Operation
23	Table 17: Footswitch Override Display - Operation
24	Table 18: Reset Defaults Display - Operation
24	Table 19: Information Screen Display - Operation
25	Table 20: Oscillation Mode - Operation
26	Table 21: Two Handpiece Buttons - Operation
27	Table 22: Accessory Handpiece Speed Limits and Delta
28	Table 23: Handpiece Accuracy
30	Table 24: Toggle Function
37	Table 25: Troubleshooting: Faults, their Causes and Solutions
39	Table 26: Guidance and Manufacturer's Declaration – Electromagnetic Emissions
39	Table 27: System Cables
40	Table 28: Guidance and Manufacturer's Statement - Electromagnetic Immunity
42	Table 29: Guidance and Manufacturer's Statement – Recommend Separation Distance between portable and mobile RF communications equipment and the Model AR-8305

General Warnings and Safety Notices - Read This First

Important Symbols and Conventions

It is imperative that the symbols and conventions listed below be clearly understood. The Synergy^{Resection} system *User Manual* identifies critical, important, and useful information using these symbols and conventions.

Users of this device are encouraged to contact their Arthrex representatives if they require a more comprehensive surgical technique.

 **WARNING!** The WARNING! is the most important safety symbol. It identifies **critical** information that must be followed precisely to avoid injury or death.

1.  Caution: Federal law restricts this device to sale by or on the order of a physician.
2.  This device is intended to be used by a trained medical professional. Use this device only for purposes described in the User Manual, under the supervision of a trained and licensed physician.
3. No modifications of the console (AR-8305) or accessories are allowed.
4. Use this device only under the supervision of a trained and licensed physician. This device should not be used by untrained personnel or used for indications other than those described in this User Manual.
5.  Do not open or attempt to service this system, as this may void your warranty. There are no user-serviceable parts inside.
 Removing the cover may introduce an electric shock hazard by exposing you to dangerous high voltages or other risks. If the system malfunctions, return it for service immediately.

6. To avoid the RISK of electric shock, this equipment must only be connected to a SUPPLY MAIN with protective earth ground.
7. Do not have device in direct contact with patient if high-frequency devices are in use, or if the patient requires defibrillation.
8. DO NOT stack or place equipment adjacent to the AR-8305 console, if possible. If such a configuration is necessary, carefully observe the configuration in question to ensure that electromagnetic interference does not degrade performance.
9. Use only Arthrex approved accessories. Other accessories may result in increased emissions or decreased immunity of the system. Contact your Arthrex representative for a complete list of accessories. **DO NOT** modify any accessory. Failure to comply may result in patient and/or operating room staff injury.
10. When not in use, place active handpieces in a clean, dry, and highly visible area not in contact with the patient. Inadvertent contact with the cutting window of the attachment while in operation outside the intended articular cavities may cause physical/ thermal injuries to the user/patient, compromised sterility, and/or other ignition hazards within the environment of use.
11. When using a footswitch, ensure that the footswitch is not accidentally activated when a handpiece is not intended to be used (i.e., outside of a joint filled with fluid). Inadvertently depressing a footswitch if outside of a joint space, run dry, or resting on surgical drapes could result in overheating, smoke, and potentially a fire.

12. Arthrex warns against the use of the device without irrigation for any amount of time to ensure safety for patients and device users. Misuse of the device for only a few seconds with a complete lack of fluid will likely result in damage/adverse events such as fire/ smoke/burn/overheating/mechanical failure.
13. This equipment is **NOT** suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen rich or nitrous oxide environment.
14. After autoclaving, the accessory devices are VERY HOT. Handle with care to avoid burns.
15. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
16. Serious incidents should be reported to Arthrex, Inc. or an in-country representative, and to the health authority where the incident occurred.

The PRECAUTION! symbol identifies methods and procedures that must be followed to avoid damaging the device or causing it to malfunction.

1. **DO NOT** – under any condition or for any reason – open the console (AR-8305) or any other shaver accessory.
2. **Only use** replacement power cords that meet medical grade standards according to International Electrotechnical Commission (IEC) 60320-1, detachable power supply cords or electrical standards for the designated country where the AR-8305 console is being used. Contact your Arthrex representative for further information.
3. **Avoid** positioning the console so that it is difficult to disconnect coupler or plug from supply main.
4. *Connecting an extension cord to the AR-8305 may result in a reduced level of safety.*
5. **Always** use fuses with the correct values to avoid allowing overcurrent to enter the system.
6. An incorrect fuse may increase the risk of electrical shock or fire hazard.
7. This device has passed testing for Electromagnetic Interference (EMI) / Radio Frequency Interference (RFI) radiation and susceptibility and Electromagnetic compatibility (EMC). This device may cause interference to other devices in the near vicinity if not set up and used as Arthrex instructs.
8. **Do not** attach handpieces or footswitches during the Self Test or the Programming Modes.
9. **Only use** footswitches developed by Arthrex specifically for the AR-8305 Synergy^{Resection} console.
10. **Only use** shaver handpieces and cables that have been developed by Arthrex specifically for the Synergy^{Resection} console.
11. **Only use** handpieces and cables that have been developed by Arthrex specifically for the Synergy^{Resection} system. The footswitch cable connects and locks to the console to prevent accidental separation during use. To avoid damage, disconnect the footswitch by pulling on the cable connector shell (plug) only.
12. The accessory handpiece cable connects and locks to the console to prevent accidental separation during use. To avoid damage, disconnect the footswitch by pulling on the cable connector shell (plug) only.
13. **Always** comply with the instructions issued by the manufacturer of the **cleaning disinfectant regarding concentration, exposure times, temperature and material compatibility.**
14. **NEVER** clean the console with liquid cleaners. If there is dust, remove it with dry compressed air.

17. **Do NOT** clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.
18. The footswitch is **NOT** suitable to be cleaned and disinfected in a thermo washer disinfectant.
19. **NEVER** allow the console receptacles to have any contact with liquids. If there is dust or moisture on the receptacles, remove it with dry compressed air. **ONLY** dry connectors should be plugged into the console.
20. Refer to the Directions for Use (DFU-0154-XX) for detailed shaver handpiece cleaning and sterilization instructions. This can be obtained from the eDFU website at arthrex.com/eDFU, or by contacting your local Arthrex representative.
21. **NEVER** place the shaver handpiece in Cidex® or other aldehyde disinfectant solutions.
22. After sterilization in the autoclave, let the accessory device cool down slowly. **NEVER** use cold water to cool the handpieces. This will damage the electronic components and seals.
23. Liquid on the cable connector of the accessory device can damage the device. Before connecting the cable, ensure the receptacles are clean and dry.
24. Open the suction control completely to ensure proper circulation, if applicable.
25. Open the chuck completely to ensure proper circulation, if applicable.
26. Remove all accessories, including shaver blades, drill bits and saw blades before sterilizing a handpiece.
27. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or contact your Arthrex representative for an onsite demonstration.

NOTE: *This identifies information that can simplify the setup and operation of this device.*

1. *Read this User Manual thoroughly before operating the device and save it for future reference.*
2. *Extension cords must meet local electrical standard. Extension cords are not recommended to be used with the AR-8305.*
3. *The AR-8305 console incorporates a universal alternating current (AC) input power supply. A voltage selection switch is not required.*
4. AR-8305 has been designed to accept EMC from other devices within the limitations as described in the **Electromagnetic Emissions** section.
5. *The setup for a footswitch is the same for all models. The console detects which version is attached and allows the appropriate functions.*
6. While operating the handpieces, the settings can be changed without risk of damaging the motor or console.

In CE Accepting Countries: Procedures carried out using these devices may be used on the general population.

In CE Accepting Countries: The clinical benefits associated with the use of these devices outweigh the known clinical risks.

In CE Accepting Countries: There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

Symbol Definitions

All symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.



Handpiece connection

Shipping, Unpacking, and Warranty Information

 Carefully unpack and inspect all components for shipping damage. Any damage could compromise patient safety and should be reported immediately to Arthrex or any authorized Arthrex distributor. The warranty could be voided if shipping or first-installation damage is not reported within 7 business days of receiving the device. Refer also to our General Terms of Business.

The warranty is not valid if modifications are made to the product or repairs are completed outside of Arthrex or an authorized Arthrex distributor. Arthrex will answer any questions referring to the quality, reliability and/or shelf life of any product identified in this *User Manual*.

Product Description

Product Description and Intended Use

The AR-8305 Synergy^{Resection} system is a motorized, suction/cutting device used for the resection and aspiration of soft tissue, cartilage and bone during arthroscopic surgical procedures and spinal surgical procedures.

The console is microprocessor controlled and features two fully functional channels that allow optional simultaneous handpiece operation. Both channels can be operated with a single foot switch or individually when a second foot switch is attached to the console. The Synergy^{Resection} system has a maximum speed of 8,000 revolutions per minute (RPM) providing precise burr control during surgical procedures.

The Synergy^{Resection} console is designed for continuous operation. The handpieces are individually protected by a resetting thermal fuse. In the event of a motor overload fault (e.g., excess current) the fuse heats up and the device is deactivated. When the fault is removed and the fuse cools down (typically about 30 seconds) operation resumes at the pre-failure settings.

The cable and all handpieces can be sterilized, so handpieces can be changed during the surgical procedure without endangering surgical sterility.

The AR-8305 Synergy^{Resection} console and included Applied Parts and Accessories:

- › AR-8305 console with power cord
- › User Manual optional

Applied Parts

- › Footswitch controlled (FC) shaver handpieces
- › Handpiece controlled (HC) shaver handpieces
- › NanoResection™ small hub shaver handpiece

Other Optional Accessories:

- › Low-profile footswitch
- › Corded multi-function footswitch
- › Wireless multi-function footswitch
- › Single pedal footswitch

Save the packaging for later transport of the device.

 **WARNING!** Use this device only under the supervision of a trained and licensed physician. This device should not be used by untrained personnel or used for indications other than those described in this *User Manual*.

Product Features

AR-8305 Console - Front View

Figure 1 shows the front panel of the AR-8305 console. Features and symbols are identified in Table 1 below.




Figure 1: Front Panel of Console

Table 1: Front Panel Elements	
1	On / Off Power Switch
2	Handpiece Connectors
3	Synergy Product Logo
4	Handpiece Symbols
5	Touch Panel Visual Display
6	Footswitch Symbols
7	Footswitch Connectors
8	Arthrex Logo
9	Type body floating (BF) Symbols (Electric Shock Protection)
10	Power Switch IEC 60417-5010 symbol (On / Off)

AR-8305 Console - Rear View

Figure 2 shows the rear panel of the console. Features and symbols are identified in Table 2 below.




Figure 2: Rear Panel of Console

Table 2: Rear Panel Elements	
1	Main power input Plug
2	Equipotential Bonding Pin with Stamped Equipotential Bonding Symbol
3	Address
4	Fan
5	USB port (For use ONLY with thumb drive)
6	Serial ports for Arthrex integration
7	Serial number label
8	Model number label
9	Fuse holder for the Power Entry Module
10	Fuse label

AR-8305 Display Messages

The console's touch panel visual display (TPVD) Figure 1 [5] provides information about the status of the AR-8305 settings in real time. A list of informational and error messages is shown in Table 3 and 4 below. Messages that represent the active motor states are displayed in green on the TPVD.

Table 3: AR-8305 Display Messages

Message	Explanation
EEPROM Reset	Temporarily displayed first time micro-controller is programmed, or if internal memory has faulted and is reset.
Motor Fault	Temporarily displayed when motor fault condition occurs.
Overload	Temporarily displayed when motor over-current condition occurs.
Self Test Fail	Console has failed a critical power on test.
CRUISE NOT ALLOWED	Temporarily displayed when the cruise button is pressed with tool that doesn't allow cruise control.
FS PEDAL STUCK Remove FS, Repower System	Either a footswitch pedal is being pressed or is stuck while connecting to the console. Communication error results. Remove the footswitch and repower the system to reset the console.
FWD or REV ONLY	Temporarily displayed for tools that do not support oscillation (OSC) mode.
HP BUTTON STUCK Remove HP, Repower System	Either a handpiece button is being pressed or is stuck while connecting to console. Communication error results. Remove the handpiece and repower the system to reset the console.
TOGGLE NOT ALLOWED	Temporarily displayed when the toggle button is pressed with tool that doesn't allow oscillation mode.
TOOL FAILURE ERROR, XXX (Where "XXX" is an error code as defined in Table 4.)	Critical tool failure caused either by a faulty sensor or an over-current condition. Remove and replace the tool or repower the console to reset.
Unable to Install	Handpiece is not recognized.
[CC] ¹	Speed locked on cruise control.
CC	Cruise control enabled, about to lock.
AGG	Current mode: Aggressive oscillation mode (stopped).
Communication Error	Internal data communication error.
Critical Error	Critical error, unit will shut down.
DFS	Digital Footswitch is connected.
Drill	Drill is connected.
EFF	Current mode: Efficient oscillation mode (stopped).
FC Shaver	Footswitch controlling shaver is connected.
Fwd	Current mode: Forward (stopped).
Fwd ¹	Current status: Motor is moving forward.
Gas Pedal Footswitch	Temporarily displayed when the gas pedal footswitch is connected.
GFS	Gas pedal footswitch is connected.
HC Shaver	Handpiece controlling shaver is connected.
No Footswitch	Temporarily displayed when the footswitch is disconnected.
Osc ¹	Current status: Motor is oscillating.
Power Failure	Power supply error.
Rev	Current mode: Reverse (stopped).
Rev ¹	Current status: Motor is moving in reverse.
Run ¹	Current status: Motor is running (for saw).
Safe Mode	Excessive HP communication errors have occurred.
Sag Saw	Sagittal Saw is connected.
SFS	Standard Footswitch is connected.
SJ Shaver	NanoResection™ Small Hub Shaver Handpiece is connected.
Standard Footswitch	Temporarily displayed when a standard footswitch is connected.
STD	Current mode: Standard oscillation mode (stopped).

¹Displayed in green, denoting an active tool.

Table 4: *TOOL FAILURE* Error Code Definitions

Message	Explanation
# \$ @ Example: 022 Sensor 3 = Stuck Low Sensor 2 = No Errors Found Sensor 1 = No Errors Found	Represents Handpiece Sensor Failures # = Sensor 3, \$ = Sensor 2, @ = Sensor 1 0 = Stuck Low, 1 = Stuck High, 2 = No Errors Found.
H13	Multiple sensors shorted
H14	Multiple sensors stuck high/low
H15	Invalid Hall sensor states detected
C31	Over current 7.5A during retries
C32	Over current 5.0A during retries
C33	Over 5.0A for 6 seconds
C34	Over 7.0A for 3 seconds
C35	Over 7.5A for 1 second
C36	Over 8.0A for 250ms
C37	Over 8.1A in INTO
C38	Over 8.1A in A/D ISR
C39	Over 8.0A in FWD/REV ramp up
S51	Invalid Gear Ratio
S52	Over restart tries
S54	General system failure

Applied Parts and Accessory Features

Standard Footswitch (AR-8310)

Figure 3 shows the standard footswitch. Features and symbols are identified in Table 5 below.



Figure 3: Standard Footswitch - Description

Table 5: Elements of the Standard Footswitch

1	Cable
2	Speed Increase button
3	Toggle button
4	Forward pedal
5	Oscillation pedal
6	Reverse pedal

Single Pedal Footswitch (AR-S8311)

Figure 4 shows the single pedal footswitch. Features and symbols are identified in Table 6 below.




Figure 4: Single Pedal Footswitch – Description

Table 6: Elements of the Single Pedal Footswitch

1	Pedal – Activates handpiece in the selected mode and direction
2	Direction Button – Changes the direction of the handpiece between FORWARD, REVERSE, AND OSCILLATE*
3	GFS / SFS Mode Button – Switches between GFS (gas pedal footswitch) mode and SFS (standard footswitch) mode

- › After pressing the Direction Button, the direction setting indicated on the console display will not update until the pedal is activated. The user should activate the handpiece in a safe location to confirm the direction setting.
- › The Cruise Control and Toggle functions described later in this manual are not present on the Single Pedal Footswitch.

Technical Specifications

Console

Table 7: AR-8305 Console Specifications

Width	40.64 cm (16.00 inches)
Height	13.335 cm (5.250 inches)
Depth	30.686 cm (12.081 inches)
Weight	6.8 kg (15 lbs.)
Water protection	IP22
Mains cable	10 A/250 V
Power entry module	IEC 320/C13
Fuse value	6.3 A, 250 V~, 2.0 cm (0.75 inches) Type T
AC input	100-240V~, 50 - 60 Hz, 6.3A
Applied Part Type	BF
Cleaning, disinfecting and sterilization	Refer to the Cleaning and Disinfecting , and Sterilization sections

Standard Footswitch

Table 8: Standard Footswitch Specifications

Functions	Reverse, Forward, Oscillation, Speed Increase, and Toggle
Width	23.6 cm (9.3 inches)
Height	2.5 cm (1.0 inches)
Depth	20.8 cm (8.2 inches)
Weight	2.2 kg (5 lbs.)
Water protection	IP68
Cleaning, disinfecting, and sterilization	Refer to the Cleaning and Disinfecting , and Sterilization sections

Single Pedal Footswitch

Table 9: Single Pedal Footswitch Specifications

Functions	Reverse, Forward, Oscillation, and GFS/SFS Mode Change
Width	14.0 cm (5.5 inches)
Height	6.1 cm (2.4 inches)
Depth	21.3 cm (8.4 inches)
Weight	1.4 kg (3 lbs.)
Water protection	IP68
Cleaning, disinfecting, and sterilization	Refer to the Cleaning and Disinfecting , and Sterilization sections

Ambient Conditions for Operation

Table 10: AR-8305 Ambient Conditions for Operation

Temperature	10 °C to 40 °C (50 °F to 104 °F)
Relative humidity	30 % to 75 %
Air Pressure	700 hPa to 1060 hPa (21 in. Hg to 31.3 in. Hg)


Ambient Conditions for Storage (in shipping packaging)

Table 11: AR-8305 Ambient Conditions for Storage

Temperature	-30 °C to 70 °C (-22 °F to 158 °F)
Relative humidity	10 % to 90 %, non-condensing
Air Pressure	500 hPa to 1060 hPa (15 in. Hg to 31.3 in. Hg)

Safety, Electromagnetic Compatibility (EMC), and Regulatory Requirements

Table 12: Safety, EMC, and Regulatory Requirements

Parameter	Parameter Value	
System Classification	FDA Class	Class 1
	EU Class	Class IIA
	Health Canada Class	Class 2
Safety Certifications	Domestic Certification	UL 60601-1:2003
	Canadian Certification	CAN/CSA C22.2 No. 60601-1-2008-02
	EU Certification	EN 60601-1:2006 IEC 60601-1 3rd Ed. + AM1 (2012)
EMC Certifications	CISPR 11 EMC Class	Class A
	CISPR 11 EMC Group	Group 1
	EMC Certification	IEC 60601-1-2:2014 Class A
Safety Certification Marking		
CE Markings	MDD Classification	CE Marking for MDD93/43/EEC Class IIA
	MDD Marking	Rule 9

For all other accessories refer to their accompanying DFUs for more information.
Refer to the **Electromagnetic Emissions** section for further details on EMC Certification.

Setup

How to Set Up the Console

Users are encouraged to contact their Arthrex representative if they require a more comprehensive surgical technique.

AC Power Safety Considerations

The AR-8305 is powered by a medically rated universal AC input switching power supply. This power supply allows users to connect the console to any local AC mains outlet. Please use the appropriate plug and a reliable ground conductor.

Arthrex supplies separate power cords for the U.S. and Europe CEE 7/7 with the AR-8305. Contact your Arthrex representative if you need a power cord that must meet the electrical standards of another country.

⚠️ Only use replacement power cords that comply with medical grade standards, IEC 60320-1, Detachable Power Supply Cords or electrical standards for the designated country where the AR-8305 is being used. Contact your Arthrex representative for further information.

⚠️ Avoid positioning the console so that it is difficult to disconnect coupler or plug from supply main.

NOTE: Extension cords must meet the local electrical standard. Extension cords are not recommended to be used with the AR-8305.

⚠️ Connecting an extension cord to the AR-8305 may result in a reduced level of safety.

The console is designed to meet power-saving guidelines. The console has an AC mains switch on the front panel. When the AC mains switch is OFF, no electrical power is drawn by the console.

When the AC mains switch is ON, the console executes a series of self-diagnostic tests. Upon successful completion of these tests, the console displays the name and model number. If the tests discover a problem, the problem is shown on the display. Refer to Tables 3 and 4 for a listing of display messages.

In the event of an AC power interruption, the console can run continuously without a fault for up to 10 milliseconds. If an AC power failure lasts for longer than 10 milliseconds, the system will return to the **default settings** when AC power is restored.

⚠️ WARNING! To avoid the RISK of electric shock, this equipment must only be connected to a SUPPLY MAIN with protective earth ground.

⚠️ WARNING! Do not have the device in direct contact with patient if high-frequency devices are in use, or if the patient requires defibrillation.

Replacing the Fuses

The mains fuse is replaced with T6.3AL250v as follows:

1. Disconnect the device from the AC mains.
2. Open the fuse tray in the AC inlet, by pulling out on the tabs.
3. Replace the fuses with T6.3AL250v line fuses, as noted on the rear panel.
4. Push the fuse holder back into the AC inlet.
5. Ensure that the fuse holder is fully seated and that the tabs snap back.

⚠️ Always use fuses with the correct values to avoid allowing overcurrent to enter the system.

⚠️ An incorrect fuse may increase the risk of electrical shock or fire hazard.

NOTE: The AR-8305 console incorporates a universal AC input power supply. A voltage selection switch is not required.

Electromagnetic Compatibility

 **This device has passed testing for EMI / RFI radiation and susceptibility and EMC compatibility. This device may cause interference to other devices in the near vicinity, if not set up and used as Arthrex instructs.**

AR-8305 has been designed to accept EMC from other devices within the limitations as described in the **Electronic Emissions** section.

To determine if the AR-8305 is causing interference to other devices, power OFF the AC main power switch and then ON again.

Try to correct the interference by following one or more of these measures:

1. Reorient or relocate the receiving device.
2. Increase the separation between devices.
3. Connect the device to an outlet on a different circuit than the other device(s) are connected.
4. Consult the manufacturer or field service technician for receiving devices for guidance.

Basic Setup Procedure for the AR-8305 Console

*NOTE: The **Operation and Frequently Used Functions** section explains how to use the console.*

1. Place the AR-8305 on a flat, dry surface, such as the AR-6481 Arthrex arthroscopy pump cart.
2. Connect the receiver end of the power cord for the AR-8305 into the AC socket and the plug end to the facility AC mains supply.
3. Power ON the console.
4. Allow it to fully initialize.
5. Attach the shaver handpiece.
6. Attach a footswitch, if applicable.

 **WARNING! DO NOT** stack or place equipment adjacent to the AR-8305 console, if possible. If such a configuration is necessary, carefully observe the configuration in question to ensure that the electromagnetic interference does not degrade the performance.

 **Do not attach handpieces or footswitches during the Self Test or the Programming Modes.**

 **WARNING!** Use only Arthrex approved accessories. Other accessories may result in increased emissions or decreased immunity of the system. Contact your Arthrex representative for a complete list of accessories. **DO NOT** modify any accessory. Failure to comply may result in patient and/or operating room staff injury.

 **WARNING!** This equipment is NOT suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen rich or nitrous oxide environment.

Setting up a Footswitch

 **Only use footswitches developed by Arthrex specifically for the AR-8305 Synergy^{Resection} console.**

NOTE: Setup for a footswitch is the same for all models. The console detects which version is attached and allows the appropriate functions.

Insert the footswitch connector into the console's footswitch receptacle. Align the red dots on the connector and receptacle so they engage easily.

Setting up the Shaver Handpieces

 **Only use shaver handpieces and cables that have been developed by Arthrex specifically for the Synergy^{Resection} console.**

Insert the handpiece connector into the console's handpiece receptacle. Align the red dots on the connector and receptacle so they engage easily.

Setting up Other Accessory Handpieces

 **Only use handpieces and cables that have been developed by Arthrex specifically for the Synergy^{Resection} system.**

Other accessory handpieces are available with the Synergy^{Resection} system and may require an accessory cable. Specific instructions for use of an accessory are packaged with that handpiece.

Insert one end of the accessory cable connector into the cable receptacle of the accessory handpiece. Align the red dots on the connector and receptacle so they engage easily.

Connect the other end of the accessory cable to the handpiece receptacle of the console. Align the red dots on the connector and receptacle so they engage easily.

Shutdown Procedure

The AR-8305 can be safely shut down at any time by powering down the console.

Operation and Frequently Used Functions

Main Screen




Figure 5: Main Screen Display - One Tool Attached




Figure 6: Main Screen Display - Two Tools Attached

Table 13: Main Screen Display Elements

1	The main screen , divided into two sections when two tools are attached
2	Settings button (on the lower right corner in the single screen mode)
3	Speed Increase button
4	Speed Decrease button
5	Tool name
6	Footswitch , displayed in green if Footswitch Override is disabled (SFS = Standard Footswitch, GFS – Gas Footswitch)
7	Speed setting
8	Direction or Mode indicator, displayed in green if active
9	Direction and Oscillation buttons

- > The main screen is divided horizontally into two sections, one for each channel, when tools are connected to both channels.
- > Each section contains a display area for messages and five buttons.
- > The Settings button on the right configures both main screens. It is on the lower right corner in the single screen mode or between the two sections' **Speed Increase** and **Speed Decrease** buttons in dual screen mode.
- > The **tool name** is positioned as depicted in the above diagrams.
- > **Footswitch data** is positioned below the tool name.
- > The **speed setting** is displayed in large numeric digits. They are visible from a distance by the surgeon in the operating room (OR).
- > **Direction** or **mode indicators** are below the speed setting in large characters. They are visible from a distance by the surgeon in the OR.
- > **Direction** and **oscillation** buttons.
- > **Cruise control status** is positioned next to the footswitch data.

Settings Menu Main Screen

The **Settings Menu** provides access to settings and information screens. There are six buttons described in Table 14 below.




Figure 7: Settings Menu Display - Operation

Table 14: Settings Menu Display - Operation		
1	Select Language button	When the Select Language button is pressed, the software displays the Language Selection screen, as shown in Figure 7.
2	Select Oscillation Mode button	When the Select Oscillation Mode button is pressed, the software displays the Oscillation Mode screen, as shown in Figure 8.
3	Footswitch Override button	When the Footswitch Override button is pressed, the software displays the Footswitch Override screen, as shown in Figure 9.
4	Reset Defaults button	When the Reset Defaults button is pressed, the software displays the Reset Defaults screen which asks for confirmation before resetting all settings to the factory defaults, as shown in Figure 10.
5	Information button	When the Information button is pressed, the software displays the Information screen, as shown in Figure 11.
6	Home button	When the Home button is pressed, the software displays the Main screen, as shown in Figure 4 or Figure 5.

Language Selection Screen

The **Language Selection** screen displays a list of available languages with radio buttons. The software supports the following languages: English, German, French, Italian and Spanish. It also includes three navigation buttons: Menu, Home and Cancel as described in Table 15 below.




Figure 8: Language Selection Display - Operation

Table 15: Language Menu Display - Operation		
1	Select Language radio buttons	Select the language of choice.
2	Menu button	When pressed, the software stores the selected language and returns to the Menu screen.
3	Home button	When pressed, the software stores the selected language and returns to the Main screen.
4	Cancel button	When pressed, the software returns to the Menu screen without making any changes.

Oscillation Mode Selection Screen

The **Oscillation Mode** selection screen has radio buttons for available oscillation modes for each channel. This screen also includes Menu, Home and Cancel buttons as described in Table 16 below:




Figure 9: Oscillation Mode Selection Screen - Operation

Table 16: Oscillation Mode Display - Operation		
1	Channel 1	The software supports each handpiece channel with its own oscillation mode.
2	Channel 2	
3	Oscillation Mode radio buttons	Select the oscillation mode: Standard, Efficient or Aggressive.
4	Menu button	When pressed, the software stores the selected mode and returns to the Menu screen.
5	Home button	When pressed, the software stores the selected mode and returns to the Main screen.
6	Cancel button	When pressed, the software returns to the Menu screen without making any changes.

Footswitch Override Screen

The **Footswitch Override** screen has radio buttons to select the footswitch override mode for each channel.

At the end user's discretion, the default setting can be changed so that the shaver handpiece can be controlled by both the footswitch and/or the hand controls as described below. If both the footswitch and the hand controls are enabled and activated simultaneously, the footswitch control will be dominant.

In the default (enabled) mode, the buttons on the hand-control shaver handpiece will not function when a footswitch is attached to the console. When disabled, both buttons on the handpiece and the footswitch will function. This screen also includes Menu, Home and Cancel buttons as described in Table 17 below:




Figure 10: Footswitch Override Selection Screen - Operation

Table 17: Footswitch Override Display - Operation		
1	Channel 1	The software displays the current settings for footswitch override for each channel with two sets of radio buttons.
2	Channel 2	
3	Footswitch Override radio buttons	Each channel has an Enabled and Disabled radio button, with the factory default set to Enabled.
4	Menu button	When pressed, the software stores the selected mode and returns to the Menu screen.
5	Home button	When pressed, the software stores the selected mode and returns to the Main screen.
6	Cancel button	When pressed, the software returns to the Menu screen without making any changes.

Reset Defaults Screen

The **Reset Defaults** screen prompts the user to confirm resetting to factory defaults. This screen is displayed when the Reset Defaults button is pressed on the **Settings** screen.




Figure 11: Reset Defaults Screen - Operation

The software displays two buttons that allows the user to confirm resetting of factory defaults as described in Table 18 below.

Table 18: Reset Defaults Display - Operation

1	Prompt to confirm resetting to factory defaults.
2	Yes to reset to factory defaults.
3	No to exit without resetting.

Information Screen

The **Information** screen displays the model, software versions and builds for each controller, along with Menu and Home buttons as described in Table 19 below:




Figure 12: Information Screen - Operation

Table 19: Information Screen Display - Operation

1	Model, Software Versions & Builds for each controller.	
2	Menu button	Returns to the menu screen.
3	Home button	Returns to the main screen

Buttons

Most buttons send a message to the controller when activated and/or released. The **Speed Increase** and **Speed Decrease** buttons auto repeat when pressure is maintained. The selected speed is not sent to the controller until the button is released.

Buttons - Main Screen

The Main Screen is displayed in one of two configurations depending on how many handpieces are connected to the system. A single handpiece screen is displayed with a single screen while a two handpiece screen has two equally sized sections, top and bottom. A single handpiece screen has the same content as a dual handpiece screen section, with the exception of the OSC Mode button on a single handpiece screen.



Figure 13: Single Handpiece Screen Oscillation Mode Location

Table 20: Oscillation Mode - Operation	
1	Oscillation Mode button The OSC Mode button is only displayed when a single handpiece is installed and the device is in OSC mode. When displayed, it can be pressed to toggle between the different OSC modes. This button is not displayed on a dual handpiece screen.




Figure 14: Two Handpiece Buttons - Operation

Table 21: Two Handpiece Buttons - Operation

1	Channel 1 Reverse button	When the Reverse button is pressed, the handpiece motor moves in a counterclockwise direction when observed from the tip of the shaver.
2	Channel 2 Reverse button	
3	Channel 1 Oscillation button	When the Oscillation button is pressed, the handpiece motor moves in an alternating clockwise/ counterclockwise direction.
4	Channel 2 Oscillation button	
5	Channel 1 Forward button	When the Forward button is pressed, the handpiece motor moves in a clockwise direction when observed from the tip of the shaver.
6	Channel 2 Forward button	
7	Channel 1 Speed Increase button	When the Speed Increase button is pressed, the speed of the motor increases. The increase of the motor speed depends on the type of handpiece.
10	Channel 2 Speed Increase button	When the Speed Increase button is initially pressed, during the first second, the speed changes by one increment. If the Speed Increase button is still pressed after the initial second, the speed continues to increase until the limit is reached.
9	Settings button	The Settings button is not displayed if either of the two motors are in motion. When the motors are not in motion, the Settings button is displayed between the Speed Increase and Speed Decrease buttons for both sections in a dual screen mode, or in the lower right corner in a single screen mode. When the Settings button is pressed for 3 seconds, the software switches the controller into speed test mode. See Speed Test Mode.
8	Channel 1 Speed Decrease button	When the Speed Decrease button is pressed, the speed of the motor decreases. The decrease of the motor speed depends on the type of handpiece.
11	Channel 2 Speed Decrease button	When the Speed Decrease button is initially pressed, during the first second, the speed changes by one increment. If the Speed Decrease button is still pressed after the initial second, the speed continues to decrease until the limit is reached.

Speed Test Mode

The software enters the Speed Test mode when the Menu button is pressed and held for 3 seconds until a triple beep is sounded.

Speed Adjustment

The **Speed Increase** and **Speed Decrease** buttons auto-repeat when pressure is maintained. The selected speed is not sent to the controller until the button is released.

Oscillation Mode Adjustment

The **OSC Mode** button is only displayed when a single handpiece is installed and the device is in OSC mode.

When displayed, it can be pressed to toggle between the different OSC modes. This button is not displayed on a dual handpiece screen.

Accessory Handpiece Speed Limits and Deltas

The console supports individual handpiece motor-speed limits, defaults, and speed increment and decrement “deltas”. These are defined below for the supported handpieces.

Table 22: Accessory Handpiece Speed Limits and Delta

Handpiece	Forward / Reverse Speeds (RPM)			Oscillation Speeds (RPM)		
	Min.	Max.	Speed Delta	Min.	Max.	Speed Delta
Shaver	500	8000	300	500	3000	250
Drill	100	1400	See Drill Speeds Section	N/A	N/A	N/A
Sag Saw	18000	18000	N/A	N/A	N/A	N/A
SJ/NanoResection	500	6200	300	500	3000	250

Drill Speeds

The drill supports only 100, 300, 500, 900, and 1400 (RPM). It does not cycle through regular increments.

NOTE: While operating handpieces, the settings can be changed without risk of damaging the motor or console.

 **WARNING!** Do not have the device in direct contact with patient if high-frequency devices are in use or if the patient requires defibrillation.

Console Speed Controls

The Speed Increase and Speed Decrease buttons adjust the speed. Each time a button is pressed, the speed is adjusted as shown in Table 22. The commands are ignored if multiple buttons are pressed simultaneously.

To make the handpiece speed increase or decrease automatically, press the Speed Increase or Speed Decrease button for more than 1 second. The values will not exceed the minimum and maximum values defined for the attachment.

Table 23: Handpiece Accuracy

Handpiece	ACCURACY
All Shavers	±5 % OR ±100 RPM (greater of the two)
Drill	±5 % OR ±20 RPM (greater of the two)
Sag Saw	±6 %

Attachment Motor Control

Definition of Direction

The definition of direction as viewed looking from the proximal to the distal (top of the tip) of the shaver handpiece.

1	Reverse	Counterclockwise (CCW)
2	Forward	Clockwise (CW)
3	Oscillation	CW then CCW
4	Stop	No rotation

Forward, Reverse and Oscillation Modes

The buttons and pedals for the Reverse, Forward, and Oscillation modes select which direction the motor rotates. If the Reverse and Forward buttons or pedals are detected simultaneously then the motor defaults to Oscillation mode. During Oscillation mode, the motor will rotate forward, then backward. Any other combination of simultaneous button and pedal activations are ignored by the device. This operation applies to both the console front-panel buttons and the supported footswitch pedals.

Cruise Control

Use the gas pedal footswitch to engage the cruise control function. The gas pedal footswitch is available as an optional system accessory. To operate the cruise control, you must:

> Cruise Set

When you select the Cruise Control Button, the Cruise function becomes ARMED (signaled by a single audible 1 second alarm) but the LOCK IN period does not start until one of the three pedals is pressed. If a pedal is not pressed within 15 seconds then the Cruise function becomes DISARMED (signaled by a single audible 1 second alarm).

If the Cruise button is pressed on an attachment that does not allow cruise control, an error alarm occurs (3 fast tones), and displays the message “*CRUISE NOT ALLOWED*”.

> Cruise Lock

The LOCK IN period lasts for up to 15 seconds. During the LOCK IN period the operator adjusts the speed and must hold the pedal steady for 3 seconds in order for the LOCKED state to occur. Steady is defined as ± 500 RPM for forward or reverse, and ± 400 RPM for oscillation. If the pedal is held steady, the shaver console emits an audible double alarm (2 brief tones) that indicates the LOCKED state has been entered, and “*Lock*” is displayed. If the pedal was not held steady within the 15 second LOCK IN period, then the Cruise function becomes DISARMED (signaled by a single audible 1 second alarm).

> Remove Foot

The operator can move the same pedal for the next 3 seconds without leaving the LOCKED state. This is necessary, since the act of taking the foot off of the pedal moves the pedal position. Once the pedal is fully released, or after 3 seconds have elapsed, the system again responds to this pedal and disables the Cruise function.

> Disable Cruise Lock

If any other footswitch or console panel button is activated while in the LOCKED state, the Cruise function becomes DISARMED (with the exception noted in item Remove Foot).

> Disable Cruise Set

If the Toggle, Cruise Control, or Front Panel buttons are pressed at any time after the Cruise Control button is selected (i.e. in the ARMED or LOCKED state) then the Cruise function becomes DISARMED.

Toggle Function

A Toggle switch is available on both the standard and the gas pedal footswitches. The toggle function is performed by pressing a button on the footswitch. The toggle function is available when one footswitch and two handpieces are connected to the Synergy^{Resection} console. The toggle function allows the operator to control the handpieces in channel 1 and channel 2, without moving the footswitch connection.

At power on, or at first connection, the footswitch controls the attachment which is on the same channel as the one the footswitch is connected to. When the user presses the Toggle switch, a beep indicates the footswitch has been toggled to the attachment on the opposite channel. The handpiece on the same channel as the footswitch is controlled by the front panel buttons only.

When the user presses Toggle again, the footswitch returns control to the handpiece on the same channel as the footswitch. The channel on the opposite side as the footswitch is controlled by front panel buttons only.

When controlled by the footswitch, the **HC shaver** is controlled by the Footswitch Override setting. When not controlled by the footswitch, the handpiece buttons and control panel are enabled.

If there is a connected attachment for each installed footswitch, and the attachment supports oscillation mode, the toggle button switches to the next oscillation mode in this sequence: efficient mode, aggressive mode, standard mode and then back to efficient mode.

If two footswitches are installed with only one connected attachment, the footswitch without a connected attachment is inactive.

If the Toggle switch is activated when only one attachment is connected and this attachment does not support oscillation modes, an error beep sounds and “TOGGLE NOT ALLOWED” message is temporarily displayed.

Table 24: Toggle Function

Installed Handpieces	Installed Footswitches	Action	Beep type
0	1	None	None
0	2	None	None
1	1	Change oscillation mode	Single
1	2	Change oscillation mode (footswitch on installed HP channel)	Single
		None (footswitch with no installed HP)	None
2	1	Move control to next channel	Single
2	2	Change oscillation mode	Single

Footswitches

 **The footswitch cable connects and locks to the console to prevent accidental separation during use. To avoid damage, disconnect the footswitch by pulling on the cable connector shell (plug) only.**

Standard Footswitch

The Speed Increase button adjusts the speed in the positive direction only. Each time the button is pressed, the speed adjusts as shown in Table 22. When the device is at maximum speed setting, the user can return to minimum speeds by pressing the Speed Increase button one more time. The commands are ignored if multiple buttons are pressed simultaneously.

If the Speed Increase button is continuously pressed for more than 1 second then the speed auto-increments accordingly. The speed setting will not wrap when auto-incrementing. Release the Speed Increase button and then press it again to get the speed setting to cycle to the minimum setting.

Gas Pedal Footswitch

The Forward, Reverse, and Oscillation pedals linearly adjust the speed. When the pedal is pressed, it returns an analog signal that is converted to a corresponding RPM between the minimum and maximum speed of the attachment connected to the footswitch.

Multi-Function Footswitches

Wireless and corded multi-function footswitches can operate as either a standard footswitch or as a gas pedal footswitch, as described in the preceding sections. Function selection is made by pressing the button on the multi-function footswitch labeled “GFS” or “SFS”. These footswitches are available as optional accessories. For more information on the corded footswitches, please see DFU-0203-XX. For more information on the Synergy^{Resection} system wireless footswitch (AR-8315W), please see DFU-0320-XX. For the Synergy^{Resection} system wireless footswitch Quick Start Guide, please see DFU-0323-SUB.

Accessory Handpieces

 **The accessory handpiece cable connects and locks to the console to prevent accidental separation during use. To avoid damage, disconnect the footswitch by pulling on the cable connector shell (plug) only.**

Cleaning and Disinfecting

Console AR-8305



The AR-8305 console is provided non-sterile and must not be sterilized.

The AR-8305 console can be cleaned/disinfected using a cloth and commercially available surfactants or surface disinfectants only. Use a commercially available cleaning disinfectant for surfaces. Alkaline agents with a pH value > 8 are to be avoided. Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.

 The AR-8305 must not be submersed in any liquid.

Always place the main power switch to “Off (O), position” and disconnect the power before cleaning the AR-8305 console.

 **Always comply with the instructions issued by the manufacturer of the cleaning disinfectant regarding concentration, exposure times, temperature and material compatibility.**

 **NEVER clean the console with liquid cleaners. If there is dust, remove it with dry compressed air.**

 **Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.**

Footswitches



The footswitch is provided non-sterile and must not be sterilized.

To clean and disinfect, unplug the cord from the console and add the protective cap to the connector. Use a commercially available cleaning disinfectant for surfaces. Alkaline agents with a pH value > 8 are to be avoided. Clean the footswitch with an enzymatic cleaner without subsequent acid neutralization.

For additional information, please refer to DFU-0203-XX.

Accessory Cable

Clean the cable with an enzymatic cleaner without subsequent acid neutralization. Use a commercially available cleaning disinfectant for surfaces. Alkaline agents with a pH value > 8 are to be avoided. Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.

Rinse the cable thoroughly after cleaning.

After cleaning, disinfect the cable with a commercially available surface disinfectant.

Thoroughly rinse the cable under lukewarm distilled water or preferably demineralized sterile water. See the **Sterilization** section for more information.

 **Always comply with the instructions issued by the manufacturer of the disinfectant.**

 **NEVER allow the console plug pins to have any contact with liquids. If there is dust or moisture on the receptacles, remove it with dry compressed air. ONLY dry connectors may be plugged into the console.**

Accessory Handpieces

 **Refer to the Directions for Use (DFU-0154-XX) for detailed shaver handpiece cleaning and sterilization instructions. This can be obtained from the eDFU website at arthrex.com/eDFU, or by contacting your local Arthrex representative.**

Remove any accessories in the device. This includes shaver blades, saw blades or drill bits.

If the handpiece has suction control, make sure it is fully open.

Clean the handpiece with an enzymatic cleaner without subsequent acid neutralization.

If the handpiece has suction control, clean the connection with a cleaning brush. Rinse the handpiece thoroughly after cleaning.

Disinfect the handpiece with a commercially available disinfectant.

Thoroughly rinse the handpiece under lukewarm distilled water or preferably demineralized sterile water.

See the **Sterilization** section for more information.

 **Always comply with the instructions issued by the manufacturer of the disinfectant.**

 **NEVER place the shaver handpiece in Cidex® or other aldehyde disinfectant solutions.**

 **NEVER allow the console plug pins to have any contact with liquids. Remove dust or moisture with dry compressed air. Only DRY connectors may be plugged into the console.**

Sterilization

Sterilization capabilities, cleaning, disinfecting, handling, and storage of instrumentation are the responsibility of qualified facility and/or user personnel. Qualified personnel must still properly clean and disinfect the instruments prior to sterilization.

 **Refer to the Directions for Use (DFU-0154-XX) for detailed shaver handpiece cleaning and sterilization instructions. This can be obtained from the eDFU website at arthrex.com/eDFU, or by contacting your local Arthrex representative.**

 **After sterilization in the autoclave, let the accessory device cool down slowly. NEVER use cold water to cool the handpieces. This will damage the electronic components and seals.**

 **Liquid on the cable connector of the accessory device can damage the device. Before connecting the cable, ensure the receptacles are clean and dry.**

 **WARNING!** After autoclaving, the accessory devices are VERY HOT. Handle with care to avoid burns.

The following special instructions are to be observed when sterilizing accessory handpieces.

 **Open or remove the suction control completely to ensure proper circulation, if applicable.**

 **Open the chuck completely to ensure proper circulation, if applicable.**

 **Remove all accessories, including shaver blades, drill bits and saw blades before sterilizing a handpiece.**

Transmissible Spongiform Encephalopathy (TSE) Agents

It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy (TSE) Agents.

The agents for transmission of Creutzfeldt-Jakob disease are believed to be resistant to normal processes of disinfection and sterilization. Therefore, the normal processing methods of decontamination and sterilization as described above may not be appropriate where CJD transmission is a risk.

In general, the tissues that come into contact with orthopedic surgical instruments are those of low TSE infectivity. However, take particular precautions when handling instruments that have been used on known, suspected, or at risk patients.

References:

ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in health care facilities

Maintenance

Regular and proper maintenance of your Synergy^{Resection} system is the best way to protect your investment and avoid non-warranty repairs.

Recommended care and handling of the Synergy^{Resection} system includes proper day- to-day operation, cleaning, and sterilization which are extremely important to ensure safe and efficient operation. It is important to visually inspect the foot pedal and handpiece, cable and connectors before each use.

Your authorized Arthrex service department is the most knowledgeable about the Arthrex medical resection systems and/ foot pedal and handpiece, and will provide competent and efficient service. Any services and/or repairs done by any unauthorized repair facility may result in reduced performance of the instruments or instrument failure.

Life Expectancy

Life expectancy for the console is approximately 5 years under normal use and standard of care.

Periodic Maintenance

The product should be inspected prior to and after each use to ensure that the foot pedal, handpiece, cable, strain relief, overmold, or connector contacts are not damaged or worn. If it becomes necessary to return the foot pedal and / or handpiece to Arthrex for service, please disinfect the foot pedal and/or sterilize the handpiece before shipping.

Technical Support

For assistance in using the products identified in this *User Manual*, contact an Arthrex representative or contact the **Arthrex Technical Support Hotline** at ☎+1 888 420 9393, Monday through Friday from 8:00 AM to 8:00 PM EST; at ☎ +49 89 909005-0 or techsupport@arthrex.de from 8:00 AM to 5:00 PM CET

How to Display the Software Version

Technical Support may request the software version of the AR-8305. Follow these instructions to display the software version.

1. Power On the AC mains power switch Figure 1 [3] on the AR-8305.
2. The software version is displayed on the TPVD during the power-up sequence.

Additional Technical Information

Contact your Arthrex representative if require more comprehensive technical information. Circuit diagrams, component part lists, descriptions, calibration instructions, or other information will be provided upon request by Arthrex APPROVED SERVICE PERSONNEL.

Troubleshooting

Refer to Table 25 for device troubleshooting if problems occur after cleaning, transporting or changing operating staff.

Table 25: Troubleshooting: Faults, their Causes and Solutions

Fault	Cause	Solution
Unit does not power on.	1. Not receiving power from the wall receptacle or the power strip. 2. Blown fuses. 3. Console has an internal power failure.	1. Check the main power plug and wall receptacle. 2. Check the fuses. 3. Send in for repairs.
Handpiece does not work.	1. Handpiece is not connected to the correct channel. 2. Not fully or properly connected or connector is damaged. 3. Liquid in the receptacle or connector. 4. Connector is damaged. 5. Footswitch is defective. 6. Footswitch is functional. 7. Handpiece is too hot after sterilization. 8. Handpiece over-current is detected. 9. Defective handpiece.	1. Check the connector is plugged into the correct channel receptacle. 2. Check both connectors are plugged in properly. 3. Dry thoroughly. 4. Replace the cable. 5. Operate the device using the console unit. 6. Replace the handpiece. 7. Allow it to cool down as indicated in the WARNING in the Sterilization section. 8. Allow 30 seconds for the console to reset the affected channel. 9. Send in for repairs.
Motor runs but shaver blades/burrs do not work.	1. Shaver blade/burr incorrectly inserted. 2. Damaged hub. 3. Bent sleeve.	1. Insert the shaver blade so that it engages. Fully close the locking mechanism. 2. Replace the shaver blade/burr. 3. Replace the shaver blade/burr.
No suction.	1. Clogged suction connection. 2. Clogged shaver blade or burr. 3. Closed suction control. 4. Damaged suction control.	1. Rinse with a syringe. 2. Clean with a brush. 3. Open suction control. 4. Send in for repairs.
System runs with the control unit but not with the footswitch.	1. Not fully, properly connected or damaged connector. 2. Liquid in receptacle or connector. 3. Defective footswitch.	1. Check the connector from the footswitch to the control unit. 2. Dry thoroughly. 3. Send in for repairs.
Connected handpiece not shown on display.	1. Not fully or properly connected or connector is damaged. 2. Defective handpiece or connector. 3. Connector is damaged. 4. Defective handpiece.	1. Check that both connectors are plugged in properly. 2. Check the connector with another handpiece. 3. Replace. 4. Send in for repairs.

If the problems persist, disinfect the Synergy^{Resection} system and send in to Arthrex using the original packaging. Always send the corresponding handpiece together with the console, footswitch and cable. Please enclose a brief explanation of the detected malfunction. Refer to the **Repair Policy** section for more information.

Troubleshooting Interference with Other Devices

Try one or more of the following to correct interference:

- › Reorient or relocate the receiving device.
- › Increase the distance between devices.
- › Connect the device to an outlet on a different circuit than the other device(s).
- › Consult the manufacturer or field service technician for the receiving device for assistance.

Repair Policy

Contact Arthrex for a Return Authorization Number and instructions prior to returning the device.

End of Life, Environmental Directives



WEEE Directive **[2012/19/EU]** on Waste Electrical and Electronic Equipment

The Directive on Waste Electrical and Electronic Equipment obliges manufacturers, importers, and/or distributors of electronic equipment to provide for recycling of the electronic equipment at the end of its useful life.

Do not dispose of WEEE in unsorted municipal waste.

The WEEE symbol on the product or its packaging indicates that this product must not be disposed of with other waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of Waste Electrical and Electronic Equipment. The separate collection and recycling of your waste equipment at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your medical electronic equipment at the end of its useful life for recycling, please contact Arthrex Customer Service Department.

Electromagnetic Emissions

The AR-8305 Synergy^{Resection} system is intended for use in the electromagnetic environment specified below. The customer or the user of the AR-8305 Synergy^{Resection} system should assure that it is used in such an environment.

Table 26: Guidance and Manufacturer's Declaration – Electromagnetic Emissions

Emissions test	Compliance	Electromagnetic environment – guidance
Radio Frequency (RF) emissions CISPR 11	Group 1	AR-8305 Synergy ^{Resection} system uses RF energy only for its internal function. Therefore its RF emissions are very low and not likely to cause any interference in nearby equipment. The AR-8305 Synergy ^{Resection} system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 27: System Cables

Type	Use	Shielded	Ferrite	Maximum Length
Power Cords	Supply Line Power to the Console	No	No	3.048 m (10 feet)

The AR-8305 Synergy^{Resection} system is intended for use in the electromagnetic environment specified below. The customer or the user of the AR-8305 Synergy^{Resection} system should ensure that it is used in such an environment.

Table 28: Guidance and Manufacturer's Statement - Electromagnetic Immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines 100kHz Cycling Frequency 100	± 2 kV for power supply lines ± 1 kV for input/output lines 100kHz Cycling Frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	U _T = 0%, 0.5 cycle (0, 45, 90, 135, 180, 225, 270, and 315°) U _T = 0%; 1 cycle U _T = 70%; 25/30 cycles (@0 degrees) U _T = 70%; 250/300 cycles	U _T = 0%, 0.5 cycle (0, 45, 90, 135, 180, 225, 270, and 315°) U _T = 0%; 1 cycle U _T = 70%; 25/30 cycles (@0 degrees) U _T = 70%; 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the AR-8305 Synergy ^{Resection} system requires continued operation during power mains interruptions, it is recommended that the AR-8305 Synergy ^{Resection} system be powered from an uninterruptible power supply.
Power frequency (50 - 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m @ 50 - 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the AC mains voltage prior to application of the test level.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Model AR-8305, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3.5 / V1] \sqrt{P} = 1.17 \sqrt{P}$ $d = [3.5 / V1] \sqrt{P} = 1.17 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = [7 / E1] \sqrt{P} = 2.33 \sqrt{P} \text{ 800 MHz to 2.7 GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model AR-8305 is used exceeds the applicable RF compliance level above, the Model AR-8305 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model AR-8305.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Model AR-8305

The Model AR-8305 Synergy^{Resection} system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model AR-8305 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model AR-8305 as recommended below, according to the maximum output power of the communications equipment.

Table 29: Guidance and Manufacturer's Statement – Recommend Separation Distance between portable and mobile RF communications equipment and the Model AR-8305

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Contact Information

Any serious incident that occurs in relation to the device should be reported to the manufacturer and to the health authority where the incident occurred.

 Arthrex, Inc.	1370 Creekside Blvd. Naples, FL 34108, USA ☎+1 800 934-4404 arthrex.com arthrex.com/eDFU
Customer Service	☎+1 866 267 9138
Technical Support	☎+1 888 420 9393
 Arthrex GmbH	Erwin-Hielscher-Strasse 9 81249 München, Germany ☎+49 89 90 90 05-0 arthrex.de info@arthrex.de
 Arthrex Distribution Hub EMEA B.V.	Ampèrestraat 9 5928 PE Venlo, Netherlands ☎+31 88 712 9800 arthrex.nl info@arthrex.nl
 Arthrex Ltd.,	Unit 1 Bessemer Park Shepcote Lane Sheffield S9 1DZ United Kingdom ☎+44 0 114 232 9180
 confinis ch-rep ag,	Hauptstrasse 16, 3186 Düringen, Switzerland ☎+41 26 494 8 494
Manufacturer's Australian Sponsor Arthrex Australia Pty Ltd	Suite 501, 20 Rodborough Road Frenchs Forest, NSW, 2086 Australia ☎+1 800 950 637 arthrex.com.au csAU@arthrex.com

Synergy^{Resection}™ is a trademark of Arthrex, Inc. Other trademarks are owned by their respective owners.

This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the "Arthrex U.S. Product Warranty" section of the Arthrex, Inc. website, found at www.arthrex.com whose provisions are incorporated herein by reference.



R_X ONLY

