# Certofix<sup>®</sup> Duo HF

<b>GB</b>	Central venous catheter set for Seldinger technique.	<b>10</b>
<b>DE</b>	Zentralvenöses Katheterbesteck in Seldingertechnik.	<b>12</b>
<b>BG</b>	Набор за централна венозна катетеризация по метода на Селдингер.	<b>14</b>
<b>CN</b>	用于 Seldinger 技术的中心静脉导管。	<b>16</b>
<b>CZ</b>	Souprava s katétrem ke katetrizaci centrální žíly Seldingerovou metodou.	<b>18</b>
<b>DK</b>	Centralvenøst katetersæt i Seldingertechnik.	<b>20</b>
<b>EE</b>	Komplekt tsentraalvenoosseks kateetri paigalduseks (Seldingeri tehnika).	<b>22</b>
<b>ES</b>	Set del catéter venoso central en técnica Seldinger.	<b>24</b>
<b>FI</b>	Seldingerin tekniikan keskilaskimokatetrivälineet.	<b>26</b>
<b>FR</b>	Set de cathétérisme veineux central pour méthode de Seldinger.	<b>28</b>
<b>GR</b>	Κεντροφλεβικό σετ καθετήρ α σε τεχνική Seldinger.	<b>30</b>
<b>HR</b>	Set za centralni venski kateter za uvođenje Seldingerovom metodom.	<b>32</b>
<b>HU</b>	Centrális vénás katéterkészlet Seldinger-módszer szerinti bevezetéshez.	<b>34</b>
<b>ID</b>	Set kateter vena pusat untuk teknik Seldinger.	<b>36</b>
<b>IT</b>	Set catetere venoso centrale conformemente alla tecnica Seldinger.	<b>38</b>
<b>KZ</b>	Сельдингер әдісі бойынша орталық қүретамырлардың катетеризациясына арналған жинақ.	<b>40</b>
<b>LT</b>	Centrinės venos kateterio rinkinys Seldingerio metodu.	<b>42</b>
<b>LV</b>	Centrālais venozais katetrs, iestatīts Seldinger metodei.	<b>44</b>
<b>NL</b>	Centraalveneus katheterbestek met Seldingertechniek.	<b>46</b>
<b>NO</b>	Sentralvenøst katetersett for Seldingerteknikk.	<b>48</b>
<b>PL</b>	Zestaw do cewnikowania żył centralnych metodą Seldingera.	<b>50</b>
<b>PT</b>	Conjunto de cateter venoso central para técnica de Seldinger.	<b>52</b>
<b>RO</b>	Set cateter venos central pentru tehnica Seldinger.	<b>54</b>
<b>RS</b>	Set centralnih venskih katetera za Seldingerovu tehniku.	<b>56</b>
<b>RU</b>	Набор с центральным венозным катетером для постановки по методу Сельдингера.	<b>58</b>
<b>SE</b>	Centralvenös katetersats i Seldingerteknik.	<b>60</b>
<b>SI</b>	Centralnovenski kateter (set) v Seldingerjevi tehniki.	<b>62</b>
<b>SK</b>	Súprava s katétrom na katetrizáciu centrálnej žily Seldingerovou metódou.	<b>64</b>
<b>TH</b>	ชุดสายสวนหลอดเลือดดำส่วนกลางสำหรับเทคนิค Seldinger	<b>66</b>
<b>TR</b>	Seldinger tekniği ile kullanım için santral venöz kateter seti.	<b>68</b>
<b>UA</b>	Набір для катетеризації центральних вен за методом Сельдингера.	<b>70</b>
<b>VN</b>	Bộ ống thông tĩnh mạch trung tâm cho kỹ thuật Seldinger.	<b>72</b>

**B | BRAUN**

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







Double-lumen cava catheter set for high flow rates and acute dialysis for catheterisation of the superior vena cava using the Seldinger technique. Opaque catheter made of polyurethane and guide wire with flexible J-tip. Includes connecting cable for intra-atrial ECG lead via the guide wire for simultaneous verification of catheter position.

#### Flow rates

Catheter size	Distal (ml/min)	Proximal (ml/min)
9F/ 20 cm 13/13G	95	120
12F/ 15cm 11/11G	260	260
12F/ 20 cm 11/11G	230	230

## GB Instructions for use

### Contents

1. Valve needle (A)
2. Guide wire (B) with length markings and flexible J-tip (r=3mm)/straight tip, in dispenser
3. Scalpel (I)
4. Dilator (C)
5. Double-lumen, opaque catheter made of polyurethane with soft tip (see packaging for length) (D):
  - With Safsite® valves
  - Catheter markings to verify position
  - Colour-coded Luer Lock connections
  - Fixation wing (G) at the catheter fork to secure the catheter
  - Attachable, adjustable fixation wing (E) to secure the catheter by suture at the insertion site
6. Connecting cable for intra-atrial ECG lead (J)
7. Attachable clip (F)
8. Omnifix® syringe, 5 ml (H)

### Materials used

ABS, Brass, EP, HD-PE, Ink, IR, LD-PE, Nitinol, PA, PA6 dyed, PC, PE, PE-HD, PP, PS, PUR, PVC-white, SBS, SI, Silicon oil, SIR, Stainless steel

### Indications

Catheterisation of the superior vena cava using the Seldinger technique for short-term (<=30 days) infusion and volume therapy or parenteral nutrition, for administration of highly osmolar or very veinirritating solutions, for continuous or intermittent monitoring of the central venous pressure, for blood sampling, or when peripheral venous puncture is not possible in state of shock, in patients with injured extremities or no detectable peripheral veins. Acute haemofiltration, detoxication procedures, other pheresis procedures.

### Contraindications

Inflammation of the skin at the site of puncture; clotting disorders e.g. in therapy with anticoagulants; anatomical anomalies, e.g. enlarged struma, tumours in the neck regions, very severe pulmonary emphysema and postoperative changes at the site of puncture.

**ECG lead** Do not record ECG lead during defibrillation, cardioversion or HF surgery.

### Risks

Haematoma at the site of puncture, catheter sepsis, pneumothorax, haemothorax, infusion hydrothorax or chylothorax due to incorrect puncture and incorrect catheter position, cardiac arrhythmias due to incorrect intracardiac placement of the catheter, atrial rupture, endocarditis due to mechanical irritation, arterial injuries due to incorrect punctures, catheter-induced thrombosis and thrombophlebitis of the superior vena cava, thromboembolism, thoracic duct injuries, damage to the brachial plexus, damage to the phrenic nerve.

### Warnings

- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Do not subject the Certofix catheter to pressures in excess the operating pressure of 1.2 bar or the maximum pressure limit of 2 bar in case of emergency.
- Use strictly aseptic techniques.
- To avoid air embolism make puncture lege artis with the patient in head-down position.
- Use X-ray to check the position of the catheter tip. If in doubt

about the position, inject contrast medium or record intra-atrial ECG lead.

### - When recording an ECG lead, please note:

- Use only in rooms that are protected against electrostatic charge.
- Wear antistatic operating theatre shoes.
- Observe relevant safety regulations (e.g. VDE 0750, VDE 0107 or IEC-publications) and national specifications and deviations.
- Use an ECG monitor with floating input (CF type) only.
- Earth the monitor for additional equipotential bonding.
- Ensure that all connections are securely fastened.
- Ensure that the guide has no contact to external conductive parts.

- Secure the catheter adequately.
- Daily strictly aseptic catheter care.

- Change contaminated or soaked dressings when required.

- Use different sites for blood sampling and administration of infusions or medication to reduce the risk of infection.

- Rinse catheter thoroughly with physiological saline after transfusions or blood sampling. Fill unused lumens with heparin saline solution.

- To prevent air embolisms, fasten the connection to the catheter securely and use only infusion sets with Luer Locks.

- Do not withdraw guide wire against the needle level to reduce risk of damages or possible shear of guide wire. If insertion is difficult, perform further actions in consideration of the given situation with regards to the clinical risks and benefits for patient.

- Alternatively, first remove the puncture needle and then insert the dilator in its place. This considerably reduces the risk of damaging the guide wire by the dilator when the guide wire is pulled back. If it is difficult to remove the guide wire from the catheter, remove both the guide wire and the catheter.

- If the vein is known to be constricted, it can be of help to turn the guide wire in the needle or remove the guide wire and insert the straight guide wire tip.

- Do not bring the PUR catheter into contact with acetone as acetone dissolves the material causing the catheter to become porous and to leak.

- Cardiac arrhythmias may occur during an MR examination with an inserted catheter. There should be no electrolyte solutions in the catheter during an MR examination.

### Acute dialysis:

Do not use Safsite® valves for haemofiltration/haemodialysis.

- Strict indication in patients with increased tendency to bleed and risk of bleeding.

- Rinse catheter thoroughly with physiological electrolyte solution after transfusions or blood sampling.

- Ensure that catheter is properly connected as otherwise there is a risk of air embolism. Use only infusion sets with Luer Locks.

Do not resterilise.

Use only if the packaging is intact.

### Directions

#### CAUTION:

If an ECG lead is to be used to check the position of the catheter and sinus rhythm can be recorded, first connect the patient to a CF type monitor with approval to record intra-atrial leads (see below) (type: CF) using a Certodyn® Universal Adapter or a switchable ECG patient cable.

1. When making a puncture, always wear sterile clothing with face mask, cap and gloves. After thoroughly cleaning and disinfecting the skin (use alcohol or iodine based cleaning agents, e.g. Braunol® or Softasept® N), cover puncture site with a sterile fenestrated drape.

Then fill the proximal channel of the catheter with physiological saline to verify the patency of the lumens.

2. Puncture the vein using the valve needle included and with a syringe attached for aspiration.

#### CAUTION:

Do not bend the introducer needle, because it may cause difficulty in advancing the guide wire through the needle or in removing the needle from the guide wire. When the needle is already bended, stop using it.

3. Guide wire in dispenser. Remove the protective cap from the dispenser. Then attach the dispenser with the guide wire to the side of the Seldinger needle. Use your thumb to advance the guide wire into the vein. Use the length markings to check the insertion depth.

4. When the guide wire is in the desired position, remove the dispenser and the cannula, retain-

- ing the position of the guide wire.
5. Dilate the puncture track. This is done by holding the extracorporeal section of the guide wire and advancing the short dilator – rotating it slightly – over the extracorporeal part of the guide wire and into the puncture site. Then remove the dilator from the puncture track, retaining the position of the guide wire.
  6. Advance catheter over the extracorporeal part of the guide wire up to the puncture site. Then hold the end of the guide wire and advance the catheter up to the desired position, turning it slightly
  7. The intravascular catheter length can be read from the permanent length markings.
  8. To verify the position of the catheter tip by ECG, first position the tip of the catheter and the guide wire at the same position. (The first wide marking on the guide wire protrudes out of the catheter connection. Now secure the clamp of the ECG connecting cable to the guide wire directly behind the catheter connection. Insert the connector of connecting cable into the CERTODYN® UNIVERSAL ADAPTER and switch over to the heart symbol (see separate instructions for connection of CERTODYN® UNIVERSAL ADAPTER). Then advance the catheter with the guide wire towards the right atrium under ECG control. If the ECG shows a higher than normal P wave, it can be assumed that the catheter tip is in an intra-atrial position (intra-atrial ECG lead). If the catheter with the guide wire is pulled back, the P wave returns to normal (→ transition point). The correct catheter position in the superior vena cava is reached when the catheter and the guide wire are pulled back a further 2–3 cm. Then remove the guide wire.
- CAUTION:**  
With a 15 cm catheter the catheter tip usually does not reach the atrium. In this case advance the guide wire under ECG control until a higher than normal P-wave appears. Thus, the correct direction of the catheter is at least checked.
9. Secure the catheter to the skin using the integrated fixation wing (not included in Econoline sets). It can also be secured directly at the puncture site with the moveable fixation wing. Clip can be used to secure the slit wing on the catheter.
  10. There are pressure clamps on the extension lines. These are used to close the lumens when infusions or transfusions are changed. If SafSite® is connected to the catheter, it is not necessary to close the catheter with the clamps as SafSite® automatically closes the infusion line if disconnected. However, ensure that SafSite® is securely connected to the catheter.

**CAUTION:**

Always use fixation wings to secure the catheter tubing as the catheter may otherwise be torn apart and flushed into the vein.

**Haemofiltration/haemodialysis:**

Disconnect all SafSite® valves from the central venous catheter.

Connect the arterial tube segment of the dialysis unit to the proximal channel of the catheter (red clamp). Connect the venous tube segment of the dialysis unit to the medial channel of the catheter (blue clamp). Then perform extracorporeal filtration as usual. After dialysis the lumens of the catheter can be rinsed and the catheter used again as an infusion catheter.

**CAUTION:**

Always keep unused channels of the catheter filled with heparin saline solution. Fill the channel directly after use with heparin saline solution and replace solution frequently if necessary.

**CAUTION:**

When sampling blood through the proximal lumen, ensure that the other lumens are kept closed to avoid incorrect laboratory results due to addition of infusions.

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B. Braun Melsungen AG