Medtronic

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High Range (HR-ACT) Low Range (LR-ACT) Recalcified (RACT)

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High Range (HR-ACT) Low Range (LR-ACT) Recalcificeret (RACT)

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 **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

Trademarks may be registered and are the property of their respective owners.
Les marques commerciales mentionnées peuvent être déposées et appartiennent à leurs propriétaires respectifs.
Marken sind eventuell eingetragen und sind das Eigentum ihrer jeweiligen Inhaber.
Las marcas comerciales pueden estar registradas y pertenecen a sus respectivos propietarios.
Handelsmerken kunnen zijn geregistreerd en zijn het eigendom van de desbetreffende eigenaars.
I marchi di fabbrica possono essere registrati e sono di proprietà dei rispettivi titolari.
Varemerkene kan være registrert og tilhører de respektive eierne.
Tavaramerkit voivat olla rekisteröityjä, ja ne ovat omistajiensa omaisuutta.
Varumärken kan vara registrerade och tillhör respektive ägare.
Τα εμπορικά σήματα μπορεί να είναι καταχωρημένα και είναι ιδιοκτησία των αντίστοιχων κατόχων τους.
Varemærker kan være registrerede og tilhører deres respektive ejere.
As marcas comerciais podem ser registadas e são propriedade dos respectivos detentores.
Ticari markalar tescilli olabilir ve ilgili sahiplerinin mülkiyetindedir.

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf der Verpackung / Explicación de los símbolos que aparecen en la etiqueta del envase / Verklaring van de symbolen op het verpakkingsetiket / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkausetiketin symboleiden selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στη σήμανση της συσκευασίας / Forklaring af symboler på emballagens mærkater / Explicação dos símbolos na documentação da embalagem / Ambalaj belge ve etiketlerindeki sembollerin açıklamaları

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1 Intended Use

For use in determining the activated clotting time of whole blood.

For in vitro diagnostic use.

2 Summary

The activated clotting time test was originally described by Hattersley in 1966.¹ The test and its variations have been widely used to monitor the anticoagulant effect of heparin.^{2,3,5,6,7,12} Medtronic Activated Clotting Time Cartridges utilize either fresh whole blood or citrated whole blood. Three types of Medtronic Activated Clotting Time Cartridges are available. The High Range Activated Clotting Time Cartridge (HR-ACT) and Low Range Activated Clotting Time Cartridge (LR-ACT) are for use with fresh drawn whole blood. The Recalcified Activated Clotting Time Cartridge (RACT) is for use with citrated whole blood.

3 Principles

The activated clotting time test measures the clotting time of fresh whole blood. Hattersley's original procedure used diatomaceous earth as the activating agent. Medtronic Activated Clotting Time Cartridges use kaolin suspended in a buffer as an activating agent.

The end point of the test is the detection of clot formation. The Medtronic ACT Plus Automated Coagulation Timer detects clot formation by measuring the rate of fall of the plunger contained in each cartridge channel. The plunger assembly falls rapidly through an unclotted sample. A fibrin web formed during clotting impedes the rate of plunger assembly descent. Clot formation is detected by a photo optical sensor located in the actuator assembly of the instrument. The ACT Plus instrument performs simultaneous duplicate tests using a two-channel cartridge. The clotting times for each channel and the average are displayed in seconds following completion of the test.

4 Reagents

Table 1 describes the cartridge reagents (100 µL per channel).

Table 1.

Cartridge Type	Reagent
HR-ACT	12% Kaolin ^c 0.05M CaCl ₂ HEPES buffer ^a Sodium Azide ^b
LR-ACT	0.75% Kaolin ^c 0.0025M CaCl ₂ HEPES buffer Sodium Azide
RACT	2.2% Kaolin ^c 0.05M CaCl ₂ HEPES buffer Sodium Azide

^a HEPES (hydroxyethyl-piperazine-ethanesulfonic acid) buffer

^b Sodium Azide is a bacteriostatic agent.

^c Kaolin concentration may vary slightly.

5 Storage

Room Temperature

ACT cartridges are designed for room temperature storage from 2°C to 25°C (35.6°F to 77°F) until the expiration date on the cartridge tray label. Cartridges should be stored in the original tray to reference the expiration date. Cartridges should be disposed of if the expiration date is unknown. Do not use if the cartridge appears discolored, cracked or disfigured; the reagent appears to have evaporated or to have been contaminated; or the Use-By Date has passed.

Note: Do not use if there is evidence of freezing.

6 Instrument

Medtronic Activated Clotting Time Cartridges are designed specifically for use in the ACT Plus instrument. For operating instructions, consult the ACT Plus instrument operator's manual.

7 Cartridges

The following Medtronic Activated Clotting Time Cartridges are available for use with the ACT Plus instrument (Table 2).

Table 2.

Cartridge	Sample Type	Sample Volume per Channel	Intended Use
HR-ACT	Fresh Whole Blood	400 µL	Cardiovascular OR Cardiac catheterization lab Heparin levels of 1.0 u/mL or higher
LR-ACT	Fresh Whole Blood	200 µL	Therapeutic heparinization Heparin levels of 0.0 to 1.5 u/mL
RACT	Citrated Whole Blood	200 µL	Therapeutic heparinization Heparin levels of 0.0 to 1.5 u/mL

8 Specimen Collection and Preparation

Blood samples that are improperly drawn and handled can yield erroneous coagulation test results; samples must be collected carefully, and handling techniques must be consistent.

If blood is obtained by venipuncture, the venipuncture should be clean and the first several milliliters of blood discarded to avoid contamination of the sample with tissue thromboplastin. Blood should not be taken from heparinized indwelling catheters or other anticoagulated lines without thoroughly flushing the line.

HR-ACT and LR-ACT

When testing fresh whole blood, samples should be tested as quickly as possible after collection. Baseline samples, which contain no anticoagulant, should be tested within 60 seconds of collection. Heparinized samples should be tested within 2 minutes of collection.

RACT

When testing citrated whole blood, collect samples in 3.2% sodium citrate using a whole blood-to-anticoagulant ratio of 9:1. Citrated samples should be tested within 6 hours when stored at room temperature.

Test stored samples from the same patient within the same timeframe. For example, if a patient sample is tested 15 minutes following collection, then subsequent samples for that patient should also be tested 15 minutes following collection.

9 Procedure

This product is intended for use only by trained individuals.

Materials supplied:

- High Range Activated Clotting Time Cartridge (HR-ACT)
- Low Range Activated Clotting Time Cartridge (LR-ACT)
- Recalcified Activated Clotting Time Cartridge (RACT)

Materials required, not supplied:

- ACT Plus instrument
- Syringe with 19-ga, 1 7/16-in (36-mm), sterile, single use, blunt needle or ACT Plus Easy Fill Accessory
- Fresh whole blood or citrated whole blood sample, depending on the test performed

Steps:

1. From the Main Menu, select "Cartridge Type" and press until the correct cartridge is displayed. Press "Enter" to confirm selection.
Note: Lot numbers and expiration dates for cartridges and controls must be entered prior to performing a test.
2. From the Main Menu, select "Enter ID". Enter the **Patient ID** and **User ID** numbers. Press "Enter" to confirm.
3. Prewarm the cartridge before collecting the test sample.
 - Minimum prewarm time: 3 minutes
 - Maximum prewarm time: 12 hours
4. Tap the cartridge to resuspend the kaolin reagent.
5. Collect the appropriate volume of sample using the sample collection guidelines. Refer to Section 8. Specimen Collection and Preparation.
6. Fill each cartridge channel (Figure 1) with the appropriate sample volume. Refer to Table 2.

Fill line indicators are located on each cartridge channel.

Note: Do not allow blood to drip onto the cartridge flag or cartridge walls.


Figure 1. Filling the Cartridge

7. Insert the cartridge into the ACT Plus instrument and close the actuator heat block to start the test.
8. Clot formation is signaled by an audible tone, the actuator heat block opens, and the results are displayed.

10 Results

The response of the activated clotting time to heparin varies considerably from individual to individual.^{2,4,6} Various drugs affect the activated clotting time.¹³ A number of other factors can affect the response of an individual's activated clotting time to heparin (for example, antithrombin III levels, heparin potency, coagulation factor deficiencies, sample activation, consumptive coagulopathies, excessive sample dilution, sample temperature). These factors should be taken into account when interpreting the results of activated clotting time tests.^{8,9,11}

Duplicate channel results should fall within 10% of each other for baseline (unheparinized) samples and within 12% of each other for extended (anticoagulated) samples.

Result calculation example:

Channel 1 clotting time: 210 seconds

Channel 2 clotting time: 200 seconds

Mean clotting time: 205 seconds

Difference: 10 seconds

12% of the mean: 25 seconds

These results are acceptable since the difference of 10 seconds is less than 12% of the mean (25 seconds).

11 Limitations of the Procedure

Medtronic ACT cartridges are for in vitro diagnostic use only and are to be used only with the Medtronic ACT Plus instrument.

Cartridges that have been used for testing should be considered hazardous. It is the responsibility of the user to dispose of the cartridge in accordance with local regulations and institution procedures.

To obtain accurate and consistent activated clotting times, the following guidelines must be observed:

- Both the instrument heat block and cartridge temperature should be at $37^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ ($98.6^{\circ}\text{F} \pm 0.9^{\circ}\text{F}$).
- It is critical to clean the instrument as soon as possible if blood contaminates the light path.
- The kaolin activator must be thoroughly resuspended prior to introduction of the test sample.
- The test sample must be free of tissue thromboplastin and tested as soon as possible following collection. Refer to Section 8, Specimen Collection and Preparation.
- Patient diagnosis and medications should be noted. Medications can alter clotting times.¹³ Test samples that give unexplained values should be recollected and retested. If necessary, additional test samples should be collected and sent to the laboratory for analysis.
- To optimize precision, all technique variables should be held constant from test to test.
- A prolonged activated clotting time is not specific for heparin. Prolonged ACTs may be due to heparin, other anticoagulant drugs (warfarin, direct thrombin inhibitors), hemodilution, administration of large volumes of citrated blood products, some antiplatelet drugs, hypothermia, severe thrombocytopenia, markedly abnormal platelet function, or other coagulopathies.¹⁰

- Excessive blood dilution during bypass can lead to postsurgically hemodiluted patients. This dilution can be great enough to extend the activated clotting time.¹⁰ In severe instances, the dilution may lead to an unmeasurable clotting time.

Warning: Do not use the HR-ACT cartridge to monitor regional anticoagulation with citrate.

For ACT Plus instrument precautions and limitations, consult the ACT Plus instrument operator's manual.

12 Expected Values

Response to an increased concentration of heparin may vary between individuals.

Clotting times using High Range Activated Clotting Time Cartridge (HR-ACT), the Low Range Activated Clotting Time Cartridge (LR-ACT), and the Recalcified Activated Clotting Time Cartridge (RACT) were determined on 3 donors. Table 3 is an example of the response to heparin using those 3 donors. The HR-ACT and the LR-ACT cartridges were tested using fresh drawn whole blood, and the RACT cartridges were tested using citrated whole blood (3.2% sodium citrate). Response to an increased concentration of heparin may vary between individuals.

Note: Each laboratory/facility should establish its own expected (reference) range based on its local population.

Table 3.

High Range ACT Cartridge		Low Range ACT Cartridge		Recalcified ACT Cartridge	
Heparin Level ^a	Avg. Clot Time (sec)	Heparin Level ^a	Avg. Clot Time (sec)	Heparin Level ^a	Avg. Clot Time (sec)
Donor 1					
0.0	156	0.0	133	0.0	151
1.0	255	1.0	489	1.0	324
2.0	350	2.0	670	2.0	532
3.0	468	3.0	801	3.0	791
4.0	574	4.0	>999	4.0	>999
5.0	723	5.0	>999	5.0	>999
6.0	925	6.0	>999	6.0	>999
Donor 2					
0.0	130	0.0	108	0.0	124
1.0	198	1.0	429	1.0	275
2.0	289	2.0	641	2.0	396
3.0	350	3.0	769	3.0	525
4.0	425	4.0	>999	4.0	651
5.0	517	5.0	>999	5.0	736
6.0	635	6.0	>999	6.0	>999
Donor 3					
0.0	117	0.0	102	0.0	114
1.0	180	1.0	218	1.0	241
2.0	241	2.0	427	2.0	405
3.0	292	3.0	655	3.0	549
4.0	360	4.0	757	4.0	684
5.0	408	5.0	>999	5.0	925
6.0	466	6.0	>999	6.0	965

^a Heparin level in units of heparin per mL of whole blood. Heparin was added to each blood sample immediately after sample collection.

13 Specific Performance Characteristics

Medtronic Activated Clotting Time Cartridges in the clotting time range of 0 to 600 seconds typically do not exceed a variation of 12% between duplicate channels for an individual. In subjects with extreme sensitivity to heparin, an activated clotting time in excess of 600 seconds is generally not considered adequately reliable to predict heparin effect.

The HR-ACT cartridge is intended for monitoring heparin levels >1.0 u/mL using fresh drawn whole blood.

The LR-ACT cartridge is intended for therapeutic heparin monitoring (<1.5 u/mL) using fresh drawn whole blood.

The RACT cartridge is intended for therapeutic heparin monitoring (<1.5 u/mL) using citrated whole blood.

14 Quality Control

ACT liquid controls (Table 4) are used to verify instrument performance, cartridge performance, and user technique. Results should fall within the specified ranges listed by the manufacturer or ranges determined by the user.

Note: When using only liquid controls, 2 levels must be performed every 8 hours of patient testing.

Table 4.

Cartridge Type	CLOTtrac Controls Used
High Range Activated Clotting Time (HR-ACT)	HR-ACT Control (Normal) HR Abnormal Control
Low Range Activated Clotting Time (LR-ACT) Note: CaCl ₂ required	CWB Control (Normal) LR Abnormal Control
Recalcified Activated Clotting Time (RACT)	CWB Control (Normal) RACT Abnormal Control

The ACTtrac Electronic Quality Control is used to supplement liquid controls and does not replace them. Refer to the current regulatory agency guidelines on the acceptability and use of electronic controls.

15 Disclaimers of Warranty

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

DISCLAIMER OF WARRANTY

ALTHOUGH THE ACTIVATED CLOTTING TIME CARTRIDGES, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

DISCLAIMER OF WARRANTY

ALTHOUGH THE ACTIVATED CLOTTING TIME CARTRIDGES, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Tabela 3.

Cartucho ACT - valores altos		Cartucho ACT - valores baixos		Cartucho ACT recalcificado	
Nível de heparina ^a	Tempo médio de coagulação (s)	Nível de heparina ^a	Tempo médio de coagulação (s)	Nível de heparina ^a	Tempo médio de coagulação (s)
Doador 1					
0,0	156	0,0	133	0,0	151
1,0	255	1,0	489	1,0	324
2,0	350	2,0	670	2,0	532
3,0	468	3,0	801	3,0	791
4,0	574	4,0	>999	4,0	>999
5,0	723	5,0	>999	5,0	>999
6,0	925	6,0	>999	6,0	>999
Doador 2					
0,0	130	0,0	108	0,0	124
1,0	198	1,0	429	1,0	275
2,0	289	2,0	641	2,0	396
3,0	350	3,0	769	3,0	525
4,0	425	4,0	>999	4,0	651
5,0	517	5,0	>999	5,0	736
6,0	635	6,0	>999	6,0	>999
Doador 3					
0,0	117	0,0	102	0,0	114
1,0	180	1,0	218	1,0	241
2,0	241	2,0	427	2,0	405

Tabela 3.

Cartucho ACT - valores altos		Cartucho ACT - valores baixos		Cartucho ACT recalcificado	
Nível de heparina ^a	Tempo médio de coagulação (s)	Nível de heparina ^a	Tempo médio de coagulação (s)	Nível de heparina ^a	Tempo médio de coagulação (s)
3,0	292	3,0	655	3,0	549
4,0	360	4,0	757	4,0	684
5,0	408	5,0	>999	5,0	925
6,0	466	6,0	>999	6,0	965

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



1. Hattersley PG. Activated coagulation time of whole blood. *JAMA*. 1966;196:436.
2. Bull BS, Huse MH, Brauer FS, Korpman RA. Heparin therapy during extra-corporeal circulation. II. The use of a doseresponse curve to individualize heparin and protamine dosage. *J Thorac Cardiovasc Surg*. 1975;69:685-689.
3. Dauchot PJ, Berzina-Moettus L, Rabinovitch A, Ankeney JL. Activated coagulation and activated partial thromboplastin times in assessment and reversal of heparin-induced anticoagulation for cardiopulmonary bypass. *Anesth Analg*. 1983;62:710-719.
4. Effenev DJ, Goldstone J, Chin D, Krupski WC, Ellis RJ. Interoperative anticoagulation in cardiovascular surgery. *Surgery*. 1981;90:1068-1074.
5. Esposito RA, Culliford AT, Colvin SB, Thomas SJ, Lackner H, Spencer FC. The role of the activated clotting time in heparin administration and neutralization for cardiopulmonary bypass. *J Thorac Cardiovasc Surg*. 1983;85:174-185.
6. Bull BS, Korpman RA, Huse WM, Briggs BD. Heparin therapy during extra-corporeal circulation. I. Problems inherent in existing heparin protocol. *J Thorac Cardiovasc Surg*. 1975;69:674-684.
7. Hattersley PG, Mitsuoka JC, King JH. The use of the activated coagulation time (ACT) as a monitor of heparin therapy. In: Lundblad RL, Brown WV, Mann KG, Roberts HR, eds. *Chemistry and Biology of Heparin*. North Holland: Elsevier; 1981;659.
8. Henry JB. *Clinical Diagnosis and Management by Laboratory Methods*. Vol 1. 16th ed. Philadelphia, PA: WB Saunders Company; 1979.
9. Moorehead MT, Westengard JC, Bull BS. Platelet involvement in the activated coagulation time of heparinized blood. *Anesth Analg*. 1984;63:394-398.
10. Niinikoski J, Laato M, Laaksonen V, Jalonen J, Ingerg MV. Use of activated clotting time to monitor anticoagulation during cardiac surgery. *Scand J Thor Cardiovasc Surg*. 1984;18:57-61.
11. Reno WJ, Rotman M, Grumbine FC, Dennis LH, Mohler ER. Evaluation of the BART test (a modification of the wholeblood activated recalcification time test) as a means of monitoring heparin therapy. *Am J Clin Path*. 1974;61:78-84.
12. Roth JA, Cukingnan RA, Scott CR. Use of activated coagulation time to monitor heparin during cardiac surgery. *Ann Thorac Surg*. 1978;28:69-72.
13. Young DS, Pastaner LC, Gibberman V. Effects of drugs on clinical laboratory tests. *Clin Chem*. 1975;21:1D.

ACT_Perf_bcv.fm 5/9/14 1:36 pm
6 x 8 inches (152.4 mm x 203.2 mm)

Medtronic Confidential
CS0032

M956357A001 Rev. 1A

Medtronic

Europe

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH - 1131 Tolochenaz
Switzerland
Internet: www.medtronic.co.uk
Tel. 41-21-802-7000
Fax 41-21-802-7900

EC	REP
----	-----

Authorized Representative in the European Community

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Asia-Pacific

Japan

Medtronic Japan
Comodio Shiodome 5F
2-14-1 Higashi-Shimbashi, Minato-ku
Tokyo 105-0021
Japan
Tel. 81-3-6430-2011
Fax 81-3-6430-7140

Australia

Medtronic Australasia Pty. Ltd.
97 Waterloo Road
North Ryde NSW 2113
Australia
Tel. 61-2-9857-9000
Fax 61-2-9878-5100

Asia

Medtronic International Ltd.
49 Changi South Avenue 2
Singapore 486056
Tel. 65-6436-5000
Fax 65-6776-6355

Americas

Latin America

Medtronic Latin America
3750 NW 87th Avenue
Suite 700
Miami, FL 33178
USA
Tel. 305-500-9328
Fax 786-709-4244

Canada

Medtronic of Canada Ltd.
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
Tel. 905-460-3800
Fax 905-826-6620
Toll-free: 1-800-268-5346

United States

Manufacturer:
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
Internet: www.medtronic.com
Tel. 763-514-4000
Fax 763-391-9100
Toll-free: 1-800-328-2518
(24-hour consultation service)

Manufactured In:

Medtronic Perfusion Systems
7611 Northland Drive
Minneapolis, MN 55428
Tel. 763-391-9000
Fax 763-391-9100
Customer Service and Product Orders
Toll-free: 1-800-854-3570
Visit Medtronic at
www.perfusion.medtronic.com

* M 9 5 6 3 5 7 A 0 0 1 *

© 2004, 2007, 2012, 2013 Medtronic, Inc.

M956357A001 Rev. 1A