EverCross™

0.035 in OTW PTA Dilatation Catheter

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Instructions for use

Device description

Dilatation balloon catheters are used to exert radial force to dilate narrow vessel segments. The EverCross™ 0.889 mm (0.035 in) OTW PTA dilatation catheter is an over-the-wire (OTW) dual lumen catheter (1). A semi-compliant inflatable balloon is mounted at the distal end (2). The distal atraumatic tip is tapered (3). The distal portion of the catheter has a hydrophilic coating (4).

The luer marked **THRU** (6) on the manifold (5) is the proximal opening for the central lumen. The central lumen ends at the distal tip of the catheter. The catheter passes over a 0.889 mm (0.035 in) maximum diameter guidewire. The luer marked **BALLOON** (7) is used to inflate and deflate the dilatation balloon, using a mixture of contrast media and saline solution. The balloon has 2 radiopaque marker bands (8) for positioning the balloon relative to the stenosis. The dilating or working section of the balloon (9) is marked by the radiopaque marker bands.


Figure 1. EverCross 0.889 mm (0.035 in) OTW PTA dilatation catheter

The EverCross catheter is available in multiple balloon sizes. Nominal balloon diameter and length are printed on the strain relief (10). Refer to the package label for information about catheter length and sheath compatibility.

All balloons distend to sizes above the nominal size at pressures greater than the nominal pressure. Consult the balloon compliance chart that is packaged with the device for balloon diameters at given pressures.

How supplied

- The EverCross catheter is supplied sterile for one-time use only.
- The EverCross catheter is packaged with a protective sheath positioned over the balloon.

Storage

Store the EverCross catheter in a cool, dry, and dark place, at room temperature. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

Indications for use

The EverCross 0.889 mm (0.035 in) OTW PTA dilatation catheter is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Contraindications

The EverCross 0.889 mm (0.035 in) OTW PTA dilatation catheter is contraindicated for use in coronary arteries or the neurovasculature, or when unable to cross the target lesion with a guidewire.

Warnings

- This device is supplied sterile for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.
- The EverCross catheter should only be used by physicians who have experience with and who understand the clinical and technical aspects of percutaneous transluminal angioplasty.
- To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Never use air or any gaseous medium to inflate the balloon.
- Use direct fluoroscopic observation while manipulating the EverCross catheter in the vascular system.
- Do not manipulate the balloon when it is inflated.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the lumen. If resistance is felt, carefully withdraw the dilatation catheter.
- Consider precautions to prevent or reduce clotting. Physician experience and discretion determines the appropriate anticoagulation therapy for each patient.
- Use the EverCross catheter with caution for procedures involving calcified lesions or synthetic vascular grafts.
- Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.

Precautions

- Federal law restricts this device to sale by or on the order of a physician.
- Balloons should not be inflated in excess of the rated burst pressure.
- Do not use contrast media that are contraindicated for intravascular use with this device.
- Inspect the packaging and the device before use. Do not use the device if the packaging or the device is damaged.
- Do not attempt to pass the EverCross catheter through a sheath introducer that is smaller than the size indicated on the packaging label. Refer to the product labeling for sizing information.
- Ensure that the guidewire is in place before changing the position of the balloon catheter.
- Larger and longer models of the EverCross catheter may exhibit longer deflation times, particularly on long catheter shafts.
- Do not retract the balloon catheter after stent postdilatation unless the balloon is free from the stent and fully deflated under vacuum.
- To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Do not wipe the catheter surface with dry gauze.
- This device may be a biohazard after use. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Potential complications /adverse effects

The following complications may result from balloon dilatation procedures. Complications or adverse effects may not be limited to this list:

- Allergic reaction to device materials or procedure medications
- Aneurysm
- Arrhythmia
- Arteriovenous fistula
- Artery perforation or rupture
- Bleeding requiring transfusion
- Contrast media reaction/renal failure
- Death
- Embolism
- Endocarditis
- Fever
- Hematoma
- Hemorrhage
- Hypertension / hypotension
- Infection
- Inflammation
- Intraluminal thrombus
- Myocardial infarction
- Pseudoaneurysm
- Sepsis
- Shock
- Stroke
- Transient ischemic attack
- Venous thromboembolism

Device selection

- The diameter of the expanded balloon should not exceed the diameter of the artery immediately distal to or proximal to the stenosis.

Warning: To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.

- Verify that the selected accessories can accommodate the balloon catheter based on sizing information on the package labeling. Refer to the product labeling for size information.
Caution: Inspect the packaging and the device before use. Do not use the device if the packaging or the device is damaged.

Recommended items

Prepare the following items using a sterile technique:

- 10 cc syringe filled with sterile heparinized saline
- Three-way stopcock
- Contrast media (The standard solution of inflation media is a 1:1 mixture of contrast media and normal saline.)
Caution: Do not use contrast media with this device that is contraindicated for intravascular use with this device.
- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

Preparation

Note: As supplied, the balloon and the guidewire lumen of the EverCross catheter contains air. This air must be displaced to make certain that only liquid fills the balloon while the catheter is in the bloodstream. Refer to the **Displacing air** section for air removal instructions.

1. Carefully remove the catheter from the inner package.
2. Remove and discard the protective tubing.

Activating the hydrophilic coating

Activate the hydrophilic coating before inserting the catheter. To activate the hydrophilic coating, immerse the catheter in normal saline for approximately 30 to 60 seconds, or wipe the catheter shaft with a saturated gauze sponge.

Caution: Do not wipe the catheter surface with dry gauze.

Displacing air

1. Flush out the lumen for the guidewire with heparinized saline solution through the wire luer marked **THRU**.
2. Mix equal volumes of contrast media and normal saline.
3. Fill a 10 cc syringe with approximately 4 cc of contrast media solution.

Warning: Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.

Caution: Do not use contrast media that is contraindicated for intravascular use with this device.

4. Evacuate air from the balloon segment:
 - a. Attach syringe to the balloon luer, marked **BALLOON**. Hold the syringe pointing down with the plunger pointed upward.
 - b. Apply negative pressure and aspirate the lumen for 15 seconds. Slowly release the pressure to neutral so that contrast media fills the shaft of the dilatation catheter.
 - c. Disconnect the syringe from the **BALLOON** luer of the dilatation catheter.
 - d. Remove all air from the syringe barrel. Reconnect the syringe to the **BALLOON** luer of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
 - e. Slowly release the device pressure to neutral.

Note: On first use, avoid any positive pressure to avoid unfolding and partially inflating the balloon outside of the vessel.
5. Repeat steps 4b to 4e as necessary to remove all air from the system.
6. Replace the syringe with a manometer-controlled dilatation system. Do not introduce air into the system.

Insertion and dilatation

1. If applicable, use an appropriately sized introducer sheath to introduce the EverCross catheter percutaneously. Refer to product labeling for sizing information.

Caution: Do not attempt to pass the EverCross catheter through a sheath introducer that is smaller than the size indicated on the packaging label. Refer to product labeling for sizing information.

2. After the guidewire is in place through the lesion, place the prepared EverCross catheter over the guidewire and advance the catheter tip to the introduction site. Use a guidewire that is of suitable length to allow for maintaining control and positioning the guidewire at all times.

Warning: Use direct fluoroscopic observation while manipulating the EverCross catheter in the vascular system.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the lumen. If resistance is felt, carefully withdraw the dilatation catheter.

Caution: To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

3. Position the catheter with the center of the balloon in the stenotic area of the vessel. The radiopaque marker bands indicate the working length of the balloon.

Caution: Ensure that the guidewire is in place before changing the position of the balloon catheter.
4. When the radiopaque marker bands are appropriately positioned, inflate the balloon to dilate the target area.

Warning: Never use air or any gaseous substances to inflate the balloon.

Warning: Do not manipulate the balloon when it is inflated.

Caution: Balloons should not be inflated in excess of the rated burst pressure.

5. Using the inflation device, deflate the balloon by aspirating the catheter.

Caution: Larger and longer models of the EverCross catheter may exhibit longer deflation times, particularly on long catheter shafts.
6. Withdraw the catheter while maintaining a vacuum in the balloon.

Caution: Do not retract the balloon catheter after stent postdilatation unless the balloon is free and fully deflated under vacuum.

Note: Multiple balloon inflations and deflations may cause resistance during device withdrawal. If resistance is encountered, use a gentle counterclockwise twisting motion to ease withdrawal through the sheath. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath together, while maintaining vessel access with the guidewire.

7. Use angiography to verify dilatation results.

Disposal

Caution: This device may be a biohazard after use. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Warranty

This product is subject to standard Medtronic warranty terms.

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© 2018 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3, Inc.