
GE HealthCare

SIGNA™ Champion

Fișă tehnică

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SECRET DE AFACERI

Cuprins

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- 2D DT-EPI
- 2D DW-EPI
- 2D FGRE
- 2D FIESTA
- 2D FIESTA FastCARD
- 2D FIESTA FastCINE
- 2D FIESTA Fat Sat
- 2D FRFSE
- 2D FRFSE-XL IDEAL
- 2D FSE IDEAL
- 2D FSE
- 2D FSE Double IR
- 2D FSE-IR
- 2D FSE Triple IR
- 2D FSE-XL IDEAL
- 2D FSPGR
- 2D GRE-EPI
- 2D MDE
- 2D MFGRE
- 2D SE-EPI
- 2D SSFSE
- 2D SSFSE 3-Plane
- 2D SSFSE-IR
- 2D SSFSE MRCP
- 2D T1FLAIR
- 2D T2MAP
- 3D BRAVO
- 3D COSMIC
- 3D Cube T1
- 3D Cube T2
- 3D Cube T2FLAIR
- Cub 3D DIR
- Cub 3D PD
- 3D Delta Flow
- 3D Dual Echo
- 3D Fast TOF GRE
- 3D Fast TOF SPGR
- 3D FGRE
- 3D FGRE IDEAL
- 3D FIESTA
- 3D FIESTA-C
- 3D FRFSE
- 3D FRFSE MRCP
- 3D FSPGR
- 3D FSPGR IDEAL
- 3D Heart
- 3D LAVA
- 3D LAVA FLEX
- 3D MDE
- 3D MERGE
- 3D QuickSTEP
- 3D SWAN
- 3D TOF GRE
- 3D TOF SPGR
- 3D TRICKS
- 3D Velocity Inflow
- 3D VIBRANT
- 3D VIBRANT FLEX
- Cine IR
- eDWI
- Fast 2D Phase Contrast
- FGRE Timecourse
- IFIR
- Inhance Inflow
- PROPELLER MB
- SWAN 2,0
- PS-MDE
- BB SSFSE
- 3D PROMO
- DISCO
- DW Duo (LX DWI Propeller)
- Flex
- HyperSense
- IR & SR Prepared
- PROMO

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PRODUCT DETAILS

I. SIGNA Champion 1.5T 30.1 Winner EMEA

Qty.	Catalog Ref	Description
1	S7530HW	SIGNA™ Champion 1.5T 64-channel
1	S7530HY	SIGNA™ Champion 1.5T Magnet with DPP - US
1	S7530RF	AIR™ Recon DL 2D and 3D and AIRx™ Package
1	S7750BF	Spectro EL pack
1	M7000JC	3D ASL
1	S7529HS	HyperSense 2.0
1	M7006NA	1.5T 16-channel AIR Anterior Array (AA)
1	M7100SG	1.5T 16-channel Shoulder Coil by NeoCoil
1	M6000CB	1.5T 8ch T/R Knee Coil
1	M7006LF	1.5T 16-ch Medium Flex Extremity Coil (without Interface)
1	M7006LG	1.5T 16-ch Large Flex Extremity Coil (without Interface)
1	M7006LH	1.5T Flex Interface Module 16-channel
1	M7005BE	Positioner for 16ch Flex coils
1	M3335CE	1.5T Calibration Phantoms
1	M66022TA	MR 30.1 for SIGNA™ System Software - China
1	M66022ED	SIGNA MR30.1 R04 eDelivery item - Champion
1	M70024AB	SIGNA™ Hero Host PC and Operator Console - China
1	M71080EB	DL Gen 7 Performance ICN - China
1	M66022SP	SIGNA™ Champion Cooling Cabinet (ICC) Assembly - China
1	M60122AT	F50SH Compressor for Cooling Type 2: 380/400/415 V, 50 Hz; 460-480 V, 60 Hz - US
1	M66023PC	Wired Gating Site Kit - China
1	M66022LL	SIGNA™ Champion Low-Height Fixed Table
1	M6602MC	High Channel Professional System
1	M70072GA	Remote Site Gradient Cable Kit Placeholder – SELECTION REQUIRED – Scan Room
1	M6001AA	Vent Adapter, Standard 8" Straight Up
1	M70022MB	Main Disconnect Panel - CE 380V/400V/415V 50/60Hz with RCD
1	M6001AG	International Shipping Collector
1	M1060KN	Remote Magnet Rundown Unit Kit
1	M1000MW	Operator Console Table
1	M7007NB	Air AA MCQA test PA Phantom Positioner
1	M50002NB	Romanian Language Kit
1	R34023AC	Secure Service Access key for warranty period
1	W99991MR	Continuity Standard for MR
1	E8823NA	MRI Audio 1505 Complete system
1	E8823NE	MRI Audio Over ear headphone
1	E80331BE	Physician Mobile Arm Chair
1	M3335CA	Calibration Kit Phantom Holder Cart
1	E8812CF	MR CCTV System with 17 inch LCD Monitor
1	E8811EJ	Small Bellows Assembly for MRI Respiratory/Cardiac Gating
1	M81511KA	AW VS7 workstation
1	M80281AA	AW VolumeShare 7 Monitors
1	M81501PG	AW HP Workstation Power Cords Kit for EMEA and ASIA
1	M80501DV	Field Engineer Letter - Operator Manual on paper for AW

CONFIDENTIAL

SECRET DE AFACERI

1	M30321CJ	ReadyView NRR
1	M30321BX	Brain View
12	W50092MR	24 Months warranty in total
4	A82000MR	16 Get Started Education Credits MR
1	A82160MR	160-Training Credits Education Package for MR
1	A33331MR	Digital Academy in MR

CONFIDENTIAL

SECRET DE AFACERI

LONG DESCRIPTIONS

All pictures featured are for illustration purposes only, and do not necessarily relate to products or services mentioned in the text.

S7530HW - SIGNA™ Champion 1.5T 64-channel

The SIGNA™ Champion from GE HealthCare, fueled by the MR 30 for SIGNA™ productivity platform, is designed to deliver uncompromised clinical performance and true patient comfort. The SIGNA™ Champion s-catalog comprises the system electronics, computing platform and operating software which includes:

- TDI 2.0 RF System
- TDI Coil Suite
- Quiet Technology

TDI 2.0 RF System:

TDI 2.0 is GE new generation digital-in-digital-out RF platform, which re-optimizes the Total Digital Imaging technology. TDI 2.0 brings higher precision and greater image quality to the scanning. It combines the power of Direct Digital Transmit (DDT) for 100% better stability and 200% higher accuracy, Direct Digital Receive (DDR) to digitize inputs from each RF channel for an SNR increase of up to 30%.

Number of available RF receive Channels *: 64 channel

Maximum number of channels per Field of View without table movement each generating an independent partial image: 64

*Independent number of digitizers

- Receiver sampling per channel: 80 MHz @ 16 bits per channel
- ADC sampling resolution: 16 bit 48
- Digital quadrature demodulation
- Receive dynamic range: >165 dB/Hz
- Receive signal resolution: up to 32 bits

TDI eXpress Coil Suite:

The Total Digital Imaging Suite of coils is designed to enhance patient comfort and image quality while simplifying workflow. The Coil Package includes:

- Integrated T/R Body Coil
- TDI Posterior Array (PA)
- TDI eXpress Head-Neck Array (HNA)

The TDI eXpress24-element Posterior Array (PA) ensures high repeatability and reproducibility during patient positioning.

- Spine imaging - CT spine (with Head-Neck Array) and TL spine.
- Body and pelvic imaging - in conjunction with the TDI eXpress Anterior Array (AA).
- Mechanically and electrically robust design allows placement of other surface coils, such as breast array, knee array, quad extremity, wrist array and Flex Arrays over the embedded array.
- Significantly improves workflow, for not only the head-neck-spine, but also abdominal and pelvic imaging.
- The scannable SI coverage is 91 cm, making it ideally suited to whole body imaging.
- Optimized for parallel imaging.

CONFIDENTIAL
SECRET DE AFACERI

The 16-element TDI eXpress Head Neck Array (HNA) is designed for high resolution and high sensitivity imaging of head, neck and C-spine.

- Superb high-SNR brain applications.
- Extensive SI coverage from the aortic arch to the circle of Willis.
- Imaging of head, neck and c-spine with high element-density.
- Offers 16 channel imaging in the Head Neck Array mode.
- Provides excellent productivity during sessions that need seamless coverage of the Central Nervous system and abdomen.
- The open-face anterior attachment of the array is removable. This helps in achieving higher patient comfort while the TL spine is being scanned with the TDI PA or abdomen is being scanned with the 9 Element TDI eXpress Body Array.
- Removable lookout mirror.
- Optimized for parallel imaging.

To improve the patient experience and provide high image quality, no other component of an MRI system has greater impact than the magnet. The SIGNA Champion 1.5T system features a high-performance magnet that delivers a large field of view and a robust fat saturation required for abdominal, breast and off-centered FOV musculoskeletal imaging. The magnet geometry has been optimized to reduce patient anxiety by providing more space in the bore and more exams with the patient's head outside of the magnet. The 50 x 50 x 50 cm field of view provides uniform image quality and can reduce exam times since fewer acquisitions may be necessary to cover large areas of anatomy. Complemented by GE's active shielding technology, the SIGNA Champion has very flexible installation specifications to provide easy siting. And with zero-boil-off magnet technology, helium refills are effectively eliminated even during installation, thus reducing operating costs and maximizing uptime.

Quiet Technology:

GE has implemented Quiet Technology on critical components of the SIGNA MR system to reduce acoustic noise and improve the patient environment. This technology enables full use of the UHE Gradient Platform for excellent image quality, while maintaining a safe environment for the patient. The technology encompasses the gradient coil, RF body coil, and magnet mounting.

- > 50 cm spacious vertical space between the table-top and the gantry inside the bore
- In-bore lighting and ventilation
- Two-way In-bore Intercom system
- Feet-first positioning
- Laser alignments for axial and sagittal reference planes
- Dual sided controls

CONFIDENTIAL
SECRET DE AFACERI

In-Room Displays (IRD):

By consolidating all controls into one place, the In-Room Displays (IRD) provides real-time feedback to the operator to improve exam room efficiency. With controls available at either side of the magnet as standard, the technologist always has all the control he needs at his fingertips, irrespective of which side he is operating from. The display provides real time interaction with the scanner and the host computer. The user has direct control or selection of the following:

- Display of patient name, ID, study description
- Display and entry of patient weight
- Display and entry of patient orientation and patient position
- Respiratory waveform display

AIR Touch™ activates an optimized set of coils that is selected based on the patient's anatomy. This advanced technology selects from unlimited coil combinations to efficiently set up patients.

- AutoStart – initiate the scanner to automatically acquire, process, and network images
- Display connected coils and coil status
- Display of table location and scan time remaining
- Screen saver
- Control multiple levels of in-bore ventilation and lighting

MR 30 for SIGNA™ ADVANCED APPLICATIONS

The MR 30 for SIGNA™ clinical imaging tools are organized and optimized to address six clinical work areas: NeuroWorks, OrthoWorks, BodyWorks, OncoWorks, CVWorks and PaedWorks. Each clinical toolkit comprises pre-programmed protocols, clinical applications and visualization tools designed for the challenges of each imaging area. In addition, MR 30 for SIGNA™ provides advanced applications that extend and enhance the clinical capability and performance of the MR 30 toolkits (quoted and described separately).

Advanced Application for NeuroWorks

- eDWI enhanced diffusion with Multi-B value and SmartNEX
- DTI diffusion tensor imaging
- FiberTrak post-processing for diffusion tensor to display white matter tracking
- fMRI - Functional imaging or BOLD provides fMRI analysis using the correlation coefficient algorithm to analyse an image set.
- 3D SWAN 2.0 GRE-based multi-echo susceptibility imaging including phase image
- PROBE PRESS SV brain spectroscopy
- Inhance 2.0 non-contrast MRA suite (3D velocity, 2D inflow, inflow IR, and Deltaflow)

Advanced Applications OrthoWorks

- MAVRIC SL 3D FSE-based spectral imaging for MR-Conditional implants with T1, PD, T2 and STIR
- CartiGram T2 cartilage mapping

Advanced Applications for BodyWorks

- 3D LAVA GRE 2-point Dixon fat-water separation for dynamic or single-phase imaging (breath-hold or free-breathing)
- IDEAL FSE 3-point Dixon fat-water separation
- Flex 2-point Dixon fat-water separation for 2D FSE, 3D Cube and GRE
- Inhance 2.0 non-contrast MRA suite with 3D velocity, 2D inflow, inflow IR and Deltaflow
- StarMap iron assessment for liver and heart (acquisition)

Advanced Applications for OncoWorks

- eDWI enhanced diffusion with Multi-B value and SmartNEX
- 3D LAVA GRE 2-point Dixon fat-water separation for dynamic or single-phase imaging (breath-hold or free-breathing)

Advanced Applications for CVWorks

- Cine IR fast gradient echo with IR-prep pulse
- 2D MDE IR-prep and gated, fast gradient echo imaging with wide bandwidth suppression and single-shot
- 2D PS MDE phase sensitive tissue characterization with wide bandwidth suppression and single-shot

CONFIDENTIAL
SECRET DE AFACERI

S7530HY - SIGNA™ Champion 1.5T Magnet with DPP - US

1.5T Magnet:

- Manufactured by GE Healthcare.
- Operating field strength 1.5T (63.86 MHz).
- Active magnet shielding
- Zero boil-off Cryogenics.
- Magnet length 174cm.
- Magnet Weight 7,275 lbs (3,300 kg).
- Patient Bore Diameter 70cm.
- Maximum Field of View (x,y,z) 50 cm x 50 cm x 50 cm.

Magnet Homogeneity: Typical ppm and Guaranteed ppm shown.

- 10cm DSV 0.007 and 0.02.
- 20cm DSV 0.035 and 0.06.
- 30cm DSV 0.10 and 0.15.
- 40cm DSV 0.33 and 0.43.
- 45cm DSV 0.88 and 1.0.
- 48cm DSV 1.75 and 2.0.
- 50cm DSV 2.8 and 3.3.

DSV = Diameter Spherical Volume.

Fringe field (axial x radial):

- 5 Gauss = 4.0 m x 2.5 m.
- 1 Gauss = 5.8 m x 3.2 m.

UHE Gradient Technology:

The gradients of an MR system play a crucial role when it comes to imaging performance, throughput, and consistency during clinical practice. Gradient speed, accuracy, and reproducibility often determine the success of demanding acquisitions like diffusion imaging and Fiesta. SIGNA™ Champion introduces Ultra High Efficiency (UHE) gradient system that includes Intelligent Gradient Control technology. This novel technology enables the SIGNA™ Champion to deliver excellent TR and TE values that enable a superior clinical performance.

Gradient performance

- 35 mT/m - 140 T/m/s

M70024AB - SIGNA™ Hero Host PC and Operator Console - China

The Host Dell T5820 GOC Upgrade delivers the next generation reconstruction engine, core pulse sequences, applications and imaging options as well as workflow enhancements and visualization tools to enable high productivity with exceptional image quality through optimized toolkits.

M70022MB - Main Disconnect Panel - CE 380V/400V/415V 50/60Hz with RCD

The Main Disconnect Panel safeguards the MR system's critical electrical components, by providing complete power distribution and emergency-off control.

M1060KN - Remote Magnet Rundown Unit Kit

The Remote Magnet Rundown Unit supplements the magnet rundown unit provided with the MR system for installations with more than one magnet.

M1000MW - Operator Console Table

The Operator Console Table is designed specifically for the color LCD monitor and keyboard.

R34023AC - Standard Service License

The Standard Service License provides access to service tools used to perform basic level service on the Equipment and is included at no charge for the warranty period.

S7530RF - AIR™ Recon DL 2D and 3D and AIRx™ Package

Comprehensive AIR™ Recon DL package for 60cm 1.5T MR Systems

AIR™ Recon DL is a pioneering, deep-learning based reconstruction algorithm that improves SNR and sharpness of 2D, DWI, PROPELLER and 3D images. This propriety technique improves image quality at the foundational level by removing image noise and ringing artifacts while enabling shorter scan times. With AIR™ Recon DL, customers will be able to:

- Remove noise in the images through trained deep learning algorithms.
- Increase productivity by enabling shorter scan times
- Eliminate Gibbs and truncation artifacts with intelligent ringing suppression
- Delivers sharper, clearer, beautiful MR images
- Apply a tailored level of AIR™ Recon DL based on preference
- Enable the most commonly applied 2D, DWI, PROPELLER and 3D sequences without anatomical limitations
- Visualize AIR™ Recon DL images directly at the MR console without reconstruction delays

AIR x™ Auto Graphic Prescription package

AIR x™ Auto Graphic Prescription replaces traditional atlas-based methods with deep learning algorithms to automatically identify anatomical structures and prescribe slices locations for brain and knee exams. As a result of the deep learning algorithms, AIR x™ automatically adapts slice prescriptions for normal patient variants in various patient anatomies and structures to enable consistency and productivity for slice positioning from technologist to technologist, patient to patient and the same patient overtime. The AIR x™ Auto Graphic Prescription package provides solutions for two high-volume exams, brain, and knee.

- AIR x™ auto graphic slice prescription for brain exams
- AIR x™ auto graphic slice prescription for knee exams

NOTE:

AIR™ Recon DL requires GEN 7 DL ICN, and AIR™ Recon DL PROPELLER and 3D licenses also require AIR™ Recon DL 2D license. This Application bundle is only available for SIGNA™ Explorer and Creator systems when sold with MR30 Application software.


CONFIDENTIAL
SECRET DE AFACERI

S7750BF - SPECTRO ELITE PACK


PROBE-PRESS and STEAM

PROBE-PRESS and STEAM Single-Voxel Spectroscopy allows the user to evaluate non-invasively and in vivo the relative concentrations of metabolites. It acquires and displays volume-localized, water-suppressed ^1H spectra in single-voxel mode. This package includes PROBE-P (PRESS) and PROBE-S (STEAM) pulse sequences, as well as automated reconstruction, acquisition set-up and graphic prescription of spectroscopic volumes.

PROBE 2D-CSI

PROBE 2D CSI expands proton brain spectroscopy capability enabling simultaneous acquisition of multiple in-plane voxels. PROBE 2D CSI uses the PRESS pulse sequence to acquire and display volume-localized, water suppressed ^1H spectra in a multi-voxel mode for the non-invasive assessment of in vivo metabolites. Metabolite maps are automatically generated in FuncTool on the operator console.

PROBE 3D-CSI

This option extends your PROBE-P 2D CSI spectroscopic capabilities by allowing the user to perform 3-dimensional multi-voxel acquisitions. Post-processing, including the creation of metabolite maps, is automatically generated with the FuncTool Performance Package (included as part of Scan Tools).

M7000JC - 3D ASL

3D ASL utilizes water in arterial blood as an endogenous contrast media to help visualize tissue perfusion and provide quantitative assessment of cerebral blood flow (CBF) in ml/100 g/min. The quantitative CBF maps can be generated and stored in DICOM format. 3D ASL deploys stacked spiral FSE readout with modulated flip angle to acquire 3D data with increased SNR and less image distortion compared to conventional 2D EPI-based ASL techniques. A pulsed-continuous labeling is applied to label arterial blood close to the imaging volume thus improving conspicuity of flowing blood. Selective, interwoven pulses are then used to saturate and invert the imaging volume, in order to achieve better background suppression, and reduce sensitivity to motion.

The isotropic 3D volume data can be reformatted to axial, sagittal, coronal or oblique planes. 3D ASL helps generate robust, reproducible images and perfusion maps with high SNR, reduced motion artifacts and less distortion in high magnetic susceptibility regions.

S7529HS - HyperSense 2.0

HyperSense 2.0 package is GE's compressed sensing application providing scan time reduction technique while maintaining SNR through an iterative reconstruction algorithm. This application can deliver higher spatial resolution images or reduced scan times, enabling faster imaging without the penalties commonly found with conventional parallel imaging.

HyperSense has been expanded to include T1 acquisitions including MP-RAGE & BRAVO for neuro imaging and LAVA, LAVA-Flex, DISCO and DISCO-Flex for body applications, and Vibrant for breast applications. In addition, HyperSense is now compatible with other 3D gradient echo sequences, such as MERGE, FIESTA and COSMIC.

CONFIDENTIAL

SECRET DE AFACERI

GE HealthCare

M7006NA - 1.5T 16-channel AIR Anterior Array (AA)

The 16-channel AIR Anterior Array (AA) is the next generation anterior array coil that allows flexibility in any direction to conform to the patient's anatomy. Based on the innovative AIR™ Coil technologies, the 1.5T 16ch AIR AA provides excellent image quality and acceleration performance, while improving the overall patient and user experience. The coil has been designed to adapt to various patient shapes and sizes, expanding positioning versatility.

M7100SG - 1.5T 16-channel Shoulder Coil by NeoCoil

The 1.5T Shoulder Coil by NeoCoil consists of a soft and light anterior array paired with a formed posterior array that together are designed to aid flexible patient positioning and heightened comfort. The coil is a phased array design with 16-channel receive and parallel imaging compatibility to also deliver enhanced SNR and speed for shoulder imaging at 1.5T.

- Elements: 16
- Dimensions: 28.2 cm x 35.3 cm x 18.2 cm
- Weight: 3.0 kg
- Parallel imaging compatible

M6000CB - 1.5T 8ch T/R Knee Coil

The 1.5T T/R Knee Coil is designed for high-definition MR imaging of the knee. This coil uses unique hybrid technology with separate birdcage coils for transmit and receive functions. Designed uniquely for GE, the 8-element receive coil delivers 30% to 100% more SNR than the standard extremity coil. The array is compatible with PURE for uniform signal intensity and ASSET for accelerated imaging speed.

M7006LF - 1.5T 16-ch Medium Flex Extremity Coil (without Interface)

The high density 16-channel receive coil is designed to give high quality images in a wide range of applications. The high degree of Flexibility was achieved by removing all non-essential electronics to an external interface assembly, ensuring reduced weight on the patient and better conformance to the anatomy. The high degree of Flexibility is particularly advantageous when imaging patients that do not fit the constraints of rigid coils, improving patient and technologist experience, and enabling most exams to be completed with the same level of image quality expected from dedicated rigid coils. This coil covers a broad range of muscular skeletal applications, including hand, wrist, elbow, shoulder, hip (unilateral and bilateral), knee, ankle, and foot. In addition, the coil's versatility has been shown in a range of general purpose applications that include head, neck, and spine exams.

Requires Interface Module (sold separately).

M7006LG - 1.5T 16-ch Large Flex Extremity Coil (without Interface)

The high density 16-channel receive coil is designed to give high quality images in a wide range of applications. The high degree of Flexibility was achieved by removing all non-essential electronics to an external interface assembly, ensuring reduced weight on the patient and better conformance to the anatomy. The high degree of Flexibility is particularly advantageous when imaging patients that do not fit the constraints of rigid coils, improving patient and technologist experience, and enabling most exams to be completed with the same level of image quality expected from dedicated rigid coils. This coil covers a broad range of muscular skeletal applications, including hand, wrist, elbow, shoulder, hip (unilateral and bilateral), knee, ankle, and foot. In addition, the coil's versatility has been shown in a range of general purpose applications that include head, neck, and spine exams.

Requires Interface Module (sold separately).

CONFIDENTIAL
SECRET DE AFACERI

M7005BE - Positioner for 16ch Flex coils

Flex Array Positioner

The Flex Array Positioner is a multipurpose support for a broad range of exams including foot, ankle, forefoot, knee, and head. The Flex Array Positioner is initially dedicated for use with the 16ch Flex Coils. A dedicated forefoot attachment allows the flex array elements to be wrapped tightly around the foot, yielding improved image quality. A repositionable support pad in the foot and ankle attachment allows for selection of a 90-degree position, or a relaxed position of the ankle. The pads and straps included with the stabilizer facilitate rapid setup and allow for flexibility in how the anatomy is secured leading to an increased reliability.

W99991MR - Continuity Standard for MR

Continuity Standard proactively provides critical and non-critical security patches to the Equipment's operating system, as well as software updates consisting of error corrections or modifications to the Equipment to the extent and in the manner further described here below. Continuity Standard shall be provided during Warranty and after Warranty, only for systems under a valid GEHC Service contract ("Continuity Standard Term").

GE Healthcare reserves the right to determine which software updates are eligible for applicable Equipment.

Failure of Customer to maintain connectivity may eliminate Customer's eligibility to receive updates.

1. Continuity Standard

- 1.1. **OS Patches.** During Continuity Standard Term, GE Healthcare will provide Critical and Non-Critical security patches ("Patches") to the Equipment's operating system **when eligible and available** that have been released by the operating system's OEM and validated by GE Healthcare to be compatible with the Equipment, but only if the operating system is supported by the operating system's OEM. "Critical Patches" address operating system vulnerabilities determined by GE Healthcare to be critically exploitable and/or have a potential critical impact if exploited. "Non-Critical Patches" address operating system vulnerabilities determined by GE Healthcare to be non-critical. GE Healthcare's installation of Non-Critical Patches is limited to one (1) per Agreement year. Installation of additional Non-Critical Patches in an Agreement year will be at GE Healthcare's then-current rates. Installation of Patches will be performed by GE Healthcare during the next scheduled Equipment Planned Maintenance ("PM") if GE Healthcare is performing PM on the Equipment under this Agreement, or at a mutually agreed upon time.
- 1.2. **Software Updates.** During Continuity Standard Term, GE Healthcare will provide "**Software Updates**" which consist of any error correction or modification to Equipment that (a) address cyber vulnerability, or (b) maintain existing Software features and functionality developed for GE Healthcare's installed customer base, when available and eligible to the covered product. Software Updates are distinct from OS Patches which are identified above. Installation of Software Updates will be performed by GE Healthcare through eDelivery when available, during the next scheduled Equipment PM, if GE Healthcare is performing PM on the Equipment under this Agreement, during the installation of an OS Patch, or at another mutually agreed upon time. Additional hardware and/or software (including upgrades to third party software or operating system software) required for Software Updates, training, project management and integration services are excluded.
2. **Price and Payment.** Pricing for Continuity Standard is set out in the new equipment order as a separate line-item. No part of the Continuity Standard price is refundable or subject to reduction, and payments are not contingent on GE Healthcare's delivery, or Customer's acceptance, of any particular OS Patch, Software Updates, upgrades and enabling hardware, base-system software and/or hardware upgrade, or any other products or service.
3. **Term.** The start date of the Term will commence on the installation completion date and continue for the duration of the Warranty and/or the Service contract.
4. **Exclusions.** Except as identified above, Continuity Standard excludes: (i) Product replacements; (ii) Product upgrades; (iii) accessories, supplies and consumable items; (iv) any separately licensed software modules which provide additional functionality related to an application or feature for the hardware or Software; (v) advanced operating system applications or features; (vi) physicist testing and calibration; and (vii) training. Products are excluded from coverage under this Agreement and Customer is not entitled to any remedy if GE Healthcare's failure to perform hereunder is due to: (a) Customer cancellation, rescheduling, or inability of GE Healthcare to access the Product; (b) Customer's default; (c) improper care of the Product; or (d) any cause beyond GE Healthcare's control. GE Healthcare is not responsible for

CONFIDENTIAL

SECRET DE AFACERI

GE HealthCare

11/18

providing system database maintenance for Customer, including but not limited to, activities related to backup, new users, user privileges, physician list updates, and archive/data entry.

5. Miscellaneous.

- 5.1. Customer is responsible for: (i) Customer will provide GE Healthcare with, and maintain, a GE Healthcare-validated remote access connection if a Product has remote access capability; Failure of Customer to maintain connectivity may eliminate Customer's eligibility to receive some updates. This remote access and collection of machine data (e.g., temperature, helium level) will continue after the end of this Agreement unless Customer requests in writing that GE Healthcare disable it. (ii) site preparations, construction and rigging that may be required for Continuity Standard; (iii) ensuring that all data is appropriately backed up prior to installation of OS Patches, Software Updates; and (iv) purchasing any additional hardware and/or software (including upgrades to third party software or operating system software), training, project management and integration services if required for delivering Software Updates under Continuity Standard.
- 5.2. PATCHES, OS UPDATES, AND/OR SOFTWARE UPDATES, PROVIDED UNDER THIS AGREEMENT ARE "AS IS" AND "AS AVAILABLE" WITH NO WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF NON-INFRINGEMENT, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GE HEALTHCARE MAKES NO REPRESENTATIONS, WARRANTIES OR CONDITIONS THAT SUCH PATCHES, OS UPDATES, UPGRADES AND ENABLING HARDWARE, SOFTWARE UPDATES, OR EQUIPMENT WILL BE UNINTERRUPTED, TIMELY, SECURE, ERROR-FREE OR MEET CUSTOMER REQUIREMENTS OR ANY NATIONAL GUIDELINE OR INDUSTRY STANDARD.
- 5.3. CYBER-RELATED SECURITY IS NOT AN EXACT SCIENCE AND IS EVER EVOLVING IN THE TYPES, VULNERABILITIES, RISKS AND THREATS THAT ARE PRESENT. CONTINUITY DOES NOT PROVIDE COMPLETE OR COMPREHENSIVE PROTECTION AGAINST ALL POSSIBLE SECURITY VULNERABILITIES, RISKS, THREATS, BREACHES, SECURITY EVENTS, SECURITY INCIDENTS OR UNAUTHORIZED INTRUSIONS, AND IS NOT INTENDED TO BE A COMMITMENT, GUARANTEE, REPRESENTATION OR WARRANTY THAT CUSTOMER'S EQUIPMENT, NETWORK OR IT INFRASTRUCTURE WILL BE PROTECTED, AGAINST ALL POSSIBLE SECURITY VULNERABILITIES, RISKS, THREATS, BREACHES, SECURITY EVENTS, SECURITY INCIDENTS OR UNAUTHORIZED INTRUSIONS.
- 5.4. Products, equipment, components, software and/or hardware replaced or removed by GE Healthcare under this Agreement will become GE Healthcare's property.
- 5.5. Equipment that is declared end of life/support by GE Healthcare is not eligible for Continuity Standard.

E8823NA - MRI Audio 1505 Complete music system for Performance & Premium MRI systems

The MRI Audio premium sound system is designed for comfort and allows the patient to listen to music while being scanned in an MRI. The technologist is in full control of the system headphones, microphone, sound source and volume controls. Standard 3.5 mm plug for music source allows any compatible music player, tablet or phone. In-ear headphones work with any head coil.

Package includes:

- Digital amplifier
- iPad Mini
- iPad Mini mount with lock
- 3G transducer
- In-ear headphones, 29dB noise reduction
- Disposable ear tips (300 pairs)
- Technologist's speakers
- 6 ft RCA 3.5 mm cable
- Auto-voice/MIC adapter

1.

Compatible with GE MR systems:

- Discovery MR750
- Discovery MR750w

CONFIDENTIAL
SECRET DE AFACERI

- Optima MR450w
- Signa PET/MR
- Signa HDxt
- Signa Architect
- Signa Artist
- Signa Pioneer (PX26.1 and beyond)
- Signa Voyager

E8823NE - MRI Audio over-ear headphone (29dB NRR)

Premium over ear headphone with a one size fits all adjustable headband. Works for all MRI procedures except head coil exams. Rigorously tested to a 29dB noise reduction rating (NRR) and provides excellent sound quality when paired with the MRI Audio music system.

E80331BE - Physician Mobile Arm Chair

Overview

Our mobile chair is an excellent addition to your medical suite.

Specifications

- Armrests
- Castors
- Color: Charcoal gray
- Adjustable height: 19.7 in/50 cm to 24.8 in/63 cm
- Chair contains some ferro-magnetic material
- Not safe for use in an MRI Suite's Zone IV, use in control room only

M3335CA - Calibration Kit Phantom Holder Cart

Calibration Kit Phantom Holder Cart

E8812CF - MR CCTV System Single

MR CCTV System with 17 inch LCD Monitor.

This MRI closed-circuit TV camera package allows for visual monitoring of your patients during their scans. Ceiling or wall-mounted fixed camera is unobtrusive and transmits high-quality color images in real time to the operator's console where a 17" TFT LCD monitor is located.

Tested and approved by GE on a 3-Tesla magnet (camera and mount).

Includes:

- 17" TFT LCD color monitor
- Compact color camera with lens
- Signal interface box
- Exterior-mount power supply
- Mounting brackets for camera

CONFIDENTIAL
SECRET DE AFACERI

- One set of cables: 10', 40', 100'

E8811EJ - Small Bellows Assembly for MRI Respiratory/Cardiac Gating

Overview :

- Bellows assembly for respiratory gating and respiratory compensation packages
- Non-magnetic
- Sold individually; choose adult or small size
- Adult size also available separately (E8811ED)

Note: This product contains natural rubber latex which may cause allergic reactions.

- 600 pieces per case.
- Packs of 4.

A82000MR - 16 Get Started Education Credits MR

Customised Training:

Get the best out of your GE HealthCare equipment with LEVEL UP! A 3-step education programme tailored to your team needs.

1. **Get Ready:** Together we define your site application requirements, advise on pre-training via Digital Academy (LMS) and plan your Get Started sessions.
2. **Get Started:** Training and handover of your system by a member of our expert Clinical Education team.
3. **Equipment Lifecycle:** We provide you with privileged access to WeConnect, our user community and central hub for educational resources. We also offer opportunities to access our long-term learning solutions via LEVEL UP! Credits.

Continuous Education:

LEVEL UP! Credits may be used to access any of our educational solutions. These include (but are not limited to):

- On-site training at Your facility (1 day = 8 credits)
- Remote training sessions (1 hour = 1 credit)
- Classroom Session at a GE HealthCare Academy (1 day/1 attendee = 4 credits)
- Immersion training at a Partner Site (1 day/1 attendee – 4-8 credits)
- 12 months access to Digital Academy (10 learners + 1 manager = 8 credits)

Credits are valid for 12 months from the Product warranty start date when associated with a product purchase; or the effective date of the signed Agreement if purchased as standalone. Unused credits within this timeframe will expire without refund.

A33331MR - Digital Academy in MR

Your dedicated online learning platform

Maximize the use of your system with the Digital Academy online learning platform.
Develop your skills with diverse clinical and product training across all modalities.

Benefit from a year of access for 10 learners and 1 admin, enabling training in total autonomy & flexibility.
Assign courses adapted to the team's learning plan, track progress and earn certificate of completion for each course completed.

Target attendees: Radiographers, Technologists, Radiologists, Physicians, Cardiologists.

Certification: GE HealthCare Education Certificate.

Language: English and local Language when available.

Number of Participants: 10 learners and 1 admin (additional learners can be added for 8 education credits per 10 learners).

Expiry date: Access granted for 1 year – Possibility to renew every year.

M81511KA - AW VS7 workstation

AW VolumeShare 7 is a multi-modality image review, comparison and post processing workstation built with simplicity and power at its core. Powerful software is optimized to take advantage of state of the art 64 bit technology and multiple cores to ensure leading edge performance.

AW VolumeShare 7 features include:

Hardware:

- HP Z4G4 Workstation
- CPU: Intel Xeon W-2135 Six-Core @ 3.7 GHz with 8.25 MB L3 Shared Cache
- RAM: 64GB DDR4 2666 MHz ECC Registered DIMM
- Upgradeable to 64GB (8x8GB)
- Graphics: NVIDIA Quadro NVS P620 with 2 GB Video cards (optionally upgradeable with certain applications)
- 1x 256GB Solid State Drive for OS and Apps
- 2x 512GB Solid State Drive in RAID-0 for image cache

Software:

- GE Healthcare HELIOS 6 operating system
- Volume Viewer for advanced post-processing
- Demo Exams for training and exploration
- Fast access to information you need through optional RIS integration & priors post-fetch
- Efficient workflow through dynamic load, end review and Key Image Notes features
- Productivity package to pre-process exams and allow up to 8 simultaneous sessions
- Applications usage monitor to track and view usage of your system
- Smart layouts with Volume Viewer General review protocol that optimizes comparison and single exam layouts
- Enhanced multi-modality contouring tool with support for PET SUVs
- Support for external DICOM USB media and preference management tool to exchange preferences across users
- Support for optional, broad suite of multi-modality advanced applications

M80281AA - AW VolumeShare 7 Monitors

CONFIDENTIAL
SECRET DE AFACERI

VolumeShare 7 Monitors are two high-quality monitors offering bright and high contrast imagery suited to the display of medical images per the AW VolumeShare Indications for Use. Each provides a 19" 1280x1024 (5:4 aspect ratio) display that complies with international medical and patient safety standards and offers the following specifications:

- Maximum luminance (panel typical) : 330 nit
- DICOM Part 14 calibrated luminance: 215 nit
- Contrast ratio (panel typical) : 900:1
- An ambient light sensor
- Brightness non-uniformity (measured as per DIN6868-157) : +/-25%

M80501DV - Field Engineer Letter

Field Engineer Letter - Operator Manual on paper for AW

M30311CTED - Ready View Complete

Includes:

READY View

READY View AD/RD

READY View NRR

READY View R2 Error

READY View is an application designed to improve multi-parametric exams by enabling the analysis of MR generated data sets containing multiple images for each scan location. The MR data sets may be any of the following:

- o Time series.
- o Diffusion weighted scan.
- o Diffusion tensor scan.
- o Variable echo imaging.
- o Blood oxygen level dependent imaging.
- o Spectroscopy-Single Voxel and 2D or 3D CSI.

The READY View application provides a combination of protocols, applications and tools that enables a fast, easy and quantified analysis of the multiple data sets. The combination of acquired images, reconstructed images, calculated parametric images, tissue segmentation, annotations and measurement performed by the clinician enables a multi parametric analysis and may provide clinically relevant information for diagnosis.

READY View is available on the Advantage Workstation (AW) and the AW Server Enterprise Gen2. Through the AW server capabilities, READY View will be available throughout the enterprise and at any remote reading location for viewing and processing images. READY View offers the following:

Personalized workflow:

- o Provides guided workflows to help process MR data.
- o Provides custom workflows and applications with adjustable layouts, personalized parameter settings and custom review steps.

o Optimizes workflows by enabling access to READY View from the location of your choice (ie PC, laptop or PACS/RIS workstation).

High-Tech framework:

- o Provides additional clinical information for diagnosis by creating graphs (time curves) and color parametric images.
- o Offers a set of protocols for processing multi-parametric data.
- o Enables the fusion of color parametric images with anatomical 2D or 3D images by a simple "drag and drop."
- o Enables MR to MR image registration to reduce the effect of a patients' movement between series.
- o Provides tools and workflows for processing 2D and 3D CSI data.

Improved productivity:

- o Provides a simple workflow to review and manipulate acquired images (Maximum Intensity Projection (MIP), Multi Planar Reconstruction (MPR), 3D rendering) and to process and fuse functional data.
- o Automatically selects the most relevant protocols.
- o Enables Save State to restore the state of previously processed data.
- o Provides easy to use slide bars that enable real time parametric image calculation and/or segmentation.

M30321BX - Brain View

Brain View is a post processing image analysis software package that provides advanced techniques to aid in the diagnosis of neurological and oncological diseases. Brain View Plus is integrated with the READY View platform and benefits from the generic READY View tools and applications. In addition, it offers four advanced protocols:

FiberTrak: FiberTrak is an advanced post processing software program that reconstructs fiber pathways from the diffusion tensor series data, based on algorithmically predicted water molecule direction. FiberTrak provides:

- Advanced parameters to stop the tracking of fibers.
- A choice of color scales including a color orientation and a user-defined color scale.
- A palette of color to select the color for each individual tract.
- A review step to set symmetry axis of the brain to mirror a bundle of fibers to the contralateral hemisphere.
- A screen to display and export the statistics of all created tracts.
- Ability to segment all DTI maps (Fractional Anisotropy, Color Orientation, ADC, Volume Ratio Anisotropy, Exponential Attenuation, T2 weighted Trace, Isotropic Image, and Anisotropy Index) based on the fiber tract coordinates.
- The fiber maps may be saved as DICOM, secondary screen capture or JPEG image format.

Arterial Spin Labeling (ASL): ASL is an automated post processing software option that is designed to process an Arterial Spin Labeling series of MR images acquired in the brain.

- The CBF maps are automatically calculated and displayed when a 3D ASL series is selected from the browser and READY View is launched.
- Real time noise thresholding of color parametric images (CBF).
- Advanced parameters to apply correction to the CBF map (e.g. "Labeling Efficiency" to account for incomplete inversion during tagging).
- Review step to set symmetry axis of the brain to mirror 3D ROI to the contralateral hemisphere.

BrainStat and BrainStat AIF: BrainStat and BrainStat AIF are designed to process a time series of MR images acquired in the brain. They can be used to represent parametric images that are calculated from the variations in image intensity over time.

- BrainStat automatically generates processed maps that are based on a Gamma Variant Fitting (GVF) of the concentration curve.

- These processed maps include: relative Cerebral Blood Volume (rCBV), relative Cerebral Blood Flow (rCBF), Mean Transit Time (MTT), Time to Peak (TTP).
- BrainStat AIF automatically generates processed maps that are based on Arterial Input Functions. There are three modes to detect artery pixels.
- Auto Vessel Selection: arterial pixels are automatically detected and displayed.
- Semi-Auto Vessel Selection: the user selects an ROI and arterial pixels are automatically detected and displayed for the specified region.
- Manual Vessel Selection: the user defines an arterial pixel.
- These processed maps include: relative Cerebral Blood Volume (rCBV), relative Cerebral Blood Flow (rCBF), Mean Transit Time (MTT), Time to Peak (TTP), Bolus Arrival Time (BAT) and Time to maximum value (Tmax).

System Requirements:

- Brain View installation requires READY View as prerequisite.

Note: All software are Non-Transferable to other hardware and are Non-Returnable.

CONFIDENTIAL
SECRET DE AFACERI

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©2024 GE HealthCare. SIGNA, AIR, AIR x și AIR Touch sunt mărci comerciale ale GE HealthCare. Intel și Xeon sunt mărci comerciale ale Intel Corporation sau ale filialelor sale din SUA și/sau din alte țări. NVidia și Quadro sunt mărci comerciale ale NVIDIA Corporation. GE este o marcă comercială a General Electric Company, utilizată sub licență de marcă comercială.